# ARJOHUNTLEIGH

GETINGE GROUP

## ALPHA ACTIVE 3

Instructions for Use

Brugervejledning • Instrucciones de uso

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**English**

## Contents

<b>General Safety .....</b>	<b>iii</b>
<b>Introduction .....</b>	<b>1</b>
About this Manual .....	1
Intended Use .....	1
About Alpha Active 3 .....	1
Alpha Active 3 Pump .....	1
Alpha Active 3 Mattress Overlay .....	2
<b>Clinical Applications .....</b>	<b>4</b>
Indications .....	4
Contraindications .....	4
Cautions .....	4
<b>Installation .....</b>	<b>5</b>
Preparing the system for use.....	5
Installing the Mattress .....	5
<b>Controls, Alarms and Indicators .....</b>	<b>7</b>
Controls .....	7
Alarms and Indicators .....	8
<b>Mattress - Pump Operation .....</b>	<b>9</b>
Quick Start.....	9
Support Setting Procedure .....	9
Changes in Patient Position.....	10
To Deflate and Store the Alpha Active 3 Mattress.....	11
CPR Control.....	12
To activate CPR.....	12
To reset CPR.....	12
<b>Decontamination .....</b>	<b>13</b>
<b>Routine Maintenance .....</b>	<b>15</b>
Alpha Active 3 System.....	15
Alpha Active 3 Pump .....	15
Alpha Active 3 Mattress Overlay .....	15
Serial Labels .....	15
<b>Troubleshooting .....</b>	<b>16</b>
<b>Technical Description .....</b>	<b>17</b>
Pump .....	17
Pump Symbols.....	18
Mattress .....	18
Mattress Size Information .....	19
Cover Specification .....	19

Cleaning Symbols . . . . .	19
----------------------------	----

## **GENERAL SAFETY**

Before you connect the system pump to a mains socket, read carefully all the installation instructions contained within this manual.

The system has been designed to comply with regulatory safety standards including:

- EN60601-1:1990/A13:1996 and IEC 60601-1:1988/A2:1995
- UL60601-1, UL2601-1 and CAN/CSA C22.2 No. 601.1-M90
- EN60601-1:2006 and IEC 60601-1:2005
- AAMI/ANSI ES60601-1:2006 and CAN/CSA C22.2 No.60601.1(2008)

### **Safety Warnings**

- **It is the responsibility of the care giver to ensure that the user can use this product safely.**
- Whilst the patient is unattended, safety sides should be used based on clinical assessment and in line with local policy.
- Alignment of the bed frame, safety sides and the mattress should leave no gap wide enough to entrap a patient's head or body, or to allow egress to occur in a hazardous manner where entanglement with the mains power cable and tubeset or air hoses may result. Care should be exercised to prevent occurrence of gaps by compression or movement of the mattress. Death or serious injury may occur.
- Make sure that the mains power cable and tubeset or air hoses are positioned to avoid causing a trip or other hazard, and are clear of moving bed mechanisms or other possible entrapment areas. Where cable management flaps are provided along the sides of the mattress, these should be used to cover the mains power cable.
- Electrical equipment may be hazardous if misused. There are no user-serviceable parts inside the pump. The pump's case must only be removed by authorised technical personnel. No modification of this equipment is allowed.
- The mains power socket/plug must be accessible at all times. To disconnect the pump completely from the electricity supply, remove the plug from the mains power socket.
- The CPR control and/or the CPR indicator tag must be visible and accessible at all times.
- Disconnect the pump from the mains power socket before cleaning and inspecting.
- Keep the pump away from sources of liquids and do not immerse in water.
- Do not use the pump in the presence of uncontaminated flammable liquids or gasses.
- The cover of this product is vapour permeable but not air permeable and may present a suffocation risk.
- Only the pump and mattress combination as indicated by ArjoHuntleigh should be used. The correct function of the product cannot be guaranteed if incorrect pump and mattress combinations are used.

## **Precautions**

For your own safety and the safety of the equipment, always take the following precautions:

- Placing extra layers between the patient and the mattress potentially reduces the benefits provided by the mattress and should be avoided or kept to a minimum. As part of sensible pressure area care, it is advisable to avoid wearing clothing which may cause areas of localised high pressure due to creases, seams, etc. Placing objects in pockets should be avoided for the same reason.
- Do not expose the system, especially the mattress, to naked flames, such as cigarettes, etc.
- In the event of a fire, a leak in the seat or mattress could propagate the fire.
- Do not store the system in direct sunlight.
- Do not use phenol-based solutions to clean the system.
- Make sure the system is clean and dry prior to use or storage.
- Never use sharp objects or electrically heated under blankets on or under the system.
- Store the pump and mattress in the protective bags supplied.
- Pets and children must be supervised in the vicinity of the system.

## **Electromagnetic Compatibility (EMC)**

This product complies with the requirements of applicable EMC Standards. Medical electrical equipment needs special precautions regarding EMC and needs to be installed in accordance with the following instructions:

- The use of accessories not specified by the manufacturer may result in increased emissions by, or decreased immunity of, the equipment, affecting its performance.
- Portable and mobile radio frequency (RF) communications equipment (e.g. mobile/cell phones) can affect medical electrical equipment.
- If this equipment needs to be used adjacent to other electrical equipment, normal operation must be checked before use.
- For detailed EMC information contact ArjoHuntleigh service personnel.

## **Environmental Protection**

Incorrect disposal of this equipment and its component parts, particularly batteries or other electrical components, may produce substances that are hazardous to the environment. To minimise these hazards, contact ArjoHuntleigh for information on correct disposal.

## **Expected Service Life**

The **Alpha Active®** pump has an expected service life of seven years. To maintain the condition of the pump, have the pump serviced regularly according to the schedule recommended by ArjoHuntleigh.

Do NOT use unapproved accessories or attempt to modify, disassemble or otherwise misuse the **Alpha Active** system. Failure to observe this caution could result in injury, or in extreme cases, death.

## **Design Policy and Copyright**

® and ™ are trademarks belonging to the ArjoHuntleigh group of companies. As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice. The content of this publication may not be copied either whole or in part without the consent of ArjoHuntleigh.

© ArjoHuntleigh 2012.

# 1. Introduction

**About this Manual** This manual is your introduction to the **Alpha Active® 3**.

You must read and fully understand this manual before using the system.

Use this manual to initially set up the system, and keep it as a reference for day-to-day routines and as a guide to maintenance.


If you have any difficulties in setting-up or using the **Alpha Active 3** system, contact your local ArjoHuntleigh sales office, listed at the end of this manual.

**Intended Use** The intended use of this product is to prevent and/or manage pressure ulcers for patients up to 135 kg (200 lb).

The **Alpha Active 3** system should be used as part of a prescribed plan of care (refer to Page 4 “Indications”).

**About Alpha Active 3** The **Alpha Active 3** systems comprise of a mattress overlay and pump. The support system can be used on hospital and domestic beds.


**Alpha Active 3 Pump** The **Alpha Active 3** pump comprises of a moulded case with non slip feet on the base and integral hanging brackets.

The controls are situated on the front of the pump. An alarm system differentiates between normal operation and genuine system faults. If an alarm situation is detected an indicator will illuminate on the front of the pump and an audible warning will sound.

### **Alpha Active 3 Mattress Overlay**

The **Alpha Active 3** mattress overlay comprises the following components:

- Detachable Cover** The standard cover comprises of a 2-way stretch polyurethane (PU) coated knitted fabric zipped to a durable nylon base. The zips are protected by flaps to prevent ingress of contaminants and allow easy removal of the cover for cleaning.
- Cells** The mattress comprises of 16 polyurethane (PU) cells, 13 providing support to the user in either Alternating (Active) or Static (Reactive) mode and 3 Static head cells.
- CPR function** A CPR (Cardio-Pulmonary Resuscitation) control is positioned at the head end of the mattress to allow rapid deflation of the mattress overlay.
- Tube-set** The tube-set has a 2-way pneumatic connection which incorporates a flexible, compact anti-kink tube that is resistant to crushing and any subsequent obstruction of air flow.

When disconnecting the tube-set, place the attached cover over the end to place the mattress in transport mode.

***Overlay Base Cover***

The base cover for the mattress overlay is PVC coated nylon (a non PVC version is also available) on the underside. Four corner retention straps are incorporated, which slide under the corners of the base mattress.

A full technical description of the **Alpha Active 3** system can be found in the Service Manual, part number SER0017, available from your ArjoHuntleigh sales office.

## 2. Clinical Applications

---

**Indications** The **Alpha Active 3** system is indicated for the prevention and/or management of all categories<sup>1</sup> of pressure ulcer, when combined with an individualised, comprehensive pressure ulcer protocol: for example, repositioning, nutritional support, skin care. Selection should be based upon a holistic assessment of the patient's individual care needs.

The systems represent one aspect of a pressure ulcer management protocol; all other aspects of care should be considered by the prescribing clinician.

If existing wounds do not improve or the patient's condition changes the overall therapy regimen should be reviewed by the prescribing clinician.

The above are guidelines only and should not replace clinical judgement.

The **Alpha Active 3** mattress is designed for patients weighing up to 135 kg (300 lb).

**Contraindications** Do not use **Alpha Active 3** system for patients with unstable spinal fractures.

**Cautions** If patients have other unstable fractures, or conditions which may be complicated by a soft or moving surface, advice should be sought from an appropriate clinician before use.

While the **Alpha Active 3** system has been designed to manage patients up to the weight limits indicated above, those approaching this upper limit are likely to have additional care and mobility needs and may be better suited to a specialist bariatric system.

 *The above are guidelines only and should not replace clinical judgement or experience.*

1. NPUAP/EPUAP International Pressure Ulcer Guideline, 2009.

### 3. Installation

#### Preparing the system for use

Remove the system from the packaging. You should have the following items:

- **Alpha Active 3** pump including mains power cord and hanging brackets.
- **Alpha Active 3** mattress overlay with integral tube-set.
- Cover.

#### Installing the Mattress

##### Caution

- **Do not use the mattress overlay directly on the bed frame.**

#### Mattress Overlay

The **Alpha Active 3** mattress overlay system should be installed as follows:

1. Place the overlay on top of the base mattress, with the tube-set located near the foot end of the bed and the CPR at the head end. The cells of the mattress must be uppermost.
2. Secure the overlay to the base mattress by placing the four long straps under the corners of the base mattress.

#### To Complete the Mattress Installation

Complete the installation of the mattress overlay as follows:

1. If not already fitted, place the protective cover over the mattress. Ensure that the logo is uppermost and at the foot end of the mattress.
2. Zip the cover onto the mattress starting from the head end and taking care not to trap any material in the zip.
3. Ensure that the CPR unit is secured in it's closed position.

*The CPR must be accessible at all times.*

## **WARNING**

**Make sure that the mains power cable is positioned to avoid causing a hazard and is clear of moving bed mechanisms or other possible entrapment areas.**


***Installing the Pump*** The pump should be installed as follows:

1. Position the pump, feet down, on any convenient horizontal surface or alternatively suspend from the bed foot rail by means of the integral hanging brackets.
2. Ensure that the mattress tube-set is not "kinked" or twisted and connect it to the pump until it clicks into place. Ensure that the tube-set is securely connected to the pump.
3. Insert the mains power plug into a suitable mains power socket.

***System Operation*** The system is now ready for use. Refer to Page 9 “Mattress - Pump Operation” for day-to-day operating instructions.

## 4. Controls, Alarms and Indicators

### Controls

**On/Off Button** Pressing the On/Off button will activate the pump. The switch will illuminate when the pump is on.

#### Alarm Mute

During an alarm condition the sound of the alarm can be muted by pressing this button. The indicator will then remain on but the alarm will be muted.


#### Mode Selector

The mode selector allows the mode of operation to be selected.

Two modes of operation are available, Alternating or Static (non-alternating). The indicator next to the respective icon indicates the currently selected mode.

#### Pressure Control

This is a rotary action control used to manually set the pressure inside the mattress overlay and therefore the amount of support the patient receives.

## Alarms and Indicators

### Low Pressure Indicator

The **Low Pressure** indicator is illuminated whenever the pump detects low pressure within the mattress overlay. An audible alarm will sound unless cancelled by the **mute** button.

The indicator will extinguish once normal pressure is reached.

 *The Low Pressure alarm is inactive for the first 30 minutes of the pump being switched on.*

 *See Page 16 “Troubleshooting” for possible causes of Low Pressure.*

### Service Indicator/Pump Fault

The **Service/pump fault** indicator will illuminate and remain on if the pump has detected an internal fault. A Service Engineer should be called.

### Power Fail Indicator

The **Power Fail** indicator will illuminate when a mains power failure has been detected. An audible alarm will sound until power is resumed or the pump is switched off using the on/off button.

 *If the operation of the pump changes during use, refer to Page 16 “Troubleshooting” of this IFU before calling a service engineer or contacting your local ArjoHuntleigh sales office.*

## 5. Mattress - Pump Operation

These instructions cover day-to-day operation of the system. Other operations, such as maintenance and repair, should only be carried out by suitably qualified personnel.

 Refer to Page 7 “Controls, Alarms and Indicators” for a comprehensive description of the controls and indicators on the pump.

### WARNING

**Do not place the patient on the mattress until it is fully inflated.**

**Quick Start** Before using the **Alpha Active 3** mattress overlay make sure it has been installed correctly in accordance with Page 5 “Installation” and ensure that the CPR unit on the mattress is clicked into the closed position.

1. Connect the pump to the mains power supply using the supplied cable and switch on the pump.
2. Press the On/Off button on the side of the pump.
3. Allow approximately 30 minutes for the mattress overlay to inflate fully.
4. Place a bed sheet over the mattress and tuck in loosely. Ensure that the CPR unit is clearly visible at the head end of the bed.

**Support Setting Procedure**

It is important to follow the correct support setting procedure to ensure the patient receives adequate support while achieving maximum pressure redistribution and comfort.

1. Lie or sit the patient on the mattress.
2. Set pressure dial on pump according to patient weight. This should serve as an approximate guide only. An independent clinical determination needs to confirm that the patient is properly supported.
3. Wait 10 minutes while the pump adjusts the pressures.
4. Ensure that the patient is not ‘bottoming out’ by unfastening the cover and sliding a outstretched hand (palm up) underneath the deflated cells below the part of the body where the pressure ulcer exists, or the area at risk for a pressure ulcer.

5. If the caregiver feels less than an inch of support material, the patient has bottomed out and the support pressures should be adjusted accordingly.
6. Bottoming out should be checked at various anatomical sites and while the patient assumes various body positions.

## **Changes in Patient Position**

When a patient is in the lying or supine position, their body weight is dispersed over a large area. When in the sitting position, their body weight is concentrated within a much smaller area and therefore will require more support than in the lying position.

Therefore, when the patient changes position, it may be necessary, in order to maximise the benefit of the support surface, to make adjustments to the setting on the pressure dial.

**From Lying to Sitting** - Increase pressure control.

**From Sitting to Lying** - Decrease pressure control.

This adjustment should be in conjunction with independent clinical determination of appropriate support.

**Static** Provides a stable, non-moving support surface for instances where a active therapy surface is contraindicated e.g. to carry out nursing procedures or for patients unable to tolerate a moving surface. In **Static** mode the support surface remains constant (all cells are equally inflated). Additional nursing assessment must be undertaken in order to direct an individualised repositioning programme.

When operating the system in **Static** mode it may be necessary, where possible, to reduce the pressure setting to increase patient comfort and safety.

## **Power Fail Condition**

If a **Power Fail** condition arises disconnect the tubeset from the pump and place the attached cap over the end of the tubeset to put the mattress into transport mode. Transport mode is non-therapeutic offering support only for up to 12 hours. It is recommended that when in transport mode the patient is frequently monitored. Once power is resumed, re-connect the tubeset to the pump to continue therapy.

<b>To Disconnect the Tubeset</b>	To disconnect the tube-set at any time, depress the buttons on the top and bottom of the tubeset connector pull the tube-set connector away from the pump.
	To deflate the mattress Refer to Page 11 “To Deflate and Store the Alpha Active 3 Mattress”.
<b>Transport Mode</b>	To transport a patient using the <b>Alpha Active 3</b> mattress overlay, disconnect the tube-set from the pump and place the attached cap over the end of the tubeset to put the mattress into transport mode. This will automatically switch the mattress into transport mode.
	The patient will remain supported by the mattress for up to 12 hours.
	To resume normal operation, simply reconnect the tube-set and run the pump.

### **Caution**

**Transport mode is non-therapeutic offering support only for up to 12 hours. It is recommended that when in transport mode the patient is frequently monitored.**

## **To Deflate and Store the Alpha Active 3 Mattress**

- To deflate the mattress:**
1. Disconnect the tube-set from the pump.
  2. Activate the CPR control to deflate the mattress.

- To store the mattress** Following deflation:



1. Bring the tubeset over the mattress to lie parallel to the foot end of the mattress.
2. Roll the mattress from the foot end toward the CPR connector at the head end of the mattress.

## CPR Control

### IMPORTANT IN THE EVENT OF CARDIAC ARREST

In the event of a patient suffering cardiac arrest and CPR needing to be administered:

- To activate CPR** Located at the head end of the mattress overlay (on the same side as the tubeset) is a red strap labelled CPR. In the event of a cardiac arrest pull this from the mattress to deflate.

- To reset CPR** To re-inflate the mattress , simply replace the stopper securely into the manifold.

## **6. Decontamination**

The following processes are recommended, but should be adapted to comply with the local or national guidelines (Decontamination of Medical Devices) which may apply within the Healthcare Facility or the country of use. If you are uncertain, you should seek advice from your local Infection Control Specialist.

The **Alpha Active 3** system should be routinely decontaminated between patients and at regular intervals while in use; as is good practice for all reusable medical devices.

### **WARNING**

**Remove the electrical supply to the pump by disconnecting the mains power cord from the mains power supply before cleaning.  
Protective clothing should always be worn when carrying out decontamination procedures.**

### **Caution**

**Do not use Phenol-based solutions or abrasive compounds or pads during the decontamination process as these will damage the surface coating. Do not boil or autoclave the cover.  
Avoid immersing electrical parts in water during the cleaning process. Do not spray cleaning solutions directly onto the pump.**

**To clean** Clean all exposed surfaces and remove any organic debris by wiping with a cloth moistened with a simple (neutral) detergent and water. Dry thoroughly.

Do not allow water or cleaning solutions to collect on the surface of the pump.

**Chemical Disinfection** To protect the integrity of the cover we recommend a chlorine-releasing agent, such as sodium hypochlorite, at a strength of 1,000ppm available chlorine (this may vary from 250ppm to 10,000ppm depending on local policy and contamination status).

Wipe all cleaned surfaces with the solution, rinse and dry thoroughly.

Alcohol based disinfectants (strength 70%) may be used as an alternative.

Ensure the product is dry before storage.

If an alternative disinfectant is selected from the wide variety available we recommend that suitability for use is confirmed with the chemical supplier prior to use.

**DO NOT WRING/MANGLE, AUTOCLAVE OR USE PHENOLIC BASED SOLUTIONS.**

- Thermal Disinfection*** For information for the mattress top cover, including laundering guidelines, refer to Page 19 “Cover Specification”.

## 7. Routine Maintenance

---

### Alpha Active 3 System

- Maintenance** The equipment has been designed to be virtually maintenance-free between service periods.
- Servicing** **ArjoHuntleigh** will make available on request service manuals, component parts lists and other information necessary for **ArjoHuntleigh** trained personnel to repair the system.
- Service Period** **ArjoHuntleigh** recommend that the **Alpha Active 3** system should be serviced by an **ArjoHuntleigh** authorised service agent, after 12 months running time has elapsed.

The Service symbol illuminates to indicate that the pump is ready for a service (refer to Page 8 “Service Indicator/Pump Fault”).

### Alpha Active 3 Pump

#### **General Care, Maintenance and Inspection**

Check all electrical connections and power cord for signs of excessive wear.

In the event of the pump being subjected to abnormal treatment, e.g. immersed in water or dropped, the unit must be returned to an authorised service centre.

### Alpha Active 3 Mattress Overlay

- General Care** Remove the top cover and inspect for signs of wear or any tears.
- Check all zips are secure.
- Check integrity of all connectors, including cell to manifold connections.
- Ensure all cell fasteners are correctly connected to the mattress base sheet and are not loose or damaged.

### Serial Labels

- Pump** The serial number for the pump is on the label on the back of the pump case.
- Mattress** The mattress serial label can be found just inside the base cover above the tubeset, refer to the illustration on page Page 2 “Alpha Active 3 Mattress Overlay”.

## 8. Troubleshooting

The following table provides a troubleshooting guide for the **Alpha Active 3** system in the event of malfunction.

Indicator	Possible Cause	Remedy
<b>LOW PRESSURE</b> 	The tube-set is not connected properly. CPR not fully closed. There is a leak in the system.	Check the tube-set connector and ensure it is securely fitted to the pump. Close CPR unit. Call service engineer.
<b>POWER FAILURE</b> 	Power has been removed from the pump.	Re-apply power or switch the pump off.
<b>SERVICE</b> 	Pump has detected an internal fault.	Switch the pump off and call service engineer.

## **9. Technical Description**

<b>PUMP</b>		
Model:	Alpha Active 3	
Part Numbers:	648301	UK
	648309AU	AUSTRALIA
Supply Voltage:	230V	
Supply Frequency:	50Hz	
Power Input:	1A	
Size:	(L)280mm x (W)205mm x (H)112mm	
Weight:	2.5kg	
Case Material:	ABS Plastic	
Plug Fuse Rating:	5A to BS1362 (UK ONLY)	
Pump Fuse Rating:	2 x T1AL 250V	
Degree of protection against electric shock:	Class II Type BF	
Degree of protection against liquid ingress:	IPX0	
Mode of operation:	Continuous	
Cycle Times:	12 mins	
	Inflate - 5.5 mins	
	Crossover - 30 secs	
	Deflate - 5.5 mins	
	Crossover - 30 secs	

## PUMP SYMBOLS

	The operator must read this document (Instructions for Use) before use. Note: This symbol is blue on the product label.	O (Off)	Power Disconnects from the mains supply	I (On)	Power Connects to the mains supply
 25EA CAN/CSA-C22.2 No 60601-1 (2008)	With respect to electric shock, fire and mechanical hazards only in accordance with CAN/CSA-C22.2 No. 60601.1 (2008). MEDICAL EQUIPMENT		Double Insulated		Type BF
	Refer to this document (Instructions for Use) for a description of the product classification (3rd Edition).	SN:	Serial Number	Ref:	Model number
	Dangerous voltage		Do not dispose of in domestic refuse		Manufacturer: This symbol is accompanied by the name and the address of the manufacturer.
	Refer to this document (Instructions for Use) for a description of the product classification (2nd Edition).				

## PUMP ENVIRONMENTAL INFORMATION

Condition	Temperature Range	Relative Humidity	Atmospheric Pressure
Operating	+10°C to +40°C (+50°F to +104°F)	30% to 75% (non-condensing)	700hPa to 1060 hPa
Storage (Long Term)	+10°C to +40°C (+50°F to +104°F)	20% to 95% (non-condensing)	700 hPa to 1060 hPa
Storage (Short Term)	-20°C to +50°C (-4°F to +122°F)	20% to 95% (non-condensing)	500 hPa to 1060 hPa
 If the pump is stored in conditions outside the “Operating” ranges, allow time for its temperature to stabilise to normal, before use.			

## MATTRESS

Description	Cell Material	Base Pad Material
ALPHA ACTIVE 3 OVERLAY 90	Polyurethane	PVC Nylon Coated
ALPHA ACTIVE 3 OVERLAY 85	Polyurethane	PVC Nylon Coated
ALPHA ACTIVE 3 OVERLAY 85 NON-PVC	Polyurethane	Nylon PU Coated

## MATTRESS SIZE INFORMATION


Part No.	Description	Spare Cover	Length mm	Width mm	Height mm
648321	ALPHA ACTIVE 3 OVERLAY 90	648430	1950 (77")	900 (35")	125 (5")
648323	ALPHA ACTIVE 3 OVERLAY 85	648460	1950 (77")	850 (33")	125 (5")
648325	ALPHA ACTIVE 3 OVERLAY 85 PU	648460	1950 (77")	850 (33")	125 (5")

## COVER SPECIFICATION

Feature	Standard Cover (Dartex) <sup>®</sup>
<b>Removable Cover</b>	Yes
<b>Moisture Vapour Permeable</b>	Yes
<b>Air Permeable</b>	No
<b>Low Friction</b>	Yes
<b>Water Resistant / Repellent</b>	Yes
<b>Infection Control</b>	Material coating is Bacteriostatic, fungistatic, antimicrobial
<b>Fire Retardant</b>	BS 7175: 0,1 & 5
<b>2-Way Stretch</b>	Yes
<b>Washing Conditions</b>	MAX 95°C (203°F) for 15 mins <sup>1</sup>
<b>Drying Conditions</b>	Tumble Dry up to 130°C (266°F) or Air Dry
<b>Life Span</b>	50 Wash Cycles (minimum)
<b>Application Area</b>	Acute and Homecare

1. Check your local policy to determine the time/temperature ratio required to achieve thermal disinfection.

## CLEANING SYMBOLS

	Wash at 80°C (176°F)		Do not tumble dry above 130°C.
	Do not iron		Do Not Use Phenol-based cleaning Solutions
	Wipe surface with damp cloth	 1000ppm NaOCl NaDCC	Use solution diluted to 1000 ppm of Available Chlorine

## Guidance and manufacturer's declaration - electromagnetic emissions

The pump is intended for use in the electromagnetic environment specified below. The customer or the user of the pump should assure that it is used in such an environment.

Emissions Test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR - 11	Group 1	The pump uses RF energy only for its internal function. therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR - 11	Class A	The pump is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-2	Complies	

## **Recommended separation distances between portable and mobile RF communications equipment and the pump**

The pump is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the pump can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the pump as recommended below, according to the maximum output power of the communications equipment.

<b>Rated maximum output power of transmitter W</b>	<b>Separation distance according to frequency of transmitter m</b>		
	<b>150 kHz to 80 MHz</b> $d = 1.2\sqrt{P}$	<b>80 MHz to 800 MHz</b> $d = 1.2\sqrt{P}$	<b>800 MHz to 2.5 GHz</b> $d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	2.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

**Note:** At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

**Note:** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

## Guidance and manufacturer's declaration - electromagnetic immunity

The pump is intended for use in the electromagnetic environment specified below. The customer or the user of the pump should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 V	<p>Portable and mobile RF communications equipment should be used no closer to any part of the pump, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p><b>Recommended separation distance</b></p> <p><math>d = 1.2\sqrt{P}</math></p> <p><math>d = 1.2\sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}</math></p> <p><math>d = 2.3\sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}</math></p> <p>where <math>P</math> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and <math>d</math> is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey <sup>a</sup>, should be less than the compliance level in each frequency range.<sup>b</sup></p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5GHz	3 V/m	

**Note:** At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

**Note:** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

<sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the pump is used exceeds the applicable RF compliance level above, the pump should be observed to verify normal operation. If abnormal operation is observed, additional measures may be necessary, such as reorientating or relocating the pump.

<sup>b</sup> Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

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**CANADA**

ArjoHuntleigh Canada Inc.  
1575 South Gateway Road  
Unit "C"  
MISSISSAUGA, ON, L4W 5J1  
Tel/Tél: +1 905 238 7880  
Free: +1 800 665 4831 Institutional  
Free: +1 800 868 0441 Home Care  
Fax: +1 905 238 7881  
E-mail: info.canada@arjohuntleigh.com

**ČESKÁ REPUBLIKA**

ARJO Hospital Equipment s.r.o.  
Hlinky 118  
CZ- 603 00 BRNO  
Tel: +420 549 254 252  
Fax: +420 541 213 550

**DANMARK**

ArjoHuntleigh A/S  
Vassingerødvej 52  
DK-3540 LYNGE  
Tel: +45 49 13 84 86  
Fax: +45 49 13 84 87  
E-mail: info.dk@arjohuntleigh.com

**DEUTSCHLAND**

ArjoHuntleigh GmbH  
Peter-Sander-Strasse 10  
D-55252 MAINZ-KASTEL  
Tel: +49 (0) 6134 186 0  
Fax: +49 (0) 6134 186 160  
E-mail: info-de@arjohuntleigh.com

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C. Psimitis Co Ltd  
Dimitriou Andr. 59  
GR-16121 KAISSARIANI ATTIKIS  
Τηλ: 21 0724 36 68  
Φάξ: 21 0721 55 53

**ESPAÑA**

ArjoHuntleigh Ibérica S.L.  
Ctra. de Rubí, 88 1<sup>a</sup> planta - A1  
08173 Sant Cugat del Vallés  
ES- BARCELONA 08173  
Tel: +34 93 583 11 20  
Fax: +34 93 583 11 22  
E-mail: info.es@arjohuntleigh.com

**FAR EAST**

ARJO Far East Limited  
Unit 3A, 4/F., Block B Hoi Luen  
Industrial Centre  
55 Hoi Yuen Road,  
Kwun Tong, Kowloon  
HONG KONG  
Tel: +852 2508 9553  
Fax: +852 2508 1416

**FRANCE**

ArjoHuntleigh SAS  
2 Avenue Alcide de Gasperi  
CS 70133  
59436 RONCQ CEDEX  
Tél: +33 (0) 3 20 28 13 13  
Fax: +33 (0) 3 20 28 13 14  
E-mail: info.france@arjohuntleigh.com

**INTERNATIONAL**

ArjoHuntleigh International Ltd.  
ArjoHuntleigh House  
Houghton Hall Park  
Houghton Regis  
UK-DUNSTABLE LU5 5XF  
Tel: +44 (0) 1582 745 800  
Fax: +44 (0) 1582 745 866  
E-mail:  
international@ArjoHuntleigh.com

**ITALIA**

ArjoHuntleigh S.p.A.  
Via di Tor Vergata 432  
00133 ROMA - ITALIA  
Tel: +39 (0) 6 87426211  
Fax: +39 (0) 6 87426222  
E-mail: Italy.promo@arjohuntleigh.com

**MURANAKA MEDICAL INSTRUMENTS CO. LTD**

2-8-2 Ayumino  
Izumi-shi  
OSAKA, 594-1157  
Tel: +81 (0)725-53-5546  
www.muranaka.co.jp

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**GETINGE GROUP** is a leading global provider of products and systems that contribute to quality enhancement and cost efficiency within healthcare and life sciences. We operate under the three brands of **ArjoHuntleigh**, **GETINGE** and **MAQUET**. **ArjoHuntleigh** focuses on patient mobility and wound management solutions. **GETINGE** provides solutions for infection control within healthcare and contamination prevention within life sciences. **MAQUET** specializes in solutions, therapies and products for surgical interventions and intensive care.

648934ML\_01:12/2012

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ArjoHuntleigh Inc.
ArjoHuntleigh AB