
PROVOX[®] HME System[™]

Patient's and Physician's Manual

U.S. federal law restricts this device to sale by or
on the order of a physician.

Read
instructions
before use

For
single
use only

ATOS
MEDICAL

Caution

United States Federal law restricts devices identified herein with the symbol "Rx" to sale, distribution and use by or on the order of a licensed practitioner trained in post-laryngectomy vocal and pulmonary rehabilitation. The availability of this product without prescription outside the United States may vary from country to country.

Disclaimer

Atos Medical makes no warranty, either expressed or implied, as to the lifetime of the product delivered to the purchaser hereunder, which may vary with individual use and biological conditions. Further Atos Medical makes no warranty of merchantability or fitness of the product for any particular purpose.

Patents and Copyrights

Remove and Skin-Prep are registered trademarks of Smith & Nephew United Inc., USA.

Provox[®] is a registered trademark owned by Atos Medical AB, Sweden.

Provox[®] HME has US patents (5.738.095, 5.042.468), other patents and patents pending.

Provox[®] LaryButton[™] and Provovox[®] LaryClip[™] are trademarks owned by Atos Medical AB.

CONTENTS

ENGLISH.....	4
DEUTSCH	11
NEDERLANDS.....	19
FRANÇAIS.....	26
ITALIANO.....	34
ESPAÑOL	42
PORTUGUÊS.....	50
SVENSKA	58
DANSK.....	65
NORSK.....	72
SUOMI.....	79
LATVIEŠU.....	86
ΕΛΛΗΝΙΚΑ.....	93
TÜRKÇE.....	102
ORDER INFORMATION.....	110

Descriptive Information

Indications for Use

The Provox® HME Cassette is a specialized device for laryngectomized patients. It is a heat and moisture exchanger (HME) that heats and humidifies inhaled air by retaining heat and moist from exhaled air in the device. It partially restores lost breathing resistance. For patients with a voice prosthesis or surgical speech fistula it may also facilitate voicing.

Description of the Device

The Provox® HME Cassettes are part of the Provox® HME System, which consists of HME Cassettes and Attachment Devices.

Fig. 1

Provox® HME Cassette

The HME Cassette is a single use device that features calcium chloride treated foam sponge in a plastic housing. This housing has a top lid that can be pushed down with a finger during speech to close the Cassette airtight. *Fig. 2*

After releasing the finger, the top lid will get back automatically to its rest position (spring mechanism).

The Provox® HME Cassettes are available in two versions: The **Normal** HME Cassette, which should be worn day and night under normal physical effort and the **HiFlow** HME Cassette with a lower breathing resistance intended for use during physical activities. The latter can also be used in a two-step approach in order to get adapted to a more normal breathing resistance. In that case, start with the HiFlow HME Cassette for the first one or two weeks after which you should switch to the Normal HME Cassette.

Technical Data:

<i>Pressure drop at 30 l/min:</i>	Normal: 89 Pa (0.908 cmH ₂ O) HiFlow: 66 Pa (0,673 cmH ₂ O)
<i>Pressure drop at 60 l/min:</i>	Normal: 207 Pa (2.11 cmH ₂ O) HiFlow: 172 Pa (1.75 cmH ₂ O)
<i>Moisture loss:</i>	Normal: 23.7 mg/l, HiFlow: 25.4 mg/l

(according to ISO 9360-2; 2001)

Expectations of the device

The Provox® HME Cassette is a device for pulmonary rehabilitation. When continuously used, the pulmonary function is likely to improve and respiratory problems, e.g. coughing and mucus production, to decrease. Clinical studies have shown also positive psychological effects.

Especially in the beginning, the slightly increased breathing resistance may be experienced as discomfort. Starting with HiFlow Cassettes may therefore be advisable (see above ‘Description of the Device’).

During the first weeks of use the mucus production may appear increased due to thinning of the mucus with retained water.

PRECAUTIONS:

Always test the spring mechanism of the Provox® HME Cassette prior to use. The top lid should immediately open after releasing the finger.

Do not disassemble the HME Cassette since this may interfere with its proper function.

Do not reuse the HME Cassette or clean it with water.

The function of the heat and moisture exchanger will be substantially reduced if the calcium chloride is washed out of the foam. Also, the risk of potential infections may increase with the time of use, due to bacterial colonization in the foam.

Checkout procedures

Check the proper function of the spring mechanism in the HME Cassette prior to use.

Instructions for use

Insert the HME Cassette. The HME Cassette can easily be inserted and removed manually from the adhesive as needed. *Fig. 3*

Closing the tracheostoma for voicing

In order to obtain an airtight seal for speech with a voice prosthesis or surgical speech fistula, press the top lid of the HME Cassette down with a finger. *Fig. 2*

You can also close the stoma, if a pressure increase of the air in the lungs is needed (e.g. coughing).

Coughing

The HME Cassette can be removed during or after coughing if the stoma needs to be cleaned from mucus. Generally, the Adhesive base plate can stay in place during cleaning of the stoma or Cassette replacement.

Replacement of the HME Cassette

The Cassette can be replaced as needed. To ensure its proper function, the same Cassette should not be used for more than 24 hours.

Storage

Store in a clean, dry and dark area at room temperature (15°C–25°C/60F–80F).

Re-order information

For order codes, please see the table on the last page.

Instructions on related, additional devices

The HME Cassette can also be used in combination with the Provox® LaryTube System (Laryngectomy tube) *Fig. 1E* or with the Provox® LaryButton *Fig. 1D* with/without using Provox® Lary-Clip, or HME Cassette Adaptor for tracheostomy tubes with ISO 15 mm connection (the latter is not for sale in the USA).

Other information

Information about clinical studies, local representatives and further information are available from our customer service department.

Date of Printing

See version number on the back cover of this Manual.

Attachment devices

The Provox Micron HME and other Provox HME Cassettes are intended to be used with the following attachment devices and accessories in Provox HME System.

For more information read Instructions for use for each of the products.

Provox® Adhesive

Provox® Adhesive base plates are single use devices that are intended to hold Provox® Micron HME and other Provox® HME Cassettes in front of the tracheostoma and to provide an airtight seal. The self-adhesive tape discs are supplied on a peel-away paper liner and come in two different shapes, i.e. oval (anatomical) and round, and three different types of glue: **OptiDerm**, **FlexiDerm** and **Regular**.

A special version, **Provox® XtraBase**, is a Flexi-Derm type adhesive that has a conically shaped base. It is designed to facilitate and improve skin attachment for deep stomas and / or if an automatic tracheostoma valve Provox® FreeHands HME is used.

OptiDerm (*fig 1a*) is intended for patients with a sensitive skin or if the skin is not yet healed properly, e.g. after surgery and/or irradiation. It absorbs moisture, is gentle to the skin and needs less frequent replacement. It can be applied gently on the 1st postoperative day. Average durability 36 hours.

FlexiDerm (*fig 1b*) is very flexible and has the strongest adhesive properties. It is especially suitable for patients with a deep and/or irregular stoma or who prefer a very soft and flexible adhesive. Due to the type of glue used, a suitable skin cleanser is recommended to clean the skin of any residual glue. Average durability 48 hours.

Regular (*fig 1c*) is a transparent, perforated adhesive. It is less flexible and sticky than FlexiDerm, but leaves less residual glue and is therefore easier and more -comfortable to use. This adhesive is mainly intended for persons with a regular stoma and skin. Average durability 24 hours.

WARNINGS

- The adhesive glue may irritate the skin. Stop using the adhesive if skin irritation develops and consult your physician.
- The Provox® HME Adhesives should not be used during radiotherapy of the skin around the stoma or for 2 weeks after radiotherapy.
- Do not use the OptiDerm adhesive as attachment for the Provox® ShowerAid. The glue may dissolve in contact with water and the Adhesive may become loose, allowing water to enter the trachea (*fig 4*).

Operating Information

Set-up Instructions:

1. Clean the stoma and the skin around it with Provox® Cleaning Towel and/ or soap and water, and then dry the area carefully.
2. You may apply skin preparation agents, like Skin-Prep™, or, if needed, Provox® Silicone Glue, according to their instruction for use.
3. If OptiDerm is used: Please carefully pre-warm the Adhesive base plate between your hands for improved attachment.
4. Remove the protective paper from the adhesive. Stretch the skin gently away from the stoma, and then apply the adhesive, in order to get optimal attachment to the skin, even in skin folds.
5. Massage the adhesive carefully gently for approximately one minute in order to improve skin attachment.

Replacement of the Adhesive base plate

The Adhesive base plate can stay in place as long as it is glued properly to the skin. To facilitate removal of partially attached Adhesives, the use of adhesive remover, e.g. Remove™ may be helpful. After Remove™ is used the skin should always be cleaned with a cleaning towel or soap and water. Read instructions. Before applying a new Adhesive base plate, the skin should be carefully cleaned and dried.

CAUTION: When cleaning the skin from residual glue, prevent particles/fluids from entering the stoma. In case of using OptiDerm in the postoperative period or sensitive skin, the Adhesive should only be replaced when completely loose.

Storage

Store in a clean, dry and dark area at room temperature (15°C–25°C/60F–80F).

Provox® LaryButton™

Soft silicone laryngectomy stoma buttons (fig 1d).

- with or without using Provox® LaryClip or Provox® TubeHolder.
Read separate instructions.

Provox® LaryTube™

Soft silicone laryngectomy tubes (fig 1e).

- with or without using Provox® LaryClip or Provox® TubeHolder.
Read separate instructions.

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Atos Medical AB, P.O. Box 183 SE 242 22 Hörby, Sweden
© Atos Medical AB, 200712B, REF 90719