DorsoFX Rückenorthese

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Back Brace - Technical description | Orthèse dorsale - Description technique
Órtesis para espalda - Descripción técnica | Ortesi dorsale - Descrizione tecnica

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Information for the patient

Functionally and effect:

- Immobilization of thoracic spine/lumbar spine
- Prevents detrimental movements
- Individual adaptation and subsequent adjustment possible
- Adjustment possibilities (Velcro closures and multivariable adjustment mechanism)
- Pain reduction at the affected spine segment
- Relief of stress on the affected spine segment
- Reduced risk of falling thanks to improved body statics

Instructions for use for the patient:


Fig. 1

1. First, put on the provided support. The support is closed at the front. If you put it on correctly, the support is wrapped around your body at waist height (**see Fig. 1**).


Fig. 2

2. The **BORT DorsoFX back brace** which was adapted to your needs by your orthopedic technician, is put on from the back to the front. Ensure that the upper braces are positioned just below the armpit and the straps are not twisted (**see Fig. 2**).

Tip: Wear a tight-fitting shirt or T-shirt under your orthopedic aid. This increases wearing comfort.


Fig. 3

3. First, close the lower tensioning strap at the hip to secure the corset and prevent it from slipping out of place.


Fig. 4

4. Close the chest strap and check again for correct positioning (braces below the armpits). (**see Fig. 4**).

Tip: Ensure that no creases form between your shirt and the corset to prevent pressure sores.


Fig. 5

5. Close the Velcro closure in the abdominal area and re-tension the bottom Velcro closure in the hip area (**see Fig. 5**). Please ensure that the support is crease-free and comfortable and can be smoothly attached to the corset.

Tip: Make sure that neither the support nor the corset put pressure on the groin area when sitting! If there is pressure, put your orthopedic aid on again.

Special notes:

Please ensure that the correct fit and function of your BORT DorsoFX back brace are regularly checked by your orthopedic technician. This is particularly important at the beginning of your treatment.

The correct fit of your back brace is decisive for the success of your treatment.

It is very important that you practice with your physician/physiotherapist/orthopedic technician how to put on and take off your support corset to gain the required routine in using your orthopedic aid.

Application period:

The application period of your BORT DorsoFX back brace should be determined by your physician in a treatment schedule taking into account the indication and your personal situation.

Please observe the daily minimum period of application set by your physician. This is crucial to support the healing process.

For fresh fractures in particular, a long period of application can be expected.

Indications:

Conservative treatment of fresh, traumatic fractures of the thoracic spine (Th 8-12) and lumbar spine (L 1-5), conservative treatment of osteoporotic compression fractures of the thoracic spine (T 8 – 12) and lumbar spine (L 1-5), immobilization in the case of acute and/or chronic spinal diseases (in the above mentioned section), if immobilization is indicated according to the physician's instructions.

Contraindications:

Peripheral arterial disease (PAD), lymphatic drainage disorders and indistinct soft tissue swellings remote from the applied aid, sensory loss and blood flow disorders in the affected body part, skin disorders in the treated part of the body, acute disk prolaps with neurological symptoms, unstable fractures requiring surgery, tumors requiring surgery, instabilities, insufficient compliance.

Important note:

Not reusable – this orthopedic aid is intended for use by one patient only. If the back brace is to be worn at night, make sure it does not hamper your blood circulation. If you feel any numbness, loosen the corset or take it off if necessary. Contact your physician if symptoms persist. Please do not use creams or ointments as they may damage the material.

The product contains latex and may cause allergic reactions.

Information for the technician:

Fitting of the BORT DorsoFX back brace:

The fitting of the BORT DorsoFX back brace may only be performed by trained professionals!

Ideally, you attend to the patient in interdisciplinary cooperation with the prescribing physician and the competent physiotherapist.

If you attend to a patient who is in a hospital at the time, it is recommended to inform the nursing staff about the BORT DorsoFX back brace.

Important note:

If the measured dimensions are directly at the transition to the next size, choose the larger size. The material is easy to process and allows for easy modification in the workshop.

Most modifications will concern the back plate. You can shorten the upper end (cranial) by 3 cm and the lower end (sagittal) by 2 cm.

Fitting of the BORT DorsoFX back brace on the patient

According to the indication and following consultation with the prescribing physician, the corset is put on while the patient is **standing** or **lying down**!


Fig. 1

1. Fit the support onto the patient. You have to wrap it around the patient's body at waist level and close it. On the back of the support, you find a freely positionable, double sided Velcro closure. It connects the support to the back brace (Fig. 1).


Fig. 2

2. Put the corset on the patient from the back to the front (dorsal-ventral). **Ensure that the position of the reclination bars does not put any pressure on the armpits of the patient (Fig. 2).**


Fig. 3

3. First, close the Velcro closure of the lower pelvis buckle (Fig.3). This ensures that the corset does not slide down. Thanks to the Velcro closures in the ventral area of the back brace, the connection between support and corset provides additional support.


Fig. 4

4. Next, close the Velcro closure of the chest buckle (reclination bar). **(Fig. 4).**


Fig. 5

5. Close the upper closing strap of the pelvic buckle **(Fig. 5).**
Please check again for correct positioning of all closures of the corset and of the support.
 If the corset was put on for the first time, ask the patient to stand up (Caution: keep diagnosis in mind).


Fig. 6

6. Check the corset for correct fit. With the help of the provided screw driver, you can continuously adjust the height of the thoracic and pelvic buckles individually **(Fig. 6).**
 Depending on the mobility of the patient, this step may also be performed while the patient is wearing the corset. If the patient has very limited mobility, ask them to lie down and adjust the buckles while the back brace is off. Put the corset on the patient and check again for correct fit of the orthopedic aid.


Fig. 7

7. Finally, adjust the closures. To do this, open the crocodile clip at the end of each strap and shorten the strap using conventional scissors so that the strap is aligned to the crocodile clip approx. 2 cm before the end **(Fig. 7).**

Tip: We recommend the patient wears a long, tight-fitting brace shirt or T-shirt underneath the orthopedic aid. Ensure that no creases form underneath the corset since this may cause pressure sores.

Important note:

In the case of any problems with the fit, you as a technician can modify the DorsoFX back brace at any time. The pelvic and thoracic buckles as well as the back plate are made of thermoplastically deformable polypropylene and can therefore be reworked using a hot air gun. This step can also be performed directly in the hospital.


Fig. 8

To do this, remove the covers attached by Velcro strips to enable smooth reworking using a hot air gun (**Fig. 8**).

Immediately after the material has cooled down, the covers can be put back to their original positions and reattached.

Material labeling (according to TextilKennzG [German Textile Labeling Act]):

Covers (chest-hip):	52 % PA microfibre, 31 % PA, 17 % spandex
Cover of back plate:	70 % PA Velcro type fabric, 30 % PUR
Support:	54 % Polyamide, 36 % monofil, 10 % spandex

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BORT GmbH | Postfach 1330 | D-71367 Weinstadt | www.bort.com