**Bedienungsanleitung · Working instructions
Mode d'emploi · Instrucciones de trabajo
Istruzioni d'uso**

deutsch

english

français

español

italiano

Drufosmart

Basis Druckformgerät 2,5 bar

Basic pressure thermoforming unit 36 psi

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Contents (GB)

1.	Description of unit	23
2.	Technical data	23
3.	Declaration of conformity	24
4.	Safety advice	25
5.	Initiation – scope of delivery	26
6.	Functional components	27
7.	Operation	28
7.1	Adjustment heating time	28
7.2	Start of programme	29
7.3	Thermoforming	30
7.4	Opening of the pressure chamber	31
8.	Deviation to the rule	32
8.1	Modification of given times	32
8.2	Plasticizing of blanks with unknown heating times	32
8.3	Abort / termination of programme	32
9.	Additional optional functions	33
9.1	Fabrication of cast form parts	33
10.	Exchange of fuses	34
11.	General advice	35
12.	Care and maintenance	35
13.	Guarantee	36
14.	Trouble shooting	37
15.	Spare part list	39
16.	Explosion sketch	91
17.	Operating diagram	92
18.	Pneumatic plan	93

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1. Description of unit

The Drufosmart is a basic pressure thermoforming unit as an alternative to common vacuum thermoforming units for daily use in the dental laboratory and dental surgery.

Simplest handling as well as continuous digital and acoustic control of all working steps ensure problem-free working, even with poor precognition.

A forming pressure of 2.5 bars in connection with a vertically forming blank guarantee for very precise results.

2. Technical data

W x H x D	320 x 420 x 290 mm
Weight	15 kg
Supply voltage	230 V/50-60 Hz, 115 V/50-60 Hz
Max. power consumption	315 W
Fuse of unit	T 3.15 A
Operating pressure	2.5 bars
Min. operating pressure	2 bars
Max. pressure in pipe	10 bars

3. Declaration of Conformity

According to the guideline machines (2006/42/EC), annex II A, according to EC-guideline for electromagnetic compatibility EMC (2004/108/EC) and guideline for electrical appliances for low voltage (2006/95/EC).

Herewith we declare that the machine described below is in its conception and design and in the shape delivered by us in accordance with the fundamental requirement for safety and health as prescribed in the applicable EC-guidelines. In the case of any change or modification of the machine not authorised by us this declaration becomes invalid.

Designation of machine: Drufosmart
D-3200/D-3200-A/D-32001/D-32001-A

Type of machine: Pressure thermoforming unit

Applicable EC-guidelines:
2006/42/EC EC-guideline machines
2006/95/EC EC-guideline for electrical appliances for low voltage
2004/108/EC EC-guideline for electromagnetic compatibility EMC

Applied harmonised standards:

- DIN EN ISO 12100-1:2004-04 Safety of machinery; basic concepts, general principles for design; part 1: basic methodology (ISO 12100-1:2003)
- DIN EN ISO 12100-1/A1 2009-10
- DIN EN ISO 12100-2:2004-04 Safety of machinery; basic concepts, general principles for design; part 2: technical principles and specifications (ISO 12100-2:2003)
- DIN EN ISO 12100-2/A1 2009-10
- EN 13732-1:2006 Ergonomics of the thermal environment – Methods for the assessment of human responses to contact with surfaces; part 1: Hot surfaces (ISO 13732-1:2006)
- EN 60204-1:2006 Safety of machinery: electric equipment of machinery; part 1: General requirements (IEC 60204-1:2005)
- EN 61010-1:2001 Safety requirements for electrical equipment for measurement, control and laboratory use; part 1: general requirements
- EN 61326-1:2006 Electrical equipment for measurement, control and laboratory use – EMC requirements; part 1: General requirements (IEC 61326-1:2005); emission according to living area, ranch of commercial and small enterprises; interference resistance according to industrial area

Dreve Dentamid GmbH
Max-Planck-Straße 31 · 59423 Unna/Germany
Tel.: +49 2303 8807-40 · Fax: +49 2303 8807-55
E-mail: info@dreve.de · www.dreve.com

4. Safety advice

Attention! Read the following notes before connecting and putting the unit into service. The operating service and the function of the unit can only be guaranteed for, if the general safety and accident prevention instructions of the legislation as well as the safety instructions of the working instructions are observed.

1. This unit must only be used according to the present working instructions. We cannot be held liable for damages which have been caused by improper treatment resp. faulty manners of operation.
2. Put the unit on a stable and even basis which is suitable for the operation weight of approx. 16 kg of the unit.
3. In order to avoid the intrusion of water into the unit (e. g. splash water) the unit is to be placed into a dry surrounding.
4. Do not store any flammable substances in the close surrounding of the unit.
5. The indicated voltage on the type plate has to meet the voltage of the power source.
6. The unit must only be connected with a socket with protective conductor. Never touch the plug with wet hands.
7. Authorised operators: The operating company of the machine has to hold the working instructions at the operators' disposal and has to make sure that authorised operators have read

and understood the manual. Only in that case the operators may operate the unit.

8. The unit has to be examined before operation in view of proper condition and operating safety. If the unit is not in proper condition, it must not be used and has to be marked accordingly.
9. The pressure of the compressed air line must not exceed 10 bars.
10. Do not intrude any objects into the unit.
11. Do not lock-up the operation elements of the unit.
12. Do not block the downwards-moving piston.
13. Do not touch the downwards-moving piston.
14. The swivelling housing of the radiator reaches high temperatures and must therefore not be touched – not even for short moments.
15. If the unit is not used, please disconnect the mains.

16. Plates and labels have to be maintained in a well legible condition and must not be removed.

17. Before cleaning and maintenance of the unit or the exchange of parts the plug has to be disconnected.

18. The opening of the unit and repairs are only to be carried out by accordingly trained experts.

19. Only those appliances and spare parts may be used which are released by the manufacturer. We are not responsible for damages caused by the use of foreign parts.

5. Initiation

Before initiation make sure that the voltage indicated on the type plate meets the voltage of the power source.

Connect the unit to an air compressor. For that purpose push the pressure hose that is included in delivery onto the hose clip (9) of the unit and fix it by means of a hose clamp. Assemble the fast coupling at the other end of the tube and connect it with the pressure line.

The pressure of the compressed air line must not exceed 10 bars.

Put the power cord into the female connector (33) and create a connection with the Euro/USA plug.

20. Arbitrary rebuilding or modifications are not permitted due to safety reasons.

21. The required operation and maintenance conditions mentioned in these working instructions are to be complied with imperatively. When working with the Drufosmart the general accident prevention regulations have to be observed.

Important: Before the first initiation please read the working instructions thoroughly. In case of comprehension problems contact your responsible dealer or Dreve Dentamid GmbH, if necessary.

Scope of delivery


- 1 Drufosmart
- 1 drawer for granules
- 1 insert bowl and nub plate
- 1 tension ring with exhaust valve
- 1 base range of blanks
- 1 power cord
- 1 air connection hose
- incl. hose clamps and fast coupling
- 1 working instructions in 5 languages
- 1 information folder thermoforming foils

6. Functional components

No. Description

- 5 Heating handle
- 6 Nub plate and insert bowl
- 7 Drawer for granules
- 9 Compressed air supply
- 10 Handle plate reception
- 11 Plate reception
- 12 Ventilation valve
- 13 Tension ring
- 20 Heating
- 27 Display
- 33 Mains supply/fuse carrier
- 34 Activation lever
- 43 Main switch
- 44 Turnable pressure controller

7. Operation

Turn the unit on by means of the main switch (43) at the back of the unit. By this the display (27) at the front of the

unit will be activated and shows 0:00. An acoustic signal shows that the unit is ready for operation.

7.1 Adjustment heating time

By means of the turnable pressure controller (44) every time between 0:05 min. and 9:55 min. can be adjusted by turning the controller clockwise. Adjustment is effected in 5 sec. steps.

Please refer to the enclosed information folder in order to find out the optimum times for Dreve thermoforming foils.

7.2 Start of programme

Before starting the adjusted heating time put a corresponding blank into the plate reception (11) and fix it by means of the tension ring (13).

The model can either be put onto the nub plate (6) or fixed by taking out the nub plate and embedding it into granules directly in the insert bowl.

Info

For an optimum adaptation of the blank towards the model always place the models in a way that the anterior teeth are directed towards the centre of the blank. When using hard blanks it is recommended to embed the models into granules. Models that are placed onto the nub plate should not exceed a max. height of 23 mm.

Swivel the plate reception inwards at the handle (10) just up to the mechanical interlock in the centre of the unit (magnetic interlocking). Also swivel in the heating element (20) at the cor-

responding handle (5) just up to the mechanical interlock in the centre of the unit. Close the ventilation valve (12) of the tension ring (13). In order to start the heating process press the turnable pressure controller (44). The display shows a countdown towards 0:00 of the adjusted time.

Info

For safety reasons the heating is equipped with a final switch and only works in swivelled in position after having started a programme. If the heating carrier (20) is swivelled outwards at the handle (5) after having started a programme, the heating process will be interrupted automatically. The indication in the display will stop and show the remaining time. By swivelling the heating carrier (20) back in the programme will be re-activated.

Attention! The swivelling housing of the heating element reaches high temperatures. Therefore it must not be touched – even not for short moments.

After the adjusted time has elapsed, there will be an acoustic signal, the indicator in the display will blink for 4 seconds and the cooling time will be counted up again automatically.

7.3 Thermoforming

After the heating time has elapsed there will be an acoustic signal; the heating time in the display will be 0:00. The heating will be switched off automatically. Pull the heating carrier (20) with the right hand at the handle (5) completely out of the centre of the unit. Then push the activation lever (34) with the left hand downwards as far as it will go.

By pulling the heating handle (5) again for approx. 3 seconds as far as it will go pressurising of the piston will be carried out automatically. This is confirmed by an acoustic signal.

Info

For safety reasons the pressure chamber can only be pressurised, if it has been moved downwards completely manually. Always make sure that no objects are in the contact area between piston and insert bowl while the piston is closing! In order to protect the operating personnel securely against bruises of the fingers, pressurising can only be effected, if the unit is operated with two hands at the same time (two-hand control).

The forming pressure will be verified continuously and regulated, if necessary. The display shows the duration since turning off the heating (cooling time). The maximum time that can be shown is 9:55 min. Please refer to the enclosed information folder in order to find out the optimum cooling times for the used Dreve thermoforming foil.

7.4 Opening of the pressure chamber

In order to attain a complete cool down of the thermoforming blank the valve at the tension ring (12) is to be opened slightly after 2/3 of the adjusted cooling time. Through this the warm air inside the piston will be released and replaced by cool air due to the automatic pressure regulation. After the cooling time has elapsed pull the heating handle (5) as far as it will go and maintain it in this position (approx. 3 seconds) until an acoustic signal can be heard. Release the heating handle and the pressure will be relieved automatically.

Then swivel the pressure chamber into its original position by leading back the activation lever (34).

Swivel out the plate reception (11) at the corresponding handle (10) and remove tension ring and thermoformed result.

8. Deviation to the rule

8.1 Modification of given times

The specified times may be prolonged or shortened individually, if required.

Info

The performance of the heating might decrease slightly after 1000 operating hours. By prolonging the heating times by 5-10 seconds opti-

mum thermoforming results can be attained again. The given times are average times and designed for general use. In case of extreme undercut areas or very high models a prolongation of the heating time can optimise the results.

8.2 Plasticizing of blanks with unknown heating times

Create a heating programme corresponding to the foil thickness in question with a sufficient prolongation of the heating time (see point 7.1). Start the programme (see point 7.2). You may test the plasticization degree of the blank by means of an edgeless metal instrument (e. g. wax knife).

Info

The complete deformability of the blank can only be tested at the edges, as that is where warmth will leak through the metal con-

tacts and the necessary flexibility will be attained at last!

After the blank has attained the necessary flexibility initiate the thermoforming process (see point 7.3). The cooling time can be terminated at any time by pushing the heating handle (5). However, a sufficient cooling down time is to be observed, as warm foils may deform during removing them from the model and therefore lead to bad resp. useless thermoforming results.

8.3 Abort/termination of programme


It is possible to abort/terminate all actions of the Drufosmart by pushing

the turnable pressure controller (44) for at least 3 seconds.

9 Optional additional functions

9.1 Fabrication of cast form parts

With the K-7 insert it is possible to thermoform up to 7 resin bases for the crown and bridge technique at the same time. Fix the single dies by means of a flexible blocking out material (e. g. Fillin) in the holes of the K-7 insert.


All further working steps as described from point 7.1 on.

10. Exchange of fuses

Disconnect the unit from mains supply. The fuse carrier is situated at the rear of the unit, integrated in the mains supply.


In order to control or exchange the fuses push the push catch of the carrier, keep it pressed and remove the carrier from the unit.

Pull the defective fuse out of the clamp carrier applying slight traction and replace it by a new, equivalent one.

Important: Only use the fuse strength indicated on the rear.

After exchanging the fuse insert the fuse carrier into the housing applying slight pressure until the fixation latches noticeably.

11. General advice

Repairs on the Drufosmart units may only be executed by experts. Disconnect the unit from mains supply at any rate. Dreve company is only responsible for the safety, reliability and performance of the unit, if:

- the unit is exclusively used for the corresponding field of application.
- operation is only carried out by accordingly trained personnel.

12. Care and maintenance

Generally disconnect the unit from mains supply before executing care and maintenance.

Ideally cleaning is to be effected dryly with a soft piece of cloth, if necessary

- the unit is operated in accordance with these working instructions.
- extensions, new adjustments, changes or repairs are executed by persons authorised by the company.

As we are currently developing our products, we reserve the right for technical changes.

with a slightly humid sponge and a mild cleaning liquid. Water and cleaning liquids must not get into the unit.

The O-rings should be greased regularly with silicone grease.

13. Guarantee

The guarantee time amounts to two years after the date of delivery. Guarantee is only accepted in case of material or fabrication defects. The guarantee obligation does not refer to natural wear and tear or to damage occurring after the transfer of risk due to improper or unsuitable handling, excessive stressing, unsuitable operating material, and chemical, electrochemical or electrical influences of a nature not provided for according to the contract. On improper alterations, repairs or maintenance or breach of seals carried out by the purchaser or third parties our warranty obligations and liability for consequences caused thereby are rendered void.

Liability

The right of the purchaser to bring actions arising out of defects shall in all cases be in lapse one year after the date of delivery. The warranty obligation does not refer to natural wear and tear or to damage occurring after the transfer of risk due to improper or unsuitable handling, excessive stressing, unsuitable operating material, and chemical, electrochemical or electrical influences of a nature not provided for according to the contract. On improper alterations, repairs or maintenance or breach of seals carried out by the purchaser or third parties our warranty obligations and liability for consequences caused thereby are rendered void.

14. Trouble shooting

Problem	Cause	Aid
Unit does not work	Power line is not correct	Control power supply
	Main breaker is off	Switch on main breaker
	Miniature fuse is defect	Check miniature fuse and exchange it, if necessary
Main breaker on, electronics show no function	Electronics are out of order	Electronics have to be exchanged by an expert
Piston cannot move downwards	Heating is not swivelled outwards	Swivel heating outwards
	Minimum pressure of 2 bars has not been reached	Hold escapement for longer, see working instructions
	Electronics are out of order	Electronics have to be checked by an expert
Piston has moved downwards, unit does not take over self-maintenance	Unit is not connected to compressed air	Connect unit to compressed air
	Wrong operation of push button	Pull push button for 3 sec.
	Contact pressure has not been relieved	Relieve contact pressure, see working instructions
Piston does not move upwards	Final switch is defect	Final switch has to be exchanged by an expert
	Magnet valve or electronics are defect	Unit has to be controlled by an expert

Problem	Cause	Aid
Electronics can be set, but radiator does not heat	Heating element is not swivelled-in completely Micro-switch defect	Swivel the heating element in correctly Micro-switch has to be exchanged by an expert
Heating time elapses, but radiator does not heat	Quartz radiator defect Temperature monitor was activated	Exchange quartz radiator Exchange temperature monitor
Heating time elapsed, blank is not plasticized sufficiently	Lifetime of quartz radiator exceeded Heating time was too short	Exchange quartz radiator Prolong heating time by 5-10 sec.
Unit does not reach the pre-adjusted pressure	Line pressure too low	Increase line pressure
Unit makes fizzy sounds during thermoforming	Ventilation valve not closed Tension ring has been placed wrongly into plate reception Connection hose leaky between valve and cylinder	Close ventilation valve Put tension ring into correct position Unit has to be controlled by an expert, connection hose has to be exchanged, if necessary

15. Spare part list

No.	Description	Item-no.
5	Heating handle	D-6505519
6	Insert bowl with nub plate	55591/50592
7	Drawer for granules	5565
9	Pressure reduction valve	55637
10	Handle plate reception	D-6505519
11	Plate reception	D-6505512
12	Ventilation valve	50761
13	Tension ring	50038
17	Housing red silver	55500rot 55500silber
19	Micro-switch	55790
20	Heating housing	D-6505513
21	Quartz-infrared-radiator	D-6505542
22	Temperature monitor	5542
23	Pillar heating element	55773
25	Magnetic valve	52021-24
26	Cylinder	52020
27	Display foil, adhesive	55820
28	Limit plate reception	xxxxx
29	Support plate reception	D-6505514
30	Pillar plate reception	55772
31	Bottom plate	5511
33	Mains supply/fuse carrier	51208
34	Activation lever	D-65885
35	Connector lever/axis	xxxxx
36	Bushings	55708
37	Final switch activation lever	5521
38	Final switch heating	50037
39	Recuperating spring	55861/55862
40	Lever axis	55847
43	Main switch	50841
without drawing no.		
	Power cord Europe	51021
	Power cord England	51213+51021
	Power cord America	51207+51163
	Electronics incl. display, key pad, connection cable	5540
	Fuse F3.15 AT	51321
	O-rings 112 x 3 mm	51077
	Cable harness complete	55520

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Dreve Dentamid GmbH · Max-Planck-Straße 31 · 59423 Unna/Germany