*TeKne Dental*

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**January 2009**

Gennaio 2009

**TN250 ELECTRONIC BOARD**  
SCHEMA ELETTRONICA TN250

*TeKne Dental* s.r.l.

Via del Pescinale 77 - 50041 Calenzano (FI) - ITALY  
info@teknedental.com      www.teknedental.com

**DESCRIPTION**

The sophisticated TN250 electronic board (REF 386.00) allows to control operating power of any piezo-electric scaler working at  $28 \div 31$  KHz frequency.

The board automatically adjusts operating frequency of vibrations to the real resonance load, which changes with the quantity of water delivered and the type of insert used.

Once wired to the TN250 electronic board, either the inserts and the scaler handle are galvanically insulated, thus assuring absolute safety to the user.

**CONTENTS**

TN250 electronic board, potentiometer for power adjustment, potentiometer dial, two plastic screws for board assembly.

**ASSEMBLY**

Install and fix TN250 heat sink on the frame of the Unit using the two supplied screws and leaving the board sufficiently far away (at least 10 mm) from other components and/or the Unit frame.

Connection with TITANUS®E scaler (REF 400.00 or REF 450.00): connect red and black wires coming from the supply silicone hose (REF 330.1x) to the TN250 connector, as shown on attached wiring diagrams.

The board is activated by closing the START-STOP switch. The board cannot be activated if the handpiece is not connected to the supply silicone hose.

If the scaler is to be used for Endodontic procedures, it is necessary to limit output power: this is accomplished by closing a switch from the connector no.4 to the common ground (connector no.6).

A 24 V solenoid valve for water flow control can be connected as shown on the attached wiring diagrams.

To protect the board it is necessary to connect an external fuse (F 2 A) as shown on the attached wiring diagrams.

**POWER ADJUSTMENT**

Power can be adjusted by means of the supplied potentiometer (4,7 K $\Omega$  linear) or by an external continuous voltage (0  $\div$  5 Vdc). Please see attached wiring diagrams.

**SERVICE**

In case of malfunction and for any overhaul and repair work, please contact an authorized technician or directly to TKD.

**REFUSE DISPOSAL**

This device needs to be recycled. Electrical and electronic equipment may contain dangerous substances which constitute health and environmental hazards. The user must return the equipment to its dealer or establish direct contact with an approved body able to process and derive value from this type of equipment (European Directive 2002/96/EEC).

**DECLARATION OF RESPONSIBILITY**

The manufacturer will not be held responsible in case of:

- use on a electric connection not meeting the requirements of CEI EN 60601-1 Standard
- operations or repairs carried out by technicians not authorized by the manufacturer

## TN250 SPECIFICATIONS

The device:

- does not produce electromagnetic interference and meets the requirements of CEI EN 60601-1-2 Standard
- is not suitable for use in the presence of a flammable anaesthetic mixture with air or with oxygen or nitrous oxide
- must be connected to a power transformer having a minimum insulation strength of 500 V between secondary winding and earth and a minimum insulation strength of 1500 V between secondary winding and mains

Classification:	device with type B applied part according to CEI EN 60601-1
Power supply:	24 Vac (50/60 Hz) or 32 Vdc
Maximum consumption:	21 VA at 24 Vac, 850 mA at 32 Vdc
External protection fuse:	F 2 A
Operating frequency:	28 ÷ 31 KHz
Maximum Output voltage:	125 V <sub>RMS</sub>
Maximum Output power:	8 W
Absolute maximum Input voltage:	40 Vdc
Sizes:	60 x 50 x 30 mm
Working temperature:	0 ÷ 40 °C (32 ÷ 104 °F)
Storage temperature:	-20 ÷ 100 °C (-4 ÷ 212 °F)
Working humidity:	30 ÷ 90 %
Storage humidity:	10 ÷ 90 %
Guarantee:	2 years

Protection: main oscillator is shut off in the following events:

- 1) Breaking of one supply wire
- 2) Short circuit between handpiece supply wires
- 3) General malfunction of the electronic control

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