Komet Dental

Gebr. Brasseler GmbH & Co. KG

Trophagener Weg 25

32657 Lemgo · Germany

Telefon +49 (0) 5261 701-0

Telefax +49 (0) 5261 701-289

info@kometdental.de

www.kometdental.de

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EN

Ultrasonic tips for prophylaxis and periodontal treatments

1. Indication

- Ultrasonic scaler tips are used in prophylactic treatments. They are intended for the sub and supragingival scaling of calculus, removal of soft and hard concretions and in some cases, for the removal of residual cement.

- Ultrasonic periodontal tips are used for the removal of soft and hard dental plaque from deeper periodontal pockets, root furcations and concavities.

- The tip PL9 can be used in the subgingival and supragingival region.

2. Contra-indications

Cardiac pacemakers and defibrillators: Patients wearing pacemakers or defibrillators are excluded from treatments with ultrasonic tips because the ultrasonic vibrations of the instruments and /or hand pieces may impair the proper function of these devices.

3. Warnings

The useful life of ultrasonic tips is limited. They therefore have to be used with the utmost care.

Attention:

- Risk of fracture: Spontaneous and premature fracture must be expected. The ultrasonic tips may fracture unexpectedly and without warning. In order to reduce the risk of fracture, make sure to stick to the correct power setting (see instructions for use) and observe all recommendations for use. Never allow the instruments to jam. Do not use instruments for leverage. Avoid excessive contact pressure. Use ultrasonic tips with the utmost care as even minor mechanical damage can provoke the fracture of the instruments.
- Risk of ingestion: Fracture of the ultrasonic tip leads to the risk of the tip being accidentally swallowed. Take all appropriate steps to prevent this.
- Premature wear: Due to the nature of their range of application, the functional parts of these sonic tips will always be subject to moderate wear, even when the tips are applied correctly and according to their intended use. A test card (411698) can be used to determine the extent of wear of certain tips (A; P; PS; PL3, PL9). Excessively worn tips have to be exchanged.
- Work on hard dental substance can cause premature wear of the diamond coating. Excessive contact pressure can provoke damage to the working part or the diamond coating. The instruments therefore have to be checked for signs of wear or damage after every use and replaced if necessary.
- Excessive heat generation: Excessive contact pressure and /or worn or blunt instruments can cause hazardous heat generation which increases the risk of thermal necrosis (damage to tissue). Excessive contact pressure therefore has to be avoided. Damaged ultrasonic tips (changed shape, corroded surfaces or defect thread etc.) have to be replaced.
- Excessive heat generation: The use of a high volume evacuator might lead to an excessive evacuation of the cooling liquid, in which case the dental surface or the instrument will not be cooled sufficiently. The use of a high volume evacuator should therefore be avoided.
- Excessive heat generation: The tip can heat up very strongly and quickly in case of an insufficient flow rate of the cooling liquid through the ultrasonic drive. To prevent this, always observe the flow rate indicated for each instrument. Do not under any circumstances work without cooling.
- Non-sterile instruments: Instruments that have not been sterilized before the first use or not correctly reprocessed after each subsequent use can cause infection. The instruments therefore have to be cleaned and sterilized prior to each use in compliance with the reprocessing directions indicated in the present instructions for use.
- Damage to hard tissues: Excessive contact pressure might damage the hard dental substance and /or root surfaces. Excessive contact pressure therefore has to be avoided. This applies in particular to the periodontal tips which may only be used with low contact pressure.
- Damage to hard tissues: Ultrasonic tips perform longitudinal oscillations relative to the longitudinal axis of the instrument. In order to prevent damage to the hard dental substance, avoid frontal positioning of the instrument on the dental surface.
- Ultrasonic tips may not be used on metal or ceramic restorations or dentures as these might get damaged by the use these instruments.
- Do not use the hand piece in the vicinity of flammable anesthetics or gases.
- The ultrasonic tips may only be used by properly trained and qualified personnel.
- Always wear protective gloves when checking, inserting or removing the ultrasonic tips to avoid risks of injury and infection. Exchange tips regularly.
- Incorrect use can lead to the above detailed undesirable effects as well as to injury to the user, the patient or third persons.

4. Appropriate use

Ultrasonic tips delivered in a non-sterile condition have to be cleaned, disinfected and sterilized prior to first use. The ultrasonic tips can be used in the following hand-pieces:

- Piezon® Master 700, MiniMaster®, MiniMaster® LED and MiniPiezon® of the company EMS
- Titanus® E of the company TEKNE DENTAL

To insert or remove the ultrasonic tips, use the enclosed torque wrench only. Screw the chosen tip into the hand piece, turning the torque wrench in a clockwise direction. To fasten the tip correctly, turn the torque wrench another quarter-turn after the tip has been completely screwed in.

To remove the tip, place the torque wrench onto the tip and unscrew in an anticlockwise direction.

5. Recommendations for use

Start up ultrasonic tip prior to application to the tooth, supplying an adequate amount of spray cooling. All ultrasonic tips described in the present instructions require plenty of cooling liquid, at least 50 ml/min.

Recommended use of each individual ultrasonic tip:

Group of tips used for scaling:

A.EM1.: Supragingival scaling in all quadrants. In case of hard concretions or an initial treatment, choose a high power setting. Patients with a low pain tolerance can be treated at a medium power setting.

B.EM1.: Removal of supragingival concretions from lingual surfaces in all quadrants and removal of residual cement. Thin deposits on occlusal surfaces can also be removed, provided that a low power setting is chosen. Patients that are sensitive to pain should be treated at a low power setting, whereas normal treatments can be carried out at normal setting. To remove particularly hard concretions, a high power setting may be chosen.

C.EM1.: Removal of large supragingival concretions from front teeth and removal of dental cement. Thanks to its flat tip, this instrument is particularly suitable for work on outer dental surfaces. For normal treatments, choose a medium setting, whereas patients with a lower pain tolerance should be treated at a low setting. To remove particularly hard concretions, a high power setting may be chosen.

P.EM1.: Supragingival scaling or removal of subgingival concretions in all quadrants, in particular in the interproximal cavities and in the region of the sulcus. For normal treatments, choose a medium setting, whereas for patients that are sensitive to pain and follow-up treatments, preference should be given to a low power setting. To remove particularly hard concretions, a high power setting may be chosen.

PS.EM1.: Removal of subgingival deposits from root surfaces and for rinsing periodontal pockets as well as periodontal follow-up treatments. Use at a low power setting. Even particularly hard concretions may only be removed at a medium power setting.

Group of tips used for periodontal treatments:

PL1.EM1. + PL2.EM1.: Periodontal debridement on hard-to-access interproximal surfaces. Always use the tip at a low or medium setting, whatever the application. For safety reasons, use the outer lateral surface of the tips only and work with low pressure.

PL3.EM1.: For cleaning and disinfecting periodontal pockets.

PL4.EM1. + PL5.EM1.: Periodontal debridement on furcations and concavities. Always use the tip at a low or medium setting, whatever the application. All lateral sides of the instruments may be used with low pressure.

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Instructions for use

Please read carefully and retain for future reference.

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PL9.EM1 : The scope of application of this tip is situated between scaling and periodontology, which means that it is suitable for removing concretions in the subgingival and supragingival region of all quadrants. It can be used at a high power setting, however, patients that are sensitive to pain can also be treated at medium setting.

DPL3.EM1. : Polishing of dental surfaces after cleaning. Work at a low or medium power setting.

HPL3.EM1. : Thorough scaling of the root under direct vision. Also suitable for smoothing overhanging restoration material and enlarging furcation roofs. Work at a low or medium power setting.

Failure to follow the above instructions can lead to the undesirable effects described in chapter 3 "warnings".

6. Maintenance/care

Place ultrasonic tips in a cleaning/disinfection tank (Fräsator) filled with a suitable detergent/disinfectant (e.g. Komet® DC1*) immediately after use in the mouth to prevent drying of residues on the instruments (protein fixation). Hold ultrasonic tips at an angle when placing them into the liquid to ensure that all hollow spaces are reached. It is recommended to reprocess the instruments within one hour of use at the very latest. The instruments should be in the cleaning/disinfection tank (Fräsator) when transported to the site where the reprocessing is to take place.

7. Reprocessing

Manual pre-cleaning (prior to mechanical and manual reprocessing):

- Remove ultrasonic tip from the cleaning/disinfection tank (Fräsator) immediately before reprocessing.
- Check patency of the cooling passage by guiding a fine needle or wire through the passage. Replace instrument if the passage is obstructed.

Mechanical reprocessing (validated):

- Rinse ultrasonic tip thoroughly under running water to prevent any residues of the detergent/disinfectant from getting into the washer/disinfect.
- Unscrew the nozzle on the injector bar of the washer/disinfect. Then screw the external thread of the rinse adapter SF1977 in the silicon holders on the injector bar. Now screw the ultrasonic tip in the internal thread at the top of the rinse adapter. Observe the directions in the instructions for use provided with the SF1977 (fig. 2).
- Put chemical detergent into the washer/disinfect, following the indications on the label and the instructions of the manufacturer of the washer/disinfect.
- Start the Vario TD program including thermal disinfection. Thermal disinfection takes place allowing for the AO value and observing national provisions (prEN/ISO 15883).
- On completion of the cycle remove ultrasonic tips from the washer/disinfecter and dry (preferably with compressed air according to the recommendations of the Commission for Hospital Hygiene and Infectious Disease Prevention of the Robert Koch Institute). When drying the inside of the instruments, please make sure to hold the compressed air gun close to the cooling passage, in order to guarantee that sufficient air flows through the ultrasonic tip.
- Visual examination to ensure that the instrument is clean and undamaged. If after mechanical reprocessing there are still visible residues of contamination, repeat the cleaning and disinfecting process until no visible contamination is left.

Manual reprocessing (alternative):

- Place ultrasonic tip onto a 10 ml syringe and rinse with detergent/disinfectant. In case of visible contamination emerging from the perforations, rinse again with detergent/disinfectant.
- Remove stubborn contamination with a nylon brush under running water, turning the ultrasonic tip constantly.
- Rinse ultrasonic tip thoroughly with running water.
- Visual examination to ensure that the instrument is clean. If there are still visible residues of contamination, repeat the cleaning process until no visible contamination is left.
- Place ultrasonic tip in a suitable container into the ultrasonic device or instrument bath filled with detergent/disinfectant.

- During chemical disinfection in the ultrasonic device or instrument bath, observe the instructions of the manufacturer regarding concentration and immersion time. Please note that, due to their internal orifice, the chemical disinfection of sonic and ultrasonic tips (instruments with internal cooling) in an ultrasonic bath using a 2% solution always takes 10 minutes. Be sure to observe the full correct immersion time which does not start until the last tip has been immersed. Attention: do not exceed 45 °C (risk of protein coagulation)!

- On completion of the immersion time, rinse ultrasonic tip thoroughly with suitable water (preferably with demineralized water to avoid residues of lime). Rinse the perforation of the ultrasonic tip with at least 10 ml of demineralized water by means of a cannula, to ensure that no residues of the detergent are left in the cooling passage.
- Dry ultrasonic tip (preferably with compressed air according to the recommendations of the Commission for Hospital Hygiene and Infectious Disease Prevention of the Robert Koch Institute). When drying the inside of the instruments, please make sure to hold the compressed air gun close to the cooling passage, in order to guarantee that sufficient air flows through the ultrasonic tip.
- Visual examination to ensure that the instrument is clean and undamaged. If there are still visible residues of contamination, repeat the cleaning and chemical disinfecting process until no visible contamination is left.

According to the recommendations of the Commission for Hospital Hygiene and Infectious Disease Prevention of the Robert Koch Institute, reprocessing should preferably be carried out mechanically.

Information on the reprocessing of instruments can be downloaded from our website www.kometdental.de or requested directly from the manufacturer Gebr. Brasseler.

Sterilization in the autoclave:

Make sure that the packaging is suitable for ultrasonic tips classed as risk group critical B and the chosen method of sterilization.

Single pack: The packaging must be large enough to ensure that there is no pressure on the seal.

Multipack: Place ultrasonic tips inside the torque wrench into a suitable sterilization container (e.g. 97507). For sterilization, the container has to be sealed in suitable packaging.

For this procedure, the following instructions have to be observed:

Steam sterilization using a vacuum process at 134 °C in a device that complies with the provisions of DIN EN 13060; with validated processes.

- Fractionated pre-vacuum (type B)
- Sterilization temperature: 134 °C
- Hold time: at least 5 minutes (full cycle)
- Drying time: at least 10 minutes

In order to prevent staining and corrosion, the steam must be free of particles. The recommended limits for particle contents in feed water and condensed steam are defined by standard DIN EN 13060. Make sure not to exceed the maximum capacity of the sterilizer when sterilizing several instruments. Follow the instructions of the device manufacturer.

Universally valid notes:

Observe the legal provisions regarding the reprocessing of medical products valid in your country (e.g. www.rki.de). The manufacturer confirms that the above detailed reprocessing methods are suitable for preparing the above named instrument group to enable their reuse. The person carrying out the reprocessing is responsible for ensuring that the applied method is carried out with appropriate equipment, materials and trained personnel at the reprocessing site and that it actually achieves the desired result. To guarantee this, routine controls of the validated mechanical preparation methods are necessary. Any deviation from the above detailed process must be carefully checked by the operator to ensure effectiveness and to avoid possible adverse consequences.

8. Liability:

It is the responsibility of the user to check the products prior to use to ensure that they are suited for the intended purpose. In case of contributory negligence by the user, Gebr. Brasseler partially or totally declines liability for all resulting damages, particularly due to non-observance of our recommendations for use or warnings as well as inadvertent misuse by the user.

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© 07/2015, Gebr. Brasseler GmbH & Co. KG