

EURODENT


ISOTRON


ISO+plus

ISO+plus lite

isoconcept

ISTRUZIONI PER L'INSTALLAZIONE
INSTALLATION INSTRUCTIONS
NOTICE D'INSTALLATION
MONTAGEANLEITUNG
GUÍA DE INSTALACIÓN

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INDEX

Page

SERVICE REQUIREMENT	26
DECLARATION OF RESPONSIBILITY	26
ARRANGEMENT	27
A PACKAGE OF THE ISORAMA AND ISO CHAIRS	28
TRANSPORT AND INSTALLATION	29
B INSTALLING THE ISORAMA CHAIR	30
PREPARATION	30
RELEASING THE BACKREST	30
FASTENING THE BACKREST	30
OPERATING TEST	31
BACKREST PAD	31
SEAT PAD	32
C INSTALLING THE ISO CHAIR	33
BACKREST	33
BACKREST UPHOLSTERY	33
SEAT UPHOLSTERY	34
FITTING OF THE ARMRESTS	34
D HEADREST FOR THE ISORAMA AND ISO CHAIRS	36
E INSTALLING THE ISOTRON, ISOPLUS, ISOPLUS LITE, ISOCONCEPT DENTAL UNITS	37
F ELECTRICAL AND ELECTROPNEUMATICAL CONNECTIONS	42
ISOTRON, ISOPLUS, ISOPLUS LITE	42
ISOCONCEPT	44
G CONNECTING THE SUCTION PUMP	46
H ARRANGEMENT OF ISOSUN LAMP FOR ISOTRON, ISOPLUS, ISOPLUS LITE	46
I ISOCONCEPT ADJUSTEMENTS	47
PRESSURE ADJUSTMENT TO INSTRUMENTS	47
TRIMMER FUNCTIONS	47

SERVICE REQUIREMENT

The unit to be installed is CLASS I, with attached parts of TYPE B , with electronic bistoury or electrobistoury with attached part of TYPE BF  (ISOTRON, ISOPLUS, ISOPLUS LITE only)

 **0051** The apparatus comply with the fundamental requirements of EC Directive 93/42.

Classification headings (EN60 601-1):

- 5.3 Common equipments (equipments with no waterproof casing);
- 5.5 Equipments which cannot be used in presence of an anaesthetic mixture inflammable by air or by oxygen or by nitrous oxide;
- 5.6 Equipments for continuous working with intermittent load.

Protection earth

Warning, see enclosed documents

Power: 220 V / 50 Hz / 1350 VA (550 VA only unit) for ISOTRON, ISOPLUS, ISOPLUS LITE.
220 V / 50 Hz / 1000 VA (200 VA only unit) for ISOCONCEPT.

Operating system: continuous, with intermitting load

Max. lifting capacity: 120 kg (dental unit and lamp excluded)

The floor must be in accordance with DIN 1055 B1.3 and DIN 18560 T1 regulations.

It is necessary to instal an outer bipolar main switch to wall, in compliance with the law and with the following electric characteristics: 250 V / 10 A for disconnecting the unit from the main feed in case of intervention or inspection in it.

The electrical system must be in accordance with IEC regulations 64.4, and made up of difference protective conductors (with index < 30 mA) at 16A/250V.

Moreover, earthing must conform to the art. 3.2.01 of the regulations mentioned above.

Water pressure at the unit inlet must range from 1,8 to 3,0 bar. The pipes must be capable to supply at least 5 l/minute of water flow rate.

Water characteristics must conform to DIN 1988 rules, with a hardness gradient $\leq 8^\circ$ dH (corresponding to 1,432 mmol/lit). We recommend installing a pressure reducing tap at the water hose outlet (preferably of 15 mm \varnothing).

The air pressure provided by the compressor must range from 5,0 to 7,0 bar. Use a dry compressor, with air output of 60 lit/min. at 5 bar. We recommend installing a tap at the hose outlet (preferably of 15 mm \varnothing).

The diameter of the emptying hose must be 40 mm, and 32 mm at its inlet; moreover the descent gradient must be at least 1 cm/mt, and the bend angle should not be less than 135° .

The diameter of the suction pipe must be approx. 32/40 mm. We recommend following the instructions provided by the suction pump's manufacturer.



All arrangements must be carried out according to the indications provided on the following page and consulting the full scale installation plan (code 633580010) attached to each dental unit.

Please refer to Spare Parts catalogue and apply to EURODENT for further information about what is not included in it.

EURODENT accepts responsibility for the equipment safety and reliability if and only if:

- installation, assembly, extension, adjustment and repair operations have been carried out by authorized people.
- the general and electric installation of the room complies with the above written requirements.
- the equipment has been installed in compliance with the "INSTALLATION INSTRUCTIONS".
- the equipment is used in compliance with the "OPERATOR'S HANDBOOK".

ARRANGEMENT

- A.** Water: 3/8" G
- B.** Air: 3/8" G
- C.** Outlet: 32 Ø
- D.** Suction: 32 Ø
- ⊙ Pump starting: 2 x 1 mm²
- E.** ⊙ Bell: 2 x 0.5 mm²
- ⊙ Net: 3 x 2.5 mm²

A - PACKAGE OF THE ISORAMA AND ISO CHAIRS


The special packing system of these chairs has been thought in order to reduce its overall dimension.

For such a reason, the body of the chair is placed into the packing case ① with the backrest bent over the seat (see figure A0), whereas upholstery, seat, backrest and headrest are contained into the packing case ②.

When removing the chair from the packing case, cut the box ① along the dotted line.

A0

ISORAMA


ISORAMA

Net weight 131
Gross weight 152
..... 0,6 m³
cm. 125 x 69 x 70

ISO

Net weight 116
Gross weight 137
..... 0,6 m³
cm. 125 x 69 x 70

ISO

ISORAMA

Net weight 13
Gross weight 18
..... 0,2 m³
cm. 115 x 65 x 25

ISO

Net weight 10
Gross weight 16
..... 0,12 m³
cm. 93 x 51 x 25

TRANSPORT AND INSTALLATION

Remove the chair from its packing and, before moving it, take the covering guard of the connecting box off, so as to make the seizing of the chair and transport easier. Do not lift the chair by seizing the pantograph head or the legrest.

Position the chair so as to place all couplings inside the connecting box (placed at the bottom of the chair parallelogram). Make sure that you do not inadvertently squash or tear the electrical cables of the chair when slipping it on the floor. Pay special care to the cable of the foot-control.

Check the level evenness of the floor with respect to the bearing surface of the base. Fasten the chair by fitting the expansion stoppers type Fischer (supplied with the accessories) into the two housings (see figure A1).




M6 expansion bolts and proper spanner are supplied with the chair.

The hole in the floor must be at least 115-mm deep.


As regards installation, we advise connecting **temporarily** the chair to the main supply in order to have it lifted simply by operating the foot-control.

BE CAREFUL! Connect the backrest to the motor and to the Trendelenburg device before making it perform any movement. When the chair reaches the height right for operating comfortably, cut out voltage.

A1


A2

IMPORTANT: if assembly does **NOT** allow connection inside the chair base box, you can use the front base frame as outlet window, by removing the baffle plugs with a cutter.


B - INSTALLING THE ISORAMA CHAIR

PREPARATION

Position properly the chair and remove the armrest by unscrewing the two M6 nuts with a 10-mm spanner.


Cut the 2 plastic clamps which hold the worm screw of the backrest motor and the tie rod for the Trendelenburg movement.

RELEASING THE BACKREST

By means of a 6-mm Allen wrench unloose the M4 dowel and remove the backrest pin from its housing using a 6-mm Allen wrench.

Note: Pay attention to the shim washer.


FASTENING THE BACKREST

Lift the backrest so as to make the hitch point match with the worm-screw eye and insert the backrest pin into its housing.

Tighten the pin firmly by means of an Allen wrench, then unloose the pin about 30° (less than one fourth turn of the wrench) in order to make it run freely and with minimum play. Finally, fasten it by tighten the M4 dowel.


Remove the polystyrene cube which protects the headrest.

B4

OPERATING TEST


Now, proceed to fit the tie rod operating the Trendelenburg movement.
 By using a 13-mm spanner remove the nut, then insert the tie rod eye and lock the rod again by the proper nut.
 Supply voltage to the chair and perform some trial movements in order to ascertain that the stop microswitches work properly.
 Finally, mount upholstery.

B5

BACKREST PAD

Position upholstery higher than backrest. Then, slide it downwards, in order to secure it properly and firmly to the three fastening screws, as shown in figure.

B6


Fasten upholstery by the T-clamp, and press it against the two upholstery tags, as shown in figure.

Relocate the armrest.


SEAT PAD

Connect the earth cable as shown in figure.

If not previously fitted, located the dowels in their seats (the long one goes in the external rear hole), and fasten them by means of the proper lock nut.
Fit the two front dowels as well.


B7


B8

Then, as shown in the figure, position upholstery taking care of inserting the dowels into the proper holes.


B9

Now, secure the pad by tighten the long cylindrical nuts of the dowels.
To carry out this operation, use the long spanner supplied.

B10


In order to fix the front side, use M4 x 40 screws.

C - INSTALLING THE ISO CHAIR**C1****BACKREST**


Remove the chair from the packing by cutting along the dash-line, and lift the backrest until the hitch point matches the push rod.

Carefully insert the shim in its housing, and after fitting the screw in its housing, tighten it with the suitable wrench.

Now the chair can operate and some test movements should be performed to ensure that stop microswitches operate properly.

**C2**


Then proceed by mounting upholstery.

**C3****BACKREST UPHOLSTERY**

Position upholstery slightly over the backrest and make it slip down so that it fits in the two fastening screws.

When upholstery is properly positioned, fit the screw into the suitable hole on the back side and tighten the M4 dowel by using the 2 mm Allen wrench.


If the pad does not match properly with the backrest, adjust the two flathead screws as described for the headrest pad.

C4

SEAT UPHOLSTERY


Position the pad on the seat and secure it with the three M4 screws by tightening the lock nuts, so as to make the parts join perfectly.

C5


FITTING OF THE ARMRESTS

First slip the shaped $\varnothing 8$ peg into its backrest housing, keeping the flat part of the peg towards the fore face of the backrest (that is, the armrest and the backrest must determine a 90° angle).


C6

Then rotate the armrest in order to lock the peg into its housing, now you can fix the lower part of the armrest.


C7

For the left armrest this part is fixed with its bolt, whereas for the right one, it shall be driven in its sliding block.


C8

The right armrest is sliding: to declutch it just pull upward the fore part.

C9

Slide the armrest towards the chair front, until it reaches the backstop onto the sliding carriage.


C10

Moreover, the right armrest may be bent to reach its negative angle: also in this position it follows the tilting movements of the backrest.

It can be easily removed by putting the movable tab of the extractable pin in the proper extraction position (tab aligned at 180° with the pin axis) and then drawing the backrest pin following, contrariwise, the aforesaid instructions.


D - HEADREST FOR THE ISORAMA AND ISO CHAIRS

Assembling is easy and intuitive:

Position the headrest. Once the screw heads are placed in their seats, push down and insert the two dowels in the respective holes.

Tighten the fastening dowels by pressing onto both parts and using a 2 mm Allen key.

Proper matching of pad and backrest can be obtained by loosening the nuts that block the two flathead screws, adjusting the latter ones in depth so as to accommodate properly the headrest pad, and then tightening again the nuts.

E - INSTALLING THE ISOTRON, ISOPLUS, ISOPLUS LITE, ISOCONCEPT DENTAL UNITS

ISOTRON ISOPLUS ISOPLUS LITE	ISOCONCEPT
------------------------------------	------------

"BASE"

Net weight 71	Net weight 46
Gross weight ... 88	Gross weight ... 65
..... 0,555 m ³ 0,555 m ³
cm. 96 x 89 x 65	cm. 96 x 89 x 65


"INSTRUMENT TABLE"

Net weight 22	Net weight 18
Gross weight ... 33	Gross weight ... 30
..... 0,428 m ³ 0,428 m ³
cm. 88 x 76 x 64	cm. 88 x 76 x 64

The dental unit is supplied into two packing cases: the first one labelled "BASE" contains the body and the assistant arm; the second one labelled "INSTRUMENTTABLE" contains the operator retraction arm, the modules and the tray holder.

Provide air and water hose couplings with taps. It is absolutely necessary to clean out air and water pipes before connecting the dental unit to them.


E0

E1

In addition to the guard of the connecting box, we recommend that you remove even the lid of the chair parallelogram, in order to operate easily when introducing hoses and cables into the dental unit. Remove such mentioned lid by unscrewing the two long brass nuts which are placed under it and which fasten it to the upper rods of the parallelogram.

The chair shall be fastened to floor by inserting expansion bolts into the two holes made on the front side of the base footboard.

E2


Now, mount the body of the dental unit on the chair support by means of the stud bolts and the socket head screws supplied (placed inside the "ACCESSORY" box, into the packing case labelled "BASE").

ISOTRON, ISOPLUS, ISOPLUS LITE

Before fastening screws and nuts, level the body properly horizontally and vertically. Please, note that stud bolts shall be inserted into the support lower holes, whereas screws into the higher ones.


ISOCONCEPT

To make such operation easier, fit the additional Allen-head stud bolt supplied into the second of the two lower holes. Level the body properly (horizontally and vertically) before tightening the two screws and the middle nut; remove the Allen-head stud bolt from the lower hole and replace it with the third screw supplied.

E3

Into the shaped hole on the parallelogram arm of the chair first insert the suction hose, then the electric cables wrapped with tape to the Rilsan hoses (green and white), finally insert the drain hose. Lead all hoses and cables up to the connecting box.


E4

ISOTRON, ISOPLUS, ISOPLUS LITE

The accessories boxes contain a foot-control, the suction tubings with their proper hoses, an assistant syringe (on request) with a clip, four screws, a grey cord to connect the foot-control cable to the connecting box, main solenoid valves for water and air and safety electrical connectors.

The cord has a squared-plug terminal and a circular-base terminal: push the first terminal mentioned into the hose which, inside the base footboard, joints the outside to the connecting box. Then, pull it until the frame of the circular base adheres to the frame of the hose. Tighten the four self-tapping screws in the holes of the frame, so that they penetrate completely the polyurethane and fasten properly the cord base.


E5

ISOCONCEPT

The accessories boxes contain the foot-control, the suction tips with their proper hoses, the nurse syringe (on request), four screws, 2 couplings, the water and air main valve, safety electric connectors and bridge-p.c.b.

Slip electric cables and hoses of the foot-control through the sleeve which, inside the base footboard, leads to the connection box, and push the flange of the foot-control cord against the sleeve frame. Tighten the four small self-tapping screws in the holes of the frame, until they have entered completely the polyurethane and lock the flange.

E6

Now, install the instrument table arm of the dental unit (however first make the cable and hoses run through the body) and block the travel of the arm by inserting the screw and the nut (supplied with accessories) into the pin hole.


ISOTRON, ISOPLUS, ISOPLUS LITE

Plug the arm cord into CN1 (on left top side), on the low-tension rod. Then, connect the white hose to the quick fitting which provides modules with air, placed on the right side of the air reduction units. Moreover, connect the red hose (ISOTRON only) to the fitting placed on the lid of the amalgam header, and the green one to the water distribution unit (for tumbler, assistant syringe and modules), placed on the chair side.

ISOCONCEPT


Plug the arm cord into X1 and X2 (on left top side), on the p.c.b. Then, connect the white 3/5 hose to the union for permanent air, placed on the right side of the air reduction unit. Moreover, connect the green one to the water distribution unit. The other pipes should be connected to the proper unions placed on the water unit.

IMPORTANT: When assembly, always respect the correspondences of the pipe colours.

E7


Fasten the tray holder under the module-holder table (to make it easier, use a 17 socket spanner).

E8

Now, insert the cable and the pin of the assistant arm into the proper housing (at the bottom of the body, on the chair side), and lock the pin with the snap ring and its lock space washer. Make the cable run through the body and plug it in the C5 socket (for ISOTRON) or in the C2 (for ISOPLUS and ISOPLUS LITE), placed on the low-tension p.c.b., or in the X5 socket of the body p.c.b. (for ISOCONCEPT).

NOTE: ensure that the cable does not hang under the body). Insert the friction unit into the hole beside the pin housing and tighten the screw properly.

E9


ISOTRON

If the dental unit is equipped with assistant syringe, follow the steps below:

Insert the syringe hose into its housing on the tubings carrier, and make it run upwards through the hole placed between the suction filter and the tank for tubings cleaning. Then, lock the end bush with the clip supplied.

Unscrew the screws which lock the lid of the box of the 6-function syringe rod located on the suction cleaning tank and connect the hoses and wires of the syringe hose. Finally, replace the lid.

Now, insert the suction pipes into their housings on the tubing carrier and connect the couplings to the proper fastening bushes.

E10


ISOPLUS, ISOPLUS LITE

If the dental unit is equipped with assistant syringe, follow the steps below:

Insert the syringe hose into its housing on the tubings carrier, and make it run upwards through the hole placed between the suction filter and the tank for tubings cleaning. Then, lock the end bush with the clip supplied.

Unscrew the screws which lock the lid of the box of the 6-function syringe rod located on the suction cleaning tank and connect the hoses and wires of the syringe hose. Finally, replace the lid.

Now, insert the suction pipes into their housings on the tubing carrier and connect the couplings to the proper fastening bushes.


E11

ISOCONCEPT


If the dental unit is equipped with assistant syringe, follow the steps below:

Insert the syringe hose through the hole at the left side of the suction filter. Then, lock the end bush with the clip supplied. Connect the water hose to the corresponding coupling placed near the suction filter.


Connect the air hose to the air reduction unit for the syringe.

**E12**

Insert the lamp-stand into its housing and lock it by tightening the two screws already fitted on the lamp (body inner side). Make a test chain run through the lamp-stand and bind it with adhesive tape to the cable coming out the pin of the lamp arm. By pulling the test chain, insert the pin into the stand, then connect the lamp cable to the 3-4 terminals of the M4 socket on the AT/BT rod (on the chair side). (For the ISOCONCEPT model, the terminal in question is the XT1 on the body p.c.b.).

**E13**

Position the bowl in its housing.

**E14**

Connecting box: fix the solenoid valve unit for air and water on the left side, by using the two screws supplied. Pay attention to the flow direction inside the valves (shown by an arrow on the brass body) when joining the air and water hoses. The electrical connectors unit shall be fastened on the box right side, by means of the screws and nuts provided.

F - ELECTRICAL AND ELECTROPNEUMATICAL CONNECTIONS



ISOTRON, ISOPLUS, ISOPLUS LITE

The box receives several cables, i.e. the grey cable which connects the foot-control cable to the box, the black four-wire cable for serial connection with the dental unit, and the white three-wire cable (brown, blue, and yellow/green) for chair mains supply.

The following ones are coming from the dental unit:

one three-wire (blue, brown, yellow-green) black cable, one three-wire (blue, brown, yellow-green) white cable, one more white (two-wire) and one more black one (two-wire), a grey cable for dental unit connection to foot-control extension, a four-wire black cable for serial connection to the chair, and one more black cable ending by two faston connections.

The safety connectors (a black and a white one) have ten terminals.

F1

Plug both fast-on wires to the main solenoid valve for air and water (24 V).

Connect the foot-control extension cable (1) to the foot-control cable (2).

Connect the black four-wire cable (3) for chair serial connection to the black four-wire cable (4) from the dental unit.


The wires of the white cable (to switch the suction pump control circuit) must be connected to the first two brown contacts, whereas the wires of the black cable must be connected to the other two brown contacts (to switch the circuit of nurse bell or any other service).

In the same way, the three wires of the white thick cable shall be connected to the terminals of the white connector (blue to N terminal, brown to L terminal and yellow/green to earth terminal). To the terminals mentioned the wires of the white cable for feeding the chair shall be connected according to their colours and L, N and earth signs.

Leave the black cable from the chair footcontrol connected to the black cable for chair movements control.

Finally, the wires of the black feeding cable coming from the dental unit shall be connected to the three terminals of the black connector (blue to N, brown to L and yellow/green to earth). To the terminals mentioned shall be connected the wires of the cable coming from power mains (220 VAC). Take care to protect all of the hanging connections by the plastic shield supplied among the accessories.

F2



Replace the parallelogram lid and fasten it by means of the two brass nuts/tie rods, then close the connecting box by using the proper guard. Mount the two lateral panels of the dental unit body and fasten the cord provided with a cable terminal to the screw placed on the external panel.

ISOCONCEPT

Electric connections in the connecting box: the box receives several cables: the black cable of the chair foot-control (male flat plug connected to the black cable for chair control), the black cable for chair movement control (female flat plug connected to the chair foot-control black cable), the black four-wire cable for serial connection to the dental unit (not to be used), and the white three-wire cable (brown, blue, yellow-green) for chair mains supply.

Coming from the dental unit:

one black three-wire mains supply cable (blue, brown, yellow-green), one white three-wire cable (blue, brown, yellow-green) another black cable with two black and two white wires, one black cable with nine-wire female flat plug for parallel connection with the chair, one black cable of the dental unit foot-control and one black cable with two fast-on plugs. There are four sockets on the bridge-p.c.b., while the safety connectors (a black and a white one) have ten terminals.



F3

Plug both fast-on wires to the main solenoid valve for air and water (24V).

The dental unit foot-control cable and the nine-wire parallel cable must be plugged into the relevant connectors. Disconnect the seven-wire parallel cable of the chair from the one of the chair foot-control and plug them in the relevant connectors of the bridge-p.c.b. The bridge-p.c.b. must not be used if the dental unit is not there or it consists of a simple Isoconcept cuspidor; in this case leave the seven-wire parallel cable of the chair and the one of the chair foot-control one plugged into the other. The white wires (for closing the suction pump control circuit) must be connected to the first two brown terminals, while the black wires must be connected to the other two (to close the circuit for nurse bell or any other service circuit).

In the same way, the three wires of the white thick cable shall be connected to the terminals of the white connector (blue to N terminal, brown to L terminal and yellow/green to earth terminal). To the terminals mentioned the wires of the white cable for feeding the chair shall be connected according to their colours and L, N and earth signs. Leave the black cable from the chair footcontrol connected to the black cable for chair movements control. Finally, the wires of the black feeding cable coming from the dental unit shall be connected to the three terminals of the black connector (blue to N, brown to L and yellow/green to earth). To the terminals mentioned shall be connected the wires of the cable coming from power mains (220 VAC). Take care to protect all of the hanging connections by the plastic shield supplied among the accessories.


F5

Connect the foot-control hoses colour-to-colour by the relevant couplings, keeping the correspondences according to the colours of the hoses.

Replace the parallelogram lid and fasten it by means of the two brass nuts/tie rods, then close the connecting box by using the proper guard. Mount the two lateral panels of the dental unit body and fasten the cord provided with a cable terminal to the screw placed on the external panel.

G - CONNECTING THE SUCTION PUMP

G1


When picking up any of the two suction tubings, the dental unit provides a contact between the two white wires from the black cable (i.e. a switch which conducts when a tubing is up and cuts off when tubings are at rest). Therefore, the pump must be provided with control switchboard.



In case of pumps as Dürr or Cattani equipped with control switchboard, which require a simple switch between two preset contacts (generally 1 and 3), connect such contacts to the two brown ones **6** and **7** in front of the white wires, as the figure shows.

H - ARRANGEMENT OF ISOSUN LAMP FOR ISOTRON, ISOPLUS, ISOPLUS LITE

The dental unit shall be set according to the type of Isosun lamp installed, i.e. with photo cell or without photo cell. In order to carry out such programming, switch on the dental unit by keeping the keys  and  pressed: the display will show one of the two following arrangements:


 control for Isosun lamp without photo cell

 control for Isosun lamp with photo cell

By keeping the keys  and  pressed, press the key  to set the  arrangement (Isosun with photo cell), or the key  to store  (Isosun without photo cell).

Release the keys to confirm the selected arrangement and to restore the usual operating order.


I - ISOCONCEPT ADJUSTMENTS

PRESSURE ADJUSTMENT TO INSTRUMENTS

The pressure of the air feeding instruments is adjusted with screw **A** placed on the proper unit. Turning rightwards the pressure decreases, leftwards the pressure increases.

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