Photometer-System CheckitDirect

CSB • COD • DCO • COD • DOC

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Seite 2–9

(GB) Instruction Manual

Page 10–17

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Page 18–25

(IT) Istruzioni d'uso

Pagina 26–33

(ES) Instrucciones

Página 34–41

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Display backlight

Press the [!] key to turn the display backlight on or off. The backlight is switched off automatically during the measurement.

Recall of stored data

If the instrument is switched on, press the [!] key for more than 4 seconds to access the recall menu.

Measuring the absorbance

AbS

This Mode function allows measurement of the absorbance at the specified wavelength. The result appears in the display in mAbs.
1000 mAbs = 1 A (Absorbance unit)

Abs1: $\lambda = 430 \text{ nm}$

Abs2: $\lambda = 605 \text{ nm}$

Method notes

- Prior to measurement ensure that the sample is suitable for analysis (no major interferences) and does not require any preparation i.e. pH adjustment, filtration etc.
- Reagents are designed for use in chemical analysis only and should be kept well out of the reach of children.
- Ensure proper disposal of reagent solutions.
- Material Safety Data Sheets: available on request.

Chemical method notes:

Method:

The organic material present in the sample is oxidised by a standard amount of a potassium dichromate oxidising mixture. After oxidation is complete, the excess of this reagent is measured photometrically.

Application:

Samples can be measured if the chloride content does not exceed 1000 mg/l (LR/MR) or 10 000 mg/l (HR).

In exceptional cases, compounds contained in the water cannot be oxidized adequate. This results in minimum findings, compared with the reference method.

Different methods of sampling, preparation of the sample itself and the time elapsed between taking the sample and analysis can all affect the obtained results.

CAUTION

The accuracy of the instrument is only valid if the instrument is used in an environment with controlled electromagnetic disturbances according to DIN 61326. Wireless devices, e.g. wireless phones, must not be used near the instrument.

Photometer COD vario

Select the appropriate vial for the desired range:

LR-Range: 0 – 150 mg/l, Order code 2420720

MR-Range: 0 – 1500 mg/l, Order code 2420721

HR-Range: 0 – 15000 mg/l, Order code 2420722

Preparing the Sample


Open a reaction vial with a white cap and add the specified volume.
(Ensure that appropriate personnel safety equipment is used.)

LR-/MR: 2 ml water sample

HR: 0.2 ml water sample

Prepare a **blank** (Note 1) by using deionised water (TOC-free) instead of the sample (LR/MR: 2 ml, Hr: 0.2 ml).

Replace the cap **tightly**. Invert the vial gently several times to mix the contents (**The vial will become hot during mixing!**) and digest the vials for **120 minutes** in the reactor at a temperature of **150°C**. Remove the vials from the reactor and allow them to cool down to 60°C or less. Mix the contents by inverting each vial several times while still warm. Then allow the vials to cool to ambient temperature before measuring.

Lr

RANGE

RANGE

0.0.0

Measurement

Fix the adapter for 16 mm vials on the sample chamber.

Switch the unit on using the [ON/OFF] key.

The display shows the following:

Select the required test using the [MODE] key:

LR → MR → HR → AbS1 → AbS2 → LR →..... (Scroll)

The display shows the following:

Place the **blank** in the adapter (Note 1–4) making sure that the marks Δ are aligned. Place the cover on the adapter.
Blanks are specially prepared for each individual test range.

Press the [ZERO/TEST] key.

The "Range" symbol flashes for approx. 8 seconds.

The display shows the following:

After zeroing remove the vial from the adapter.

Place the **sample** in the adapter (Note 2–4) making sure that the marks Δ are aligned. Place the cover on the adapter.

Press the [ZERO/TEST] key.

The "Range" symbol flashes for approx. 3 seconds.

The result appears in the display.

LR-/MR-range: in mg/l

HR-range: in g/l

The result is saved automatically.

Tolerance: ± 3.5 % (full scale)

Repeating the test:

Press the [ZERO/TEST] key again.

New zero calibration:

Press the [MODE] key until the desired range symbol appears in the display again.

Guidelines for photometric measurements

1. Run samples and blanks with the same batch of vials.
The blank is stable when stored in the dark and can be used for further measurements with vials from the same batch.
2. Don't place hot vials in the adapter. Allow the vials to cool to room temperature for minimum 45 minutes. It is recommended to leave the vials to cool over night.
3. Suspended solids in the vial lead to incorrect measurements. For this reason it is important to place the vials carefully in the adapter. The precipitant at the bottom of the sample should not be suspended.
4. Clean the outside of the vials with a towel. Finger-prints or other marks will be removed.
5. Avoid spillage of water or reagent solution into the sample chamber because this can lead to incorrect test results.
6. Contamination of the optics (LED and photo-sensor) in the sample chamber will result in incorrect measurements.
The windows of the sample chamber must be checked at regular intervals and cleaned as necessary. Use a moistened cloth and cotton buds for cleaning purposes.
7. Large temperature differences between the photometer and ambient conditions can lead to incorrect measurements, or build-up of condensate around the optics of the vial.
8. To avoid errors caused by stray light do not use the instrument in bright sunlight.

Mode

On
Off

!

Mode

Menu selections

Press the [MODE] key and **hold**.

Switch the unit on using the [ON/OFF] key.

Allow the 3 decimal points to be displayed before releasing the [MODE] key.

The [!] key allows for selection of the following menu points:

- ▲ recall stored data
- ▲ ▽ setting the date and time
- ▽ user calibration

The selected menu is indicated by an arrow in the display.

Confirm the selection with the [MODE] key.

▲ Recall of stored data

The photometer shows the last 16 data sets in the following format (automatically proceeds every 3 seconds until result is displayed):

Number n xx (xx: 16...1)
Year YYYY (i.e. 2008)
Date mm.dd (monthmonth:dayday)
Time hh:mm (hourhour:minute) minute)
Test Method
Result x,xx

The [ZERO/TEST] key repeats the current data set.

The [MODE] key scrolls through all stored data sets.

Quit the menu by pressing [!] key.

Zero
Test

Mode

!

Mode

SET

DATE

YYYY

(2 sec.)

Mode

Zero
Test

!

2 3 Setting date and time (24-hour-format)

After confirming the selection with the [MODE] key the value to be edited will be shown for 2 sec.

The setting starts with the year (YYYY) followed by the actual value to be edited. Same applies for month (mm), day (dd), hour (hh) and minutes (mm). Set the minutes first in steps of 10, press the [!] key to continue setting the minutes in steps of 1.

Increase the value by pressing the [MODE] key.

Decrease the value by pressing [ZERO/TEST] key.

Proceed to the next value to be edited by pressing [!] key.

After setting the minutes and pressing the [!] key the display will show "IS SET" and instrument returns to the measurement mode.

ATTENTION: If the battery is removed for more than one minute the date and time menu starts automatically when the photometer is switched on the next time.

4 User calibration

Note:

user calibration (Display in calibration mode)

factory calibration (Display in calibration mode)

cAL

CAL

CAL

Lr

Zero
Test

RANGE

0.0.0

CAL

Zero
Test

RANGE

RESULT

CAL

Mode

Zero
Test

CAL

RESULT + x

On
Off

:

After confirming the selection with the [MODE] key the instrument will show CAL/Lr.

Scroll through ranges using the [MODE] key.

Place the blank in the adapter making sure that the marks Δ are aligned.

Press the [ZERO/TEST] key.

The "Range" symbol flashes for approx. 8 seconds.

The display shows the following in alternating mode:

Perform calibration with a standard of known concentration (see "Operation").

Press the [ZERO/TEST] key.

The "Range" symbol flashes for approx. 3 seconds.

The result is shown in the display, alternating with CAL.

If the reading corresponds with the value of the calibration standard (within the specified tolerance), exit calibration mode by pressing the [ON/OFF] key.

Changing the displayed value:

Pressing the [MODE] key once increases the displayed value by 1 digit.

Pressing the [ZERO/TEST] key once decreases the displayed value by 1 digit.

Press the corresponding key until the reading equals the value of the calibration standard.

By pressing the [ON/OFF] key, the new correction factor is calculated and stored in the user calibration software.

Confirmation of calibration (3 seconds).

Recommended calibration values

LR: 100 mg/l O₂

MR: 500 mg/l O₂

HR: The instrument calibration in the range MR is automatically taken also for HR.

Factory calibration reset

Resetting the user calibration to the original factory calibration will reset all methods and ranges.

A user calibrated method is indicated by an arrow while the test result is displayed.

To reset the calibration press both the [MODE] and [ZERO/TEST] key and **hold**.

Switch the unit on using the [ON/OFF] key.

Release the [MODE] and [ZERO/TEST] keys after approx. 1 second.

The following messages will appear in turn on the display:

SEL

CAL

The factory setting is active.
(SEL stands for Select)

or:

SEL

cAL

Calibration has been set by the user.

(If the user calibration is to be retained, switch the unit off using the [ON/OFF] key).

Mode

SEL

CAL

Calibration is reset to the factory setting by pressing the [MODE] key.

The following messages will appear in turn on the display:

On
Off

Switch the unit off using the [ON/OFF] key.

Technical Data

Light source:	Lr-Range: LED, Filter ($\lambda = 430 \text{ nm}$) Mr-/Hr-Range: LED ($\lambda = 605 \text{ nm}$)
Battery:	9 V-block battery (Life approx. 600 tests, without display light)
Auto-OFF:	Automatic switch off 7 minutes after last keypress
Ambient conditions:	5–40°C 30–90% rel. humidity (non-condensing).
CE:	Certificate for Declaration of CE-Conformity at www.lovibond.com

To ensure maximum accuracy of test results, always use the reagent systems supplied by the instrument manufacturer.

Operating messages

Hi	Measuring range exceeded or excessive turbidity.
Lo	Result below the lowest limit of the measuring range.
	Replace 9 V battery, no further tests possible.
btLo	Battery capacity is too low for the display backlight; measurement is still possible.

Error codes

E 01	Light absorption too great. Reasons: e.g. dirty optics.
E 10/E 11	Calibration factor "out of range"
E 20/E 21	Too much light reaching the detector
E 22	Battery capacity was too low during measurement. Change battery.
E 70	Lr: Factory calibration incorrect / erase
E 71	Lr: User calibration incorrect / erase
E 72	Mr: Factory calibration incorrect / erase
E 73	Mr: User calibration incorrect / erase

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are Trademarks of the
Tintometer AG
5212 Hausen AG