GARANTIEZERTIFIKAT

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GUARANTEE CERTIFICATE

**GB** Faro offers the final customer a **24 month** guarantee starting from the date of purchase. Repairs under guarantee must be performed at FARO; expenses and transport risks are at the risk of the purchaser.  
**Repair under guarantee is considered valid only when:**  
- **all sections of the certificate have been filled in and sent in advance to FARO by Fax (039.6010540).**

The guarantee covers faults due to the bad quality of the material or manufacturing defects; in the case of valid claims, the guarantee covers free repair or replacement. **Claims for damages and/or interest are excluded.** The guarantee is not considered valid, at the sole discretion of FARO, if the fault is due to tampering, damage, incorrect use, improper maintenance and normal wear and tear.

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## CONTENTS

|                                 |            |
|---------------------------------|------------|
| Safety regulations.....         | page 23-24 |
| Characteristics.....            | page 25-26 |
| Description of the parts* ..... | page 27    |
| Installation.....               | page 28    |
| Instructions for use.....       | page 29    |
| Maintenance.....                | page 35    |
| Guide to problems.....          | page 41    |
| Technical specifications .....  | page 42    |

\* The pictures are shown on the last few pages

## SYMBOLS

### **DANGER**

Paragraphs marked with this symbol contain instructions that must be followed carefully to prevent damage occurring to the device, the operator and possibly the patient.

### **WARNING**

These instructions warn that great care must be taken to avoid situations that could damage the device.

### **PROHIBITION**

This icon highlights what must be avoided to prevent damage occurring to the device.

### **SUGGESTIONS**

With this icon, information is supplied that enables the device to be used in a more efficient way.

Disposal container. Dispose the device observing the standards about the differentiated collection of electrical material.

### **PLEASE NOTE**

## SAFETY REGULATIONS

Follow all instructions given in this manual.

- The autoclave should be used by staff with specific training in the process of sterilization within doctor's surgeries, podiatrist's studios and beauty centres, under the responsibility of the health director.

- The device is to be used to sterilize, with steam at a temperature of 121°C or 134°C, tools, rotating instruments, reusable medical devices and porous materials (e.g. fabrics).

- The autoclave must be installed in an especially hygienic environment.

- The device must be connected to the mains.

- The electrical system must comply with IEC 364-1 standard and the "Italian regulations for installation of electrical systems in areas assigned for medical use".

- Check that the supply voltage, given on the ratings plate on the rear panel of the autoclave, corresponds to that of your mains. If in doubt, contact your retailer.

- The autoclave plug is earthed. Ensure that you use a socket with an earth connection and position the device so that it can be easily disconnected.

- When using an extension cable, ensure that the type of cable used is suitable for the absorption of the equipment.

- Before using the device for the first time, clean the chamber (as described on page 42).

- Do not rest the autoclave on unstable tables or surfaces.

- Do not attempt to perform maintenance on the autoclave when the power is on; disconnect the power supply from the mains before intervening on the machine.

- If, during a sterilization cycle, an error message should appear, **the cycle must be repeated.**

- Remove the trays with the pincers provided to avoid being burnt.

- Accidentally dropping the device may cause it to be damaged and not work correctly; if this occurs, it should be examined by technical assistance.

- The autoclave works with distilled water (see technical specifications on page 42).

**Any other substance used could compromise the device's operation.**

**If an inappropriate liquid is introduced, consult technical assistance.**

- The device in its packaging can withstand the following environmental conditions for a period no greater than 15 weeks:

Ambient temperature from -20°C to +70°C

Relative humidity from 10% to 90%

Atmospheric pressure from 500 to 1060 mbar.

## **SAFETY REGULATIONS**

**The equipment must be used at the following environmental conditions:**

- indoor use
- altitudes up to 3000 m
- temperature from 5° to 40°C
- relative humidity max 85%
- max variation of mains voltage  $\pm 10\%$
- installation category (overvoltage category) II
- pollution level I

**Do not perform any maintenance operations other than those described in the manual. Any operation not indicated therein could compromise the safety aspects provided by the machine. For all operations not specified, consult the technical assistance department.**

## CHARACTERISTICS

### THE DEVICE IS DESIGNED TO STERILIZE USING WET SATURATED STEAM.

The vacuum pump enables an optimal process of sterilization to be achieved with a correct interaction between the steam and the processed load, even in the cavities.

### STERILIZATION PROGRAMMES

- |                  |                              |
|------------------|------------------------------|
| - 121°C or 134°C | <b>cycle type B standard</b> |
| - 134°C          | <b>cycle type B Prion</b>    |
| - 134°C B Flash  | <b>cycle type B fast</b>     |

### TEST PROGRAMMES

- |                                 |                                            |
|---------------------------------|--------------------------------------------|
| - <b>Vacuum leak test:</b>      | tests the ability to maintain the vacuum   |
| - <b>Bowie &amp; Dick test:</b> | tests steam penetration in the porous load |
| - <b>Helix test 121°C</b>       |                                            |
| - <b>Helix test 134°C:</b>      | tests steam penetration in hollow loads    |

### AUTOMATIC LAST ADVANCED CYCLE OPERATION

At the end of the cycle, with the door closed, the autoclave goes into standby after 30 minutes. Furthermore, when the temperature in the chamber reaches 60° C, the drying phase is repeated.

### WATER FILLING SYSTEMS

**Automatic:** with the aid of a pump found inside the autoclave.

**Manual:** alternatively the device is equipped with a manual filling system.

Possibility of connecting a demineralizing system to increase the length of operating time.

### AUTOMATIC CONTROL OF THE STERILIZATION PARAMETERS

The parameters of pressure, temperature and time are continuously monitored during the sterilization phase by means of an automatic control system.

### BIO-X FILTER

Inputs sterile air into the chamber during the drying phase and in the pressure balance phase.

### USB CONNECTION

A computer can be connected to monitor machine operation and set the parameters (this must be done by an authorized technician).

### LCD

A menu system ensures complete management of the autoclave.

## CHARACTERISTICS

### TEST METHODS PERFORMED ON THE SK07 AUTOCLAVE

Test methods in compliance with standard EN13060 B type cycles

|                                               |
|-----------------------------------------------|
| Seal                                          |
| Dynamic pressure of the sterilization chamber |
| Vacuum chamber                                |
| Solid load                                    |
| Reduced porous load                           |
| Full porous load                              |
| Hollow load A                                 |
| Multiple packaged                             |
| Dryness, solid load                           |
| Dryness, porous load                          |

## DESCRIPTION OF THE PARTS

### SEE FIG. 1

1. On/off switch
2. Handle for opening-closing door
3. LCD
4. Pushbutton (\\)
5. Pushbutton (←)
6. Pushbutton (-)
7. Pushbutton (+)
8. USB socket
9. Bio-X filter
10. Integrated printer compartment/memory card slot/compartment for external printer

### SEE FIG. 2

11. Tank discharge
12. Waste water tank filling
18. Handles

### SEE FIG. 3

13. Clean water tank discharge
14. Manual water filling
15. Demineralizing water-fill coupling
16. Water discharge to the external can

### SEE FIG. 4

17. Door adjustment push-button
19. Inspection window

### SEE FIG. 5

Minimum overall dimension of the housing (mm)

### SEE FIG. 6

Dimensions required when opening the door 3 (mm)

## INSTALLATION

**THE DEVICE MUST BE INSTALLED BY SPECIALIZED TECHNICIANS**

On the base of the autoclave packaging are holding points.

Brief handling operations of the equipment can be carried out by two people; the holding points for lifting the autoclave are positioned on the upper part of the equipment (18 fig. 2). For longer handling operations, the pallet supplied should be inserted under the equipment.

- Take the autoclave out of its packaging.
- Position the autoclave onto its support (table or unit) with the aid of the handles.
- Remove the handles and insert the cover plugs. Keep the handles and their screws for future handling.
- Open the door by turning the knob 2 in an anti-clockwise direction (fig. 1) and remove the accessories supplied from the chamber.

**Check that the equipment includes:**

- Documentation (instruction manual)
- Tray support
- 3 trays
- 1 tray extraction pincer
- 1 filling pipe
- 1 draining pipe
- 1 funnel
- 4 plugs

**Check with a spirit level that the resting surface is level in both directions; if not, level the autoclave using the adjustable front feet.**

**If the autoclave is to be built in to the unit, the dimensions in fig. 5 must be observed.**

**Check that the power plug used is earthed.**

Plug in the autoclave plug taking care that the cable does not remain squashed and that it can easily be removed from the socket.

## INSTRUCTIONS FOR USE

To use the autoclave, proceed in the following way:

- Open the door by turning the knob (2 fig.1) in an anti-clockwise direction.
- Press the main luminous switch (1 fig. 1).

### USING THE MENU

The main menu consists of 3 submenus:

- 1 - Programs
- 2 - Test
- 3 - Options

To access the various submenus, move the selection indicator by pressing the (+) or (-) key and then confirming the selection with the  key.

Press the (\) key to cancel the operation.

### SWITCHING ON FOR THE FIRST TIME

When the device is switched on for the first time, the message “Add H<sub>2</sub>O” will appear on the display. **In order to carry out the sterilization cycles, the clean water tank must be filled as described below.**

### FILLING THE CLEAN WATER TANK

- Use only distilled water (see technical specifications on page 42).
- When the message “Add H<sub>2</sub>O” appears on the LCD, fill the clean water tank in one of the following ways:

#### AUTOMATIC FILLING

- Insert the part of the pipe with the quick coupling into the water-fill coupling (12 fig. 2) and place the other end in the can.
- Press the (+) button to start the water pump.
- Press the (-) or (\) key to stop the operation.
- Filling automatically stops once the maximum level of the tank has been reached.

#### MANUAL FILLING

- Insert the part of the pipe with the quick coupling into the hole (14 fig. 3).
- Place the funnel into the other end of the pipe and pour max. 2 liters of distilled water.

### INTRODUCING MATERIAL TO BE STERILIZED INTO THE CHAMBER

Package the load very carefully ensuring that the packaging chosen is compatible with the medical devices to be sterilized.

**The material used must be:**

- Compatible with the sterilization process.
- Compatible with the labeling system.
- The packaging material must comply with standard EN 868.

**Place the thoroughly cleaned material on the trays, distributing it uniformly.**

**The chamber can contain a maximum of:**

- 2 Kg** of packaged material
- 4 Kg** of unpackaged material
- 1.3 Kg** of porous load

The packaged items must be placed on the trays with the transparent part of the packet facing downwards.

## SELECTING THE STERILIZATION PROGRAMMES

The preheating time from cold is 10 minutes.

By selecting the submenu **“Programs”**, the sterilization cycles that the autoclave is able to carry out can be accessed:

- 1 - 134°C Standard type B cycle**
- 2 - 134°C Prion type B cycle**
- 3 - 121°C Standard type B cycle**
- 4 - 134°C B Flash type B cycle\***

The selection indicator is positioned on “134°C Standard B” cycle.

To scroll the cycles and move the selection indicator, press the (+) or (-) keys, then confirm with the  key.

Press the “\” key to cancel the operation.

## STANDARD CYCLES

- Cycle at 134°C : 4 minutes of sterilization plus 14 minutes of drying
- Cycle at 134°C B Flash\* : 3,5 minutes of sterilization plus 8 minutes of drying
- Cycle at 121°C : 16 minutes of sterilization plus 14 minutes of drying

\* Maximum load: 0,5 Kg of packaged material – 0,2 Kg of porous material..

NOTE: For each available sterilization cycle must be done by manufacturer the Works Test as specified by EN 13060 (par. 7.3.1). The available cycles for that the manufacturer hasn't done the Works Test must be performed by Installation Tests following the EN 13060 par. 7.4.

## PRION CYCLE

- Cycle at 134°C : 18 minutes of sterilization plus 14 minutes of drying

## “LAST ADVANCED CYCLE” FUNCTION




The “Last Advanced Cycle” function is enabled for all programs.

30 minutes after the end of sterilization cycle, the heating elements switch off and a **second drying cycle starts automatically** when the temperature inside the chamber reaches 60°C.





## MODIFYING STERILIZATION TIMES

The sterilization times can be increased; this operation requires a specialized technician.

## CYCLES DIAGRAMS

## CYCLES DIAGRAMS

## EXAMPLE OF PRINTED VALUES

### PRINTING SEQUENCE:

**Logo Faro** or **FARO S.p.A.**

Ornago (Mi)

**SNXX STXXXX**

← Autoclave serial no.

**Program:**

← type of cycle or type of test

**Cycle nr XXXX**

← cycle no.

**Start:**

← date and time of cycle start

DD/MM/YY HH:MM:SS

**PRESS TEMP TIME**  
(kPa) (°C) (H:MM:SS)

← unit of measurement  
(time calculated from the start of the cycle)

### Fractional Vacuum

XXXX XXXX X:XX:XX  
+++

← Pressure, temperature and time of 1<sup>st</sup> vacuum phase

← Top-ups during pressure increase

XXXX XXXX X:XX:XX  
+++

← Pressure, temperature and time of 2<sup>nd</sup> vacuum phase

← Top-ups during pressure increase

XXXX XXXX X:XX:XX  
+++

← Pressure, temperature and time of 3<sup>rd</sup> vacuum phase

← Top-ups during pressure increase

### Sterilization

XXXX XXXX X:XX:XX

← Pressure, temperature and start time each min of sterilization

XXXX XXXX X:XX:XX

← Pressure, temperature and start time each min of sterilization

XXXX XXXX X:XX:XX

← Pressure, temperature and start time each min of sterilization

XXXX XXXX X:XX:XX

← Pressure, temperature and start time each min of sterilization

XXXX XXXX X:XX:XX

← Pressure, temperature and start time at end of sterilization

+++

← Top-ups in sterilization

**P<sub>MAX</sub>=XXXXkPa T<sub>MAX</sub>=XXXX°C** ← Maximum pressure and temperature during sterilization

**P<sub>MIN</sub>=XXXXkPa T<sub>MIN</sub>=XXXX°C** ← Minimum pressure and temperature during sterilization

### Discharge

X:XX:XX

← Discharge phase start time

### Dry

X:XX:XX

← Drying phase start time

### Press. Balance

X:XX:XX

← Pressure balance phase start time

### End:

X:XX:XX

← Total time of cycle

dd/mm/yy hh:mm:ss

← Date and time cycle end

**OK**

← Formatting bold and expanded type, larger typeface

**Operator signature**

## **TEST PROGRAM SELECTION**

By selecting the submenu “test”, the list of test cycles that the autoclave performs can be accessed:

- 1 - Vacuum Leak Test
- 2 - Helix Test 134°C
- 3 - Helix Test 121°C
- 4 - Bowie&Dick Test
- 5 - Quit

The selection indicator is positioned on “Vacuum Leak Test”.

Select the desired program, moving using the (+) and (-) keys and confirm with .

### **VACUUM LEAK TEST**

In order to perform this function, the chamber must be cold and without instruments.

### **BOWIE & DICK TEST**

In order to perform this function, the chamber must be cold and without instruments.

The Bowie & Dick “pack” for checking the porous load must be placed on the lower tray.

### **HELIX TEST**

In order to perform this function, the chamber must be cold and without instruments.

The control tool for examining the **“Hollow A” load must be placed on the lower tray.**

## MAINTENANCE

- **Periodic service:** this activity shall be performed by the end user.  
*Frequency:* daily, weekly or monthly.
- **Yearly service:** this activity shall be performed by a qualified technician from FARO After-Sales Technical Support Centre.  
*Frequency:* yearly or on expiry of a certain number of cycles.

### PERIODIC SERVICE

- Vacuum Test
- Penetration Test (Bowie & Dick / Helix Test)
- Cleaning of the chamber and accessories (tray holder and trays)
- porthole cleaning
- External frontal panel cleaning
- porthole gasket cleaning
- Chamber filter cleaning
- porthole regulation
- Bacteriological filter replacement
- porthole gasket replacement

#### Vacuum Test

*Check:* daily.

#### Penetration Test - Bowie & Dick Test

*Check:* weekly (suggested). Refer to national regulations and instructions.

#### Penetration Test - Helix Test

*Check:* weekly (suggested). Refer to national regulations and instructions.

 **Before proceeding with the operations described below, always turn off the main electrical power supply switch (1 Fig. 1) and ensure that the chamber is cold.**

#### Cleaning of the sterilisation chamber and of internal accessories

*Check:* monthly or every 50 cycles.

Use a soft cloth soaked in water to avoid damaging the protective surface coat of the chamber. Rinse thoroughly with water and, on completion of the operation, dry.

However, clean the chamber only in case of extreme need.

If there are white calcium deposits, immediately check the reason for their presence in to chamber. Formation of white patches at the base of the chamber reveals the use of either poor quality water or of water that lacks the required characteristics.

#### Porthole cleaning

*Check:* daily or every 5 cycles.

Carry out the procedures described below with a soft cloth soaked in water to avoid damaging the protective coat of the hatch. If necessary, moisten the cloth with white vinegar to remove calcium scale.

Rinse thoroughly with water and, on completion of the operation, dry.

During the cleaning process, DO NOT USE either metallic or highly abrasive cloths as the latter could cause abrasion of the external porthole surface.

## External frontal panel cleaning

*Check:* monthly or every 50 cycles.

Use a soft cloth soaked in alcohol, taking care don't let alcohol penetrate into the parts inside the autoclave.

Any ammonia-based glass cleaner is tolerated as an alternative to alcohol, but it is absolutely forbidden to use detergents based on trichloroethylene, benzene, turpentine or solvents in a broad sense.

## Porthole gasket cleaning

*Check:* daily or every 5 cycles.

Clean the gasket with a soft cloth moistened with either water or white vinegar to remove calcium scale.

This procedure must be performed to remove impurities that can cause the sterilisation chamber to lose pressure and can even produce cuts in the gasket.

Avoid external elements (hair, dust or the like) from remaining deposited on the contact surface of the porthole gasket.

## Chamber filter cleaning

*Check:* monthly or every 50 cycles

- Extract the tray holders with the trays from the chamber.
- Rotate the filter upward inside the chamber, and draw it towards you to extract it.
- Immerse in a container filled with water and vinegar, wait a few minutes and then use a small brush to clean the filter area.
- Rinse the filter unit with water.
- Reassemble by applying the same procedure in reverse order.
- Reposition the filter at the bottom of the chamber.

In case of failure to restore initial conditions, call Technical Support for a replacement.

## Porthole door regulation

*Check:* when necessary. Generally not before six months after the porthole gasket has been changed.

The autoclave closure system requires occasional adjustments due to the normal misalignment of mechanical parts and of the porthole gasket.

Imperfect sealing on the part of the gasket could prevent correct performance of sterilisation cycles with impairment of the expected outcome.

## Activating always when the machine is turned off and cold.

*Instructions:*

- Keep button 17 pressed. (Fig. 4)
- Rotate porthole 19 (Fig. 4) in an anticlockwise direction for a quarter of a circle.
- Release the button and keep rotating porthole 19 (Fig. 4) until it locks.
- Ensure that the door closes smoothly. If the handle is too hard to close, slightly rotate it in the opposite direction, keeping the button pressed.
- Run a test cycle to check efficacy of the implemented adjustments.

This procedure must usually be performed at the first slight leak or sealing problem of the gasket. Before tightening the porthole closure system, clean the hatch gasket as described above, and repeat the sterilisation cycle.

## Bacteriological or Bio X filter replacement

*Check:* quarterly or every 500 cycles.

*Instructions:*

- Open the door of the printer compartment/CF Card of the autoclave.
- Unscrew filter 9 (Fig. 1) by rotating it anticlockwise.
- Screw it clockwise, ensuring it is fully screwed on.

**Porthole gasket replacement**

*Check:* yearly or every 1000 cycles

*Instructions:*

- Open the sterilising chamber door.
- Manually remove the gasket by grasping the lip with two fingers (observe how the lip of the seal is fitted).
- Clean the site of the gasket with a cloth moistened with alcohol or with a cotton-bud.
- Moisten the new gasket with water or alcohol.
- Fit the gasket in 6 points and then use one finger to arrange the parts of the gasket inside the site.
- Insert the new gasket in its site with lips turned outward, arranging it evenly around the circumference, using the fingers with constant pressure along the entire perimeter of the gasket.
- On completion of the procedure, check that no points have been incorrectly inserted or that the gasket presents no protrusion. In case of protrusions, remove the gasket and perform the insertion process once again.

**The gasket must be correctly positioned to avoid leakage. It must be inserted along the entire circumference without either protrusions or deformations.**

**YEARLY SERVICE**

After 1000 cycles or one year after installation of the autoclave.

On reaching the 1000th cycle, a memo message will appear on performing service. Service includes replacing, cleaning and checking components, as specified in the table below.

Yearly service shall only be performed by a qualified FARO technician from a FARO After-Sales Technical Service Network.

The intervention is already defined, and includes:

**Replacement:**

- porthole gasket
- bacteriological Bio X filter
- chamber filter
- water filter
- vacuum pump membrane kit (if necessary)

**Checks:**

- safe door closure system
- porthole door regulation
- safety valve
- water loading/unloading tank
- hydraulic connections
- electrical connections
- pressure and temperature calibration with SW
- electrical test protection
- effective vacuum pump function.

**On completion of the intervention, the qualified technician will issue a certificate stating that autoclave yearly service was performed.**

## SELECTING THE OPTIONS:

By selecting the “Options” submenu, the list of submenus for configuring the machine can be accessed.

- 1 - Configuration
- 2 - Water filling
- 3 - Cycle printing
- 4 - Memory deletion
- 5 - Quit

### Submenu Configuration

By selecting the “Configuration” submenu, the list of configuration options for the machine can be accessed:

- 1 - Basic
- 2 - Advanced
- 3 - Special
- 4 - Quit

### Basic selection:

By selecting the “Basic” submenu, the list of machine configuration can be accessed:

- 1 - Language
- 2 - Date
- 3 - Time
- 4 - Quit

### Advanced selection:

Access to this menu is only permitted by an authorized technician.

By selecting the “Advanced” submenu, the list of machine configurations can be accessed:

- 1 - Water filling options
- 2 - Storage
- 3 - Quit

### Special selection:

Access to this menu is only permitted by an authorized technician.

By selecting the “Special” submenu, the list of machine configuration can be accessed:

- 1 - Ambient pressure
- 2 - LCD Contrast
- 3 - Menag Password
- 4 - Quit

## MAINTENANCE:

Maintenance must only be performed by authorized technicians.

## DISPLAY SYMBOLS

Maximum level  
of discharge tank

Minimum level  
of filling tank

Alarm

Sterilization  
OK

Door lock on

Cycle in  
progress

Heaters  
active

Cycle end  
OK

## SETTING THE CLOCK

### **Adjusting the Date and Time:**

From the main menu select “Options”.

From the options menu select “Configuration”.

From the configuration menu select “Basic”.

### **To change the date:**

From the basic menu select “Date”:

Press the “+” key  
to increase the date.

Press the “-” key  
to decrease the date.

Press the “Enter” key  
to move and select.

Select “Save” and press the “Reset”  
key to store the settings.

### **To change the time:**

From the basic menu select “Time”:


Press the “+” key  
to increase the date.

Press the “-” key  
to decrease the date.

Press the “Enter” key  
to move and select.





Select “Save” and press the “Reset”  
key to store the settings.

## CHANGING THE PAPER ROLL

To change the paper roll proceed as follows:

1. Open the cover of the printer by pressing the OPEN button as shown (**fig. 1**).
2. Position the paper roll inside the printer, making sure the roll rotates in the direction indicated (**fig. 2**).
3. Pull the paper until it protrudes out of the printer and close the cover (**fig. 3**). The cover will lock into place.
4. Tear off any excess paper along the serrated edge (**fig. 4**).
5. The printer is ready to print.

## GUIDE TO PROBLEMS

### ERRORS DURING THE CYCLES

| CODE  | MESSAGE                     | DESCRIPTION                                                                                      |
|-------|-----------------------------|--------------------------------------------------------------------------------------------------|
| ER001 | <b>Safety Cut Off</b>       | the safety temperature of the chamber or band has been exceeded                                  |
| ER002 | <b>Top-up</b>               | the maximum number of top-ups permitted has been exceeded                                        |
| ER003 | <b>Pressure balance</b>     | the pressure balance has not been reached                                                        |
| ER004 | <b>No Ready</b>             | the ready temperature has not been reached                                                       |
| ER005 | <b>Timeout 1</b>            | failure to reach final pressure during thermodynamic discharge                                   |
| ER006 | <b>Vacuum</b>               | the second vacuum threshold has not been reached                                                 |
| ER007 | <b>Timeout 2</b>            | the sterilization conditions have not been reached                                               |
| ER008 | <b>chamber T High</b>       | chamber temperature is above the limit allowed                                                   |
| ER009 | <b>chamber T Low</b>        | chamber temperature is below the limit allowed                                                   |
| ER010 | <b>High Pressure</b>        | the pressure is above the limit allowed                                                          |
| ER011 | <b>Low Pressure</b>         | the pressure is below the limit allowed                                                          |
| ER012 | <b>High theoretic T</b>     | the theoretical temperature is above the limit allowed                                           |
| ER013 | <b>Low theoretic T</b>      | temperature below the limit allowed                                                              |
| ER014 | <b>Delta Temperature</b>    | the difference between chamber temperature and theoretic temperature is beyond the limit allowed |
| ER015 | <b>Vacuum Test: vacuum</b>  | vacuum threshold not reached during vacuum leak test                                             |
| ER016 | <b>EN 13060</b>             | parameters not respected in first phase of vacuum leak test                                      |
| ER017 | <b>EN 13060</b>             | parameters not respected in second phase of vacuum leak test                                     |
| ER018 | <b>EN 13060</b>             | failure to observe temperature limits during vacuum test                                         |
| ER019 | <b>Door open</b>            | Door not closed                                                                                  |
| ER020 | <b>Timeout water charge</b> | failed to fill the charge tank                                                                   |
| ER021 | <b>Emptying water tank</b>  | failed to empty the discharge tank                                                               |

### MACHINE FAILURES

| CODE  | MESSAGE                     | DESCRIPTION                                   |
|-------|-----------------------------|-----------------------------------------------|
| AL001 | <b>Chamber Probe Fault</b>  | chamber probe open                            |
| AL002 | <b>Chamber Probe Fault</b>  | chamber probe short-circuited                 |
| AL003 | <b>Band Probe Fault</b>     | band probe open                               |
| AL004 | <b>Band Probe Fault</b>     | band probe short-circuited                    |
| AL005 | <b>Timeout water charge</b> | exceeding the maximum time for charging water |
| AL006 | <b>Pressure Transducer</b>  | pressure transducer fault                     |
| AL007 | <b>Door Micro</b>           | unforeseen door micro contact opening         |
| AL008 | <b>Safety Micro</b>         | unforeseen door safety micro contact opening  |
| AL009 | <b>Blackout</b>             | blackout occurred                             |
| AL010 | <b>Safety Pressure</b>      | pressure above safety limit                   |
| AL011 | <b>Corrupted Memory</b>     | memory corrupted                              |
| AL012 | <b>Safety Micro</b>         | Door safety off                               |
| AL013 | <b>PUSH Error</b>           | Software error                                |
| AL014 | <b>POP Error</b>            | Software error                                |
| AL015 | <b>Reset Pressed</b>        | Reset button pressed                          |
| AL016 | <b>Safety Micro-switch</b>  | Door safety on                                |

Machine release code:

Press the keys “←”, “+” and “-” simultaneously, and wait for the end of the procedure. If the machine has broken down, contact Faro.

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