• LTM630

NEOPRO 10000 MANUALE D'ISTRUZIONE

NEOPRO 10000 INSTRUCTION MANUAL

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The Serie Neo is in conformity with the Medical Devices Directive (93/42/EEC) and has been subject to the conformity assurance procedures laid down in the Council Directive.

Codes LTM630

Thank you for purchasing the mattress antidecubitus and pump **LEVITAS** by Moretti, which have been manufactured according to the highest quality and safety standards (93/42/CEE).

IMPORTANT SAFEGUARDS READ ALL INSTRUCTIONS BEFORE USING

Use: for antidecubitus treatment.

1.1 General Information

Neo-series systems use most advanced technologies and mattress materials to deliver the best quality, comfort and whisper-quiet operations to you and your customers. Most systems in Neo-series are recommended for prevention and treatment of pressure ulcers up to and including stage IV pressure ulcers.

The Neo-series range of systems have been independently tested and successfully approved to the following standards:

EN 60601-1/ UL 2601
EN 60555-2 Class A
IEC 1000-4-4

EN 60601-2
EN 60555-3
IEC 1000-4-5

EN 55011 Class B
EN 50141
IEC 1000-4-2

1.2 Intended Use


This product is an advanced alternating pressure relief replacement system suitable for high-risk patients. This low air-loss system is intended to help and reduce the incidence of pressure ulcers while optimizing patient comfort. It is recommended for prevention and treatment of pressure ulcers up to stage IV.

2. The Product

The digital controlled pump unit, which provides adjustable pressure and cycle time, is compact and effective in pressure control and management. The replacement mattress, which offers a low air loss therapy, provides patients a comfortable pressure support suiting for each individual. The patented cell-in-cell design allows the inflation and deflation of air cells to be achieved in the shortest amount of time while providing the good pressure dispersion. The inner cell also offers a sufficient weight support during a power failure (the cell-in-cell bladder design is only available in models of 10" and 8" mattress).

2.1 Product Description

1. Mattress
2. Pump unit
3. Air tube connector
4. Control panel

3. Installation

3.1 Unpacking

The pump unit and mattress are packaged in a separate box to secure its contents inside. Unpacking these boxes to remove the pump unit and mattress and check for any damage, which may have occurred in shipping. If there are damages, please contact your dealer immediately.

3.2 Setting Up

1. Place the mattress on top of the bed frame. Please note for the foot end. Mattress can be secured firmly by fixing the elastic straps at the bottom bed frame.

Warning! Please cover the mattress with a cotton sheet to avoid direct skin contact and for patient's comfort.

2. Hang the pump onto bed rail (foot-end), and adjust a suitable wideness of hangers.

3. Connect air tube connector from air mattress to the pump unit. When "click" is felt or heard, the connection is completed and Secured.

Warning! Check and ensure the air tubes are not kinked or tucked under mattress.


4. Plug into mains electrical outlet.

Warning! Before inserting the plug into the mains, make sure it is suitable for the local main voltage.

Warning! Put the plug into socket of grounded AC power.

5. Then turn the main power switch ON. The power switch is located on the right side panel of the pump unit.

Warning! Ensure the CPR connectors are sealed completely.

Grounding

Before any connection to the output connectors is made, the unit shall be connected to grounding conductor via the three-core main cable; the grounding is only a "functional grounding".

Installation tips are listed below

After installation, the extra wire of the power cord, if any, should be neatly arranged to avoid any tripping accidents. The equipment should be firmly placed at position where users/doctors can easily access.

4. Operation

4.1. Function Description

Control Button *Function Description*

Power

To turn on control unit. It takes approximately 15-min to inflate the mattress. The default settings are 5-min alternating cycle time, level 1 comfort level, 10-min flat level pause interval, and 4-min turn level pause interval.

Alarm/Mute

Enable and disable the audible alarm feature. Press this button to silent the audible alarm. The alarm buzzer will alert again after 3-min if the pressure has not yet returned to normal condition.

Panel Lock

To lock the current setting and prevent other persons from changing them. To unlock the panel, press this button and hold it for 3 seconds.

Low Pressure Indicator

3 Conditions:

1. No light - normal condition; mattress is operated accordingly.
2. Yellow light - pressure has not reached the preset setting.
3. Flashing yellow light - low pressure condition.

Autofirm

For maximum inflation of mattress. All other functions are disabled unless this button is pressed again. This mode will last for 20-min and then it will return to default setting automatically.

Static

For non-alternating mode of mattress. To go back to alternating mode, simply press this button again (such as regular air mattress).

Seat Inflate

The Body zone's pressure will be increased higher than the preset setting to give a proper support for purposes of user's sitting up and easiness of getting in and out of the mattress. This function should be disable when user is in lying position to prevent too high pressure support at body zone (only available in 10" and 8" mattress models).

Comfort Range

6 Comfort ranges of mattress. Use UP/DOWN button to adjust the comfort levels. 1 is the softest & 6 is the firmest.

Cycle time

4 Alternating cycle time. Use UP/DOWN button to adjust the cycle time of 5, 10, 15, 20 minutes.

4.2 Operating Instruction

- 1 Press the POWER button on the display panel to start the system.
- 2 For faster inflation, the AUTOFIRM function can be activated. This mode will last for 20 minutes and then it will return to default setting automatically.
- 3 After the mattress is fully inflated, user can change the setting according to their preference or recommended by physician.

4.3 CPR

When there is an emergency to perform CPR on the patient, quickly pull the CPR valves to release air from mattress. The CPR valves are located at the head-end, right-hand-side of the mattress. The air tube connector (located on the pump unit indicated with "PULL CPR" sign) can be disconnected for even faster deflation.

4.4 Alternating Function Setup

User can adjust pressure of air mattress by adjusting the comfortable range from level 1 to 6. The system also provides various alternating cycle time mode, 5, 10, 15, 20 minutes. Please consult your physician for a suitable setting. When the Low Pressure indicator is off, it indicates the mattress has reached to its preset softness. Patient can lie on the mattress.

Warning! Check to see if the suitable comfort range (or pressure setting) is selected by sliding one hand between the deflated air cell and the patient to feel the patient buttock. Users should be able to feel the minimal contact.

4.5 Intubation Function

To activate the emergency intubation, simply disconnect the connector located at the head section to quickly deflate the first 4 air cells. The whole mattress will still function normally except at this region. To make mattress back to normal, simply secure the connector and wait for a few minutes for inflation. (only available in models of 10" and 8" mattress).

4.6 Seat Inflate Function

When it is necessary for patient to sit on the bed (or for easiness of getting in and out of mattress), the seat inflate function provide extra pressure compensation in body region's cells for more suitable surface support. (only available in models of 10" and 8" mattress) This function should be disable when user is in lying position to prevent too high pressure support at body zone.

4.7 Heel Protection Function

The last 5 air cells of mattress are interchangeable with a set of 5" and 18" air cells to provide a proper heel protection. The 5" air cell provides the heel protection and 18" air cell provides the foot stopper. These integrated air cells help patient to have proper position while reducing the interface pressure at the heel region (only available in 10" mattress model).

4.8 Low Pressure and Alarm Function

When an abnormal low pressure is lasted for a long period, the Low Pressure Indicator (flashing yellow LED) will light up with audible alarm. Check if all the connections are securely made and that they are correctly installed as per installation instructions.

Warning! If the pressure level is consistently low, check for any leakage at air cells or connecting

tubes. If necessary, replace any damaged cells or tubes or contact local qualified dealer for repair. If the control unit is equipped without audible alarm feature, the Low Pressure Indicator will flash until the low-pressure fault condition is solved. If the control unit is equipped with audible alarm feature, the Low Pressure indicator will flash and the alarm will sound when air pressure is below normal (preset pressure). To silent the alarm, simply press the alarm mute button, and it will be silent for 3 minutes; however, the Low Pressure Indicator still remains flashing. After 3 minutes, alarm buzzer will alert again if the pressure has not yet returned to normal.

5. Cleaning

Warning! Ensure the pump unit is disconnected from the mains before cleaning.

5.1 Pump Unit

Wipe with a damp cloth and a mild detergent and keep it away from dust. If other detergent is used, choose one that will have no chemical effects on the surface of the plastic case of the pump unit. All parts should be air dried thoroughly before use.

Warning! Do not immerse or soak pump unit.

5.2 Mattress

Wipe down with warm water containing a mild detergent. The cover may also be cleaned using sodium hypochlorite diluted in water. All parts should be air dried thoroughly before use. Do not use Phenolic based products for cleaning.

For suitable and locally available cleaner please ask your distributor.

Warning! Dry the mattress on SUNLESS capitalized area after cleaning.

6. Storage

- Pull the CPR tags and make sure all cells are deflated.
- Lay the mattress out flat and upsides down.
- Protect the air tube coupler by putting it in the inner side of mattress.
- Roll the mattress from the head end towards the foot end.
- Foot-end straps can then be stretched around the rolled mattress to prevent unrolling.

7. Maintenance


7.1 General

- Check main power cord and plug if there is abrasions or excessive wear.
- Check mattress' cover for signs of wear or damage. Ensure mattress cover and tubes are stubbed together correctly.
- Check the air tube connections for any damage.
- Report any problems or malfunctions to local distributors.

7.2 Air Filter Replacement

The air filter should be checked and dusted periodically.

- Snap out the plastic cover, located on the back of the pump unit.
- The filter is reusable and can be washed gently with a mild detergent and water.
- Suggest to replace air filter regularly if the environment is dirty. This will help to improve the airflow. Spare filters are available from the local distributors.

7.3 Fuse Replacement

1. Disconnect the plug from mains power when a fuse is to be replaced.
2. Remove the cover of the fuse holder by means of a small screwdriver.
3. Insert a new fuse of the correct rating in, and replace the cover of the fuse holder back. Check the technical specification for correct fuse rating.

8. Troubleshooting

Problem*

Power is not ON

Alarm is on (audible & visual)

Patient is bottoming out

Mattress becomes too firm and pressure can't be lower down

Mattress form is loose

Air tube connector can't fit into the pump unit

No air produced from some air outlets of the air tube connector

Solution

- Check if the plug is connected to mains.
- Check if there is any blown fuse.
- Check if the connection between air tube connector to pump unit is tightly secured.
- Check if all tubing connections along mattress are secured.
- Check if the CPR is sealed.
- Check if there are any leakage on air cells.
- Pressure setting might be inadequate for the patient, adjust comfort range 1 to 2 levels higher and wait for a few minutes.
- The AUTOFIRM is activated for too long, release some air by disconnecting the air tube connector and change to your preferred setting.
- Check if all the snap buttons and straps of mattress on all air cells are all securely fastened.
- Check if the mattress is fixed to the bed frame by elastic straps.
- Press handles on both sides of air tube connector slightly before push the connector into the pump unit.
- Make sure the connector is tightly secured by hearing "click" while pushing it into pump unit.
- This is normal since there are alternating and rotating mode. Air outlets take turn to produce air during their preset cycle time.

* If your questions can't be answered with above information, please contact your local agent directly. They might require a technician to take care the problem.

9. Technical Description

Pump	Specification
Power supply	AC100-120V,50/60Hz AC220-240,50/60Hz
Fuse Rating	1A T1A
Dimension (LxWxH)	37 x 27x 13 (cm)or 14.4" x 10.5" x 5.0"
Weight	6.0 kg or13.2 lb
Environment	<ul style="list-style-type: none"> • Temperature Operation: 10℃~50℃ Storage:-10℃~70℃ Shipping: -10℃~70℃ • Humidity Operation: 30%~ 75% non-condensing Storage: 10%~100% non-condensing • Atmosphere 700 hPa ~1060 hPa
Classification:	<ul style="list-style-type: none"> • Class II, Type BF Applied Parts Mattress • Not suitable for use in the presence of a flammable anesthetic mixture • IPX0 Enclosed equipment without protection against ingress of water • Continuous operation
Mattress	Specification
• Dimension (inflated, LxWxH)*	200 x 90x 25.4(cm) 78.7" x 35.4"x 10"
• Weight*	10 kg or22 lb

* The weight and dimension of mattress is measured without foam.

Warning! The specifications, also apply to those areas operating with the same power supply.

"BF" symbol, indicates this product is according to the degree of protecting against electric shock for type BF equipment.

The Symbol for class II equipment, 2 concentric squares indicating double insulation.

Attention, should read the instructions.

Enclosed equipment without protection against ingress of water.

Disposal of Electrical & Electronic Equipment (WEEE): This product should be handed over to an applicable collection point for the recycling of electrical and electronic equipment. For more detailed information about the recycling of this product, please contact your local city office, household waste disposal service or the retail store where you purchased this product.

Important safeguard

Warning! To reduce the risk of electrocution:

- Always unplug this product immediately after using.
- Do not use while bathing.
- Do not place or store product where it can fall or be pulled into a tub or sink.
- Do not place in or drop into water or other liquid.
- Do not reach for a product that has fallen into water. Unplug immediately.

Warning! To reduce the risk of burns, electrocution, fire, or injury to persons:

- This product should never be left unattended when plugged in.
- Close supervision is necessary when this product is used by, on, or near children or invalids.
- Use this product only for its intended use as described in this manual.
- Do not use attachments not recommended by the manufacturer.
- Never operate this product if it has a damaged cord or plug, if it is not working properly, if it has been dropped or damaged, or dropped into water. Return the product to a service center for examination and repair.
- Keep the cord away from heated surfaces.
- Never block the air openings of this product or place it on a soft surface, such as a bed or couch, where their openings may be blocked. Keep the air opening free of lint, hair, and other similar particles.
- Never drop or insert any object into any opening or hose.
- Connect this product to a properly grounded outlet only.
- This device is Class II equipment with functional earth, no Protective Earth is provided.

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
WARRANTY

All kit mattress antidecubitus and pump Serie Neo sets are guaranteed for **24 months** against any possible manufacturing fault. The warranty covers neither damages, perforations or bad functioning caused by mishandling or faulty installation nor parts subject to wear and tear. Postage must be paid and the appliance must be sent together with the present guarantee coupon duly filled in.

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Acquistato in data / Purchasing date _____

Rivenditore / Authorized dealer _____

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Venduto a / Purchased by _____

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Distribuito da/Distributed by: **Moretti S.p.A.**