

R E S C U E

everytime
needs... someone

:: kong **LECCO**
870.00

KONG s.p.a.
via XXV Aprile, 4 – (zona Industriale)
23804 MONTE MARENZO (LC) - ITALY
TEL. + 39 0341 630506
FAX. + 39 0341 641550
www.kong.it

INFORMATIONS GIVEN BY THE MANUFACTURER

LECCO STRETCHER

1 - INFORMATIVE NOTE on CE certified Kong products conforming to 89/686/EEC dir.

Users must carefully read and understand the informative note and instructions for use before using the product.

Warning: these manufacturer's informations are not exhaustive and do not pretend to teach neither the practise of use and the procedures.

Practical use and procedures, above all in natural environments, must be learned/thought by qualified, trained and skilled instructors.

Warning: rock and ice-climbing, double-cord descent, klettersteig, speleology, alpine skiing, canoeing, exploration, rescue work, tree climbing and working at heights are all highly dangerous with possibilities of even lethal accidents. Train yourself to use this article and make sure you fully understand how it works, if you have any doubts, don't run any risks and ask.

Remember that:

- this article shall be used by a trained and or otherwise competent person or the user should be under the direct supervision of such a person,
- you are personally responsible to know and understand this article, how to use it and safety measures to be taken, and that you alone are solely and totally liable for all damages, injuries or death eventually incurred by you or third parties from using all KONG Spa products, of any whatsoever kind.

Avoid using this article if you are not able to undertake this responsibility and take these risks.

2 - GENERAL INFORMATION

2.1 - Carry out all the checks described in the specific "instructions for use" for all the products, before use and be particularly sure that the article is:

- in top condition and working correctly,
- suited for the use you intend to make of it: use the techniques shown without crosses only, all other usage is absolutely forbidden: risk of death!

2.2 - Immediately replace the product if you have any doubts as to its safety and efficiency. Do not use the product after falling as its strength could be considerably impaired by eventual invisible breakage or deformations inside.

Incorrect use, mechanical deformation, equipment accidentally dropped from a height, wear, chemical contamination, exposure to heat above normal climatic conditions (max 100°C for purely metallic articles and max. 50°C for articles with parts in fabric), are just some examples of other reasons possibly reducing, limiting and even terminating the service life of the article.

2.3 - Users are responsible for using this article correctly.

2.4 - This article may be used in combination with individual protective devices conforming to the 89/686/EEC Directive, compatibly with relevant instructions for use.

2.5 - The resistance of natural or not natural anchorages in the rocks is not *a priori* guaranteed, it is therefore essential that users carefully judge for themselves beforehand to ensure adequate protection.

2.6 - For safety's sake,

- users shall evaluate all the possible risks and be sure that the whole safety system (anchor point, attachments, safety line, harness, etc. etc.) of which this equipment is a component, is reliable and safe, in addition to this, with systems blocking a fall check:
- whenever possible the anchorage point is above the user,
- use a full body harness,
- organise the work with the purpose of reducing to a minimum any potential fall and its heights.

2.7 - The anchorage position is essential for safely blocking a fall: carefully judge the height of the fall, length of cord and "pendulum" effect to avoid all possible obstacles (e.g. ground, articles scraping against the rocks, etc.).

2.8 - Your life depends on the continual efficiency of your equipment (we thoroughly recommend personal use of equipment) and how it is kept (use, storage, controls, etc.). We strongly recommend having pre and post-use controls carried out by a qualified person. An accurate overhaul by a manufacturer's qualified technician at least once a year is also recommended. The user is responsible for recording and conserving the information given in the "Record" table together with the instructions.

2.9 - KONG Spa shall not be held liable in any way for damages, injuries or death caused by: improper use, changes made to product, repairs by non-authorized persons or use of non-original spare parts.

2.10 - No special precautions should be taken for transportation, do however avoid contact with chemical reagents or other corrosive substances and protect all the pointed, sharp parts adequately. Warning: never leave your equipment in cars parked in the sun!

2.11 – Warning: it is compulsory for dealers to have these instructions for use translated for products sold in countries other than those selected by KONG Spa.

3 - MAINTENANCE AND STORAGE

No particular maintenance has to be done by the user, the article just has to be cleaned and lubricated, as explained below.

3.1 - Cleaning: frequently rinse the product under lukewarm fresh water (max. 40°C), eventually adding a delicate detergent (neutral soap). Leave it to dry naturally away from direct sources of heat.

3.2 - Lubrication: frequently lubricate the moving parts with silicone based oil. Avoid contact between oil and fabric parts. This operation must be done after cleaning and thoroughly drying.

3.3 - Disinfection: submerge the product for an hour in lukewarm water adding disinfectant containing salts of quaternary ammonium, then rinse under fresh water, dry and lubricate.

3.4 - Storage: after cleaning, drying and lubricating store the equipment loosely in a dry, cool, dark (avoid U.V. rays), chemically neutral (absolutely avoid saline environments) place, away from sharp corners, sources of heat, humidity, corrosive substances or other potentially harmful conditions. Do not store when wet!

4 - TESTING AND CERTIFICATION

This product has been certified by the notified register organization under n. 0123 - TUV Product Service GMBH, Ridlerstraße 65, 80339 Munchen, Germany in conformity to the norms given in the marking.

All KONG products are tested/inspected piece by piece in conformity to the certified UNI EN ISO 9001 standard Quality System procedures. The notified register organization n. 0426 - ITALCERT, v.le Sarca 336, 20126 Milan, Italy, supervises the manufacture of class III individual protection devices, in conformity to article 11B, 89/686/EEC Directive.

Warning: laboratory tests, testing and inspection, instructions for use and standards are not always able to reproduce the case in practice, so that results achieved with the article under effective conditions of use in natural environments may sometimes differ to a considerable degree. The best instructions can be had from continuous use under the supervision of qualified skilled instructors.

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N.B.: The cloth bag accepts vacuum-operated mattress or any other blocking system for an injured person.

5.1 CLOTH BAG (INTERNAL HOOKING BELTS):

If a rescue board or a vacuum-operated mattress is not used to block the injured person, the cloth bag is equipped with several internal hooking belts (fig. 01 - 02) for immobilization of injured.

Special Velcro® strips allow the correct fixation of the rescue board or of the vacuum mattress inside the cloth bag.

The head-rest pillow is fixed by Velcro® strips (fig. 03), and removable if necessary (i.e.: when a rescue board is used).

5.2 CLOTH BAG (EXTERNAL HOOKING BELTS):

The external hooking belts assure the optimal fixing of injured (fig. 04) and avoid wind effects on the bag (i.e.: rescue by helicopter).

Two double zippers allow the complete opening of the cloth bag or only a partial lateral control (fig. 05) without completely uncover the injured.

WARNING!

Even if the zippers are not self-locking, do not open the cloth bag by pulling the cover; use the pullers on the zipper cursors (fig. 06).

5.3 HAND SHAFTS:

The two pairs of shafts (two front shafts and two rear shafts) are fitted with a comfortable padding:

- o The front shafts can be inserted in straight position for normal carrying (fig. 07) or in opposite position for extremely steep slopes or overhangs, giving a better control of stretcher drive (fig. 08).
- o The shape of the rear shafts allows a good visibility of the ground (fig. 09).

5.4 WORKING LOAD:

The stretcher is designed and tested to support a max weight evenly distributed of 720 Kg.

The safety working load (SWL) of the cloth bag is calculated with a safety factor 1:4 to support a weight evenly distributed of 180 Kg. **The load shall never be applied in one only point.**

The hanging kits (for helicopters or cableways) have a SWL of 500 Kg (Safety factor 1:4 – breaking load 2000 Kg).

The use of Kong Stretcher in hard conditions may cause an adjustment of the connections between heads and longitudinal bars with a play between the metallic parts: this is not prejudicial to the functionality of the stretcher provide it is correctly assembled with the lateral and longitudinal webbings correctly stretched (fig. 23 – 27).

5.5 LIFE TIME OF THE STRETCHER:

- o The life time of the stretcher frame is technically unlimited, provided it is correctly stored and that after an intensive use (we suggest at least every two years) the stretcher is revised by Kong s.p.a. or persons authorized by Kong s.p.a.
- o The life time of the cloth bag, if correctly stored, is of ten years, provided that after an intensive use (we suggest at least every two years) the cloth is revised by Kong s.p.a. or persons authorized by Kong s.p.a.
- o Hanging kits “Orion” – “Gemini” – “Evolution” and “Virgo”, should be replaced every three years (even if not used).

6. ASSEMBLING INSTRUCTIONS:

6.1 HOW TO ASSEMBLY THE METAL FRAME

The longitudinal bars with the same coloured adhesives shall be inserted each into the other (fig. 10); the click of safety pins (fig. 11) will show the correct assembling of the two components.

- o lay on the ground the two longitudinal bars assembled (with three slots on the same side) and the head (follow the markers: red for the head part – close to the three slots; blue for the feet part) (fig. 12);
- o insert a longitudinal bar in the head, checking the right marking and pushing the safety pin (fig. 13). The click of safety pin (fig. 14) will show the correct assembling.
- o insert one of the three spacing bars in the holes of longitudinal bars; insert the longitudinal bar in the head and check the click of safety pin (fig. 15);
- o insert the other two spacing bars in the holes of longitudinal bars and connect the other head (feet side - blue markers) (fig. 16). Check the click of safety pins in the holes of longitudinal bars (fig. 17);

IMPORTANT WARNING!

- o **The safety pins are made only to assembly the metal frame; the strength of the stretcher is ONLY guaranteed by the correct assembly and stretching of the hooking belts (lateral and longitudinal) of the cloth bag.**

6.2 HOW TO FIX THE CLOTH BAG ON THE METAL FRAME:

- o Unfold the cloth upside-down with the lower part upwards and the hooking belts lying on the ground (fig. 18);
- o Put the metal frame on it, with the three hooking belts in correspondence to the three slots of the longitudinal bars (fig. 19);
- o Insert the hooking belts through the slots of the longitudinal bars (from external to internal); the buckle and ring shall be inserted through the slots one after the other (fig.20) (for safety reasons slots are narrow in order not to allow the simultaneous passage of buckle and ring) (fig.29);
- o Fasten without stretching the four transversal hooking belts (details of the buckle clasping system in fig. 21 - 22 – 23);
- o The longitudinal hooking belts overturn the transversal bar of the heads; fasten the first buckle (fig. 24). Stretch the webbings paying attention that the angle between longitudinal and lateral belts shall be 90° (fig. 25);
- o stretch definitely the transversal hooking belts (fig. 23) fixing the end part around the ring (fig. 26) in order to avoid accidental opening;
- o stretch definitely the longitudinal hooking belts, fastening the second safety buckle, then fix the end part as done for transversal hooking belts (fig. 27);

WARNING!

- o The buckles of transversal hooking belts are alternated to grant a connection between cloth and frame in case of accidental release (fig. 28);
- o Due to the thickness of the buckle and ring together compared to slot width, is not possible an accidental complete release (fig. 29);
- o Don't knot the end of hooking belts but fix them in the way shown in figure 26.

6.3 HOW TO ASSEMBLY THE SHAFTS:

- o to fix the shafts (fitted with removable washable padding), insert the end into the hole of head, pulling the safety pin (fig. 30);
- o the shafts will be correctly inserted when you heard the clutch of the locking pin after releasing – check also visually;

6.4 SHAFTS ASSEMBLING POSSIBILITIES IN DESCENDING TRANSPORT:

Normally an injured should be carried with the head higher than the feet, (but in certain cases doctors may decide the opposite).

- o for normal down carrying, the two front shafts (feet side) can be inserted in straight position (fig. 31). The two rear shafts (head side) shall be inserted in the position shown in fig. 32;
- o for extremely steep slopes or overhang, the front shafts (feet side) can be inserted in the position shown in fig. 33. The rear shafts (head side) shall be inserted in the same position as before (fig. 32);
- o the normal position of shafts (for descent - with the head of injured higher than feet) is indicated by numbers marked on the heads and shafts (fig. 34);

6.5 SHAFTS ASSEMBLING POSSIBILITIES IN ASCENDING TRANSPORT:

- o shafts shall be inserted in opposite heads (fig. 35) to keep the head of injured higher than feet.

WARNING!

- o Shafts are not conceived for hard lateral strains: therefore, in the case of such stresses, they can have a permanent deformation – in such cases do not attempt to make them straight but change them.

7. ACCESSORIES:

7.1 RESCUE BOARD (# 855.05):

Doctors say that the complete immobilization of a person with spine injuries is possible only using a spine rescue board.

KONG S.P.A. has developed a rescue board foldable in three pieces, made of a special treated wood washable and disinfected, easy to transport, with a special pillow to block the head of injured.

The rescue board is radio-transparent for a first X-ray diagnosis without need to move the injured (fig. 36 - 37).

7.2 "TAYLAN" BOARD BINDING (# 855.02):

The complete immobilization of an injured is possible with the special adjustable board binding (fig. 38).

7.3 SAFETY WEBBING "VIRGO" (# 871.02):

For lowering or lifting the stretcher on steep slopes or overhang, it is compulsory the use of the safety webbing "Virgo" to be inserted in and around the frame (better with empty stretcher) and then firmly connected to the belaying ropes.

HOW TO ASSEMBLY "VIRGO":

- o Insert one Virgo's end into the head of the stretcher (head side) (fig. 39), then pull it between the spacing bars of the frame and the cloth bag (fig. 40);
- o Pull out Virgo from the feet head as shown in fig. 41;
- o insert Virgo under longitudinal hooking belts (feet side)(fig. 42);
- o repeat the same operations through the other side of the frame (fig. 43).
- o pull Virgo until the two ends are aligned (fig. 44);
- o Connect both belaying ropes to both Virgo's end loops (fig. 45).

WARNING - DANGER!

Virgo shall go all around all the stretcher: don't pass Virgo only in the rear head part (head side) otherwise the strength of the structure is not guaranteed.

7.4 HEAD-REST PILLOW WITH LEXAN® VISOR (# 871.30)

Designed to improve protection of injured in case of rain, hail, debris and recovery with helicopter.

HOW TO ASSEMBLY LEXAN VISOR:

Insert the shaped visor into the lateral pockets of the pillow.

Fix it with the two lateral Velcro® strips (fig. 46).

The pillow is equipped with a special Velcro® front strip to retain the head of injured (fig. 47).

7.5 BOTTOM PROTECTION (# 871.40)

Designed to protect the back of the injured and to make easier the sliding/pulling of the stretcher in particular conditions (canyoning, rough, stones, etc...).

HOW TO ASSEMBLY THE BOTTOM PROTECTION:

The assembling of the two sheets shall be made with the stretcher empty.

- o Overturn the stretcher and place sheet n°1 (feet side) (fig. 48). Push it towards the feet head until the two blocks ("L" shaped) engage the spacing bars (fig. 49);
- o Lay it on the central spacing bar and pull it back towards the stretcher centre to engage the central "L" shaped blocks (fig. 50 e 51)
- o Lay sheet n°2 (head side) - it shall be overlapping sheet n°1 - push it towards the centre to engage the blocks on the central spacing bar (fig. 52); pull it back until the two blocks engage the external spacing bar (fig. 53);
- o Now the four "U" shaped blocks of the two sheets shall be in correspondence to the vertical lateral bars of the heads;
- o Insert the two long "U" bars on to the blocks (fig. 54 e 55);
- o Secure the "U" long bars inserting the split pins supplied (fig. 56).

7.6 HAND SHAFTS (# 871.50 – kit of 2 pairs)

Hand shafts allow two rescuer to drive Kong Stretcher when fitted with a single and/or double wheel (fig. 57).

This shafts can be used also for shoulder transport by four people.

7.7 SINGLE WHEEL (# 871.20) AND DOUBLE WHEEL (# 871.22)

Designed to make easier the injured transport on flat and not rough route. (fig. 58 e 59).

HOW TO ASSEMBLY WHEEL(S)

- o Insert the hand shafts (# 871.50) into the heads (procedure is the same as shoulder shafts);
- o Lift the stretcher carrying it close to the head (fig. 60) and insert the profiles under the longitudinal bars;

N.B. it might be necessary to slightly enlarge or to narrow the two longitudinal bars

- o It is possible to fix the wheel(s) in any position to find the correct barycentre (it is normally located where shown in fig. 61);
- o Screw or unscrew the threaded hook to find the better snap position (fig. 62);
- o The high of the wheel(s) can be varied from mm 460 up to mm 570 unblocking the spring pins on the fork sides (fig. 63);

WARNING!

O A stretcher fitted with one wheel shall be driven and controlled by two rescuers.

It is also possible a combined installation of a double wheel (preferably on the back side) and of one single wheel (preferably front side).

N.B. to make a bent, it is necessary to lift the single wheel

7.8 HOW TO USE HANGING KIT "ORION" (# 871.00):

Hanging kit Orion has been studied to lower or lift Kong Stretcher by winches, fishing poles, tripods, helicopters, etc...

- Fix hanging kit to stretcher by inserting the four screw sleeves connectors stitched on the four arms; insert the connectors from outside towards inside first onto the head and then into the webbing loop stitched on the cloth bag (fig.64). Follow reference colours (red: head side -

blue: feet side).

Fully close the connector screw sleeves

- During the manoeuvres the rescuer shall hook himself to the suspension ring above the central bar and not over the webbing (fig.65).

7.9 HOW TO USE HANGING KIT "EVOLUTION" (# 871.04):

Hanging kit Evolution, made of two separate arms, has been conceived to make easier the use of Kong Stretcher when the injured is suspended in vertical or over hanging position. Better than Orion in case of a sturdy injured.

The fixing of Evolution arms is identical to that of Orion (red with red – blue with blue – connectors from outside to inside, etc...).

To lift (or to lower) the stretcher use a twist lock connector (# 962) to connect the two central rings. It is advisable to interpose between the hook and the Evolution a "Ciapin" (# 274S – multi-anchor static webbing) that minimizes the risks caused by rigid connections and rationalize the positioning of different loads (stretcher, first rescuer, second rescuer, safety longes, etc...) (fig. 67).

If the "Ciapin" is not used it is advisable to insert both rings into the lifting hook – in this last case the twist lock connector it's only for safety (fig. 66).

7.10 HOW TO USE "GEMINI" HANGING KIT FOR CABLEWAYS (# 871.01)

Hanging kit Gemini has been studied to keep the stretcher parallel to the bearing cable.

The fixing of Gemini arms is identical to that of Orion (red with red – blue with blue – connectors from outside to inside).

- The best sliding of Kong Stretcher fitted with Gemini is obtained by using two double pulleys Pamir (# 826)(fig. 68).
- Otherwise two "Turbo Roll" pulleys (# 821) may be used as an alternative, connecting them to the two rings of the Gemini via a screw sleeve connector Oval Alu Classic (# 730) as shown in fig. 69
- In case of "critical" environment, it is advisable to install two static bearing ropes and the use of two pulleys "Heavy duty" (# 819) inserted onto both ropes (fig. 70).

WARNING!

Hanging kit Gemini shall not be used for vertical lowering or lifting – in that occasions use only either Orion or Evolution.

7.11 BLOCKING KIT HELICOPTER'S CABIN FLOOR (# 871.03)

The kit is composed of four adjustable webbings fitted with a quick link and a stainless steel quick release connector. It allows to fix the stretcher to the cabin floor of helicopters or other vehicles (fig. 71). The pulling of the red webbing makes possible the quick release of the connectors (fig. 72).

WARNING!

BEFORE THE FIRST USE OF THIS PRODUCT IT IS COMPULSARY TO READ AND UNDERSTAND INSTRUCTIONS AND FILL THIS CONTROL SHEET

KONG LECCO			
Series number:		Year of manufacture:	
User's name:		Buying place:	
Date of first use:		Buying date:	
Control date	Remarks on correction operations	OK (Y/N)	Signature

KIT DIMENSIONS




A comfortable bag includes all the kit!

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