

757L20

| | |
|--|----|
| DE Gebrauchsanweisung | 3 |
| EN Instructions for use | 6 |
| FR Instructions d'utilisation | 9 |
| IT Istruzioni per l'uso | 12 |
| ES Instrucciones de uso | 15 |
| PT Manual de utilização | 19 |
| NL Gebruiksaanwijzing | 22 |
| SV Bruksanvisning | 25 |
| DA Brugsanvisning | 28 |
| PL Instrukcja użytkowania | 31 |
| HU Használati utasítás | 35 |
| CS Návod k použití | 38 |
| EL Οδηγίες χρήσης | 41 |
| ZH 使用说明书 | 45 |

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English

These instructions for use were last updated on 2015-03-19.

Before operating the Li-Ion Battery Charger, carefully read these operating instructions as well as the operating instructions for the Ottobock Energy-Pack 757B20/757B21. Please pay particular attention to the Safety Precautions!

1 Application

The Li-Ion Battery Charger 757L20 is to be used **solely** for charging one or two Ottobock EnergyPacks 757B20/757B21.

2 Safety precautions

Attention!

- Use the Charging Station only in combination with the Power Supply Unit 757L16-* or with the Charging Cable 4X74 12V.
- Use the battery charger under dry conditions only.
- Do not use the battery charger if any damage is visible. If you discover any damage, send the battery charger immediately to an Ottobock service facility.
- Do not let foreign particles or liquids get into the battery charger.
- Do not expose the battery charger to excessive heat (avoid direct sunlight), excessive smoke, mechanical vibrations or shock.

These products may not be disposed of with household waste in some jurisdictions. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the information provided by the responsible authorities in your country regarding return and collection processes.

3 Scope of delivery

- 1 pce. Charging Station
- 1 pce. Power Supply Unit with EU plug
- 1 pce. US plug
- 1 pce. Operating Manual

4 Accessories

For using the battery charger with the 12V cigar lighter of a car, the Charging Cable 4X74 is optionally available.

5 Description and function

Unlike NiCad batteries, the rechargeable li-ion battery has the benefits of a lower self-discharge, a higher cell potential and a higher energy density. Li-ion technology offers higher capacities, has no memory effect and requires considerably reduced charging times.

The Li-Ion Battery Charger 757L20 consists of a Charging Station (fig. 1, pos. 1), a Power Supply Unit 757L16-* with EU plug (fig. 1, pos. 2), and an additional US plug (fig. 1, pos. 3). For exchanging the plug, please refer to figure 4.

The battery charger has an operating indicator (green LED, fig. 2, pos. 3) as well as four other LED indicators:

- Green LEDs (fig. 2, pos. 1 and 1a): Indicate the charging status of the EnergyPacks
- Yellow/red LEDs (fig. 2, pos. 2 and 2a): Indicate malfunctions of the Energy Packs or battery charger

Before charging, the battery charger conducts a self-test.

In a test procedure, short-circuited or damaged cells are identified and prepared for the CCCV process (constant current - constant voltage) during the conditioning phase of the charging process.

During the CCCV process, the EnergyPacks are charged with a continuously high current until a certain voltage level is reached. The current then automatically decreases to a lower level until it is switched off completely.

If the voltage of the EnergyPacks falls below a certain level, the battery charger automatically switches to trickle charge.

The charging process will also be interrupted when the temperature exceeds 50 °C/122 °F or falls below 0 °C/32 °F. This is indicated by the yellow LED. The charging process will resume as soon as the EnergyPack has reached room temperature.

6 Function

6.1 Self-test

Connect the Charging Station to the Power Supply Unit 757L16-* (see fig. 2, pos. 4) and plug into an electrical outlet. The operating indicator lights up (see fig. 2, pos. 3) showing the battery charger is ready for use.

Next the battery charger conducts a self-test, during which all other indicator lights will light up briefly.

6.2 Charging

Insert one or two EnergyPacks into the charging station (fig. 3).

The charging process begins when the indicator lights1 and/or 1a light up.

The charging time for a depleted EnergyPack is approximately 3-1/2 hours. As soon as the EnergyPack is fully charged, the 1/1a indicator light(s) flash and the charging current is cut off.

If the EnergyPack voltage falls below a certain value due to self-discharge, the Li-Ion Battery Charger automatically switches to trickle charging. The fully charged battery may be left in the battery charger without becoming damaged.

Notice:

To avoid malfunctions, the electrical contacts of the EnergyPack and of the battery charger must be kept clean and free of grease. For cleaning please use a dry cotton pad stick. Clean the electric contacts with a cotton bud and some alcohol on a regular basis.

6.3 Trouble-shooting tips

- The li-ion battery charger is plugged into an electrical outlet, but the operating indicator is not lit:
→ **Re-test charger function at a different electrical outlet.**
 - If the operating indicator still does not light up, the li-ion battery charger is defective and must be sent to an Ottobock service facility.
 - If the indicator lights 2 and 2a alternately flash red, the battery charger is defective and must be sent to an Ottobock service facility.
 - The li-ion battery charger is plugged into an electrical outlet, the operating indicator is lit, and two EnergyPacks are inserted, but only one of the indicators 1 or 1a is illuminated:
→ **Switch the locations of the EnergyPacks:**
- a) If the same LED does not light up, the battery charger is defective and must be sent to an Ottobock service facility for repair.
- b) However, if the other LED now fails to light, the corresponding EnergyPack is defective and must be sent to an Ottobock service facility for repair.

Other messages:

| Battery Compartment | Indicator 1/1a green | Indicator 2/2a yellow | Indicator red | Function | Remarks |
|---------------------|----------------------|-----------------------|---------------|----------|--|
| empty | off | off | off | okay | stand-by |
| empty | off | off | flashing | error | self-test – return charger |
| with EnergyPack | off | off | flashing | error | self-test – return charger |
| with EnergyPack | on | off | off | okay | charging |
| with EnergyPack | blinking | off | off | okay | EP is fully charged |
| with EnergyPack | off | on | off | error | temperature of EP > 50 °C: let cool down (see Section 5) |
| with EnergyPack | off | off | on | error | EP is defective – return EP |
| empty | off | on | off | error | charger defective – return charger |
| empty | off | off | on | error | charger defective – return charger |
| empty | on | off | off | error | charger defective – return charger |

7 Technical data

Charging Time..... approx. 3-1/2 h
 Power supply..... 100 - 240 V
 Power supply frequency47 - 63 Hz
 Operating temperature0 - 50 °C (32-122°F)

8 Liability

Otto Bock Healthcare Products GmbH, hereafter referred to as manufacturer, assumes liability only if the user complies with the processing, operating and maintenance instructions as well as the service intervals. The manufacturer explicitly states that this device may only be used in combination with components that were authorised by the manufacturer (see instructions for use and catalogues). The manufacturer does not assume liability for damage caused by component combinations which it did not authorise.

The device may only be opened and repaired by authorised Ottobock technicians.

9 CE conformity

The product meets the requirements of guideline 2004/108/EC for electromagnetic compatibility. The declaration of conformity was therefore created by Ottobock according to appendix II of the guidelines.

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