**MissAPPY****AIRTHERAPY**

Aerosol a pistone

MissAPPY

Istruzioni d'uso

Artsana S.p.A. - Via Saldarini Catelli 1 - 22070 Grandate (Co) Italy - Tel 031-382111
www.picindolor.com

Aerosol a pistone Miss APPY

Istruzioni d'uso

Numero verde consumatori Italia 800 188 898

pag. 2-9

Miss APPY piston aerosol

Operating instructions**pag. 10-16**

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Mode d'emploi**pag. 17-24**

Aerosol de pistón Miss APPY

Instrucciones de uso**pag. 25-32**

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Dear client

Thank you for choosing the PiC Indolor Miss APPY, the aerosol therapy device designed and produced using the most advanced technology. For an even easier use, please follow the instructions on preparing and using the Miss APPY aerosol unit. However, before proceeding, we strongly recommend you to read the short list of warnings in this manual to make sure you have fully understood how the device works. Miss APPY is a class IIa medical device, which turns medicines in a liquid form into an aerosol spray, to be directly administered to the respiratory tract. With the Sidestream nebulizer cup, medicines penetrate the respiratory tracts more easily, thanks to the optimal size of the aerosol particles. Sidestream significantly increases particles flow and reduces therapy time, which has a critical impact on the treatment's efficacy and its acceptance, especially with children.

The device, which is particularly attractive for younger users, has been specifically designed to make treatment more pleasant for children, whilst ensuring all the quality and safety standards of PiC Indolor products.

Note: despite its appearance, this device is not a toy. Keep out of reach of children and do not leave children alone during treatment.

The Miss APPY pack contains all accessories for the most suitable type of treatment, with a mouthpiece, a nose fork and a mask (both for adult and paediatric use). Accessories are in a practical apple-shape kit.

Miss APPY is reliable and easy to use, plus no special maintenance is required. The special air inlet filter retains dust and can be replaced at regular intervals. Thank you for choosing the new Miss APPY aerosol unit.

We are sure you will be extremely pleased with and will benefit from your purchase.

DEVICE COMPONENTS

- a-On/off button
- b-Compressed air outlet
- c-Filter
- d-Nebulizer cup base
- e-Plug
- f- Handle
- g-Air vents

ACCESSORIES

- 7- Adult mask
- 8- Paediatric mask
- 9- Nose fork
- 10- Mouthpiece
- 11- Connection tube
- 12- Sidestream nebulizer cup

GENERAL WARNINGS

- Before use, carefully read the information in this manual and keep for future reference.
- Device for aerosol therapy. Follow the instructions of a medical practitioner as to the type of medicine to use, the dose, frequency and duration of inhalations. All other uses apart from the intended use of the device are improper and therefore dangerous. Artsana cannot be held liable for any damage caused by incorrect or unreasonable use or if the device is used in electrical systems that do not conform to safety regulations in effect.
- After removing the device from its packaging, check that it is intact, without any visible damage that could have been incurred during transport. If in doubt do not use the device, but contact a qualified technician.
- Packaging materials (bags, box, etc.) must not be left within reach of children as they are potential hazards

- Before connecting the device, check that your mains voltage corresponds to that one shown on the rating label. The rating label is situated on the bottom of the device.
- If the plug of the device is not suitable for the mains socket, use an adaptor approved for use in the country where you are using the device, or have the plug replaced by a qualified technician authorised by Artsana SpA.
- To prevent overheating and compressor damage, turn the device off after 20 minutes of use and leave off for 40 minutes.
- Make sure filters are dry to ensure the device works properly.
- The device must not be used in the presence of nitrous oxide, oxygen or anaesthetic mixtures flammable with air.
- Turn off the device and disconnect from the mains after use and before adding medicine. Do not overfill the chamber.
- Keep the device and power lead away from hot surfaces.
- Do not use the device while having a bath or shower, or when in a damp place, or close to bath-tubs, sinks, washbasins, etc.
- Never touch the device with wet or damp hands.
- Do not drop or put the device into water or other liquids. If the device is dropped or falls into water or other liquids, remove the plug immediately, do not use the device and contact a qualified technician.
- Do not obstruct air vents during use.
- Do not use the device if you feel drowsy.
- Children and physically and mentally challenged users must be closely supervised when using the device.
- Do not use connectors or accessories not recommended by the manufacturer.
- Make sure the device is placed on level, stable surfaces when used (and stored) in order to prevent spillage.
- Make sure the surface the device rests on when used is free of any items that may

obstruct air flow.

- Never leave the device plugged into the power supply when it is not operating or when it is unattended.
- Disconnect the device from the mains before cleaning and/or maintenance operations.
- If you decide to no longer use the device, disconnect from the mains and cut off the plug so the device cannot work. We also recommend eliminating any parts that are potential hazards, particularly for children.
- To dispose of accessories subjected to wear, please refer to applicable legislation. To dispose of the device, please refer to the EU directive 2002/96/EC.
- This device must not be used if any part has been damaged or weakened, after it has fallen or been dropped. If in doubt contact a qualified technician.
- In the case of failure or malfunctioning of the device, switch it off by removing the plug from the power supply. Do not try to open or tamper with the device.

SIDESTREAM® NEBULIZER CUP

Effective aerosol treatment depends on how well the medicine circulates in the respiratory tree. The nebulizer cup - and its capacity to create particles of an ideal size and distribute them in the respiratory tract - is essential, in order to maximise breathable flow and reduce treatment times.

The Sidestream nebulizer cup (12) has been developed basing on numerous clinical studies on aerosol therapy.

Sidestream ensures a faster and more regular therapy, with a greater quantity of nebulised medicine and particles at the ideal size for treatment of the entire respiratory tract.

Sidestream uses the exclusive Diamond Jet technology, guaranteeing a more uni-

form nebulisation: in fact, the medicine crosses through 4 channels independent of the air flow channel. This avoids premature nebulisation within the tubes [13].

The Sidestream nebulizer cup features the unique Venturi Double Effect, where air enters the nebulizer cup from the bottom (compressed air) - like conventional nebulizer cups, as well as from the top. The additional air flow from outside the nebulizer cup combined with the air flow from the compressor conveys even more particles to the lungs [14].


At the nebulizer cup summit, the medicine is nebulised by the air produced by the compressor and air drawn from outside. A central deflector creates two currents of nebulised air, which hit the nebulizer cup walls. This impact causes the largest particles to drop to the nebulizer cup floor and reinitiate the cycle, while those with a suitable dimension leave the nebulizer cup and are breathed in by the patient.

This particular technology significantly increases particles flow and reduces therapy time, which has a critical impact on the treatment's efficiency and its acceptance, especially with children.

NEBULIZER CUP EFFECTIVENESS

The special configuration of Sidestream enables medicine to penetrate the respiratory tracts more easily, thanks to the optimal size of nebulised aerosol particles. With the Sidestream nebulizer cup, 80% of particles are smaller than 5 micron (from 1 to 5 micron) and can reach and deposit in the lower airways, while particle size in traditional nebulizer cups varies, so only 50% of particles reaches the lower airways.


T. Loffert, D.Ikle, H.S. Nelson / Chest, 1994

NEBULISATION RATE

The design of the Sidestream nebulizer cup also reduces treatment length, with medicine administered quickly and effectively. As the table shows, Sidestream nebulises the same amount of medicine as a traditional chamber in half the time.

PK Hurley, SW Smye/Technology and Health Care, 1994

NEBULISATION QUALITY

Sidestream produces same-size, optimal particles to maximise breathable flow. Breathable flow with the Sidestream nebulizer cup is increased by 140% compared to a conventional nebulizer cup.

PREPARING AND USING THE DEVICE

- Remove the device and accessories from the pack and take the Sidestream nebulizer cup (12)
- Unscrew the top of the Sidestream nebulizer cup in an anticlockwise direction (4)
- Add the medicine and/or saline solution to the nebulizer cup tank. Use doses recommended by a medical practitioner or pharmacist.
- Close the nebulizer cup, tightening the top in a clockwise direction
- Connect the tube (11) to the nebulizer cup (5) and compressed air outlet nozzle (b)
- Choose one of the following accessories for treatment: nose fork (9), mouthpiece (10), adult mask (7) or paediatric mask (8). If in doubt as to which accessory to use, consult a medical practitioner.
- Connect the devices to the mains
- Turn on the device, put the ON switch (a) to position I and start treatment.
- Start inhalation. During treatment, the nebulizer cup can be placed on its base (d)
- Carry out the treatment preferably when seated and in a comfortable position
- When the aerosol flow becomes intermittent, stop the treatment for a few seconds and let the droplets of medicine deposit on the walls of the nebulizer cup. Start treatment again and end the session when all nebulised medicine has been expelled out of the nebulizer cup
- After treatment, turn off the device, remove the plug (e) from the socket, remove the used accessory and clean the device and accessories as indicated in the section Cleaning and Maintenance

CLEANING AND MAINTENANCE

Before cleaning or maintenance, always turn off the device and remove the plug from the socket.

Thoroughly clean the device at regular intervals to ensure it stays in good working order.

Cleaning and maintenance of the compressor

Caution! Do not wet or put the device into water or other liquids. Use only a soft, dry cloth to clean the device.

Cleaning and maintenance of the nebulizer cup and accessories

Remove and separate all nebulizer cup components and accessories, immerse in boiling water for at least 5 minutes. Rinse the accessories under running water. After cooling, disinfect parts that come into contact with the user with denatured alcohol.

Carefully dry all components and the device before reassembling.

Store in a cool, dry place, away from light and heat.

Never use benzene, thinners or other inflammable chemical substances for cleaning.

For greater hygienic safety, we recommend you to avoid using the same accessories for more than one person and purchase a kit for each user (Aerosol kit cod. 00 038405 000 000).

CHECKING AND REPLACING FILTERS

We recommend checking filters at regular intervals. The filter protects the com-

pressor. Proper maintenance of the filter prolongs the life of the aerosol unit. Filter replacement times depend on the external conditions of use of the device. We recommend replacing filters after around 20 hours of use or if they are visibly dirty.

The filter is on the front of the aerosol unit (c)

To replace the filter :

- Open the filter cover
- Remove the filter from its housing
- Carefully fit the new filter
- Re-fit the filter cover

TECHNICAL DATA

Nominal voltage 230 V

Frequency 50Hz

Power: 140VA

Fuse: T 1,6A 250~

Maximum flow rate : 14 l/min

Operating flow : 7 l/min

Maximum pressure : 3 atm

Working pressure : 0.6 atm

Protection against liquid seepage : IPX0

Weight: 2Kg

Conditions for using the device: 20 ON

– 40 OFF

MMAD: 2,63 µm

% of particles below 5µ /

Breathable part: 77%

Nebulizer cup capacity : 5 ml

Noise level at 50 cm* : 65 dB

*Data on noise levels were measured on a new device. Values may vary with use.

The information on performance provided by the manufacturer in accordance with standard EN 13544-1 may not be applicable to medicines sold as suspensions or with viscosity.

EN 60601-1 Medical electrical equipment
– Part 1 : general requirements for safety

EN 60601-1-2 Medical electrical equip-

ment – Part 1 : general requirements for safety – Collateral standard : electromagnetic compatibility – Requirements and tests

EN 13544-1 Respiratory therapy equipment – Part one: Nebulizing systems and their components

LEGEND OF SYMBOLS:

Caution! Read instructions

Type B device

Double insulation

Complies with directive
MDD 93/42/EEC

Alternate current

Conditions for storage :

Relative humidity: 15% / 85%

Temperatures: -10°C / +55°C

“The medical device requires special precautions regarding electromagnetic compatibility and shall be installed and put into service according to the information provided in the user manual”

Guidance and Manufacturer's Declaration – Electromagnetic Emissions

The Miss APPY aerosol therapy unit is intended for use in the electromagnetic environment specified below. The client or user of the Miss APPY aerosol therapy unit shall make sure that the device is used in this environment

Emissions test	Compliance	Electromagnetic Environment - Guidance
Radiated/Conducted Emissions CISPR11	Group 1	The Miss APPY aerosol therapy unit only uses RF energy for its internal function. Therefore its RF emissions are very low and do not cause any interference in nearby electronic equipment.
Radiated/Conducted Emissions CISPR11	Class [B]	The Miss APPY aerosol therapy unit is suitable for use in all environments, including domestic ones and those directly connected to the public power supply network that supplies buildings used for domestic purposes.
Harmonic Emissions IEC/EN 61000-3-2	Class [A]	
Voltage Fluctuations/Flicker Emissions IEC/EN 61000-3-3	[Satisfies]	

Guidance and Manufacturer's Declaration – Electromagnetic Immunity

The Miss APPY aerosol therapy unit is intended for use in the electromagnetic environment specified below. The client or user of the Miss APPY aerosol therapy unit shall make sure that the device is used in this environment

Emissions test	Compliance	Electromagnetic Environment - Guidance
Electrostatic discharge [ESD] IEC/EN 61000-4-2	±6kV on contact ±8kV in the air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should not exceed 30%.
Electrical fast transient/burst IEC/EN 61000-4-4	2 kV power supply	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC/EN 61000-4-5	1 kV differential mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations IEC/EN 61000-4-11	<5% Ut for 0.5 cycles 40% Ut for 5 cycles 70% Ut for 25 cycles <5% Ut for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Miss APPY aerosol therapy unit requires continued operation, it is recommended that the unit is powered from an uninterruptible power supply.
Magnetic field IEC/EN 61000-4-8	3A/m	The magnetic field should be that of a typical commercial or hospital environment.
Immunity to conducted disturbances IEC/EN 61000-4-6	3V/m 80MHz to 2.5GHz (for non-life-supporting equipment)	
Irradiated immunity IEC/EN 61000-4-3	3V/m 80MHz to 2.5GHz (for non-life-supporting equipment)	

NoTe: Ut is the mains voltage

**This product complies with
EU Directive 2002/96/EC.**

The crossed bin symbol on the appliance indicates that the product, at the end of its life, must be disposed of separately from domestic waste, either by taking it to a separate waste disposal site for electric and electronic appliances or by returning it to your dealer when you buy another similar appliance. The user is responsible for taking the appliance to a special waste disposal site at the end of its life. If the disused appliance is collected correctly as separate waste, it can be recycled, treated and disposed of ecologically; this avoids a negative impact on both the environment and health, and contributes towards the recycling of the product's materials. For further information regarding the waste disposal services available, contact your local waste disposal agency or the shop where you bought the appliance.

WARRANTY

Miss APPY is covered by a warranty in the case of manufacturing defects in materials or workmanship. The warranty is for 2 years. Faulty components will only be replaced free of charge if the device has not been tampered with and the enclosed instructions for use have been scrupulously followed.

Any tampering/repairs carried out on the device by unauthorised persons will immediately invalidate the warranty.

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Sidestream è un marchio registrato di Respirationics (UK) Ltd. - Sidestream is a registered trademark of Respirationics (UK) Ltd.