IgG-IEF

Instructions For Use

REF 102100 (Double Antibody)

REF 102200 (Single Antibody)

IEF-IgG kit

Fiche technique

REF 102100 (Double anticorps)

REF 102200 (Simple anticorps)

IgG IEF Kit

Anleitung

REF 102100 (Double Antikörper)

REF 102200 (Single Antikörper)

IgG IEF KIT

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Kit IgG-IEF

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REF 102100 (anticuerpo doble)

REF 102200 (anticuerpo simple)

Contents

English	I
Français	8
Deutsch	15
Italiano	23
Español	30

INTENDED PURPOSE

The IgG IEF Kit is intended for the identification of IgG-specific oligoclonal banding in serum and CSF by agarose gel isoelectric focusing and immunoblotting. The Helena BioSciences IgG-IEF Kit is based upon the method originally described by Keir, G. *et al*.

The diagnosis of multiple sclerosis is ultimately a clinical decision, although the analysis of cerebrospinal fluid gives evidence of an inflammatory response within the central nervous system, showing oligoclonal bands on the immunoglobulin electrophoresis pattern. In 1994, the Committee of the European Concerted Action for Multiple Sclerosis recommended isoelectric focusing as the most sensitive method for the detection of oligoclonal bands in serum and CSF¹.

The IgG IEF Kit separates serum or CSF proteins according to isoelectric point in an agarose IEF gel, pH 3-10. The proteins are then transferred to nitrocellulose, and the nitrocellulose immunofixed to visualise IgG specific bands. The patterns are interpreted qualitatively by comparing the presence or absence of oligoclonal bands in serum and/or CSF.

WARNINGS AND PRECAUTIONS

All reagents are for in-vitro diagnostic use only. Do not ingest or pipette by mouth any kit component. Refer to product safety data sheet for risk and safety phrases and disposal information.

COMPOSITION

1. **IgG IEF Gel ($\times 10$)**
Contains agarose and separators, pH 3-10 with thiomersal as preservative. The gel is ready for use as packaged.
2. **IgG IEF Anode Solution (1 $\times$ 25ml)**
Contains 0.3M acetic acid. The solution is ready for use as packaged.
3. **IgG IEF Cathode Solution (1 $\times$ 25ml)**
Contains 1M NaOH. The solution is ready for use as packaged.
4. **IgG IEF Blocking Agent (1 $\times$ 10g)**
Contains bovine milk protein solids. IgG-IEF Blocking Agent is used to prepare 2 working solutions: Blocking Solution A: Dissolve 1g powder in 50ml of saline solution. Blocking Solution B: Dilute 5ml of blocking solution A to 50ml with saline solution.
5. **IgG IEF Acetate Buffer Concentrate (1 $\times$ 25ml)**
Contains 0.2M acetate buffer pH 5.1. Dilute 2.5ml of Acetate Buffer concentrate to 25ml with purified water prior to use and use immediately.
6. **IgG IEF Chromagen (1 $\times$ 0.05g)**
Contains 3-amino-9-ethyl-carbazole powder. Dissolve the contents of the vial in 50ml of pure methanol. Store 5ml aliquots at -20°C in the dark. Use glassware dedicated for this purpose to prevent protein contamination of the chromagen.
7. **Antibodies for 102100 Kit:**
IgG IEF Primary Antibody (1 $\times$ 240 μ l)
Contains sheep anti-human IgG antiserum and is ready for use as packaged.
IgG IEF Secondary Antibody (1 $\times$ 200 μ l)
Contains anti-sheep IgG peroxidase conjugate and is ready for use as packaged.
8. **Antibodies for 102200 Kit:**
IgG IEF Single Antibody (1 $\times$ 200 μ l)
Contains sheep anti-human IgG peroxidase conjugate and is ready for use as packaged.

9. Other Kit Components

Each kit contains Instructions For Use and sufficient Blotting Membranes, Blotters C and D, Electrode Wicks, Wick Blotters and Soak-away Blotters to complete 10 gels.

STORAGE AND SHELF-LIFE

1. IgG IEF Gel

Gels should be stored at 15...30°C and are stable until the expiry date indicated on the package. DO NOT REFRIGERATE OR FREEZE. Deterioration of the gel may be indicated by 1) crystalline appearance indicating the gel has been frozen, 2) cracking and peeling indicating drying of the gel or 3) visible contamination of the agarose from bacterial or fungal sources.

2. IgG IEF Anode Solution

The anode solution should be stored at 15...30°C and is stable until the expiry date indicated on the label. Cloudiness or particulate contamination may indicate deterioration.

3. IgG IEF Cathode Solution

The cathode solution should be stored at 15...30°C and is stable until the expiry date indicated on the label. Cloudiness or particulate contamination may indicate deterioration.

4. IgG IEF Blocking Agent

The blocking agent should be stored at 15...30°C and is stable until the expiry date indicated on the label. The IgG-IEF Blocking Agent should be a free-flowing white or off-white lumpy powder, any changes to this appearance, or unusual colouration may indicate deterioration.

5. IgG IEF Acetate Buffer Concentrate

The buffer concentrate should be stored at 15...30°C and is stable until the expiry date indicated on the label. Cloudiness or particulate contamination may indicate deterioration.

6. IgG IEF Chromagen

The chromagen should be stored at 2...6°C and is stable until the expiry date indicated on the label. Frozen aliquots are stable for 6 months. The Chromagen solution should be a red-brown colour, changes from this appearance or failure of the chromagen to perform in the test may indicate deterioration.

7. Antibodies for 102100 Kit:

IgG IEF Primary Antibody

The primary antibody solution should be stored at 2...6°C and is stable until the expiry date indicated on the label. Particulate contamination or cloudiness may indicate deterioration.

IgG IEF Secondary Antibody

The secondary antibody solution should be stored at 2...6°C and is stable until the expiry date indicated on the label. Particulate contamination or cloudiness may indicate deterioration.

8. Antibodies for 102200 Kit:

IgG IEF Single Antibody

The single antibody solution should be stored at 2...6°C and is stable until the expiry date indicated on the label. Particulate contamination or cloudiness may indicate deterioration.

ITEMS REQUIRED BUT NOT PROVIDED

Cat. No. 4065 Helena IEF Chamber
 Cat. No. 1525 EPS 600 Power Supply
 Cat. No. 1349/1111 REP IEF/SPIFE IEF Electrode Assembly
 Cat. No. 5014 Development Weight
 Cat. No. 3039 Chamber Cooling Device
 30 volume hydrogen peroxide
 Methanol
 10% methanol solution
 Drying oven with forced air capable of 60...70°C
 Saline Solution (0.85% NaCl)
 Purified water

SAMPLE COLLECTION AND PREPARATION

Freshly collected serum or CSF is the specimen of choice. Samples can be stored refrigerated at 2...6°C for up to 3 days or 2 weeks at -20°C. CSF samples should be applied neat. Serum samples should be diluted in physiological saline to an IgG concentration equivalent to the matched CSF sample.

STEP-BY-STEP PROCEDURE

1. Remove the gel from the packaging and place on a paper towel. Remove the overlay and blot the gel surface with a blotter C, discard the blotter.
2. Soak 1 electrode wick in the anode solution and blot almost dry between two wick blotters. Position this wick along the anode (+) edge of the gel. Repeat this for the cathode solution and position along the cathode (-) edge of the gel. **NOTE:** Take care not to mix the anode and cathode solutions.
3. Apply 5µl of sample to each sample well.

NOTE: Follow procedure according to instrument.

4a. Running in a REP unit:

- i. Use a REP IEF electrode assembly. Dispense 500µL of di-water onto the chamber. Place the gel into the chamber by positioning the anode wick under the central anode electrode of the REP IEF Electrode Assembly and the cathode wick under the outside cathode electrode. Position one soak-away blotter on the outside of the cathode electrode and one soak-away blotter on the outside of the anode electrode to act as excess fluid blotters. Ensure the soak-aways are in contact with the gel.
- ii. Program the following parameters into the REP (leaving all other sections of the program screen blank):

Electrophoresis Voltage	600-800 volts
Electrophoresis Time	45:00 mm:ss
Electrophoresis Temperature	10...20°C
Standby Temperature	16°C

- iii. Slide the chamber lid to close and press F1 to commence the run.

4b. Running in a SPIFE / SAS-3 unit:

- Use a SPIFE / SAS-3 IEF electrode assembly. Position gel as described in 4a (i).
- Program the following parameters into the SPIFE / SAS-3:

Electrophoresis	45:00 mm:ss	10°C	800 volts
End of test			

- Close the chamber lid and press START to commence electrophoresis.

4c. Running in the SAS-I Plus unit:

- Dispense 400 μ L of di-water onto the heat sink. Place the gel onto the heat sink, agarose side up, aligning the positive and negative sides with the corresponding electrode posts, taking care to avoid air bubbles under the gel.
- Attach the electrodes onto the top side of the electrode posts so that they are in contact with the buffer wicks.
- Perform the IgG-IEF electrophoresis:

Number of applications	0
Electrophoresis Voltage	600 volts
Electrophoresis Temperature	15°C
Electrophoresis Time	60:00 mm:ss
Incubation	OFF

4d. Running in conventional IEF systems:

- The following parameters should be used if using alternative IEF systems:

Power	10-20 watts
Volts Limit	600-1200 volts
Current Limit	50mA
Run Time	30:00 mm:ss
Run Temperature	10...20°C

4e. Running in the Helena IEF Chamber:

- The following parameters should be used with the Helena IEF chamber:

Electrophoresis Voltage	400 volts
Electrophoresis Time	60 mm:ss

NOTE: The voltage, time and temperature chosen for the focusing should be determined by each individual laboratory according to pattern length required. The temperature should also be chosen to minimise condensation on the gel surface which can lead to smearing of the patterns. Avoid overfocusing of the gel, as cathode drift of the gradient can result in the loss of proteins into the cathode wick, as well as the introduction of breaks in the gradient. When using Helena IEF chamber ensure the cooling device has been cooled to 5°C before use.

5. Following focusing, place the gel onto a blotter D. Remove the electrode wicks.
6. Pre-wet a sheet of transfer membrane in purified water. Remove excess water by gently blotting between two paper towels or tissues. Carefully place this sheet of transfer membrane, red dot down, onto the gel surface. Allow the membrane to take up moisture as it is placed on the gel surface to prevent the production of air pockets. Leave for 10 seconds and discard.
7. Pre-wet a second sheet of transfer membrane in purified water or 10% methanol solution. Allow excess fluid to drain from the membrane by holding the membrane vertical for 10-20 seconds. Carefully lay this transfer membrane, red dot down, onto the gel surface, and rub the surface of the membrane to ensure good gel contact and the removal of air pockets between the gel and the membrane (avoid excessive pressure which could damage the gel). Place a blotter C, three blotter D's and a development weight on top. Press for 30 minutes.
8. Discard the blotters and place the transfer membrane, red dot up, into a staining dish containing 20 ml of blocking solution A so that the whole surface of the membrane is covered. Gently agitate for 30 minutes.
9. Wash the membrane in purified water several times.

10a.DOUBLE ANTIBODY PROCEDURE

Place the membrane in 15ml of blocking solution B containing 15 μ l of primary antibody for 30 minutes with gentle agitation. Wash the membrane in purified water several times and, finally, in saline solution for 5 minutes. Place the membrane in 15ml of blocking solution B containing 15 μ l of secondary antibody for 30 minutes with gentle agitation.

10b.SINGLE ANTIBODY PROCEDURE

Place the membrane in 15ml of blocking solution B containing 15 μ l of single step antibody for 30 minutes with gentle agitation.

11. Wash the membrane in purified water several times and, finally, in saline solution for 5 minutes. Rinse in purified water. Place the membrane in a solution of 25ml acetate buffer (diluted as specified), 5ml chromagen solution and 25 μ l of hydrogen peroxide. Allow colour to develop for 20 minutes with gentle agitation. It is normal for precipitation of chromagen to occur during this step.
12. Wash the membrane in purified water several times and allow to air dry.

INTERPRETATION OF RESULTS

For a full review of the interpretation of IgG-IEF patterns, see Andersson, M. *et al*². The IgG bands become visible as a red/red brown colour on the membrane. Comparison of paired serum/CSF samples should always be the standard procedure. Various disease states show monoclonal or oligoclonal IgG bands. Multiple Sclerosis (MS) patients typically show oligoclonal IgG in the CSF without bands in the serum. The interpretation of the patterns obtained should be performed with the aid of the clinical history of the patient.

QUALITY CONTROL

The use of CSF from a donor positive for oligoclonal bands is recommended for use on each gel as a control sample to confirm all stages of the procedure in the absence of a commercially available control material.

TROUBLESHOOTING

The following table details possible reasons and solutions to some of the most commonly encountered problems associated with isoelectric focusing followed by immunoblotting.

Problem Observed	Possible Reasons and Solutions
No staining of sample lanes.	a) No samples applied. b) Primary and/or secondary antibody omitted. c) No hydrogen peroxide added.
Background heavily and uniformly stained.	a) Blocking step omitted. b) Inadequate washing of the transfer membrane during the immunofixation and colouration steps.
Background heavily but unevenly stained.	a) Omitted or inadequate pre-blot step
Sample lanes stained but faint.	a) Insufficient sample applied. b) Insufficient staining time. c) Old or contaminated hydrogen peroxide. d) Decay of antibodies due to incorrect storage. e) Oxidised chromagen (will appear black). f) pH of acetate buffer incorrect (should be 5.0 - 5.2).
Sample lanes intensely stained.	a) Excessive sample applied losing resolution. b) Overdeveloped (rare).
Sample lanes distorted/smeared.	a) No soak-aways used. b) Condensation on gel surface due to excessive cooling of the gel.
All staining in the sample well.	a) No power applied during electrophoresis.
White spots in the stained sample.	a) Air bubbles between gel and membrane at press stage. b) Ensure a pre-wet membrane is used and that bubbles are rubbed out from between the membrane and gel.

PERFORMANCE CHARACTERISTICS

Sensitivity

A discrete band was evident on the membrane at a concentration of IgG equivalent to 0.28mg/L in the CSF.

Reproducibility

a. Within-gel

A single CSF sample showed the same pattern of bands when tested in all 10 wells of a single gel.

b. Between-gel

The same sample was tested on 3 different gels and showed the same pattern of bands on each of the gels.

Specificity

The antibodies used in the kit were shown to be specific for IgG, and did not cross-react with human IgA or IgM.

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