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
Figure 6 - X-Ray image/Image radiographique

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Integra NPH™ Low Flow Valve

DESCRIPTION

The Integra NPH™ Low Flow Valve is an implantable hydrocephalus valve system for controlled cerebrospinal fluid (CSF) drainage from the ventricles to the peritoneal cavity or other appropriate drainage site such as the heart's right atrium. Unlike conventional valves, the Integra NPH™ Low Flow Valve is a variable resistance valve that maintains a drainage at a constant rate, between 8 and 17 ml/hr, within the physiological range of intracranial pressure (ICP). The mechanism incorporates a safety pressure relief mode to prevent accidental intracranial hypertension.

Caution: The Integra NPH™ Low Flow Valve may appear similar to the Integra OSV II®. However it features a black radiopaque flow direction arrow to differentiate both products, visually during implantation, and by X-ray, after implantation. X-ray image of Integra NPH™ Low Flow Valve is included as Figure 6.

WARNING:

The Integra NPH™ Low Flow Valve should only be used in patients who need a low drainage rate (8-17 ml/hr). For patients who need a higher drainage rate (18-30 ml/hr), the Integra OSV II® should be used.

Principles of Operation

A variable flow restrictor (**Figure 1**) consists of a silicone elastomer diaphragm (4), synthetic ruby seat (5), and notched pin (6). The diaphragm reacts to differential pressure (DP) variations. The seat is inserted into the center of the diaphragm. As pressure varies, the clearance between the seat and the pin increases or decreases depending on seat movement along the pin. A typical flow/pressure curve is shown in **Figure 2**.

Note - Differential pressure/flow characteristics for each valve are verified to be within specification at time of manufacture.

The three stages of operation of the Integra NPH™ Low Flow Valve are defined below:

Stage I - Low Differential Pressure

This stage begins when the flow rate through the valve reaches 5ml/hr (DP will be between 30 and 80 mm H₂O). The valve remains in Stage I with CSF flow rates up to 8 ml/hr (DP will be between 40 and 120 mm H₂O).

Stage II - Flow Regulation

When the DP increases, the valve operates as a variable resistance flow regulator. At DP ranges between 120 and 300 mm H₂O, the valve restricts flow between 8 and 17 ml/hr.

Stage III - Pressure Relief Mode

Should the intraventricular pressure (IVP) elevate abruptly, the shunt operates in a rapid flow rate mode to facilitate IVP normalization. The valve then reverts to Stage II or I, depending upon conditions.

Configurations

The system is available in different configurations; a system generally consists of a ventricular (proximal) catheter, a valve unit, and a drainage (distal) catheter. Some configurations contain an antechamber or a burr hole cap. Configurations consisting of one-piece (integral), and two-piece (separate ventricular catheter) systems are available. Included with the valve systems that have a ventricular catheter are one introducing rod or one ventricular catheter introducer, one Luer connector, one straight connector and one Right-Angle Guide (RAG). Configurations with a burr hole cap contain burr hole reservoirs. Configurations without a ventricular catheter do not contain a ventricular catheter introducer (or rod), or Right-Angle Guide. A malleable tunneler is supplied with some configurations for subcutaneous tunneling of the distal tubing. Marked ventricular catheters feature radiopaque length markings at 2 cm intervals from the tip. Connectors are designed for use with 1.1 to 1.4 mm inner diameter silicone elastomer tubings.

Implanted Materials

The implanted product may contain one or more of the following: Silicone elastomer with or without barium sulfate, polysulfone, synthetic ruby, silicone adhesive, tantalum, polypropylene and epoxy.

Compatibility with other Medical Procedures

The Integra NPH™ Low Flow Valve Systems contain no ferromagnetic parts, and are safe and compatible with diagnostic procedures such as Magnetic Resonance Imaging (MRI) and CT scans.

INDICATIONS

The Integra NPH™ Low Flow Valve is an implantable system used in the treatment of patients with hydrocephalus, to shunt CSF from the ventricles to the peritoneal cavity or other appropriate drainage site such as the heart's right atrium.

CONTRAINDICATIONS

This valve system should not be implanted when an infection along the shunt pathway (e.g. meningitis, ventriculitis, peritonitis, septicemia, bacteremia) is suspected. It is advisable to postpone valve implantation if an infection is present anywhere in the body. Atrial shunting is not advised for patients with congenital heart disease or other serious cardiopulmonary abnormalities.

Integra NPH™ Low Flow Valve should not be implanted in patients with untreated choroid plexus tumors. Such tumors produce CSF at rates in excess of the specification of the flow regulation Stage II, and the valve would underdrain under these conditions.

Integra NPH™ Low Flow Valve should not be used for drainage of extraventricular fluid collections such as hygromas or cysts; such conditions are typically treated with very low differential pressure valves.

SIDE EFFECTS

- In addition to the risks associated with shunt system or component implantation, major complications include mechanical failure, such as tubing fracture, material deterioration with time, or shunt pathway obstruction, infection, reaction to materials, and CSF leakage along the shunt pathway. System disconnection can lead to catheter migration into the atrium, peritoneum or lateral ventricles.
- Ventricular shunting has been associated with epilepsy and shunt dependency. Intracranial hypertension syndrome is a specific result of obstruction. Overdrainage may result in subdural hematomas, craniostenosis, intracranial hypotension syndrome or sunken fontanelles (in infants).
- Shunt system failure should be suspected if any of these is noted: Continuing symptoms of hydrocephalus; CSF leakage; redness, tenderness or erosion of the skin along the shunt pathway. These complications require prompt replacement of the shunt system or affected component(s).
- Low-grade bacterial colonization may cause recurrent fevers, anemias, splenomegaly, shunt nephritis or pulmonary hypertension.
- With atrial catheterization, the possibility of clot formation around the atrial catheter tip could lead to pulmonary tree embolization.
- Perforation of an abdominal viscus is a rare complication of ventriculoperitoneal shunts.

PRECAUTIONS

- Closely observe patients with an implanted valve system for shunt failure symptoms (see "Side Effects").
- During implantation, have duplicates of each component available.
- Handle valve tubing carefully. Cover forceps and clamp jaws with silicone elastomer tubing.
- Do not lubricate valves. They will be lubricated adequately by CSF.
- Silicone elastomer requires special handling to prevent contamination. Avoid contact with glove talc or skin oils.
- The drainage catheter may be trimmed to the desired length during implantation. Other modifications are not recommended.
- In children, insert sufficient drainage catheter length into the abdominal cavity to allow for growth.
- To minimize the possibility of underdrainage in patients confined to a supine position (i.e., patients under 3 months of age and bedridden patients), these patients should be elevated at a 45° angle for several hours a day (using a pillow, etc.) until they can sit or stand independently.

Caution - Do not overpressurize the valve system. Excessive flushing pressures may lead to valve damage.

Caution - Do not use in conjunction with other devices, such as antisiphon devices or distal slit valves. Addition of such components would change the flow/pressure characteristics of the valve.

Storage and Transport: The valve systems in Integra packaging may be stored at room temperature. They will withstand the conditions of normal transport. This product is recommended for single use only. Do not implant explanted systems or components in another patient. *This product is sterilized with ethylene oxide. Do not use if the package is open or damaged. Use the device prior to the "Use by" date on the package label.*

Caution - Do not resterilize. Integra will not be liable for any direct, indirect, incidental, consequential or punitive damages resulting from or related to resterilization.

Integra NPH™ Low Flow Valve PRESSURE / FLOW TEST

Each valve is tested at the end of the manufacturing process. Integra does not recommend to test the valve in the operating room, before implantation, for obvious reasons of possible contamination, and because the conditions required for this testing are difficult to meet in a clinical environment.

PREIMPLANT TESTING FOR PATENCY

Should verification of valve patency be desired prior to implantation, place the inlet of the valve system in a container of sterile apyrogenic saline solution. Connect a syringe, with the Luer connector attached to the drainage catheter, and gently draw the solution through the system to displace air (**Figure 3**). Remove syringe and Luer connector. Holding the valve's inlet end in the sterile solution, allow the drainage catheter to hang freely (**Figure 3**). The valve is patent if solution drips from the drainage catheter.

Note - The observed outflow is lower than that observed with conventional DP shunts because of the flow restriction mode of the valve.

Note - This procedure only tests the valve patency. It does not allow verification of the 3 different stages of operation of the valve.

RECOMMENDED IMPLANTATION PROCEDURE

Before implantation, the valve serial number should be recorded in the patient's file.

Note - After implantation, complete the patient identification card for the patient or family.

Clinical Judgment

Incision size, size and type of nonabsorbable suture material, etc. should be selected according to the surgeon's experience and preference. Procedure may be varied in accordance with the surgeon's clinical judgment.

Ventriculoperitoneal Approach

The described procedure is based on the use of a separate ventricular catheter and a specific type of peritoneal introducer. If a valve with an integral ventricular catheter is used, follow the procedure "Introduction of Integral Ventricular Catheter".

1. Place the patient in the supine position. Administer anesthesia. To facilitate tunneling, the head is turned to the contralateral side, placing the abdomen, thorax, neck and skull of the occipital region on a straight plane. Padding may be placed under the neck to facilitate alignment.

Note - The skull incision location depends on the ventricular approach (frontal, temporal or occipital).

2. Expose the skull by making a 25 mm linear incision. Using a perforator, make a burr hole 3.0 mm or larger. Using an electrocoagulator, make an opening in the dura just large enough for ventricular catheter passage. For systems with burr hole cap: Using a perforator, make a burr hole larger than 6.5 mm. The technique may be varied according to the surgeon's experience.

3. Using blunt dissection, create a small pocket under the skin to avoid unnecessary pull strength on the drainage catheter when positioning the shunt. Pass the Integra disposable, malleable tunneler from the subcutaneous pocket to the abdominal area, directing the tunneler tip by trapping it between the thumb and the index finger. The tunneler may be deformed to follow body contours.

Note - An alternate method may be necessary when the tunneler is not long enough or when the patient's position precludes a single pass without an interim incision.

Make the abdominal incision. Remove the handle of the tunneler. Attach the drainage catheter to the connector of the internal tunneler wire, and pass the catheter from the head to the abdominal incisions. Remove the metal outer shaft. Pull on the distal end of the drainage catheter. Position the valve under the skin.

4. Place the introducing rod into the ventricular catheter. If the Right Angle Guide (RAG) is used, introduce the ventricular catheter into the center of the RAG and insert the catheter into the ventricle. Holding the catheter with covered forceps, remove the introducing rod. Verify adequate CSF flow.

For systems with burr hole cap: Trim the ventricular catheter to the desired length. Insert the ventricular catheter into the burr hole reservoir and tie with non absorbable suture. Using the introducing rod, insert the catheter into the ventricle. Verify adequate flow.

Connect a syringe with the provided Luer connector to the ventricular catheter and gently aspirate 2 to 3 ml of CSF to eliminate possible debris from the CSF.

Note - If desired, sample CSF for analysis or inject antibiotics at this time.

5. Place the RAG into the burr hole. The channel should point in the direction of the valve. Bend the exposed catheter into the channel of the RAG until the catheter snaps in place.

Note - The use of the RAG in premature infants or neonates is at the discretion of the surgeon. If the RAG is not used, the surgeon should trim the rim of the burr hole where the catheter emerges to provide a smooth angular transition for the catheter.

Trim the ventricular catheter to length (approximately 2 cm from the burr hole).

6. Carefully insert the integral inlet straight connector of the valve into the ventricular catheter. Do not dislodge the catheter. Tie the catheter onto the connector and to the underlying fascia, using appropriate sutures.

Note - For configurations with burr hole cap, attach the cap to the reservoir.

7. **Connect a syringe with the provided Luer connector to the drainage catheter and rinse the system by gently aspirating 2 to 3 ml of CSF to eliminate possible debris and purge air.**

8. If desired, trim the drainage catheter.

9. Introduce the drainage catheter in the peritoneal cavity using a trocar or standard laparotomy technique consistent with the surgeon's experience.

10. **Caution - Prior to the introduction of the drainage catheter into the abdominal cavity, it is advisable to confirm the patency of the system by observing the formation of CSF droplets at the end of the drainage catheter.**

Note - The CSF droplets formation rate may be lower than that observed with conventional DP shunts because of the flow restriction mode of the Integra NPH™ Low Flow Valve.

11. Close incisions.

Note - X-ray the complete system just after implantation for future reference to determine whether system components have shifted.

Introduction of Integral Ventricular Catheter (Alternative for Steps 4-5)

If the RAG is used, slide the catheter through the center hole of the RAG. Position the RAG as far along the catheter as possible. Lock the stylet of the ventricular catheter introducer into its cannula as shown in **Figure 4**. Insert the tip of the introducer into the hole nearest the tip of the ventricular catheter. Align the catheter with the introducer cannula and snap it into the adjustable plastic collar.

Caution - To avoid improper placement, do not overstretch the ventricular catheter when positioning it on the introducer. Puncture the ventricle with this assembly. Carefully withdraw the stylet from the introducer cannula, causing the ventricular catheter tip to disengage from the stylet.

Note - Sample CSF at this time if desired by withdrawing fluid from the cannula.

Disengage the catheter from the collar on the cannula, ensuring that the catheter's position is maintained. Remove the introducer cannula from the burr hole. For instructions on the use of the RAG, refer to Steps 4 and 5 of the "Ventriculoperitoneal Approach" procedure.

Ventriculoatrial Approach

Caution - The procedure is based on the use of a separate ventricular catheter. If a valve with an integral ventricular catheter is used, follow the procedure "Introduction of Integral Ventricular Catheter".

1. Position the patient and open the skull as described in Steps 1 and 2 of the "Ventriculoperitoneal Approach" procedure above.
Using blunt dissection, create a small pocket under the skin to avoid unnecessary pull strength on the drainage catheter when positioning the shunt.
2. Make a neck incision over the point where the angle of the mandible crosses the anterior edge of the sternocleidomastoid muscle. If the internal jugular vein is used, facilitate entry by mobilizing the common facial vein at its junction with the internal jugular vein.
3. Tunnel the drainage catheter from the skull subcutaneous pocket to the neck incision. Position the valve under the skin.
4. Insert the ventricular catheter following the Steps 4, 5, and 6 of the "Ventriculoperitoneal Approach" procedure.
5. Trim the drainage catheter at the neck incision. Fill the distal part with sterile apyrogenic saline solution. Clamp the proximal end. Introduce the drainage catheter into the vein. Position the distal end of the drainage catheter in the heart's right atrium at the level of the 6th or 7th thoracic vertebra. Determine the exact location of the catheter by X-ray, radioscapy, ECG, pressure analysis or echography.
Note - To enhance X-ray visualization, the atrial catheter may be filled with contrast.
6. Carefully maintaining the catheter in place, trim its proximal part, and the valve outlet tubing to length. Connect the two together using a straight connector. Tie securely using appropriate sutures.
7. Close incisions.
Note - X-ray the complete system just after implantation for future reference to determine whether system components have shifted.

Valve Pumping after Implantation

Percutaneous depression of the antechamber with the finger will force CSF out of the valve in both the proximal and distal directions. Occlude either the distal or proximal side of the antechamber by finger pressure (according to **Figure 5**), then depress the antechamber. This will push the fluid in the direction opposite to the occlusion.

Warning - Shunt obstruction may occur in any component of a shunt system and should be diagnosed by clinical symptoms and X-rays. Valve pumping testing may not be adequate to diagnose occlusion of catheters.

Injection into the Valve System

The antechamber can be used to access CSF for injections. **Use of 25G or smaller beveled needles is recommended.** Care must be taken during insertion and removal of the needle to avoid tearing the silicone material. It is also advised to avoid puncturing several times at the same place. Taking into account the above recommendations, the antechamber is qualified to sustain up to 40 punctures without leaking.

Caution - Do not overpressurize the valve system. Excessive flushing pressures may lead to valve damage.

HOW SUPPLIED

Integra Valves Systems are supplied sterile and non-pyrogenic in a double-wrap packaging.

PRODUCT INFORMATION DISCLOSURE

INTEGRA NEUROSCIENCES HAS EXERCISED REASONABLE CARE IN THE SELECTION OF MATERIALS AND THE MANUFACTURE OF THESE PRODUCTS. INTEGRA NEUROSCIENCES EXCLUDES ALL WARRANTIES, WHETHER EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO, ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. INTEGRA NEUROSCIENCES SHALL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL LOSS, DAMAGE, OR EXPENSE, DIRECTLY OR INDIRECTLY ARISING FROM USE OF THIS PRODUCT. INTEGRA NEUROSCIENCES NEITHER ASSUMES NOR AUTHORIZES ANY PERSON TO ASSUME FOR IT ANY OTHER OR ADDITIONAL LIABILITY OR RESPONSIBILITY IN CONNECTION WITH THESE PRODUCTS.

Returned Goods Policy

Products must be returned in unopened packages with manufacturer's seals intact to be accepted for replacement or credit unless returned due to a complaint of product defect.

Determination of a product defect will be made by Integra NeuroSciences.

Products will not be accepted for replacement if they have been in the possession of the customer for more than 90 days.

Product Order Information

All products can be ordered through your Integra NeuroSpecialist or customer service representative or by contacting:

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OSV II is a registered trademark of Integra LifeSciences Corporation.
Integra NPH Low Flow Valve and the Integra wave logo are trademarks of Integra LifeSciences Corporation.