
homecare

**JOYCE vented**

**JOYCE vented 40 hPa**

**JOYCE non vented**

**JOYCE non vented 40 hPa**

**JOYCE GEL vented**

**JOYCE GEL non vented**

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---

### Legend

**(Description of individual parts: see enclosed leaflet)**

- 1 Headgear
- 2 Forehead cushion
- 3 Fine adjustment component of forehead support
- 4 Coarse adjustment component of forehead support
- 5 Forehead support
- 6 Headgear clips
- 7 Pressure measurement ports
- 8 Port cap
- 9 Rotating sleeve (non vented: blue coloured)
- 10 Elbow (non vented: blue coloured)
- 11 Retaining ring (non vented: blue coloured)
- 12 Mask frame
- 13 Mask cushion

### Intended use

The nasal mask is an accessory and acts as a seal between the patient and the therapy device. It is used for the treatment of sleep apnoea and for the non-invasive ventilation of patients with ventilatory insufficiency.

### Before use by the patient

Patients must be instructed in the use of the mask by qualified personnel before therapy begins.

## 2 Safety instructions

---

### Safety regulations

Following the requirements of Richtlinie 93/42 EWG please observe the following points:

- Please read these instructions for use very carefully. They are a part of the product and should be available at all times.

- Use the mask only in the therapy pressure ranges required for the individual variants:
  - 4-30 hPa: JOYCE vented  
JOYCE non vented
  - 4-40 hPa: JOYCE vented 40 hPa  
JOYCE non vented 40 hPa  
JOYCE GEL vented  
JOYCE GEL non vented
- If third-party items are used, functional failures may occur and fitness for use may be restricted. Biocompatibility requirements may also not be met. Please note that in such cases, any claim under warranty and liability will be voided if neither the accessories nor genuine replacement parts recommended in the instructions for use are used.
- Patients with limited spontaneous respiration must be monitored at all times. The appropriate negative pressure/leak alarm on the ventilation device must be activated.
- Patients with obstructive and restrictive lung diseases should be examined on an individual basis during use of the JOYCE masks (for example with regard to re-inhalation of CO<sub>2</sub>).
- Do not close the openings in this mask except in situations described in these instructions for use. Otherwise there is a risk of re-inhaling CO<sub>2</sub> and of suffocation.
- Only put the nasal mask over your face for a prolonged period when the therapy device is switched on. Otherwise there is a risk of re-inhaling CO<sub>2</sub>.
- Do not use a breathing hose that is electrically conductive (exception: breathing tubes which can be heated electrically).
- Follow the instructions in section 4 "Hygiene treatment" to avoid infection, bacterial contamination and poor functionality of the mask.

### Side-effects

The following side-effects may occur:

nasal congestion, dry nose, dry mouth in the morning, sinus pressure, irritation of the conjunctiva, irritating noises during

exhalation, pressure sores on the face, reddening of the skin.

## 3 Operation

---

**Fig. A-E: see enclosed leaflet.**

### 3.1 Adjust forehead support (Fig. A)

#### Fine adjustment

Squeeze the upper side tabs of the coarse adjustment component inwards. Fit the fine adjustment component into one of the three notches.

To change the height of the forehead cushion, i.e. the fine adjustment (by approx. 1 cm), turn the fine adjustment component through 180°. Then position it as described above.

#### Coarse adjustment

Squeeze the lower side tabs of the coarse adjustment component inwards. Place it in one of the slots in the mask frame.

### 3.2 Putting on the mask

1. Unhook one of the lower headgear clips from the mask.
2. Guide the headgear over your head and hook the headgear clips back onto the strap.
3. Adjust the headgear with the aid of the hook-and-loop closures so that the mask seals on the face but is not too tight.

#### **WEINMANNstrap only:**

The WEINMANNstrap headgear provides the following additional options (see also "Overview" on the enclosed leaflet):

- hook-and-loop closures on both sides.
- hook-and-loop closure on top of the head: release the hook-and-loop closure, pull gently on the ends of the headgear with both hands and push the hook-and-loop closure back on.

You can tell whether the headgear is correctly attached by the hook-and-loop tabs. The hook-and-loop tabs on the forehead are the same color as the inside of the headgear.

#### **HEADstrap only:**

The HEADstrap headgear provides the following options for you to do this (see also "Overview" on the attached leaflet):

- coarse adjustments: slot "a" or "b" of the headgear.
  - fine adjustments: four hook-and-loop straps."
4. If necessary, re-adjust the forehead support.

### 3.3 Exhaled air system (Fig. E)

#### **vented variants**

vented variants (all mask components are transparent) have a so-called **integrated exhaled air system**: the retaining ring and mask frame are shaped in such a way that a gap is formed between these parts. Through this gap the exhaled air can escape. The special geometry of the parts guarantees that when putting on the retaining ring the gap cannot be closed. **Do not** connect an external exhaled air system, otherwise the therapy pressure will not be reached.

#### **non vented variants**

non vented variants (retaining ring, elbow and rotating sleeve in blue) have no exhaled air system. It is essential to connect an external exhaled air system. Otherwise there is a risk of reinhaling CO<sub>2</sub> and of asphyxia.

Follow the instructions for use of your exhaled air system.

### 3.4 Connecting the mask to the therapy device (Fig. E)

Push the hose of your therapy device onto the rotating sleeve (vented variants) or onto your external exhaled air system (non vented

variants). Ensure that the connection to the hose is secure and airtight. The connection must not be loosened by the therapy pressure. If you want to detach the mask temporarily from the hose, pull the rotating sleeve off the elbow.

### 3.5 Gel-filled mask components

On JOYCE GEL vented and JOYCE GEL non vented, the forehead cushion and the mask cushion are gel-filled as standard. Gel forehead cushions and gel mask cushions are also available as replacement parts for other JOYCE variants.

#### Caution!

Gel-filled mask components are highly sensitive components. Always handle them with care to prevent deformation and damage.

### 3.6 Therapy pressure ranges

The individual JOYCE mask variants are each available for a particular therapy pressure range (see "7. Scope of supply" on page 17):

- 4-30 hPa: JOYCE vented  
JOYCE non vented
- 4-40 hPa: JOYCE vented 40 hPa  
JOYCE non vented 40 hPa  
JOYCE GEL vented  
JOYCE GEL non vented

The JOYCE masks for the therapy pressure range 4 - 30 hPa have a soft mask cushion, the JOYCE masks for the therapy range 4 - 40 hPa have a harder mask cushion or a gel mask cushion (gel-filled variants).

The harder mask cushion for the therapy pressure range 4 - 40 hPa is marked with the number "40".

The soft (up to 30 hPa) and harder (up to 40 hPa) mask cushion and the gel-filled mask cushion (gel variants) are also available as individual replacement parts to enable you to adapt your mask to changed therapy requirements.

### 3.7 Ports for pressure measurement (Fig. C)

JOYCE masks are fitted with two ports that can be used to measure pressure or to feed in oxygen.

#### Caution!

- Oxygen therapy can have side-effects if incorrectly administered or used in the wrong dosage. Only carry out this sort of treatment if prescribed by a physician.
- When carrying out oxygen therapy, strictly follow the instructions for use of your O<sub>2</sub> system and your therapy device.

When the pressure measurement ports are not in use, always close them with the cap provided. Otherwise air will escape from the ports and the therapy pressure will not be reached.

### 3.8 Compatible devices

#### Therapy devices

JOYCE masks can be combined with all devices designed for carrying out the therapies described under "Intended Use".

Note that if unsuitable device combinations are used, the actual pressure in the mask may not correspond to the therapy pressure you have been prescribed.

This may also happen even if the correct pressure is indicated on the therapy device. Have your device combination adjusted by a physician or distributor so that the pressure in the mask corresponds to your therapy pressure.

#### Air humidifiers

JOYCE masks can be used with either cold-air or warm-air humidifiers.

Never set the humidification level so high that condensation forms in the air hose. Follow the instructions for use of your humidifier and therapy device.

### 3.9 After use

1. Unhook one of the lower headgear clips from the mask and remove the mask.

2. Loosen the headgear clips and remove the headgear from the mask.
3. Proceed as described in 4 "Hygiene treatment".

### 3.10 Disassembling/reassembling the mask

#### Disassembling the mask (Fig. D - A)

1. Detach the mask cushion from the mask frame.
2. Remove the cap from the pressure measurement ports.
3. Rotate the retaining ring counter-clockwise. Now pull the retaining ring and the elbow away from the mask frame.
4. Press the elbow through and out of the retaining ring and pull the rotating sleeve off the elbow.
5. Detach the coarse adjustment component from the mask frame.
6. Detach the fine adjustment component from the coarse adjustment component. Then detach the forehead cushion from the fine adjustment component.

#### Reassembling the mask (Fig. A - D)

1. Attach the forehead cushion to the fine adjustment component. Then attach the fine adjustment component to the coarse adjustment component.
2. Attach the coarse adjustment component to the mask frame.
3. Push the rotating sleeve onto the elbow until you hear it click into place.
4. Press the elbow into the retaining ring until you hear it click into place.
5. Attach the retaining ring and elbow onto the mask frame. Turn the retaining ring clockwise until you hear it click into place.
6. Push the cap onto the pressure measurement ports.

7. Attach the mask cushion to the mask frame.
8. Fasten the headgear to the mask with the clips.

## 4 Hygiene treatment

---

### 4.1 Schedule

#### Daily

Clean the mask in warm water to remove all residues. You can disassemble the mask for cleaning. Use a mild cleaning agent (e.g. dishwasher detergent). After cleaning, rinse all parts carefully with clean water.

#### Caution!

- For patients with a weakened immune system or a special medical history, **daily** disinfecting or sterilization of the mask components may be necessary. In this case, please consult your doctor.
- Please take great care when cleaning and disinfecting the contact surfaces between the mask frame and the retaining ring. Residues on these surfaces can reduce the success of the therapy and the functioning of the integrated exhaled air system (vented variants).

#### Weekly

Disassemble the mask and subject the parts of the mask to a hygiene treatment (see "4.2. Cleaning, disinfection, sterilisation" on page 15).

### 4.2 Cleaning, disinfection, sterilisation

#### Permissible procedures

#### Caution!

Always clean the components of the mask thoroughly before disinfecting or sterilising them as indicated in the above Table.

Procedure				
Components	Dishwasher	Hand wash	Disinfection*	Sterilisation (alternative to disinfection)
Mask cushion	•	•	•	•
Mask cushion GEL	•	•	•	
Forehead cushion	•	•	•	
Forehead cushion GEL		•	•	
Port cap	•	•	•	•
Mask frame	•	•	•	
Coarse adjustment component	•	•	•	
Fine adjustment component	•	•	•	
Elbow	•	•	•	
Rotating sleeve	•	•	•	
Retaining ring	•	•	•	
Headgear clips	•	•	•	
Headgear WEINMANNstrap		•		
Headgear HEADstrap		•		
Headgear JOYCEstrap		•		

Gel-filled mask components are highly sensitive components. Always handle them with care to prevent deformation and damage.

### Description of procedures

Dishwasher	Rinse cycle at 65 °C. Then rinse all parts carefully with clean water.
Hand wash	Hand wash in warm water. Use a mild detergent. Then rinse all parts in clean water and allow to air-dry.
Disinfection*	Clean mechanically, e.g. with a brush, in dilute solution. After disinfection rinse all parts thoroughly with distilled water and dry thoroughly.
Sterilisation (alternative to disinfection)	Steam sterilisation in equipment rated according to EN 285. Temperature 134 °C, minimum sterilisation time 3 minutes.

\*We recommend the disinfectants gigasept® FF (new), Cidex® OPA and Aniosyme DD1. Evidence of tolerance to these disinfectants has been provided. It is

essential to follow the manufacturer's instructions.

Notes on the hygiene treatment and on the hygiene treatment for a change of patient can be found in a brochure from our information portal at [www.weinmann.de](http://www.weinmann.de). If you do not have Internet access, contact Weinmann.

#### Note

Wear suitable gloves (e.g. disposable gloves) during the hygiene treatment.

Wash the HEADstrap headgear before using for the first time, as the colour might run.

The headgear may only be washed by hand!

Do not iron the headgear, otherwise the hook-and-loop closures will no longer work!

Do not dry the headgear in a dryer!

### 4.3 Following treatment

Thoroughly dry all parts following the hygiene treatment. Check the components for residual contamination or damage. Repeat the hygiene treatment if required. Replace damaged parts. Reassemble the mask.

### 4.4 Multiple patient use

If the mask is used by a different patient, it must be subjected to a hygiene treatment first.

Disassemble the mask. Disinfect or sterilize all parts as described in section 4.2 "Cleaning, disinfection, sterilisation". Use a new headgear.

## 5 Functional checks

Carry out a visual check after every hygiene treatment.

If the mask cushion is damaged and is no longer airtight, replace it with a new one.

Check the functionality and air-tightness of the ball-and-socket joint on the elbow. If it no longer moves easily or is no longer airtight, replace the elbow and retaining ring.

Parts that are no longer usable can be disposed of as normal household waste.

The forehead cushion and the mask cushion may become discolored, but this does not impair function.

## 6 Troubleshooting

Problem	Cause	Solution
Painful pressure on the face.	Mask is too tight.	Loosen headgear slightly. Adjust position of forehead support.
Draught in the eyes.	Mask is too loose.	Tighten headgear slightly. Adjust position of forehead support.
	Mask does not fit.	Use a different size of mask; contact your distributor.
	High therapy pressure.	Use firmer mask cushion (for pressures up to 40 hPa). Available from your distributor.
Therapy pressure is not reached.	Mask not adjusted correctly.	Adjust headgear. Adjust position of forehead support.
	Mask cushion is damaged.	Replace mask cushion.
	Hose system is leaking.	Check push-fit connector and fit of hose sleeves.
	Air is escaping from pressure measurement ports.	Cover pressure measurement ports with cap (included in package).

## 7 Scope of supply

### Note

All masks are supplied fully assembled and with headgear.

### vented variants

Article	Article no. with WEINMANNstrap/HEADstrap	Article no. with JOYCEstrap
JOYCE vented, Size S	WM 26110	WM 26110/1
JOYCE vented, Size M	WM 26120	WM 26120/1
JOYCE vented, Size L	WM 26130	WM 26130/1
JOYCE vented, Size XL	WM 26140	WM 26140/1
JOYCE vented 40 hPa, Size S	WM 26111	WM 26111/1
JOYCE vented 40 hPa, Size M	WM 26121	WM 26121/1
JOYCE vented 40 hPa, Size L	WM 26131	WM 26131/1
JOYCE vented 40 hPa, Size XL	WM 26141	WM 26141/1
JOYCE GEL vented, Size S	WM 26112	WM 26112/1
JOYCE GEL vented, Size. M	WM 26122	WM 26122/1
JOYCE GEL vented, Size. L	WM 26132	WM 26132/1

### non vented variants

Article	Article no. with WEINMANNstrap/HEADstrap
JOYCE non vented, Size S	WM 26160
JOYCE non vented, Size M	WM 26170
JOYCE non vented, Size L	WM 26180
JOYCE non vented, Size XL	WM 26190
JOYCE non vented 40 hPa, Size S	WM 26161
JOYCE non vented 40 hPa, Size M	WM 26171
JOYCE non vented 40 hPa, Size L	WM 26181
JOYCE non vented 40 hPa, Size XL	WM 26191
JOYCE GEL non vented, Size S	WM 26162
JOYCE GEL non vented, Size M	WM 26172
JOYCE GEL non vented, Size L	WM 26182

## Replacement parts

Article	Article no.
Ball-and-socket joint set: retaining ring, elbow, rotating sleeve	WM 26206
Ball-and-socket joint set non vented: retaining ring, elbow, rotating sleeve	WM 26207
Rotating sleeve vented	WM 26254
Forehead support: Forehead cushion, forehead support coarse adjustment component, forehead support fine adjustment component	WM 26201
Forehead support GEL: Forehead cushion GEL, forehead support coarse adjustment component, forehead support fine adjustment component	WM 26215
Forehead cushion	WM 26200
Forehead cushion GEL	WM 26209
Mask cushion 30 hPa, Size S	WM 26210
Mask cushion 30 hPa, Size M	WM 26220
Mask cushion 30 hPa, Size L	WM 26230
Mask cushion 30 hPa, Size XL	WM 26240
Mask cushion 40 hPa, Size S	WM 26211
Mask cushion 40 hPa, Size M	WM 26221
Mask cushion 40 hPa, Size L	WM 26231
Mask cushion 40 hPa, Size XL	WM 26241
Mask cushion GEL, Size S	WM 26212
Mask cushion GEL, Size M	WM 26222
Mask cushion GEL, Size L	WM 26232
WEINMANN <i>strap</i> headgear, incl. headgear clips	WM 26350
Headgear HEAD <i>strap</i>	WM 26360
JOYCE <i>strap</i> headgear, blue, incl. headgear clips	WM 26290
Headgear clips (4 pcs.)	WM 26205
Port cap	WM 26208
Instructions for use DE, EN, FR, NL, IT, ES, PT, TR, EL, PL, HR RU (pdf, Internet: www.Weinmann.de)	WM 16950 WM 16814

## Accessories

Article	Article no.
Exhaled air system: Silentflow 2 Noise reduction system	WM 23600 WM 23685

## 8 Specifications

### Technical data

JOYCE	Value
Product category according to 93/42/EEC:	II a
Dimensions (WxHxD):	80 x 130 x 90 mm <sup>1</sup>
Weight:	c. 85 grams <sup>1</sup>
Therapy pressure range: mask cushion soft: mask cushion hard: mask cushion GEL:	4 to 30 hPa 4 to 40 hPa 4 to 40 hPa
Hose connector: cone complying with EN ISO 5356-1 vented: non vented:	Ø 22 mm (male) Ø 22 mm (female)
Pressure measurement ports:	Ø 4 mm
Width of headgear clips:	max. 20 mm
Temperature range Operation: Storage:	+ 5 °C to + 40 °C - 20 °C to + 70 °C
Flow resistance vented: at 50 l/min at 100 l/min non vented: at 50 l/min at 100 l/min	0.02 hPa 0.11 hPa 0.09 hPa 0.37 hPa
Quoted two-figure noise emission value to ISO 4871: - sound pressure level - sound power level (uncertainty factor: 3 dB(A))	15 dB(A) 23 dB(A)
Working life	up to 12 months <sup>2</sup>
Standards applied	EN ISO 17510-2:2009

<sup>1</sup> depending on size and design

<sup>2</sup> If the materials used to produce masks are exposed to aggressive cleaning agents, for example, they age. In individual cases, therefore, you may need to replace your mask sooner (see "5. Functional checks" on page 17).

**CE 0197**

Subject to design changes


## Materials

Article	Material
Mask cushion, port cap	Silicone
Mask cushion GEL	Silicone, silicone gel
Forehead cushion, retaining ring	TPE, Polypropylene
Forehead cushion GEL	Silicone, silicone gel, polycarbonate
Mask frame, elbow	Polycarbonate
Headgear WEINMANN <i>strap</i>	Elasthane, Polyester, Polyurethane, Nylon, Cotton, Polyamide
Headgear JOYCE <i>strap</i>	Cotton, polyamide, elasthane
Headgear HEAD <i>strap</i>	Nylon/Spandex Polyurethane Polyester/Lycra
Rotating sleeve, forehead support coarse adjustment component, forehead support fine adjustment component	Polypropylene
Headgear clips	POM

All mask components are latex-free.

### Pressure/flow characteristic curve

The pressure/flow characteristic curve shows outlet flow as a function of therapy pressure.

## 9 Warranty

Weinmann gives the customer a limited manufacturer warranty on new genuine Weinmann products and any replacement part fitted by Weinmann in accordance with the warranty conditions applicable to the product in question and in accordance with the warranty periods from date of purchase as listed below. The warranty conditions can be downloaded from [www.weinmann.de](http://www.weinmann.de) on the Internet. We can also send you the warranty conditions on request. In the event of a claim under warranty, please contact your specialist dealer.

Product	Period of guarantee
Weinmann devices including accessories (except masks) for sleep diagnosis, sleep therapy, home ventilation, oxygen medicine and emergency medicine	2 years
Masks including accessories, rechargeable batteries, batteries*, sensors, tube systems	6 months
Disposable products	None

\*unless there is information to the contrary in the technical documentation.

## 10 Declaration of conformity

Weinmann Geräte für Medizin GmbH + Co. KG declares herewith that the product complies fully with the respective regulations of the Medical Device Directive 93/42/EEC.

The unabridged text of the Declaration of Conformity can be found on our website at [www.weinmann.de](http://www.weinmann.de)

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