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Mod. Econeb-jet

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
1 - Adult Mask
2 - Children mask
3 - Nasal Prong
4 - Mouthpiece
5 - Connector accessories
6 - Nebulizer Tank
7 - Pisper
8 - Lower chamber
9 - Air Tube

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Specifiche Tecniche/Technical specifications as provided by EN13544-1

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Results of measurements series with NEBJET nebulizer in compliance with EN 13544-1

Operating Flow l/min	Fill Volume ml	MMAD µm	GSD
5	2,5ml 2,5%NaF	2,61	2,75
7,5	2,5ml 2,5%NaF	2,03	2,97

MMAD = Mass Median Aerodynamic Diameter

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Results of INHALATION SIMULATOR measurements with NEBJET nebulizer in compliance with EN 13544-1

Operating Flow l/min	Fill Volume ml	Output µl	Output Rate µl/min
5	2ml 1%NaF	143	87
7,5	2ml 1%NaF	187	120

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Storage and transport conditions

Storage and transportation temperature: - 20°C (-4°F) ÷ +70°C (158°F)
Storage and transportation humidity percentage: 5% ÷ 95% RH

Operating data

Operating Minimum flow 5 l/min +/- 10%
Operating pressurewith nebulizer 60 kPa (0.6 bar, 8.7psi) +/- 10%
Operating Max flow 7.5 l/min +/- 10%
Operating Max pressure 150 kPa (1.5 bar, 21.75psi) +/- 10%
Min. capacity 2 ml
Max. capacity 12 ml

Read carefully instructions for use

0434 NEBJET is a Medical Device in conformity with 93/42/EEC Directive

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Mod. Econeb-jet

Use the unit as described in this manual, i.e. as an aerosol therapy system, following the instructions of your physician. Use for any purpose other than that for which is designed is to be considered dangerous; the manufacturer and distributor cannot be held liable for damage caused by improper, incorrect and unreasonable use or if the unit is connected to an electrical system that does not comply with current safety standards.

Econeb-jet is a unit for the treatment of respiratory apparatus. It is designed principally for domestic use.

Econeb-jet is equipped with a powerful piston compressor that does not need lubrication. The unit includes an air filter to be periodically replaced (every 70 applications). It is built in compliance with the current European Manufacturing and Safety Standards for electromedical equipments (EN 60601-1; EN13544-1).

Econeb-jet is supplied with the following Realcheck accessories: Nebjet nebulizer with mouthpiece, nasal prong, children mask, adult mask, air connection tube, filters (3pcs),

SAFEGUARDS

Econeb-jet is a medical device and has to be used only upon medical prescription and in accordance with this manual and the instructions of medication used. It is important that the user reads and understands the information provided for use of the device. Please contact the Technical service of the distributor.

Keep the unit and its accessories out of the reach of children, since some components are small enough to be swallowed and thus represents a choking hazard. Children and disabled persons should always use the unit under supervision of an adult who had read this manual.

- Never use an extension cord or electrical adapter and during use place unit on a stable, horizontal surface. Keep the cord unwind and away from heated surfaces. Never obstruct the air intakes.
- Do not handle the plug with wet hands, and not use the product while bathing. Never submerge the unit in water; should this occur, unplug immediately. Before touching a product fallen into water, unplug it and do not use it again.

Always unplug the power cord immediately after use.

Service should be performed only by specialized personnel. Unauthorised repairs invalidates the warranty.

At the end of its useful life, dispose of it in accordance with current regulations. (see the relevant chapter).

INSTRUCTIONS FOR USE AND CLEANING OF UNIT

Plug the power cord into the electrical outlet corresponding to the voltage of the unit as indicated in the data label under the unit.

To use the nebulizer and its accessories, see the relevant chapter.

When the therapy has finished, turn the unit off and carefully clean the accessories. Do not wash the unit under water or do not submerge it in water; wipe it only with a damp cloth using mild detergents, not abrasive and suitable for plastics.

Replacing the spare filter: lift the filter support lid (3) by levering with a flat-ended screwdriver as illustrated in figure 1a; take off the filter to be replaced and insert the spare filter as illustrated in figure 1b.

Finally, insert the filter support lid (3) into its place, making it adhere well to the surface. In compliance with EMC standard (EN60601-1-2).

Unit case not water – proof.

The unit should not be used in the presence of anaesthetic mixture inflammable with air, oxygen or nitrous oxide.

MICROBIC CONTAMINATION

A personal use of the accessories and the nebulizer is suggested, particularly in the presence of pathologies with the risk of infection and microbic contamination, (always consult your physician).

Designed and created in accordance with current standards, NEBJET is an Medical Device (Risk Class IIa - MDD 93/42/EEC) device suitable for aerosol therapy.

SAFETY WARNINGS

• A single-patient use is recommended. The device should be used only with nebulizers for aerosoltherapy that comply with standards in force and whose use is indicated in the instruction manual • Device not suitable for anaesthesia and lung ventilation
• The device should be used only with original accessories shown in this manual • Information concerning connection and use with devices for aerosol therapy are shown in the instruction manual of the devices.
• Use the device and its accessories as per physician recommendations • Use medications prescribed by physician according to his indications • Do not use the device for other use than Nebulizer for Aerosol therapy. Manufacturer will not be held responsible for inappropriate use of the device • The device does not come in sterile pack; always carry out cleaning and sterilization operations before and after each treatment. When using the device follow ordinary hygiene and cleaning precautions, (wash your hands and carry out cleaning operations as described in the specific section of this manual) • The device has small components which might be removed and easily swallowed. Use by minors and disabled people require presence of an adult with his faculties. Do not leave the device unattended in places easily accessible by minors and disabled people • Always store the device in a dry and clean place far from light, heat and weather • Dispose of device as provided by standards in force • In presence of pathologies with microbial contamination and infection hazard a single-patient use of accessories and nebulizer is recommended (always consult your doctor)

HOW TO PREPARE AND TO OPERATE THE DEVICE

The device is not sterile. Before use carry out cleaning and disinfection operations as described in the specific section

To pour the medication rotate the nebulizer top counter clockwise and remove it. Pour the quantity of medication prescribed by physician into the nebulizer tank. For proper atomization and to reduce the loss of Active Principle to minimum, the drug should diluted with a few ml of sterile physiological solution. Follow the instruction in the medication data sheet and the physician advice. IMPORTANT: If you overfill, empty the nebulizer, clean it thoroughly and proceed with new filling. When medication has been poured, connect the nebulizer top onto the tank, rotate clockwise and insert the accessory prescribed by your doctor into the nebulizer top. Firmly connect one end of the air tube to the nebulizer and the other end to the air outlet on the compressor nebulizer for aerosol therapy. Start to begin. During final nebulization phase we recommend shaking the nebulizer gently to favour the descent of non-nebulised particles from the walls to ptimise medication use.

IMPORTANT: NEVER INHALE IN HORIZONTAL POSITION. NEVER BEND THE NEBULIZER OVER 60°.

CLEANING INSTRUCTIONS

After each treatment clean thoroughly each component of the nebulizer removing medication residuals and possible impurities. Clean each component as described here. To clean the **nasal prong and the mouthpiece**, carefully wash under running water, keeping in mind that they can be disinfected by boiling (5-10 minutes) or using chemical agents available in pharmacies; following the instructions provided with product. Rinse thoroughly making sure that all deposits are washed away and let dry. To clean the **nebulizer**, disassemble in each part, wash with warm water and soap and rinse thoroughly. Remove excess water with a towel and allow to dry completely on a suitable surface. To clean **mouthpiece and mask**, wash them with warm water or boil them (5-10 min.) and allow to dry. Keep the accessories in a well-ventilated place protected from heat and sunlight.

NEBULIZER LIFE

A properly kept nebulizer will last for: 5 YEARS

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SSL HEALTHCARE ITALIA SPA-VIA M.E. LEPIDO 178/5 - 40132 BOLOGNA ITALY

WARRANTY CONDITIONS:

The unit is covered by a 2 year warranty against defects in materials and workmanship. The warranty doesn't extend to units that have undergone unauthorized repairs or been tampered with. The warranty period starts on the purchase date. It is mandatory to attach the sales receipt or invoice as proof of the purchase date to the warranty certificate below, which must be filled out in every part. Under the warranty, any components with manufacturing defects will be replaced and/or repaired free of charge. The purchaser will bear all expenses for the delivery to the distributor.

The warranty doesn't cover:

The accessories supplied with the unit or parts subject to normal wear. Repair of unsubstantiated defects. Repairs of units that have failed or been rendered defective as a result of improper use, carelessness or negligence or if the damage is not referable to the manufacturer (accidental drops, careless transport, etc.)

The warranty doesn't cover any direct or indirect damage or injury caused to persons, animals or property as a result of improper use of the product or during its period of inefficiency.

WHEN SUBMIT A WARRANTY CLAIM, IN OBSERVANCE OF THE ABOVE CONDITIONS, SEND THE UNIT, COMPLETE AND CLEAN TOGETHER WITH THE CERTIFICATE BELOW AND THE PROOF OF PURCHASE TO:

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The household rubbish box symbol crossed showed on pack means that the product, at the end of his life, must be collected separately from the other rubbish, therefore the costumer should dispose them to the specific collection points provided for this product; or given to the retailer in the same moment you buy another similar product, in reason of one to one. The suitable waste recycling, the treatment and disposal compatibly with environment disposal, avoid negative effects on environment, on health and encourage the recycling of raw materials. The non-accomplishment of these rules and illegal disposal could be subject of administrative sanctions by competent authorities according to this law.

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