INNO-LiPA HLA-B Update Plus

IVD

Manufactured by:

INNOGENETICS N.V.
Technologiepark 6
9052 Ghent
Belgium
☎+32-9 329 13 29

Distributed by:

INNOGENETICS GmbH
Lembecker Straße 19
46359 Heiden (Westfalen)
Germany
☎+49-2867 99 07 0

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English

Intended use

The INNO-LiPA HLA-B Update Plus is a line probe assay, for *in vitro* use, designed for the molecular typing of the human leukocyte antigen (HLA) B alleles at the allele group level.

Test Principle

The INNO-LiPA HLA typing tests are based on the reverse hybridization principle. Amplified biotinylated DNA material is chemically denatured, and the separated strands are hybridized with specific oligonucleotide probes immobilized as parallel lines on membrane-based strips. This is followed by a stringent wash step to remove any mismatched amplified material. After the stringent wash, streptavidin conjugated with alkaline phosphatase is added and bound to any biotinylated hybrid previously formed. Incubation with a substrate solution containing a chromogen results in a purple/brown precipitate. The reaction is stopped by a wash step, and the reactivity pattern of the probes is recorded.

With the INNO-LiPA HLA-B Update Plus, an amplification kit (INNO-LiPA HLA-B Multiplex Plus) is available for standardized preparation of biotinylated amplified material. The amplification kit is based on the polymerase chain reaction (PCR*).

Amplification products are subsequently hybridized using 2 typing strips on which 66 sequence-specific probe lines and 2 control lines are fixed (Figure 1).

* Use of this product is covered by a license from F. Hoffmann - La Roche Ltd. and Roche Molecular Systems, Inc.

INNO-LiPA HLA-B Update Plus: steps involved

- step 1 Multiplex amplification of exons 2 to 4 of the HLA-B alleles.
 step 2 Hybridization and stringent wash with 67 probes immobilized on two INNO-LiPA HLA-B Update Plus strips (56°C).
 step 3 Color development.
 step 4 Interpretation of the probe reactivity pattern.
 step 5 If required, the LiRAS™ interpretation software advises for a further amplification and hybridization step.

Using the INNO-LiPA HLA-B Update Plus kit, it is advised to start with (a) and proceed to (b), if ambiguities involving Bw4/Bw6 alleles occur. This will be advised by the LiRAS™ Interpretation Software.

- (a) the amplification of exons 2, 3 and 4 of the HLA-B alleles with the HLA-B Multiplex primer solution prior to analysis with the INNO-LiPA HLA-B Update Plus strips 1 and 2.
 (b) group-specific amplification of exon 2 of the HLA-Bw4 alleles with primer solution HLA-Bw4 prior to analysis with the INNO-LiPA HLA-B Update Plus strip 1.
 The HLA-Bw4 primer solution provided, amplifies exon 2 of the HLA-Bw4 alleles. Exon 2 of other HLA-B alleles will not be amplified, except for the alleles B*4411, B*4703 and B*5708.

The INNO-LiPA HLA-B Update Plus is designed to give the best possible resolution, using the reverse hybridization assay format, at the allele group level (this means the first two digits after the asterisk in an allele name when following standard HLA nomenclature e.g. HLA-B*07).

In addition to the allele group interpretation, the software will give all possible allele combinations. This information however, should be considered as extra and not as decisive.

Reagents**Description, preparation for use and recommended storage conditions**

- If kept at 2 - 8°C, opened and unopened, and stored in the original vials, the reagents are stable until the expiry date of the kit. Do not freeze reagents. Do not use the reagents beyond the expiry date.
- The kit should be stored isolated from any source of contaminating DNA, especially amplified products.
- All reagents should be brought to room temperature (20 - 25°C) approximately 60 minutes before use and should be returned to the refrigerator immediately after use.
- Alterations in physical appearance of the kit components may indicate instability or deterioration.
- To minimize the possibility that strips curl before use, it is recommended to store the tube horizontally.

Reagents supplied:

Component	Quantity	Ref.	Description
Strips 1	1 x 20	58188	Containing 20 strips for INNO-LiPA HLA-B Update Plus marked with a brown marker line.
Strips 2	1 x 20	58189	Containing 20 strips for INNO-LiPA HLA-B Update Plus marked with a chrome colored marker line.
LiPA Control	1 x 0.05 ml	58190	Containing ethylenediaminetetraacetic acid (EDTA) and 0.01% methylisothiazolon (MIT) 0.1% chloroacetamide (CAA) as preservative.

Denaturation Solution	1 x 1 ml	56718	Alkaline solution containing EDTA. The vial should be closed immediately after use; prolonged exposure of this solution to air leads to a rapid deterioration of the denaturing strength.
Hybridization Solution	2 x 80 ml	56719	Saline sodium phosphate EDTA (SSPE) buffer containing 0.5% sodium dodecyl sulphate (SDS). The Hybridization Solution should be pre-warmed to a temperature of at least 37°C and must not exceed 56°C (all crystals should be dissolved before use).
Stringent Wash Solution	2 x 200 ml	56720	SSPE buffer containing 0.1% SDS. The Stringent Wash Solution should be pre-warmed to a temperature of at least 37°C and must not exceed 56°C (all crystals should be dissolved before use).
Conjugate Diluent	1 x 150 ml	55671	Phosphate buffer containing NaCl, Triton®, protein stabilizers, and 0.01% MIT/0.1% CAA as preservative.
Conjugate 100x	1 x 1.5 ml	56716	Streptavidin labeled with alkaline phosphatase in Tris buffer containing protein stabilizers and 0.01% MIT/0.098% CAA as preservative, to be diluted 1/100 in Conjugate Diluent before use. Prepare 2 ml Conjugate working solution for each test trough + 2 ml in excess (Conjugate working solution can be prepared during stringent wash) for manual testing. For <i>Auto</i> -LiPA, make 10 ml in excess. The Conjugate working solution is stable for 24 hours at room temperature (20 - 25°C) if stored in the dark.
Substrate Buffer	1 x 235 ml	55732	Tris buffer containing NaCl, MgCl ₂ and 0.01% MIT/0.1% CAA as preservative.
Substrate BCIP/NBT 100x	1 x 1.5 ml	55723	5-Bromo-4-chloro-3-indolyl phosphate p-toluidine salt (BCIP) and 4-nitroblue tetrazolium (NBT) in dimethylformamide (DMF), to be diluted 1/100 in Substrate buffer before use. Prepare 2 ml Substrate working solution for each test trough + 2 ml in excess (Substrate working solution can be prepared during conjugate incubation) for manual testing. For <i>Auto</i> -LiPA, make 10 ml in excess. The Substrate working solution is stable for 24 hours at room temperature (20 - 25°C) if stored in the dark.
Rinse Solution 5x	2 x 80 ml	56721	Phosphate buffer containing NaCl, Triton®, and 0.05% MIT/0.48% CAA as preservative, to be diluted 1/5 (1 part + 4 parts) in distilled or deionized water before use. Prepare 8 ml Rinse working solution for each test trough + 10 ml in excess. For <i>Auto</i> -LiPA, make 20 ml in excess. The Rinse working solution is stable for 2 weeks at 2 - 8°C.
Incubation tray	5	-	Containing 8 troughs each.
Reading card	1	-	For identification of positive probes.
Data reporting sheet	2	-	For storage of developed strips.

Materials required but not provided

- Water bath with shaking platform (80 rpm; with inclined lid; temperature adjustable to 56°C ± 0.5°C).
- Aspiration apparatus.
- Calibrated thermometer.
- Orbital, reciprocal or rocking platform shaker.

Recommendations:**For an orbital shaker:**

- the diameter of the circular motion should be equal or superior to 13 mm
- recommended speed for a 13 mm circular motion is 160 rpm.

For a reciprocal shaker:

- recommended speed for the to and from motion is 80 movements per minute.

For a rocking platform shaker:

- the shaking angle should not exceed 13° to avoid spilling of liquid
- recommended speed is 50 rpm.

- Vortex mixer or equivalent.
- Graduated cylinders (10, 25, 50, and 100 ml).
- Distilled or deionized water.
- Disposable gloves.
- Disposable sterile pipette tips (preferably cotton-plugged).
- Forceps for strip handling.
- Adjustable pipettes to deliver 1 - 20 µl, 20 - 200 µl, and 200 - 1000 µl.
- Dispensing Multipipette (Eppendorf, optional).
- Timer, 2 hours (± 1 minute).

Safety and environment

- **Please refer to the manufacturer's safety data sheet and product labeling for information on potentially hazardous components.**

R20/21, R36, R45, R61, S36/37, S45-53

Toxic! (T) Harmful by inhalation and in contact with skin. Irritating to eyes. May cause cancer. May cause harm to the unborn child. Wear suitable protective clothing and gloves. In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible). Avoid exposure - obtain special instructions for use. Restricted to professional users.

Contains Dimethylformamide, 5-Bromo-4-chloro-3-indolyl phosphate p-Toluidine salt:
Substrate BCIP/NBT 100x.

R43, S24-37

Irritant! (Xi) Avoid contact with skin. May cause sensitization by skin contact. Wear suitable gloves.

Contains 2-Chloroacetamide: Rinse Solution, Substrate Buffer, Conjugate Diluent and LiPA Control.

R36/38, S23-24-26

Irritant! (Xi) Irritating to eyes and skin. Do not breathe vapour. Avoid contact with skin. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.
Contains sodium hydroxide: Denaturation Solution.

- Specimens should be handled as potentially infectious. Therefore, all blood components and biological materials should be considered as being potentially infectious and should be handled as such. Only adequately trained personnel should be permitted to perform the test procedure. All blood components and biological materials should be disposed of in accordance with established safety procedures.
 - Autoclave for at least 15 minutes at 121°C.
 - Incinerate disposable material.
 - Mix liquid waste with sodium hypochlorite so that the final concentration is $\pm 1\%$ sodium hypochlorite. Allow to stand overnight before disposal.
CAUTION: Neutralize liquid waste that contains acid before adding sodium hypochlorite.
- Use of personal protective equipment is necessary: gloves and safety spectacles when manipulating dangerous or infectious agents.
- Waste should be handled according to the institution's waste disposal guidelines. All federal, state, and local environmental regulations should also be observed.

Specimens

Since the LiPA test utilizes biotinylated amplified DNA material as specimen, an amplification kit, INNO-LiPA HLA-B Multiplex Plus, is available as an accompanying tool.

Manipulations procedures

Strip handling

- The strips are designed to be used only once!
- Do not touch the strips with bare hands; use clean forceps.
- Use a **pencil** for identification of the test strips. Do not use ballpoints, etc. Write the ID above the marker line on the strips.
- Throughout the different incubation steps, test strips should always remain in the same trough.
- Unused or developed strips should be kept away from intense light and heat.
- Allow the developed strips to dry completely before interpretation, covering, and storing.
- Developed dry strips should be stored preferably in the dark at 20 - 25°C.
- Do not reuse the troughs.

Directions for manual incubation

- Incubation at 56°C \pm 0.5°C during hybridization and stringent wash are the most critical steps in avoiding false positive (temperature too low) or false negative/very weak signals (temperature too high). A shaking water bath with inclined lid allows a good control of temperature variations. Strict temperature control with a calibrated thermometer is necessary.
- Always **close** the lid of the water bath during incubations in order to avoid false positive signals.
- **Do not use a hot air shaker** for the hybridization and stringent wash.
- For hybridization and stringent wash, the troughs should be placed on the shaking platform of the water bath. Adjust the water level between 1/3 and 1/2 of the height of the trough. Make sure that the troughs do not float on the water. The water should be in direct contact with the troughs.
- Incubation steps for the color development should be between 20 and 25°C. If the temperature is below 20°C, weaker results may be obtained. If the temperature is above 25°C, high background and/or false positive signals may be obtained.
- Always incubate exactly for the duration as mentioned in the protocol.
- The amplitude of the motion generated by both the shaking water bath (hybridization procedure) and the shaker (color development procedure) is critical in achieving maximum sensitivity and homogeneous staining. The strip surface should be completely submerged.

The amplitude should be as high as possible. **However, spilling of liquid over the edges of the troughs should be avoided!** This can lead to cross-contamination and invalid results.

- Shaking during incubation of the strips should be performed in such a way that both the liquid and the test strips move back and forth in the trough, without liquid being spilled over the edge of the troughs.
- Do not cover the tray. During hybridization and stringent wash incubations, the troughs can be left uncovered in the water bath. Covering the troughs with microplate sealers may cause cross-contamination.

Directions for manual changing liquid in the troughs

- The liquid is aspirated from the trough with a pipette, preferably attached to a vacuum aspirator. The tray is held at an angle to allow all liquid to flow to one end of the trough.
- Add 2 ml of the appropriate solution to each trough and follow the protocol. A dispensing Multipipette (Eppendorf) is useful for this purpose.
- Repeat this step as many times as indicated in the test procedure.

NOTE:

- Do not allow the strips to dry between two steps.
- Make sure not to damage the surface of the strips when aspirating. Aspirate the liquid at the top of the strip above the marker line.
- Make sure all liquid is aspirated.
- Make sure the whole strip is thoroughly washed by complete submersion in the solution.
- Adapt the speed of the shaker when necessary.

Remarks and precautions

- For professional use only.
- In order to avoid DNA contamination, a maximum physical separation between the pre- and post-amplification steps is recommended: separate rooms, separate pipettes and other lab material, separate lab coats and gloves (and their stock) are minimum precautions for a good laboratory practice.
- Avoid any return from the post-amplification room to the pre-amplification room.
- The use of autoclaved disposable lab material is recommended.
- Do not reuse disposable lab material.
- Use a new sterile pipette tip for each aliquoted specimen.
- Do not use the reagents beyond the expiry date.
- Do not mix reagents between kits, unless the components have identical lot numbers.
- Avoid microbial contamination of reagents.

Manual test procedure

NOTE:

- Throughout the different incubation steps, the test strips should always remain in the same trough.

Before incubation, check the temperature of the water bath using a calibrated thermometer, and adjust the temperature if necessary before placing the tray in the water bath. Always close the lid.

The LiPA Control should be included in each test run.

Samples

1. HLA-B amplified product (use 10 μ l). For sample preparation: see INNO-LiPA HLA-B Multiplex Plus instructions for use.
NOTE:
 - Make sure exactly 10 μ l of amplified sample is added. Too much or too little sample may lead to a wrong typing result.
2. LiPA Control sample for INNO-LiPA HLA-B Update Plus (use 10 μ l; **amplification is not required!**).
3. Blank amplified control sample (negative control; use 10 μ l).

Denaturation and hybridization

1. Heat a shaking water bath to $56^{\circ}\text{C} \pm 0.5^{\circ}\text{C}$. Check the temperature using a calibrated thermometer, and adjust the temperature if necessary. Do not exceed the indicated temperature. Pre-warm the Hybridization Solution and Stringent Wash Solution in a water bath of at least 37°C and do not exceed 56°C . Mix before use. All crystals should be dissolved.
2. Using forceps, remove the required number of INNO-LiPA HLA-B Update Plus Strip 1 and Strip 2 from the tube (one strip 1 and one strip 2 per test sample for the INNO-LiPA HLA-B Multiplex Plus, one strip 1 for the INNO-LiPA HLA-Bw4 amplicon). Include both strips for the LiPA Control sample for INNO-LiPA HLA-B Update Plus and for the blank amplified control sample. Pencil an identification number above the brown/chrome colored marker line on the strip.
3. Take the required number of test troughs (two troughs per test sample for HLA-B Multiplex Plus and one trough per test sample for HLA-Bw4) and place them into the tray.
4. Pipette 10 μ l Denaturation Solution into the upper corner of each trough.
NOTE:
 - Close the vial immediately after use.
5. Add 10 μ l sample (see Samples; see Instructions for use of INNO-LiPA HLA-B Multiplex Plus) and carefully mix by pipetting up and down. Always use sterile pipette tips. Allow denaturation to proceed for **5 minutes** at $20 - 25^{\circ}\text{C}$.
6. Shake the pre-warmed Hybridization Solution and **gently** add 2 ml to the denatured amplified product into each trough. Mix by gentle shaking. Take care not to contaminate neighboring troughs during pipetting.
7. Immediately place the strip into the trough. The strips should be completely submerged in the solution.
NOTE:
 - Wear disposable gloves and use forceps.
8. Place the tray into the $56^{\circ}\text{C} \pm 0.5^{\circ}\text{C}$ shaking water bath (80 rpm; see Directions for manual incubation), close the lid, and incubate for **30 minutes**.
NOTE:
 - Avoid splashing water from the water bath into the trough. Adjust the water level between 1/3 and 1/2 of the height of the trough. To prevent the tray from sliding, immobilize the tray between two heavy weights.

Stringent wash

1. After hybridization, remove the tray from the water bath.
2. Hold the tray at a low angle and aspirate the liquid from the trough with a pipette, preferably attached to a vacuum aspirator. Add 2 ml pre-warmed Stringent Wash Solution into each trough, place the tray into the $56^{\circ}\text{C} \pm 0.5^{\circ}\text{C}$ shaking water bath, close the lid and incubate for **3 minutes**.
3. After incubation, remove the tray from the water bath and aspirate the liquid from the trough.
4. Repeat the **3 minutes** incubation step and the aspiration step once.

- At last, add 2 ml pre-warmed Stringent Wash Solution into each trough and place the tray into the $56^{\circ}\text{C} \pm 0.5^{\circ}\text{C}$ shaking water bath, close the lid and incubate for **4 minutes**.

Before incubation, check the temperature of the water bath using a calibrated thermometer, and adjust the temperature if necessary. Always close the lid.

NOTE:

- Dilute the concentrated Rinse Solution 5x and Conjugate 100x during stringent wash. See Reagents.

Color development

All subsequent incubations are carried out at **20 - 25°C on a shaker**.

During the incubations, the liquid and test strips should move back and forth in the trough for homogeneous staining.

- Wash each strip twice for 1 minute using 2 ml of the Rinse working solution (see Directions for changing liquid in the troughs).
- Add 2 ml of the Conjugate working solution to each trough and incubate for **30 minutes** while agitating the tray on the shaker.

NOTE:

- Dilute the Substrate BCIP/NBT 100x solution about 10 minutes prior to the end of the conjugate incubation. See Reagents.
- Wash each strip twice for **1 minute** using 2 ml of the Rinse working solution and wash once more using 2 ml Substrate Buffer.
 - Add 2 ml of the Substrate working solution to each trough and incubate for **30 minutes** while agitating the tray on the shaker.

CAUTION:

- Wear gloves and protective goggles.
- Stop the color development by washing the strips twice in 2 ml distilled water while agitating the tray on the shaker for at least **3 minutes**.
 - Using forceps, remove the strips from the troughs and place them on absorbent paper. Let the strips dry completely before reading the results. Store the developed and dried strips in the dark.

Automated test procedure: *Auto*-LiPA

The LiPA test procedure is extremely well-suited for automation. Therefore, the *Auto*-LiPA was designed to fully handle hybridization, stringent wash and color development steps. The *Auto*-LiPA is featured as a walk-away system with automated heating and cooling, and with automated aspiration and pipetting.

The standard *Auto*-LiPA 48 protocol describes the testing of up to 24 samples when using a single strip per trough. However, this number can be increased to a maximum of 48. To perform this, two strips from the same locus are tested in the same trough, while keeping all other protocol parameters the same. Note also that weaker probe reactivities can occasionally be observed. Therefore, those customers who intend to increase the number of samples are advised to validate this '2 strips in one trough' procedure in his/her laboratory.

For more information and specific protocols on the *Auto*-LiPA, please contact your local distributor.

Results

Reading

Figure 1 illustrates the position of the different oligonucleotide probes on the INNO-LiPA HLA-B Update Plus strips. A line is considered positive when a clear purple/brown band appears at the end of the test procedure.


Figure 1: Location of the marker line (brown line on Strip 1 and chrome line on Strip 2), the conjugate control line (conj. control), the HLA-B Update Plus control line (HLA-B control), and the 66 sequence-specific DNA probes on the INNO-LiPA HLA-B Update Plus strips.

Validation

- Include one negative control and the LiPA Control (for INNO-LiPA HLA-B Update Plus) each time a test is performed. As with any new test procedure, the inclusion of additional positive and negative controls is to be considered until a high degree of confidence is reached in the ability to correctly perform the test procedure. It is suggested that reference DNA should be used to validate that the proper test conditions have been established.
- The uppermost lines on the INNO-LiPA HLA-B Update Plus strips (Figure 1) are the marker lines (brown and chrome colored line). These lines allow correct orientation of the strips.
- The next line controls the addition of reactive Conjugate and Substrate working solution during the detection procedure. This line should be lined up with the conjugate control line on the plastic reading card. This line should always be positive (if not: no interpretation should be done!) and should have approximately the same intensity on each strip in the same test run.
- The next line is the HLA-B Update Plus internal hybridization control line (HLA-B control on the reading card) which indicates whether or not an appropriate amount of HLA-B amplified material was added for hybridization. This line should always be positive, except for the negative control of the HLA-B multiplex amplification, but its intensity may vary between samples.
- The assay result of each negative control should give no apparent signal for any of the lines on the strip, except for the HLA-Bw4 amplification where the internal hybridization control line (HLA-B Control) will be positive as well.
- The **LiPA Control** sample for INNO-LiPA HLA-B Update Plus should yield the following pattern of positive probes: strip 1: conjugate control, HLA-B Control, 22 and 35; strip 2: conjugate control, HLA-B control, 53, 56 and 57. Every other probe should be negative. This control sample is composed of biotinylated oligonucleotides, complementary to the mentioned hybridization probes.

Positive hybridization on these probes thus demonstrates satisfactory performance of the assay including hybridization, stringent washing, color development, and detection. A different pattern may indicate assay performance problems such as incorrect temperature during hybridization or deviations from the prescribed incubation times or temperatures.

- Color intensities of probes on a strip may differ from one line to another.

Interpretation of results

- Check for the correct reactivity pattern of the LiPA Control sample.
- The conjugate control line on the strip should be lined up with the conjugate control line on the plastic reading card.
- First check the positivity of the control lines (first two lines) in order to validate each individual strip (see Validation).
- Typing results are based on the reactivity of the probes in the kit. A list of probe specificities is available.
- Identify all probe numbers that are positive on the INNO-LiPA HLA-B Update Plus strips and deduce the HLA-B type by using the INNO-LiPA HLA-B Update Plus typing table or a version of the LiPA interpretation software (see Interpretation software: LiRAS™).

The software often gives more results than the intended use suggests: interpretation at allele group is given, as well as all possible allele combinations. This information however, should be considered as extra, not as decisive.

Interpretation software: LiRAS™

The LiRAS™ software for INNO-LiPA HLA is designed to assist with the interpretation of the LiPA results.

Please contact your local distributor to obtain the latest updated version.

Limitations of the procedure

- Use of this product should be limited to personnel trained in the techniques of hybridization.
- Only good laboratory practice and careful performance of the specified procedures, will allow specific hybridization and correct typing of target DNA.
- The INNO-LiPA HLA-B Update Plus strips have been designed to give resolution at allele group level (B*07 to B*83). However, a combination of certain alleles could lead to a non-unique probe combination and thus to an ambiguous answer of two or more possible allelic combinations.
- New alleles might have polymorphisms outside probe regions, this should be taken into consideration when interpreting the results.
- The HLA-B Update Plus control lines are not only specific for HLA-B, but also for some HLA-A and some HLA-C alleles. Therefore, the HLA-B control lines might also react with HLA-A and HLA-C amplified products.

Test performance

Test performance using DNA amplified with the HLA-B Multiplex Primer Solution

The performance of INNO-LiPA HLA-B Update was evaluated in 4 European tissue typing laboratories on routine typing samples. A total of 347 samples were analyzed. DNA was extracted from fresh or frozen EDTA-anticoagulated or citrate-anticoagulated blood. DNA extraction methods included salting out, QIAamp 96 DNA blood kit (Qiagen), Genomic DNA Purification kit (Promega) and Puregene (Gentra Systems). Fresh or frozen DNA was amplified using the Perkin Elmer 9600 or 9700 thermal cycler. Amplicons were stored at 2 - 8°C or -15/-25°C before they were applied on the strips.

No amplification failures were encountered, although repeat testing was needed for 2 samples from one center.

The samples were processed using the *Auto*-LiPA instrument and interpreted with a clinical trial version of the LiRAS™ interpretation software.

Accuracy

An HLA-B typing result was obtained for 342 of the 347 samples. Repeat testing was needed for 4 samples from one center. One sample resulted in 'no typing deducible' but according to the investigator, this was probably due to contamination of the DNA. For 4 samples from one center, no typing result was provided by the investigator, although options were available via editing for one of them.

Accuracy at allele group level was calculated as a percentage of observed concordance with alternative DNA assays (Biotech ELPHA HLA-AB LowRes Typing kit, Dynal SSP, LiPA HLA-B, sequencing, Micro SSP™, OLERUP SSP™, Biotech HLA-B SSP). The INNO-LiPA HLA-B Update results are considered concordant with the reference assay result if the two results were identical, if the INNO-LiPA HLA-B Update results were more specific than the reference result, or if the INNO-LiPA HLA-B Update results showed a lower resolution but correct result by comparison to the reference methods. Accuracy was 100% (342/342).

Resolution

(typing table version v1.0/2001-03-31)

The INNO-LiPA HLA-B Update has been designed to resolve HLA-B alleles to the allele group level, i.e. both at serological broad and split level. 'Allele group' is equivalent to the designation of the first two digits following the asterisk in standard HLA nomenclature (e.g. HLA B*07). Furthermore, the kit is designed to detect the following null alleles: 0808N (exon 3), 1526N (exon 3) and 5111N (exon 4). The resolution at allele group level observed during this study is 92.4% (316/342).

The ambiguities observed were B*08xB*51 or B*08xB*78 (8 cases), B*15xB*35 or B*15xB*15 (2 cases), B*35xB*38 or B*39xB*53 (2 cases), B*08xB*38 or B*08xB*39 (2 cases), B*15xB*51 or B*15xB*52 or B*15xB*78 (2 cases), B*35xB*51 or B*53xB*78 (2 cases), B*35xB*38 or B*35xB*53 (1 case), B*35xB*49 or B*50xB*53 (1 case), B*35xB*47 or B*53xB*47 (1 case), B*40xB*78 or B*40xB*52 (1 case), B*35xB*55 or B*35xB*56 (1 case), B*15xB*49 or B*15xB*50 (1 case), B*49xB*78 or B*50xB*51 (1 case), B*15xB*51 or B*52xB*78 (1 case). For six of these cases, both the reference assay and INNO-LiPA HLA-B Update had an ambiguous result. Two ambiguities (6 samples) observed with the reference method were reduced to a single clear-cut typing result using the INNO-LiPA HLA-B Update assay.

Probe reactivity

Four probes reacted false positive and nine probes reacted weak to false negative at one or more occasions during this clinical trial. This information, if not already present, has been implemented in the software and/or guidelines. The false reactivities were all identified as such and concordant typing at allele group level was achieved when the LiRAS™ interpretation software was used.

Consult the guidelines for specific information on probe reactivity.

A complete list of probe and primer specificities is available. All probes are designed to hybridize specifically with their complementary sequence at the level of one mismatch, unless stated otherwise in the list.

Test sensitivity

Results obtained during an internal evaluation of a dilution series of 3 DNA samples, varying from 0.01 µg/µl to 2 µg/µl, recommended the use of DNA with a minimum concentration of 0.1 µg/µl.

The external performance evaluation included 347 DNA samples with a concentration between 20 ng/µl and 1500 ng/µl and with a purity (A_{260nm}/A_{280nm}) ≥ 1.5 .

The DNA was diluted to a concentration between 20 ng/µl and 250 ng/µl prior to

amplification. No amplification failures were encountered, although repeat testing was needed for 2 samples from one center.

Precision

Internal data:

Three samples were tested in quadruplicate, by two different persons using the manual method, or by one person on two different *Auto*-LiPA instruments. Two samples were analyzed on three different lots in two different runs. Interpretation of the results was the same for all samples. Probe reactivities were comparable in all cases.

External validation:

- Inter-lot variability was assessed by analyzing a proficiency panel (5 samples with known reactivity and interpretation) with two different batches of product. For each sample the same typing result was obtained, independent of the batch used.
- Inter-lab variability was assessed by having the same proficiency panel tested by the four centers. Each center obtained the same typing results for the five samples.

Test performance using DNA Amplified with the HLA-Bw4 Primer Solution

The performance of the HLA-Bw4 Primer Solution of INNO-LiPA HLA-B Multiplex Plus was evaluated in-house and in two Belgian tissue typing laboratories. A total of 82 anonymous samples was analyzed (50 samples externally; 32 in-house). All samples had previously been typed with INNO-LiPA HLA-B Update. The samples included 31 ambiguous (18 different ambiguities) and 51 unambiguous samples. All ambiguous samples were further typed to a unique result with one or more alternative DNA-based typing assays.

All hybridizations were performed using the *Auto*-LiPA instrument, using program HLA56v3, and the interpretation was done with a clinical trial version of the LiRAS™ for LiPA HLA v3.00 interpretation software.

Accuracy

Calculations of diagnostic accuracy were only carried out on samples that were successfully amplified and for which a unique reference result was available. As a consequence, two samples were excluded from this calculation. A unique HLA-B typing result was obtained for the remaining 80 samples.

Diagnostic accuracy of the assay at allele group level was calculated as the number of correctly obtained typing results, compared to the unique typing result obtained with INNO-LiPA HLA-B Update (Innogenetics) or an alternative typing assay (Olerup SSP (GenoVision), SBT (Abbott), Biotest HLA-B SSP (Biotest), Allset+ HLA-B Low Resolution (Dynal), and sequencing (in-house method)).

The typing results of all 80 samples were concordant with the reference assay, resulting in a diagnostic accuracy of 100% (80/80; 95% CI [95.5%; 100%]).

Resolution (typing table v.1.4/031010)

The HLA-B Multiplex Primer Solution, used in combination with the INNO-LiPA HLA-B Update strips, is designed for typing HLA-B alleles at allele group level.

The HLA-Bw4 Primer Solution, used in combination with the INNO-LiPA HLA-B Update strip 1, is designed to resolve the most frequent ambiguities, ie, ambiguities involving Bw4/Bw6 alleles.

All ambiguities, at allele group level, observed after testing with the HLA-B Multiplex Primer Solution, were reduced to an unambiguous typing result using the HLA-Bw4 Primer Solution: B*08, B*38 or B*08, B*39 (2 cases); B*08, B*51 or B*08, B*78 (4 cases); B*13, B*49 or B*13, B*50 (2 cases); B*15, B*38 or B*15, B*39 (1 case); B*15, B*51 or B*15, B*52 or B*15, B*78 (1 case); B*15, B*51 or B*15, B*52 or B*15, B*78 or B*35, B*52 (1 case); B*35, B*38 or B*39, B*53 (4 cases); B*35, B*40 or B*40, B*53

(1 case); B*35, B*47 or B*47, B*53 (1 case); B*35, B*49 or B*50, B*53 (2 cases); B*35, B*52 or B*53, B*78 (2 cases); B*38, B*44 or B*39, B*44 (1 case); B*38, B*50 or B*39, B*49 (1 case); B*38, B*78 or B*39, B*51 (1 case); B*40, B*49 or B*40, B*50 (2 cases); B*41, B*53 or B*42, B*53 (1 case); B*44, B*49 or B*44, B*50 (2 cases).

For all ambiguous results, the LiRAS™ interpretation software advised Bw4 amplification followed by hybridization with the INNO-LiPA HLA-B Update strip 1.

Probe reactivity

On one or more occasions during this study, some probes gave weak reactions. Concordant typing at allele group level was achieved when the LiRAS™ interpretation software was used.

Test sensitivity

Test sensitivity of the assay at allele group level was calculated as the number of successful typing results over the total number of samples tested.

Test sensitivity in this evaluation, after one repeat amplification testing of four samples, was 98.8% (81/82; 95% CI [93.4%; 100%]). The sample for which amplification failed after retesting at the center, was successfully amplified in-house and a correct typing result was obtained.

Precision

Repeatability (intra-assay), inter-assay and inter-person variation were evaluated in-house. One sample, amplified in duplicate, using three different cyclers and processing two runs per cycler, was tested using the *Auto*-LiPA instrument. This protocol was performed in duplicate by two different persons.

No significant variation was observed and an identical typing result was obtained in all cases without the need for editing.

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