
ALARIS[®]
MEDICAL SYSTEMS

IVAC[®] P7000 Mk II

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English Getting Started

Page

Introduction	2
Features of the IVAC® P7000 Syringe Pump	2
Controls and Indicators	3
Front Panel and Main Display	3
Operating Precautions	4 - 5
Installation	6
Pole Clamp	6
Replacing the AC Fuses	6
Functional Earth	6
Loading a Syringe	7
Starting the Pump	8
Power ON/OFF	9
Purge	9
Bolus Infusions	9
Review	10
Pressure Level with Pressure Set fitted	10
Pressure Level without Pressure Set fitted	10
Drugs and Dosing - Drug Name Only	11
Drugs and Dosing - Using Pre-configured Drug Dosing Protocol	11
Drugs and Dosing - Using User-programmed Drug Dosing	11
Volume to be Infused (VTBI)	12
Set VTBI Over Time	12
Set by Doserate / Set by ml/h	13
Selecting the Set by Doserate Option	13
Selecting the Set by ml/h Option	13
Clear Volume	14
Induction	14
Multidose	15
Clear multidose	15
24 H Log	16
Event Log	17
Rate Lock	17
Disable Rate Lock	17
Re-Enable Rate Lock	17
Alarm Procedures	18 - 19

Configured Options

Drug Set-up	20 - 21
General Options	22 - 23
Clock Set	24
Hospital Name	25
Enable Syringes	26
Enable Units	27

Specifications

RS232 / Nursecall Feature	85
Self Test Routine	87 - 90
Configuration Record	91
Symbol Definition and Equipment Classifications	93
Specifications	94 - 99
Routine Maintenance Procedures	100
Disposal	100
Cleaning and Storage	100
Battery Operation	104
Occlusion Pressure Limits for IVAC® 50ml Syringes	106 - 107
Spare Parts	108
Service Equipment	109
Trumpet and Start-up Curves	110 - 113
Technical Description	114
Service Contacts	116
Warranty	117

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
Getting Started

Introduction

The IVAC® P7000 is a fully featured high end variable pressure syringe pump suitable for critical care applications.

The IVAC® P7000 variable pressure pump is one of the IVAC® family of "P" Series syringe pumps. Many aspects of general operation, performance, accuracy and safety checking systems designed into the P7000 unit are common throughout the range.

The IVAC® P7000 functions with a range of standard, single use, disposable Luer-lock syringes of various types and accepts sizes from 5ml to 100ml.

Features of IVAC® P7000 Mk II

- ◆ **Rate range 0.1 to 1200ml/hr.**
- ◆ **In-line pressure monitoring 0 to 750mmHg.**
- ◆ **Optional auto setting of pressure alarm.**
- ◆ **Large graphics format display includes pressure trending.**
- ◆ **24 hour logging of volume and mass.**
- ◆ **Event logging records operation of the P7000.**
- ◆ **Optional auto bolus reduction following occlusion alarm.**
- ◆ **Optional fast start facility reduces infusion start-up delays.**
- ◆ **Communications and nurse call interfaces.**
- ◆ **Multi-position pole clamp.**
- ◆ **Configurable drug protocols for simplified drug dosing.**
- ◆ **Programmable multidose and induction protocols.**

Getting Started

Controls and Indicators

START Button - Press to start the infusion.

STOP Button - Press to stop the infusion. The Amber light will flash to indicate an alarm.

PRESSURE TRANSDUCER - Detects if an infusion line with a pressure disc is fitted. The pressure transducer will measure positive infusion line pressures.

MAIN DISPLAY - For screen information see Front Panel and Main Display.

Press **PURGE/BOLUS** to purge the extension line during set up while the pump is stopped or to bolus at an accelerated rate while an infusion is running.

OPTIONS button - Press the **OPTIONS** button to access drug name and other optional features.

AC POWER - When illuminated the unit is connected to AC power and the internal battery is being charged.

BATTERY - When illuminated the pump is running on the internal battery. When flashing the battery power is low, less than 30 minutes of use remains.


PRESSURE - Use the **PRESSURE** button to display the pumping pressure.

Use "Blank Soft Keys" in conjunction with the prompts shown on the display.

ON/OFF - Press once to switch the pump ON. Press and hold down for 3 seconds to switch the pump OFF.

Use "Double Chevron" for fast increase / decrease and the "Single Chevron" for slower increase / decrease of values shown on the display in set-up and configuration. Are also referred to as the rate keys.

Front Panel and Main Display

Getting Started

Operating Precautions

This IVAC® pump has been calibrated for use with single use disposable syringes. To ensure correct and accurate operation, only use Luer-Lock versions of the syringe make and size specified on the pump or described in the directions for use. Use of non-specified syringes or administration sets may impair the operation of the pump and the accuracy of the infusion.

Uncontrolled flow or syphoning may result if the syringe is located on the pump without its finger grips and plunger correctly located in the slots provided, or if it is removed from the pump before the extension line is properly isolated from the patient. Isolation may include closing a tap in the patient line or activating a flow stop clamp.






When combining several apparatus and/or instruments with administration sets and other tubing, for example via a 3 way tap, the performance of the pump may be impacted and should be monitored closely.


Do not mount the pump in a vertical position with the AC power inlet or the syringe pointing upwards. This could affect electrical safety in the event of a fluid spill over the unit or lead to an infusion of air which may be in the syringe. To protect against the introduction of air the user should regularly monitor the progress of the infusion, syringe, extension line and patient connections and follow the priming procedure specified in the purge procedure in Getting Started.

This is a positive pressure device designed to achieve very accurate fluid administration by automatically compensating for resistance encountered in the infusion system.


The pumping pressure alarm system is not designed to provide protection against, or detection of, infiltration conditions which can occur at low pressures.

Several alarm conditions detected by this pump will stop the infusion and generate audible alarms. Users must perform regular checks to ensure that the infusion is progressing correctly and no alarms are operating.


This instrument is protected against the effects of external interference, including high energy radio frequency emissions, magnetic fields and electrostatic discharge (for example, as generated by electrosurgical and cauterising equipment, large motors, portable radios, cellular telephones etc.) and is designed to fail safe if unreasonable levels of interference are encountered.

In some circumstances the unit may be affected by an electrostatic discharge of above 8kV. At test levels lower than these values the unit will operate normally. In rare circumstances the unit may be affected by radio frequency radiation above a level of 3V/m. If the unit is affected by this external interference the unit will fail safe or reset, (a call back alarm will occur after 2 minutes). Should false alarm conditions be encountered either, remove the source of the interference, or regulate the infusion by another appropriate means.

This unit emits a certain level of electromagnetic radiation which is within the levels specified by IEC60601-2-24 and IEC60601-1-2. If however the unit interacts with other equipment measures should be taken to minimise the effects, for instance by repositioning or relocation.

If this instrument is dropped, subjected to excessive moisture, humidity or high temperature, or otherwise suspected to have been damaged, remove it from service for inspection by a qualified service engineer.

An explosion hazard exists if the instrument is used in the presence of flammable anaesthetics. Exercise care to locate the unit away from any such hazardous sources. An electrical shock hazard exists if the units casing is opened or removed. Refer all servicing to qualified service personnel.

A comprehensive service manual containing circuit descriptions, servicing and testing information is available for this unit. It can be ordered from your ALARIS Medical Systems® authorised distributor (Technical Service Manual Part Number 6000PB00001).

Getting Started

Installation

Check that the pump is complete, undamaged and that the voltage rating specified on the base plate is compatible with your AC power supply. Items supplied with this ALARIS Medical Systems® syringe pump are;

- ◆ **IVAC® P7000 SYRINGE PUMP**
- ◆ **MULTI-POSITION POLE CLAMP**
- ◆ **DIRECTIONS FOR USE**
- ◆ **AC POWER CABLE (AS REQUESTED)**
- ◆ **PROTECTIVE PACKAGING**

Connect the unit to the AC power supply for 24 hours to ensure that the internal battery is fully charged.

Should the pump fail to perform correctly, replace it in its original protective packaging and contact a qualified service engineer for investigation.

Pole Clamp

The pole clamp is supplied fitted to the rear of the unit and will provide secure fixing to standard I.V. poles of a diameter of up to 40mm.

The pole clamp can also be fitted in a choice of 4 fixing positions allowing the unit to be mounted to vertical and horizontal poles, equipment rails and hospital furniture in a variety of convenient operating orientations.

The pole clamp may be adjusted for use with horizontal fittings by using the existing fixings screws with the alternative fixing holes in the pole clamp.

The multi position pole clamp hinge will support the unit at a range of angles on an I. V. pole. To adjust the angle the unit needs to be secured to a pole; using a hand at either end of the instrument, change the screen viewing and syringe access angle of the instrument.

Important:

Do not mount the unit with the AC power inlet or the syringe pointing upwards. This could affect the electrical safety in the event of a fluid spill or lead to the infusion of air which may be in the syringe.

Replacing the AC Fuses

If the pump continually illuminates the battery symbol and the AC power indicator light does not illuminate when the pump is connected to the AC power supply and switched ON, suspect that either, the power supply fuse in the AC power plug, or, the internal AC fuses have blown.

First check the power supply fuse in the AC mains plug. If the AC power indicator light does not illuminate remove the pump from service.

It is recommended that only a qualified service engineer replaces the AC fuses. For further information regarding the replacement of the internal AC fuses refer to the technical service manual.

Functional Earth

It should be noted that the PSU fitted to this pump is functionally earthed. The functional earth feature has been designed to prevent electrical interference to other medical or electrical equipment.

The AC power supply applied to this pump must have an earth conductor.

Important:

This provides a functional earth, not a protective earth.


Getting Started

Loading a Syringe

Place the pump on a stable horizontal surface or secure using the pole clamp fitted. Prepare, load and prime the single use disposable syringe and extension line using standard aseptic techniques.


Important:

Only use a syringe of the type stated on the pump labelling.

1. Squeeze the finger grips on the plunger holder and slide the mechanism to the left. Lift the syringe clamp and rotate to the left.


2. Insert the syringe into the slots on the plunger holder.

3. Squeeze the finger grips on the plunger holder and slide the mechanism to the right until the syringe finger flanges locate in the V slot.


Important:

Advance the syringe until the finger flanges touch the front of the V slot closest to the syringe clamp. This is important to prevent delay at the start of the infusion.

4. Release the finger grips. Apply gentle pressure on the plunger holder to ensure that the drive is engaged. Rotate the syringe clamp forward until it locks onto the syringe barrel.

5. Check that the syringe plunger and finger flanges are correctly located in their slots.

Important:

Only use a syringe of the type and size indicated on the pump or in this manual. Using an incorrect syringe could adversely affect the accuracy of the infusion and the performance of the pump.






When initially loading the syringe, allow for the volume of fluid contained in the extension line and retained in the syringe at the end of infusion as this "dead-space" will not be infused.

Getting Started

Starting the Pump

1. Connect the pump to an AC power supply using the AC power cable. Press the ON/OFF button
2. NEW PATIENT? - Answering **NO** will retain all previous patient history. **YES** will automatically reset the patient history to zero, press the **CONFIRM** softkey to confirm selection. Check time and date is correct.
3. CONFIRM SYRINGE - Check that the syringe type and size being used matches the display. If required, the type of syringe can be changed by pressing the **TYPE** button.
4. Press **CONFIRM** when the correct type and size are shown.
5. SET A RATE - Check the rate shown if old patient data has been retained and change the rate if necessary using the **ADJUST** buttons.
6. PURGE (if required) - Press the **PURGE/BOLUS** button then the two **PURGE** softkeys together until fluid flows and priming of the syringe extension line is complete. The audible alarm will operate during use of the **PURGE** softkeys and the volume used during priming will be shown on the display.
7. CONNECT PATIENT - Connect the extension line to the patient access device, and insert the pressure disc into the pressure transducer.
8. START - Press **START** to commence operation. **INFUSING** will be displayed.

The **AMBER STOP** light will be replaced by the flashing **GREEN START** light to indicate that the pump is operating.


Important:

The unit will automatically operate from its internal battery if the pump is switched on without being connected to the AC power supply.

Each time the unit is switched ON, check that the alarm beeps twice and that all the segments of the display, the green and amber lights are illuminated during the self test routine.

Getting Started

Power ON/OFF


To power up the unit press the power on/off button once.

To power off the unit press and hold down the power ON/OFF button for 3 seconds. The display will count down and switch the unit OFF.

Purge

The **PURGE** button allows the delivery of a limited volume of fluid in order to purge the extension line prior to being connected to a patient or after changing a syringe. The purge function cannot be used if the rate lock is enabled.

1. To use the purge option stop the pump and press the **PURGE** button.

2. The display will change and show two purge softkeys. Press the two **PURGE** softkeys simultaneously.

The purged volume is not added to the volume infused.


3. When the purging function has been completed press the **QUIT** softkey to exit the menu.

Bolus Infusions

The bolus feature can be used during induction as well as during infusion.

If the volume of the bolus reaches the bolus volume limit the bolus will stop and the pump will automatically revert to infuse at the set rate. If the volume to be infused is reached during a bolus, the volume to be infused complete alarm will operate and the unit will revert to its previous state. Press **MUTE** to stop the alarm or **CANCEL** to acknowledge the alarm.

1. During infusion press the **BOLUS** button once to display the bolus screen.

2. Use the chevron keys to set the bolus dose required.

3. To deliver the bolus press the two **BOLUS** softkeys simultaneously. During the bolus the unit will display the volume being infused.

4. When the desired bolus has been delivered, release the two bolus softkeys. The bolus volume will be added to the total volume infused. Press **QUIT** to exit the bolus screen.

Important:


Alarms are not disabled during the operation of the bolus feature. The pressure alarms are temporarily increased to their maximum levels during the operation of the bolus feature.

A bolus cannot be administered if the rate lock is enabled or if a multidose set-up is in use.

Getting Started

REVIEW Review

1. To review the protocol data press the **REVIEW** softkey.
2. The protocol data will be displayed in the main display.
3. Press the **REVIEW** softkey to toggle between the pressure setting and the VTBI status.
4. Repeat until the protocol data required is displayed in the main display.



**REVIEW
DATA**


Pressure level with pressure set fitted (dedicated)

1. To check and adjust the pressure level press the **PRESSURE** button. The display will change and a bar graph will show the pressure alarm level and current pressure level.
2. To adjust the alarm level press the rate keys to increase or decrease the level. The new level will be indicated on the display.
3. When the pressure level has been selected press **QUIT** to exit the screen.
4. If the Auto Pressure Option has been enabled (see General Options) an automatic pressure alarm level can be calculated and set by pressing the **AUTO** softkey.

Pressure level without pressure set fitted (non-dedicated)


1. To check and adjust the pressure level press the **PRESSURE** button. The display will change and a bar graph will show the pressure alarm level and current pressure level.
2. To adjust the alarm level press the rate keys to increase or decrease the level. The new level will be indicated on the display.
3. When the pressure level has been selected press **QUIT** to exit the screen.

Getting Started

These options enable the unit to be set-up for use with a specific drug and/or dosing protocol. Drugs are pre-configured (see configured options) to enable rapid selection of the drug name, dosing units and default rate. For increased security using a configured drug maximum and minimum safety limits are programmable for concentration and dose rates.

Drugs and Dosing - Drug name only ?



1. Press the ? button to access the options menu.
2. Select the **DRUGS AND DOSING** option using the chevron keys
3. Press the **OK** softkey indicated on the screen to confirm the selection.
4. To select a drug name press the **YES** softkey indicated on the display.
5. Press the **NO** softkey to not select dosing or drug protocols.

Caution :

The mass infused display is a sum of all drug masses infused which may be of different concentrations and even drug type. Therefore it should be noted that the relationship between currently displayed volume infused and mass infused may not directly relate to the current concentration.



Drugs and Dosing - Using pre-configured drug dosing protocol ?

When a protocol is used the drug name will be followed by a *. The set by ml/h / set by doserate option is now available.

1. Press the ? button to access the options menu.
2. Select the **DRUGS AND DOSING** option using the chevron keys
3. Press the **OK** softkey indicated on the screen to confirm the selection.
4. Select the drug name from the list displayed using the chevron keys, press the **YES** softkey indicated on the display to confirm the selection.
5. Press the **YES** softkey indicated on the screen to select DOSING.
6. Press the **YES** softkey indicated on the screen to select PROTOCOL. This will select the predefined protocol for the drug name selected.
7. Enter the dosing information prompted on the screen for the drug selected using the **OK** softkey.
8. Press the **CONFIRM** softkey to enter the drug name, dosing information and the protocol selected.

Drugs and Dosing - Using user-programmed drug dosing ?


1. Press the ? button to access the options menu.
2. Select the **DRUGS AND DOSING** option using the chevron keys
3. Press the **OK** softkey indicated on the screen to confirm the selection.
4. Select the drug name from the list displayed using the chevron keys, press the **YES** softkey indicated on the display to confirm the selection.
5. Press the **YES** softkey indicated on the screen to select DOSING.
6. Press the **NO** softkey indicated on the screen to not select PROTOCOL.
7. Enter the required dosing information via the screen prompts, using the chevron keys and the **OK** softkey.
8. Press the **CONFIRM** softkey to enter the drug name, dosing information and the protocol entered.

? Volume To Be Infused (VTBI)

This option allows you to set a specific volume to infuse, and rate at the end of the volume to be infused, from stop, 1ml/h, 2ml/h and continuous infusion at the set rate.

1. Press the “?” softkey to access the options menu.
2. Select the **SET VTBI** option using the chevron keys and press the **OK** softkey indicated on the screen.
3. Enter the volume to be infused using the chevron keys. When the volume to be infused has been entered, press the **OK** softkey.
4. Select the rate at the end of the VTBI using the rate selection chevrons. The default is STOP.
5. Press the **OK** softkey to exit the VTBI menu.
6. When the pump has delivered the set volume it will alarm. Press the **CANCEL** softkey to clear the alarm.
7. Press the **CLEAR** softkey to turn the VTBI function off, or set a new VTBI using the chevron keys.


? Set VTBI Over Time

Set 'VTBI Over Time' must have been activated in the Configured Options.

This option allows a VTBI and delivery time to be specified. The rate necessary to deliver the required volume within the required time will be calculated and displayed. Infusion must be stopped to select this option.

1. Press the “?” button to access the options menu.
2. Select the **SET VTBI OVER TIME** option using the chevron keys and press the **OK** softkey indicated on the screen.
3. Adjust the volume to be infused using the chevron keys, when the volume to be infused has been entered press the **OK** softkey.
4. Enter the time over which the volume is to be infused. The infusion rate will automatically be calculated. Press the **OK** softkey to enter the value.
5. Select the rate at VTBI end and press the **OK** softkey. The default is STOP.

Getting Started


Set by Doserate / Set by ml/h ?

To set rates accurately in doserate or flowrate increments it may be necessary to switch between the rate adjust options SET BY DOSERATE and SET BY ml/h. An arrow to the left of the rate display indicates the rate that will change when the chevron buttons are used to increase/decrease the infusion rate.

To set a doserate precisely the arrow must be pointing to the doserate (mg/kg/h); the flowrate will be calculated from the doserate.

To accurately set a flowrate the arrow must be pointing to flowrate (ml/h); the doserate will be calculated from the flowrate.


Selecting the Set By Doserate Option

1. Press the ? button to access the options menu.
2. Select the **SET BY DOSERATE** option using the chevron keys.
3. Press the **OK** softkey. This will select the set by doserate option and an arrow will automatically point to the doserate on the display. If necessary the doserate can be increased or decreased using the chevron keys.

INFUSING	DRUG NAME	
56.0 ml/h	VOLUME	WEIGHT
8.00 mg/kg/h	19.0ml	70.0kg
190mg		
+ ADJUST -		REVIEW

Selecting the Set By ml/h Option


1. Press the ? button to access the options menu.
2. Select the **SET BY ml/h** option using the chevron keys.
3. Press the **OK** softkey. This will select the set by flowrate option and an arrow will automatically point to the flowrate on the display. If necessary the flowrate can be increased or decreased using the chevron keys.

INFUSING	DRUG NAME	
56.0 ml/h	VOLUME	WEIGHT
8.00 mg/kg/h	19.0ml	70.0kg
190mg		
+ ADJUST -		REVIEW

Getting Started

? Clear Volume


1. Press the ? button to access the options menu.
2. Select the **CLEAR VOLUME** option using the chevron keys and press the **OK** softkey indicated on the screen.
3. The volume infused will be displayed.
4. Press the **YES** softkey to clear the volume infused displayed. Press the **NO** softkey to retain the volume.


? Induction

The induction option is enabled within the general options. See General Options in this user manual.


1. Press the ? button to access the options menu.
2. Select the **INDUCTION** option using the chevron keys
3. Press the **OK** softkey indicated on the screen to confirm the selection.
4. Adjust the induction **RATE** using the chevron keys, then press the **OK** softkey.
5. Select the induction **VOLUME** using the chevron keys, then press the **OK** softkey. If an incorrect value has been entered, press the **BACK** softkey to return to the previous stage.
6. Enter the **MAINTENANCE** infusion rate using the chevron keys, then press the **OK** softkey.
7. Review the induction data on screen, then press **CONFIRM**. If necessary use the **BACK** softkey to return to the data. Once confirmed, the pump will return to the main display.
8. The induction can be cleared. Stop the infusion then press the ? button. Select **CLEAR INDUCTION** with the chevron keys and press the **OK** softkey to confirm.


The multidose function is enabled within the general options. See General Options in this user manual.

1. Press the **?** button to access the options menu.
2. Select the **MULTIDOSE** option using the chevron keys
3. Press the **OK** softkey indicated on the screen to confirm the selection.
4. Adjust the multidose **RATE** using the chevron keys, then press the **OK** softkey.
5. Select the multidose **VOLUME** to be infused using the chevron keys, then press the **OK** softkey. If an incorrect value has been entered, use the **BACK** softkey to return to the previous stage.
6. Enter the **MAINTENANCE** infusion rate using the chevron keys, then press the **OK** softkey.
7. Select the **FREQUENCY** (hr:mins) of each multidose using the chevron keys, then press the **OK** softkey.
8. Review the multidose data on screen, then press **CONFIRM**. If necessary use the **BACK** softkey to return to the data. Once confirmed, the pump will return to the main display. The review section of the display will show either multidose volume remaining or time between multidoses remaining.

Clear Multidose ?


1. Stop the infusion then press the **?** button.
2. Select **CLEAR MULTIDOSE** with the chevron keys and press the **OK** softkey to confirm.

Getting Started


? 24H Log

This option allows the 24H log of volume infused to be reviewed.


1. Press the "?" button to access the options menu.
2. Select the **24H LOG** option using the chevron keys and press the **OK** softkey indicated on the screen.
3. Press the **NEXT** softkey indicated on the screen to access the hourly volume infused log.
4. Press the **QUIT** softkey indicated on the screen to exit the log.

2a 2b

3


4

? Event Log

1. Press the ? button to access the options menu.
2. Select the **EVENT LOG** option using the chevron keys.
3. Press the **OK** softkey indicated on the screen.
4. Press the **QUIT** softkey indicated on the screen to exit the menu.

2 3


4

The rate lock feature can be enabled in the general options of the pump.


When the infusion rate has been set and the infusion started, the rate lock prompt will appear in the main menu.

When the rate lock is enabled the rate cannot be altered if the pump is stopped, the infusion can be started and stopped in the normal way. The rate lock disables the bolus and purge operations.


1. To select the rate lock function press the **YES** softkey. Press the **NO** softkey if the rate lock is not required.

Disable Rate Lock ?

1. Press the "?" button to access the options menu.
2. Select the **UNLOCK RATE** option using the chevron keys and press the **OK** softkey indicated on the screen.

Re-Enable Rate Lock ?


1. Press the "?" button to access the options menu.
2. Select the **RATE LOCK** option using the chevron keys and press the **OK** softkey indicated on the screen.

Getting Started

Alarm Procedures

Alarms are indicated by a combination of an audible alarm, flashing amber **STOP** light and a descriptive message in the display.

1. A continuous audible alarm indicates that the infusion has stopped. First press **MUTE** to silence the alarm for a maximum of 2 minutes. Then check the display for a alarm message. Press **CANCEL** to cancel the alarm message.
2. When the cause of the alarm has been rectified, press the **START** button to resume the infusion.

Description

Display

PUMP DRIVE DISENGAGED - The drive system has been disengaged during the infusion. Check the finger grip and the position of the syringe.

DRIVE DECLUTCHED

EXCESSIVE DRIVE PRESSURE - Pumping pressure has reached the alarm limit. Following an occlusion the drive will backoff (if enabled - see general options) to reduce excess pressure; alternatively squeeze finger grips on the plunger holder to release the drive mechanism and relieve any excessive pressure in the syringe and patient line. Identify and remove the cause of the blockage in the administration system before restarting the infusion.

OCCLUSION

SYRINGE CLAMP LOCATION - Incorrect size of syringe has been fitted, the syringe clamp has not been positioned correctly, or the syringe has been disturbed during operation. Check syringe location, position of syringe clamp, and syringe.

CHECK SYRINGE LOCATION

PLUNGER NOT DETECTED -The syringe plunger is not correctly fitted in the plunger holder. Check the syringe plunger location.

CHECK PLUNGER LOCATION

BATTERY CHARGE LOW WARNING - Battery charge low with up to 30 minutes operation remaining. Battery indicator will flash and after 30 minutes a continuous audible alarm will indicate that the battery is exhausted. Reconnect to AC power to continue operation and charge the internal battery.

BATTERY LOW

Display	Description
BATTERY EMPTY	BATTERY EMPTY - Internal battery empty. To silence the alarm switch the pump OFF and reconnect the pump to AC power. Restart operation on AC power to charge the internal battery.
NEAR END OF SYRINGE	NEAR END OF SYRINGE WARNING (NEOI) - This is a timed event and can be set from 1 to 15 minutes. To set / change this value see the general options section in this user manual.
SYRINGE EMPTY	SYRINGE EMPTY - END OF INFUSION (EOI) -The pump has reached the end of the infusion. About 1% (can be set in configuration) of the syringe volume will remain in the syringe to minimise the risk of the infusion of air bubbles into the set.
AC POWER FAIL	AC POWER SUPPLY DISCONNECTED WARNING AC Power has been disconnected and the pump is operating on battery power. Reconnect AC power or press START to silence the alarm and continue battery operation. The display will flash AC POWER FAIL. The alarm will automatically cancel if the AC power is reconnected.
ERROR	INTERNAL MALFUNCTION - The alarm system has detected an internal malfunction. Note the malfunction code. Remove the unit from service for examination by a qualified service engineer.
BATTERY WAS LOW AT LAST POWER OFF. HAS PUMP BEEN RECHARGED	BATTERY MAY BE LOW/EMPTY - When the pump was last turned off the battery was low. If the pump has not been charged since then operate on AC power only.




Alarm Procedures - Alarms without screen prompts

	Description
3 BEEPS	CALL BACK WARNING - Unit left switched ON for over 2 minutes without starting operation. Press STOP or any of the control buttons to silence the alarm for a further 2 minutes. Alternatively, press and hold down STOP for 2 seconds to delay call back for 15 minutes.

Configured Options

Drug Set-up

1. Press the ON/OFF power button OFF.
2. Hold down the **START** button and press the power button ON.
3. The main display will show "000". Enter the Access Code "251" using the chevron keys.
4. Enter the first digit and press **NEXT** indicated on the screen.
5. When the complete code has been entered press **OK** softkey indicated on the screen.
6. Select the **DRUG SETUP** option using the chevron keys and press the **OK** softkey indicated on the screen.
7. Select the required drug with the chevron keys and press the **OK** softkey.
8. To use a drug it must be enabled. Press the **YES** softkey indicated on the screen to enable the selected drug.
9. To add or change a drug name use the chevron keys to toggle through the alphabet. To select a letter press the **NEXT** softkey. On completion press the **OK** softkey indicated on the screen.
10. Select the Dose Units required for maintenance doses. Use the chevron keys to select the dose units required. When the dose units have been selected press the **OK** softkey.
11. Select the Default Doserate using the chevron keys to increase / decrease the doserate displayed. When the required doserate has been selected press the **OK** softkey.
12. Select the Minimum Doserate using the chevron keys to increase / decrease the doserate displayed. When the minimum doserate has been selected press the **OK** softkey.
13. Select the Maximum Doserate using the chevron keys to increase / decrease the doserate displayed. When the maximum doserate has been selected press the **OK** softkey.

ng/ml ug/ml mg/ml	DRUG 01 ← CONC UNITS
↑ SELECT ↓ OK	BACK

40.0 ug/ml	DRUG 01 ← DEFAULT CONC
+ ADJUST - OK	BACK

40.0 ug/ml	DRUG 01 ← MIN CONC
+ ADJUST - OK	BACK

100 ug/ml	DRUG 01 ← MAX CONC
+ ADJUST - OK	BACK

900 ml/h	DRUG 01 ← BOLUS RATE
+ ADJUST - OK	BACK

30 ml	DRUG 01 ← MAXIMUM BOLUS
+ ADJUST - OK	BACK

650 mmHg	DRUG 01 ← PRESSURE ALARM
+ ADJUST - OK	BACK

DOSE 10.0ug/ml (8.0 . . 20.0) CONC. 40.00ug/ml (40.0 . . 100) BOLUS 600 ml/hr (50ml max) PRESSURE ALARM ▲ 650mmHg	DRUG 01 ? CONFIRM
OK	BACK

14. Select the Concentration Units using the chevron keys, press the **OK** softkey to select the concentration units.

15. Enter the Default Concentration of the drug selected, using the chevron keys to increase / decrease the default concentration shown on the screen. When the default concentration has been selected press the **OK** softkey.

16. Select the Minimum Concentration of the drug selected, using the chevron keys to increase / decrease the minimum concentration shown on the screen. When the minimum concentration has been selected press the **OK** softkey.

17. Select the Maximum Concentration of the drug selected, using the chevron keys to increase / decrease the maximum concentration shown on the screen. When the maximum concentration has been selected press the **OK** softkey.

18. Select the Bolus Rate required, using the chevron keys to increase / decrease the bolus rate. When the rate has been selected press the **OK** softkey.

19. Select the Maximum Bolus using the chevron keys to increase / decrease the volume. When the volume has been selected press the **OK** softkey.




20. Select the Pressure Alarm level using the chevron keys to increase / decrease the level. When the level has been selected press the **OK** softkey.

21. Review the drug setup data, press the **OK** softkey to confirm the data displayed. This will return the pump to the drug setup menu. Use the **BACK** softkey at any time to return to the previous stage.

Configured Options

General Options

1. Press the ON/OFF power button to the OFF position.
2. Hold down the **START** button and press the power button ON.
3. The main display will show "000". Enter the Access Code "251" using the rate selection arrows.
4. Enter the first digit and press **NEXT** indicated on the screen.
5. When the complete code has been entered press the **OK** softkey indicated on the screen.
6. Select **GENERAL OPTIONS** using the chevron keys and press the **OK** softkey indicated on the screen.
7. Select the option you wish to enable/disable or adjust and press the **MODIFY** softkey indicated on the screen.
8. When all the options that you require to enable/disable or adjust have been selected press the **OK** softkey indicated on the screen.

The general options available are as follows:


Fast Start and Backoff -	Enables drive slack fast start and occlusion auto bolus reduction.
VTBI -	Enables or disables the VTBI function.
Set VTBI over Time -	Sets a specific volume to be infused over a fixed period of time.
Induction -	Enables or disables the induction option.
Multidose -	Enables or disables the multidose option.
Rate Lock -	If enabled, the rate can be locked to prevent tampering of the set infusion rate.
Quiet Mode -	If enabled, the low priority alarms and button pushes are muted.
Pressure Alarm -	Sets the pressure alarm value. Default is 300mmHg.
Max. Pressure Alarm -	Sets the maximum pressure limit. Default is 750mmHg.
Auto Pressure Alarm -	Enables or disables the automatic pressure alarm option.
AC Fail Alarm -	Enables or disables the AC Fail Alarm.
Rate Titration -	If enabled, the rate can be changed whilst the pump is infusing. To titrate a rate while infusing, use the rate adjust keys to select the new rate and press the START button to confirm the new rate.
NEOI Warning -	Sets the NEOI warning from 1 to 15mins.
EOI Point -	Sets the end of infusion volume as a percentage of the syringe volume from 0.5% to 5.0%.
KVO at EOI -	If enabled the pump will switch to run at the KVO rate - 1.0ml/h (or the current set rate if this is lower) at the EOI point. If disabled the pump will stop at the EOI point.
Beam Alarm -	Sets up the default beam level between 0 and 7.
Weight -	Sets up the patient weight between 0.1kg and 150kg. Default is 70.0kg.
Purge Rate -	Sets the purge rate between 100 and 500ml/h. The upper limit prevents excessive loss of fluid whilst purging.
Bolus -	Enables bolus delivery.
Default Bolus -	Sets the default bolus rate between 10 and 1200ml/h.
Cap Bolus Rate -	Sets maximum permissible value for the bolus rate.
Cap Rate -	Sets maximum permissible value for the infusion rate.
Display Event Log -	If enabled the event log can be accessed via the options menu.
Log Last Patient Only -	If enabled the event log shows only the last patient (since new patient selected). If disabled the complete log can be viewed.
Nurse Call -	Enables or disables the nurse call feature.
Nurse Call Inverted -	When enabled the nurse call feature is inverted.
Comms Pump Address -	Sets the Pump Communications Address.
Comms Monitor Only -	When disabled, allows remote control of the unit from the comms link. Monitor is always possible.
Comms Odd Parity -	If enabled, communications parity is odd.
Comms ASCII -	If enabled, communications protocol is in ASCII.

Configured Options


Clock Set

Clock Set is used to set the pumps internal clock.


- 1. Press the ON/OFF power button OFF.
- 2. Hold down the **START** button and press the power button ON.
- 3. The main display will show “000”. Enter the Access Code “251” using the rate selection arrows.
- 4. Enter the first digit and press **NEXT** indicated on the screen.


- 5. When the complete code has been entered press the **OK** softkey indicated on the screen.
- 6. Select **CLOCK SET** using the chevrons and press the **OK** softkey indicated on the screen.
- 7. To adjust the value displayed use the chevrons.


- 8. When the displayed value is correct press the **NEXT** softkey to access the next field.
- 9. Repeat the process until the correct time and date is displayed. Press the **OK** softkey to store the values.


The hospital name allows the user to programme in the name of the hospital, ward or department to appear during the power up display sequence.

1. Press the ON/OFF power button OFF.
2. Hold down the **START** button and press the power button ON.
3. The main display will show "000". Enter the Access Code "**251**" using the rate selection arrows.
4. Enter the first digit and press **NEXT** indicated on the screen.

5. When the complete code has been entered press the **OK** softkey indicated on the screen.
6. Select **HOSPITAL NAME** using the chevrons and press the **OK** softkey indicated on the screen.
7. To adjust the character displayed use the chevrons.


8. When the displayed character is correct press the **NEXT** softkey to access the next field.
9. Repeat the process until the correct Hospital Name is displayed. Press the **OK** softkey to store the name.

Configured Options


Enable Syringes

Enable syringes is used to pre-configure the type and size of syringe allowed for use by an operator. Select all the syringes which may be used with the pump. Disable any syringes which will not be used.

1. Press the ON/OFF power button OFF.
2. Hold down the **START** button and press the power button ON.
3. The main display will show "000". Enter the Access Code "251" using the rate selection arrows.


4. Enter the first digit and press **NEXT** indicated on the screen.
5. When the complete code has been entered press the **OK** softkey indicated on the screen.


7. Use the chevron to toggle up and down the list of available syringe types and sizes. Press the **MODIFY** softkey to enable/disable the syringe.

8. Press the **OK** softkey to store the syringes enabled.


9. When the set-up is complete, turn the pump OFF and return the pump into service, or select the next configuration option.


Enable units is used to pre-configure the dose units allowed for use by the operator when selecting units in the Drugs and Dosing option.

1. Press the ON/OFF power button OFF.
2. Hold down the **START** button and press the power button ON.
3. The main display will show "000". Enter the Access Code "251" using the rate selection arrows.

4. Enter the first digit and press **NEXT** indicated on the screen.
5. When the complete code has been entered press the **OK** softkey indicated on the screen.
6. Select **ENABLE UNITS** using the chevron keys and press the **OK** softkey indicated on the screen.

7. Use the chevrons to toggle up and down the list of available units. Press the **MODIFY** softkey to enable/disable the units.
8. Press the **OK** softkey to store the units.
9. When the set-up is complete, turn the pump OFF and return the pump into service, or select the next configuration option.

Notes

[illegible]

Pagina

English

Getting Started	2 - 19
Configured Options	20 - 28
Specifications	85 - 116

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


The RS232 / Nursecall feature fitted to this IVAC® "P" Series Syringe Pump allows the pump to be monitored remotely and/or controlled via a suitable central monitoring or computer system.

When the pump is started by a command from the serial interface, communication must take place over the serial interface, a communication must take place every 15 seconds or the pump will alarm, display communications failure and stop infusing. This failure protects against failure of the communications, including the removal of the RS232 cable.

Important:

The nursecall interface provides a remote backup to the internal audible alarm. It should not be relied upon to replace monitoring of internal alarm.

Refer to the technical service manual for further information regarding the RS232 interface. Since it is possible to control the syringe pump using the RS232 interface at some distance from the pump and hence remote from the patient, responsibility for the control of the pump is vested in the software run on the computer control system. The assessment for the suitability of any software used in the clinical environment to control or receive data from the pump lies with the user of the equipment. This software should include detection of the disconnection or other failure of the RS232 cable. The protocol is detailed in the technical service manual and is for general information only. This relates to IVAC® "P" Series Syringe pumps fitted with the RS232 communication interface.

Any connected analogue and digital components are required to meet EN60950 for data processing and EN60601 for medical devices. Anyone connecting additional devices to the signal input or output is a system configurator and responsible for meeting the requirements of the system standard EN60601-1-1.

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
The self test routine is designed to allow confirmation of many of the pump functions, defaults and calibrations without requiring internal inspection. It does not represent a full calibration check.


1. Press the ON/OFF power button OFF.
2. Hold down the **START** button and press the power button ON.
3. The main display will show "000". Enter the Access Code "123" using the rate selection arrows.
4. Enter the first digit and press **NEXT** indicated on the screen.
5. When the complete code has been entered press **OK** indicated on the screen.
6. The pump will now proceed through a series of tests. Press the **NEXT** button to move to the next test.

Important:

Should the pump fail the test sequence, it should be taken out of service and inspected by a qualified service engineer.

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
Spezifikationen, Dati tecnici, Specifications

Configuration Record

Use the following sheet to record the configuration settings.

1	FAST START & BACK OFF	Drive slack fast start and occlusion auto bolus reduction enabled.
2	VOLUME OVER TIME	Volume to be delivered over a fixed time enabled
3	INDUCTION	Induction option enabled
4	MULTIDOSE	Multidose option enabled
5	RATE LOCK	Rate lock feature enabled
6	QUIET MODE	Quiet mode enabled
7	PRESSURE ALARM	Default pressure alarm limit
8	MAX PRESSURE	Maximum pressure alarm limit
9	AUTO PRESSURE	Automatic pressure option enabled
10	AC FAIL ALARM	AC Fail alarm enabled
11	RATE TITRATION	Rate titration feature enabled
12	NEOI WARNING	Near end of infusion warning time (minutes)
13	EOI POINT	End of infusion: percentage of volume in syringe at end of infusion
14	KVO AT EOI	KVO enabled at end of infusion
15	BEAM ALARM	Default occlusion alarm level
16	WEIGHT	Default patient weight value
17	PURGE RATE	Purge rate value
18	BOLUS	Bolus Enabled
19	DEFAULT BOLUS	Default bolus value
20	CAP BOLUS RATE	Maximum cap on bolus rate
21	CAP RATE	Maximum cap on infusion rate
22	DISPLAY EVENT LOG	Display event log enabled
23	LOG LAST PATIENT	Event log last patient or whole log
24	NURSECALL	Nurse call enabled
25	NURSECALL INVERTED	Nurse call inverted
26	COMMS PUMP ADDRESS	Pump communications address
27	COMMS MONITOR ONLY	Comms. monitor only or monitor and control
28	COMMS ODD PARITY	Comms parity odd enabled
29	COMMS ASCII	Comms ASCII protocol enabled

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
Specifications, Dati tecnici, Spezifikationen

Syringe Types / Tipi di siringa / Sprizentypen -

The unit is calibrated and labelled for use with single use disposable luer-lock syringes. Only use the size and type of syringe specified on the pump display .

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	5ml	10ml	20ml	30ml	50ml	100ml
IVAC®					✓	✓
BD Plastipak	✓	✓	✓	✓	✓	
Terumo		✓	✓	✓	✓	
B Braun Omnifix	✓	✓	✓	✓	✓	
Sherwood Monoject	✓	✓	✓	✓	✓	
Nipro	✓		✓	✓	✓	
Once					✓	
Fresenius					✓	
Zeneca			✓		✓	
B Braun Perfusor			✓		✓	
JMS						✓
BD Perfusion					✓	
BD Precise			✓		✓	

Compatible Extension Lines and Syringes / Tubi di prolunga e siringhe compatibili / Kompatible Infusionsleitungen -

The unit uses a standard, single use, disposable extension line and syringes with luer-lock connectors. Of type designed for use on syringe pumps.

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30602N IVAC® 50/60ml Luer-Lock Syringe / Siringa / Spritze

30120 IVAC® 100/120ml Luer-Lock Syringe / Siringa / Spritze

Syringe extension set with occlusion sensing disc

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Opaque syringe extension set with occlusion sensing disc

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Opaque low sorbing extension set with occlusion sensing disc

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It is recommended that the extension sets are changed according to the hospital protocols.

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Spezifikationen, Dati tecnici, Specifications

Infusion Rate / Portata infusione / Infusionsrate -

Maximum infusion rate can be set as part of the configuration.

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Maximum Bolus rates can be set as part of the configuration. Bolus rates are user adjustable.

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Continuous Mode Operation, Transportable

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Critical Volume / Volume critico / Kritisches Volumen -

The maximum-over infusion which can occur in the event of a single fault condition is 0.5ml

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Important:

(System accuracy is +/-2% typical by volume as measured using the trumpet curve test method defined in IEC601-2-24 (DRAFT) at rates of 1.0ml/h and above when the instrument is used with the recommended syringes. Differences in factors such as size and plunger force in compatible syringes can cause variations in accuracy and trumpet curves.) also see trumpet curves section.

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Rechargeable Sealed Lead Acid. Automatically charges when the pump is connected to AC power.)

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10 hours from discharge to 80% charge and 24 hours to 100% charge

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The electronic memory of the unit will be retained for more than 6 months when not powered up.

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Dimensions / Dimensioni / Abmessungen -

400mm (w) x 115mm (h) x 180mm (d). Weight: 3.5kg (excluding pole clamp and power cable).

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30% to 90%

700mbar to 1060mbar

Transport/Storage / Trasporto/Stockage /

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This instrument is designed and manufactured in the U.K. by ALARIS Medical UK Ltd under patent GB222444. ALARIS Medical UK Ltd reserves the right to alter product specifications without notice.

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


- Description
- 1 Nursecall Normally Closed
 - 2 Received data (RXD) Input
 - 3 Transmit data (TXD) Output
 - 4 Not used
 - 5 Ground (GND)
 - 6 Nurse Call Normally Open
 - 7 Request to send (RTS) Input
 - 8 Clear to send (CTS) Output
 - 9 Nursecall Common

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
Important:

The IBM connector pins grouped in pairs above should be linked at the connector.

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Routine Maintenance Procedures

To ensure that this instrument remains in good operating condition, it is important to keep it clean and carry out the routine maintenance procedures described below. All servicing should only be performed by a qualified service engineer with reference to the technical service manual (TSM) for this product (Service Manual reference: 6000PB00001).

Technical service features are available using access code "376", refer to the technical service manual.

Important:

If the pump is dropped, damaged, subjected to excessive moisture or high temperature, immediately take it out of service for examination by a qualified service engineer.

INTERVAL

As Required

12 Monthly

ROUTINE MAINTENANCE PROCEDURE

Thoroughly clean external surfaces of the pump before and after prolonged period of storage.

1. Inspect AC power supply plug and cable for damage.
2. Perform functional OPERATIONAL AND ALARM TEST - see TSM
3. Perform the SELF TEST ROUTINE.
4. Perform rate accuracy verification test - see TSM.
5. Perform pressure calibration verification test - see TSM.
6. Operate the pump on battery power until the battery low alarm then charge the battery to confirm battery operation and charging.

Disposal

The pump should be disposed of taking environmental factors into consideration. To ensure no risk or hazard remove the internal rechargeable battery and the nicad battery (component B1) from the control board and dispose of as outlined by the local country regulations. Do not send back to manufacturer. All other components can be safely disposed of in the normal manner.

Cleaning and Storage

Before the transfer of the pump to a new patient and periodically during the use, clean the pump by wiping over with a lint-free cloth lightly dampened with warm water and a standard disinfectant / detergent solution.

The syringe and extension line are disposable single use items and should be discarded after use according to their manufacturers instructions.

If the pump is to be stored for an extended period it should be first cleaned and the internal battery fully charged. Store in a clean, dry atmosphere at room temperature and, if available, employ the original packaging for protection.

Once every 3 months during storage, carry out functional tests as described in the OPERATIONAL AND ALARM TEST in the technical service manual and SELF TEST ROUTINE and ensure that the internal battery is fully charged.

Important:

Before cleaning always switch OFF and disconnect from the AC power supply. Never allow liquid to enter the casing and avoid excess fluid build up on the pump. Do not use aggressive cleaning agents as these may damage the exterior surface of the pump. Do not steam autoclave, ethylene oxide sterilise or immerse this pump in any fluid.

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Battery Operation

The internal rechargeable battery allows continued operation when the AC power is unavailable, for example during patient transfer or AC power failure. A fully charged battery will provide over 6 hours operation at typical infusion rates. From the battery low alarm it will take about 24 hours to fully recharge when reconnected to the AC power supply, whether the unit is in use or not. The battery is automatically charged during AC operation and whenever the unit is connected to the AC power supply and the AC power indicator is illuminated.

It is good practice to periodically operate the unit on battery power until the battery low alarm then charge the battery to confirm battery operation and charging. When not in use, connect the unit to the AC power supply in order to maintain the battery in the fully charged state.

The battery is a maintenance free, sealed lead acid type and requires no routine servicing. However, to achieve optimum operation, ensure that the battery is fully recharged after full discharge, before storage, and at regular 3 month intervals during storage.

The internal rechargeable battery will retain charge if utilised as described above on a regular basis. Charge retention will eventually degrade. Where retention is critical the internal battery should be replaced every 3 years.

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Specifications, Dati tecnici, Spezifikationen

Occlusion Pressure Limits for IVAC® 50ml Syringes

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The following tables show the worst case values for line pressure, time to alarm and bolus volume that can be expected in the event of an occlusion when the IVAC® 50ml syringe is selected, G40020 administration set.

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* Tests at these levels may alarm immediately - the force at these levels is commonly less than the friction in the syringe (with no additional fluid pressure). The result is that the pressure relating to the low forces will be less than the nominal quoted occlusion pressure.

Bolus volume following occlusion will be minimised by the back off feature if enabled. The back off will reduce the line pressure by removing the volume stored in the occluded line and deducting this volume from volume infused.

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A comprehensive list of spare parts for IVAC® Syringe pumps is included within the service manual. This can be ordered from ALARIS Medical Systems®, or authorised distributor. For part number please refer to summary parts list below:

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Spezifikationen, Dati tecnici, Specifications

Service Equipment, Attrezzi per l'assistenza, Wartungsgeräte

These IVAC® "P" Series Syringe pumps have been designed to allow simple and low cost servicing. Standard components are employed where possible so that no special test, calibration or tools are required. However, the following items may be useful for general servicing.

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In this instrument, as with all infusion systems, the action of the pumping mechanism and variations in individual syringes cause short-term fluctuations in rate accuracy.

The following curves show typical performance of the system in two ways: 1) the accuracy of fluid delivery over various time periods is measured (trumpet curves), and 2) the delay in onset of fluid flow when infusion commences (start-up curves).

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Trumpet Curves, Curve a tromba, Trompetenkurven

Trumpet curves are named for their characteristic shape. They display discrete data averaged over particular time periods or 'observation windows', not continuous data versus operating time. Over long observation windows, short term fluctuations have little effect on accuracy as represented by the flat part of the curve. As the observation window is reduced, short term fluctuations have greater effects as represented by the "mouth" of the trumpet.

Knowledge of system accuracy over various observation windows may be of interest when certain drugs are being administered. Short term fluctuations in rate accuracy may have clinical impact depending on the half-life of the drug being infused and the degree of inter vascular integration, the clinical effect cannot be determined from the trumpet curves alone.

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The start-up curves represent continuous flow versus operating time for two hours from the start of the infusion. They exhibit the delay in onset of delivery due to mechanical compliance and provide a visual representation of uniformity. Trumpet curves are derived from the second hour of this data. Tests performed per IEC601-2-24 (draft) standard.

Important:

Start-up and trumpet curves may not be indicative of operation under negative pressure. Differences in factors such as size and plunger force in compatible syringes produced by other manufacturers can cause variations in accuracy and trumpet curves as compared to those represented. Additional curves for compatible syringes are available upon written request. For applications where flow uniformity is a concern, rates of 1.0ml/hr or above are recommended.

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Specifications, Dati tecnici, Spezifikationen

In this instrument, as with all infusion systems, the action of the pumping mechanism and variations in individual syringes cause short-term fluctuations in rate accuracy.

The following curves show typical performance of the system in two ways: 1) the delay in onset of fluid flow when infusion commences (start-up curves), and 2) the accuracy of fluid delivery over various time periods (trumpet curves).

The start-up curves represent continuous flow versus operating time from the start of the infusion. They exhibit the delay in onset of delivery due to mechanical compliance and provide a visual representation of uniformity. Trumpet curves are derived from the second hour of this data. Tests performed per IEC60601-2-24 standard.

Trumpet curves are named for their characteristic shape. They display discrete data averaged over particular time periods or 'observation windows', not continuous data versus operating time. Over long observation windows, short term fluctuations have little effect on accuracy as represented by the flat part of the curve. As the observation window is reduced, short term fluctuations have greater effects as represented by the "mouth" of the trumpet.

Knowledge of system accuracy over various observation windows may be of interest when certain drugs are being administered. Short term fluctuations in rate accuracy may have clinical impact depending on the half-life of the drug being infused, therefore the clinical effect cannot be determined from the trumpet curves alone.

Important: Start-up and trumpet curves may not be indicative of operation under negative pressure.

Differences in factors such as size and plunger force in compatible syringes produced by other manufacturers can cause variations in accuracy and trumpet curves as compared to those represented. Additional curves for compatible syringes are available upon written request.

For applications where flow uniformity is a concern, rates of 1.0ml/hr or above are recommended.

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

Spezifikationen, Dati tecnici, Specifications

Technical Description, Descrizione tecnica, Technische Beschreibung

The following details outline the safety checks designed into IVAC® Syringe Pumps to minimise the possibility of under/over infusions.

MONITORING OF THE SYSTEM CLOCK FREQUENCY / WATCHDOG

The system clock, which is used to control the rate of the pump is derived from the microprocessor crystal oscillator. A watchdog circuit is implemented in the pump to monitor the correct time period. The watchdog circuit requires the microprocessor to send a reset pulse every 10mS to stop the counter from timing out and triggering the watchdog alarm. The reset signal has to be in a time 'window' of between 8 to 12mS. If it is either too fast or too slow the watchdog hardware will detect this and generate a alarm and disable the motor drive. Additionally, on power up the watchdog is allowed to time-out and the period is measured and tested to be within a set tolerance. This then confirms that both the microprocessor crystal and the watchdog crystal frequencies are correct.

DETECTION OF LINEAR MOVEMENT

A linear potentiometer is incorporated in the unit to detect the movement of the pumping mechanism. This movement is monitored in the electronics and software of the system. If it is detected that the mechanism is either moving too fast, too slow, or not moving at all then an error code will be displayed, the motor drive will be disabled and the pump will stop.

CONTROL OF LINEAR SPEED / INFUSION RATE

The pump mechanism is driven using a d.c. motor; feedback for the control system is provided by two opto switches. To enable the motor drive there are three transistors which need to be turned on. The correct operation of these transistors is tested on power up. The control system monitors the feedback from the opto encoders and adjusts the motor on time to maintain the required speed. If there are no encoder signals feedback, indicating that an opto has failed or the transmission has jammed, an error code will be displayed and the pump will fail-safe. If there are too many encoders detected by the encoder feedback indication that a transistor has gone short circuit an error code will be displayed and the pump will fail-safe. The pump calculates the appropriate motor control rate from the set infusion rate and syringe constant data stored within the pump software. The syringe constant data converts ml's to mm's of movement for each syringe type and size.

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Spezifikationen, Dati tecnici, Specifications

Service Contacts, Assistenza, Service Kontaktadressen

For service contact your local ALARIS Medical Systems® Affiliate Office or Distributor.

ALARIS Medical Systems® Service Centre Addresses:

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AU

ALARIS Medical Australia Pty Ltd,
3/167 Prospect Highway,
Seven Hills, NSW 2147,
Australia.
Tel: (61) 2 9838 0255
Fax: (61) 2 9674 4444

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NZ

ALARIS Medical NZ Ltd.,
Unit 14, 13 Highbrook Drive,
East Tamaki, Auckland,
New Zealand.
Tel: (64) 9 273 3901
Fax: (64) 9 273 3098

SE

ALARIS Medical Nordic, AB
Box 452,
191 24 Sollentuna
Sverige.
Tel: (46) 8 544 43 200
Fax: (46) 8 544 43 225

SG

ALARIS Medical Systems Office,
65 Chulia Street,
#40-04 OCBC Centre,
Singapore 049513.
Tel: (65) 5345351
Fax: (65) 5345516

US

ALARIS Medical Systems, Inc.
10221 Wateridge Circle
San Diego, CA 92121, USA.
Tel: (1) 800 854 7128
Fax: (1) 858 458 6179

ZA

ALARIS Medical S.A. (Pty) Ltd.,
Unit 2 Oude Molen Business Park,
Oude Molen Road, Ndbeni, Cape
Town 7405, South Africa.
Tel: (27) 0860 597 572
Fax: (27) 21 5107567

B D is a trademark of /è un marchio registrato di Becton-Dickinson, Inc.

MONOJECT is a trademark of /è un marchio registrato di Tyco/Healthcare Kendal-Monoject

INJECTOMAT is a trademark of /è un marchio registrato di Fresenius A.G.

PERFUSOR is a trademark of /è un marchio registrato di B Braun Melsungen A.G.

ALARIS Medical Systems® (hereinafter referred to as "ALARIS") warrants that:

(A) Each new instrument (pump, controller or peripheral instrument) is free from defects in material and workmanship under normal use and service for a period of one year from the date of delivery by ALARIS to the first purchaser.

(B) Each new accessory is free from defects in material and workmanship under normal use and service for a period of ninety (90) days from the date of delivery by ALARIS to the first purchaser.

If any product requires service during the applicable warranty period, the purchaser should communicate directly with the local ALARIS service centre to determine appropriate repair facility. Repair or replacement will be carried out at ALARIS's expense, subject to the terms of this warranty. The product requiring service should be returned promptly, properly packed, and postage prepaid. Loss or damage in return shipment to ALARIS shall be at purchaser's risk.

In no event shall ALARIS be liable for any incidental, indirect or consequential damages in connection with the purchase or use of any ALARIS product. This warranty shall not apply to, and ALARIS shall not be responsible for, any loss arising in connection with the purchase or use of any ALARIS product which has been repaired by anyone other than an authorised ALARIS service representative or altered in any way so as, in ALARIS's judgement, to affect its stability or reliability, or which has been subject to misuse or negligence or accident, or which has had the serial or lot number altered, effaced or removed, or which has been used otherwise than in accordance with the instructions furnished by ALARIS.

This warranty is in lieu of all other warranties, express or implied, and of all other obligations of liabilities on ALARIS's part, and ALARIS neither assumes nor authorises any representative or other person to assume for it any other liability in connection with the sale of ALARIS products. See packing inserts for international warranty.

ALARIS DISCLAIMS ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, INCLUDING ANY WARRANTY OF MERCHANTABILITY OR FOR FUNCTION OF FITNESS FOR A PARTICULAR PURPOSE OR APPLICATION.

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