INNO-LiPA HLA-DRB1 Plus

IVD

Manufactured by:

INNOGENETICS N.V.
Technologiepark 6
9052 Ghent
Belgium
☎+32-9 329 13 29

Distributed by:

INNOGENETICS GmbH
Lembecker Straße 19
46359 Heiden (Westfalen)
Germany
☎+49-2867 99 07 0

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English

Intended use

The INNO-LiPA HLA-DRB1 Plus is a line probe assay, for *in vitro* use, designed for the molecular typing of human leukocyte antigen (HLA) DRB1 alleles at the allele group level (DRB1*01 to DRB1*16).

Test principle

The INNO-LiPA HLA typing tests are based on the reverse hybridization principle. Amplified biotinylated DNA material is chemically denatured, and the separated strands are hybridized with specific oligonucleotide probes immobilized as parallel lines on membrane-based strips. This is followed by a stringent wash step to remove any mismatched amplified material. After the stringent wash, streptavidin conjugated with alkaline phosphatase is added and bound to any biotinylated hybrid previously formed. Incubation with a substrate solution containing a chromogen results in a purple/brown precipitate. The reaction is stopped by a wash step, and the reactivity pattern of the probes is recorded.

With the INNO-LiPA HLA-DRB1 Plus kit, an amplification kit (INNO-LiPA HLA-DRB1 Amp Plus) is available for standardized preparation of biotinylated amplified material. This amplification kit is based on the polymerase chain reaction (PCR*).

Amplification products are subsequently hybridized using 1 typing strip on which 37 sequence-specific DNA probes and 2 control probes are fixed (fig.1).

* *The use of this product is covered by a license of F. Hoffmann - La Roche Ltd and Roche Molecular Systems, Inc.*

INNO-LiPA HLA-DRB1 Plus: steps involved

- | | |
|--------|---|
| step 1 | Amplification of exon 2 of the HLA-DRB1 alleles. |
| step 2 | Hybridization and stringent wash with 37 probes immobilized on one INNO-LiPA HLA-DRB1 Plus strip (56°C). |
| step 3 | Color development. |
| step 4 | Interpretation of the probe reactivity pattern. |
| step 5 | If required, the LiRAS™ interpretation software advises for a further amplification and hybridization step. |

Using the INNO-LiPA HLA-DRB1 Plus kit, it is advised to start with (a) and proceed to (b), if ambiguities involving DRB1*03,11,13 or 14 alleles with other DRB1 alleles occur. This will be advised by the LiRAS™ interpretation software.

- the amplification of exon 2 of the HLA-DRB1*03,11,13,14 alleles with the HLA-DRB1 primer solution prior to analysis with the INNO-LiPA HLA-DRB1 Plus strip.
- group-specific amplification of exon 2 of the HLA-DRB1*03,11,13,14 alleles with HLA-DRB1*03,11,13,14 primer solution prior to analysis with the INNO-LiPA HLA-DRB1 Plus strip. The primer solution HLA-DRB1*03,11,13,14 provided, amplifies exon 2 of the HLA-DRB1*03,11,13,14 alleles. Exon 2 of other DRB1 alleles will not be amplified, an exception is the allele DRB1*0820.

The INNO-LiPA HLA-DRB1 Plus is designed to give the best possible resolution, using the reverse hybridization assay format, at the allele group level (this means the first two digits after the asterisk in an allele name when following standard HLA nomenclature e.g. HLA-DRB1*01).

In addition to the allele group interpretation, the software will give all possible allele combinations. This information however, should be considered as extra and not as decisive.

Reagents

Description, preparation for use and recommended storage conditions

- If kept at 2 - 8°C, opened and unopened, and stored in the original vials, the reagents are stable until the expiry date of the kit. Do not freeze reagents. Do not use the reagents beyond the expiry date.
- The kit should be stored isolated from any source of contaminating DNA, especially amplified products.
- All reagents should be brought to room temperature (20 - 25°C) approximately 60 minutes before use and should be returned to the refrigerator immediately after use.
- Alterations in physical appearance of the kit components may indicate instability or deterioration.
- To minimise the possibility that strips curl before use, it is recommended to store the tube horizontally.

Reagents supplied:

Component	Quantity	Ref.	Description
Strips	1 x 20	58191	Containing 20 strips for INNO-LiPA HLA-DRB1 Plus marked with a black marker line.
LiPA Control	1 x 0.05 ml	58192	Containing ethylenediaminetetraacetic acid (EDTA) and 0.01% methylisothiazolon (MIT)/0.1% chloroacetamide (CAA) as preservative.
Denaturation Solution	1 x 1 ml	56718	Alkaline solution containing EDTA. The vial should be closed immediately after use; prolonged exposure of this solution to air leads to a rapid deterioration of the denaturing strength.
Hybridization Solution	1 x 80 ml	56719	Saline sodium phosphate EDTA (SSPE) buffer containing 0.5% sodium dodecyl sulphate (SDS). The Hybridization Solution should be pre-warmed to a temperature of at least 37°C and must not exceed 56°C (all crystals should be dissolved before use).
Stringent Wash Solution	1 x 200 ml	56720	SSPE buffer containing 0.1% SDS. The Stringent Wash Solution should be pre-warmed to a temperature of at least 37°C and must not exceed 56°C (all crystals should be dissolved before use).
Conjugate Diluent	1 x 80 ml	56951	Phosphate buffer containing NaCl, Triton®, protein stabilizers, and 0.01% MIT/0.1% CAA as preservative.
Conjugate 100x	1 x 0.8 ml	56952	Streptavidin labeled with alkaline phosphatase in Tris buffer containing protein stabilizers and 0.01% MIT/0.098% CAA as preservative, to be diluted 1/100 in Conjugate Diluent before use. Prepare 2 ml Conjugate working solution for each test trough + 2 ml in excess (Conjugate working solution can be prepared during stringent wash) for manual testing. For <i>Auto-LiPA</i> , make 10 ml in excess. The Conjugate working solution is stable for 24 hours at room temperature (20 - 25°C) if stored in the dark.
Substrate Buffer	1 x 180 ml	56953	Tris buffer containing NaCl, MgCl ₂ and 0.01% MIT/0.1% CAA as preservative.

Substrate BCIP/NBT 100x1 x 0.8 ml	56954	5-Bromo-4-chloro-3-indolyl phosphate p-toluidine salt (BCIP) and 4-nitroblue tetrazolium (NBT) in dimethylformamide (DMF), to be diluted 1/100 in Substrate buffer before use. Prepare 2 ml Substrate working solution for each test trough + 2 ml in excess (Substrate working solution can be prepared during conjugate incubation) for manual testing. For <i>Auto</i> -LiPA, make 10 ml in excess. The Substrate working solution is stable for 24 hours at room temperature (20 - 25°C) if stored in the dark.
Rinse Solution 5x 1 x 80 ml	56721	Phosphate buffer containing NaCl, Triton®, and 0.05% MIT/0.48% CAA as preservative, to be diluted 1/5 (1 part + 4 parts) in distilled or deionized water before use. Prepare 8 ml Rinse working solution for each test trough + 10 ml in excess. For <i>Auto</i> -LiPA, make 10 ml in excess. The Rinse working solution is stable for 2 weeks at 2 - 8°C.
Incubation tray 3	-	Containing 8 troughs each.
Reading card 1	-	For identification of positive probes.
Data reporting sheet 2	-	For storage of developed strips.

Materials required but not provided

- Water bath with shaking platform (80 rpm; with inclined lid; temperature adjustable to 56°C ± 0.5°C).
- Aspiration apparatus.
- Calibrated thermometer.
- Orbital, reciprocal or rocking platform shaker.

Recommendations:

For an orbital shaker:

- the diameter of the circular motion should be equal or superior to 13 mm
- recommended speed for a 13 mm circular motion is 160 rpm.

For a reciprocal shaker:

- recommended speed for the to and from motion is 80 movements per minute.

For a rocking platform shaker:

- the shaking angle should not exceed 13° to avoid spilling of liquid
- recommended speed is 50 rpm.

- Vortex mixer or equivalent.
- Graduated cylinders (10, 25, 50, and 100 ml).
- Distilled or deionized water.
- Disposable gloves.
- Disposable sterile pipette tips (preferably cotton-plugged).
- Forceps for strip handling.
- Adjustable pipettes to deliver 1 - 20 µl, 20 - 200 µl, and 200 - 1000 µl.
- Dispensing Multipipette (Eppendorf, optional).
- Timer, 2 hours (± 1 minute).

Safety and environment

- **Please refer to the manufacturer's safety data sheet and product labeling for information on potentially hazardous components.**

R20/21, R36, R45, R61, S36/37, S45 53

Toxic! (T) Harmful by inhalation and in contact with skin. Irritating to eyes. May cause cancer. May cause harm to the unborn child. Wear suitable protective clothing and gloves. In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible). Avoid exposure - obtain special instructions for use.

Restricted to professional users. **Contains Dimethylformamide, 5-Bromo-4-chloro-3-indolyl phosphate p-Toluidine salt:** Substrate BCIP/NBT 100x.

R43, S24-37

Irritant! (Xi) Avoid contact with skin. May cause sensitization by skin contact. Wear suitable gloves. **Contains 2-Chloroacetamide:** Rinse Solution, Substrate Buffer, Conjugate Diluent and LiPA Control.

R36/38, S23-24-26

Irritant! (Xi) Irritating to eyes and skin. Do not breathe vapour. Avoid contact with skin. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. **Contains sodium hydroxide:** Denaturation Solution.

- Specimens should be handled as potentially infectious. Therefore, all blood components and biological materials should be considered as being potentially infectious and should be handled as such. Only adequately trained personnel should be permitted to perform the test procedure. All blood components and biological materials should be disposed of in accordance with established safety procedures.
 - Autoclave for at least 15 minutes at 121°C.
 - Incinerate disposable material.
 - Mix liquid waste with sodium hypochlorite so that the final concentration is $\pm 1\%$ sodium hypochlorite. Allow to stand overnight before disposal. CAUTION: Neutralize liquid waste that contains acid before adding sodium hypochlorite.
- Use of personal protective equipment is necessary: gloves and safety spectacles when manipulating dangerous or infectious agents.
- Waste should be handled according to the institution's waste disposal guidelines. All federal, state, and local environmental regulations should also be observed.

Specimens

Since the LiPA test utilizes biotinylated amplified DNA material as specimen, an amplification kit, INNO-LiPA DRB1 Amp Plus, is available as an accompanying tool.

Manipulations procedures

Strip handling

- The strips are designed to be used only once!
- Do not touch the strips with bare hands; use clean forceps.
- Use a **pencil** for identification of the test strips. Do not use ballpoints, etc. Write the ID above the marker line on the strips.
- Throughout the different incubation steps, test strips should always remain in the same trough.
- Unused or developed strips should be kept away from intense light and heat.
- Allow the developed strips to dry completely before interpretation, covering, and storing.
- Developed dry strips should be stored preferably in the dark at 20 - 25°C.

- Do not reuse the troughs.

Directions for manual incubation

- Incubation at $56^{\circ}\text{C} \pm 0.5^{\circ}\text{C}$ during hybridization and stringent wash are the most critical steps in avoiding false positive (temperature too low) or false negative/very weak signals (temperature too high). A shaking water bath with inclined lid allows a good control of temperature variations. Strict temperature control with a calibrated thermometer is necessary.
- Always **close** the lid of the water bath during incubations in order to avoid false positive signals.
- **Do not use a hot air shaker** for the hybridization and stringent wash.
- For hybridization and stringent wash, the troughs should be placed on the shaking platform of the water bath. Adjust the water level between 1/3 and 1/2 of the height of the trough. Make sure that the troughs do not float on the water. The water should be in direct contact with the troughs.
- Incubation steps for the color development should be between 20 and 25°C . If the temperature is below 20°C , weaker results may be obtained. If the temperature is above 25°C , high background and/or false positive signals may be obtained.
- Always incubate exactly for the duration as mentioned in the protocol.
- The amplitude of the motion generated by both the shaking water bath (hybridization procedure) and the shaker (color development procedure) is critical in achieving maximum sensitivity and homogeneous staining. The strip surface should be completely submerged. The amplitude should be as high as possible. **However, spilling of liquid over the edges of the troughs should be avoided!** This can lead to cross-contamination and invalid results.
- Shaking during incubation of the strips should be performed in such a way that both the liquid and the test strips move back and forth in the trough, without liquid being spilled over the edge of the troughs.
- Do not cover the tray. During hybridization and stringent wash incubations, the troughs can be left uncovered in the water bath. Covering the troughs with microplate sealers may cause cross-contamination.

Directions for manual changing liquid in the troughs

- The liquid is aspirated from the trough with a pipette, preferably attached to a vacuum aspirator. The tray is held at an angle to allow all liquid to flow to one end of the trough.
- Add 2 ml of the appropriate solution to each trough and follow the protocol. A dispensing Multipipette (Eppendorf) is useful for this purpose.
- Repeat this step as many times as indicated in the test procedure.

NOTE:

- Do not allow the strips to dry between two steps.
- Make sure not to damage the surface of the strips when aspirating. Aspirate the liquid at the top of the strip above the marker line.
- Make sure all liquid is aspirated.
- Make sure the whole strip is thoroughly washed by complete submersion in the solution.
- Adapt the speed of the shaker when necessary.

Remarks and precautions

- For professional use only.
- In order to avoid DNA contamination, a maximum physical separation between the pre- and post-amplification steps is recommended: separate rooms, separate pipettes and other lab material, separate lab coats and gloves (and their stock) are minimum precautions for a good laboratory practice.
- Avoid any return from the post-amplification room to the pre-amplification room.
- The use of autoclaved disposable lab material is recommended.
- Do not reuse disposable lab material.

- Use a new sterile pipette tip for each aliquoted specimen.
- Do not use the reagents beyond the expiry date.
- Do not mix reagents between kits, unless the components have identical lot numbers.
- Avoid microbial contamination of reagents.

Manual test procedure

NOTE:

- Throughout the different incubation steps, the test strips should always remain in the same trough.

Before incubation, check the temperature of the water bath using a calibrated thermometer, and adjust the temperature if necessary before placing the tray in the water bath. Always close the lid. The LiPA Control should be included in each test run.

Samples

1. DRB1 amplified product (use 10 µl). For sample preparation: see INNO-LiPA HLA-DRB1 Amp Plus instructions for use.
2. LiPA Control sample for INNO-LiPA HLA-DRB1 Plus (use 10 µl; **amplification is not required!**).
3. Blank amplified control sample (negative control; use 10 µl).

Denaturation and hybridization

1. Heat a shaking water bath to $56^{\circ}\text{C} \pm 0.5^{\circ}\text{C}$. Check the temperature using a calibrated thermometer, and adjust the temperature if necessary. Do not exceed the indicated temperatures. Pre-warm the Hybridization Solution and Stringent Wash Solution in a water bath of at least 37°C and do not exceed 56°C . Mix before use. All crystals should be dissolved.
2. Using forceps, remove the required number of INNO-LiPA HLA-DRB1 Plus strips from the tube (one strip per test sample). Include a strip for the LiPA Control sample for INNO-LiPA HLA-DRB1 Plus and for the blank amplified control sample. Pencil an identification number above the black marker line on the strip.
3. Take the required number of test troughs (one trough per test sample) and place them into the tray.
4. Pipette 10 µl Denaturation Solution into the upper corner of each trough.

NOTE:

- Close the vial immediately after use.
5. Add 10 µl sample (see Instructions for use of INNO-LiPA HLA-DRB1 Amp Plus) and carefully mix by pipetting up and down. Always use sterile pipette tips. Allow denaturation to proceed for **5 minutes** at $20 - 25^{\circ}\text{C}$.
 6. Shake the prewarmed Hybridization Solution and **gently** add 2 ml to the denatured amplified product into each trough. Mix by gentle shaking. Take care not to contaminate neighbouring troughs during pipetting.
 7. Immediately place the strip with the marked side of the membrane up into the trough. The strips should be completely submerged in the solution.
- #### NOTE:
- Wear disposable gloves and use forceps.
8. Place the tray into the $56^{\circ}\text{C} \pm 0.5^{\circ}\text{C}$ shaking water bath (80 rpm; see Directions for manual incubation), close the lid, and incubate for 30 minutes.

NOTE:

- Avoid splashing water from the water bath into the trough. Adjust the water level between 1/3 and 1/2 of the height of the trough. To prevent the tray from sliding, immobilize the tray between two heavy weights.

Stringent wash

1. After hybridization, remove the tray from the water bath.
2. Hold the tray at a low angle and aspirate the liquid from the trough with a pipette, preferably attached to a vacuum aspirator. Add 2 ml pre-warmed Stringent Wash Solution into each trough and rinse by rocking the tray briefly (10 - 20 seconds) at 20 - 25°C. Aspirate the solution from each trough. Repeat this brief washing step once.
3. At last, aspirate solution and incubate each strip in 2 ml pre-warmed Stringent Wash Solution in the shaking water bath at 56°C ± 0.5°C for **10 minutes**. Close the lid of the water bath.

Before incubation, check the temperature of the water bath using a calibrated thermometer, and adjust the temperature if necessary. Always close the lid.

NOTE:

- Dilute the concentrated Rinse Solution 5x and Conjugate 100x during stringent wash. See Reagents.

Color development

All subsequent incubations are carried out at **20 - 25°C on a shaker**. During the incubations, the liquid and test strips should move back and forth in the trough for homogeneous staining.

1. Wash each strip twice for **1 minute** using 2 ml of the Rinse working solution.
2. Add 2 ml of the Conjugate working solution to each trough and incubate for 30 minutes while agitating the tray on the shaker.

NOTE:

- Dilute the Substrate BCIP/NBT 100x solution about 10 minutes prior to the end of the conjugate incubation. See Reagents.
3. Wash each strip twice for **1 minute** using 2 ml of the Rinse working solution and wash once more using 2 ml Substrate Buffer.
 4. Add 2 ml of the diluted Substrate working solution to each trough and incubate for **30 minutes** while agitating the tray on the shaker.

CAUTION:

- Wear gloves and protective goggles.
5. Stop the color development by washing the strips twice in 2 ml distilled water while agitating the tray on the shaker for at least **3 minutes**.
 6. Using forceps, remove the strips from the troughs and place them on absorbent paper. Let the strips dry completely before reading the results. Store the developed and dried strips in the dark.

Automated test procedure: Auto-LiPA

The LiPA test procedure is extremely well-suited for automation. Therefore, the *Auto-LiPA* was designed to fully handle hybridization, stringent wash and color development steps. The *Auto-LiPA* is featured as a walk-away system with automated heating and cooling, and with automated aspiration and pipetting.

For more information and specific protocols on the *Auto-LiPA*, please contact your local distributor.

Results

Reading

A line is considered positive when a clear purple/brown band appears at the end of the test procedure.


Figure 1: Location of the marker line (black), the conjugate control line (conj.control), the HLA-DRB1 Plus control line (HLA-DRB. Control) and the 37 sequence-specific DNA probes on the INNO-LiPA HLA-DRB1 Plus strip.

Validation

- Include one negative control and the LiPA Control for INNO-LiPA HLA-DRB1 Plus each time a test is performed. As with any new test procedure, the inclusion of additional positive and negative controls is to be considered until a high degree of confidence is reached in the ability to correctly perform the test procedure. It is suggested that reference DNA should be used to validate that the proper test conditions have been established.
- The uppermost line on the INNO-LiPA HLA-DRB1 Plus strip is the marker line (black line). This line allows correct orientation of the strip.
- The next line controls the addition of reactive Conjugate and Substrate Solution during the detection procedure. This line should be lined up with the Conjugate control line on the plastic reading card. This line should always be positive and should have approximately the same intensity on each strip in the same test run.
- The next line is the HLA-DRB positive control line (HLA-DRB control on the reading card) which indicates whether or not an appropriate amount of HLA-DRB1 amplified material was added for hybridization. This line should always be positive, except for the negative control, but its intensity may vary between samples.
- The assay result of each negative control should give no apparent signal for any of the lines on the strip, except for the conjugate control line.
- The LiPA Control sample for INNO-LiPA HLA-DRB1 Plus should yield the following pattern of positive probes: conjugate control, HLA-DRB control, 2, 8, 23, 25, 26, 30, 33, 37. This control sample is composed of biotinylated oligonucleotides, complementary to the mentioned hybridization probes. Positive hybridization on these probes thus demonstrates satisfactory performance of the assay including hybridization, stringent washing, color development, and detection. A different pattern may indicate assay performance problems such as incorrect temperature during hybridization or deviations from the prescribed incubation times or temperatures.
- Color intensities of probes on a strip may differ from one line to the other.

Interpretation of results

- Check for the correct reactivity pattern of the LiPA Control sample.
- The conjugate control line on the strip should be lined up with the conjugate control line on the plastic reading card.
- Check the positivity of the control lines (first two lines) in order to validate each individual strip (see Validation).
- Typing results are based on the reactivity of the probes in the kit. A list of probe specificities is available on request.
- Identify all probe numbers that are positive on the INNO-LiPA HLA-DRB1 Plus strip and deduce the DRB1 type by using the INNO-LiPA HLA-DRB1 Plus typing table or a version of the LiPA interpretation software (LiRAS™). The software often gives more results than the intended use suggests: interpretation at allele

group is given, as well as all possible allele combinations.

This information however, should be considered as extra, not as decisive.

Interpretation software: LiRAS™

The LiRAS™ software is designed to assist with the interpretation of the LiPA results. Please contact your local distributor to obtain the latest updated version.

Limitations of the procedure

- Use of this product should be limited to personnel trained in the techniques of hybridization.
- Only good laboratory practice and careful performance of the procedures specified, will allow specific hybridization and correct typing of target DNA.
- The INNO-LiPA HLA-DRB1 Plus is designed to give resolution at the allele group level of HLA-DRB1 (this means the first two digits after the asterisk in an allele name when following standard HLA nomenclature e.g. HLA-DRB1*01). However, combinations of certain alleles could lead to a non-unique probe combination and thus to an ambiguous answer of two or more possible allelic combinations. Knowledge of the allele frequencies in the population under study will indicate the most probable DRB1 type. In case a higher resolution typing is required, use the INNO-LiPA HLA-DRB decoder kit.
- New alleles might have polymorphisms outside the probe regions, this should be taken into consideration when interpreting results.
- The sequence at the primer binding site is not known for all alleles. The risk of mismatch at the primer sites is considered to be very low as primer failure has never been observed and more than 80 alleles spanning all groups have been tested. However, the possibility of allelic drop-out should be considered in the event of a homozygous result.

Test performance

Test performance using DNA obtained with DRB1 Primer Solution in combination with PCR protocol II (DRB1; DRB1*3+4+5; 86G; 86V; DQB1)

INNO-LiPA HLA-DRB1 has been evaluated in three European and two North American tissue typing laboratories on routine typing samples.

In the study a total of 243 samples were tested with INNO-LiPAHLA-DRB1.

The results were compared to established classical typing techniques (DNA-based assays (INNO-LiPA DRB *key/decoder*, SSP) and/or serology).

There were no INNO-LiPA HLA-DRB1 amplification failures during the study.

Probe reactivity

In total 74 false positive bands were observed with INNO-LiPA HLA-DRB1. No false negative bands were reported during the study. Consult the guidelines for specific information about probe reactivity. A complete list of probe and primer specificities is available. All probes are designed to hybridize specifically with their complementary sequence at the level of one mismatch, unless stated otherwise in the list.

A report on the quality control of the reactivities is available for accreditation purposes on request.

Application range

The lower limit of the application range was determined using a DNA titration series. As little as 0.05 µg extracted DNA ($A_{260\text{ nm}}/A_{280\text{ nm}} > 1.5$) can be used in the accompanying amplification assay, yielding more than sufficient DNA for detection in the LiPA.

Resolution

(Typing table used: K-1067/v.5.x/1998-02-20)

The INNO-LiPA HLA-DRB1 strips, with 37 DNA probes immobilized as parallel lines, allow DRB1 allele-group typing. The theoretical resolution at broad level (DRB1*01 to

DRB1*16) is calculated as 89.3%. During the study the overall allele-group resolution achieved was 86.8% (211/243). A total of 32 ambiguous results were described. These arose from seven probe combinations and were reported with the following allele combinations:

DRB1*03 x *11/*13 (8 samples)
DRB1*04 x *11/*13 (3 samples)
DRB1*01 x *13/*14 (3 samples)
DRB1*04 x *13/*14 (6 samples)
and DRB1*03 x *13/*15 (12 samples)

One antigen was ambiguous for **DRB1*11/*13** in 11/83 (13.3%) samples where either DRB1*11 or DRB1*13 was present. These ambiguity causing alleles were reported by the LiPA-Expert software as follows:- DRB1*0308 x *1324 (two samples), DRB1*0415 x *1322/*1323 (three samples) and DRB1*0308 x *1305/*1307/*1311/*1314c75L/*1314c75V (six samples).

In the presence of DRB1*03, where DRB1*13 is given as the alternative allele-group to DRB1*11, the specific DRB1*03 allele is, in each case, DRB1*0308. This allele has been described as having a very similar sequence to DRB1*03011 except at codon 58 where the similarity is with the DRB1*11 allele group. In addition, both the DRB1*0308 and the specific DRB1*13 alleles (DRB1*1305, *1307, *1311, *1314 and *1324), with which it was ambiguously reported, also have an expected low frequency.

In the presence of DRB1*04, either of two rare alleles of DRB1*13 (*1322, *1323) is given as the alternative to DRB1*11 allele-group typing and the only DRB1*04 allele they were reported with is DRB1*0415. There is documented sequence similarities between the specific DRB1*13 alleles involved and the DRB1*11 allelegroup. The frequencies of the DRB1*13 alleles causing ambiguity and of the DRB1*0415 allele are also expected to be low.

One antigen was ambiguous for **DRB1*13/*14** in 9/66 (13.6%) samples where either DRB1*13 or DRB1*14 was present. These ambiguity causing alleles were reported by the LiPA-Expert software as follows:- DRB1*0103 x *1417/*1430 (three samples), DRB1*0402/*0414c13H x *1421 (five samples) and DRB1*0402/ *0414 x *1417/*1430 (one sample).

In the presence of DRB1*04, rare and recently described alleles of DRB1*14 (*1417, *1421 and *1430) are given as the alternative to DRB1*13 allele group typing. The specific DRB1*14 alternatives (*1417 and *1430 - see above) to the DRB1*13 alleles in the presence of DRB1*01 are both expected to occur rarely.

One antigen was ambiguous for **DRB1*13/*15** in 12/97 (12.4%) samples where either DRB1*13 or DRB1*15 was present. These ambiguity causing alleles were reported by the LiPA-Expert software as DRB1*0304 x *1309 (all twelve samples).

The DRB1*1309 alternative to DRB1*15 in the presence of DRB1*03 is a rare allele of Hispanic origin. The DRB1*0304 allele it is reported with is also rare having first been detected in three unrelated Caucasians. This allele has sequence similarities with DRB1*0301, which was given as the first alternative.

Reproducibility

Repeated typing, between differing batches of INNO-LiPA HLA-DRB1, using the DNA derived from clinical samples during the study, resulted in reproducible probe reactivities in every case.

Accuracy

During the external evaluation INNO-LiPA HLA-DRB1 showed complete concordance with the reference methods.

Summary of sample numbers

243	Total number of amplification attempts (no failures)
211	Total number of unambiguous results at allele group level (32 ambiguities)
203/243	Number of samples compared with INNO-LiPA DRB key (used to measure concordancy)
40/243	Number of samples tested by PCR-SSP (also used to measure concordancy).

Test performance using DNA obtained with DRB1 Primer Solution in combination with PCR protocol I (A; B; Bw4; C; DRB1; DRB1*03,11,13,14)

In order to evaluate the equivalence between PCR protocols I and II, 42 samples were tested in-house using both methods. The PCR products of both protocols gave identical typing results for all samples.

Test performance using DNA obtained with DRB1*03,11,13,14 Primer Solution in combination with PCR protocol I (A; B; Bw4; C; DRB1; DRB1*03,11,13,14)

The performance of the DRB1*03,11,13,14 Primer Solution of INNO-LiPA HLA-DRB1 Plus was evaluated in-house and in two Belgian tissue typing laboratories. A total of 89 anonymous samples was analyzed (50 samples externally; 39 in-house). All samples had previously been typed with INNO-LiPA HLA-DRB1. The samples included 14 ambiguous (six different ambiguities) and 75 unambiguous samples. All ambiguous samples were further typed to a unique result with one or more alternative DNA-based typing assays.

All hybridizations were performed using the *Auto*-LiPA instrument, using program HLA56v3, and the interpretation was done with a clinical trial version of the LiRAS™ for LiPA HLA v3.00 interpretation software.

Test sensitivity

Test sensitivity of the assay at allele group level was calculated as the number of successful typing results over the total number of samples tested.

Test sensitivity in this evaluation was 100% (89/89; 95% CI [95.9%; 100%]).

Probe reactivity

On one or more occasions during this study, some probes gave weak reactions. Concordant typing at allele group level was achieved when the LiRAS™ interpretation software was used.

Resolution

(typing table v5.8/021112)

The DRB1 Primer Solution, used in combination with the INNO-LiPA HLA-DRB1 strip, is designed for typing HLA-DRB1 alleles at allele group level.

The DRB1*03,11,13,14 Primer Solution, used in combination with the INNO-LiPA HLA-DRB1 strip, is designed to resolve the most frequent ambiguities, i.e., ambiguities involving DRB1*03, 11, 13, and/or 14 alleles.

All ambiguities, at allele group level, observed after testing with the DRB1 Primer Solution, were reduced to an unambiguous typing result using the DRB1*03,11,13,14 Primer Solution: DRB1*01, DRB1*13 or DRB1*01, DRB1*14 (3 cases); DRB1*04, DRB1*13 or DRB1*04, DRB1*14 (6 cases); DRB1*08, DRB1*11 or DRB1*11, DRB1*13 (1 case); DRB1*11, DRB1*12 or DRB1*12, DRB1*13 (2 cases).

For all ambiguous results, the LiRAS™ interpretation software advised DRB1*03,11,13,14 amplification followed by hybridization with the INNO-LiPA HLA-DRB1 strip.

Precision

Repeatability (intra-assay), inter-assay and inter-person variation were evaluated in-house. One sample, amplified in duplicate, using three different cyclers and processing two runs per cycler, was tested using the *Auto-LiPA* instrument. This protocol was performed in duplicate by two different persons.

No significant variation was observed and an identical typing result was obtained in all cases without the need for editing.

Accuracy

Calculations of diagnostic accuracy were only carried out on samples for which a unique reference result was available. As a consequence, two samples were excluded from this calculation.

A unique HLA-DRB1 typing result was obtained for the remaining 87 samples.

Diagnostic accuracy of the assay at allele group level was calculated as the number of correctly obtained typing results, compared to the unique typing result obtained with INNO-LiPA HLA-DRB1 (Innogenetics) or an alternative typing assay (Olerup SSP (GenoVision) or Allset⁺ DR Low Resolution (Dynal)).

The typing results of all 87 samples were concordant with the reference assay, resulting in a diagnostic accuracy of 100% (87/87; 95% CI [95.9%; 100%]).

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