Vitrification Freeze Kit (Vit Kit® - Freeze)

Catalog No. 90133-HSV Includes:

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
IrvineScientific®

Figure 1:

H = Modified HTF with HEPES (90126) + protein
ES = Equilibration Solution (90131)


Figure 2:


Figure 3:

Prepare CryoTip for loading


Figure 4:

Prepare HSV straw for Loading


Figure 5:


Figure 6a:


Figure 6b:



Figure 7:

INTENDED USE

Vit Kit® - Freeze is intended for use in assisted reproductive procedures for vitrification and storage of human oocytes (MI) and embryos (zygote to blastocyst). This kit is designed for use with Irvine Scientific's CryoTip® (Catalog No. 40709) or HSV Straw (Catalog No. 019921), and Vitrification Thaw Kit (Vit Kit® - Thaw) for optimal recovery of specimens.

PRODUCT DESCRIPTION

Equilibration Solution-ES is a HEPES buffered solution of Medium-199 containing gentamicin sulfate (35 µg/mL), 7.5% (v/v) of each DMSO and ethylene glycol and 20% (v/v) Dextran Serum Supplement.

Vitrification Solution-VS is a HEPES buffered solution of Medium-199 containing gentamicin sulfate (35 µg/mL), 15% (v/v) of each DMSO and ethylene glycol, 20% (v/v) Dextran Serum Supplement and 0.5 M sucrose.

These two solutions are to be used in sequence according to the step-wise microdrop vitrification protocol.

QUALITY ASSURANCE

The solutions in Vit Kit-Freeze are membrane filtered and aseptically processed according to manufacturing procedures which have been validated to meet a sterility assurance level (SAL) of 10⁻³.

Each lot of Vit Kit-Freeze receives the following tests: Solutions and CryoTips.

- Endotoxin by Limulus Amebocyte Lysate (LAL) methodology
- Biocompatibility by Mouse Embryo Assay (one-cell)
- Sterility by the current USP Sterility Test <71>

All results are reported on a lot specific Certificate of Analysis which is available upon request.

MATERIALS REQUIRED BUT NOT INCLUDED

- Irvine Scientific CryoTip® (Catalog No. 40709) or HSV Straw (Catalog No. 0019921)
- Irvine Scientific Connector (Catalog No. 40736) or other adaptor
- Sterile Petri Dishes (50 X 9 mm, Falcon 351006 or equivalent)
- Cryotubes (4.5 mL) or goblets or cryocanes
- Modified HTF - HEPES (Catalog #90126) culture medium supplemented with protein
- Disposable gloves
- Hamilton GASTIGHT® Syringe (50 µL, Catalog #80901) or other aspiration tool
- Transfer pipettes (pulled glass pipettes or micropipette tips with an inner tip diameter of ~200µm)
- Tweezers
- Impulse Heat Sealer
- SYMS Sealer for HSV Straw (Catalog No. 016296)
- Stopwatch or timer
- Liquid Nitrogen Reservoir (Dewar or Styrofoam container with lid, 1-2 L volume)
- Liquid Nitrogen (sufficient volume to achieve 4 inch depth in reservoir)

DIRECTIONS FOR USE

- Vit Kit-Freeze component requirements (per application):
- Equilibration Solution (ES):
60 µL for **Oocyte Vitrification Protocol**
Or
50 µL for **Embryo Vitrification Protocol**
 - Vitrification Solution (VS):
50 µL for Vitrification Protocol
 - 1 CryoTip or HSV Straw (**stores up to 2 specimens**)
 - 1 Connector

VITRIFICATION PROTOCOL:

NOTE: Procedures are to be done at room temperature (20-27° C). DO NOT use heated microscope stage for the following procedures. CAUTION: Minimize exposure of specimen to light during equilibration in ES and VS solutions.

- Bring the quantity to be used of ES and VS to room temperature (20-27° C). NOTE: Avoid bringing the entire vials of ES and VS to room temperature repeatedly when a partial of the solution is needed each time. It is better to aliquot the quantity to be used and return the vials to 2-8°C right after aliquoting. Modified HTF (HEPES) with protein is also required for oocyte vitrification protocol.
- Fill the liquid nitrogen reservoir with liquid nitrogen (sufficient to achieve a depth of 4 inches or to completely submerge cryotube on cane) and place close to microscope. Attach a cryotube or goblet (uncapped) to the bottom clamp of a cryocane and submerge in the liquid nitrogen in preparation for storage of the vitrified specimens.
- Determine the number of specimens to be vitrified.
- Label each sterile petri dish (or lid) and Cryo storage device with necessary information.
- Gently invert each vial of ES and VS twice to mix contents before use.
- Prepare dish with droplets of solutions for Vitrification Procedure as follows:
 - OOCYTE (MI) Vitrification Protocol:** NOTE: Refer to Section B for embryo vitrification protocol.
 - Aseptically dispense 20 µL drop of culture medium, Modified HTF - HEPES with protein, and ES in close proximity on an inverted lid of sterile petri dish as shown in **Figure 1**, and place the dish on the microscope stage:
 - one- 20 µL drop of Modified HTF (HEPES with protein)
 - three- 20 µL drops (60 µL total) of ES (ES1, ES2, ES3)
 - Remove the culture dish containing MI oocytes from the incubator and check the quality of the specimens under microscope. **Where possible, select only the best quality MI stage oocyte(s).** **CAUTION: Minimize the exposure of the specimen(s) to light during equilibration in the H, ES and VS drops.**
 - Transfer the oocyte (up to 2 at a time) with minimal volume of medium from the culture dish (in incubator) into the 20 µL drop of H for one minute.
 - Merge the drop of H to ES1 (See Fig. 1, arrow 1) with the tip of the transfer pipette and allow spontaneous mixing of the two solutions to occur for **2 minutes**.
 - Then merge the drop of ES2 (arrow 2) to the previously merged drops and leave for **2 minutes**.
 - Transfer the oocyte(s) with minimal volume of solution from merged drop to ES3 drop for **6-10 minutes**. **Note: equilibration of oocyte(s) in ES3 is complete when the thickness of the zona pellucida and perivitelline space is equal. The oocyte(s) will settle to the bottom of the drop within 3 minutes.**
 - During the equilibration time in ES3: Aseptically dispense one (1) 50 µL drop of VS 2 minutes prior to complete equilibration and prepare the cryotip (fig. 3) or HSV straw (fig. 4) for loading:**
 - CryoTip:** connecting to the Hamilton syringe or appropriate aspiration tool using a connector or adaptor to assure a tight seal. NOTE: Keep the metal cover sleeve over the fine pulled tip to protect it until ready to load specimens.
 - HSV Straw:** connecting the longer end of the blue plastic insertion device to the colored end of the handling rod.
 - The following steps (9-13) should be completed in 90-110 seconds. CAUTION: Exposure of specimens to VS should be limited to prevent cytotoxicity. Specimen(s) tend to float in VS, so adjust the focus through the microscope to maintain continuous visualization during exposure and keep the tip of the transfer pipette nearby to assure rapid transfer between VS drops. Refer to Figure 5.**
- After equilibration in ES is complete, draw up some ES into the transfer pipet and transfer the specimen(s) with minimal volume from the drop of ES into the **drop of VS for 30 seconds**.
- Load and heat seal the CryoTip as follows (See Figure 6a):**
 - Slide the metal cover sleeve up along the CryoTip to expose the fine fragile tip end.
 - Handling the CryoTip and Hamilton syringe while observing under the microscope, carefully aspirate a small volume of VS to the **Mark #1** on the CryoTip.
 - Continue observation under the microscope and gently aspirate the specimen with VS to the **Mark #2** on the CryoTip.
 - Now observe the CryoTip directly and aspirate more VS to the **Mark #3**.
 - Specimen must be located between **Mark #2 and Mark #3**.
 - Heat seal (**Seal #1**) the CryoTip on (or just below) **Mark #1**, and slide the cover sleeve back down to cover and protect the fine fragile tip.
 - Carefully remove the CryoTip from the aspiration tool and adaptor and then heat seal (**Seal #2**) at the thick end of the CryoTip above the **Mark #4**.
 - Plunge the covered CryoTip directly into liquid nitrogen (cooling at a rate of ~12,000° C/min) (See Figure 6b).

Load and seal the HSV straw as follows:

- Using a micropipette, carefully deposit the specimen(s) into the gutter of the capillary rod at 1 mm from the end. The drop holding the specimen(s) must be under 0.5 µL. Maximum of 2 oocytes or embryos for each capillary rod.
- Immediately place the capillary rod and handler into the straw and push until the rectangular portion of the handler comes in contact with the flared end of the straw.
- Slightly pinch the straw between your thumb and finger and remove the insertion device.
- While still holding the straw in place, seal the open end using a SYMS sealer.
- Hold the straw using tweezers in the area of the handling rod.
- Quickly plunge the entire straw into LN₂ vertically. Gently stir the straw in LN₂ for a few seconds so as to avoid formation of an isolating air bubble layer around the straw.
- Place the vitrified CryoTip or HSV straw into the submerged LN₂ filled cryotube or goblet (on the cryocane). Cap the cryotube (or goblet) or attach up-side-down with another uncapped cryotube in order to secure the vitrified CryoTip or HSV straw in liquid nitrogen.
- Move the LN₂ reservoir close to the LN₂ cryo-freezer and transfer the cryocane with contents to the cryo-freezer for long-term storage.

B. EMBRYOS (PN to Blastocyst)

Vitrification Protocol:

- Aseptically dispense one- 50 µL drop of ES on an inverted lid of Petri dish.
- Remove the culture dish with embryo(s) from the incubator and check the quality of the specimen(s) under microscope. **Where possible, select only the best quality embryo(s) for vitrification.**
- Carefully transfer the specimen (up to two at one time) with a minimal volume of medium from the culture dish to the drop of ES and start the timer. Embryos should equilibrate in the ES drop slowly by free-fall for 6-10 minutes. **Note: The specimen will shrink and then gradually return to its original size, which indicates that the equilibration is complete. CAUTION: Minimize the exposure of specimen(s) to light during equilibration in ES and VS drops.**
- During this equilibration time in ES:**
 - set up one-50 µL drop of VS solution as shown in Fig 7 and prepare the CryoTip or HSV Straw for loading.

Follow protocol as written above (Section A - Oocyte [MI] Vitrification Protocol) from steps 9 through 12 for exposure to VS solutions, loading of CryoTip or HSV Straw, plunging in LN₂ and long term storage.

For additional details on the use of these products, each laboratory should consult its own laboratory procedures and protocols which have been specifically developed and optimized for your individual medical program.

STORAGE INSTRUCTIONS AND STABILITY

Store the unopened vials refrigerated at 2°C to 8°C. When stored as directed, Vitrification Freeze Kit Solutions are stable until the expiration date shown on the vial labels.

Do not use media for more than eight (8) weeks once containers have been opened.

As human source material is present in the product it may develop some particulate matter during storage. This type of particulate matter is not known to have an effect on product performance.

PRECAUTIONS AND WARNINGS

This device is intended to be used by staff trained in assisted reproductive procedures that include the indicated application for which the device is intended.

As an added precaution during the preparation procedure, we recommend that each Cryotip be carefully examined when taken out of the package. Prior to use, the CryoTips should be examined under a suitable magnification (40x power) for possible damage (such as tip breakages or cracks) that may have occurred during transport.

Do not use any vial of solution that shows evidence of particulate matter or cloudiness.

To avoid problems with contamination, handle using aseptic techniques.

Vitrification Freeze Kit Solutions contain the antibiotic gentamicin sulfate. Appropriate precautions should be taken to ensure that the patient is not sensitized to this antibiotic.

Currently, research literature indicates the long-term effects of vitrification on embryos remains unknown.

***Human blood derivatives used in the manufacture of this product has been tested by FDA licensed kits, and found to be non-reactive for the Hepatitis B surface antigen (HBsAg), antibodies to Hepatitis C (HCV) and antibodies to Human Immunodeficiency Virus (HIV). However, no test method offers complete assurance that products derived from human sources are noninfectious. Handle all Human blood derivatives as if it were capable of transmitting infection, using universal precautions. Donors of the source material have also been screened for risk of exposure to Cruetzfeldt-Jakob Disease (CJD).**

Do not use any bottle in which the sterile packaging has been compromised.

Caution: Federal (U.S.) law restricts this device to sale by or on the order of a physician.

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
REFERENCES

1. Kuwayama M, Vajta G, Kato O, Liebo SP. Highly efficient vitrification of human oocytes. *RBM Online* 11:300-308, 2005.
2. Stehlik E, Stehlik J, Katayama KP, Kuwayama M, Jambor V, Brohammer R, Kato O. Vitrification demonstrates significant improvement versus slow freezing of human blastocysts. *RBM Online* 11:53-57, 2005.
3. Takahashi K, Mukaida T, Goto T, Oka C. Perinatal outcome of blastocyst transfer with vitrification using cryoloop: a 4-year follow-up study. *Fertil Steril* 84:88-92, 2005.
4. Kuwayama M, Vajta G, Ieda S, Kato O. Comparison of open and closed methods for vitrification of human embryos and the elimination of potential contamination. *RBM Online* 11:608-614, 2005.
5. Cremades N, Sousa M, Silva J, Viana P, Sousa S, Oliveira C, Teixeira de Silva J, Barros A. Experimental vitrification of human compacted morulae and early blastocysts using fine diameter plastic micropipettes. *Human Reproduction* 19, 300-305, 2004.
6. Mukaida T, Nakamura S, Tomiyama T, Wasa S, Oka C, Kasai M, Takahashi K. Vitrification of human blastocysts using cryoloops: clinical outcome of 223 cycles. *Hum Reprod* 18:384-391, 2003.
7. Yoon TK, Kim TJ, Park SE, Hong SW, Ko JJ, Chung HM, Cha KY. Live births after vitrification of oocytes in a stimulated in vitro fertilization- embryo transfer program. *Fertil Steril* 79:1323-1326, 2003.
8. Mukaida T, Takahashi K, Kasai M. Blastocyst cryopreservation: ultrarapid vitrification using cryoloop technique. *RBM Online* 6: 221-225, 2002.
9. Yokota Y, Sato S, Yokota M, Ishikawa Y, Makita M, Asada T, and Araki Y. Successful pregnancy following blastocyst vitrification. *Hum. Reprod.* 15:1802-1803, 2000.
10. Yoon TK, Chung HM, Lim JM, Han SE, Ko JJ, and Cha KY. Pregnancy and delivery of healthy infants developed from vitrified oocytes in a stimulated in vitro fertilization-embryo transfer program. *Fertil. Steril.* 74:180-181, 2000.
11. Chen S-U, Lien Y-R, Cheng Y-Y, Chen H-F, Ho H-N and Yang Y-S. Vitrification of mouse oocytes using closed pulled straws (CPS) achieves a high survival and preserves good patterns of meiotic spindles, compared with conventional straws, open pulled straws (OPS) and grids. *Hum. Reprod.* 16: 2350-2356, 2001.
12. Isachenko V, Selman H, Isachenko E, Montag M, El-Danasouri I, and Nawroth F. Modified vitrification of human pronuclear oocytes: efficacy and effect on ultrastructure. *Reproductive BioMedicine Online* 7,211-216, 2003.
13. Van den Abbeel E, Camus M, Van Waesberghe L, Devroey P, Van Steirteghem AC. A randomized comparison of the cryopreservation of one-cell human embryos with a slow controlled-rate cooling procedure or a rapid cooling procedure by direct plunging into liquid nitrogen. *Hum. Reprod.* 12:1554-1560, 1997.
14. Walker DL, Tummon IS, Hammitt DG, Session DR, Dumesic DA, Thornhill AR. Vitrification versus programmable rate freezing of late stage murine embryos: a randomized comparison prior to application in clinical IVF. *Reprod. BioMed Online* 8: 558-568, 2004.
15. Kasai M. Vitrification: a refined strategy for the preservation of mammalian embryos. [Review] *J. Mamm. Ova Res.*, 14:17-28, 1997.
16. Ohta N, Nohara M and Kojimahara T. Ultrarapid freezing of human embryos by vitrification method, A case of delivery. *Jpn. J. Fertil. Steril.* 41:276-279, 1996.
17. Feichtinger W, Hochfellner C, and Ferst C. Clinical experience with ultra-rapid freezing of embryos. *Hum. Reprod.*, 6:735-736, 1991.
18. Rall WF. Factors affecting the survival of mouse embryos cryopreserved by vitrification. *Cryobiology* 24:387-402, 1987.
19. Rall WF, Fahy GM. Ice-free cryopreservation of mouse embryos at -196°C by vitrification. *Nature* 313:573-575, 1985.

GASTIGHT® is a Registered Trademark of Hamilton Co.

IrvineScientific®

2511 Daimler Street, Santa Ana, California 92705-5588

Telephone: 1 949 261 7800 • 1 800 437 5706

Fax: 1 949 261 6522 • www.irvinesci.com

PN 40695 Rev.17