**Manuale di Istruzione Uso e Manutenzione**  
**User and Maintenance manual**  
**Manuel d'Utilisation et d'Entretien**  
**Gebrauchs- und Wartungsanleitung**  
**Manual de Instrucciones de Uso y Mantenimiento**

## **Intra Max Plus**

**carlo de giorgi** s.r.l.

---

**Dati modello****Model data / Donnés du modèle / Modelldaten / Datos del modelo**

---

<b>Modello</b> <b>Model</b> <b>Modèle</b> <b>Modell</b> <b>Modelo</b>	<b>Alimentazione</b> <b>Power supply</b> <b>Alimentation</b> <b>Betriebs-spannung</b> <b>Alimentación</b>	<b>Anno di costruzione</b> <b>Year of manufacture</b> <b>Année de fabrication</b> <b>Baujahr</b> <b>Año de fabricación</b>
559/00	<input type="checkbox"/> 230 V 50 Hz <input type="checkbox"/> 120 V 60 Hz	<input type="checkbox"/> 2005 <input type="checkbox"/> 2006 <input type="checkbox"/> 2007 <input type="checkbox"/> 2008 <input type="checkbox"/> 2009

**carlo de giorgi s.r.l.**

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---

## CONTENTS

1 - STANDARDS AND GENERAL PRECAUTIONS .....	3
1.1 - INTRODUCTION .....	3
1.2 - DESCRIPTION OF THE APPLIANCE .....	3
1.3 - REFERENCE STANDARDS .....	3
1.3.1 - Obligatory compliance .....	3
1.3.2 - Voluntary compliance .....	4
1.4 - RESPONSIBILITIES OF THE CUSTOMER.....	4
1.5 - REQUESTS FOR SERVICE AND SPARE PARTS.....	4
2- TECHNICAL FEATURES .....	4
3 - COMPONENTS.....	5
4 - OPERATION AND USE .....	5
4.1 - OPERATION .....	5
4.2 - OPERATOR .....	8
4.3 – INTENDED USE .....	8
4.4 - SAFETY WARNINGS.....	8
5 - OPERATOR INSTRUCTIONS .....	9
5.1 - CONTROLS.....	9
5.2 - INSTRUCTIONS FOR USE .....	11
5.2.1 - Setting into operation.....	12
5.2.2 - Normal deactivation and emergency deactivation procedures.....	12
5.2.3 – Editing a program .....	12
6 - MAINTENANCE .....	13
6.1 - CLEANING .....	14
6.2 - STERILIZATION.....	14
7 - DIAGNOSTICS.....	15
8 - ACCESSORIES.....	16
8.1 - MOTOR .....	16
8.2 – TUBE FOR IRRIGATION .....	16
8.3 - UNION .....	16

# 1 - STANDARDS AND GENERAL PRECAUTIONS

## 1.1 - INTRODUCTION

This manual is the property of CARLO DE GIORGI S.r.l. The contents of this manual cannot be reproduced or disclosed to other parties. All rights reserved.

## 1.2 - DESCRIPTION OF THE APPLIANCE

The appliance is a micromotor for use in the dental sector for surgery and implantology, only when equipped with an implantology handpiece.

This device must be used exclusively by suitably trained persons who are authorised to perform the practices in question and who have accumulated sufficient experience in the use of appliances of this type.

Before using the micromotor it is mandatory to read the contents of this manual in their entirety, learning about the position of and function performed by all the controls and paying due attention to all the warnings and precautions that are prescribed. Failure to observe the above obligatory instruction or improper use of the micromotor may cause injury to the operator, to the patient, or to other persons present in the area.

IT IS STRICTLY PROHIBITED TO ATTEMPT TO OPEN OR REPAIR THE APPLIANCE;  
IN THE EVENT OF FAULTS OR ANOMALIES ALWAYS CONTACT CARLO DE GIORGI.

## 1.3 - REFERENCE STANDARDS

### 1.3.1 - Obligatory compliance

- EEC Directive n° 93/42 – Medical Devices (Italian Legislative Decree D.Lgs. n° 46/1997)
- EEC Directive n° 98/37 – Machinery Directive (Italian Presidential Decree D.P.R. n° 459/1996)
- EEC directive n° 73/23 – Low voltage Directive (Italian legislation: Law n° 791/1977, Legislative decrees D.Lgs n° 626/1996 and D.Lgs. n° 277/97)
- EEC directive n° 89/336 relative to Electromagnetic Compatibility (EMC) - (Italian legislation: Legislative decree D.Lgs n° 615/1996)
- DPR 27.4.1955 n° 547, "Norme per la prevenzione degli infortuni sul lavoro" (work accident prevention regulations)
- DPR 27.4.1956 n° 303, "Norme generali per l'igiene del lavoro" (general work hygiene regulations)
- Italian Legislative Decree 12 August 1991 n° 277 implementing EEC directives n° 80/605; n° 83/477; n° 86/188; n° 88/642 concerning the protection of workers against risks deriving from exposure to chemical, physical and biological agents in the workplace, as defined by art. 7 of Italian Law 30 July 1990, n° 212.

### 1.3.2 - Voluntary compliance

- EN ISO 12100 (2005) Machine safety – Basic concepts; general design principles – Part 1 – Basic terminology and methods (ex EN 292-1) /Part 2 – Technical principles (ex EN 292-2)

### 1.4 - RESPONSIBILITIES OF THE CUSTOMER

The user must install the appliance in a suitable place equipped with an electrical system that complies with established legislation.

The appliance must be installed in a dry environment that is illuminated in accordance with the requirements of established legislation.

NOTE: The expression "established legislation" refers to the laws in force in the user's country.

### 1.5 - REQUESTS FOR SERVICE AND SPARE PARTS

For all electrical maintenance requirements, contact CARLO DE GIORGI of Baranzate di Bollate, Via Tonale, 1.

In the case of faults, notify CARLO DE GIORGI technical personnel of the exact nature of the problem, specifying also the lot number shown on the label affixed to the rear of the appliance.

## 2- TECHNICAL FEATURES

		559/00
Voltage	V	230/120
Frequency	Hz	50/60
Max. power consumption	W	200
Nominal output voltage	V	30
Delayed acting fuse T	A	2x2.5 (230V) 2x5 (120V)
Mass of the appliance	kg	9

Class II app.

type BF

The appliance features double electrical insulation (Class II).

### 3 - COMPONENTS

The package contains:

- a) power supply unit;
- b) power cord;
- c) micromotor;
- d) support for micromotor;
- e) footswitch;
- f) disposable irrigator;
- g) saline solution pipeline connection;
- h) stand for saline solution flask;
- i) flask of saline solution with plastic cradle;
- j) user instructions;
- k) motor endcap;

The appliance can be moved by hand without requiring lifting devices.

### 4 - OPERATION AND USE

#### 4.1 - OPERATION

The power supply unit must be placed on a surface next to the operator in such a way that:

- it is in a stable position without risk of falling;
- the micromotor connection lines are not subjected to excessive strain - risk of breakage of the wires or dragging of the unit;
- the functional status of the appliance can be checked easily.

Referring to Fig. 1, insert micromotor supports in the lateral sockets (1); connect the micromotor to the socket (2-3) by inserting the plug and then twisting the locking ring without exerting undue force;


Fig. 1

With reference to Fig. 2, connect the footswitch to socket (4) on the rear panel by inserting the plug and twisting the locking ring; Insert the saline solution stand into upper socket (8) and then hang the saline solution bottle by the plastic cradle.


Fig. 2

With reference to Fig. 3 open the side door of peristaltic pump (1) and remove the tube (keeping it safe for future uses with irrigators other than the one supplied); insert the larger section tube (2) (included with the disposable irrigator supplied) into the peristaltic pump, paying attention to:

- the positioning of supports (3), since, as can be observed in the figure, the supports perform a restraining function on the exterior of the pump;
- the connection of the "IN" and "OUT" sides (4) respectively to the saline solution bottle and the surgical handpiece (not supplied).






Fig. 3

If the handpiece is of the dual irrigation type (see Fig. 4), connect the double outlet union (Fig. 5-1) to the irrigator supplied; in the event of single irrigation types, connect the single outlet union (Fig. 5-2). Fig. 5 shows the two types of unions supplied with the irrigator in the package.

**Fig. 4****Fig. 5****Fig. 5-1****Fig. 5-2**

Finally, plug the power plug into the mains socket outlet, connecting the appliance side (see Fig. 2) to socket (5) located on the rear of the unit; the same figure shows the main power switch (6) and the two fuses (7). At this point attach the required surgical handpiece (not included) to the micromotor.

**Always switch OFF the appliance after use by setting the main power switch to "O".**

## 4.2 - OPERATOR

After selecting the required program the operator starts the appliance by means of voluntary activation of the footswitch thereby assuming direct control of its operation.

**NOTE: For production reasons the 32 preset programs contain standard parameters that are not of specific relevance for the implantologist; the programs must therefore be personalised in relation to individual requirements.**

## 4.3 – INTENDED USE

The implantological micromotor is designed for use in dental and veterinary surgeries for surgical and implantological operations.

The appliance is not designed for any uses other than those specified by the manufacturer.

## 4.4 - SAFETY WARNINGS

- The appliance must be used exclusively by skilled personnel who are aware of the potential hazards;
- Do not use the appliance in the vicinity of inflammable materials or gases;
- Do not open or disassemble the appliance;
- Wear suitable work clothing and observe the safety regulations in force in the workplace;
- Pay attention to the use of flammable liquids or gases such as disinfectants, anaesthetics, solvents, etc.; the appliance is not designed for use in atmospheres subject to explosion hazards so ensure that any such substances have been completely eliminated or evaporated and dispersed before switching on the appliance (ensure there are no residual substances in the area);
- For the above reason, also pay attention to the presence of endogenous gases and materials impregnated with oxygen (e.g. cotton wadding, gauzes, etc.).

## 5 - OPERATOR INSTRUCTIONS

### 5.1 - CONTROLS

The appliance features the following control devices on the front panel (Fig. 6):

- water pump on/off button "H<sub>2</sub>O" (1);
  - "reverse" button (2) serves to reverse the direction of rotation of the tools;
  - the two commands described above are replicated also on the footswitch (Fig. 7), where the green left-hand button (3) controls starting and stopping of the water pump, and the red right-hand button (4) serves to reverse the direction of rotation of the tools;
  - "↑" (5) and "↓" keys (6) to edit the field currently selected on the display (8);
  - "shift" key (7) to select the field to edit on the display (8);
  - display (8) shows the editable parameters on the two outer lines; the parameters (graphically) set and the relative numeric intervals according to the handpiece used;
  - the editable parameters can be selected by pressing the "shift" key or the footswitch
  - the "shift" key can therefore be used (the cursor flashes on the selected field) to edit the values using the "↑" (5) and "↓" (6) keys.
  - the footswitch (Fig. 7) is equipped with a yellow button (9) that serves to change the preset program and vary the editable parameters when used in conjunction with the green (3) or red (4) buttons;
- it is also equipped with a lever (10) that can be moved to select the parameters to be edited or start the motor and control the speed up to the maximum value previously entered in the program.

To change the program or the editable parameters with the footswitch first press yellow button (9) and then

Position lever (10) on the parameter to be edited and set the required value by means of the green (3) and red (4) buttons:

The green button (3) corresponds to ↑,

and the red button (4) corresponds in this function to ↓.

When you press the yellow button a continuous audible signal sounds to indicate that the green (3) and red (4) buttons are active only in programming mode and can therefore not perform their primary function of, respectively, starting/stopping the water pump and reversing the direction of rotation of tools. When you press yellow button (9) again the audible signal stops and the green (3) and red (4) buttons are once more available to perform their primary functions.


Fig. 6


Fig. 7

Figure 8 shows the data that appear on the display as described below:

- 1) "Step" indicates the sequential program number (from 00 to 31);
- 2) "Torque" (Ncm) - "xxx": describes the torque (power) available at the tool tip;
- 3) "Pump": indicates the quantity of water delivered;
- 4) indicates the direction of rotation of the micromotor selected;
- 5) "MM1/MM2": indicates which of the two micromotors is selected (MM1 or MM2);
- 6) "PUMP\_\_\_\_\_%": graphic display of the quantity of water delivered in percentage;
- 7) "TORQUE\_\_\_\_\_%": graphic display of the percentage value of the torque (power) available at the tool tip;
- 8) "xx ÷ xxx": displays the torque (power) numeric interval that can be used with the handpiece selected;
- 9) "SPEED xxx : indicates the numeric value of the speed at the tool tip;
- 10) "\_\_\_\_\_%": graphic display of the percentage value of the tool tip speed;
- 11) "xx ÷ xxx" : displays the numeric interval of the speed that can be used with the handpiece selected;

- 12) "xx/x": transmission ratio between motor and speed reducer or speed multiplier (this parameter is set in accordance with the handpiece installed); reducers are identified by the first digit, multipliers by the second: e.g. 32/1 indicates that speed reducer 32 has been selected, while 1/5 corresponds to speed multiplier 5; when the reducers or multipliers are changed the torque (2 and 8) and speeds (11 and 14) values are updated automatically on the display since these values depend on the currently selected reducer or multiplier;
- 13) "SP-SW%": graphic display of footswitch percentage speed
- 14) "Speed (rpm) "xxxxx": shows the speed set
  - for the motor if no reducers or multipliers are selected;
  - for the spindle if a reducer or multiplier has been selected.


Fig. 8

## **5.2 - INSTRUCTIONS FOR USE**

### **5.2.1 - Setting into operation**

After performing the installation procedures the appliance can be set into operation:

- 1) equip the micromotor with the required handpiece and tool;
- 2) switch on the appliance by means of the main power switch (6, Fig. 2);
- 3) select the program number by positioning the cursor using the "shift" key (7, Fig. 6) in field 1 shown in figure 8, and then using the "↑" (5 in Fig. 6) and "↓" (6 in Fig. 6) keys until obtaining the program required;
- 4) if there is no program saved in the memory that is suitable in relation to the operation you are about to perform, use the edit procedure to change one of the existing programs appropriately (see § 5.2.3);
- 5) push the footswitch lever (10 in Fig. 7) to start the motor and change the speed; use the previously described control devices on the footswitch (see § 5.1) to start or stop the water pump and reverse the direction of run of the motor as required.

### **5.2.2 - Normal deactivation and emergency deactivation procedures**

The appliance can be deactivated as follows:

- by releasing the footswitch (Fig. 7);
- by setting the main power switch to OFF (6 in Fig. 2);
- by unplugging the power cable from the electrical socket outlet.

### **5.2.3 – Editing a program**

The programs can be edited by using the commands on the front of the machine, as well as by using the footswitch.

To edit a program using the commands on the front (see Figs. 6 and 8): position the cursor on the "Step" field (1 in Fig. 8) and use the "↑" (5 in Fig. 6) and "↓" (6 in Fig. 6) keys to select the number of the first free or unused program.

Once you have selected a program to edit, move around the fields using the "shift" key (7 in Fig. 6) and change the relative parameter values with the "↑" (5 in Fig. 6) and "↓" (6 in Fig. 6) keys; when you have finished the edit procedure to save the program return to the "program n°" field and then quit the program by going to the next or the previous program by pressing "↑" (5 in Fig. 6) and "↓" (6 in Fig. 6). If the appliance is switched off before this final step the changes to the parameter values will not be saved in the memory and the new program settings will be lost.

To edit a program using the footswitch (see Fig. 7 and 8), press the yellow button (9 in Fig. 7), bring the flashing cursor on the "Step" field (1 in Fig. 8) by means of lever (10 in Fig. 7) and select the first free or unused program using the green (3 in Fig. 7) and red (4 in Fig. 7) buttons.

Once you have found the program to be edited, move around the field using lever (10 in Fig. 7), editing the parameter values using the green (3 in Fig. 7) and red (4 in Fig. 7) buttons; when the setting is complete, to save the program, return to the "Step" field (1 in Fig. 8) and quit the program, by moving to another program by pressing the green (3 in Fig. 7) and red (4 in Fig. 7) buttons. If the appliance is switched off before this final step the changes to the parameter values will not be saved in the memory and the new program settings will be lost.

A maximum of 32 programs can be saved in memory.

**Note:** The parameter on which the cursor is positioned continues to flash. ALWAYS SWITCH OFF THE APPLIANCE WHEN YOU FINISH WORKING WITH IT AND FOLLOWING ALL STOPPAGES, IRRESPECTIVE OF THE MOTIVE. IF YOU NOTICE MALFUNCTIONS, STOP THE APPLIANCE IMMEDIATELY AND DISCONNECT IT FROM THE ELECTRICAL POWER SUPPLY. BEFORE RECONNECTING THE POWER SUPPLY FIND THE CAUSE OF THE FAULT AND REMEDY IT. IN THE CASE OF PROBLEMS, CONSULT CARLO DE GIORGI.

## 6 - MAINTENANCE

Periodically check the condition and operation of the accessories and parts subject to wear and tear; in particular check the condition of the burs and the micromotors in order to identify any possible damage.

Periodically check the condition of cable insulation. The appliance cables must be replaced exclusively by suitably skilled and authorised personnel.

## 6.1 - CLEANING

The appliance does not require any particular maintenance operations, with the exception of regular cleaning of the working areas: before cleaning the appliance always disconnect it from the power supply.

The micromotor supports, stand for saline solution bottle and footswitch IPX7 supplied with the appliance can be cleaned with alcohol;

Any flammable substances utilized for cleaning or disinfecting, or any adhesive solvents, must be allowed to evaporate completely before the appliance can be used.

If you use the tube supplied with the peristaltic pump, flush it with clean water to remove any traces of sodium chloride before stowing the appliance away after use.

## 6.2 - STERILIZATION


The tube supplied in the peristaltic pump can be cleaned with alcohol or non-toxic sterilizing liquids.

Particular attention must be paid to cleaning and sterilizing the micromotor complete with hose and support.

The micromotor can be sterilized in an autoclave at 135 °C, while the support can be sterilized by autoclaving at 123 °C.

### **IMPORTANT:**

- before placing the micromotor in the autoclave fit the specially designed endcap in the front of the motor (Fig. 9-1 and Fig. 9-2):

**Fig. 9-1**

**Fig. 9-2**

- do not sterilize the micromotor with liquids that conduct electricity.

## 7 - DIAGNOSTICS

If the appliance fails to operate first check the fuses at the rear (7, Fig. 2). If the fuses must be changed refer to the specifications in the table in chapter 2 of this manual. Before checking the fuses, set the main switch to 0 (OFF) and unplug the mains power cable; check the condition of the fuses and replace if necessary. Once the fuse caps have been refitted you can reconnect the mains power cable.

**If malfunctions and faulty operation persist, unplug the appliance and consult CARLO DE GIORGI.**

### **WARNING!**

**ALL WORK ON THE APPLIANCE MUST BE PERFORMED WITH  
THE POWER SUPPLY DISCONNECTED**

## 8 - ACCESSORIES

### 8.1 - MOTOR


Autoclavable motor complete with support

### 8.2 – TUBE FOR IRRIGATION

Alternative tube for the connection between irrigator and handpiece (different with respect to the tube supplied as standard) (See Fig. 10 - A).

### 8.3 - UNION

Alternative union for the connection between irrigator and handpiece (different with respect to the union supplied as standard) (See Fig. 10 - B).

**Fig. 10**

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