Medtronic

XPS® Bone Bridge™ Drill Bit

Mèche XPS® Bone Bridge™

XPS® Bone Bridge™ punta da trapano

XPS® Bone Bridge™ Bohraufsatz

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Product Information and Instructions

FR Informations et instructions concernant le produit **IT** Informazioni sul prodotto ed istruzioni per l'uso **DE** Produktinformationen und Gebrauchsanweisung **ES** Instrucciones e información sobre el producto

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EN**Device description**

The XPS Bone Bridge Drill Bit is an interchangeable drill bit intended for use with the Medtronic Xomed XPS System and with Bone Bridge Drill Guide. The reusable drill is packaged singly and provided sterile for first time convenience.

Indications for use

The XPS Bone Bridge Drill Bit is indicated for use in incising or removing bone and tissues in creating a bone tunnel fixation point for suture in a browlift surgery.

Contraindications

There are no known contraindications.

Warnings

- Do not operate drill continuously. Use irrigation. Drill is intended for use with the 3747101 Bone Bridge Drill Guide.
- Excessive noise from the drill or eccentricity of the drill may cause hearing damage when drilling close to the cochlea or ossicular chain.
- The use of a drill may cause an inordinate amount of heat buildup resulting in thermal injury to tissue. This heat buildup may be accentuated by eccentricity of the drill.
- Excessive force may lead to drill bending, breakage or excess vibration. If the drill fractures during use, the resultant rough surfaces may cause lacerations of the body tissue. A drill which fractures inside a patient may also cause tissue injury and necessitate emergency procedures for extraction.
- Eccentricity of the drill can cause drill vibration and may result in excess tissue and bone destruction.
- Always keep the cutting tip of the drill away from fingers and loose clothing.
- Bending or prying may cause the drill to break causing harm to patient or staff.

Precautions

- Confirm depth of penetration prior to use. CT Scan is an important reference for use to avoid penetration.
- Always ensure that the drill is securely engaged into the drill handpiece prior to operating the system.
- Always examine operation of each drill in a handpiece before use.

Instructions for use

Refer to the XPS System User's Guide for detailed set up instructions of the system.

- Apply irrigation to both Drill Guide and drill.
- Operate drill in forward briefly, remove and irrigate again for reduced risk of thermal injury.

Reprocessing Instructions for Reusable Instruments

When applying dry times to Medtronic instrument trays and their accessories, dry times outside the standard healthcare prevacuum parameters may be required. The current recommended dry times may range from a standard 15 minutes to an extended time of 40 minutes. This is especially important for polymer-based (plastic) instrument trays. The user should have adequate means to verify drying time, such as visual inspection.

Reprocessing Instructions (Per ISO17664)

Cleaning and Sterilization instructions are subject to change without notice. Up to date instructions are available online at manuals.medtronic.com.
68E3979 B

Reuse

The XPS Bone Bridge Drill Bit is provided sterile for first time convenience. Maintaining a quantity of sterile provided drills is recommended. All drills must be cleaned and sterilized using the following process before each reuse.

Warnings and Precautions

- Before use, carefully inspect the device for any irregularities or eccentricities. Damaged or worn product should be replaced.
- Do not cold soak in glutaraldehyde, chlorine, or ammonium solutions, or dry sterilize, as damage to the instrument finish may occur.

Limitations

- Discard any device that shows signs of damage or wear.
- After cleaning and sterilization, verify functionality prior to re-use.

Point of Use

- This product must be cleaned and sterilized before any reuse.
- Rinse instrument immediately after use in distilled or demineralized water. Keep instrument moist (for example, cover with a wet drape or saturate with tap water - do not use saline) after use and prior to decontamination, cleaning and sterilization to ensure adequate cleaning.
- Instruments with cleaning ports and/or suction tubes shall be flushed with neutral pH enzymatic detergent, followed by a tap water rinse.

Containment and Transportation	It is recommended that devices are reprocessed as soon as is practical following use.				
Preparation for decontamination	<ul style="list-style-type: none"> Promptly and thoroughly rinse with deionized water after each use. Fully disassemble modular designed instruments for effective cleaning. Place instruments so they do not touch each other. <p>Note: Flushing in running water is essential between decontamination and cleaning to prevent any risk of reaction between the two solutions.</p>				
Cleaning: Automated	<ul style="list-style-type: none"> All instruments shall be thoroughly rinsed with tap water prior to placement in the automatic washer. Remove instruments and equipment from any sterilization trays before placing into washer baskets. Orient devices following recommendations of washer/disinfector manufacturers. Use mildly alkaline or neutral pH detergent (pH 7.0 - 8.5) recommended by washer/disinfector or detergent manufacturers. These products have been validated for effective cleaning using an automatic washer/disinfector cycle consisting of the following: 				
Phase	Recirculation Time	Water temperature	Detergent type		
Pre-wash 1	3 minutes	cold tap water	not applicable		
Enzyme Wash	5 minutes	hot tap water	neutral pH enzymatic detergent		
Wash 1	5 minutes	66°C set point	neutral pH detergent		
Rinse 1	3 minutes	hot tap water	not applicable		
Thoroughly examine instruments for any residual soil.					
Cleaning: Manual	<ul style="list-style-type: none"> Soak in lukewarm (less than 43°C), mild (pH 7.0 - 8.5), enzymatic detergent and deionized water for a minimum of 2 minutes. Use a soft instrument brush to scrub instruments while submerged in cleaning solution to remove organic matter. Rinse with deionized water, then clean ultrasonically in a lukewarm (less than 43°C), mild (pH 7.0 - 8.5), detergent and deionized water for 10 minutes. To prevent stains, use distilled or demineralized water as a final rinse prior to sterilization. 				
Disinfection	Note: Do not cold soak in glutaraldehyde, chlorine, or ammonium solutions, or dry heat sterilize, as damage to the instrument finish may occur.				
Packaging	<ul style="list-style-type: none"> A standard sterilization wrap may be used. In the US, an FDA approved surgical wrap is required. In sets: Instruments may be loaded into dedicated Medtronic Xomed instrument trays. Ensure that cutting edges are protected. Wrap trays using appropriate method. 				
Sterilization (Temperatures are minimum required; times are minimum required)	<ul style="list-style-type: none"> Check the cleanliness and operation of the instrument. Clean again if debris is present and remove from use any damaged instrument. Close instruments with catches and racks on the first notch. If applicable, arrange the instruments in the dedicated instrument tray. Follow the appropriate cycle listed in the table below. All steam cycles have been validated in the wrapped configuration and can be sterilized wrapped or unwrapped. These devices have only been validated for steam sterilization methods. 				
	Instruments Only		Instruments Only or Instruments in Medtronic Xomed Instrument Trays		
Cycle	Gravity	Gravity	Pre-vac	Pre-vac (FR/WHO)	Pre-vac (UK)
Temperature	121°C	132°C	132°C	134°C	134°C
Time	30 minutes	10 minutes	4 minutes	18 minutes	3 minutes
Drying	15 – 40 minutes, or until visibly dry				
Maintenance, Inspection and Testing	Inspect components for any damage before and after each use. If damage is observed do not use the instrument until it is repaired or replaced. After cleaning and sterilization, verify functionality prior to re-use.				

Storage	Store instruments in a clean, dry area.
Additional Information	None
<p>Note: The instructions provided above have been validated by the manufacturer as being CAPABLE of preparing the product for re-use. It remains the responsibility of the processor to ensure that the reprocessing is performed using validated equipment to achieve the desired result. This normally requires validation and routine monitoring of the process. Some devices have specific assembly instructions. In this case, refer to the assembly insert provided with the device for additional instructions.</p> <p>All validations performed per current AAMI TIR12, "Designing, testing, and labeling reusable medical devices for reprocessing in health care facilities: A guide for medical device manufacturers."</p> <p>Medtronic recommends incineration of devices that have directly contacted patients suspected or confirmed with Transmissible Spongiform Encephalopathy (TSE)/ CJD diagnosis. NHS Estates HTM 2010 Parts 4 & 6: Appendix 2, Items contaminated with TSE Agents and WHO Infection Control Guidelines for Transmissible Spongiform Encephalopathies refers to a TSE decontamination cycle using a steam autoclave at a temperature of 134-137°C for a single cycle of 18 minutes or repeated for a total of six 3-minute cycles.</p> <p>Note: In France, soaking the device in one mole of Soda (NaOH) per liter of solution is the recommendation of French Instruction DGS/RI3/2011/449 in order to prevent Creutzfeldt-Jakob Disease transmission.</p>	

Customer service information

For further information regarding the use of this product or to report any problems, please contact Medtronic Xomed using the appropriate information provided on the blue and white contact information card packaged with each device; or contact your local distributor.

Limited warranty

- A. This LIMITED WARRANTY provides assurance for the customer who purchases a Medtronic Xomed Product (hereinafter the "Product") that should the Product fail to function to Medtronic Xomed's published specifications during the term of this LIMITED WARRANTY (one year from the date of shipment for new Product, 90 days from date of shipment for refurbished or used Product), Medtronic Xomed will either replace, repair, or issue a credit (adjusted to reflect the age of the Product) for the Product or any portion thereof. This LIMITED WARRANTY is extended only to the buyer purchasing the Product directly from Medtronic Xomed or from its affiliate or its authorized distributor or representative.
- B. To qualify for this LIMITED WARRANTY, the following conditions must be met:
1. The Product must be used on or before its "Use By" or "Use Before" date, if applicable.
 2. The Product must be used in accordance with its labeling and may not be altered or subjected to misuse, abuse, accident or improper handling.
 3. Medtronic Xomed must be notified in writing within thirty (30) days following discovery of a defect.
 4. The Product must be returned to Medtronic Xomed within thirty (30) days of Medtronic Xomed receiving notice as provided for in (3) above.
 5. Upon examination of the Product by Medtronic Xomed, Medtronic Xomed shall have determined that: (i) the Product was not repaired or altered by anyone other than Medtronic Xomed or its authorized representative, (ii) the Product was not operated under conditions other than normal use, and (iii) the prescribed periodic maintenance and services have been performed on the Product.
- C. This LIMITED WARRANTY is limited to its express terms. THIS LIMITED WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED WHETHER STATUTORY OR OTHERWISE, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. In no event shall Medtronic Xomed be liable for any consequential, incidental, prospective or other similar damage resulting from a defect, failure, or malfunction of the Product, whether a claim for such damage is based upon the warranty, contract, negligence or otherwise.
- D. The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. Users may benefit from statutory warranty rights under legislation governing the sale of consumer goods. If any part or term of this LIMITED WARRANTY is held by any court of competent jurisdiction to be illegal, unenforceable, or in conflict with applicable law, the validity of the remaining portion of the LIMITED WARRANTY shall not be affected, and all rights and obligations shall be construed and enforced as if this LIMITED WARRANTY did not contain the particular part or term held to be invalid.

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© 2014 Medtronic, Inc.
Medtronic Xomed, Inc.
Medtronic B.V.