Instructions for Use

Part 1

Instructions for Use – GB
Gebrauchsanweisung – DE
Instructions d'utilisation – FR
Bruksanvisning – SE
Istruzioni per l'uso – IT
Οδηγίες Χρήσης – GR
Manual de instrucciones – ES
Brugsanvisning – DK

PC-based Audiometer Equinox²

***Valid from Equinox 2.0 SoftwareSuite version 2.00
80702307 – 05/2012***

Instructions for Use - GB

PC-based Audiometer Equinox²

Table of Contents

1.	INTRODUCTION	1
1.1	About this Manual	1
1.2	Intended Use.....	1
1.3	Product Description.....	1
1.4	Warnings.....	3
2.	UNPACKING AND INSTALLATION	5
2.1	Unpacking and Inspection.....	5
2.2	Marking	6
2.3	Hardware installation	6
2.4	Connection Panel Dictionary.....	8
2.5	Software installation.....	9
2.6	Driver Installation	20
2.7	How to Install a Shortcut in the Noah3 Toolbar for direct access to the Software Suite	23
2.8	License.....	23
3.	OPERATING INSTRUCTIONS	25
3.1	Using the Tone Screen	26
3.2	Using the Speech Screen	31
3.2.1	Speech Audiometry in Graph Mode	33
3.2.2	Speech Audiometry in Table Mode	34
3.3	Using the Print Wizard	36
4.	MAINTENANCE	39
4.1	General Maintenance Procedures	39
4.2	How to clean Interacoustics Products.....	39
4.3	Concerning Repair	40
4.4	Warranty.....	40
5.	GENERAL TECHNICAL SPECIFICATIONS	43
5.1	Reference Equivalent Threshold Values for transducers	44
5.2	Maximum hearing level settings provided at each test frequency	44
5.3	Pin Assignments	44
5.4	Electromagnetic Compatibility (EMC)	44

1. Introduction

1.2 About this Manual

This manual is valid for the Equinox². These products are manufactured by:

Interacoustics A/S

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DK 5610 Assens

Denmark

Tel.: +45 6371 3555

Fax: +45 6371 3522

E-mail: info@interacoustics.com

Web: www.interacoustics.com

1.2 Intended Use

The audiometry software, AC440, for the Equinox² is a fully NOAH compatible air, bone, speech audiometer that can be operated either from a PC keyboard or a dedicated audiometer keyboard.

The AC440 is a device for diagnosing a hearing loss. Output and specificity of this type of device are based on the test characteristics defined by the user, and may vary depending on environmental and operating conditions. The diagnosing of hearing loss using this kind of audiometer depends on the interaction with the patient. However, for patients not responding well possibilities of various tests allow the tester of having at least some evaluative result. Thus, a "normal hearing" result should not allow for ignoring other contra indications in this case. A full audiologic evaluation should be administered if concerns about hearing sensitivity persist.

It is intended to be used by hearing aid dispensers and hearing healthcare professionals.

1.3 Product Description

Equinox² is a PC based audiometer with integrated audiologic software modules on a PC.

The systems consist of the following included and optional parts:

AC440
Included parts: Affinity2.0 AC440 CD OtoAccess™ database CD TDH39 Audiometric headset MTH400 Headset EMS400 Talk back microphone B71 Bone conductor APS3 Patient response button Standard USB cable Power cable 120 or 230V Mouse pad Operation and multilingual CE manual
Optional parts: DD45 Audiometric headset DAK70 Audiometer keyboard with live voice mic. EARTONE 3A Audiometric headset (5As may be

substituted)

ACC60 Affinity2 carrying case

CIR22 Insert masking

earphones

Audiocup enclosures

Peltor noise excluding headset

HDA200 Audiometric headset

HDA280 Audiometric headset

KOSS R80 high frequency

headset

AP12 Power amplifier 2x12

Watt

AP70 Power amplifier 2x70

Watt

ALS7 Loudspeaker

AFC8 Sound cabin installation

panel

UCO15 Optical USB extension

cable

Optional special tests:

High Frequency audiometry

(HF440)

Masking Level Difference

(MLD440)

Multi Frequency module

(MF440)

Speech from hard-drive

(SFH440)

SISI test

Master Hearing Aid (MHA440),

Hearing Loss Simulator

(HLS440)

Loudness Scaling (LS440)

1.4 Warnings

Throughout this manual the following meaning of warnings, cautions and notices are used:

	WARNING indicates a hazardous situation which, if not avoided, could result in death or serious injury.
	CAUTION , used with the safety alert symbol, indicates a hazardous situation which, if not avoided, could result in minor or moderate injury.
	NOTICE is used to address practices not related to personal injury.

2. Unpacking and Installation

2.2 Unpacking and Inspection

Check box and contents for damage

When the instrument is received please check the shipping box for rough handling and damage. If the box is damaged it should be kept until the contents of the shipment have been checked mechanically and electrically. If the instrument is faulty please contact your local distributor. Keep the shipping material for the carrier's inspection and insurance claim.

Keep carton for future shipment

Equinox 2 comes in its own shipping carton, which is specially designed for the Equinox 2. Please keep this carton. It will be needed if the instrument has to be returned for service. If service is required please contact your local distributor.

Reporting Imperfections

Inspect before connection

Prior to connecting the product it should once more be inspected for damage. All of the cabinet and the accessories should be checked visually for scratches and missing parts.

Report immediately any faults

Any missing part or malfunction should be reported immediately to the supplier of the instrument together with the invoice, serial number, and a detailed report of the problem. In the back of this manual you will find a "Return Report" where you can describe the problem.

Please use "Return Report"

Please realise that if the service engineer does not know what problem to look for he may not find it, so using the Return Report will be of great help to us and is your best guarantee that the correction of the problem will be to your satisfaction.


Storage

If you need to store the Equinox 2 for a period, please ensure it is stored under the following conditions:

Temperature;	0-50°C
Relative Humidity:	10-95% Non condensing

2.2 Marking

The following marking can be found on the instrument:

Symbol	Explanation
	Type B applied parts. Patient applied parts that are not conductive and can be immediately released from the patient.
 or 	Refer to instruction manual
	WEEE (EU-directive) This symbol indicates that when the end-user wishes to discard this product, it must be sent to separate collection facilities for recovery and recycling. Failing to do so may endanger the environment.
	The CE-mark indicates that Interacoustics A/S meets the requirements of Annex II of the Medical Device Directive 93/42/EEC. TÜV Product Service, Identification No. 0123, has approved the quality system.
	Year of manufacture

2.3 Hardware installation

NOTICE

- DO NOT connect the Equinox² hardware to the computer before the software has been installed!
- The instrument is intended for continued use. However, there is a risk of damaging the transducers if they are operated at the highest intensities for longer periods.
- To avoid overheating, please observe that cooling vents are unobstructed and that the air can flow freely around the instrument.
- Within the European Union it is illegal to dispose of electrical and electronic waste as unsorted municipal waste. Electrical and electronic waste may contain hazardous substances and therefore have to be disposed of separately. Such products will be marked with the crossed-out wheeled-bin image shown to the left. User cooperation is important in order to ensure a high level of reuse and recycling of electrical and electronic waste. Failure to recycle such waste products in an appropriate way may endanger the environment and consequently the health of human beings.

No installation precautions are necessary to avoid unwanted sound radiation from the audiometer.


When connecting the Equinox² to the mains supply and to a computer the following warnings must be observed:

1. The Equinox² is class I medical product and they must only be connected to a power outlet with a functional protective ground.
2. When connecting this device to a computer it is essential to insert a galvanic separation between the device and the computer, unless the computer is battery operated or supplied by a medical approved power supply.
3. If a mains powered computer is used without galvanic separation, it must be supplied from a medical isolation transformer fulfilling the requirements of IEC/ES 60601-1.
4. Although the instrument fulfils the relevant EMC requirements precautions should be taken to avoid unnecessary exposure to electromagnetic fields, e.g. from mobile phones etc. If the device is used adjacent to other equipment it must be observed that no mutual disturbance appears, e.g. unwanted noise in headphones. If disturbances are observed, try to separate the Equinox² from the disturbing device.
5. If this apparatus is connected to one or more other devices with medical CE marking, to make up a system or pack, the CE marking is only valid also for the combination if the supplier has issued a declaration stating that the requirements in the Medical Device Directive article 12 are fulfilled for the combination.
6. If connection is made to standard equipment like an active speaker, special precautions must be taken in order to maintain medical safety. If used without galvanic separation, it must be supplied from a medical isolation transformer fulfilling the requirements of IEC/ES 60601-1.
7. Do not disassemble or modify the product as this may impact the safety and/or performance of the device.

Headphones, patient response switches and other accessories shall be plugged into the appropriate connectors as indicated on the rear of the instrument and the overview below:

2.4 Connection Panel Dictionary

Position:	Symbol:	Function:
1	FF1	Connection of FF1
2	FF2	Connection of FF2
3	Left	Plug for left AC headphone
4	Right	Plug for right AC headphone
5	Ins. Left	Plug for left insert earphone
6	Ins. Right	Plug for right insert earphone
7	Bone	Plug for bone conductor
8	Ins. Mask.	Plug for insert earphone for masking
9	HF/HLS	Plug for High Frequency headphone/Hearing Loss Simulator
10	Talk Back	Plug for talk back microphone
11	Mic. 1/TF	Plug for microphone / talk forward
12	Mic. 2	Plug for microphone
13	Ass. Mon.	Plug for assistant's headset
14	Monitor	Plug for monitor headset
15	Pat. Resp. L	Plug for left patient response button
16	Pat. Resp. R	Plug for right patient response button
17	Inp. Aux. 1	Plug for input aux. 1
18	Inp. Aux. 2	Plug for input aux. 2
19	Batt. Sim.	Plug for battery simulator
20	TB Lsp.	Plug for test box loudspeaker
21	TB Loop	Plug for test box loop
22	FF Loop	Plug for free field loop
23	TB Coupler	Plug for test box coupler
24	TB Ref.	Plug for test box reference microphone
25		Box for couplers
26		Ground
27	Sp. 1-4 Power Out	Plug for loudspeaker 1-4 power out
28	FF1	Connection of power amplifier FF1
29	FF2	Connection of power amplifier FF2
30	Sp 1	Connection of loudspeaker 1
31	Sp 2	Connection of loudspeaker 2
32	Sp 3	Connection of loudspeaker 3
33	Sp 4	Connection of loudspeaker 4
34	CD1	Input plug for CD 1
35	CD2	Input plug for CD 2
36	Insitu L.	Connection of Insitu Headset Left
37	Insitu R.	Connection of Insitu Headset Right
38	Keyb.	Connection of Keyboard
39	DC	Plug for power supply for optical USB extension cable
40	USB/PC	Plug for USB cable or PC
41	USB	Plug for USB cable
42	Ref.Mic./Ext.	Connection of external reference microphone
43	Coupler/Ext.	Connection of external coupler
44	Mains	Plug for mains cable
45	Power	Turns power on/off.

2.5 Software installation

To know before you Start Installation

You must have administrative rights to the computer on which you are installing the Affinity²/ Equinox² Suite.

NOTICE

1. DO NOT connect the Equinox² hardware to the computer before the software has been installed!
2. Interacoustics will not make any guaranty to the functionality of the system if any other software is installed, with exception of the Interacoustics measurement (AC440/REM440) modules and AuditBase system4, Otoaccess™ or a Noah3.7. compatible Office Systems or later releases.

What you will need:

1. Equinox² Suite Installation DVD
2. USB Cable
3. Equinox 2 Hardware

Supported Noah3 Office Systems

- Noah 3.7 (from HIMSA)
- AuditBase System 4
- Mirage
- Practice Navigator
- Power Office
- AkuWin
- He@r-O
- Entendre 2

To use the software in conjunction with a database (e.g. Noah3 or OtoAccess™), make sure the database is installed prior to the Equinox² Suite installation. Follow the manufacturer's installation instructions provided to install the relevant database.

Note that if you are using AuditBase System 4 you must ensure to launch this office system prior to this installation of the Equinox² Suite.

Installation on various Windows® versions


Installation on Windows® XP (SP2 or later), Windows Vista and Windows 7 (32 and 64 bit) systems are supported.

Note that with Windows® "N" versions, Media Player needs to be installed manually for the software suite to run.


Software Installation on Windows7

Insert the installation DVD and follow the steps below to install the Equinox² Suite software. If the installation procedure does not start automatically, click "Start", then go to "My Computer" and double click the DVD/CD-RW drive to view the contents of the installation DVD. Double click the "setup.exe" file to initiate the installation.

1. If the dialog below is shown, click Accept.


2. If the dialog below is shown, click Install.


3. Wait for the dialog shown below to appear, then click "Next".


4. Click "Next", to start the installation.


5. Select whether the software is to be used with Affinity or Equinox hardware and click "Next"

6. Choose an installation folder and select whether the Affinity/Equinox Suite should be installed for you exclusively or for anyone else using the PC. Then press "Next"

7. If the dialog box below appears, select "Install this driver software anyway".


8. If the dialog below is shown, click Install (this will install the driver for VIOT).

9. If the dialog box below appears, select "Install this driver software anyway".


10. When the installation is complete, the dialog box below is shown. Click "Close" to finish the installation. The Affinity/Equinox Suite is now installed.


Software Installation on Windows XP

Insert the installation DVD and follow the steps below to install the Equinox² Suite software. If the installation procedure does not start automatically, click “Start”, then go to “My Computer” and double click the DVD/CD-RW drive to view the contents of the installation DVD. Double click the “setup.exe” file to initiate the installation.


1. If the dialog below is shown, click install.


2. If the dialog below is shown, click install


- Wait for the dialog shown below to appear, then click "Next".


- Click "Next", to start the installation.


5. Select whether the software is to be used with Affinity or Equinox hardware and click "Next"

6. Choose an installation folder and select whether the Affinity/Equinox Suite should be installed for you exclusively or for anyone else using the PC. Then press "Next"


7. Please wait while the Affinity/Equinox Suite installs


8. If the dialog box below appears, select "Continue Anyway".

9. If the dialogs below appear, click “Yes” to confirm file replace.

10. When the installation is complete, the dialog box below is shown. Click “Close” to finish the installation. The Affinity/Equinox Suite is now installed.

2.6 Driver Installation

Now that the Equinox² Suite software is installed, you must install the driver for the hardware.

For Windows 7:


1. Connect the Equinox² hardware to the PC via the USB connection.
2. The system will now automatically detect the hardware and a pop-up will appear on the task bar near the clock indicating that the driver is installed and the hardware is ready for use.

For Windows XP:


1. Connect the Equinox² hardware to the PC via the USB connection.
2. A pop-up indicating "New Hardware Found" will appear near the clock on the task bar.
3. Wait for the dialog shown below to appear. Select "No, not this time" and click "Next".


4. Select "Install the software automatically (Recommended)".

5. Click "Next" and wait while the software installs

- When the installation is complete, the dialog box below is shown. Click "Finish". The Equinox² driver is now installed.

- To finalize the installation process, start the Affinity/Equinox Suite by clicking Start | All Programs | Interacoustics | Suite for Affinity and Equinox | Suite Launcher. Select the desired Regional settings and Language when the pop-up below appears.

2.7 How to Install a Shortcut in the Noah3 Toolbar for direct access to the Software Suite

If you are using HIMSA'S Noah3 Office System and you want a button shortcut for the installed module, please follow these steps.

1. Start Noah3
2. Go to the menu item "Tools"
3. Select "Module buttons"
4. Choose one of the two available buttons and select the software name in the drop down list.
5. Press "OK" and your module will be launched when the appropriate button gets pushed.

How to install a shortcut for launching without Noah3

If you do not have Noah3 on your computer you can make a shortcut for direct launching of the software suite as a stand-alone module.

You can find the executable file (.exe) for this purpose through the following path start/All Programs/Interacoustics/Suite for Affinity/Equinox Suite Launcher. Right click on the last item with the mouse and select "Create shortcut".

However, you will not be able to save your recordings when using this way of working.

Working with OtoAccess™

If you want to launch the software suite from OtoAccess™ this is also possible.

1. Start the database by clicking on the OtoAccess™ shortcut or select the OtoAccess™ from the "Start" menu.
2. Select a patient from the list in the "Client" tab
3. Double click on the Affinity/ Equinox Suite icon under "Select Instrument"

The screenshot shows the OtoAccess™ software interface. The top menu bar includes 'File', 'Edit', 'View', 'Help', and 'Login'. The main window is divided into several sections:

- Client Information:** Fields for First name (Ejvind), Last name (Christensen), Initials, Person ID (Diagnostic Suite), Birth date (19-04-1949), Gender (Male), Phone, Mobile, Fax, Email, Address, Address2, Zipcode, City, Country, State, and County. There is a checkbox for 'With schedule today'.
- Select Instrument:** A list of instruments including Affinity Suite, Diagnostic Suite, Titan Suite, Callisto Suite, and VIOT Suite.
- Categories:** Fields for Category 1 and Category 2.
- Table:** A table with columns: Person ID, First name, Last name, Telephone, Mobile, and City. The table contains several rows of data, including 'Titan', 'Diagnostic Suite', 'Affinity', 'Callisto', and 'Affinity (REM 2.0.5)'.

The status bar at the bottom shows 'Login: admin', 'Created: 07-05-2010 10:51:34', and 'Modified: 16-05-2011 08:17:18'.

For further instructions about working with the database, please see the operation manual for OtoAccess™.

2.8 License

When you receive the product it already contains the licenses to access the ordered software modules. If you would like to add additional modules, please contact your dealer.

3. Operating instructions

The instrument is switched by the switch on the rear. When operating the instrument, please observe the following general precautions:

CAUTION


1. The intended operators of the instrument are ENT doctors, Audiologists and other professionals with similar knowledge. Using the instrument without adequate knowledge may lead to erroneous results and may endanger the patients hearing.
2. Only recorded speech material with a stated relationship with the calibration signal should be used. In the calibration of the instrument it is assumed that the calibration signal level is equal to the average level for the speech material. If this is not the case, the calibration of the sound pressure levels will be invalid and the instrument needs recalibration.

It is recommended that the disposable foam eartips supplied with the optional E'A'R Tone 3A or E'A'R Tone 5A insert transducers are replaced after each client tested. Disposable plugs also insure that sanitary conditions exist for each of your clients, and that periodic cleaning of a headband or cushion is no longer required.

3. The instrument must warm up for at least 3 minutes in room temperature before use.
4. Be sure to use only stimulation intensities, which will be acceptable for the patient.
4. The transducers (headphones, bone conductor, etc.) supplied with the instrument are calibrated to this instrument - exchange of transducers requires a new calibration.
5. It is recommended that parts which are in direct contact with the patient (e.g. earphone cushions) are subjected to standard disinfecting procedure between patients. This includes physically cleaning and use of a recognized disinfectant. Individual manufacturer's instruction should be followed for use of this disinfecting agent to provide an appropriated level of cleanliness.
6. To establish conformity with the IEC 60645-2 standard, it is important that the speech input level is adjusted to 0VU. It is equally important that any free field installation is calibrated at the site where it is used and under the conditions that exists during normal operation.
7. For maximum electrical safety, remove the USB cable when it is left unused.

3.1 Using the Tone Screen

The following section describes the elements of the tone screen.

Menu provides access to File, Edit, View, Tests Setup, and Help (see section **Error! Reference source not found.** for information about the menu items).

Print allows for printing the sessions acquired data (see section **Error! Reference source not found.** for information about the Print Wizard).

Save & New Session saves the current session in Noah3 or OtoAccess™ and opens a new one.

Save & Exit saves the current session in Noah3 or OtoAccess™ and exits the Suite.

Collapse the left side panel.

Go to Tone Audiometry activates the tone screen when in another test.

Go to Speech Audiometry activates the speech screen when in another test.

Extended Range +20 dB extends the testing range and can be activated when the testing dial setting gets within 55 dB of the maximum level of the transducer.


Note that the extended range button will flash when it needs activation for reaching higher intensities.

Fold an area so that it only shows the label or the buttons of that area.

Unfold an area so that all buttons and labels are visible


Show/hide areas can be found by right mouse clicking on one of the areas. The visibility of the different areas as well as the space that they take on the screen is locally saved to the examiner.

List of Defined Protocols allows for selecting a test protocol for the current test session. Please refer to section **Error! Reference source not found.** for more information about protocols. Right mouse click on a protocol allows the current examiner to set or deselect a default startup protocol.


Temporary Setup allows for making temporary changes to the selected protocol. The changes will be valid for the current session only. After making the changes and returning to the main screen, the name of the protocol will be followed by an asterisk (*).

List of historical sessions accesses historical sessions for comparison purposes. The audiogram of the selected session, indicated by the orange background, is shown in colours as defined by the used symbol set. All other audiograms that are selected by check marks show on screen in the colours as indicated by the text colour of the date and time stamp. Note that this listing can be resized by dragging the double lines up or down.


Go to Current Session brings you back to the current session.

High Frequency shows frequencies on the audiogram (up to 20 kHz for the Affinity²/Equinox²). However, you will only be able to test in the frequency range the selected headset is calibrated for.


High Frequency Zoom¹ activates high frequency testing and zooms in on the high frequency range. Please refer to section **Error! Reference source not found.** for more information on high frequency testing.

Single audiogram toggles between viewing the information of both ears in a single graph and two separate graphs.

Multi frequencies² activate testing with frequencies in between the standard audiogram points. The frequency resolution can be adjusted in the AC440 setup.

Synchronize channels locks the two attenuators together. This function may be used to perform synchronous masking.

Edit Mode button activates the editing function. Left clicking on the graph will add/move a point to the position of the cursor. By right clicking on a point it is possible to **Delete** the recorded point or the entire curve. Furthermore, right mouse click provides the option to **Add no response**, **Add masked threshold**, **Add**


¹ HF requires an additional license for the AC440. If not purchased, the button is grayed out.

² MF requires additional license for the AC440. If not purchased, the button is grayed out.

masked-no-response threshold and **Hide unmasked thresholds** where masked exist.


Mouse controlled audiometry enables you to do the audiometry using the mouse only. Frequencies and intensities are changed by mouse movements. The stimulus is presented with the left mouse button, and the threshold is stored with the right mouse button.

The **dB step size** button indicates to which dB step size the system is currently set. It rotates from 1 to 2 to 5 dB.


The **hide unmasked threshold** will hide those unmasked thresholds where masked thresholds exist.


The **Patient monitor** opens an always-on-top window with the tone audiograms and all its counselling overlays shown. The size and position of the patient monitor gets saved for each examiner individually.

The **Phonemes** counselling overlay shows phonemes as it is set up in the protocol that is currently in use.

The **Sound examples** counselling overlay shows pictures (png-files) as they are set up in the protocol that is currently in use.


The **Speech banana** counselling overlay shows the speech area as it is set up in the protocol that is currently in use.

The **Severity** counselling overlay shows the degrees of hearing loss as it is set up in the protocol that is currently in use.

The **Max. testable values** shows the area beyond the maximum intensity the system allows. This is a reflection of the transducer calibration and depends on the extended range being activated.


Talk Forward activates the Talk Forward microphone. The arrow keys can be used to set the talk forward level through the currently selected transducers. The level will be accurate when VUE meter indicates to be at zero dB.

Selecting the **Monitor Ch1** and/or **Ch2** check boxes allow you to monitor one or both channels through an external loudspeaker/headset connected to the monitor input. The monitor intensity is adjusted by the arrow keys.


The **Talk back** check box enables you to listen to the patient. Note that you need to be equipped with a microphone connected to the talk back input and an external loudspeaker/headset connected to the monitor input.


Selecting **HL**, **MCL**, **UCL** or **Tinnitus** sets the symbol types that are currently in use by the audiogram. HL stands for hearing level, MCL stands for most comfortable level and UCL stands for uncomfortable level. Note that these buttons show the unmasked right and left symbols of the currently selected symbol set.

Each type of measurement is saved as a separate curve.

Binaural and **Aided** selection is available when free field 1 or free field 2 is selected as output. The function allows for indicating if the test is performed with hearing aids or binaurally with corresponding symbols. The measurements will be saved as separate curves.

In the **Comments** section you can type comments related to any audiometric test. The used space by the comments area can be set by dragging the double line with your mouse. Pressing the **Report editor** button opens a separate window for adding notes to the current session. The report editor and comment box contain the same text. In case the formatting of the text is important, this can only be set within the report editor.

Note that after saving the session no changes can be added to the report. See section **Error! Reference source not found.** for more information about the Report Editor.

Output	Input
Phone right	Tone
Phone left	Warble
Bone right	NB
Bone left	WN
Free field 1	
Free field 2	
Insert right	
Insert left	

The **Output** list for channel 1 provides the option to test through head phones, bone conductor, free field speakers or insert phones. Note that the system only shows the calibrated transducers.

The **Input** list for channel 1 provides the option to select pure tone, warble tone, narrow band noise (NB) and white noise (WN).

Note that the background shading is according to the side that is selected, red for right and blue for left.

Input	Output
Tone	Phone right
Warble	Phone left
NB	Free field 1
WN	Free field 2
TEN	Insert right
	Insert left
	Insert mask
	Off

The **Output** list for channel 2 provides the option to test through head phones, free field speakers, insert phones or insert masking phone. Note that the system only shows the calibrated transducers.

The **Input** list for channel 2 provides the option to select pure tone, warble tone, narrow band noise (NB), white noise (WN) and TEN noise³.

Note that the background shading is according to the side that is selected, red for right, blue for left, and white when off.

Pulsation allows for single and continuous pulsating presentation. The duration of the stimulus can be adjusted in the AC440 setup.

Sim/Alt allows toggling between Simultaneous and Alternate presentation. Ch1 and Ch2 will present the stimulus simultaneously when Sim is selected. When Alt is selected, the stimulus will alternate between Ch1 and Ch2.

Masking indicates if channel 2 is currently in use as a masking channel and in that way makes sure masking symbols are used in the audiogram. For example in paediatric testing through free field speakers, channel 2 can be set as a second testing channel. Note that a separate store function for channel 2 is available when channel 2 is not used for masking.


³ TENs test requires an additional license for the AC440. If not purchased, the stimulus is grayed out.

dB HL Increase and **Decrease** buttons allow in- and decreasing the intensities of channel 1 and 2.

The arrow keys on the PC keyboard can be used for in-/decreasing channel 1 intensities.

PgUp and PgDn on the PC keyboard can be used for in-/decreasing channel 2 intensities.


Stimuli or **attenuator** buttons will light up when the mouse goes over and indicates the presentation of a stimulus.

A right mouse click in the Stimuli area will store a no response threshold. A left mouse click in the Stimuli area will store the threshold at the current position.

Channel 1 stimulation can also be obtained by pressing the space bar or left Ctrl key on the PC keyboard.

Channel 2 stimulation can also be obtained by pressing the right Ctrl key on the PC keyboard.

Mouse movements in the Stimuli area for both channel 1 and channel 2 can be ignored depending on the setup.

Frequency and Intensity display area shows what is currently presented.

To the left the dB HL value for channel 1 is shown and to the right for channel 2. In the centre the frequency is displayed.

Notice that the dB dial setting will flash when trying to go louder than the maximum available intensity.

Frequency increase/decrease increases and decreases the frequency respectively. This can also be obtained using the left and right arrow keys on the PC keyboard.

No visual


Storing thresholds for channel 1 is done by pressing **S** or by a left mouse click in the attenuator of channel 1. Storing a no response threshold can be done by pressing **N** or by a right mouse click in the attenuator of channel 1.

No visual

Storing thresholds for channel 2 is available when channel 2 is not the masking channel. It is done by pressing **<Shift> S** or by a left mouse click in the attenuator of channel 2. Storing a no response threshold can be done by pressing **<Shift> N** or by a right mouse click in the attenuator of channel 2.

The hardware indication picture indicates whether the hardware is connected. **Simulation mode** is indicated when operating the software without hardware.

When opening the Suite the system will automatically search for the hardware. If it does not detect the hardware a dialog box pops up and asks if you want to *continue in simulation mode*.


The **Examiner** indicates the current clinician who is testing the patient. The examiner is saved with a session and can be printed with the results.

For each examiner is logged how the suite is set up with regards to the use of space in the screen. The examiner will find that the suite starts up looking the same as the last time they used the software. Also can an examiner select which protocol must be selected at start up (by right mouse click on the protocol selection list).

3.2 Using the Speech Screen

The following section describes the elements of the speech screen in addition to the tone screen:

Input Levels	
Mic1	28
Mic2	27
CD1	26
CD2	26

Input levels sliders allow for adjusting the input level to 0 VU for the selected input. This ensures that correct calibration is obtained for Mic1, Mic2, CD1, and CD2⁴.

WR1
WR3
WR2

WR1, WR2 and WR3 (Word Recognition) allows selecting different speech list setups as defined by the selected protocol. The labels of these lists which go along with these buttons can also be customized in the protocol setup. (See section **Error! Reference source not found.**).

Output	Input
Phone right	WN
Phone left	Mic 1
Bone right	Mic 2
Bone left	Cd 1
Free field 1	Cd 2
Free field 2	SN
Insert right	Wavefile 1
Insert left	Wavefile 2

The **Output** list for channel 1 provides the option to test through head phones, bone conductor, free field speakers or insert phones. Note that the system only shows the calibrated transducers.

The **Input** list for channel 1 provides the option to select white noise (WN), speech noise (SN), microphone 1 or 2 (Mic1 and Mic2), CD1, CD2 and wave file.

Note that the background shading is according to the side that is selected, red for right and blue for left.

⁴ Mic2 and speech audiometry using a CD player is only available on the Affinity^{2.0}/Equinox^{2.0}.

Equinox² Instruction for Use - English

Date: 2012-03-03

Page 32/44


Input	Output
WN	Phone right
Mic 1	Phone left
Mic 2	Free field 1
Cd 1	Free field 2
Cd 2	Insert right
SN	Insert left
Wavefile 1	Insert mask
Wavefile 2	Off

The **Output** list for channel 2 provides the option to test through head phones, free field speakers, insert phones or insert masking phone. Note that the system only shows the calibrated transducers.


The **Input** list for channel 2 provides the option to select white noise (WN), speech noise (SN), microphone (Mic1 and Mic2), CD1, CD2 and wave file.

Note that the background shading is according to the side that is selected, red for right, blue for left, and white when off.

Speech Scoring:

Phoneme scoring:


- a) **Correct:** A mouse click on this button will store the word as correctly repeated. Quick keys to score correct are the **Up** arrow key and the **B**.
- b) **Incorrect:** A mouse click on this button will store the word as incorrectly repeated. Quick keys to score incorrect are the **Down** arrow key and also **X**, **C** and **V**.
- c) **Store:** A mouse click on this button will store the speech threshold in the speech graph. A point can also be stored by pressing **S**.
- a) **Phoneme scoring:** If phoneme scoring is selected in the AC440 setup, mouse click on the corresponding number to indicate phoneme score. The quick keys to score phonemes are **X**, **C**, **V** and **B** for 0, 1, 2 and 3 respectively.
- b) **Store:** A mouse click on this button will store the speech threshold in the speech graph. A point can also be stored by pressing **S**.

Frequency and Intensity display shows what is currently presented. On the left the dB value for channel 1 is shown and on the right side for channel 2.

In the centre of the current *Speech Score* in % and the *Word Counter* monitors the number of words presented during the test.


3.2.1 Speech Audiometry in Graph Mode

Graph mode presentation settings in the lower left corner and in the presentation options (Ch1 and Ch2) in the upper part of the screen you can adjust the test parameters during the test.

- 1) **The graph:** The curves of the recorded speech graph will be displayed on your screen. The x-axis shows the intensity of the speech signal and the y-axis shows the score in percent. The score is also displayed in the black display in the upper part of the screen together with a word counter.
- 2) **The norm curves** illustrate norm values for **S** (Single syllabic) and **M** (Multi syllabic) speech material respectively. The curves can be edited according to individual preferences in the AC440 setup (see section **Error! Reference source not found.**)
- 3) **The shaded area** illustrates how high in intensity the system will allow. The *Extended Range +20 dB* button can be pressed to go higher. The maximal loudness is determined by the transducer calibration.

3.2.2 Speech Audiometry in Table Mode

The AC440 Table Mode consists of two tables:

- 1) The **SRT** (Speech Reception Threshold) table. When the SRT test is active, it is indicated in orange **SRT**
- 2) The **WR** (Word Recognition) table. When WR1, WR2, or WR3 is active the corresponding label will be orange **WR1**

The SRT table

The SRT table (Speech Reception Threshold table) allows for measuring multiple SRTs using different test parameters, e.g. *Transducer*, *Test Type*, *Intensity*, *Masking*, and *Aided*.

Upon changing *Transducer*, *Masking*, and/or *Aided* and re-testing, an additional SRT entry will appear in the SRT table. This allows for multiple SRT measurements to be shown in the SRT table.

Right		SRT	Left	
SRT	SRT		SRT	SRT
Phone	Phone	Transducer Intensity Masking Test Type Aided Wordlist	Phone	Phone
30	10		10	30
15	15		15	15
HL	HL		HL	HL
	x		x	
Spondee A	Spondee B		Spondee A	Spondee B

The WR Table

The word recognition (WR) table allows for measuring multiple WR scores using different parameters (e.g. *Transducer*, *Test Type*, *Intensity*, *Masking*, and *Aided*).

Upon changing Transducer, Masking, and/or Aided re-testing an additional WR entry will appear in the WR table. This allows for multiple WR measurements to be shown in the WR table.


Right		WR1	WR2	WR3	Left	
WR1	WR1			WR1	WR2	
Phone	FF1	Transducer		Phone	FF2	
55	55	Intensity		55	30	
		Masking				
85	95	Score		90	100	
	x	Aided				
NU-6 LIST 1A	NU-6 LIST 3A	Wordlist		NU-6 LIST 1A	Spondee A	

3.3 Using the Print Wizard

In the Print Wizard you have the option to create customized print templates which can be linked to individual protocols for quick printing. The Print Wizard can be reached in two ways.

- If you want to make a template for general use, or select an existing one for printing: Go to **Menu/ File/Print Layout...** in any of the Affinity²/Equinox² or Callisto Suite tabs (AUD, REM or HIT)
- If you want to make a template or select an existing one to link to a specific protocol: Go to Module tab (AUD, REM, or HIT) relating to the specific protocol and select **Menu/Setup/AC440 setup, Menu/Setup/REM440 setup, or Menu/Setup HIT440 setup**. Select the specific protocol from the drop down menu and select **Print Wizard** at the bottom of the window.


Now the **Print Wizard** window opens and shows the following information and functionalities:

- Underneath **Categories** you can select
 - Templates** to show all available templates
 - Factory defaults** to show only standard templates
 - User defined** to show only custom templates
 - Hidden** to show hidden templates
 - My favorites** to show only templates marked as a favorite
- Available templates from the selected category are shown in the **Templates** viewing area.
- Factory default templates are recognized by the lock icon. They ensure that you always have a standard template and do not need to create a customized one. However, they cannot be edited according to personal preferences without resaving with a new name. **User defined**/created templates can be set to **Read-only** (showing the lock icon), by right clicking on the template and selecting **Read-only** from the drop down list. **Read-only** status can also be removed from **User defined** templates by following the same steps.
- Templates added to **My favorites** are marked with a star. Adding templates to **My favorites** allows quick viewing of your most commonly used templates.

5. The template that is attached to the selected protocol when entering the print wizard via the **AC440** or **REM440** window is recognized by a checkmark.
6. Press the **New Template** button to open a new empty template.
7. Select one of the existing templates and press the **Edit Template** button to modify the selected layout.
8. Select one of the existing templates and press the **Delete Template** button to delete the selected template. You will be prompted to confirm that you want to delete the template.
9. Select one of the existing templates and press the **Hide Template** button to hide the selected template. The template will now be visible only when **Hidden** is selected under **Categories**. To unhide the template, select **Hidden** under **Categories**, right click on the desired template and select **View/Show**.
10. Select one of the existing templates and press the **My Favorites** button to mark the template as a favorite. The template can now be quickly found when **My Favorites** is selected under **Categories**. To remove a template marked with a star from My Favorites, select the template and press the **My Favorites** button.
11. Select one of the templates and press the **Preview** button to print preview the template on screen.
12. Depending how you reached the Print Wizard, you will have the option to press
 - a. **Print** for using the selected template for printing or press
 - b. **Select** for dedicating the selected template to the protocol from which you got into the Print Wizard.
13. To leave the Print Wizard without selecting or changing a template press **Cancel**.

Right clicking on a specific template provides a drop down menu offering an alternative method for performing the options as described above:

4. Maintenance

4.1 General Maintenance Procedures

The performance and safety of the instrument will be kept if the following recommendations for care and maintenance are observed:

- The instrument must go through at least one annual overhaul, to ensure that the acoustical, electrical and mechanical properties are correct. This should be made by an authorised workshop in order to guaranty proper service and repair as Interacoustics provides the necessary circuit diagrams etc. to these repair shops.
- To ensure that the reliability of the instrument is kept, it is recommended that the operator at short intervals, for instance once a day, perform a test on a person with known data. This person could be the operator him/herself.
- After each examination of a patient, it should be ensured that there is no contamination on the parts in connection with the patient. General precautions must be observed in order to avoid that disease from one patient is conducted to others. If ear cushions or ear tips are contaminated, it is strongly recommended to remove them from the transducer before they are cleaned. By frequent cleaning water should be used, but by severe contamination it may be necessary to use a disinfectant. The use of organic solvents and aromatic oils must be avoided.

NOTICE

1. Great care should be exercised by the handling of earphones and other transducers, as mechanical shock may cause change of calibration.

4.2 How to clean Interacoustics Products

If the surface of the instrument or parts of it are contaminated, it can be cleaned using a soft cloth moistened with a mild solution of water and dish washing cleaner or similar. The use of organic solvents and aromatic oils must be avoided. Always disconnect the USB cable during the cleaning process, and be careful that no fluid is entering the inside of the instrument or the accessories.

CAUTION

- Before cleaning always switch off and disconnect from power
- Use a soft cloth lightly dampened with cleaning solution to clean all exposed surfaces
- Do not allow liquid to come in contact with the metal parts inside the earphones / headphones
- Do not autoclave, sterilize or immerse the instrument or accessory in any fluid
- Do not use hard or pointed objects to clean any part of the instrument or accessory
- Do not let parts that have been in contact with fluids dry before cleaning
- Rubber ear-tips or foam ear-tips are single use components
- Ensure isopropyl alcohol does not come into contact with any screens on the instruments

Recommended cleaning and disinfection solutions:

- Warm water with mild, nonabrasive cleaning solution (soap)
- 70% isopropyl alcohol

Procedure:

- Clean the instrument by wiping outer case with a lint free cloth lightly dampened in cleaning solution.
- Clean cushions and patient hand switch and other parts with a lint free cloth lightly dampened in cleaning solution.

- Make sure not to get moisture in the speaker portion of the earphones and similar parts

4.3 Concerning Repair

Interacoustics is only considered to be responsible for the validity of the CE marking, effects on safety, reliability and performance of the equipment if:

1. assembly operations, extensions, readjustments, modifications or repairs are carried out by authorised persons,
2. a 1 year service interval is maintained
3. the electrical installation of the relevant room complies with the appropriate requirements, and
4. the equipment is used by authorised personnel in accordance with the documentation supplied by Interacoustics.

It is important that the customer (agent) fills out the RETURN REPORT every time a problem arises and sends it to Interacoustics, Drejervaenget 8, DK-5610 Assens, Denmark. This should also be done every time an instrument is returned to Interacoustics. (This of course also applies in the unthinkable worst case of death or serious deterioration to patient or user).

4.4 Warranty

INTERACOUSTICS warrants that:

- The Equinox² is free from defects in material and workmanship under normal use and service for a period of 24 months from the date of delivery by Interacoustics to the first purchaser.
- Accessories are free from defects in material and workmanship under normal use and service for a period of ninety (90) days from the date of delivery by Interacoustics to the first purchaser.

If any product requires service during the applicable warranty period, the purchaser should communicate directly with the local Interacoustics service centre to determine the appropriate repair facility. Repair or replacement will be carried out at Interacoustics' expense, subject to the terms of this warranty. The product requiring service should be returned promptly, properly packed, and postage prepaid. Loss or damage in return shipment to Interacoustics shall be at purchaser's risk.

In no event shall Interacoustics be liable for any incidental, indirect or consequential damages in connection with the purchase or use of any Interacoustics product.

This shall apply solely to the original purchaser. This warranty shall not apply to any subsequent owner or holder of the product. Furthermore, this warranty shall not apply to, and Interacoustics shall not be responsible for, any loss arising in connection with the purchase or use of any Interacoustics product that has been:

- repaired by anyone other than an authorized Interacoustics service representative;
- altered in any way so as, in Interacoustics judgement, to affect its stability or reliability;
- subject to misuse or negligence or accident, or which has had the serial or lot number altered, effaced or removed; or
- improperly maintained or used in any manner other than in accordance with the instructions furnished by Interacoustics.

This warranty is in lieu of all other warranties, express or implied, and of all other obligations or liabilities of Interacoustics, and Interacoustics does not give or grant, directly or indirectly, the authority to any representative or other person to assume on behalf of Interacoustics any other liability in connection with the sale of Interacoustics products.

INTERACOUSTICS DISCLAIMS ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING ANY WARRANTY OF MERCHANTABILITY OR FOR FUNCTION OF FITNESS FOR A PARTICULAR PURPOSE OR APPLICATION.

5. General Technical Specifications

Technical Specifications of the AC440 Software

Medical CE-mark:	The CE-mark indicates that Interacoustics A/S meets the requirements of Annex II of the Medical Device Directive 93/42/EEC. Approval of the quality system is made by TÜV – identification no. 0123.	
Audiometer Standards:	Tone: IEC60645-1/ANSI S3.6 Type 1 Speech: IEC60645-2/ANSI S3.6 Type A or A-E	
Transducers & Calibration:	Calibration information and instructions are located in the Service manual. Check the accompanying Appendix for RETSPL levels for transducers	
Air Conduction		
DD45	PTB/DTU report 2009	Headband Static Force 4.5N ±0.5N
TDH39	ISO 389-1 1998, ANSI S3.6-2010	Headband Static Force 4.5N ±0.5N
HDA200	ISO 389-8 2006, ANSI S3.6-2010	Headband Static Force 4.5N ±0.5N
HDA280	PTB report 2004	Headband Static Force 5N ±0.5N
E.A.R Tone 3A/5A	ISO 389-2 1994, ANSI S3.6-2010	
CIR 33	ISO 389-2	
Bone Conduction	Placemenet: Mastoid	
B71	ISO 389-3 1994, ANSI S3.6-2010	Headband Static Force 5.4N ±0.5N
Free Field	ISO 389-7 2005, ANSI S3.6-2010	
High Frequency	ISO 389-5 2004, ANSI S3.6-2010	
Effective masking	ISO 389-4 1994, ANSI S3.6-2010	
Patient Response switch:	Hand held push button.	
Patient communication:	Talk Forward and Talk Back.	
Monitor:	Output through external earphone or speaker.	
Stimuli:		
Tone	125-20000Hz separated in two ranges 125-8000Hz and 8000-20000Hz. Resolution 1/2-1/24 octave.	
Warble Tone	1-10 Hz sine +/- 5% modulation	
Wave file	44100Hz sampling, 16 bits, 2 channels	
Masking	Automatic selection of narrow band noise (or white noise) for tone presentation and speech noise for speech presentation.	
Narrow band noise:	IEC 60645-1:2001, 5/12 Octave filter with the same centre frequency resolution as pure Tone.	
White noise:	80-20000Hz measured with constant bandwidth	
Speech Noise.	IEC 60645-2:1993 125-6000Hz falling 12dB/octave above 1KHz +/-5dB	
Presentation	Manual or Reverse. Single or multiple pulses.	
Intensity	Check the accompanying Appendix for maximum output levels	
Steps	Available Intensity Steps is 1, 2 or 5dB	
Accuracy	Sound pressure levels: ± 2 dB. Vibration force levels: ± 5 dB.	
Extended range function	If not activated, the Air Conduction output will be limited to 20 dB below maximum output.	
Frequency	Range: 125Hz to 8kHz (Optional High Frequency: 8 kHz to 20 kHz) Accuracy: Better than ± 1 %	
Distortion (THD)	Sound pressure levels: below 1.5 % Vibration force levels: below 3 %.	
Signal Indicator(VU)	Time weighting: 350mS Dynamic range: -20dB to +3dB Rectifier characteristics: RMS Selectable inputs are provide with an attnuator by which the level can be adjusted to the indicator reference position(0dB)	
Storing capability:	Tone audiogram: dB HL, MCL, UCL Speech Audiogram: SDS1, SDS2, SDS3, MCL, UCL, Aided, Unaided.	
Compatible Software:	Noah 3.0, OtoAccess™ and XML compatible	

5.1 Reference Equivalent Threshold Values for transducers

See Appendix in English in the bag of the manual.

5.2 Maximum hearing level settings provided at each test frequency

See Appendix in English in the bag of the manual.

5.3 Pin Assignments

See Appendix in English in the bag of the manual.

5.4 Electromagnetic Compatibility (EMC)

See Appendix in English in the bag of the manual.

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You can find the executable file (.exe) for this purpose through the following path start/All Programs/Interacoustics/Suite for / Equinox² Suite Launcher. Right click on the last item with the mouse and select "Create shortcut".

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1. Start the database by clicking on the OtoAccess™ shortcut or select the OtoAccess™ from the "Start" menu.
2. Select a patient from the list in the "Client" tab
3. Double click on the Equinox² Suite icon under "Select Instrument"

Person ID	First name	Last name	Telephone	Mobile	City
0101013333	Demo	Demo	+45 6371 3555		Assens
Titan	Ejvind	Christensen			
Diagnostic Suite	Ejvind	Christensen			
Affinity	Ejvind	Christensen			
Callisto	Ejvind	Christensen			
Callisto 1	Ejvind	Christensen			
Callisto 2	Ejvind	Christensen			
AffinityREM 2.0.5	Ejvind	Christensen			
Affinity/AUD 2.0.5	Ejvind	Christensen			Assens
D51	E	C			
D53	E	Christensen			
D54	E	C			
D52	E	C			

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The screenshot shows the OtoAccess™ software interface. The top menu bar includes File, Edit, View, Help, and Login. The main window is divided into several sections. On the left, there's a 'Clients' tab with a search bar and a list of clients. The central area displays the details for a selected client, Ejvind Christensen, including fields for First name, Last name, Initials, Person ID, Birth date, Gender, Phone, Mobile, Fax, Email, Address, Address2, Zipcode, City, Country, State, and County. There's also a 'Picture' field with a placeholder image. On the right, there's a 'Select Instrument' section with a list of instruments: Affinity Suite, Diagnostic Suite, Titan Suite, Callisto Suite, and VIOT Suite. Below this, there's a 'Categories' section with two dropdown menus for Category 1 and Category 2. At the bottom, there's a table with columns for Person ID, First name, Last name, Telephone, Mobile, and City. The table contains several rows of data, including clients like Demo, Ejvind Christensen, and various instrument models like D51, D53, D54, and D52.

Person ID	First name	Last name	Telephone	Mobile	City
0101013333	Demo	Demo	+45 6371 3555		Assens
Titan	Ejvind	Christensen			
Diagnostic Suite	Ejvind	Christensen			
Affinity	Ejvind	Christensen			
Callisto	Ejvind	Christensen			
Callisto 1	Ejvind	Christensen			
Callisto 2	Ejvind	Christensen			
AffinityREM 2.0.5	Ejvind	Christensen			
AffinityAUD 2.0.5	Ejvind	Christensen			Assens
D51	E	C			
D53	E	Christensen			
D54	E	C			
D52	E	C			

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The screenshot displays the OtoAccess software interface. The top menu bar includes 'File', 'Edit', 'View', 'Help', and 'Login'. The main window is divided into several sections:

- Client Information:** Fields for First name (Ejvind), Last name (Christensen), Initials, Person ID (Diagnostic Suite), Birth date (19-04-1949), Gender (Male), Phone, Mobile, Fax, Email, Address, Address2, Zipcode, City, Country, State, and County. There is also a 'Picture' field with a placeholder icon and a checkbox for 'With schedule today'.
- Select Instrument:** A list of available instruments: Affinity Suite, Diagnostic Suite, Titan Suite, Callisto Suite, and VIOT Suite.
- Categories:** Fields for Category 1 and Category 2.
- Table:** A table with columns: Person ID, First name, Last name, Telephone, Mobile, and City. It lists various instruments and their associated details.

Person ID	First name	Last name	Telephone	Mobile	City
0101013333	Demo	Demo	+45 6371 3555		Assens
Titan	Ejvind	Christensen			
Diagnostic Suite	Ejvind	Christensen			
Affinity	Ejvind	Christensen			
Callisto	Ejvind	Christensen			
Callisto 1	Ejvind	Christensen			
Callisto 2	Ejvind	Christensen			
Affinity/REM 2.0.5	Ejvind	Christensen			
Affinity/AUD 2.0.5	Ejvind	Christensen			Assens
D51	E	C			
D53	E	Christensen			
D54	E	C			
D52	E	C			

At the bottom, there is a status bar showing 'Login: admin', 'Created: 07-05-2010 10:51:34', and 'Modified: 16-05-2011 08:17:18'.

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The screenshot shows the OtoAccess software interface. The main window has a menu bar (File, Edit, View, Help, Login) and a toolbar. The 'Clients' tab is active, displaying a form for a patient named Ejvind Christensen. The form includes fields for First name, Last name, Initials, Person ID, Birth date, Gender, Phone, Mobile, Fax, Email, Address, Address2, Zipcode, City, Country, State, and County. There is also a 'Picture' field and a 'With schedule today' checkbox. On the right side, there is a 'Select Instrument' section with a list of instruments: Affinity Suite, Diagnostic Suite, Titan Suite, Callisto Suite, and VIOT Suite. Below this, there are 'Categories' for Category 1 and Category 2. At the bottom, there is a table with columns: Person ID, First name, Last name, Telephone, Mobile, and City. The table contains several rows of data, including 'Titan', 'Diagnostic Suite', 'Affinity', 'Callisto', 'Callisto 1', 'Callisto 2', 'Affinity/REM 2.0.5', 'Affinity/AUD 2.0.5', 'D51', 'D53', 'D54', and 'D52'.

Person ID	First name	Last name	Telephone	Mobile	City
0101013333	Demo	Demo	+45 6371 3555		Assens
Titan	Ejvind	Christensen			
Diagnostic Suite	Ejvind	Christensen			
Affinity	Ejvind	Christensen			
Callisto	Ejvind	Christensen			
Callisto 1	Ejvind	Christensen			
Callisto 2	Ejvind	Christensen			
Affinity/REM 2.0.5	Ejvind	Christensen			
Affinity/AUD 2.0.5	Ejvind	Christensen			Assens
D51	E	C			
D53	E	Christensen			
D54	E	C			
D52	E	C			

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The screenshot shows the OtoAccess™ software interface. The main window is titled "Ejvind Christensen" and has a menu bar with "File", "Edit", "View", "Help", and "Login". Below the menu bar are two tabs: "Clients" and "Journal". The "Clients" tab is active, showing a form for a patient record. The form has fields for "First name" (Ejvind), "Last name" (Christensen), "Initials", "Person ID" (Diagnostic Suite), "Birth date" (19-04-1949), "Gender" (Male), "Phone", "Mobile", "Fax", "Email", "Picture", "Address", "Address2", "Zipcode", "City", "Country", "State", and "County". There is also a checkbox for "With schedule today". To the right of the form is a "Select Instrument" section with a list of instruments: Affinity Suite, Diagnostic Suite, Titan Suite, Callisto Suite, and VIOT Suite. Below this is a "Categories" section with "Category 1" and "Category 2" dropdown menus. At the bottom of the window is a table with columns: "Person ID", "First name", "Last name", "Telephone", "Mobile", and "City". The table contains several rows of data, including "0101013333", "Titan", "Diagnostic Suite", "Affinity", "Callisto", "Callisto 1", "Callisto 2", "AffinityREM 2.0.5", "AffinityAUD 2.0.5", "D51", "D53", "D54", and "D52".

Person ID	First name	Last name	Telephone	Mobile	City
0101013333	Demo	Demo	+45 6371 3555		Assens
Titan	Ejvind	Christensen			
Diagnostic Suite	Ejvind	Christensen			
Affinity	Ejvind	Christensen			
Callisto	Ejvind	Christensen			
Callisto 1	Ejvind	Christensen			
Callisto 2	Ejvind	Christensen			
AffinityREM 2.0.5	Ejvind	Christensen			
AffinityAUD 2.0.5	Ejvind	Christensen			Assens
D51	E	C			
D53	E	Christensen			
D54	E	C			
D52	E	C			

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