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1 Table of Contents and Figures

1	Table of Contents and Figures	51	7	Operation	64
1.1	Contents	51	7.1	Preparing for Operation	64
1.2	Figures	51	7.1.1	Component Check	64
			7.1.2	Installing a Filter Cartridge	64
2	Overview and Intended Usage	52	7.1.3	Closing the Housing Bell	65
2.1	Overview	52	7.1.4	Monitoring and Checks	66
2.2	Intended Usage	52	7.2	Operation	66
2.3	Technical Specifications	52	7.2.1	Sterilization	67
2.3.1	Chemical Compatibilities of the Seal Materials	53	7.2.2	Filtration of Liquids or Gases	68
2.3.2	Chemical Compatibilities of the Materials	54	7.2.3	Changing Filter Elements	68
2.3.3	Manufacturer's Label	55	7.3	Cleaning the Components	68
2.3.3.1	Cartridge Housing Series 7	56	8	Servicing Cleaning	69
2.3.3.2	Cartridge Housing with up to 8 Cartridges According to Sartorius Stedim Biotech M.D.S. Software (F&B Pharmaceuticals):	56	8.1	Service Life of the Filter Cartridge Housing	69
2.3.3.3	Cartridge Housings with More than 8 Cartridges According to Sartorius Stedim Biotech M.D.S. Software (F&B)	56	8.1.1	Limit of the Service Life Due to Load Cycles	69
2.3.4	General Specifications	57	8.1.2	Service Life Limited by Wear to the Housing Bell	69
3	Safety	58	8.2	Servicing and Inspection List	70
3.1	Notes Explanations	58	8.2.1	Visual Check for Rust	71
3.1.1	Filter Cartridge Housing Identification	58	8.2.2	Checking Clamps, Flanges and Threaded Couplings	71
3.2	EC Council Directives and Standards	58	8.2.3	Checking the Filter Cartridge Housing for Leaks	71
3.3	Built-In Safety Systems	58	8.2.4	Maintenance When Used in Hazardous Areas	47
3.4	Safety Precautions (To Be Carried Out by the Operator)	59	9	Troubleshooting	72
3.5	Responsibilities of the Operator	59	9.1	Removing the Vent Valve	73
4	General Warnings	60	10	Emergencies	74
4.1	Dangers	60	11	Shutdown Disposal	74
4.2	Installation of Replacement and Wear Parts	60	1.2	Figures	
4.3	Shutdown Procedure	60	Fig. 2-1	Overview (Example of a Filter Cartridge Housing)	52
5	Installation	61	Fig. 2-2	Manufacturer's Label	55
5.1	Equipment Supplied	61	Fig. 3-1	Manufacturer's Label	58
5.2	Transportation and Packaging	61	Fig. 5-1	Height of the Filter Cartridge Housing	62
5.2.1	Equipment Supplied (Also for Spare and Replacement Parts)	61	Fig. 5-2	Legs without and with Level Adjustment	62
5.3	Interim Storage	61	Fig. 6-1	Overview (Example of a Filter Cartridge Housing)	63
5.4	Storage Packaging for Further Transportation	61	Fig. 7-1	Clamps, flanges and threaded couplings	65
5.5	Installation	62	Fig. 7-2	Clamps, flanges and threaded couplings	65
5.5.1	Assembly	62	Fig. 7-3	Schematic Diagram of Sterilization	67
5.5.2	Attaching the Filter Cartridge Housing	62	Fig. 8-1	Clips on the Clamps, Flanges and Threaded Couplings	71
5.6	Installation of Pipes and a Heater	62	Fig. 9-1	Vent Valves	73
5.6.1	Connecting the Inlet and Outlet Pipes	62			
5.6.2	Connecting a Heater (Optional)	62			
5.6.3	Assembly When Used in Hazardous Areas	38			
6	Function	63			

2 Overview and Intended Usage

2.1 Overview


Fig. 2-1 | Overview
(Example of a Filter Cartridge Housing)

Depending on the equipment and the model, some components are not included with every filter cartridge housing.

The information in the enclosed dimension sheet is binding. The order number on the dimension sheet is identical to that on the Manufacturer's Label.

The filter cartridge housing consists of the following components:

- 1 Housing bell
- 2 Clamps, flanges and threaded couplings
- 3 Inlet
- 4 Run-off valve
- 5 Housing bell connection
- 6 Vent valve
- 7 Filter cartridge (not included in equipment supplied)
- 8 Seal
- 9 Base plate (lower part)
- 10 Outlet
- 11 Legs

2.2 Intended Usage

The filter cartridge housing has been developed, designed and built exclusively for installation in a system used for commercial operations. It serves to hold filter cartridges manufactured by Sartorius Stedim Biotech. The filter cartridge housing can be used to filter various media (gases and liquids). The medium is pressed through the filter cartridge by means of a pressure differential between the inlet and the outlet. The pressure in the inlet must always be greater than the pressure in the outlet.

The filter cartridge housing is intended solely for installation in a system used for industrial or commercial operations. Any other use or modification without the

written permission of the manufacturer is considered inappropriate use according to the manufacturer's instructions. The manufacturer shall not accept any liability for resultant damage. The risk shall be borne solely by the customer | operator.

Operation of the filter cartridge housing is prohibited until it has been established that the complete system of which it is a part complies with the regulations of the applicable EC directives. The filter cartridge, its seals and clamps, flanges and threaded couplings are not part of the filter cartridge housing. Always use parts suitable for the particular application. Comply with the legal regulations that apply to the place of use.

The filter cartridge housing materials and the seal materials are described in the dimension sheet.

Media that are characteristically abrasive to the standard materials are listed in the chemical compatibility tables in Chapter 2.3. Please obtain a written release from Sartorius Stedim Biotech for any special applications or if the application is unclear.

The customer | operator shall be responsible for purchasing and employing the media to be filtered in the filter cartridge housing according to its intended use. As a result, proper handling of these media and the related dangers are the sole responsibility of the customer | operator.

The filter cartridge housings are intended for the filtration of combustible or explosive media and for use in hazardous areas, provided that the special assembly and maintenance requirements pertinent to these requirements are strictly observed as detailed in Chapters 5.6.3 "Assembly When Used in Hazardous Areas" and 8.2.4 "Maintenance When Used in Hazardous Areas."

All movable filter units are standard-equipped with antistatic and lockable castor wheels or sliding rollers. To prevent sparking, the operator must ensure that the component groups are not moved and/or cannot collide with other objects during operation.

Although these Sartorius Stedim Biotech filter cartridge housings

do not generate heat on their own, they will take on the temperature of the working medium when hot fluids are processed; this includes steam sterilization and any process malfunctions. The user must either ensure that this temperature is permitted for the hazardous area where the filter cartridge housing is to be used or implement the appropriate precautionary measures. It is the user's responsibility to protect these products against all types of foreseeable mechanical damage, such as impacts and wear that might cause leakage.

If you have any further questions about the product's suitability for use in hazardous areas, please contact your local Sartorius Stedim Biotech subsidiary or dealer.

The filter cartridge housings are designed for explosion protection category 2 and are therefore approved for use in zone 1 or 21 hazardous areas. They bear the following marking of explosion protection:

II 2 G D c

Note: The use of original Sartorius Stedim Biotech filter cartridges and seals/gaskets is recommended. If filter cartridges and seals/gaskets from other manufacturers are used, the suitability for use in hazardous areas is not affected, provided that the precautions described in Chapter 5.6.3 are also strictly observed.

If you are using a filter cartridge heater, before use, always check that it is compatible with the medium to be filtered.

2.3 Technical Specifications

The dimension sheet lists the materials used. Please refer to this list and compare it with the medium to be filtered and the compatibility tables. These show the following:

- Suitability for use
- Service life
- If there are any potential hazards, this information must be provided to the operating and maintenance personnel in the form of a hazard assessment report or operating instructions.

Please obtain a written release from Sartorius Stedim Biotech for any special applications or if the application is unclear.

2.3.1 Chemical Compatibilities of the Seal Materials

Pos.	Abrasive Agents	Concentration	Temperature	Seal material		
				Silicone	EPDM	Fluoroelastomer
1	Acetic acid	96%	20°C	⊙	○	●
2	Hydrofluoric acid	25%	20°C	●	⊙	⊙
3	Hydrofluoric acid	50%	20°C	●	⊙	⊙
4	Hexafluorosilicic acid	50%	20°C	●	⊙	○
5	Perchloric acid	25%	20°C	⊙	○	○
6	Phosphoric acid	25%	20°C	○	○	○
7	Nitric acid	25%	20°C	●	⊙	○
8	Nitric acid	65%	20°C	●	●	⊙
9	Hydrochloric acid	25%	20°C	●	⊙	○
10	Sulfuric acid	25%	20°C	⊙	○	○
11	Trichloroacetic acid 25%	25%	20°C	●	○	●
12	Ammonium hydroxide	25%	20°C	○	○	⊙
13	Potassium hydroxide	32%	20°C	⊙	○	○
14	Sodium hydroxide	32%	20°C	⊙	○	○
15	Formalin	30%	20°C	⊙	○	○
16	Hydrogen peroxide	5%	20°C	○	○	○
17	Acetone		20°C	●	○	●
18	Gasoline, naphtha		20°C	●	●	●
19	Benzene		20°C	●	●	●
20	1-Butanol		20°C	○	○	○
21	Cellosolve		20°C	●	⊙	●
22	Methanol	98%	20°C	○	○	○
23	Toluene		20°C	●	●	⊙
24	Xylene		20°C	●	●	⊙
25	Trichloroethane		20°C	●	●	⊙
26	Trichloroethylene		20°C	●	⊙	⊙
27	Isopropanol		20°C	○	○	○
28	Formamide		20°C	●	○	●
29	n-Hexane		20°C	●	●	○
30	Ferric chloride		20°C	●	○	○

Key to symbols: ○ compatible, ⊙ relatively compatible, ● incompatible, - no data

Test conditions: 7 days at 20°C

2.3.2 Chemical Compatibilities of the Materials

Pos.	Abrasive Agents	Concentration	Temperature	Stainless Steel (Material No.)		
				1.4301 1.4541	1.4404 1.4435 1.4571 316L	1.4539
1	Acetic acid	50%	20°C	○	○	-
2	Hydrofluoric acid	40%	20°C	●	●	-
3	Hydrofluoric acid	50%	20°C	●	●	-
4	Salicylic acid	10%	20°C	○	○	-
5	Perchloric acid	25%	20°C	○	○	-
6	Phosphoric acid	25%	20°C	○	○	○
7	Nitric acid	25%	20°C	○	○	○
8	Nitric acid	65%	20°C	○	○	○
9	Hydrochloric acid	25%	20°C	◎	◎	-
10	Sulfuric acid	25%	20°C	◎	◎	○
11	Trichloroacetic acid	25%	20°C	◎	◎	-
12	Ammonium hydroxide	25%	20°C	○	○	-
13	Potassium hydroxide	32%	20°C	○	○	-
14	Sodium hydroxide	32%	20°C	○	○	-
15	Formalin	30%	20°C	○	○	-
16	Hydrogen peroxide	5%	20°C	○	○	-
17	Acetone		20°C	○	○	-
18	Benzene		20°C	○	○	-
19	Gasoline, naphtha		20°C	○	○	-
20	l-Butanol		20°C	○	○	-
21	Cellosolve		20°C	◎	◎	-
22	Methanol	98%	20°C	○	○	-
23	Toluene		20°C	○	○	-
24	Xylene		20°C	○	○	-
25	Trichloroethane		20°C	-	○	-
26	Trichloroethylene		20°C	○	○	-
27	Isopropanol		20°C	○	○	-
28	Tartaric acid	10%	20°C	○	○	-
29	Fruit acids		20°C	○	○	-
30	Ferric chloride		20°C	●	◎	-

Key to symbols: ○ compatible, ◎ relatively compatible, ● incompatible, - no data
 Test conditions: 7 days at 20°C

2.3.3 Manufacturer's Label

sartorius stedim
 biotech

Modell
 Model
 Ser.-Nr.
 Ser. No.

Raum / Chamber	I	II
Volumen V Volume V	L <input type="text"/>	<input type="text"/>
Prüfdruck/Datum PT Test Pressure/Date PT	bar <input type="text"/>	<input type="text"/>
min./max. Druck PS min./max. Pressure PS	bar <input type="text"/>	<input type="text"/>
min./max. Temperatur TS min./max. Temperature TS °C	<input type="text"/>	<input type="text"/>

CE 0045

Sartorius Stedim Biotech GmbH
 37070 Goettingen, Germany

Fig. 2-2 | Manufacturer's Label

Model:

XXXXX = Order number
 See Pos. 2.331 and 2.332

Ser. No.:

AAAA | BB CCCC DD
 A = Year of manufacture, e.g. 2001
 B = Month of manufacture, e.g. August = 08
 C = Drawing number
 D = Manufacturer's ID

Room | Chamber I, II

Indicates whether there are 1 or 2 pressure area

Volume = Volume of the pressure tank

Test pressure | Date PT

Test pressure in bar and date of the test

Min. | max. pressure PS

Min. | max. pressure (PS) in bar

Min. | max. temperature TS


Min. | max. temperature (TS) in °C

Manufacturer's Label with CE Mark of Conformity


of the registered location in conjunction with Directive 97/23/EU

Pressure tanks not covered by Directive 97/23/EU (Article 3/Para. 3) are not affixed with the CE mark of conformity.


2.3.3.1 Cartridge Housing Series 7

2.3.3.2 Cartridge Housing with up to 8 Cartridges According to Sartorius Stedim Biotech M.D.S. Software (F&B | Pharmaceuticals)

2.3.3.3 Cartridge Housings with More than 8 Cartridges According to Sartorius Stedim Biotech M.D.S. Software (F&B)

2.3.4 General Specifications

Ambient conditions:
Minimum temperature - 10 °C
Maximum temperature + 80 °C

Humidity: 10 – 80 %

Air : The ambient air should have a low salinity, and no dust or oxidizing substances (rust film) should be present.

Noise level (ear protection): < 70 dB (A)

In the case of large housing bells > 10 kg, the weight in kg is indicated on the housing bell without additional devices (e.g. heater).

The customer | operator should perform a noise test after installation because any transport noise caused by the medium being pumped should be taken into account.

3 Safety

3.1 Notes | Explanations

STOP Warnings are in a shaded frame and alert you with a stop sign.

Warning Symbols are in a shaded frame and alert you with a warning triangle.

Notes are in a shaded frame and alert you with a hand symbol.

Danger of Burns signs are in a shaded frame and alert you with the symbol shown opposite.

3.1.1 Filter Cartridge Housing Identification


Fig. 3-1 | Manufacturer's Label

The information in these operating instructions specifically applies to the filter cartridge housing of the model designation listed on the dimension sheet and must correspond with that stated on the Manufacturer's Label and no other. The Manufacturer's Label [1] with the CE mark of conformity is located on the filter cartridge housing (see Fig. 3-1). The Manufacturer's Label indicates the Serial Number. Make sure that the housing bell [2] is mounted on the base plate (lower part) [3] with the same model designation.

It is important that you state the following details correctly whenever you contact us with any questions or queries:

- Model name
- Serial number

We cannot process your query properly and quickly unless this information is correct.

3.2 EC Council Directives and Standards

The filter cartridge housing was built in conformity with the German versions of the following European Council Directives:

1. Pressure Equipment Directive (97/23/EC)
- 1.1 Technical Specifications AD-2000

3.3 Built-In Safety Systems

The built-in safety systems must be checked at regular intervals and in the case of product changes (**d** = daily, **w** = weekly, **m** = monthly, **y** = yearly).

The test methods to be applied are: **V** = visual check, **F** = function test, **M** = measurement.

– **Pressure gauge (optional)**

The filter housing is equipped with a pressure gauge for pressure monitoring and checks during operation and before disassembling pressurized parts.

Test	
Interval	Method
W	F

– **Vent valve (optional)**

In the case of liquid media the trapped air is released through the vent valve. The vent valve may only be opened by approx. two turns because otherwise liquid may escape uncontrolled.

Test	
Interval	Method
W	F

– **Run-off valve (optional)**

In the case of liquid media the medium is drained off through the run-off valve. The run-off valve may only be opened by approx. two turns because otherwise liquid may escape uncontrolled.

Test	
Interval	Method
W	F

These operating instructions are an integral part of the filter cartridge housing and must be kept readily at hand for the operators at all times.

The safety information contained in it must be complied with accordingly.

The operator must ensure that the filter housing is equipped with a safety valve that prevents inadmissible overpressure.

It is strictly forbidden to render any of the safety devices inoperative or to modify their function.

3.4 Safety Precautions (To Be Carried Out by the Operator)

- We would like to point out that the operator must instruct his operating, maintenance and servicing personnel:
- in the use of the filter cartridge housing’s safety devices and,
 - ensure that the safety precautions are complied with.
 - The operator must prevent access to the filter cartridge housing danger zone by unauthorized persons (who are not operating and servicing personnel).

These operating instructions must be kept in a safe place for future reference. The frequency of inspections and checks must be maintained | complied with.

- The work described in these operating instructions is explained in such a way that
- the chapters on operation and modes of operation can be understood by a trained worker and,
 - those on transportation, installation and assembly, maintenance and troubleshooting by a **qualified technician**.

The chapters on **transportation, installation and assembly, servicing and trouble-shooting** are intended for **qualified technicians** only. Only such technicians may carry out the work described in these chapters.

Trained Workers

A person familiarized with and, if applicable, trained in his particular duties by a qualified technician and instructed about the potential dangers in the event of improper conduct and on the necessary safety devices and safety precautions.

This particularly applies to hazardous areas. Such safety training instructions must be given and documented at least once a year.

Qualified Technicians

A person who, by virtue of his technical training, know-how, experience and knowledge of relevant standards, is able to assess the work assigned to him and identify potential dangers.

Definitions adapted from EN 60204-1.

3.5 Responsibilities of the Operator

Operators within the EEC (European Economic Community) must observe and comply with the respectively valid versions of the national transposition of framework Council Directive 89/391/EEC as well as the corresponding individual directives and specifically Council Directive 89/655/EEC on the minimum requirements for protecting the safety and health of workers using tools and equipment at work. The operator must obtain any local **operating permits** and comply with the provisions contained therein.

In addition to this, he must observe local laws and regulations on

- personnel safety (accident prevention regulations),
- safety of worker equipment (safety equipment and servicing),
- disposal of products (laws on wastes),
- disposal of materials (laws on wastes),
- cleaning (cleaning agents and disposal),
- environmental protection.

4 General Warnings

4.1 Dangers

The safety systems and safety notes described in this manual should be complied with accordingly.

Keep the work area of the filter cartridge housing free of objects during operation so that unimpeded access is possible at all times.

Watch out for liquid or gas leaks. They can result in death.

Make sure the vent valve is only opened by approx. two turns. If the vent valve is opened completely, liquid may escape uncontrolled.

Before opening any clamps, flanges or threaded couplings, make sure there is no pressure in the filter cartridge housing and that no media can escape.

Ensure that the materials are compatible. Protect the seal material from corrosion and degradation. Particles from eroded materials may accumulate in the medium being filtered. This can lead to damage or destruction of the filter cartridge.

Observe the required operating conditions, especially with regard to the temperature of the medium. It is transmitted to the components. Touching components or contact with the medium can result in burns or scalding.

Wear suitable protective clothing if necessary. Don't forget that the components of the filter cartridge housing store heat and cool down slowly.

4.2 Installation of Replacement and Wear Parts

Only use replacement parts and accessories that have been supplied, tested and approved by Sartorius Stedim Biotech. There is no guarantee that parts that have not been approved by Sartorius Stedim Biotech have been designed and manufactured to be stress-resistant and safe for installation and use with the filter cartridge housing. Sartorius Stedim Biotech GmbH accepts no liability for damages arising from the use of non-original parts and accessories.

4.3 Shutdown Procedure

The following shutdown procedure must be followed before any cleaning, servicing or repair work is carried out (by qualified technicians only).

1. Turn off circulation in the filter
2. Shut off the medium inlet and outlets
 - Close shut off valve | flap
3. Drain the filter cartridge housing
 - Check to make sure the filter cartridge housing is not under pressure
 - Open the vent valve by two turns
 - Open the drain valve by two turns
 - Drain off the medium

Personnel may suffer physical injuries or fatalities if this procedure is not followed!

5 Installation

5.1 Equipment Supplied

The equipment supplied consists of:

1. Filter cartridge housing as ordered (see confirmation of order)
2. Operating Instructions
3. The dimension sheet of the filter cartridge housing may be sent to the customer | operator separately. The dimension sheet contains order-specific data.

A detailed list of the equipment supplied is included on the dimension sheet.

5.2 Transportation and Packaging

Filter cartridge housings from Sartorius Stedim Biotech GmbH are carefully inspected and packed before shipment. Nevertheless, it is still possible that they might become damaged during transit.

5.2.1 Equipment Supplied (Also for Spare and Replacement Parts)

Incoming Inspection:

- Check the shipment against the delivery note to ensure that it is complete!

If the packaging is damaged

- Check the shipment itself for damage (visual check)!

Complaints

If the shipment was damaged during transit:

- Immediately contact the last carrier!
- Keep the packaging material (in case the carrier wants to inspect it or the shipment has to be returned to us)!

Packaging for return shipment

Use the original packaging material as far as possible. Pack the filter cartridge housing in such a way that mechanical damage is prevented.

5.3 Interim Storage

The external packaging of the filter cartridge housing and its spare and replacement parts is designed for a storage period of three months from delivery.

Storage Conditions

In a closed and dry room at a room temperature of -10°C to $+25^{\circ}\text{C}$.

If the external packaging has already been removed, follow the instructions given in Chapter 5.4 "Storage | Packaging for Further Transportation"

5.4 Storage | Packaging for Further Transportation

If necessary, the components must be decontaminated before packaging. Please follow the instructions in the data sheets for all media that have come in contact with the components.

- Clean the components with a suitable liquid, e.g. soapy water.
- Rinse off the liquid completely with hot water.
- Rinse the components with distilled water.

Do not use cloths for drying the components. They may produce lint. Use dry, oil-free compressed air for drying the components.

- Dry the components with dry, oil-free compressed air.
- Do not close the inlet, outlet and housing bell connection with a sealing cap.
- Remove the legs of the housing bell. Pack the legs in such a way that they cannot hit against the housing bell. Store the legs directly next to the housing bell.
- Make sure the vent valve and run-off valve cannot be exposed to blows during transportation and storage.
- No load may be placed on the housing bell. Stacking of the housing bell is prohibited.

For further transportation

Pack the filter cartridge housing in such a way that mechanical damage is avoided.

Note for Storage of the Filter Cartridge Housing

Ground the housing bell to prevent electrostatic charging.

5.5 Installation

Prerequisites

- Make sure the load bearing capacity of the floor is sufficient for the filter cartridge housing. Remember that the medium changes the weight of the filter cartridge housing.

The weights are stated on the dimension sheet.

- There must be sufficient space around the housing bell for operation and servicing.
- A person must have free access to the filter cartridge housing from any side.
- Make sure the ceiling is adequately high. The distance between the top of the base plate (lower part) and the ceiling should be at least 2.3 times greater than the height of the housing bell.


Fig. 5-1 | Height of the Filter Cartridge Housing

5.5.1 Assembly

There are three different versions of filter cartridge housing: filter cartridge housing without legs, with legs (without level adjustment) and with legs (with level adjustment).


Fig. 5-2 | Legs without and with Level Adjustment

Assembly of Components (without Level Adjustment)

- Screw the legs into the base plate (lower part) and tighten firmly.

Assembly of Components (with Level Adjustment)

- Screw the legs into the base plate (lower part) and tighten firmly.

Level Adjustment

- Set up the filter cartridge housing so that the vent valve is located above the filter cartridge. At the same time the run-off valve must be at the lowest point. The housing bell faces up.
- Using two fork wrenches, loosen the nuts on the adjustable leg (bottom).
- Change the length of the leg so that the filter cartridge housing stands firmly and tighten the nuts with the two fork wrenches.

5.5.2 Attaching the Filter Cartridge Housing

The filter cartridge housing with legs can be freestanding and does not require additional stabilization.

The filter cartridge housing without legs should be stabilized by attaching it to the inlet and outlet pipes and be secured with an additional bracket to reduce mechanical stress.

5.6 Installation of Pipes and a Heater

5.6.1 Connecting the Inlet and Outlet Pipes

Connect the inlet and outlet ports of the filter cartridge housing to the system. Various connection options are available depending on the type of filter cartridge housing. A description of the individual connections is provided in the Appendix.

Power should be disconnected for installation.

Pipes, valves, etc. installed to the housing have to be fixed separately. Make sure that their weight and forces have no influence on the housing connection.

5.6.2 Connecting a Heater (Optional)

When using a filter cartridge housing heater, follow the instructions in the separate operating manual delivered with the heater.

Connection to AC power (mains supply) may only be carried out by a certified electrician.

5.6.3 Assembly When Used in Hazardous Areas

Only qualified personnel may install the cartridge filter housings. In addition to the assembly instructions in the previous chapters, please comply with the following additional notes and assembly instructions when assembling the filters in hazardous areas:

Static electricity can be generated particularly if fluids with low conductivity come into contact with metallic enclosures. This can lead to discharge reactions with spark-overs and cause ignition in an explosive atmosphere.

In order to essentially eliminate this problem, the filter housings and the connected pipes, valves and supporting assemblies must be electrically cross-connected with a ground strap and routed to a grounding point.

During assembly, any filter housing components that are insulated with non-conductive seals/gaskets, etc., must be electrically cross-connected and grounded. Be sure to perform a conductivity test to check that the electrical grounding function works properly.

Sartorius Stedim Biotech offers optional accessories for this purpose.

It is prohibited to install any filter heater in hazardous areas unless the heating unit is specifically designed for that purpose and is affixed with the appropriate marking of explosion protection and CE conformity.

Use only non-sparking tools in hazardous areas. Mechanical sparks are capable of triggering ignition in explosive atmospheres.

Hazardous areas are subject to specific leakage requirements. To prevent damage or quality loss to the filter cartridge housing, the operator is required to thoroughly check that all construction materials (including seals/gaskets) are compatible with the working medium under the process conditions.

6 Function


Fig. 6-1 | Overview (Example of a Filter Cartridge Housing)

Depending on the model, not all components are supplied with each filter cartridge housing.

The filter cartridge housing consists of the following components:

1. Housing bell
2. Clamps, flanges and threaded couplings
3. Inlet
4. Run-off valve
5. Housing bell connection
6. Vent valve
7. Filter cartridge (not in equipment supplied)
8. Seal
9. Base plate (lower part)
10. Outlet
11. Legs

Description of the Functions

The filter cartridge housing is supplied with the medium through the inlet. The medium rises upwards in the housing bell and surrounds the filter cartridge. The medium is pressed through the filter cartridge due to the fact that there is overpressure at the inlet and a counteracting negative pressure at the outlet. The medium drains out through the outlet.

The housing bell connection is located at the top of the housing bell. This can be used for various applications, e.g. for connecting a pressure gauge. When filtering liquid media, trapped air can be evacuated through the vent valve.

The housing bell is connected to the base plate (lower part) by clamps, flanges or threaded couplings. A seal between the housing bell and base plate (lower part) prevents the medium from escaping. When filtering liquid media, the medium can be drained off at the run-off valve. This reduces the residual volume of the medium in the filter cartridge housing to a minimum.

Three legs provide the filter cartridge housing with a stable base.

7 Operation

Selecting Seals and Filter Elements

As the operator of the filter cartridge housing, you must ensure that the materials are compatible. Protect the seal material and filter element from corrosion and degradation.

Select a seal material and filter elements that are suitable for the application. Observe the required operating conditions, especially with regard to the medium, temperature and pressure. Also, observe the operating conditions required for sterilization.

7.1 Preparing for Operation

7.1.1 Component Check

All components of every filter cartridge housing have a serial number.

- Before assembling the filter cartridge housing, check to make sure that the housing bell and the base plate (lower part) have the same serial number.
- Check that the filter cartridge housing is clean.

7.1.2 Installing a Filter Cartridge

Be aware of the weight of the housing bell, use a lifting jack if necessary (the weight is stated on the bell or listed in the dimension sheet).

Follow the instructions included with the filter cartridge.

While installing a filter cartridge, leave the plastic bag on the cartridge for protection. Remove the bag only after you have completely installed the cartridge.

- Cut open the plastic bag on the O-ring end of the filter cartridge.
- Slide the plastic bag back to expose the O-ring.
- Wet the O-ring with pyrogen-free, distilled water.
- Firmly hold the filter cartridge by the upper and lower ends.
- Slightly turn the filter cartridge to insert it into or install it onto the base plate (lower part).
- Turn the filter cartridge all the way until it stops.
- Check to make sure that the filter cartridge is seated firmly.
- Remove the plastic bag from the filter cartridge.

Pressurizing the cartridge at the downstream outlet may destroy the cartridge.

Do not use a damaged filter cartridge for filtration.

Ensure that the pressure at the outlet does not exceed the maximum allowable pressure.

7.1.3 Closing the Housing Bell

Different types of clamps, flanges or threaded couplings are available:

- Clamps, flanges or threaded couplings with a wing nut [1 + 2]
- Clamps, flanges or threaded couplings with a hex nut [3]


Fig. 7-1 | Clamps, flanges and threaded couplings

- Segment locking clamp [4]
- Eccentric locking clamp [5+6]
- The centering pin can be seen in [6] below.


Fig. 7-2 | Clamps, flanges and threaded couplings

- Check the seat of the seal on the base plate (lower part).
- Make sure the seal is not damaged or misshapen.
- Place the seal onto the base plate (lower part).
- Place the housing bell vertically onto the seal. Prevent the housing bell from coming into contact with the filter cartridge.
- Position the clamps, flanges or threaded couplings so that they connect the housing bell and the base plate (lower part) together.
- Close the clamps, flanges and threaded couplings.

Clamps, flanges and threaded couplings with wing nut:

- Tighten the wing nut manually.

Clamps, flanges and threaded couplings with hex nut:

- Tighten the hex nut with a torque wrench.

Thread	Seal material	Tightening torque in Nm
All thread sizes	Fluoroelastomer	5 to 6 Nm
	Silicone	5 to 6 Nm
	EPDM	5 to 6 Nm
	PTFE	8 Nm

- Close the vent valve by manually tightening the knurled nut slightly.
- Do not use a tool (e.g. pliers) to tighten the knurled nut.
- Close the run-off valve by manually tightening the knurled nut slightly.
- Do not use a tool (e.g. pliers) to tighten the knurled nut.

7.1.4 Monitoring and Checks

After pressurizing the filter cartridge housing, check all connections. During operation, check all connections at regular intervals.

- Inlet
- Outlet
- Run-off valve
- Vent valve
- Housing bell connection
- Housing bell | base plate (lower part)
- Following application of a load and subsequent release, all tightening torques must be checked.

7.2

Operation

Danger if the technical specifications and the recommended ambient conditions (e.g. pressure, temperature, medium) of the filter cartridge housing, seal and filter cartridge are not complied with.

The operator must comply with the technical specifications and the recommended ambient conditions of all components.

The technical specifications for the filter cartridge housing are described in Chapter 2, in the DIN EN 10204 Certificate, on the filter cartridge housing dimension sheet (replacement parts) and in the certificate on production and hydraulic pressure testing.

Danger of injury from splashing liquids or leaking gases. Is caustic to skin and eyes.

When working with the filter cartridge housing, always wear the prescribed protective equipment, e.g.

- acid-resistant protective goggles,
- acid and temperature-resistant gloves.

Warning: Hot medium. This warning applies to all components as well. Touching components or contact with the medium can result in burns or scalding.

All surfaces are as hot as the medium. If appropriate, wear the necessary protective clothing. In addition, note that filter cartridge housing components store heat and release it very slowly.

It is necessary that the operator install the following components in the system to enable monitoring of and compliance with the technical specifications:

- Pressure gauge
- Temperature indicator at the medium inlet
- Safety valve

7.2.1 Sterilization

The housing with an installed filter cartridge can be autoclaved in the place or can be steam sterilized while installed in a system. Please refer to the pertinent Instructions for Use for the filter cartridge for the sterilization temperature.

Prerequisites for Sterilization

- The filter cartridge has been installed correctly and prepared for sterilization according to the manufacturer's specifications.
- The housing bell is closed.
- All valves are closed.
- The inlet and outlet are connected.

Proceed as follows to sterilize the housing with steam:

Make sure the filter cartridge has been installed properly and all valves are closed.


Fig. 7-3 | Schematic Diagram of Sterilization

V0 to V7 Valves

P1	Pressure gauge pressure
1	
P2	Pressure gauge pressure
2	
T1	Temperature
A	Medium
B	Compressed Air
C	Steam
D	Filtrate

- Open valve V0 to drain off condensate that has collected in the steam supply line.
- Close valve V0 and open the valves V4, V5 and V7.
- Open valve V1 slowly and allow the steam to flow into the housing.
- Close valve V4 only enough to create a small jet of steam.
- Leave the valve V7 open until the steam flows through the filter cartridge and then discharges through valve V5.
- Then close valve V5 and V7 and open valve V6. To prevent a differential pressure ($P1 - P2$) of more than 0.3 bar (4 psi), the flow of steam must be regulated at valve V1 and V6. At the same time the sterilization temperature at T1, e.g. 121°C, must be maintained.
- The sterilization time begins as soon as the temperature on the sterile side of the filter has reached the recommended level.
- When the sterilization time is up, close valves V1 and V6.
- To prevent a condensation vacuum during cooling, valves V7 and V5 can remain open or be opened.
- If the sterile side remains closed, compressed air (or N₂) can be filled into the system via valve V2.
- If drainage downstream from valve V6 leads into a tank that has also been sterilized, avoid a condensation vacuum there by venting or applying supporting pressure (otherwise the filter cartridge may be damaged).
- After cooling, the housing is ready for use.

7.2.2 Filtration of Liquids or Gases

Danger of injury from leaking liquids that result in burns and reactions in the liquids.

Only open the vent valve by approx. two turns.

Filling the Housing Bell

- Feed the medium or gas that is to be filtered into the housing.
- Use a filter cartridge housing with condensate tank to filter gases.

Venting | Draining the Condensate (During Filtration)

For Liquids

- Connect a hose to the vent valve connection.
- Fasten the hose or hold the hose tightly.
- Hold a tank under the hose.
- Open the vent valve (by turning twice).
- The housing bell is full when the liquid flows out of the tube without any bubbles.
- Close the vent valve. The filter is now ready to use.

For Gases

Venting is not necessary when filtering gases.

Monitoring

- After pressurizing the filter cartridge housing, check all connections. During operation, check all connections at regular intervals. Changes in temperature and pressure can affect the seal. Check inlet, outlet, run-off valve, vent valve, housing bell connection, housing bell and base plate (lower part).

7.2.3 Changing Filter Elements

Draining the Housing Bell

- Turn off circulation in the filter.
- Make sure no new medium can enter the filter cartridge housing.
- First close the valves at the filter inlet and then at the filter outlet.
- Check the pressure gauge to ensure the filter cartridge housing is not under pressure.
- Open the vent valve by approx. two turns.
- Place an adequately large tank under the run-off valve to collect the liquid.
- Open the run-off valve by approx. two turns.
- Wait until no further liquid escapes from the valve.
- Close the run-off valve by manually tightening the knurled nuts slightly. Do not use a tool (e.g. pliers) to tighten the knurled nut.

Opening the Housing Bell

- Check the pressure gauge once more to make sure that the filter cartridge housing is not under pressure.
- Open and remove the clamps, flanges or threaded couplings.

The housing bell, seal and base plate (lower part) may be bonded together so that the housing bell cannot be pulled off straight upwards. If this is the case, firmly hold the top of the housing bell. Loosen the components by lightly tapping the lower section of the housing bell with a rubber hammer.

- Remove the housing bell by lifting it straight upwards.
- Change the filter element.
- Remove the filter element by turning it once.
- Clean the components as described in Chapter 7.3.
- Install the new filter element as described in Chapter 7.1.2.
- Close the housing bell as described in Chapter 7.1.3.
- Fill the housing bell as described in Chapter 7.2.2. before starting up operations

7.3 Cleaning the Components

- Clean all the components of the filter cartridge housing immediately after every filtration.
- Clean the components with a suitable liquid, e.g. soapy water.
- Rinse off the liquid completely with hot water.
- Rinse the components with distilled water.
- Do not use cloths for drying the components. They may produce lint. Use dry, oil-free compressed air for drying the components.

These servicing and cleaning instructions are intended for qualified technicians only. **Servicing, cleaning** and repair work may only be carried out by such technicians.

Qualified Technician

- A person who, on account of his technical training, know-how, experience and knowledge of relevant standards, is able to assess the work assigned to him and identify potential dangers. **The above definition is based on the European Standard EN 60204-1.**

For trouble free operation it is absolutely essential that the filter cartridge housing is cleaned and serviced at regular intervals.

The following shutdown procedure must be followed before any cleaning, servicing or repair work is carried out (only by qualified technicians)

1. Turn off circulation in the filter
2. Shut off the medium inlet and outlets
 - Close shut off valve | flap
3. Drain the filter cartridge housing
 - Check to make sure the filter cartridge housing is not under pressure
 - Open the vent valve by two turns
 - Open the drain valve by two turns
 - Drain off the medium

Personnel may suffer physical injuries or fatalities if this procedure is not followed!

Danger of injury from splashing liquids or leaking gases. They can result in death.

Make sure the vent valve is only opened by approx. two turns. If the vent valve is opened completely, uncontrolled release of liquid will occur.

Before opening the filter housing, make sure that the filter cartridge housing is not under pressure and that no media can escape.

Warning: Hot Medium. This warning applies to all components as well. Contact with components or the medium may cause burns and scalding.

All surfaces are as hot as the medium. If appropriate, wear the necessary protective clothing. In addition, note that filter cartridge housing components store heat and only release it very slowly.

8.1 Service Life of the Filter Cartridge Housing

8.1.1 Limit of the Service Life Due to Load Cycles

The service life of the filter cartridge housings depends on the operating conditions at the place of use. The operator must determine the service life by referring to the applicable regulations (AD 2000 standards in the German Code of Practice "AD-Merkblätter" or the European Standard EN 13445).

8.1.2 Service Life Limited by Wear to the Housing Bell

Procedure:

- Perform a visual check of the welding seams for any wear or pitting.
- Check the wall thickness at several different places on the housing bell.
- Compare the wall thickness measured with the wall thickness stated on the drawing (filter cartridge housing dimension sheet, replacement parts).

If wear is more than 25 % of the original wall thickness, the filter cartridge housing is worn out. The filter cartridge housing has reached its maximum service life. Please withdraw the filter cartridge housing from use.

8.2 Servicing and Inspection List

Servicing Intervals	Check Point Servicing Notes	Refer to Chapter
Daily	<ul style="list-style-type: none"> • Visually check clamps, flanges or threaded couplings for any damage. • Check the filter cartridge housing for leaks. • Check the threads on clamps, flanges or threaded couplings: Tighten the wing nut manually. Tighten the hex nut with torque wrench 	8.2.2 8.2.3
Whenever opening the housing bell	<ul style="list-style-type: none"> • Make sure the seal is not damaged or misshapen. • Visually check if the housing bell contains rust. • Visually check clamps, flanges or threaded couplings for damage. • Check the gaskets on the run-off valve and the vent valve. 	8.2.1 8.2.2
Weekly	<ul style="list-style-type: none"> • Check that the pressure gauge is working properly. 	
Every 6 months	<ul style="list-style-type: none"> • Check the wall thickness of the housing bell. 	8.2.4
Every 12 months	<ul style="list-style-type: none"> • Test the pressure of the filter cartridge housing with water. • The test pressure is stated on the certificate on production and hydraulic pressure testing. 	
Maximum cycle number	<ul style="list-style-type: none"> • The filter cartridge housing has reached its maximum service life. Please withdraw the filter cartridge housing from use. 	

This information provides recommendations; however, the operator must adapt servicing intervals to meet the individual requirements of his or her on-site operating conditions.

8.2.1 Visual Check for Rust

- Check to see if any rust has collected in the housing bell, paying special attention to the weld seam.
- Remove the rust with a suitable cleaner.
- If the rust cannot be removed, inform your superior.

8.2.2 Checking Clamps, Flanges and Threaded Couplings


Fig. 8-1 | Clips on the Clamps, Flanges and Threaded Couplings

- Check whether the clips [1] are bent.
- Check the threads.
- Check the hinges for axial and radial play.
- Check the hex or wing nut.

8.2.3 Checking the Filter Cartridge Housing for Leaks

- During operation, check all connections at regular intervals. Changes in temperature and pressure can affect the seal. Check inlet, outlet, run-off valve, vent valve, housing bell connection, housing bell and base plate (lower part).

8.2.4 Maintenance When Used in Hazardous Areas

Only qualified personnel may perform maintenance work on the filter cartridge housings. Conductivity testing must be performed on all electrically cross-connected and grounded filter components whenever maintenance work is performed on the filters or at the least once a year. Please adhere to the maintenance intervals prescribed by the operator.

Use only non-sparking tools in hazardous areas. Mechanical sparks are capable of provoking ignition in explosive atmospheres.

Measures must be taken to remove any static residues in the fluid on the sterile downstream side of the filter, e.g. by using conductive pipes with the proper equipotential bonding.

Precautionary measures have to be implemented to prevent any static discharge from occurring during filter replacement. This is done by allowing the housing to stand for an appropriate period of

time (which is determined by the operator and depends on the working medium and operating conditions) so that the static discharge is allowed to disperse down to a safe level or by securing the housing environment for the corresponding activities.

If the filter is used to process flammable fluids, the user must make sure that leakage of liquid is minimized, controlled or drained off to a non-hazardous area during filling, venting, converting to standard pressure, draining and when replacing filter cartridges. The user must pay particular attention that such fluids do not come into contact with hot surfaces, which might cause the fluid to ignite.

When processing flammable fluids, the user must make sure that the filter housings are vented completely during filling and subsequent operation in order to prevent potentially flammable or explosive vapors or air mixtures from forming inside the filter unit. This can be accomplished by thorough venting prior to starting up operation.

Flammable fluids escaping from the housing due to improper installation or damage to the system (including the closing gaskets) can create a source of ignition if the fluids come into contact with a hot surface. The user must make sure that the filter unit is examined at regular intervals for damage and leaks so that it can be repaired immediately, if necessary, and that the new closing gaskets are installed every time the filter is replaced.

Regular cleaning with an antistatic medium is required to prevent any dust from collecting on the filter housings.

9 Troubleshooting

The situations and information on problems, their causes and correction in this Chapter are described in such a way that they can be understood by qualified technicians (see definition in Chapter 3.5 "Safety Precautions") with training in mechanical engineering | servicing. These technicians are to be equipped with the proper tools and test materials.

The shutdown procedure (see Chapter 4.5) must be followed before any servicing or repair work is carried out. Should the measures described below not prove successful, please contact Sartorius Stedim Biotech GmbH.

Problem	Reason Source of Problem	Remedy
Vent valve does not move smoothly	Threads, needle or seal are contaminated by grease or dirt	Remove and clean the vent valve. Install a new vent valve. See Chapter 9.1.
Vent valve cannot be opened.	Threads, needle or seal are by grease or dirt.	Place an open-ended wrench at hex or pliers at the knurl and open vent valve. Remove and clean the vent valve. Install new vent valve. See Chapter 9.1.
Vent valve leaks.	Seal are contaminated by grease or dirt.	Briefly open vent valve and close again. Remove and clean the vent valve. Install new vent valve. See Chapter 9.1.
Run-off valve does not move smoothly.	Threads or seal are contaminated by grease or dirt.	Remove and clean the run-off valve. Install new run-off valve.
Run-off valve cannot be opened.	Threads or seal are contaminated by grease or dirt.	Position open-ended wrench at hex or pliers at the knurl and open run-off valve. Remove and clean the run-off valve. Install new run-off valve.
Run-off valve leaks.	Seal are contaminated by grease or dirt.	Briefly open and close the run-off valve. Remove and clean the run-off valve. Install a new run-off valve.
Filter cartridge housing wobbles.	Legs have not been aligned. Legs have not been tightened.	Level the legs. Firmly tighten legs.
Flow rate drops during filtration.	There is air in the housing bell. Filter is blocked.	Vent housing bell. Replace filter cartridge.
Seal between housing bell and base plate (lower part) leaks	Clamps, flanges or threaded couplings, seal, sealing faces are defective. Clamps, flanges or threaded couplings are fastened too tightly.	Check hex nut tightening torque. Check wing nut tightening torque. Clean sealing faces. Replace seal. Check sealing faces for damage. Check clamps, flanges or threaded couplings, see Chapter 8.4. Replace clamps, flanges or threaded couplings.
Filter cartridge broken off	When opened, the housing bell not pulled off straight upwards.	Remove the broken off filter cartridge from the base plate (lower part). The broken off section of the filter cartridge might have fallen into the outlet. Disassemble the piping and remove the broken off section of the filter cartridge.
Rattling is heard during the filtration of fluids.	Air is in the housing bell.	Vent housing bell.

9.1 Removing the Vent Valve

- Drain the housing bell as described in Chapter 7.2.3
- Refer to the DIN EN 10204 certificate (Metallic Products, Test Certificates) to find out which vent valve is installed in your filter cartridge housing.


Fig. 9-1 | Vent Valves

- 1 Pharmaceutical-grade valve G 3/8
- 2 M8 x 0.5
- 3 Pharmaceutical-grade valve
- 4 Locking screw (not standard; available as an option)

Pharmaceutical-Grade Valve G 3/8 and M8 x 0.5

- Unscrew the vent valve.

Pharmaceutical-Grade valve with Safety Screw

The vent valve is protected by a locking screw. The valve can only be removed when this screw has been removed.

- Turn the vent valve completely until it stops.
- Remove the locking screw [4] with a hexagon socket.
- Unscrew the Vent Valve.

10 Emergencies

In the event of an emergency, stop the supply of medium.

Make sure all filter cartridge housings are not under pressure.

11 Shutdown | Disposal

Sort metals, non-metals, compound materials and process materials according to type and dispose of them in compliance with environmental safety regulations.

Always observe environmental safety, health hazard protection, disposal regulations and local options for proper disposal, especially those regarding hazardous waste. More detailed information can be obtained from your local authorities.

Cleaning agents must also be disposed of in accordance with local regulations.

Collect and dispose of all residual media according to the manufacturer's instructions.

Sartorius Stedim Biotech GmbH
August-Spindler-Strasse 11
37079 Goettingen, Germany

Phone +49.551.308.0
Fax +49.551.308.3289
www.sartorius-stedim.com

Copyright by Sartorius Stedim Biotech GmbH, Goettingen, Germany. All rights reserved. No part of this publication may be reprinted or translated in any form or by any means without the prior written permission of Sartorius Stedim Biotech GmbH. The status of the information, specifications and illustrations in this manual is indicated by the date given below. Sartorius Stedim Biotech GmbH reserves the right to make changes to the technology, features, specifications and design of the equipment without notice.

Status:
May 2014,
Sartorius Stedim Biotech GmbH,
Goettingen, Germany