Medtronic

5866-40M

Lead adaptor kit

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Elektrodenadapterbausatz

Geleidingsdraad-adapterset

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Technical manual • Manuel technique

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Technische handleiding • Manuale tecnico

! USA

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

The following are trademarks of Medtronic:

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
Medtronic

Sterile package opening instructions

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Explanation of symbols on package labeling

Refer to the package label to see which symbols apply to this product.

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Conformité Européenne
(European Conformity).
This symbol means that the
device fully complies with
European Directive AIMD
90/385/EEC.

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



Description/intended use

The Model 5866-40M Lead Adaptor Kit is designed to connect a Medtronic pacing lead with a bipolar connector (3.2 mm low-profile) to a pulse generator featuring a bipolar connector block which meets the (IS-1 BI)¹ standard and has a short-pin cavity and no sealing rings.

The lead adaptor features MP35N nickel alloy conductor coils, silicone rubber insulation, and a stainless steel connector.

The Model 5866-40M is intended for single use only.

Contraindications

To date, there are no known contraindications to the use of the Model 5866-40M Lead Adaptor Kit.

¹ IS-1 refers to an International Connector Standard (ISO 5841-3) whereby pulse generators and leads so designated are assured of a basic mechanical fit.

Contents of package

1 Model 5866-40M lead adaptor

4 #2 Wrench and setscrew assemblies

1 Medical adhesive tube with applicator tip

Documentation

Instructions for use

Before the lead adaptor is installed, detach the lead from the existing pulse generator and perform the following instructions:

1. Examine the lead for signs of connector pin corrosion, pitting, or any other visible deterioration and test the lead for adequate pacing and sensing thresholds.

Caution: If the lead is in any way unreliable, do not use it.

2. Gently insert the male connector end of the pacing lead into the female receptacle end of the adaptor until seated firmly (see Figure 1). If necessary, use sterile water as an insertion lubricant. Ensure the tip of the connector pin is visible through the inner connector block hole.

Figure 1. Insert the lead connector into the adaptor receptacle

3. Use the supplied hex wrenches with setscrew attached and turn a setscrew into each threaded hole located on each side of the receptacle end of the adaptor (see Figure 2).

Caution: Do not overtighten the setscrews; this could damage the lead connector or the adaptor receptacle.


Figure 2. Turn the setscrews into the threaded holes on the adaptor

4. Remove the hex wrench from the setscrews with a gentle back and forth motion.
5. Gently pull on the leads and adaptor to test the strength of the connection.

6. Tie a nonabsorbable ligature around the adaptor boot body.
7. Perform electrical measurements to determine the integrity of the connection.
8. Attach the applicator tip to the medical adhesive tube. Fill and cover the two setscrew hole openings with medical adhesive (see Figure 3).


Figure 3. *Apply medical adhesive in the two setscrew holes*

9. Insert the male connector end of the adaptor into the connector block of the pulse generator and tighten the setscrews as described in the pulse generator technical manual.
10. Insert the pulse generator into the subcutaneous pocket by coiling excess lead and adaptor lengths within the pocket as indicated in the lead technical manual.

Warnings

An implanted lead forms a direct low-resistance current path to the myocardium. Care must be taken during surgical procedures to insulate the terminal pins of an implanted lead from leakage currents that may arise from nearby line-powered equipment.

Since the pacing lead adaptor is a functional extension of the pacing lead, it is recommended that the pacing lead manual be reviewed for all appropriate warnings, complications and precautions.

Precautions

Necessary Hospital Equipment

Defibrillating equipment should be kept nearby for immediate use during the implantation procedure and whenever arrhythmias are possible or intentionally induced.

Handling the Lead Adaptor

Lead insulators attract small particles such as lint and dust; therefore, protect the lead adaptor from materials shedding these substances.

Ensuring Package Integrity

Medtronic has sterilized the package contents with ethylene oxide before shipment. This device is for single use only and is not intended to be resterilized.

Carefully inspect the package before opening.

It is not recommended that you use the product after its expiration date.

Potential complications

The potential complications related to the use of a lead system, of which a lead adaptor is a part, include, but are not limited to, the following patient-related conditions: infection, myocarditis, constrictive pericarditis, fibrillation, hemothorax, and bleeding after placement, after repositioning, or if sutures tear out.

Other potential complications related to the lead system and the programmed parameters of a pulse generator include, but are not limited to, the following:

Complications	Possible Effects
Displacement, lead insulator or conductor fracture, poor connection to the pulse generator or lead.	Intermittent or continuous loss of sensing and pacing; intermittent or continuous loss of defibrillation or cardioversion therapy; and possible muscle stimulation in the pocket area.
Electrode fracture	Intermittent or continuous loss of pacing, sensing, defibrillation, or cardioversion therapy.
Threshold elevation	Loss of pacing, sensing, defibrillation, or cardioversion therapy.

Special notice

Medtronic Lead Adaptor Kit consist of lead adaptors and tools to connect the adaptor to implantable leads and pulse generators. Adaptors are used with leads, which are implanted in the extremely hostile environment of the human body. Leads, when used with adaptors, may not be the functional equivalent of leads that are used without adaptors. Leads or adaptors may fail to function for a variety of causes, including, but not limited to: medical complications, body rejection phenomena, allergic reaction, fibrotic tissue, or failure of leads by breakage or by breach of their insulation covering. In addition, adaptors and tools may easily be damaged by improper handling or use. For leads used with adaptors, no representation or warranty is made that failure or cessation of function will not occur, or that the body will not react adversely to the implantation of leads or that medical complications (including perforation of the heart) will not occur.

Disclaimer of warranty

The following Disclaimer of Warranty Applies to United States Customers Only:

ALTHOUGH THE MEDTRONIC MODEL 5866-40M LEAD ADAPTOR KIT HEREAFTER REFERRED TO AS “PRODUCT” HAS BEEN MANUFACTURED UNDER CAREFULLY CONTROLLED CONDITIONS, MEDTRONIC HAS NO CONTROL OVER THE CONDITIONS UNDER WHICH THIS PRODUCT IS USED. MEDTRONIC, THEREFORE DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. MEDTRONIC SHALL NOT BE LIABLE TO ANY PERSON OR ENTITY FOR ANY MEDICAL EXPENSES OR ANY DIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, FAILURE OR MALFUNCTION OF THE PRODUCT, WHETHER A CLAIM FOR SUCH DAMAGES IS BASED UPON

WARRANTY, CONTRACT, TORT OR OTHERWISE. NO PERSON HAS ANY AUTHORITY TO BIND MEDTRONIC TO ANY REPRESENTATION OR WARRANTY WITH RESPECT TO THE PRODUCT.

The exclusions and limitations set out above are not intended to, and should not be construed so as to contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held to be illegal, unenforceable or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

The following Disclaimer of Warranty Applies to Customers Outside the United States:

Although the Medtronic Model 5866-40M Lead Adaptor Kit hereafter referred to as "Product" has been carefully designed, manufactured and tested prior to sale, the Product may fail to perform its intended function satisfactorily for a variety of reasons. The warnings contained in the Product labeling provide more detailed information and are considered an integral part of this **DISCLAIMER OF WARRANTY**.

Medtronic, therefore, disclaims all warranties, both express and implied, with respect to the Product. Medtronic shall not be liable for any incidental or consequential damages caused by any use, defect or failure of the Product, whether the claim is based on warranty, contract, tort or otherwise.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this **DISCLAIMER OF WARRANTY** is held by any court of competent jurisdiction to be illegal, unenforceable or in conflict with

applicable law, the validity of the remaining portion of the DISCLAIMER OF WARRANTY shall not be affected, and all rights and obligations shall be construed and enforced as if this DISCLAIMER OF WARRANTY did not contain the particular part or term held to be invalid.

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© Medtronic, Inc. 2010
Medtronic B.V.
Medtronic Australasia Pty Ltd