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AUTO CHART PROJECTOR

Ref. **08245**

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
Ref. 08245

## WARRANTY AND SERVICE

- 1 CHECK THE CONTENTS OF ALL THE BOXES
- 2 NOTIFY CUSTOMER SERVICE OF ANY PROBLEMS WITHIN 8 DAYS

If any fault develops within the warranty period, do not open the cleaner but contact Technical Assistance immediately. All parts recognised as having manufacturing faults will be replaced under the terms of the warranty.

The warranty covers faults and breakages not caused by normal wear, provided the cleaner has not been tampered with. Consumables are not covered by the warranty.

\* Within one year of purchase

full warranty

\* from the date of invoice.

### IN CASE OF MALFUNCTIONING

**Never return the cleaner or parts without prior notice!**

Contact your nearest authorised Service Centre before returning the cleaner or any parts, or to request repairs under the warranty.

Centro Style S.p.A. Via G.D. Martinengo, 21040 Vedano Olona, Varese.

Fax: + 39 0332 270 405

E-mail: [export.dpt@centrostyle.it](mailto:export.dpt@centrostyle.it)

Call to clarify the following points:

- 1 - Cleaner code and serial number
- 1 - Type of fault
- 3 - Quick over-the-phone troubleshooting and problem solving
- 4 - Quick over-the-phone troubleshooting with delivery of spare parts
- 5 - Always state the purchase date, the number of the invoice or delivery note and your customer code.

# REMOVAL FROM SERVICE

## Decommissioning and disposal.

At the end of its life cycle, or in the event of irreparable damage, dispose of the tester immediately or render it unserviceable to prevent further use and unnecessary risk.

- Disconnect the mains power cable and cut the plug off to make the Auto Chart Projector unusable.

## Notes on the collection and recycling of different materials.

Dispose of and recycle the Auto Chart Projector in compliance with European refuse disposal standard 160200 governing "waste equipment and equipment parts".

# TECHNICAL SPECIFICATIONS

Chart	33 Selections
Chart selection speed	Average 0.3 sec.
Mask	1 open 5 horizontal lines 5 vertical lines 21 single letters
Filter	Red / Green
Projection distance	2.7 m ~ 7.0 m (5 m is standard)
Projection magnification	30X (at 5 m)
Projection size	252 mm and 330 mm (w) x 225 mm (h) (a 5 m)
Tilt angle	Ball joint ( ± 10 degree)
Lamp	6V 30W (Halogen), 2,000 hrs lifetime
Power source	AC110,220~240V 50 / 60Hz
Power consumption	50VA
Auto-OFF function	after 5 minutes
Dimensions (mm) w x p x h	Main body: 230 x 360 x 265 mm Remote controller: 64 x 195 x 20 mm
Net weight	Main body: 6.0 Kg Remote controller: 180 g
Interface	RS232
Ambient operating conditions	Temperature between 15° and 30° C. Relative humidity <80%
	The appliance complies with the following directives: 73/23 CEE Bassa Tensione 89/336 CEE Compatibilità Elettromagnetica

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## GB Description of controls

- 1 Upper case
- 2 Projection Lens
- 3 Window for remote detection
- 4 Table stand
- 5 Focusing Wheel
- 6 Pilot Lamp
- 7 Power Switch ON/OFF
- 8 RS-232C Connector
- 9 Fuse Holders
- 10 Power Connector
- 11 Infrared ray transmitter
- 12 Alphabet chart buttons
- 13 Number chart buttons
- 14 Illiterate E chart buttons
- 15 Program buttons
- 16 Lamp on/off button
- 17 Red/Green button
- 18 Children chart buttons
- 19 Vectograph chart buttons
- 20 Direction buttons
- 21 Vertical, horizontal, single letter Mask buttons
- 22 Buzzer
- 23 Reset

# MAINTENANCE AND CLEANING

-  Always disconnect the mains power cable before performing any maintenance or cleaning.
- There are no cleaning requirements other than regular office housekeeping, such as dusting.
-  Periodically check the mains power cable and the plug for wear or damage. Also check the external shell of the Auto Chart Projector.
-  Clean the external surfaces of the Auto Chart Projector with a soft cloth, using only small quantities of neutral detergents to prevent water infiltrating the circuitry.
- Do not use solvents or other aggressive substances. Do not spray liquids on to the Spectrotester. Do not use abrasive materials.

## • Cleaning the projection lens

- 1) Blow dusts with a blower brush.
- 2) If it is not clean enough, wipe with lens cleaning paper.

## • Disassembling the upper case

- 1) Turn off the power and pull out the power cord.
- 2) Unfasten the screw in counterclockwise at the back side of case ❶.
- 3) Disassembling the Upper ❶ Case by lifting it.
- 4) In case of reassembling, place the upper case ❶ on the main body and adjust the front position. Then, fasten the screw in clockwise.

## • Replacement of lamp

-  The lamp may be HOT! Do not touch the Lamp directly and allow the sufficient amount of time to cool down the heat.

- 1) Unfasten the screw at the Lamp Cover.
- 2) Disassembling the Lamp Cover by lifting it and pull out the Lamp.

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 **WARNING! Read this manual carefully before using the appliance.** Keep this manual safe for reference throughout the life of the appliance. Copyright Centro Style S.p.A.. All rights reserved. Total or partial reproduction of this manual is strictly forbidden without the prior written consent of Centro Style S.p.A.

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GB

### **Classifications under the provision of 93/42/EEC (MDD) Class I**

Auto Chart Projector is classified into Class I system

### **Protection method against electric shock Class I**

Auto Chart Projector is classified into Class I instrument.

Class I instrument in which protection against electric shock does not rely basic insulation only, but which includes an additional safety precaution with the provision of the instrument to the protective earth conductor in fixed wiring of the installation in such a way that any accessible conductive parts cannot become live in the event of a failure of the basic insulation.

### **Degree of protection against ingress of liquids IPXO**

Auto Chart Projector is the ordinary instrument (enclosed instrument without protection against ingress of liquids). Be careful not to splash water on the instrument.

### **Degree of protection against flammability**

Auto Chart Projector is classified as an equipment not suitable to be used in a Potentially flammable anesthetic mixture with air or with oxygen or nitrous oxide. Do not use near flammable material.

### **Mode of operation** Continuous operation

# SAFETY

## Improper and prohibited uses and residual risks.

- Make sure that the voltage of the main power supply corresponds to that specified on the data label underneath the Auto Chart Projector.

• The Auto Chart Projector is not impermeable to water and other liquids. Do not immerse it in water, or expose it to splashes or to ambient humidity in excess of 85%. Operate the Auto Chart Projector only in temperatures between 15° and 30° C.

- Do not attempt to repair the Auto Chart Projector yourself. Refer all repairs to a Centro Style service centre (see page 10).

• Do not damage the body of the Auto Chart Projector or its mains power cable. Risk of electric shock, accidents, malfunctions and fire! Have any damaged parts replaced immediately by an authorised service centre.

- The electrical power supply to which the Auto Chart Projector is connected must conform to applicable standards.

The Auto Chart Projector is designed and made exclusively for the purposes described in this manual as normal, specified or deducible uses. Any other use is considered improper.

Centro Style S.p.A. declines all responsibility for injury to persons or damage to property due to:

- Improper or prohibited uses of the Auto Chart Projector.
- Unauthorised modifications to the Auto Chart Projector or its accessories.
- Use of non original spare parts.

## BASIC TROUBLE SHOOTING

Always disconnect the power cord from the wall and the instrument prior to performing any of the following procedures.

### **If the instrument does not function at all:**

- 1) Check the facility power source. (Circuit Breakers)
- 2) Check the electrical connections. (Power Cord)
- 3) Check the main fuse located on the rear of the unit.

### **If the projector lamp does not light:**

- 1) Check the main ON/OFF Switch .
- 2) Replace the lamp with the spare lamp.

The lamp may be HOT! Do not touch the Lamp directly and allow the sufficient amount of time to cool down the heat.

### **If the projector turns on but does not function:**

- Check the batteries in the remote control.

- To install Programm "P2"
  - 1) Keep pressing the left button of  for 2-3 seconds and release the button upon hearing a beep sound;
  - 2) Press the charts needed to store one by one, with a limit of at most 30 charts;
  - 3) Keep pressing the right button of  for 2-3 seconds to store selected charts and release the button upon hearing a beep sound;
- Programmed Operation
  - 1) Press 

The first chart in the program will be viewed and the programmed operation started. (Press the button once again to switch over between P1 and P2)
  - 2) Press the left of right button to change charts.

### 5. Reset

Press RESET button  to re.position charts to "0.05E".

### 6. Mute

Press the button  to switch over between mute and sound states.

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## START-UP AND OPERATION

### 1. Turn the power switch of instrument on .

The projection lam light and the chart for 0.05 vision shall be viewed in 3 seconds.

### 2. Changing the chart

Press the appropriate chart button the remote controller on.  
The charts have 33 selections

### 3. Isolating the visual acuity chart

- To use Vertical line masking: Press 
  - To move the position right or left: Press  or 
  - To move the position up or down: Press  or  0.2  
0.3  
0.4
- To use Horizontal line masking: Press 
  - To move the position right or left: Press  or 
  - To move the position up or down: Press  or  0.2  
0.3  
0.4
- To use Single letter masking: Press 
  - To move the position right or left: Press  or 
  - To move the position up or down: Press  or  0.2  
0.3  
0.4
- To use Red/Green filter: Press 
- To release masking and filter: Press any chart button  
 0.2  
0.3  
0.4


### 4. Programmed operation

There are two programs for selection, which are Program "P1" and "P2". Be sure to press  to turn off the lamp before program selecton.


- To install Programm "P1"
  - 1) Keep pressing the left button of  for 2-3 seconds and release the button upon hearing a beep sound;
  - 2) Press the charts needed to store one by one, with a limit of at most 30 charts;
  - 3) Keep pressing the right button of  for 2-3 seconds to store selected charts and release the button upon hearing a beep sound;
  - 4) Press  to turn on the lamp and exit programmig.

The Auto Chart Projector should be installed according to the following procedures:

1. If using a wall mount, it is best to locate a wall stud to support the weight of instrument.
2. Determine the refracting distance (from the patient's eye to the screen) and install the projection screen.
3. Position the instrument at the same distance from the screen to the patient within the range between 2 m and 7 m. To obtain longer refracting distances in small rooms, a mirror or systems of mirrors can be used in this case, a high quality front surface mirror is required.

4. Check the three-dimensional alignment of the system. The instrument is optimized when the projection screen is angled to direct light to the patient's head. Place a mirror on the screen. The light should project where the patient's head would be.

## FOCUSING

To obtain the correct focus:

1. Measure the distance from the patient to the screen (refracting distance).
2. Turn the Auto Chart Projector on, and using the remote control, project the 0.05 "E" onto the screen.
3. Adjust the position of the instrument forward or backward for sizing.
4. Adjust for sharpness and clarity by turning the focusing wheel 5 left or right as needed.

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

# DESCRIPTION OF THE AUTO CHART PROJECTOR

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In this area of subjective refraction, the Auto Chart Projector is best used for measuring visual acuity. It is fully recommended that you carefully read, understand and follow the steps in this manual order to ensure safe operation, optimum performance and a longer service life from your new instrument. Please retain this manual for future reference.

### PURPOSE:

The machine is a professional instrument and must only be used by adequately trained personnel.

### APPLICATIONS:

The machine is designed for use in optical labs and optical outlets. The machine is a non-invasive test appliance and is not suitable for medical diagnostic purposes.

### WARNING!

Do not use this machine near by other machines of electric bulb (ex. Halogen lamp), it will improperly operate on the effect of electronic waves.

## ASSEMBLY AND INSTALLATION

Open the packaging and check that the contents are not damaged. Also make sure that all the contents listed on page 1 (Standard accessories) are present. Contact technical assistance if the Auto Chart Projector has been damaged in transit.

- Install the instrument in a stable and level place where vibration or shock does not occur. The instrument may not perform observation correctly or may malfunction.
- If using a wall mount, it is best to locate a wall stud to support the weight of instrument.
- Ambient operating conditions: Temperature between 15 and 30° C. Relative humidity: <85%. Avoid dusty environments.

- **The electrical mains supply to which the tester is connected must conform to applicable standards in the country of installation.**
- If the instrument is not to be used for a long time, disconnect the power cord from the wall outlet. Otherwise, it may cause a fire. Remove batteries from the battery case of the remote controller.

Page	1	<b>TRANSPORT AND STORAGE.</b>
Page	1	<b>STANDARD ACCESSORIES.</b>
Page	2	<b>DESCRIPTION OF THE AUTO CHART PROJECTOR</b> , proper use and applications.
Page	2	<b>ASSEMBLY AND INSTALLATION.</b>
Page	3	<b>FOCUSING.</b>
Page	4	<b>START-UP AND OPERATION.</b>
Page	5	<b>CHART DESCRIPTION.</b>
Page	6	<b>MAINTENANCE AND CLEANING.</b>
Page	8	<b>SAFETY</b> , improper and prohibited uses and residual risks.
Page	8	<b>BASIC TROUBLE SHOOTING.</b>
Page	9	<b>REMOVAL FROM SERVICE.</b>
Page	9	<b>TECHNICAL SPECIFICATIONS.</b>
Page	10	<b>WARRANTY AND SERVICE.</b>

## TRANSPORT AND STORAGE

The Auto chart Projector and the accessories supplied with it must be left in their original packaging during transport and storage.

- Do not store the instrument in a place where it may get wet or where poisonous gas or liquid is stored.
- Do not expose the Auto chart Projector to temperatures below 10° C or above 70° C, relative humidity above 85%, or particularly dusty atmospheres.
- Do not stack more than 10 boxes or weights in excess of 15 Kg. Avoid impacts and mechanical stresses.
- Atmospheric pressure range within 500 to 1060hPa

## STANDARD ACCESSORIES

- 1 x Auto chart Projector, code 08245
- 1 x Polarized Metal Screen (400 mm x 350 mm)
- 1 x Power cable with Schuko plug
- 1 x Remote Control
- 2 x Spare fuses (2A, 250V)
- 2 x AAA type batteries
- 1 x Instruction manual

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