Medtronic

Medtronic Open Pivot™

Heart Valve

Valve cardiaque

Herzklappe

Válvula cardíaca

Valvola cardiaca

hjerterklaff

Hjärtklaff

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Heart Valve

1. Device Description

The Medtronic Open Pivot Heart Valve (Models 500, 501, 503, and 505) orifice and leaflets are constructed of pyrolytic carbon. The sewing cuff is constructed of titanium or cobalt chrome, and polyester fabric. The Medtronic Open Pivot Heart Valve is available in the aortic and mitral configurations in several sewing cuff styles, the Standard, AP, APex, and AP360. The AP Series valves are supra-annular configurations of the Standard Heart Valve. The valves are available in the dimensions and sizes shown in Table 1 through Table 5.

Table 1. MEDTRONIC OPEN PIVOT STANDARD HEART VALVE SPECIFICATIONS

Valve Size (mm)	Model Number	Overall Profile Height (mm)	Inflow Profile Height (mm)	Outflow Profile Height (mm)	Inside Diameter I.D. (mm)	Tissue Annulus Diameter (mm)	Geometric Orifice Area (cm ²)
AORTIC							
19	500FA19	9.3	2.4	6.9	14.8	19.2	1.55
21	500FA21	10.3	2.9	7.4	16.8	21.2	2.02
23	500FA23	11.3	3.4	7.9	18.8	23.2	2.56
25	500FA25	12.3	3.4	9.2	20.8	25.2	3.17
27	500FA27	13.3	3.4	9.9	22.8	27.2	3.8
29	500FA29	14.3	3.9	10.4	24.8	29.2	4.59
31	500FA31	15.4	4.6	10.8	26.8	31.2	5.35
MITRAL							
19	500DM19	9.3	3.0	6.3	14.8	19.2	1.55
21	500DM21	10.3	3.0	7.3	16.8	21.2	2.02
23	500DM23	11.3	3.0	8.3	18.8	23.2	2.56
25	500DM25	12.3	3.3	9.0	20.8	25.2	3.17
27	500DM27	13.3	3.7	9.6	22.8	27.2	3.84
29	500DM29	14.3	3.7	10.6	24.8	29.2	4.59
31	500DM31	15.4	3.7	11.7	26.8	31.2	5.35
33	500DM33	15.4	3.7	11.7	26.8	33.2	5.35

Table 2. MEDTRONIC OPEN PIVOT AP360 HEART VALVE SPECIFICATIONS

Valve Size (mm)	Model Number	Overall Profile Height (mm)	Inflow Profile Height (mm)	Outflow Profile Height (mm)	Inside Diameter I.D. (mm)	Tissue Annulus Diameter (mm)	Geometric Orifice Area (cm ²)
AORTIC							
16	505DA16	9.3	1.5	7.8	14.8	16.2	1.55
18	505DA18	10.3	1.5	8.8	16.8	18.2	2.02
20	505DA20	11.3	1.5	9.8	18.8	20.2	2.56
22	505DA22	12.3	1.5	10.8	20.8	22.2	3.17
24	505DA24	13.3	1.5	11.8	22.8	24.2	3.84
26	505DA26	14.3	1.5	12.8	24.8	26.2	4.59
28	505DA28	15.4	1.5	13.9	26.8	28.2	5.35
MITRAL							
16	505DM16	9.3	4.8	4.5	14.8	16.2	1.55
18	505DM18	10.3	5.3	5.0	16.8	18.2	2.02
20	505DM20	11.3	5.8	5.5	18.8	20.2	2.56
22	505DM22	12.3	6.3	6.1	20.8	22.2	3.17
24	505DM24	13.3	6.7	6.6	22.8	24.2	3.84
26	505DM26	14.3	7.2	7.1	24.8	26.2	4.59
28	505DM28	15.4	7.7	7.7	26.8	28.2	5.35

Table 3. MEDTRONIC OPEN PIVOT AP HEART VALVE SPECIFICATIONS

Valve Size (mm)	Model Number	Overall Profile Height (mm)	Inflow Profile Height (mm)	Outflow Profile Height (mm)	Inside Diameter I.D. (mm)	Tissue Annulus Diameter (mm)	Geometric Orifice Area (cm ²)
AORTIC							
16	501DA16	9.3	1.5	7.8	14.8	16.2	1.55
18	501DA18	10.3	1.5	8.8	16.8	18.2	2.02
20	501DA20	11.3	1.5	9.8	18.8	20.2	2.56
22	501DA22	12.3	1.5	10.8	20.8	22.2	3.17
24	501DA24	13.3	1.5	11.8	22.8	24.2	3.84

Valve Size (mm)	Model Number	Overall Profile Height (mm)	Inflow Profile Height (mm)	Outflow Profile Height (mm)	Inside Diameter I.D. (mm)	Tissue Annulus Diameter (mm)	Geometric Orifice Area (cm ²)
26	501DA26	14.3	1.5	12.8	24.8	26.2	4.59
28	501DA28	15.4	1.5	13.9	26.8	28.2	5.35
MITRAL							
16	501DM16	9.3	4.8	4.5	14.8	16.2	1.55
18	501DM18	10.3	5.3	5.0	16.8	18.2	2.02
20	501DM20	11.3	5.8	5.5	18.8	20.2	2.56
22	501DM22	12.3	6.3	6.1	20.8	22.2	3.17
24	501DM24	13.3	6.7	6.6	22.8	24.2	3.84
26	501DM26	14.3	7.2	7.1	24.8	26.2	4.59
28	501DM28	15.4	7.7	7.7	26.8	28.2	5.35

Table 4. MEDTRONIC OPEN PIVOT APeX HEART VALVE SPECIFICATIONS

Valve Size (mm)	Model Number	Overall Profile Height (mm)	Inflow Profile Height (mm)	Outflow Profile Height (mm)	Inside Diameter I.D. (mm)	Tissue Annulus Diameter (mm)	Geometric Orifice Area (cm ²)
AORTIC							
16	503DA16	9.3	1.5	7.8	14.8	16.2	1.55
18	503DA18	10.3	1.5	8.8	16.8	18.2	2.02
20	503DA20	11.3	1.5	9.8	18.8	20.2	2.56
22	503DA22	12.3	1.5	10.8	20.8	22.2	3.17
24	503DA24	13.3	1.5	11.8	22.8	24.2	3.84
26	503DA26	14.3	1.5	12.8	24.8	26.2	4.59
28	503DA28	15.4	1.5	13.9	26.8	28.2	5.35

Table 5. MEDTRONIC OPEN PIVOT STANDARD HEART VALVE SPECIFICATIONS (Japan Only)

Valve Size (mm)	Model Number	Overall Profile Height (mm)	Inflow Profile Height (mm)	Outflow Profile Height (mm)	Inside Diameter I.D. (mm)	Tissue Annulus Diameter (mm)	Geometric Orifice Area (cm ²)
AORTIC							
19	500FA19	9.3	2.4	6.9	14.8	19.5	1.55
21	500FA21	10.3	2.9	7.4	16.8	21.5	2.02
23	500FA23	11.3	3.4	7.9	18.8	23.5	2.56
25	500FA25	12.3	3.4	9.2	20.8	25.5	3.17
27	500FA27	13.3	3.4	9.9	22.8	27.5	3.8
29	500FA29	14.3	3.9	10.4	24.8	29.5	4.59
31	500FA31	15.4	4.6	10.8	26.8	31.5	5.35
MITRAL							
19	500DM19	9.3	3.0	6.3	14.8	19.5	1.55
21	500DM21	10.3	3.0	7.3	16.8	21.5	2.02
23	500DM23	11.3	3.0	8.3	18.8	23.5	2.56
25	500DM25	12.3	3.3	9.0	20.8	25.5	3.17
27	500DM27	13.3	3.7	9.6	22.8	27.5	3.84
29	500DM29	14.3	3.7	10.6	24.8	29.5	4.59
31	500DM31	15.4	3.7	11.7	26.8	31.5	5.35
33	500DM33	15.4	3.7	11.7	26.8	33.5	5.35

1.1. Orifice and Leaflets

Based on its history of successful use in cardiac valve prostheses, pyrolytic carbon was chosen for the leaflets and orifice. The orifice is made of 100% pyrolytic carbon. The leaflets are made of pyrolytic carbon coated over a graphite substrate that is impregnated with tungsten (20%) for radiopacity.

Each leaflet opens to a maximum of 85° with a maximum travel arc of 60° in all sizes.

The low profile of the prosthesis results from the bileaflet design with the pivot areas located entirely within the orifice ring. This method of construction minimizes the overall height of the valve. Leaflet movement is controlled by pivot guides located on the inner circumference of the orifice ring. The leaflets move smoothly within these guides. There are no recesses or cavities in the pivot area.

The Medtronic Open Pivot Heart Valve can be rotated for optimal positioning once implanted in either the aortic or mitral position.

The strengthening band, lock rings, and lockwire surrounding the pyrolytic carbon orifice are made of titanium or a cobalt-chromium alloy (MP35N), for added strength and radiopacity.

1.2. Sewing Cuff

The sewing cuff of the Medtronic Open Pivot Heart Valve is made from double velour polyester that has a history of use in cardiac valve prostheses.

Separate sewing cuffs are available for mitral or aortic implantation (Figure 1, Figure 2, and Figure 3). Cuff markers are provided to assist in the uniform placement of sutures around the valve annulus.

- aortic cuff – 3 markers at 120° intervals
- mitral cuff – 4 markers at 90° intervals

2. Indications for Use

The Medtronic Open Pivot Heart Valve is intended for use as a replacement valve in patients with diseased, damaged, or malfunctioning heart valves. This device may also be used to replace a previously implanted prosthetic heart valve.

3. Contraindications

The Medtronic Open Pivot Heart Valve is contraindicated in patients unable to tolerate anticoagulation therapy.

4. Warnings and Precautions

4.1. Warnings

- This device was designed for single patient use only. Do not reuse, reprocess, or resterilize this product. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or create a risk of contamination of the device, which could result in patient injury, illness, or death.
- Use only Medtronic Open Pivot Heart Valve holders, rotators, and sizers.
- Leaflet motion must be tested only with the blue leaflet actuator enclosed in each valve package. Do not test leaflet motion or touch the valve with metal instruments.
- Never pass any catheter or instrument through the valve.
- Do not handle the valve except with the accessories provided. Contact with the valve by any other instruments may cause damage to the leaflets and orifice or result in imperfections on blood contact areas. Scratching the surface may lead to a loss of structural integrity.
- When securing the valve in place, suture needles should be passed through the outer half of the sewing cuff and suture ends should be cut short after the knots are tied.
- It is suggested that only taper point needles be used for suturing the cuff as taper cut or other cutting needles may cut the cuff fibers.
- Force should never be applied to the valve leaflets.
- Before closing the heart, test for leaflet motion with the blue leaflet actuator. If necessary, rotate the valve to avoid abnormal residual pathology which could interfere with leaflet motion.
- Persons allergic to cobalt, chromium, or nickel may suffer an allergic reaction to this device.

4.2. Precautions

- The outer carton label serves as a tamper-evident package seal. After removing the double barrier trays from the shipping package, examine them carefully to verify that the seals and trays are intact. If any seal is damaged or missing, do not use the valve.
- If the Use By date has expired, the device should not be used.

4.3. Patient Counseling Information

- Patients with prosthetic valves who undergo dental or potentially bacteremic procedures should receive antibiotic prophylaxis.
- All mechanical prosthetic heart valves produce sound as a function of their operation. Patients should be advised of this prior to implantation.
- Patients should be encouraged to carry the Implanted Device Identification Card, provided by Medtronic, with them at all times.

5. Potential Adverse Events

Adverse events potentially associated with the use of prosthetic aortic heart valves include, but are not limited to:

- angina
- cardiac arrhythmia
- endocarditis
- hemolysis
- hemolytic anemia
- hemorrhage, anticoagulation-related
- myocardial infarction
- leaflet entrapment (impingement)
- nonstructural dysfunction
- pannus
- perivalvular leak
- transvalvular regurgitation
- structural dysfunction
- thrombosis
- stroke
- thromboembolism

It is possible that these complications could lead to:

- reoperation
- explantation

- permanent disability
- heart failure
- death

6. Instructions for Use

6.1. Sterilization Information

The Medtronic Open Pivot Heart Valve and package contents have been steam sterilized within the trays. The valve is ready for use when received. **Do not use the valve if:**

- the Use By date has expired
- the valve container (inner packaging) is damaged
- the sterility barrier is broken

The sizers and reusable handles are supplied NONSTERILE and must be cleaned and sterilized prior to use. Refer to the Medtronic Open Pivot Accessories Instructions for Use for cleaning and sterilization instructions.

Warning: If the valve is removed from its container but not used, it must not be repackaged or resterilized.

6.2. Preparation For Implant

Circulating Nurse:

1. Remove the double barrier trays from the shipping box along with the Patient Registration Form.
2. Verify that the valve size, type, and serial number located on the tray label matches the information on the box label.
Note: Verify that the valve is attached to the valve holder BEFORE opening the outer tray. If the valve is not attached to the valve holder, DO NOT use the valve. DO NOT attempt to reassemble the valve onto the valve holder.
3. Hold the outer tray and gently peel the lid back using sterile technique until the lid has been completely removed from the tray.
4. Record the serial number of the valve in the patient's record using the stickers provided on the Patient Registration Form.

Scrub Nurse/Surgeon:

1. Remove the inner tray from the outer tray by inserting a gloved finger into the open corner and lifting the inner tray.
2. Holding the inner tray from the bottom, grasp the inner tray lid tab and peel back until it has been completely removed.
3. Pick up the retainer using the textured areas. The inner tray retainer has textured areas on both ends.
4. Discard or recycle the retainer after it is removed.
Note: DO NOT remove the valve and holder from the tray.
5. Check that the valve holder is secured to the handle/rotator. There should be no visible screw threads at the connection between the colored handle/rotator and the valve holder when secured properly. Reference Figure 5 and Figure 6 for the position of an aortic valve and a mitral valve on the holder, respectively. DO NOT use the valve if it is not properly secured to the holder.
6. Remove the blue leaflet actuator and set aside for later use.
7. Grasp the handle/rotator in the center and gently lift the valve and valve holder out of the inner tray.

Caution: Do not remove the valve holder from the valve unless directed by the implanting surgeon.

Caution: Avoid touching the valve unnecessarily to prevent structural damage and to minimize the potential for thrombus formation.

6.3. Surgical Information

Valve Sizing

- Prior to surgery, a complete set of Medtronic Open Pivot Heart Valve sizers and reusable handle must be cleaned and sterilized per the Accessories Instructions for Use and made available in the operating room (Figure 7, Figure 8, and Figure 9).
- Sizers that are damaged or cracked must not be used.
- Use only Medtronic Open Pivot Heart Valve sizers to measure the patient's tissue annulus. The Medtronic Open Pivot Heart Valve sizer has a Nitinol shaft that can be shaped allowing the same sizer to be used for mitral or aortic annulus estimations.
- The sizer with the diameter marked in millimeters has a ring which allows easy tissue annulus measurement. The sizer ring should pass easily through the patient's tissue annulus with minimal resistance. Do not force the sizer through the tissue annulus. Valve oversizing or undersizing should be avoided.
- AP Series (AP, APex, AP360) valves are intended for supra-annular placement.

Valve Holder

- The Medtronic Open Pivot Heart Valve is premounted on a holder and connected to a color-coded handle/rotator. Reference Figure 5 and Figure 6 for the position of an aortic valve and a mitral valve on the holder. The BLOOD FLOW in either the aortic or mitral position is ALWAYS INTO THE STRAIGHT EDGE of the leaflets (Figure 1, Figure 2, and Figure 3).
- The bendable, reusable handle (Figure 9) may be used if the surgeon desires a different angle for the valve holder.

Cuff Suturing

- Suture techniques may vary according to the preference of the implanting surgeon and the anatomy of the patient. The medical literature reports several preferred suturing methods.
- Sutures should be placed in the outer half of the cuff, particularly on the AP models. Deep needle placement will contact the strengthening band inside the cuff and impede smooth placement of sutures.
- After sutures are placed in the valve cuff and the valve is seated in the tissue annulus, release the valve from the valve holder. Slide a scalpel down the holder slot to cut the green suture and release the valve holder. The green suture and valve holder will remain attached as the handle/rotator is carefully withdrawn. The suture knots can be tied at this point. Long suture ends must be avoided.

Caution: Testing of leaflet motion must only be done with the enclosed blue leaflet actuator. Use of another instrument may damage the valve. If the leaflets do not move freely due to contact with tissue, the valve orifice may be rotated to a more optimal position. Rotate the orifice only after the suture knots have been tied, securing the cuff to the tissue annulus.

Rotatability

- The valve orifice and leaflets may be rotated in situ by using the handle/rotator supplied in the package. Freedom of rotation should be checked before implanting the valve by holding the valve cuff and gently turning the valve handle/holder. Each valve size and model has a specific rotator. If desired, carefully remove the rotator from the handle by turning the rotator counterclockwise. The rotator seats properly into the valve orifice when the radial ridge on the handle side of the rotator head is aligned with the straight edge of the leaflets. The rotator head is designed to contact the flat surfaces of the orifice near the pivot areas. The orifice may be rotated in a clockwise or counterclockwise direction. After orientation of the valve orifice in the patient's annulus, leaflet mobility may be determined with the blue leaflet actuator.

Caution: Use only the rotator supplied with the device. Do not use any other instrument to rotate the valve. Use of instruments other than those provided may cause damage to the orifice and leaflets and may result in valve malfunction.

Caution: Extreme care must be exercised when handling the valve holder assembly to avoid structural damage to the valve. Inspect each accessory before each use. Cracked or damaged accessories must not be used.

7. Postoperative Information

7.1. Anticoagulation

The medical literature indicates that patients implanted with mechanical heart valves should be routinely maintained on anticoagulants to reduce the risk of valve thrombosis and thromboembolism.

7.2. Diagnostic Imaging

The Medtronic Open Pivot Heart Valve provides excellent x-ray images due to the radiopaque strengthening band around the orifice and 20% tungsten in each leaflet substrate.

7.3. Magnetic Resonance Imaging (MRI) Testing

The Medtronic Open Pivot Heart Valve (Models 500, 501, 503, and 505) was determined to be MR-conditional according to the terminology specified in the American Society for Testing and Materials (ASTM) International, Designation: F2503-05. Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment. ASTM International, 100 Barr Harbor Drive, PO Box C700, West Conshohocken, Pennsylvania, 2005.

Nonclinical testing demonstrated that the Medtronic Open Pivot Heart Valve (Models 500, 501, 503, and 505) is MR Conditional. A patient with this device can be scanned safely immediately after placement under the following conditions:

7.4. Static Magnetic Field

- Static magnetic field of 3 Tesla or less
- Maximum spatial gradient magnetic field of 720 Gauss/cm or less

7.5. MRI-Related Heating

In nonclinical testing, the Medtronic Open Pivot Heart Valve (Models 500, 501, 503, and 505) produced the following temperature rise during MRI performed for 15 minutes in the 3 T (3 T/128-MHz, Excite, HDX, Software 14X.M5, General Electric Healthcare, Milwaukee, WI) MR system: Highest temperature change: +1.6°C.

Therefore, the MRI-related heating experiments for the Medtronic Open Pivot Heart Valve (Models 500, 501, 503, and 505) at 3 T using a transmit/receive RF body coil at an MR system reported whole body averaged SAR of 2.9 W/kg (ie, associated with a calorimetry measured whole body averaged value of 2.7 W/kg) indicated that the greatest amount of heating that occurred in association with these specific conditions was equal to or less than +1.6°C.

7.6. Artifact Information

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the Medtronic Open Pivot Heart Valve.

Models 500, 501, and 503:				
Pulse Sequence	T1-SE	T1-SE	GRE	GRE
Signal Void Size	1,213 mm ²	486 mm ²	2,189 mm ²	1,369 mm ²
Plane Orientation	Parallel	Perpendicular	Parallel	Perpendicular

Model 505:				
Pulse Sequence	T1-SE	T1-SE	GRE	GRE
Signal Void Size	2,257 mm ²	724 mm ²	4,244 mm ²	2,327 mm ²
Plane Orientation	Parallel	Perpendicular	Parallel	Perpendicular

Therefore, optimization of MR imaging parameters to compensate for the presence of this device may be necessary.

8. How Supplied

8.1. Packaging

The Medtronic Open Pivot Heart Valve is packaged in sterile transparent double barrier trays. The valve is premounted on a holder and attached to a handle/rotator. All package labeling and handle/rotator accessories are color-coded and marked:

- Aortic – Green
- Mitral - Red

The valve package (Figure 4) includes:

- Outer Tray
- Inner Tray
- Inner Tray Retainer
- Valve, Valve Holder, Handle/Rotator Assembly
- Instructions for Use
- Patient Registration Form
- Blue Leaflet Actuator

All packaging materials are recyclable. Manage disposal in accordance with local statutory regulations.

8.2. Storage

For maximum protection and product identification, it is recommended that the valve be stored in its original packaging. The storage environment should be clean, cool, and dry. The sterility and nonpyrogenicity of the valve is validated to remain unaffected until the Use By date printed on the product, provided the seals and containers are not opened or damaged.

8.3. Returned Product Information

For detailed information on the Medtronic Return Product Policy, please contact your local representative.

9. Accessories

A color-coded rotator is included with each valve (Figure 5 and Figure 6). Only use the supplied rotator because the rotator is matched to the valve size and position.

The Medtronic Open Pivot Heart Valve sizer sets and bendable handles are intended for use only with Medtronic Open Pivot Heart Valves (Figure 7, Figure 8, and Figure 9). Medtronic Open Pivot Accessories are provided NONSTERILE and must be cleaned and sterilized before each use.

The blue leaflet actuator supplied with each package is used to test for free leaflet motion after the valve has been implanted.

10. Patient Information

10.1. Registration Information

Note: Patient registration does not apply in countries where patient privacy laws conflict with providing patient information, including countries from the EU. A patient registration form is included in each device package. After implantation, please complete all requested information. The serial number may be found on the tray and box label. Return the original form to the Medtronic address indicated on the form, and provide the temporary identification card to the patient prior to discharge.

An Implanted Device Identification Card is provided to the patient. The card contains the name and telephone number of the patient's physician, as well as information that medical personnel would require in the event of an emergency.

11. DISCLAIMER OF WARRANTY

THE FOLLOWING DISCLAIMER OF WARRANTY APPLIES TO CUSTOMERS OUTSIDE THE UNITED STATES:

ALTHOUGH THE MEDTRONIC OPEN PIVOT HEART VALVES, HEREAFTER REFERRED TO AS "PRODUCT", HAVE BEEN CAREFULLY DESIGNED, MANUFACTURED AND TESTED PRIOR TO SALE, THE PRODUCT MAY FAIL TO PERFORM ITS INTENDED FUNCTION SATISFACTORILY FOR A VARIETY OF REASONS. THE WARNINGS CONTAINED IN THE PRODUCT LABELING PROVIDE MORE DETAILED INFORMATION AND ARE CONSIDERED AN INTEGRAL PART OF THIS DISCLAIMER OF WARRANTY. MEDTRONIC, THEREFORE, DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT. MEDTRONIC SHALL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, OR FAILURE OF THE PRODUCT, WHETHER THE CLAIM IS BASED ON WARRANTY, CONTRACT, TORT, OR OTHERWISE.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this DISCLAIMER OF WARRANTY is held by any court of competent jurisdiction to be illegal, unenforceable, or in conflict with applicable law, the validity of the remaining portion of the DISCLAIMER OF WARRANTY shall not be affected, and all rights

and obligations shall be construed and enforced as if this DISCLAIMER OF WARRANTY did not contain the particular part or term held to be invalid.

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7.3. MR-testing

Medtronic Open Pivot hjerteklaff (modell 500, 501, 503 og 505) er fastslått å være MR-betinget i henhold til terminologien angitt i American Society for Testing and Materials (ASTM) International, Designation: F2503-05. Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment. ASTM International, 100 Barr Harbor Drive, PO Box C700, West Conshohocken, Pennsylvania, 2005.

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Medtronic Open Pivot Kalp Kapağı (500, 501, 503 ve 505 Modelleri), American Society for Testing and Materials (ASTM) International tarafından şu Tayinde belirtilen terminolojiye uygun biçimde, MR koşullu olarak belirlenmiştir: F2503-05. Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment. ASTM International, 100 Barr Harbor Drive, PO Box C700, West Conshohocken, Pennsylvania, 2005.

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© 2012 Medtronic
Medtronic Open Pivot™, AP360®, APex® and AP® are trademarks of Medtronic, Inc.
Medtronic, Inc.