## MINIMAX<sup>®</sup>, MINIMAX PLUS<sup>®</sup>

Filtered Hardshell Reservoir / Réservoir rigide à filtre / Gefiltertes Hartschalen-Reservoir / Reservorio rígido filtrado / Gefilterd reservoir met harde schaal / Serbatoio rigido filtrato / Filtrerad hård reservoir / Ubøjelig beholder med filter / Filtreli Hardshell Rezervuar  
1316

Hollow Fiber Oxygenator with Plasma Resistant Fiber / Oxygénateur à fibre creuse plasma-résistante / Oxygenator mit plasmaresistenter Hohlfaser / Oxigenador de fibra hueca resistente al plasma / Hollow-fiber-oxygenator met plasmabestendige vezel (PRF) / Ossigenatore a fibre cave plasmaresistenti / Hålfiberoxygenator med plasmaresistent fiber / Hollow-fiber oxygenator med plasmaresistent fiber / Plazmaya Dayanıklı Lifli Hollow Fiber Oksijenatör  
3381/CB3381\*

Hollow Fiber Oxygenator with Plasma Resistant Fiber (PRF) with Filtered Hardshell Venous Reservoir / Oxygénateur à fibre creuse plasma-résistante muni d'un réservoir veineux rigide à filtre / Oxygenator mit plasmaresistenter Hohlfaser (PRF Plasma Resistant Fiber) mit gefiltertem venösem Hartschalen-Reservoir / Oxigenador de fibra hueca resistente al plasma con reservorio venoso rígido filtrado / Hollow-fiber-oxygenator met plasmabestendige vezel (PRF) met gefilterd veneus reservoir met harde schaal / Ossigenatore a fibre cave plasmaresistenti (PRF) con serbatoio venoso rigido dotato di filtro / Hålfiberoxygenator med plasmaresistent fiber (PRF) med filtrerad hård venreservoir / Hollow-fiber oxygenator med plasmaresistent fiber (PRF) med ubøjelig venøs beholder med filter / Filtreli Hardshell Venöz Rezervuara Sahip Plazmaya Dayanıklı Lifli (PRF) Hollow Fiber Oksijenatör  
3301

Instructions for Use • Mode d'emploi • Gebrauchsanweisung • Instrucciones de uso  
Gebruiksaanwijzing • Istruzioni per l'uso • Bruksanvisning • Brugsanvisning  
Kullanım Talimatları

 Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

- \*With Carmeda® Bioactive Surface
- \*Avec surface bioactive Carmeda®
- \*Mit Carmeda® BioActive Surface
- \*Con superficie bioactiva Carmeda®
- \*Met Carmeda® bioactieve oppervlaktelaag
- \*Con superficie bioattiva Carmeda®
- \*Med Carmeda® BioActive Surface.
- \*Med Carmeda® BioActive overflade
- \*Carmeda® Biyo Aktif Yüzeye sahiptir

Medtronic®, Minimax®, and Minimax Plus® are registered trademarks of Medtronic, Inc.  
Medtronic®, Minimax® et Minimax Plus® sont des marques déposées de Medtronic, Inc.  
Medtronic®, Minimax® und Minimax Plus® sind eingetragene Marken von Medtronic, Inc.  
Medtronic®, Minimax® y Minimax Plus® son marcas comerciales registradas de Medtronic, Inc.  
Medtronic®, Minimax® en Minimax Plus® zijn geregistreerde handelsmerken van Medtronic, Inc.  
Medtronic®, Minimax® e Minimax Plus® sono marchi registrati di Medtronic, Inc.  
Medtronic®, Minimax® och Minimax Plus® är registrerade varumärken som tillhör Medtronic, Inc.  
Medtronic®, Minimax® og Minimax Plus® er varemærker tilhørende Medtronic, Inc.  
Medtronic®, Minimax® ve Minimax Plus® Medtronic, Inc. şirketinin tescilli ticari markalarıdır.

**Explanation of symbols on package labeling / Explication des symboles des étiquettes sur l'emballage / Erläuterung der Symbole auf der Verpackung / Explicación de los símbolos que aparecen en el etiquetado del envase / Verklaring van de symbolen op de verpakkingslabels / Spiegazione dei simboli sulle etichette della confezione / Förklaring av symboler på förpackningsetiketten / Forklaring af symboler på emballagens mærkater / Ambalaj etiketi üzerindeki sembollerin açıklaması**

Refer to the outer package label to see which symbols apply to this product / Se référer à l'étiquette sur l'emballage extérieur pour savoir quels symboles s'appliquent à ce produit / Welche Symbole für dieses Produkt zutreffen, entnehmen Sie bitte dem Etikett auf der Verpackung / Consulte la etiqueta del envase exterior para comprobar qué símbolos son aplicables a este producto / Controleer het label op de buitenverpakking om te zien welke symbolen op dit product van toepassing zijn / Fare riferimento all'etichetta sulla confezione esterna per conoscere i simboli applicati a questo prodotto / Se etiketten på ytterförpackningen för de symboler som gäller denna produkt / Se mærkaten på den ydre emballage for de symboler, der gælder for dette produkt / Hangi sembollerin bu ürün için geçerli olduğunu görmek için ambalaj etiketine başvurun

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
<sup>1</sup> These results are based on ISO/AAMI standard venous inlet conditions. All information referenced in Section 6.3(t) of the ISO/AAMI guidelines is available upon request.

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### Description

#### Models: 3381, CB3381

The Minimax Plus Hollow Fiber Oxygenator with Plasma Resistant Fiber, with or without Carmeda®<sup>1</sup> BioActive Surface<sup>2</sup> (heparin), is a single-use, sterile, nonpyrogenic gas-exchange device with a self-contained, venous-side heat exchanger for regulating blood temperature. The microporous, hollow fibers are wound around a core and encased in a plastic outer shell. Oxygen flows through the hollow fibers, and blood flows from the heat exchanger around the hollow fibers. Gas exchange occurs by diffusion across the hollow fiber membrane. A gas outlet is provided for scavenging and a gas vent prevents overpressurization.

The Minimax Plus Hollow Fiber Oxygenator with Plasma Resistant Fiber with Carmeda BioActive Surface (CB3381) is coated with a nonleaching BioActive Surface, which provides thromboresistant blood contact surfaces.

#### Model: 1316

The Minimax Filtered Hardshell Reservoir is a single-use, sterile, nonpyrogenic device consisting of an outer housing with a defoamer, and a 20 micron depth filter which are used for defoaming and filtering venous and cardiotomy blood. Venous blood enters via a 1.0 cm (3/8 in) venous inlet port. Cardiotomy blood enters via the 0.6 cm (1/4 in) inlet ports on the top of the vessel. The blood is pumped to an oxygenator via the 1.0 cm (3/8 in) outlet port at the base of the reservoir.

An adapter is provided with the Minimax Filtered Hardshell Reservoir to accommodate changing the 1.0 cm (3/8 in) venous inlet and outlet ports to 0.6 cm (1/4 in) ports.

The Minimax Filtered Hardshell Reservoir may be used with the Minimax Plus Hollow Fiber Oxygenator or may be used with other membrane oxygenators.

### Specifications

#### Oxygenator Models 3381, CB3381

Static Priming Volume	Approximately 149 mL
Recommended Blood Flow Rate	0.5 to 2.3 Liters per minute
Maximum Water Pressure	40 psi/2068 mm Hg/275.8 kPa
Maximum Blood Path Pressure	6 psi/310 mm Hg/41.3 kPa
Oxygen Transfer Rate	see Figure 5
Carbon Dioxide Transfer Rate	see Figure 6
Heat Exchanger Performance	see Figure 7
Blood Side Pressure Drop	see Figure 8
Membrane	Microporous hollow fiber (0.8 m <sup>2</sup> )
Venous Inlet Port	0.6 cm (1/4 in)
Arterial Outlet Port	0.6 cm (1/4 in)
Arterial Sample Port	0.6 cm (1/4 in)
Venous Air Director	Female luer port
Gas Ports	0.6 cm (1/4 in)
Water Ports	1.3 cm (1/2 in) quick disconnect
Arterial Temperature Port	YSI Series 400

#### Hardshell Reservoir Model 1316

Recommended minimum operating level	150 mL
Maximum operating volume	2000 mL
Cardiotomy Filtration	20 micron
Suction Inlets	4 @ 0.6 cm (1/4 in)
Filtered Quick Prime	1 @ 0.6 cm (1/4 in)
Filtered Luer	1
Non-Filtered Luers	3
Venous Inlet (with Luer and Temperature Probe)	1 @ 1.0 cm (3/8 in)

<sup>1</sup> Carmeda® is a registered trademark of Carmeda AB.

<sup>2</sup> This product is coated with the Carmeda® BioActive Surface under license from Carmeda AB and is licensed for use only as part of an extracorporeal blood circulation system or circuit which includes an oxygenator or blood pump.

Outlet	1 @ 1.0 cm (3/8 in)
Recirculation Port	1 @ 0.6 cm (1/4 in)

## Indications

The Minimax Plus Hollow Fiber Oxygenator with Plasma Resistant Fiber, with or without Carmeda BioActive Surface, is indicated for use in procedures requiring the extracorporeal oxygenation of and carbon dioxide removal from blood. It is designed to operate at a maximum blood flow rate of 2.3 liters per minute for periods of up to 6 hours.

The Minimax Filtered Hardshell Reservoir is indicated for use in procedures requiring the storage and filtration of blood during extracorporeal circulation at venous blood flow rates of up to 4 liters per minute and cardiotomy filtration rates of up to 2 liters per minute.

## Contraindications

This device used for any other purposes than for the indicated intended use is the responsibility of the user.

## Warnings

Read all Warnings, Precautions, and Instructions for Use carefully prior to use. **Failure to read and follow all instructions, or failure to observe all stated warnings, could cause serious injury or death to the patient.**

- The oxygenator/hardshell reservoir should be used only by persons thoroughly trained in cardiopulmonary bypass procedures. Operation of each device requires constant supervision by qualified personnel for patient safety.
- Each oxygenator has been carefully manufactured, tested, and packaged; however, the state of the art has not developed to the point that Medtronic is able to ensure that the oxygenator will not leak, crack, or fail during use. Perfusion must be carefully and constantly monitored.
- A standby oxygenator should be available during perfusion. If an oxygenator must be replaced during perfusion, follow the instructions in the "Emergency Oxygenator Replacement" section.
- The fluid pathway is sterile and nonpyrogenic. Inspect each package and device prior to use. Do not use if the package is opened or damaged or the protective caps are not in place.
- Each device is intended for single use only. Do not reuse or resterilize. EO sterilized.
- Blood pressure must be maintained higher than gas phase pressure during perfusion to avoid the formation of gaseous microemboli.
- Do not let anesthetic agents such as isoflurane come in concentrated contact with this device. These agents may jeopardize its structural integrity.
- Do not occlude or obstruct the gas vent port.
- If the heart is vented back to the Minimax Filtered Hardshell Reservoir, care must be taken to prevent backflow through the ventricular vent.
- Do not use a vacuum source to generate suction in the reservoir.
- Ensure that the reservoir operating level is above the top of the membrane oxygenator.
- Heat exchanger water temperature should not exceed 42°C (107.6°F). High temperature may result in blood damage.
- The temperature gradient between blood and heat exchange water should not exceed 15°C (59°F).
- If air bubbles and/or leaks are observed during priming and/or operation, these conditions may result in air embolism to the patient and/or fluid loss. The extracorporeal circuit must be continually monitored. Do not use the device if these conditions are observed.
- In any extracorporeal circuit, the use of a prebypass and arterial line filter should be considered.
- During and after priming, always keep the fluid level in the venous reservoir above the oxygenator.
- All gas emboli must be cleared from the extracorporeal circuit before initiating bypass. Gas emboli are hazardous to the patient.
- Do not allow air to enter the circuit.
- Ensure that the recirculation line is primed/clamped prior to initiating bypass or reinstating bypass and that it remains clamped at all times while on bypass.
- **Caution:** Federal law (USA) restricts this device to sale by or on the order of a physician.
- **Warning:** Manufactured with HCFC-225, a substance which harms public health and environment by destroying ozone in the upper atmosphere.

## Precautions

- Refer to package labeling for storage temperature requirements.
- Use aseptic techniques in all procedures.
- Monitor oxygen concentration throughout perfusion.

- A strict anticoagulation protocol should be followed, and anticoagulation should be routinely monitored during all procedures. The benefits of extracorporeal support must be weighed against the risk of systemic anticoagulation and must be assessed by the prescribing physician.
- Additional warnings and precautions applicable to specific procedures can be found at appropriate places in the Instructions for Use.

## Instructions for Use

### Setup

#### 1. Water and Gas Lines

- a. Place the oxygenator securely in its holder. Be sure the venous reservoir is higher than the oxygenator.
 

**Note:** When the Minimax Filtered Hardshell Reservoir is used in conjunction with the Minimax Plus Hollow Fiber Oxygenator, the holder is attached to the Minimax Oxygenator holder. When used separately, the holder is attached to a holder adapter.
- b. Connect the water inlet and outlet lines to the heat exchanger water ports. Start water flow and check for leaks from the water compartment.
 

**Caution:** If leaks are observed, do not use. Set up the standby oxygenator.

**Caution:** Heat exchanger water temperature should not exceed 42°C (107.6°F). High temperature may result in blood damage.

**Caution:** The temperature gradient between blood and heat exchange water should not exceed 15°C (59°F).
- c. Connect the gas line to the gas inlet port. If a gas blender is used, connect the line from the gas blender to the gas inlet port. Refer to the gas blender manufacturer's instructions for use.
- d. Remove the cap from the gas outlet port.

#### 2. Blood Lines used with Venous Reservoir Bag

- a. Connect the venous line to the venous reservoir inlet.
- b. Connect the pump header line to the venous reservoir outlet and to the oxygenator venous inlet port.
- c. Connect the drain line between the cardiotomy reservoir and the venous reservoir bag.
- d. Connect the arterial line to the oxygenator blood outlet port.
- e. Connect the oxygenator recirculation line to the recirculation port of the venous reservoir inlet connector.
- f. Place a 3-way stopcock on the venous air director on top of the oxygenator. Close the stopcock. A male/male luer line may be attached between the stopcock and a reservoir for venting of air.
- g. Securely band all connections.

#### 3. Blood Lines used with Hardshell Venous Reservoir

- a. Remove the accessory pack from the package and examine the contents for any obvious visual defects. The contents of each pack will include:
  - **Elbow Connector (1 each)**

This elbow connector with its preattached 1.0 cm x 1.0 cm x 0.6 cm (3/8 in x 3/8 in x 1/4 in) connector with temperature probe and luer port should be placed on the 1.0 cm (3/8 in) venous inlet port of the hardshell reservoir at an orientation of choice.
  - **1.0 cm – 0.6 cm (3/8 in – 1/4 in) Adapters (2 each)**

The adapters may be used to convert the inlet or outlet ports on the reservoir from 1.0 cm (3/8 in) to 0.6 cm (1/4 in).
  - **Recirculation Line (1 each)**

This 0.6 cm (1/4 in) (ID) line may be used for recirculation between the Minimax Filtered Hardshell Reservoir and a membrane oxygenator. Follow the oxygenator manufacturer's instructions for recirculation.
  - **Sampling Manifold (1 each)**

A 3-gang manifold of stopcocks is provided. Male/male luer lines are preattached, along with a one-way duckbill valve. The one-way duckbill valve is designed to prevent the reverse flow of fluid. This manifold line may be used for the sampling of arterial and venous blood during cardiopulmonary bypass. Care should be taken to ensure that arterial blood flows in the same direction as the arrows located near the duckbill valve. The manifold line should be primed and debubbled prior to bypass.

**Note:** The venous luer port should be used specifically for the addition of any drug to the reservoir.
- b. Remove the yellow gas vent port cap located on the top of the reservoir.
- c. Attach the venous line to the 1.0 cm (3/8 in) venous inlet port or to the 0.6 cm (1/4 in) venous adapter port.

- d. Connect the pump header line to the Minimax Filtered Hardshell Reservoir 1.0 cm (3/8 in) outlet port or to the 0.6 cm (1/4 in) adapter outlet port. Connect the other end to the inlet port of the oxygenator taking care to ensure proper flow direction.
- e. Connect any suction port required.
- f. Connect the 0.6 cm (1/4 in) recirculation line to the 0.6 cm (1/4 in) recirculation port on the venous inlet connector of the reservoir. If the recirculation line is not used, a 0.6 cm (1/4 in) dead end plug must be placed on the 0.6 cm (1/4 in) port of the venous inlet connector.
- g. Follow manufacturer's directions for use when connecting tubing to the oxygenator.
- h. Securely band all connections.

## Priming

**Note:** Before priming, a CO<sub>2</sub> flush is recommended.

### 1. A/V Loop

- a. Ensure the arterial and venous lines are connected to form an A/V loop for priming and recirculation.
- b. Ensure the following lines are clamped: the venous return line at the venous reservoir; the arterial line at the oxygenator; the blood outlet line at the venous reservoir; the recirculation line; and the cardiotomy drain line.

### 2. Filling the Reservoir

- a. The Minimax Filtered Hardshell Reservoir may be primed by gravity via any of the filtered inlet ports located on top of the reservoir.
- b. Fill the reservoir with priming fluid through the quick prime port.
- c. If a venous reservoir bag is used, unclamp the priming reservoir drain line and fill the venous reservoir with the priming fluid. Purge air from the venous reservoir before recirculating.
- d. Prime the oxygenator and circuit per manufacturer's instructions for use.

**Warning:** During and after priming, always keep the fluid level in the venous reservoir above the oxygenator.

### 3. Recirculation

- a. Place the tubing in the pump head if using a roller pump.
- b. Unclamp the recirculation line between the oxygenator and the venous reservoir.
- c. Turn on the blood pump to 50 mL/min and recirculate the priming fluid through the venous reservoir, blood pump, and oxygenator.
- d. Open the venous air director line at the top of the oxygenator.
- e. Slowly increase pump speed to 500 mL/min.
- f. Close the venous air director.

### 4. Prime the A/V Loop

- a. Stop the blood pump.
- b. Clamp the recirculation line.
- c. If necessary, add more priming fluid to the venous reservoir.
- d. Unclamp the venous line.
- e. As you start the blood pump at a slow speed, unclamp the arterial line and purge the air from the loop into the venous reservoir vent.
- f. Through the venous air director, purge any air that enters the oxygenator.
- g. When all air is purged from the A/V loop, slowly increase blood pump flow to 1 liter/minute and recirculate priming fluid through the entire circuit; if desired, start the flow of ventilating gas.

**Warning:** All gas emboli must be cleared from the extracorporeal circuit before initiating bypass. Gas emboli are hazardous to the patient.

### 5. Pre-warm the Priming Fluid

- a. Start heat exchanger water flow to pre-warm the priming fluid.
- b. Add cellular prime, if desired, and recirculate.
- c. Inspect the entire system for leaks.

## Perfusion

### 1. Beginning Bypass

Gas flow rate depends on the blood flow rate and venous inlet conditions. A gas-to-blood flow ratio of 1:1 with a minimum FiO<sub>2</sub> of 80% is recommended for initiating bypass. Have a sufficient gas supply for the entire procedure.

- a. Stop the blood pump, and clamp the arterial and venous lines.
- b. To begin bypass, unclamp the venous and arterial lines and gradually increase pump speed to calculated blood flow. Ensure that the operating level in the venous reservoir is maintained above 150 mL.

**Warning:** Do not allow air to enter the circuit.

**Warning:** Ensure that the recirculation line is primed/clamped prior to initiating bypass or reinstating bypass and that it remains clamped at all times while on bypass.

## 2. Controlling Perfusion Parameters

- a. Oxygen Transfer. Transfer is controlled primarily by the concentration of O<sub>2</sub> in the ventilating gas. Increase the O<sub>2</sub> concentration to raise the patient's arterial pO<sub>2</sub> and decrease the O<sub>2</sub> concentration to lower the arterial pO<sub>2</sub>.

**Note:** The use of an in-line O<sub>2</sub> analyzer between the oxygenator and blender should be considered.

- b. Carbon Dioxide Transfer. Carbon dioxide transfer is controlled primarily by gas flow rate. Increase the gas flow rate to increase the amount of CO<sub>2</sub> removed from the blood; decrease the gas flow rate to decrease the amount of CO<sub>2</sub> removed from the blood.

## 3. Coming Off Bypass

- a. Terminate bypass as indicated by individual case and patient condition and perfusion protocol.
- b. After coming off bypass, continue to recirculate through the oxygenator at low flow rates, using either the recirculation line, a reconnected A/V loop or both.

## Emergency Oxygenator Replacement

A standby Minimax Plus Hollow Fiber Oxygenator should always be available during perfusion.

1. Turn off the blood pump and heat exchanger water supply.
2. Clamp the venous line at the venous reservoir.
3. Clamp the venous, arterial, and recirculation lines at the oxygenator.
4. Clamp the heat exchanger lines at the oxygenator and disconnect. Disconnect gas and monitor lines.
5. Disconnect the blood inlet, outlet, and recirculation lines from the oxygenator. Remove the oxygenator from its holder.
6. Place a new oxygenator in the holder. Connect all lines and securely band all connections.
7. Connect the water inlet and outlet lines to the heat exchanger water ports. Start water flow and check for leaks from the water compartment.
8. Add an additional 200 mL of fluid to the priming reservoir for oxygenator priming.
9. Remove the clamps from the priming reservoir drain line, the recirculation line, and the venous line at the oxygenator.

**Warning:** Do not remove the clamp on the venous return line or the clamp on the arterial line.

10. Turn on the blood pump and slowly fill the oxygenator. Open venous air director to purge air and close. Increase blood flow to 500 mL/min through the recirculation line. Vent air bubbles from the venous reservoir.
11. Inspect the entire system for leaks and gas emboli.
12. Remove all clamps from the venous and arterial lines, clamp the recirculation line and continue bypass.

### **Important Notice — Limited Warranty (US Customers Only)**

- A. This **LIMITED WARRANTY** provides the following assurances to the patient who receives a Medtronic MINIMAX PLUS HOLLOW FIBER OXYGENATOR and the MINIMAX FILTERED HARDSHELL RESERVOIR, hereafter referred to as "The Product".
- (1) Should the Product fail to function within normal tolerances due to a defect in materials or workmanship prior to its "Use By" date, Medtronic will, at its option: (a) issue a credit to the purchaser equal to the Purchase price, as defined in Subsection A(2), against the purchase of the replacement Product or (b) provide a functionally comparable replacement Product at no charge.
  - (2) As used herein, Purchase Price shall mean the lesser of the net invoiced price of the original, or current functionally comparable, or replacement Product.
- B. To qualify for the **LIMITED WARRANTY**, these conditions must be met:
- (1) The Product must be used prior to its "Use By" date.
  - (2) The unused portion of the Product must be returned to Medtronic and shall be the property of Medtronic.
  - (3) The Product must not have been altered or subjected to misuse, abuse or accident.
  - (4) The Product must be used in accordance with the labeling and instructions for use provided with the Product.
- C. This **LIMITED WARRANTY** is limited to its express terms. In particular:
- (1) Except as expressly provided by this **LIMITED WARRANTY**, Medtronic is not responsible for any direct, incidental or consequential damages based on any defect, failure or malfunction of the product, whether the claim is based on warranty, contract, tort or otherwise.
  - (2) **THIS LIMITED WARRANTY IS MADE ONLY TO THE PATIENT IN WHOM THE PRODUCT WAS USED. AS TO ALL OTHERS, MEDTRONIC MAKES NO WARRANTY, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY, OR FITNESS FOR A PARTICULAR PURPOSE, WHETHER ARISING FROM STATUTE, COMMON LAW, CUSTOMER OR OTHERWISE. NO EXPRESS OR IMPLIED WARRANTY TO THE PATIENT SHALL EXTEND BEYOND THE PERIOD SPECIFIED IN A(1) ABOVE. THIS LIMITED WARRANTY SHALL BE THE EXCLUSIVE REMEDY AVAILABLE TO ANY PERSON.**
  - (3) The exclusions and limitation set out above are not intended to, and should not be construed so as to contravene mandatory provisions of applicable law. If any part or term of this **LIMITED WARRANTY** is held to be illegal, unenforceable or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of the **LIMITED WARRANTY** shall not be affected, and all rights and obligations shall be construed and enforced as if this **LIMITED WARRANTY** did not contain the particular part or term held to be invalid. This **LIMITED WARRANTY** gives the patient specific legal rights. The patient may also have other rights which vary from state to state.
  - (4) No person has any authority to bind Medtronic to any representation, condition or warranty except this **LIMITED WARRANTY**.

**Note:** This **LIMITED WARRANTY** is provided by Medtronic, Inc. 710 Medtronic Parkway, Minneapolis, MN 55432-5604. It applies only in the United States. Areas outside the United States should contact their local Medtronic representative for exact terms of the **LIMITED WARRANTY**.

### **Important Notice — Limited Warranty (For Countries Outside the US)**

- A. This **LIMITED WARRANTY** provides assurance for the purchaser who receives a Medtronic MINIMAX PLUS HOLLOW FIBER OXYGENATOR and the MINIMAX FILTERED HARDSHELL RESERVOIR, hereafter referred to as the "Product", that, should the Product fail to function to specification, Medtronic will issue a credit, equal to the original Product purchase price (but not to exceed the value of the replacement Product) against the purchase of any Medtronic replacement Product used for that patient.
- The Warnings contained in the product labeling are considered an integral part of this **LIMITED WARRANTY**. Contact your local Medtronic representative to obtain information on how to process a claim under this **LIMITED WARRANTY**.
- B. To qualify for the **LIMITED WARRANTY**, these conditions must be met:
- (1) The Product must be used prior to its "Use By" date.
  - (2) The Product must be returned to Medtronic within 60 days after use and shall be the property of Medtronic.
  - (3) The Product may not have been used for any other patient.
- C. This **LIMITED WARRANTY** is limited to its express terms. In particular:
- (1) In no event shall any replacement credit be granted where there is evidence of improper handling, improper implantation or material alteration of the replaced Product.
  - (2) Medtronic is not responsible for any incidental or consequential damages based on any use, defect or failure of the Product, whether the claim is based on warranty, contract, tort or otherwise.
- D. The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this **LIMITED WARRANTY** is held by any court of competent jurisdiction to be illegal, unenforceable or in conflict with applicable law, the validity

of the remaining portion of the **LIMITED WARRANTY** shall not be affected, and all rights and obligations shall be construed and enforced as if this **LIMITED WARRANTY** did not contain the particular part or term held to be invalid.

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All Rights Reserved
© Medtronic, Inc. 2002, 2008