IMPLANTABLE PORT SYSTEMS

INSTRUCTIONS FOR USE

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ANGIODYNAMICS®

INCORPORATED

ANGIODYNAMICS®

INCORPORATED

IMPLANTABLE PORT SYSTEMS

en	INSTRUCTIONS FOR USE	1
du	INSTRUCTIES VOOR GEBRUIK	19
es	INSTRUCCIONES DE USO	39
it	ISTRUZIONI PER L'USO	59
fr	INSTRUCTIONS D'UTILISATION	79
gm	GEBRAUCHSANWEISUNG	99

INSTRUCTIONS FOR USE

Overview

Introduction	2
Product Description	2
Indications for Use	3
Contraindications	3
Potential Complications	3
Needles	4

Implantation Instructions

General Guidelines	4
Port and Catheter Preparation	5
Port Placement Considerations	5
Catheter Placement Considerations	6
Implantation of Preattached Catheter	6
Implantation of Attachable Catheter	9
Implantation of Arterial Port	11
Implantation of Peritoneal Port	12

Procedures for Use

General Guidelines	13
System Assembly	13
Bolus Injection/Continuous Infusion	14
Blood Sampling	16
System Maintenance	16
Patient Identification Card	17
Explantation of an AngioDynamics Port	17
Discontinuing System Use	17

Quick Reference

Guidelines	18
------------	----

OVERVIEW


Introduction

AngioDynamics, Inc. manufactures a complete line of implantable access ports. Ports are totally implantable devices designed to provide repeated access to the vascular system or a selected body site. These subcutaneous devices reduce the trauma associated with multiple punctures or the inconvenience of an externalized catheter.

The port is intended to facilitate frequent blood sampling or the delivery of medications, nutritions, blood products, and imaging solutions. Access is performed by percutaneous needle insertion, using an anti-coring needle (huber point).

Product Description

Ports manufactured by AngioDynamics are supplied as sterile devices, and are intended for single patient use only. They are available in a variety of port materials including polysulfone, acetal plastic, and titanium. Available catheter materials include polyurethane and silicone. Suture sites are incorporated into the port base to facilitate anchorage to the underlying fascia. Ports are available in single and dual, standard and low profile models. They are available in attachable or attached catheters. The single port has one high compression self sealing silicone septum, and the dual model has two.

Indications for Use

The AngioDynamics port line is indicated for any patient requiring repeated access of the vascular system for delivery of medications, nutritional supplementation, fluids, blood, blood products, and sampling of blood. Dual models are indicated for combination therapy, simultaneous infusions, withdrawal of body fluids, and bolus delivery during continuous infusion.

Contraindications

AngioDynamics port systems should not be implanted in the presence of known or suspected infections, bacteremia, septicemia, and peritonitis, in patients who have exhibited prior intolerance to the materials of construction, or patients whose body size or tissue is insufficient to accommodate the size of the port or catheter.

Potential Complications

Use of AngioDynamics port systems involve potential risks normally associated with the insertion or use of any implanted device or indwelling catheter including but not limited to:

Infection; occlusion; thrombophlebitis, pneumothorax; catheter malposition; migration and inadequate anchoring; hemorrhage; vessel trauma, including puncture, laceration, and erosion of vessel and the skin; catheter pinch-off (compression of the catheter between the clavicle and the first rib); hematoma; clot formation; catheter fragmentation; embolization; cardiac arrhythmia; cardiac puncture; cardiac tamponade; fibrin sheath, endocarditis; implant rejection; thoracic duct injury; thromboembolism; peritonitis; thrombosis; and drug extravasation (leakage).

Occlusion may result from clot formation inside the lumen of the catheter, precipitate formation inside the port from

incompatible drugs, or from catheter tip placement against a vein wall or valve.

Needles

Use of AngioDynamics anti-coring (19 to 22 gauge huber point) needles in all procedures is recommended. These needles have been designed and tested to ensure that septum life is preserved. Needles are for single use only.

Note: Septum Puncture Life — Under qualified testing procedures, the septum allows up to 2,000 punctures when tested at 10 PSI using a 22 gauge non-coring needle. This pressure exceeds typical levels experienced in clinical practice.

IMPLANTATION INSTRUCTIONS

General Guidelines

The following suggestions for surgical insertion are provided as an aid to facilitate safe and prolonged use of the AngioDynamics port systems. Since the reservoir body may be placed in a number of areas of the body and the outlet tube may be placed in a variety of vessels, or other selected sites, use the surgical procedure and the sterile technique which best suits your application and is appropriate for the patient. AngioDynamics recommends that the patient, as appropriate, be placed in the trendelenburg position.

PRECAUTIONS:

- Strict aseptic technique is of paramount importance when implanting any device.

- Before handling the port, ensure that fingers of surgical gloves are free of talc.
- When suturing around the silicone catheter, avoid excessive suture tightness to prevent occlusion of the catheter. Sutures should not be placed directly on polyurethane catheters.

Port and Catheter Preparation

Prime the port system prior to placement using 10 ml of normal saline or heparinized saline (100 units/ml). Attach the anti-coring needle to the syringe, penetrate the septum of the port, and flush the system.

With dual lumen systems, both lumens must be flushed with heparinized saline.

Caution: Do not exceed 40 psi pressure when administering fluid into system. The use of a 10cc or larger syringe will help prevent this from occurring.

Port Placement Considerations

- Placement needs to be supported by underlying bony structure.
- A minimum of three sutures should be used to secure port body.
- Port location should be convenient and comfortable to the patient.
- Avoid placing port system directly under port pocket incision.
- Avoid placing port too deep or too shallow (minimum 0.5 cm-maximum 2 cm under skin surface).
- Pre-operative mapping of location is recommended whenever possible.

Catheter Placement Considerations

Place catheter tip in area of high blood flow.

Warning: Avoid medial catheter placement into subclavian vein through percutaneous technique. This placement could lead to catheter occlusion, damage, rupture, shearing, or fragmentation due to compression of the catheter between the first rib and clavicle. Catheter shearing has been reported when the catheter is inserted via a more medial route in the subclavian vein.¹

1. Aiken DR, Minton JP. The “pinch-off” sign: a warning of impending problems with permanent subclavian catheters. *Am J Surgery* 1994; 148:633-636

The port catheter should be positioned at the selected site of therapy and secured by accepted surgical technique to prevent catheter dislodgement. Position should be confirmed by appropriate radiographic procedures.

Caution: Sufficient slack should be left between the catheter insertion point and the port body to preclude strain on the catheter.

Implantation of Preattached Catheter (Venous/Vascular)

Prime the port system prior to placement.

1. Select appropriate site for portal placement.
2. Measure appropriate catheter length. Provide slack from port site to allow for body movement, but not enough to allow kinking of catheter.
3. Trim excess catheter by cutting squarely across the distal end. Do not trim catheter at an angle since this could cause the catheter tip to seal off against the side of the vessel.

Percutaneous Procedure (preattached catheter)

- a. Select appropriate French-size sheath introducer.

- b. Puncture skin with introducer needle into the subclavian vein at selected venous site. Gently aspirate while inserting.
-

Warning: The use of alcohol, acetone, or solutions containing these agents may result in degradation of the plastic introducer needle hub.

- c. Remove syringe, leaving needle in place.
-

Caution: To prevent air embolism, place thumb over exposed orifice of needle.

- d. Slide the “J” guidewire straightener over the “J” tip of wire. Insert the straightened “J” tip through the percutaneous entry needle and advance the wire 5-10 cm into the vein. Verify guidewire position radiographically.
- e. Withdraw the needle and guidewire straightener, leaving the guidewire in place. Clamp guidewire with hemostat to prevent further advancement into the vascular system.
- f. Create a subcutaneous pocket for the port. An incision is made and pocket formed by either sharp or blunt dissection down to underlying fascia.
- g. Place port in pocket and pass catheter from port pocket to entry site via subcutaneous tunnel.
- h. Unclamp guidewire and advance dilator/sheath over the exposed “J” wire. Withdraw vessel dilator and “J” wire, leaving sheath in place.
-

Caution: To prevent air embolism, place thumb over exposed orifice of sheath.

- i. Insert catheter into sheath. Position the distal end of the catheter at the desired location. Peel away sheath while withdrawing it from vessel. Care should be taken not to withdraw catheter as sheath is removed. Catheter position should be confirmed radiographically. Secure catheter in place.
- j. Secure port body to underlying fascia using non-absorbable sutures and a minimum of three suture sites. Care should

be exercised so that incision does not cross septum of port after closure.

Caution: Avoid piercing catheter with suture needle.

- k.** Prior to wound closure, aspirate septum to confirm ability to withdraw blood. Flush port with 3-5 ml of 10-100 units/ml of heparinized saline. Maintain positive pressure on syringe plunger to avoid reflux of blood into catheter tip. Stabilize port while withdrawing needle.
- l.** Close incision after wound irrigation by appropriate surgical technique. Dress wound per hospital protocol.

Note: A 90 degree anti-coring needle with winged infusion set may be positioned in port septum intraoperatively for patient comfort during initial access.

Surgical Cutdown (preattached catheter)

Follow general port placement guidelines described under “Port Placement Considerations” and “Percutaneous Procedure”.

- a.** A small incision is made in deltopectoral groove to expose cephalic vein or a small transverse incision in neck to expose external jugular vein. Isolate vessel.
- b.** Introduce catheter through venotomy and advance to desired location. Confirm catheter placement by appropriate radiographic technique.
- c.** Anchor catheter at venotomy site. Avoid excessive suture tightness to prevent catheter occlusion.

Caution: Sufficient slack should be left between port and catheter insertion point to preclude strain on the catheter. When using external jugular vein, carefully position the catheter over clavicle to avoid kinking or occlusion.

Note: A 90 degree anti-coring needle with winged infusion set may be positioned in the port septum intraoperatively for patient comfort during initial access.

Implantation of Attachable Catheter (Venous/Vascular)

Percutaneous Procedure (attachable catheter)

Prime the port system prior to placement.

- a. Select appropriate French-size sheath introducer.
 - b. Puncture skin with introducer needle into the subclavian vein at selected venous site. Gently aspirate while inserting.
-

Warning: The use of alcohol, acetone, or solutions containing these agents may result in degradation of the plastic introducer needle hub.

- c. Remove syringe, leaving needle in place.
-

Caution: To prevent air embolism, place thumb over exposed orifice of needle.

- d. Slide the “J” guidewire straightener over the “J” tip of wire. Insert the straightened “J” tip through the percutaneous entry needle and advance the wire 5-10 cm into the vein. Verify guidewire position radiographically.
 - e. Withdraw the needle and guidewire straightener, leaving the guidewire in place. Clamp guidewire with hemostat to prevent further advancement into the vascular system.
 - f. Create a subcutaneous pocket for the port. An incision is made and pocket formed by either sharp or blunt dissection down to underlying fascia.
 - g. Unclamp guidewire and advance dilator/sheath over the exposed “J” wire. Withdraw vessel dilator and “J” wire, leaving sheath in place.
-

Caution: To prevent air embolism, place thumb over exposed orifice of sheath.

- h.** Insert catheter into sheath. Position the distal end of the catheter at the desired location. Peel away sheath while withdrawing it from vessel. Care should be taken not to withdraw catheter as sheath is removed. Catheter position should be confirmed radiographically. Secure catheter in place.
- i.** Trim proximal end of catheter and advance through subcutaneous tunnel to the port pocket. Attach catheter to the port body. (See System Assembly in Procedures For Use section).
- j.** Secure port body to underlying fascia using non-absorbable sutures and a minimum of three suture sites. Care should be exercised so that incision does not cross septum of port after closure.

Caution: Avoid piercing catheter with suture needle.

- k.** Prior to wound closure, aspirate septum to confirm ability to withdraw blood. Flush port with 3-5 ml of 10-100 units/ml of heparinized saline. Maintain positive pressure on syringe plunger to avoid reflux of blood into catheter tip. Stabilize port while withdrawing needle.
- l.** Close incision after wound irrigation by appropriate surgical technique. Dress wound per hospital protocol.

Note: A 90 degree anti-coring needle with winged infusion set may be positioned in the port septum intraoperatively for patient comfort during initial access.

Surgical Cutdown (attachable catheter)

Follow general port placement guidelines described under “Port Placement Considerations” and “Percutaneous Procedure”.

- a.** A small incision is made in deltopectoral groove to expose cephalic vein or a small transverse incision in neck to expose external jugular vein. Isolate vessel.

- b. Introduce catheter through venotomy and advance to desired location. Confirm catheter placement by appropriate radiographic technique. Catheter is passed to pocket site via subcutaneous tunnel.
- c. Anchor catheter at venotomy site. Avoid excessive suture tightness to prevent catheter occlusion.

Caution: Sufficient slack should be left between port and catheter insertion point to preclude strain on the catheter. When using external jugular vein, carefully position the catheter over clavicle to avoid kinking or occlusion.

Note: A 90 degree anti-coring needle with winged infusion set may be positioned in the port septum intraoperatively for patient comfort during initial access.

Implantation of Arterial Port

The common hepatic artery is the most common location for intra-arterial catheter placement, with the port body placed on either side of the lower rib cage. Other sites may be used depending on the organ targeted for drug delivery. It may be necessary to determine the catheter insertion site by direct visualization due to variations in patient anatomy.

1. Clamp gastroduodenal artery on either side of chosen entry site.
2. Trim catheter to necessary length ensuring that at least two suture beads are left on catheter.
3. Insert catheter into artery through small arteriotomy. Beaded catheter is advanced up to, but not into, lumen of common hepatic artery.
4. One fixation suture is placed around artery just behind bead at catheter tip; another suture is placed around second bead just outside arteriotomy.

5. Follow guidelines for port pocket formation as outlined previously.

Implantation of Peritoneal Port

1. Follow the guidelines for port body placement as outlined previously. The peritoneal catheter is tunneled to 3-5 cm below the umbilicus. This area is preferred due to its relative avascularity and reduced fascia resistance.
2. Make a 3-4 cm midline vertical incision. Incise fascia. Identify and separate subcutaneous tissue and pre-peritoneal fat. Grasp peritoneum with hemostat or forceps. Elevate and carefully incise under direct vision.
3. Insert index finger into peritoneal incision to check for adhesions.
4. Insert catheter into selected area, usually pelvic gutter region. Place purse-string suture around catheter to prevent leakage though not too tight as to occlude catheter. Once placement is accomplished, cuffed portion of catheter can be sutured into place above or at level of peritoneal incision.

Caution: To avoid injection into subcutaneous tissue, ensure that all catheter holes remain within the peritoneal cavity.

5. Attach an anti-coring needle to a 50 ml syringe filled with saline. Flush system with saline solution and immediately aspirate to confirm that flow is not obstructed and that no leaks exist.
6. Close fascia and skin incisions using a continuous monofilament fascia closure for watertight seal and standard skin closure. Dress wound.
7. Attach an anti-coring needle to a 10 ml syringe of normal saline. Penetrate septum and flush system.

PROCEDURES FOR USE

General Guidelines

- Each access of an AngioDynamics port should be performed using aseptic technique.
- The needle should be advanced through the septum until it contacts the base of the port body. Once positioned in the septum, the needle should not be rocked or tilted. Such movement may cause septum damage.
- At no time should the system be left open to air. Tubing clamps should be used to prevent inadvertent air embolism.
- When infusing into an AngioDynamics port system, do not exceed 40 psi.
- Do not use a syringe size smaller than 10ml. Smaller syringes may create an overpressurized system.

System Assembly

Note: Do not use surgical instruments for assembly.

1. **One Step Locking Mechanism as on Triumph-1®**–
Slide the locking mechanism over the end of the catheter, continuing to the proximal end, leaving 1 to 2 cm of catheter protruding. Advance the catheter just past the first ring on the port outlet tube. Slide the locking mechanism toward the port body, rotating clockwise, advancing the catheter into the port body.
2. **Bayonet Locking Mechanism (BLM) as on LifePort®**–
Slide the BLM over the end of the catheter. The suture tab should be in the upright position. Dry the outlet tube(s) and proximal end of the catheter. Slide the catheter approximately half way into the port outlet tube(s). Slide the BLM and catheter onto the outlet tube until the suture tab contacts

the port body. Turn Lock 90 (degrees) until the tab lays flush with the port base locking the catheter in place. With proper assembly, a short “doughnut” shaped section should be visible between the BLM and the port body.

3. **Snap-Lock™ (SL) Locking Mechanism as on Infuse-a-Port®**—Slide the SL onto the proximal end of the catheter leaving 1 to 2 cm of catheter protruding. Advance the catheter over the end of the barbed connector until the catheter begins to enter the port body. Slide the SL up to the port body and “snap” into receptacle of port.
4. **Blue Boot Locking Mechanism as on Vortex™, OmegaPort™, and TitanPort™**—Slide the boot over the end of the heparinized catheter. The tapered end of the boot should point away from the proximal end of the catheter. Dry the outlet tube and the proximal end of the catheter. Slide the trimmed end of the catheter tip onto the stem until the catheter is flush with the stem flanges. Slide the boot over the catheter and onto the stem until the boot contacts the port body.

Bolus Injection/Continuous Infusion

1. Identify the port septum by palpating outer perimeter of the port.
2. Observing aseptic technique, prepare injection site.
3. Attach syringe with normal saline to tubing and anti-coring needle.
4. Insert the anti-coring needle through the skin perpendicular to the port and advance slowly until contact with the base is made.
5. Unclamp tubing and inject 3-5 ml of normal saline to flush port catheter. Clamp tubing.

Note: If vascular access, needle placement should be confirmed by aspiration.

6. Remove syringe from tubing and attach drug syringe. Unclamp tubing and inject drug slowly.
7. For continuous infusion, connect infusion pump to extension tubing. Tighten all connections. Position and secure height adjustable wings of infusion set. Start infusion pump. Open tubing clamp.

Caution: Examine injection site closely. If patient feels an abnormal sensation or pain at injection site, it may indicate the drug has extravasated. Discontinue infusion immediately and proceed with accepted extravasation protocol. Notify physician immediately.

8. Clamp tubing and carefully disconnect syringe.
9. Reattach syringe filled with normal saline. Unclamp tubing and flush catheter.
10. If additional drug infusions are required, flush port with an adequate volume of saline between infusions and repeat steps 6 through 10.
11. Heparin Lock Procedure
 - a. Attach syringe containing 3-5 ml of heparinized saline (100 units/ml) to tubing.
 - b. Flush catheter.

Caution: Maximum flow rate of 5 ml/min is recommended for heparin lock procedure. This flow will minimize blood reflux into catheter.

- c. Maintenance of positive pressure on syringe plunger will prevent blood reflux.
12. Gently withdraw needle from port septum and apply adhesive bandage.

Caution: It is extremely important to adequately flush the port after blood withdrawal. Occlusion of the catheter can occur if blood is left in the catheter for an extended period of time.

Blood Sampling

Blood sampling may be performed as an isolated procedure, at the time of bolus injection, or during the continuous infusion process.

1. Insert the anti-coring needle into the prepared site and flush with 5-10 ml of normal saline solution.
2. Withdraw “discard sample” consisting of 5 ml of blood. Discard this sample and syringe. Perform required blood sampling.
3. Immediately flush the catheter with a minimum of 10 ml of saline followed by 3-5 ml of 10-100 units/ml heparinized saline solution to establish the heparin lock.

System Maintenance

Venous Systems

After use, the port should be flushed with at least 20 ml of normal saline. This should be followed by a “heparin lock”. If the port is not being used, a “heparin lock” should be administered every 4 weeks.

Arterial Systems

After use, the port should be flushed with at least 20 ml of normal saline. This should be followed by a “heparin lock”. Positive pressure should be maintained on the syringe plunger at all times to minimize reflux of blood into the catheter. If port is not being used, a “heparin lock” should be administered once a week.

Peritoneal Systems

The intraperitoneal port should be flushed after use, or once a month if not in use, with 10 ml of heparinized saline at a concentration of 10 units/ml.

Patient Identification Card

A Patient Identification Card is included in the Patient's Manual with accompanies each port. The completed card should be given to the patient who should be instructed to carry it at all times.

Explantation of an AngioDynamics Port

Ports and catheters which are explanted due to suspected malfunction should be returned to AngioDynamics for analysis. The AngioDynamics customer service department must be contacted to obtain a return authorization number and instructions.

AngioDynamics will provide an explant kit for use in storage and shipment of the explanted device. Hospitals must advise AngioDynamics of any infectious disease that the patient is known to have.

No returned product will be accepted without an RGA number and properly packaged in an AngioDynamics explant kit or equivalent.

Discontinuing System Use

AngioDynamics recommends that the clinician consider explantation of the system once it is determined that it is no longer required for therapy. If the clinician decides to leave the system in place, AngioDynamics recommends that periodic X-rays be taken with the patient in an upright, arms at side position. This procedure will verify catheter placement and detect problems with the system such as pinching of the catheter between the clavicle and first rib.

ANGIODYNAMICS SYSTEM CARE GUIDELINES

Site Preparation: Always access the system using aseptic technique.

Syringes: 10 ml syringes or larger are recommended for all flushing or injection procedures. Use of smaller syringes may result in system damage.

Needles: Use AngioDynamics anti-coring (huber point) needles only.

Saline Flushes: Prior to drug administration, flush the system with saline solution to remove heparin. If more than one drug is administered, flush the system with saline solution between drugs. After patient treatment is completed, always flush the system to cleanse the catheter and port chamber. With dual lumen systems, flushing must be repeated for each side.

Heparin Flush Schedule: To keep the AngioDynamics system patent, the system must be flushed with heparinized saline at regular intervals.

Heparin Concentration: (10-100 units/ml) of heparinized saline. Typical volume of 3-5 ml.

Venous Systems: “Heparin lock” once every 4 weeks.

Peritoneal Systems: “Heparin lock” once every 4 weeks.

Arterial Systems: “Heparin lock” once a week.

Note: Follow institutional guidelines for infusion set use. Center for Disease Control (CDC) recommends that I.V. tubing be changed every 48 hours.

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CAUTION: Federal (U.S.A.) Law restricts this device
AngioDynamics UK Ltd.**