Aerosol - Nebulizer **MEDIneb® F700**

FLAEM®

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Congratulations on your purchase and thank you for choosing our unit.

Visit our Internet site www.flaemnuova.it to view the whole range of Flaem products

**Carefully read these instructions and keep them in a safe place for future reference.
Use the unit only as described in this instruction manual.**

IMPORTANT SAFEGUARDS

This a medical device for home use. It is to be used only upon medical prescription.

Before using the unit for the first time, and periodically throughout its life, check the power supply cord for any damage. If there is any damage, unplug immediately and return the product to an authorized dealer for service.

- Children and disabled persons should always use the unit under strict supervision of an adult who has read this manual.
- The unit should not be used in the presence of anaesthetic substances inflammable with air, oxygen or nitrous oxide.
- Never use adapters for power supplies different from the one indicated on the bottom of the unit. Keep the cord away from heated surfaces.

Do not handle the plug with wet hands. Never use the unit while having a bath or shower. Never submerge the unit in water; should this accidentally occur, unplug immediately - do not reach for the submerged the unit, always unplug first. After removing a unit from water, do not use again, ensure it is returned to an authorized service dealer.

- Do not wash the unit under running water or do not submerge it in water.
- Unit with case not protected against ingress of liquids.
- Always unplug the unit immediately after use.
- Repairs should be performed by authorized personnel only. Unauthorized repairs will invalidate the warranty and may represent a hazard to the user.

INSTRUCTIONS FOR USE

Before each use, the nebulizer and the accessories should be cleaned according to the instructions in the section "CLEANING AND DISINFECTION".

1. Plug the power cord into a power socket corresponding to the voltage of the unit.
2. Carefully wash your hands before preparing the inhalation therapy.
3. Open the medication nebulizer by turning the upper part (9c) counterclockwise.
4. Pour the drug prescribed by the physician into the lower part (9a). Close the medication nebulizer by turning the upper part (9c) clockwise.
5. Connect the accessories as shown in the drawing on page 1. In order to avoid waste of nebulized medication, it is possible to use the drug saving device (8). Connect one end of the device to the medication nebulizer (9) and the other end to the tube (7). Close the opening with the finger to perform nebulization; reopen it to stop nebulization.
6. Start the unit by turning the switch on (3).
7. When the treatment is finished, turn off and unplug the unit.
8. Put caps A and B back in their seat in order to protect the inner part of the medication nebulizer.
9. It may happen that during the use of the nebulizer visible humidity deposits form in the connection tube (7), disconnect the tube from the nebulizer and let it dry using the compressor air flow.

CLEANING AND DISINFECTION

Accessories

Carefully wash all parts used with the medicament under running water with mild (non abrasive) detergents or in the dishwasher.

All accessories (except for the connecting tube) can be disinfected or sterilised choosing one of these methods:

- a) by boiling in a bain-marie (for max 5 min.); this method cannot be used for the mask
- b) by soaking in a solution of water and vinegar (3/4 water, 1/4 vinegar). Rinse thoroughly with drinking water.

It is recommended to replace the medication nebulizer once a year.

Unit

Wipe the case only with a damp cloth using non abrasive detergents.

The unit is equipped with an air filter which (5) should be replaced when dirty. Use only Flaem's original accessories.

IMPORTANT!

Microbial contamination

A personal use of the nebulizer and of the accessories is recommended in order to avoid any risks of infection.

Interactions

The materials used that come into contact with drugs have been tested with a wide range of drugs. However, due to the variety and continuous evolution of the drugs, the possibility of a chemical interaction cannot be excluded. Once opened, the drug should be used as soon as possible, and prolonged contact between drug and nebulizer should be avoided.

The use of ESSENTIAL OILS can cause damage to the nebulizer.

Disposal: In conformity with Directive 2002/96/EC on waste electrical and electronic equipment (WEEE), at the end of its life the unit should be taken to an electrical appliance collection point for its environmentally friendly disposal.

UNIT TECHNICAL FEATURES

Mod. MEDINEB F700

Voltage:	230V ~ 50Hz 140VA	115V ~ 60Hz 180VA	220V ~ 60Hz 160VA	100V ~ 50/60Hz
Fuse:	T 2A - 250V	T 3A - 250V	T 2A - 250V	T 5A - 250V
In compliance with:	Dir. 93/42/EEC (EN 60 601-1) 			

Operating pressure (with nebulizer): 0.8 bar approx.
Max pressure: 2.5 bar approx.
Compressor air output: 10 LPM approx.
Dimensions: 20x27x10.5 cm
Weight: 2.3 Kg
Sound level (at 1 m): 58 dB (A) approx.
Continuous use

SYMBOLS

Class II unit

Type B unit

Warning: follow the "Instructions for use"

Alternating current

Switch OFF

Switch ON

Never use the unit while having a bath or shower

Latex Free

RAPIDFLAEM 2 NEBULIZER- TECHNICAL FEATURES

Capacity max: 7 ml
***Aerosol output:** *with cap B* 0.25 ml/min approx.
without cap B 0.51 ml/min approx.

***Particle:** *with cap B* MMAD 2.8 µm
without cap B MMAD 3.1 µm

*Average values measured with 0.9% saline solution. Certified with TSI Aerosizer 3220 computerized laser system in compliance with Flaem internal procedure.

Operating requirements: temperature min. 10°C; max 40°C
air humidity min. 10%; max 95% **RH**

Storage requirements: temperature min. -25°C; max 70°C **RH**
air humidity min. 10%; max 95% **RH**

Operating/storage atmospheric pressure: min. 690 hPa; max. 1060 hPa **P**

Electromagnetic compatibility

This aerosol therapy unit is designed to meet the current requirements for electromagnetic compatibility. However if you suspect that the unit is interfering with the operation of your TV, radio or other electrical appliance, try relocating the unit until the interference ceases. Alternatively try plugging the unit into a different power socket.

The Manufacturer, the Retailer and the Importer are responsible for safety, reliability and performance of the unit only if:

- the unit is used in accordance with the instructions for use;
- the electrical system of the facility where the unit is used adheres to IEC standards;
- assembly, additions, recalibrations, modifications or repairs are carried out by authorized personnel.

The manufacturer reserves the right to introduce technical changes.

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