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General Instructions for pH Sensors of the SE 552-Memosens® Series

Caution!

Please read these instructions carefully before putting the sensor into operation, in order to ensure faultless function. Operation and servicing of the sensor must only be performed by specially trained personnel authorized by the operating company.

Warning!

After an Ex-certified sensor has been used with non-Ex components, this sensor must no longer be used in a hazardous location.

1. Introduction

The Knick Series SE 552 pH sensors with Memosens® connector system are low-maintenance sensors with pressurized electrolyte and integrated temperature detector for automatic temperature compensation of the pH signal during calibration and operation of the sensors. The SE 552 sensors are autoclavable and in-situ sterilizable.

2. Safety Information

The sensors are intended only for the use described in section 1 above. The sensors consist partly of glass. Be sure to wear safety goggles and gloves when cleaning or calibrating a sensor using acid or alkaline solutions.

The sensors may be operated and serviced only by personnel familiar with the equipment in question and who have read and understood these instructions. The sensors are designed to withstand a maximum operating pressure of 4 bars and fulfill all current PED requirements (Pressure Equipment Directive 97/23/EC, Article 3, Sec. 3). Please take extra care during installation/fitting of the sensor to ensure that neither the polymer body nor the thread (Pg 13.5) suffer any damage.

Please observe the hazardous-area concept for Memosens® sensors described in section 4 when installing the sensor in a hazardous area.

3. Product Description

The markings on each sensor or on the packaging label include the following information:

Knick	Manufacturer of sensor
SE 552X/*-MSN	Model designation
	Memosens® connector (MS) Solution Ground (A: with, N: without) Length (1: 120 mm, 2: 225 mm)

pH 0...14	pH measuring range
0...100 (135) °C	Temperature range (sterilization)

The year of production of the sensor is included in the serial number, which is printed on and also electronically stored in the sensor head. See also enclosed quality certificate.

4. Hazardous-Area Concept

The sensors are intrinsically safe equipment of Category 1. They may be installed within explosive gas atmospheres up to Zone 0.

- The sensors must only be connected to an intrinsically safe, approved cable of Type CA/MS-***X** (BVS 09 ATEX E 083 X) from Knick Elektronische Messgeräte GmbH & Co. KG.

- SE552X/*-MS*, category II 1 G, temperature codes:
T6: $-20\text{ °C} \leq T_a \leq +70\text{ °C}$
T4: $-20\text{ °C} \leq T_a \leq +120\text{ °C}$
T3: $-20\text{ °C} \leq T_a \leq +135\text{ °C}$

- The sensors must not be operated under electrostatically unfavorable process conditions. Strong vapor and dust streams acting directly on the sensors and measuring cable must be avoided. Shaft length and solution ground are not Ex-relevant.

- When using the sensors and cables, the stipulations for electrical installations in explosive gas atmospheres (EN 60079-14) must be observed.

- The stipulations/conditions for safe use of the connected intrinsically safe Memosens® sensor must be observed.

Note!

If the ambient temperatures indicated are observed, no impermissible temperatures for the temperature class in question occur at the equipment.

5. Installation and Commissioning

- On unpacking, check the sensor for mechanical damage of the pH sensitive glass and shaft. Report any damage to your Knick service team.
- Remove the watering cap and briefly rinse the sensor with de-ionized water. Carefully remove the protection droplet of silicone on the junction with the enclosed knife, (as shown on the green pictures on the packing).
- Again rinse the sensor well with de-ionized water. After rinsing, the sensor should only be dabbed dry with a tissue. Do not rub the pH-sensitive glass, since this can lead to electrostatic charging and sluggish response times.
- Check the space behind the pH-sensitive glass for the presence of any air bubbles and remove them by gently swinging the electrode in a vertical plane.
- Please refer to the instruction manual of the relative fitting in order to install the sensor correctly.
- Sensor and cable are plugged together and locked by a quarter turn of the coupling ring. See function description of Memosens® cable.
- Note:** Ex sensors are marked by an orange-red ring.

6. Operation

6.1 Calibration of Sensor and Transmitter

2-point pH-calibration is recommended for the SE 552 sensor. First remove the watering cap. Dip the sensor successively into two different buffer solutions with given pH values (e.g. pH 4.01 and pH 7.00) and calibrate it to these buffer values. Please refer to the instruction manual of the pH transmitter for further details.

Memosens® sensors should preferably be precalibrated in the lab.

6.2 Sterilization of Sensor

For application in sterile processes such as fermentation the sensors are sterilized prior to the production cycle. This sensor type allows sterilizing of the whole fermenter/reactor unit including pH electrode by means of autoclaving. Sterilization can also be effected in-situ by means of steam or superheated process medium in contact with the lower sensor parts.

6.3 Temperature Detector

As mentioned in section 1, the temperature detector integrated in the SE 552 is primarily intended for automatic compensation of the pH signal and not for any accurate and safe temperature indication or control of the process temperature.

7. Maintenance

- Carefully rinse the sensor tip and junction with de-ionized water after each production cycle. Under no circumstances must measuring solution be allowed to dry on these parts!
- When the sensor is not in operation, store it with sensor tip and junction well submerged in electrolyte (3 mol/l KCl).
- If a sensor is stored dry for a few days by mistake, let it soak in electrolyte for several hours before use.
- The Memosens® connector requires no maintenance. It is resistant against humidity and contaminations. Nevertheless observe cleanliness for hygienic reasons.

8. Troubleshooting

(sluggish response time, insufficient sensitivity etc.)
The glass membrane can be reactivated using a HF solution. Immerse the sensor tip (pH-sensitive glass part only) in the HF solution for about 5 minutes (make sure that the junction is not immersed!). Immediately afterwards, rinse the sensor thoroughly with de-ionized water and store it in electrolyte (3 mol/l KCl) for 12 hours. Possible protein contamination of the junction can be treated with pepsin/HCl.

9. Warranty

12-month warranty after date of shipment in the event of manufacturing faults.

10. Disposal


Please observe the applicable local or national regulations concerning the disposal of "waste electrical and electronic equipment".

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