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1. INTRODUCTION

This instruction manual is intended to assist the user for safe and efficient operation of the semi-automatic digital blood pressure monitor Realcheck™ BASIC 2 (hereinafter: device).

The device must be used in accordance with the procedures described in this manual. It is important to read and understand the entire manual, especially the section “Information about blood pressure measurement”. This device is intended for the non-invasive measurement of systolic and diastolic arterial blood pressure and pulse frequency. Always contact the physician if measurement is taken on children or persons with arrhythmia, as errors may occur.

**READ INSTRUCTION MANUAL CAREFULLY BEFORE USE
THE DEVICE MUST NOT BE DISASSEMBLED
FOR ANY PROBLEM PLEASE CONTACT THE TECHNICAL SERVICE OF DISTRIBUTOR
FOR SPECIFIC INFORMATION ON YOUR OWN BLOOD PRESSURE PLEASE CONTACT YOUR PHYSICIAN.**

SAFETY RULES

1. Before using the device, always check the electrical requirements on manufacturer label and verify the type of battery to use itself.
2. Particular precautions must be made concerning electromagnetic compatibility. The medical device must be set-up and used according to information supplied with the present manual;
3. Respect the safety rules and particularly:
 - Children and/or disable persons may only use the device under the supervision of an adult in possession of his/her full faculties.
 - If the device is stored at low temperature, do not use it immediately. Leave the device at room temperature for at least one hour.
 - Avoid extreme environmental conditions and the exposure to direct sunlight.
 - **NEVER DIP ANY PART OF THE DEVICE AND ACCESSORIES IN WATER, IT ISN'T WATER PROOF.**
 - The device must be used only by adult people.
 - Cleaning the device with soft and dry cloth. Do not use petrol, thinners or similar solvents. Spot on the cuff can be removed carefully with a damp cloth. Do not wash any part of the device.
 - If you do not use the device for a long time, please take the batteries off.

KEEP THE BATTERIES OUT OF THE REACH OF CHILDREN.

- It is recommended the use of alkaline battery.
 - Do not wash or iron the cuff
 - Do not squeeze the air hoses nor the cuff
 - Keep the cuff away from cutting edges and do not twist it
 - Avoid to drop the device or threat it roughly in any way.
4. **The device must not be disassembled. For any problem please contact the Technical Service of the Distributor. Failure to comply with the above invalidates the safety of the device and the guarantee claims.**
 5. **This medical device must be used exclusively for the use for which it has been designed and as described in this manual.** Any different use must be considered incorrect and therefore dangerous; the manufacturer and the distributor cannot be considered liable for damage caused to people, things or animals by improper, incorrect and/or unreasonable use.
 6. Measurement on people with serious peripheral circulation or renal problems, people who have had an ictus or if they are unconscious and the case of serious arrhythmia may be difficult and bring extra error; it is recommended to consult the physician.
 7. In case of malfunctions contact the technical service of the distributor. If by technical check the device will result opened by user, the guarantee will be expired.

PRINCIPLE OF OPERATION

This device adopts the oscillometric technology with Fuzzy Algorithm to measure the arterial blood pressure and pulse rate. The cuff is wrapped around the arm and inflated by pumping the inflation bulb. The sensor of the device relieves weak fluctuation of the pressure in the cuff produced by extension and contraction of the arm artery in correspondance to each heartbeat. The amplitude of the pressure waves is calculated and converted in millimeters of the mercury column, and it is displayed, as final result of measurement, by digital value.


Warning: This device can not provide reasonable accuracy if used or stored in the temperature or humidity beyond the range stated in the section “Technical Characteristics” of this manual.

TECHNOLOGY USED


Fuzzy Algorithm is the processing algorithm which considers the specificity of individual heartbeat, providing higher accuracy of measurement.

2. INFORMATION ABOUT BLOOD PRESSURE MEASUREMENT

1. It is necessary to know that arterial blood pressure is subjected to sharp fluctuations. The level of the arterial blood pressure depends on many factors. Generally arterial blood pressure is lower in summer and higher in winter. Arterial blood pressure changes with atmosphere pressure and it is affected considerably by many factors, e.g. physical loads, emotional excitability, stress, diet, etc. Medicines, alcohol, smoking, affect greatly the level of blood pressure. When blood pressure is measured at the hospital, the value is always higher than that at home. The reason is the stress caused by the ambient, which is known as 'White coat effect'. Blood pressure will raise in low temperature, so it is recommended to take the blood pressure at room temperature. If the device was stored at temperatures out of range recommended, it is necessary to leave it at room temperature for at least 1 hour, otherwise malfunction may occur. Blood pressure varies with age and individually, and it is recommended to write down daily the measures, to evaluate with physician the one's "normal blood pressure".

The illustration is from British Hypertension Society


2. In case of patients with cardio-vascular diseases, the measurement should be taken under physician's instruction. Not in any case it is allowed to change the dosage of the medications suggested by the physician without his consultation.
3. An accurate measurement of blood pressure may be difficult in serious arteriosclerosis, weak pulses, or in patients with manifest fluctuation of heart contraction rhythm. It is recommended to consult a physician for a correct reading of the blood pressure measurements.
4. It is necessary to keep quiet during measurement to get accurate readings. Measurement should be conducted in quiet environment at room temperature. Don't eat or smoke before a measurement. This device is supplied with the standard cuff which is fit for the arm size 22-32 cm. Care should be taken to ensure that the cuff size is appropriate for the person whose blood pressure is being taken. Children and adults with cuff size fall outside the range 22-32 cm should select special size cuffs. Please contact the dealer to get these special size cuffs.

ATTENTION: Do not use cuffs other than the original cuff contained in this kit!

3. INSTRUCTIONS

3.1 BATTERY INSERTION


1. Open the battery cover and install one 'AAA' type battery into the battery compartment as indicated. Make sure that the polarity is correct;
2. Close the battery compartment cover.
 - Replace the battery when the replacement indication  appears on the display or nothing appears after O/I button is pressed;
 - When replace battery, please use always new battery. The use of not new or bad quality battery can reduce the work capacity of the device
 - Don't use rechargeable battery;
 - If the device will be unused for long time, please take off the battery;
 - The life span of the battery included can be shorter than expected in consequence of storage period.


3.2 FITTING THE CUFF

1. Insert the edge of the cuff approximately 5 centimeters into the D-ring as shown.

2. Put the cuff on the left upper arm with the tubes toward the direction of the palm. If measurement on your left arm is difficult, you can use right arm. In this case, it is better to know that the readings may differ from the left one .


3. Wrap cuff around your upper arm with the lower edge of the cuff approximately 2-3 centimeters above the elbow. The mark **ARTERY** must be over the artery of the arm.

4. Press the cuff to make sure that it is safety closed. The cuff should not be too tight or too loose. Two fingers should be easily put in between cuff and upper arm.

5. When cuff is wrapped around your arm, the edge of the cuff must be inside the <NORMAL> area (22-32 cm). This means the cuff size is correct. If it's not in the <NORMAL> area it could be possible the needs of a special cuff.

6. Sometimes it is difficult to make a regular fit due to the shape of the user's upper arm, the cone-shape assembly of cuff is also acceptable.

7. If the clothes fits closely to the upper arm, or the sleeve has been rolled up take them off to get accurate measurement.

Once the cuff is wrapped, make sure that its position is approx at the same level of the heart, and that the forearm is extended naturally on the table.

3.3 CARRY OUT A MEASUREMENT

It is necessary to keep quiet during measurement to get accurate readings. Measurement should be conducted in quiet environment at room temperature. Don't eat or smoke before a measurement.

ATTENTION: Do not use cuffs other than the original cuff contained in this kit


Fig.1



Fig.2

1. Hold the measuring unit with the hand whose arm is assembled with cuff, and hold the inflation bulb with the other hand. Before the measurement, take 3~5 minutes of deep breathing and relax. Don't talk or move the arm.
2. Press button 'O/I', and all symbols will appear on display in 2 seconds as Fig.1. Then two short beep will sound, and the number '0' (zero) will appear on the display with upward arrow flashing to remind the user to start inflating, as Fig.2. If there is too much air in the cuff, the display will show downward arrow to remind the user to deflate the air in the cuff as Fig.3. Please push the air release valve to deflate the air until the '0' (zero) and upward arrow appear.
3. Pump the inflation bulb to increase pressure in the arm cuff. The display will show the pressure value within the cuff. Continue to pump until the pressure reaches approx 30~40 mmHg above the regular systolic pressure. If the regular systolic pressure is unknown, pump until 180 mmHg (If the pressure is not enough to get accurate measurement, the display will show upward arrow following ERR indication, to remind the user to pump to higher pressure). Then, stop inflation and the pressure begins to decrease gradually, during this time the device calculates the blood pressure and pulse rates.


4. There will be a long beep following the end of measurement. The blood pressure and pulse reading will show in the display as Fig.4.
The reading will automatically stored in the memory set.
5. Press rapid air release valve to release the remaining air in the cuff and press the button 'O/I' to turn off the device. If you'd like to carry out another measurement, please rest for at least 15 minutes. If the power supply is not switched off and the device is left unused for 3 minutes, it will automatically switch off.

3.4 FAST INTERRUPTION DURING MEASUREMENT

If you do not feel well during measurement or want to stop the measurement for any reason, press the quick air release valve to deflate the air in the cuff.

4. MEMORY FUNCTION

4.1 MEMORY RECALL

1. REALCHECK™ BASIC 2 can store 90 sets of readings and will automatically calculate the average value of the latest 3 readings. When the memory is full (90 sets of readings are stored), the oldest reading will be replaced by new one. The set of Memory keep the data until the power supply is removed;
2. After a measurement is finished or when the device is in standby, the user can press the Memory button to recall data.
 - 2a: Press once the Memory button and the display will show the average value of the latest 3 readings as Fig. 5;
 - 2b: Press again, and the display will show '01' and then turns to another screen to show the last reading as Fig. 6;
 - 2c: Press again, and the display will show '02', which means the second last reading and backwards to the first reading.

Fig.7

4.2 MEMORY CLEARANCE

After a measurement is finished or when the device is in standby, press Memory button for at least 5 seconds, the display will show 'Clr' which means that all the stored reading is removed as Fig.7.

5. ERROR MESSAGES AND LOW BATTERY INFORMATION

INDICATION	POSSIBLE REASON	CORRECTION METHODS
	The cuff is put on wrongly or the tube plug is inserted too loosely. Movement of arm/hand or talking during measurement. The cuff is not inflated to enough pressure. Arrhythmia.	Make sure that cuff is put on correctly and the tube plug is inserted tightly and repeat the measurement Repeat the measurement with following completely recommendations of manual. Repeat the measurement with pumping cuff to higher pressure. Consult the physician
	The battery is weak	Replace the battery with new one

6. MAINTENANCE AND REPAIR

1. It's necessary to protect this device against high moisture and heat, direct sunlight, shock, solvents, alcohol.
 2. Remove the battery if the device will be unused for a long time, and keep the battery out of the reach of children.
 3. Keep the cuff and inflation bulb away from sharp subjects and don't extend or twist the cuff.
 4. Use only soft and dry cloth to clean the device.
 5. The cuffs are delicate and must be handled with care. You can clean the cuff cover with damp cloth.
- WARNING: Under no circumstances any part of device and its accessories could be washed Keep the device and all of its accessories out of the reach of children**

7. TROUBLESHOOTING

SYMPTOM	CHECK POINT	REMEDY
No display when the O/I button is pressed	The battery is exhausted. The polarity of battery is wrong The contacts of battery compartment is polluted	Replace the battery with new one Install the batteries correctly Clean the battery terminals with dry cloth
The reading is extremely low or high	Is the cuff at the same level as the heart? Is the cuff wrapped right? Did you move your arm during measurement? Did you talk or move your arm (or hand) during measurement?	Make sure that your posture is right Wrap the cuff correctly Relax during measurement Keep quiet and silent during the measurement
Pulse rate is too low or too high	Did you talk or move your arm (or hand) during measurement? Did you make measurement right after exercise?	Keep quiet and silent during the measurement Take measurement again after resting for almost 15 minutes
The batteries are run down soon	Faulty batteries are used	Use only alkaline batteries of good quality
The device is automatically turned off	It is the result of automatic turn off system	This is to save the power consumption of the device, and it is not a fault

8. TECHNICAL CHARACTERISTICS

Model	REALCHECK BASIC 2
Size	Ø 65(D) x 26(H) mm
Weight	Approximately 50 g not including battery
Measuring method	Oscillometry
Measuring range	40 to 260 mmHg (pressure) 40 to 160 beats/minute (pulse rate)
Measuring accuracy	± 3 mmHg for systolic and diastolic pressure ± 5% of the reading for the pulse rate
Inflation	Manual
Rapid deflation	Manual Valve
Batteries	1“AAA” x 1.5V
Memory	90 sets of memories
Operation temperature and humidity	+10° C to + 40° C , 85% and below
Storage temperature and humidity	-20° C to + 50° C, 85% and below
Cuff size	Applicable for arm size 22-32 cm

Guidance and manufacture's declaration – Electromagnetic immunity
(as request by regulation EN 60601-1-2:2001)

The model REALCHECK BASIC 2 Blood Pressure Monitor is tended for use in the electromagnetic environment specified below. The customer or the user of the model Automatic Blood Pressure Monitor should assure that is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environmet-guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floor should be wood, concrete or ceramic tile. If floors are covered with synthetic material the relative humidity should be at least 30%.
Electrical fast transient / burst IEC 61000-4-4	Not applicable	-	-
Surge IEC 61000-4-5	Not applicable	-	-
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	Not applicable	-	-
Power frequency magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment

Guidance and manufacture's declaration – Electromagnetic immunity
(as request by regulation EN 60601-1-2:2001)

The model REALCHECK BASIC 2 Blood Pressure Monitor is intended for use in the electromagnetic environment specified below. The customer or the user of this device should assure that is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environmet-guidance
Conducted RF IEC 61000-4-6	Not applicable	-	Portable and mobile RF communications equipment should be used no closer to any part of the REALCHECK BASIC 2 Blood Pressure Monitor, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. $d = \left[\frac{3.5}{E_1} \right] \sqrt{P} \quad 80\text{MHz to } 800\text{MHz}$ $d = \left[\frac{7}{E_1} \right] \sqrt{P} \quad 800\text{MHz to } 2.5\text{GHz}$ Where P is the maximum out put power rating of the transmitter in watts (W) according to the transmitter manufacturer and is recommended separation distance in meters (m). Field strength from fixed RF transmitters, as determined by an electromagnetic site survey ((•••)) (a), should be less than the compliance level in each frequency range (b).  Interference may occur in the vicinity of equipment marked with the following symbol:
Ratiated RF IEC 61000-4-3	3V/m 80MHz-2.5GHz	3V/m	

Note 1 At 80MHz and 800MHz the higher frequency range applies.

Note 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- (a) Field strength from fixed transmitters, such as base stations for radio telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the REALCHECK BASIC 2 Blood Pressure Monitor should be observed to verify normal operation. If abnormal performance is observed additional measures may be necessary, such as re-orienting or relocating the device.
- (b) Over the frequency range 150KHz to 80MHz, field strength should be less than 3V/m

9. SYMBOLS

	Type B
	CE marking in conformity with EC directive 93/42/EEC
	Warning, consult the instruction manual
	To keep in cool and dry place
	Storage temperature
	On
	Off
	Direct current

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“Actuation of European directives 2002/95/EC, 2002/96/EC and 2003/108/EC, for reduction in use of dangerous substances in the electric and electronic device and for garbage disposal”

The symbol as over applied on the device or its packaging means that at the end of its useful life the product must not be disposed of with domestic waste. At the end of device useful, the user will must deliver it to the able collecting centres for electric and electronic garbage, or give back to the retailer in the moment of equivalent new device purchasing, one against one.

Disposing of the product separately prevents possible negative consequences for the environment and for health, deriving from inadequate disposal. It also allows the recovery of materials of witch it's made up in order to obtain an important saving of energy and resources and to avoid negative effects to the ambient and health.

In case of abusive disposal of device by user, will be applied administrative endorsements in compliance with current standard.

10. WARRANTY

- The unit is covered by a 2-year warranty against defects in materials and workmanship. The warranty does not extend to units that have undergone unauthorised repairs or been tampered with.
- The warranty period commences on the purchase date. It is mandatory to attach the sales receipt or invoice serving as proof of the purchase date to the warranty certificate below, which must be filled out in every part.
- Under the warranty, any components with manufacturing defects will be replaced and/or repaired free of charge. The Purchaser will bear all expenses for transporting, delivering and collecting the unit.

The warranty does not cover:

- The accessories supplied with the unit or parts subject to normal wear.
- Repair of unsubstantiated defects.
- Repairs of units that have failed or been rendered defective as a result of improper use, carelessness or negligence or if the damage is not attributable to the manufacturer (accidental falls, careless transport, etc.).

The warranty does not cover any direct or indirect damage or injury caused to persons, animals or property as a result of improper use of the product or during its period of inefficiency.

WHEN SUBMITTING A WARRANTY CLAIM, IN OBSERVANCE OF THE ABOVE CONDITIONS, SEND THE UNIT TOGETHER WITH THE CERTIFICATE BELOW AND PROOF OF PURCHASE TO:

SSL HEALTHCARE ITALIA S.P.A. VIA MARCO EMILIO LEPIDO, 178/5 - BOLOGNA (ITALY).

CERTIFICATE TO BE RETURNED IN CASE OF REPLACEMENT/REPAIRS

THE WARRANTY IS VALID ONLY IF THE SALES RECEIPT OR INVOICE IS ATTACHED

Unit Type/Model: _____

Lot: _____ Serial No: (SN) _____

Date of purchase: _____

PURCHASER'S DETAILS

Last Name and Name: _____

Address: _____

Telephone: _____

Description of the defect: _____

☐ I hereby authorise the use of the above information in accordance with the Law on Privacy 196/03

Signature confirming acceptance of warranty conditions _____

CE 0123

Medical Device
Classe IIa
(DIR 93/42/CEE)

Fabbricante/Manufacturer:

CA.MI snc

Via Ugo La Malfa, 31
43010 Pilastro (PR) - Italy

Distributore/Distributor

SSL Healthcare Italia SPA

Via M.E.Lepido, 178/5
40132 Bologna - Italy