PIASTRA FREDDA PF100

Codice: 23-PF100

C€ IVD

PRECAUZIONI D'USO

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
C€ IVD

USE PRECAUTIONS

Before using the instrument, read carefully the instructions and warnings contained in this manual and keep it for further reference. They supply important indications regarding the functions and safety for installing, using and maintaining the instrument.

Bio-Optica Milano S.p.A. cannot be held responsible for any damage caused by improper or incorrect use and by the non-observance of any of the prescription provided in this manual and by the safety regulations in force.

- 1. After unpacking, make sure that the instrument is complete and not damaged by transport.
- 2. Before connecting the instrument to the power supply make sure that its rating corresponds to that of the power supply.
- 3. This instrument must only be used for the purpose for which it was designed, that is, as <u>cooling plate</u> for laboratory use. Any other use is to be considered improper and therefore hazardous.
- 4. The instrument must only be used by authorized and professionally qualified technician.
- 5. The electrical safety of this instrument can be guaranteed only if it is correctly connected to an efficient earth circuit as indicated by current electrical safety regulations. It is necessary to check this fundamental safety prerequisite, and if in doubt, ask to check the circuit. The instrument is provided with a power supply cable having 2 wires + ground tap that have to be connected to the power supply socket.
- 6. Do not remove the chassis or parts of it during operation. Switch off the instrument and disconnect the power supply cable before opening it. This operation must to be effected only by authorized and professionally qualified technician.
- 7. To eliminate instrument malfunctioning risks, do not work near strong magnetic fields and do not use transmitters such as cellular phones near the instrument. In case of serious malfunctioning switch off the instrument and contact the Technical Assistance Service.
- 8. All waste material, both infectious and radioactive, deriving from the appliance working cycle must be disposed in compliance with the regulation in force.

This appliance is marked from this symbol, in compliance with EU directive 2002/96/CE regarding electric and electronic appliances waste. This mean that the instrument, at the end of its useful life, must be collected separately from other refuse. The user must deliver it to the special differentiated refuse collection centres, that are predisposed by the public authority.

- 9. The contents of this manual is subject to change without further notice.
- 10. Please find enclosed the declaration of conformity.

11. Graphic symbols indicated on the label (positioned near the instrument's power supply socket):
Symbol for CATALOGUE NUMBER:
REF
Symbol for SERIAL NUMBER:
SN
Symbol for ALTERNATING CURRENT:
~
Symbol for FUSE:
Symbol for CONSULT THE INSTRUCTIONS:
Ţ <u>i</u>
Symbol for EC MARK:
(€
Symbol for IN VITRO DIAGNOSTIC-MEDICAL DEVICE:
IVD
Symbol for DISPOSAL OF ELECTRIC AND ELECTRONIC EQUIPMENT:
Symbol for DATE OF MANUFACTURE:
Symbol for MANUFACTURER:

DESCRIPTION

Instrument used to obtain the rapid cooling of histological samples included in paraffin.

Dimensional features

Dimensions (WxDxH): 500x550x290 mm.

Cooling surface dimensions (WxD): 420x245 mm.

Weight: 30 Kg.

Electrical connections

Power supply: 230V~ 50/60Hz.

Power: 600 Watt.

Fuses: N. 2 fuses of 4 Ampere - 5x20 mm - T4AH250V.

Other connections

Water connections: Not necessary. Fumes aspiration/filtration: Not necessary.

Structural features

- Painted sheet steel chassis.
- Stainless steel cooling surface which can house up to 70 standard embedding cassettes. Provided with plastic anti-dripping edge.
- Refrigerant without CFC: R134 A, 100 grams or R413, 130 grams (starting from serial number 435NC).

Technical features

Working temperature adjustment: From +20°C to -20°C through electronic thermostat with microprocessor.

Reading precision ±2℃.

Working parameters visualization

and modification:

visualization Through control panel with digital display, bilingual (Italian and English).

Programming possibility: Instrument starting and switching off, weekly.

Installation

Position the instrument on a level and stable working bench and connect it to the power supply socket (230V~50/60Hz) using the provided cable.

Important: Do not use any extension or adapter and do not modify the provided cable. Make sure that the back grid has at least a 10 cm free space in order to allow the aeration of the cooling system.

USE

- 1) Press the green general switch to start the cooling system.
- 2) Regulate the temperature following the instructions shown on the following pages. T1 corresponds to the stainless steel plate temperature.

Backlit control panel

TIMER DELAY: Delay the automatic instrument's turn off.

TEMPERATURE SETTING: Set the desired temperature value.

MENU: Allows to access the functions.


ENTER: Accept/store values.

UP/DOWN SELECTIONS BUTTONS: Modify values or select the functions.

Display

Setting the temperature

Press to display the temperature and to change it. The value is saved automatically (after 5 seconds the main screen appears on the display) or by pressing for 4 seconds.

Setting time and date

Before setting the timer, control time and date by pressing simultaneously and . The day will flash; if necessary, modify it by pressing to view the time and, if necessary, modify it by pressing to exit.

Setting the timer

The timer allows to plan the starting and switching off of the instrument for every day of the week.

Example: From Monday to Friday: starting (ON) at 08.00 a.m. - switching off (OFF) at 18.30 p.m. On Saturday and on Sunday: switching off.

Press to view this screen:

Progr Day Work Monday TIMER

Press to set the starting or switching off for Monday. To confirm press. Repeat for every day of the week, up to Sunday.

Afterwards will be displayed this screen:

Progr Timer Start T1. 08.00

To modify the starting time, press and, to confirm,

After pressing , this screen will appear:

Progr Timer STOP 18.30

Press to modify the switching off of the resistance.

At last press for 3 seconds to exit. The main screen will appear on the display.


Timer delay


To delay the automatic switching off of the instrument, press for 4 seconds. When the switching off time will flash, it will be possible to modify it by pressing.

This operation does not modify the preset switching off time.

Timer esclusion (manual operating)

To operate the instrument manually, it is necessary to set all days to ON. (see "Setting the timer"). On the display will appear ON. To start or switching off the instrument it will be necessary to use the general switch every time.

Cleaning

Before effecting the cleaning, switch off the instrument and unplug the cable from the socket. Use only alcohol or non-aggressive detergents. Do not use abrasive products or acids in order to avoid ruining the varnished or plastic parts.

The instrument does not need to be sterilized because is not expected the treatment of fresh samples but only of histological samples fixed and included in paraffin.

Maintenance

Effect the periodic cleaning of the back grids.


Replacing the fuses

<u>If the instrument doesn't start</u>, check that the electric cable has been connected properly, there is current and the two fuses under the instrument's socket aren't burnt. If necessary, change them with fuses of identic value.

Fuses holder socket (N. 2 delayed fuses of 4 Ampere - T4AH250V).

<u>Instructions for the replacement</u>: Switch off the appliance and unplug the cable from the socket. Make a light pressure on the little carter covering the fuses (if necessary using a little screwdriver), change them, close the carter and check the instrument's ignition.

MAIN COMPONENTS


CODE	DESCRIPTION	QUANTITY
37-AR09E2F6A	6 AMPERE WIRE FILTER	1
37-B311/1N	1 TEMPERATURE THERMAL TIMER	1
37-FP606	NOISE FILTER	1
37-GL60TB	COMPRESSOR	1
37-RS208	STATIC RELAY	1
37-VI100	SCREWS FOR PLATE	6
37-129-432	GREEN WATERPROOF SWITCH	1
37-0105	FRAME	1
37-1023	RUBBER FEETS	4
37-1045	POWER SUPPLY CABLE 1	
37-4042N	POLYCARBONATE PANEL	
37-1099/1	UPPER CHASSIS 1	
37-1099/2	LOWER BASE 1	
37-1099/3	BASIN 1	
37-1099/4	COPPER PLATE 1	
MU23-PF100	USER MANUAL	1

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