Fisher
Bioblock Scientific

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The labels applied to the unit warn the user on the dangers to which he is exposed during the use or the maintenance. The labels must be left on the unit and substituted if they are no more readable.

Warn of danger

To prevent accidents which could be caused by the unit

Start the operation

Do not start to use the unit before reading directions for use and maintenance

SAFETY RULES

- 1)** *Do not turn on the instrument before placing into position the containers protecting the blades.*
- 2)** *If a fuse or a lamp is to be replaced, first of all unplug the unit.*

CLEANING

Always unplug the unit before cleaning. Use a damp cloth and not flammable, not aggressive detergents.

PERSONAL PROTECTION EQUIPMENT

The equipment used for personal protection must be compatible with the dangers due to the working materials and from the glass vessel.

MAINTENANCE

According to the law on products guarantee our instrument must be repaired at our factory, save different agreements with our agents.

GUARANTEE

Starts from the date of delivery note and is referred to the number of register of the single unit.

The manufacturer, engaged in a continuous improvement of his products reserves the right to modify any of their characteristics without notice.

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1. Forewords

The equipment is aimed to satisfy two different analytical operations:

- a) Jar tests for the optimization of dosing coagulants and polyelectrolytes in wastewaters treatment plants.
- b) Leaching tests using diluted acetic acid solutions or water saturated by carbon dioxide for the evaluation of leacheability of heavy metals by solid wastes to be sent to different kinds of dumps.

2. Instrument description

The JLT6 equipment is manufactured in a strong metallic frame with a purposely studied epoxy painting giving a high degree of protection against chemical and mechanical corrosion.

The examined samples can be lighted by a back fluorescent lamp controlled by a switch located on front panel. The stainless steel stirring shafts are settable in height by a self blocking clutch.

Rotation speed is controlled by a gear motor powered by direct current controlled electronically by a microprocessor. The programmable speed of rotation from 10 to 300 r.p.m. (within 1 r.p.m.) is continuously compared to the set value and shown by a display. The JLT6 equipment is equipped by an electronic timer allowing to program the stirring time on two different scales (minutes or hours) or to choose a continuous operation.

3. Checking the instrument as received

When received, after packing removal, check the integrity of unit.

The delivery is composed by:

- 1 Jar Test model JLT6
- 1 Connecting electric cable

4. Assembling and installing

Verify that electric power voltage corresponds to the voltage shown on the label.

The equipment is to be placed on a bench in a stable position.

5. Operation controls

Place in each stirring position an empty beaker centred in correspondance of each stirring blade.

The beakers can be placed raising the shafts and lowering them again to the bottom.

Connect the equipment to electric power by the supplied cable. Turn on the equipment by the main switch that lights showing that power is on.

The back light is operated turning on the switch located on front panel. The chosen speed is setted by the corresponding selector for each stirring place.

Programming the electronic timer

Three different modes of operation can be chosen by the key **Set** of the "Time" window:

- Continuous
- Timer controlled on a minute basis. Maximum programmable time: 999 minutes.
- Timer controlled on an hour basis. Maximum programmable time: 99 hours.

Simply pressing the key **Set** the operation mode changes from continuous to timer controlled while the display shows:

- continuous operation "..."
- timer controlled on a minute basis "000"
- timer controlled on an hour basis "h00"

The time basis is changed keeping the **Set** key for some seconds. The chosen time is set by the keys $\uparrow$ e $\downarrow$. The display shows continuously the residual time and at the same time the third decimal dot twickens showing that the counting is on. At the end of the programmed time the stirring is automatically stopped and an acoustic signal warns about the end of test, while the display shows "end".

Programming the speed of stirring

The stirring speed from 10 to 300 r.p.m., in steps of 1 r.p.m., is selected by pressing the keys $\uparrow$ $\downarrow$ of "Speed" window. Keeping one of the keys pressed the speed of selection increases. The stirring is not operative if a timer controlled operation was chosen and the residual time is 0 (zero).

Black Out

The JLT6 unit is equipped by an internal memory that retains set stirring speed and residual times allowing the return of operation at the same conditions when power is again available.

6. Start

After placing the sample containing beakers in each stirring position and turning on the main switch, set the mode of operation by pressing the key **Set** of "Time" window.

- Continuous: the display shows "..."
- Timer controlled minutes: the display shows "000"
- Timer controlled hours: the display shows "h00"

If a timer controlled operation is chosen the required time is set pressing the keys $\uparrow$ and $\downarrow$ of the "Speed" window. The lighted of sample beakers is possible at any moment by operating the light switch.

7. Operations at the end of work

Remove the sample beakers and wipe off or dry by filter paper the shafts and blades from sample residues in order to avoid drippings on the equipment

8. Accessories

1000 ml glass beaker	11.1001
Imhoff cone, plastics	11.1002
Graduated Imhoff cone, glass	11.1003
2 places Imhoff cone stand	11.1004

9. Replacement parts

	230V/60Hz F10500109	115V/60Hz F10510109
Gear motor	MOTO002	MOTO002
Fuse 0.5 A 5x20 mm	ELFU003	ELFU008
Fluorescent lamp 30 W	ELLA002	ELLA002
Electronic card	EESH080	EESH080
Driving belt 165mm	ACCT112	ACCT112
Driving belt 250mm	ACCT145	ACCT145
Driving belt 340mm	ACCT167	ACCT167
O ring 2015	ACOR001	ACOR001
Shaft gasket	ACGU001	ACGU001
Switch	ELIB003	ELIB003
Transformer	"	ELTR005

10. Technical features

Power supply	230V/50Hz	F10500109
	115V/60Hz	F10510109
Total power	50 W	
Dimensions (W x H x D)	940 x 350 x 300 mm	
Weight	17 kg	
Stirring places	N° 6	
Stirring volume (for each place)	1 l	
Single or multiple control	Single control (same speed for each place)	
Speed interval	From 10 to 300 r.p.m.	
Speed selection	1 r.p.m.	
Speed precision	± 1 r.p.m.	
Timer	From 0 to 999 minutes or from 0 to 99 hours	
Fluorescent lamp	Yes	
Safety class	Class 2 IEC 1010	
Protection degree	IP30 IEC 529	

11. Analytical methods

Coagulation and flocculation of wastewaters

Chemical coagulants are added to wastewaters for the primary purpose of removing suspended solids and also phosphorus. The addition of chemicals reduces also heavy metal concentrations and improves disinfection efficiency. Chemical coagulation can be obtained by adding to wastewaters defined amounts of lime (calcium hydroxide), alum (aluminium sulfate) or iron salts (ferric or ferrous).

The coagulation-sedimentation process typically involves:

a) injection and mixing of the coagulant that neutralize the electric charge, prevalently negative, on suspended particles. Good and rapid mixing of the coagulant and wastewater is important to ensure efficient use of the chemical. Typical detention times in plant basins for turbulent mixing are 15 to 120 seconds.

b) agglomeration of the coagulated particles into large settleable flocs. This is accomplished by stirring slowly the water to allow the solids formed by the addition of chemicals to grow in size so that they can be removed by gravity settling.

The stirring is usually performed by slowly rotating paddles during 10 to 30 minutes.

c) sedimentation of the flocculated material by gravity separation. The solids are allowed to settle by gravity from the chemically treated standing wastewater. The effluent from the chemical settling tank is very clear and does not usually require further filtration.

The quantity of chemical coagulant required to achieve good coagulation varies with time and from wastewater to wastewater.

Typical coagulant doses are:

75 to 250 g/m³ for alum,
45 to 90 g/m³ for ferric chloride,
200 to 400 g/m³ for lime.

Laboratory evaluation by jar test

The choice and dosing of chemical coagulant to be adopted for the removal of suspended solids from wastewaters are derived from the results of laboratory evaluations by the so called jar test. Multiple stirrers with reproducible stirring speeds allow to adopt standard conditions for the test, that are the basic requirement for reproducible results.

The standard conditions relate to:

- volume of wastewater in each stirring position
- dimensions and shape of containing vessel and stirring paddle
- time and speed of stirring during turbulent mixing phase
- time and speed of stirring during flocculation phase
- evaluation criteria for the obtained results.

The most diffused adopted conditions are:

- 1000 ml glass beakers (jars), tall form, 105 mm of diameter
- 600 ml wastewater samples
- stirring paddles 25 mm high, 75 mm wide, with a thickness of 1 mm
- stirring height of paddles roughly at middle height of wastewater sample
- quick turbulent stirring after chemical addition: 120 rpm for 120 seconds
- slow speed stirring during flocculation: 30 rpm for 25 minutes
- first evaluation of results after 5 minutes of sedimentation after the end of stirring.

The results can be evaluated on the basis of different criteria:

a) evaluation of floc dimensions with a numerical degree of merit (0 = no flocs, 2 = very small hardly visible, 4 = small flocs, 6 = medium size flocs, 8 = good size flocs, 10 = very large flocs).

b) time from the addition of chemicals to first appearing of flocs

c) evaluation of residual turbidity of supernatant, after a determined sedimentation time, by turbidimeter

d) measurement of electrokinetic potential of suspended particles on a sample taken immediately after the addition and mixing of chemicals. Coagulations should be easier when zeta potential of particles is very low.

e) evaluation of filterability of clarified water by standardized membrane filters under pressure. The reduction of water flow with time is related to the degree of clogging of filters due to residual unsettled suspended matter.

Other measurements relate to the pH of water after chemical treatment (lime raises the value while iron and aluminium salts lower it) and the temperature of water at which the experiments are performed.

References

- American Society for Testing Materials Norm ASTM D-2035-G4 T

- Passino R. and M. Beccari (1970).

Standardization of jar test in coagulation-flocculation processes for the removal of inorganic turbidities (in Italian). *Acqua e Aria*, Milano, 1-10, October 1970.

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- U.S. Association of Environmental Engineering Professors (1972). *Environmental Engineering Unit Operations and Unit Processes Laboratory Manual*. J.T. O'Connor Ed. 350 ppg.

Leaching tests

Are used to simulate the behaviour of a waste in a dump submitted to the leaching activity of rain water. Two different methods use leaching solutions of different type:

a) acetic acid, proposed for dumps receiving both organic and inorganic wastes;

b) carbon dioxide saturated water, for dumps receiving only inorganic wastes.

Waste sample to be examined

If a liquid phase is present it must be removed by filtration or centrifugation and preserved in a refrigerator at 2-4°C. The granulometry of the solid phase must be lower than 9.5 mm (standard sieve).

Waste extraction

20-25 g of the solid phase are weighed with a precision of ± 0.1 g.

Method a. A volume of distilled water corresponding to 16 times the weight of the sample is added.

The pH of the mixture is kept to 5 ± 0.2 using 0.5 N acetic acid. The stirring is performed for 24 hours, controlling the pH every 15 minutes during the first hour and then every hour. The pH is corrected to 5 ± 2 using again 0.5 N acetic acid. If after 24 hours of stirring the pH is higher than 5.2, a new correction is made and the stirring prolonged for other 4 hours.

Method b. The leaching solution is given by a volume corresponding to 20 times the sample weight of distilled water saturated by carbon dioxide for not less than 15 minutes. The stirring is performed during 6 hours without any pH correction.

Stirring speed

Must be able to maintain continuously mixed the solid phase with the leaching solution.

Analysis of extract



The leaching solution and the liquid phase if present in the original sample are filtered by a membrane filter of 0.45 μm porosity and submitted to analysis for the single metals by the methods commonly used for waters.

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
Electric scheme 230V/50Hz

- 1 Retarded fuse 5x20 mm 0.5A (3,15 A)
- 2 Power inlet with fuses
- 3 Lighted main switch
- 4 Lamp holder, left
- 5 Starter
- 6 Lamp holder, right
- 7 Reactor
- 8 Transformer
- 9 Electronic card
- 10 Gear motor
- 11 Speed transducer
- 12 Lamp switch

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
Electric scheme 115V/60Hz

- 1** Retarded fuse 5x20 mm 0.5A (3,15 A)
- 2** Power inlet with fuses
- 3** Lighted main switch
- 4** Lamp holder, left
- 5** Starter
- 6** Lamp holder, right
- 7** Reactor
- 8** Transformer
- 9** Electronic card
- 10** Gear motor
- 11** Speed transducer
- 12** Lamp switch
- 13** Transformer (115 V)

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Pompe

POMPE PERISTALTICHE

** distribuiti solo in Italia - distributed only in Italy*

Thanks for choosing a VELP's product !

Since 1983 Velp offers to the professional in the field a range of sophisticated and reliable equipment which make available high levels of know-how and operative capacity at competitive price.

Velp works according to **ISO 9001, ISO14001 and OHSAS 18001** Quality System Certification. Instruments are built according to the International norms **IEC 1010-1** and to the rules of **CE** mark.

We present you our product Lines:

Environment

C.O.D.
B.O.D. – B.O.D. Sensor
REFRIGERATED THERMOSTATS
INCUBATORS
JAR TEST
METALS IN SLUDGES AND SEDIMENTS
NEPHELOMETRIC TURBIDITY
NITROGEN IN DIFFERENT FORMS
PERFORMANCE OF ANAEROBIC DIGESTERS
PHENOLS
LEACHABILITY TESTS

Stirring

MAGNETIC STIRRERS
HEATING MAGNETIC STIRRERS
MULTIPLACE MAGNETIC STIRRERS
LIGHTED MAGNETIC STIRRERS
HEATING PLATES
VORTEX
STIRRERS
STIRRERS WITHOUT MOTOR

Food and feed

DIGESTERS
STEAM DISTILLING UNITS
SCRUBBER
CATALYST TABLETS
NITROGEN-FREE WEIGHING BOATS
EXTRACTORS FOR RAW FIBER
EXTRACTORS FOR DIETARY FIBER
EXTRACTORS FOR SELECTIVE
EXTRACTION OF SOLUBLE PRODUCTS BY
SOLVENTS (FOR FAT)

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Declaration of compliance

We

manufacturer VELP SCIENTIFICA s.r.l.
address Via Stazione, 16
USMATE (MI)
Italy

Under our responsibility declare that the product is in according with the following standards:
EN 61010-1/2001
EN 61326-1/2000

And satisfies the essential requirements of directives:

Machines directive 98/37/CEE

Low voltage directive 73/23/CEE

Electromagnetic compatibility directive 89/336/CEE

Plus modifications and that the documents listed in the annex V are available at Velp's offices as per the machine directive.

Velp Scientifica srl

Via Stazione, 16
20040 Usmate (Milano) Italy
Tel. +39 039 62 88 11
Fax +39 039 62 88 120

e-mail: inse@velp.it

internet: <http://www.velp.com>