# INSTRUCTIONS FOR USE

## Distal Radius Fracture Repair System

**R<sub>x</sub> ONLY** - Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

### DESCRIPTION:

The Distal Radius Fracture Repair System contains bone plates for the repair of either dorsal or volar radial fractures. Included in the instrumentation sets are K-Wires, bone screws, fixation pegs, fragment plates and specialized instruments.

The Distal Volar Radius Anatomic (DVR<sup>®</sup> Anatomic) plate utilizes 3.5mm cortical screws to affix the plate to the diaphysis and fixed angle locking pegs for the distal bone fragment(s). Alternatively, multidirectional threaded pegs (MDTP) are available to achieve trajectories other than the predefined angles.

The Dorsal Nail Plate Anatomic (DNP<sup>®</sup> Anatomic) is a bone stabilization device consisting of an intrafocal nail plate that incorporates four fixed angle locking bone pegs for the distal bone fragment(s) and up to three proximal cortical screws. The DNP Anatomic has a narrow distal plate-like section that lies on the surface of the distal fragment and a proximal nail-like section that is introduced into the diaphysis of the radius through the fracture site. Fixed angle pegs are used to stabilize the distal fragment(s) to the plate section and cortical screws are used to fix the proximal fragment to the proximal intramedullary section.

Contents of this surgical tray are covered by one or more of U.S. Pat. Nos. 6,358,250, 6,364,882, 6,730,090, 6,712,820, 6,706,046, 6,866,665, and 6,926,720, and U.S. and international patents pending.

### INDICATIONS (DVR Anatomic and DNP Anatomic Systems):

The Distal Radius Fracture Repair System is intended for the fixation of fractures and osteotomies involving the distal radius.

### INDICATIONS (Fragment Plate System):

The Fragment Plate System is intended for essentially non-load bearing stabilization and fixation of small bone fragments in fresh fractures, revision procedures, joint fusion and reconstruction of small bones of the hand, foot, wrist, ankle, humerus, scapula, finger, toe, pelvis and craniomaxillofacial skeleton.

### CONTRAINDICATIONS:

If any of the following are suspected, tests are to be performed prior to implantation. Active or latent infection. Sepsis. Insufficient quantity or quality of bone and/or soft tissue. Material sensitivity. Patients who are unwilling or incapable of following post operative care instructions.

### PRECAUTIONS:

- An implant must never be reused. Previous stresses may have created imperfections that can potentially lead to device failure. Protect implant appliances against scratching or nicking. Such stress concentration can lead to failure.
- Orthopaedic instrumentation do not have an indefinite functional life. All re-usable instruments are subjected to repeated stresses related to bone contact, impaction, routine cleaning and sterilization processes. Instruments should be carefully inspected before each use to ensure that they are fully functional. Scratches or dents can result in breakage. Dullness of cutting edges can result in poor functionality. Damaged instruments should be replaced to prevent potential patient injury such as metal fragments into the surgical site. Care should be taken to remove any debris, tissue or bone fragments that may collect on the instrument. Most instrument systems include inserts/trays and a container(s). Many instruments are intended for use with a specific implant system. It is essential that the surgeon and operating theatre staff are fully conversant with the appropriate surgical technique for the instruments and associated implant, if any.
- Do NOT open the volar wrist capsule. Doing so may cause devascularization of the fracture fragments and destabilization of the volar wrist ligaments.
- If necessary, contour the DVR Anatomic plate in small increments. Excessive contouring may weaken or fracture the plate.
- Exercise care when bending the fragment plates to avoid weakening or fracture of the plates.
- Ensure removal of all F.A.S.T. Guide™ inserts after use.
- Do NOT use fully threaded pegs (FP) with the DVR Anatomic and DNP Anatomic plates. The fully threaded pegs (FP) are designed for use with the fragment plates.
- Do NOT use peg/screw lengths that will excessively protrude through the far cortex. Protrusion through the far cortex may result in soft tissue irritation.
- SP-series screws are NOT intended to provide sub-chondral support and use should be limited to capture of remote bone fragments where partially or fully threaded pegs cannot be used.
- Do NOT permanently implant K-wires through the holes of the plate as they may back out and cause tissue damage. Use of the K-Wires allows you to provisionally secure the plates to the anatomy.
- Do NOT use the MDTPs in the distal row of the DVR Anatomic. The MDTPs are intended to be used only with the DVR Anatomic plates. Ensure the MDTPs are installed after insertion of the fixed angle pegs.

### WARNINGS:

Although the surgeon is the learned intermediary between the company and the patient, the important information conveyed in this document should be conveyed to the patient. The patient must be cautioned about the use, limitations and possible adverse effects of these implants including the potential for these devices failing as a result of loose fixation and/or loosening, stress, excessive activity, load bearing particularly when the implants experience increased loads due to a delayed union, non-union or incomplete healing. The patient must be warned that failure to follow postoperative care instructions may cause the implant or treatment to fail.

**For implantation instructions, please refer to the following Surgical Technique Guides:**

- DVR ANATOMIC
- DNP ANATOMIC
- F<sup>3</sup> FRAGMENT PLATING SYSTEM

### HOW SUPPLIED:

The required components and specialized instruments are supplied non-sterile in a container suitable for steam sterilization. Note that select items that are part of the DNP Anatomic and the Fragment Plate Systems are optional. All components and instruments may be purchased separately.

### Sterilization:

The system is provided non-sterile and is intended for steam sterilization at the healthcare institution.

The system must be cleaned prior to sterilization. Microsurgical instruments should be cleaned separately from other instruments. Disassemble instruments, as appropriate and inspect for damage. Cracked or bent instruments should be replaced. Inspect all instruments prior to sterilization or storage to ensure instruments are suitable for use.

Steam sterilization may be accomplished using the parameters shown in the table. The recommendations in ANSI/AAMI ST79, "Comprehensive Guide to Steam Sterilization and Sterility Assurance in Healthcare Facilities" or an equivalent standard should be followed with the exceptions noted in the table. The cycles are suggested based upon validation of a single, wrapped, multiple insert case, within a properly maintained autoclave. It is critical that process parameters be validated for each facility's individual type of sterilization equipment and product load configuration.

The use of flash sterilization is not recommended.

**Note:** It is NOT required to remove the F.A.S.T. Guide inserts to sterilize the plate. Validation testing has been performed to demonstrate a sterility assurance level of 10<sup>-6</sup> with the parameters in the table.

Cycle Type	Minimum Temperature	Min Exposure Time / Dry Time [Wrapped]
Prevacuum (with 4 vacuum pulses)	132 – 134°C	8 min / 20 minute dry time

### STORAGE INSTRUCTIONS:

Store in a cool dry place and keep away from direct sunlight. Prior to use, inspect product packaging for signs of tampering, damage, or water contamination.

### MRI INFORMATION

The Distal Radius Fracture Repair System has not been evaluated for safety and compatibility in the MR environment. The devices have not been tested for heating or migration in the MR environment. The risks associated with a passive implant in an MR environment have been evaluated and are known to include heating, migration, and image artifacts at or near the implant site.

### SYMBOLS:

**MATL:** Material  
**Ti:** Titanium Alloy  
**SS, SST:** Stainless Steel  
**MADE IN:** Made in <<Country>>  
**NON-STERILE:** Supplied non-sterile

Manufacturer

Authorized Representative in the European Community

**Biomet Trauma**  
 56 East Bell Drive  
 P.O. Box 587  
 Warsaw, Indiana, 46581 USA  
 Tel: 1-800-348-9500  
 Fax: 574-372-3968

**Biomet UK LTD.**  
 Waterton Industrial Estate  
 Bridgend CF31 3XA  
 UK  
 Tel: +44 (0) 1656 655221  
 Fax: +44 (0) 1656 645454

**For recognized manufacturer and model designation, please refer to the product label.**

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