INSTRUCTIONS FOR USE

Diaphyseal Plates

R_x ONLY - Caution: Federal (USA) law restricts this device to sale by or on the order of a physician

Indications

The Diaphyseal Plate is intended for the fixation of fractures, osteotomies and nonunions of the clavicle, scapula, olecranon, humerus, radius, pelvis, distal tibia, and fibula, particularly in osteopenic bone.

Description

The Diaphyseal Plate utilizes 3.8 mm multidirectional screws and 90° Locking screws with set screws are used to affix the plate to the bone.

Material Specification

The Diaphyseal Plates and implantable components are made of a Stainless Steel alloy (316L). This material meets the requirements of ASTM Designation F 138, *Standard Specification for Wrought 18 Chromium-14 Nickel-2.5 Molybdenum Stainless Steel Bar and Wire for Surgical Implants (UNS S31673)*.

Sterilization

The Diaphyseal Plate System is provided non-sterile. This system is intended for steam sterilization at the health-care institution.

Prior to sterilization the shoulder bone depth gauge (SBDG) must be disassembled and cleaned. Microsurgical instruments should be cleaned separately from other instruments. Disassemble instruments, as appropriate and inspect for damage. Cracked or bent instruments should be replaced. Inspect all instruments prior to sterilization or storage to ensure instruments are suitable for use.

Steam sterilization may be accomplished using the parameters shown in the table. The recommendations in ANSI/AAMI ST79, "Comprehensive Guide to Steam Sterilization and Sterility Assurance in Healthcare Facilities" or an equivalent standard should be followed with the exceptions noted in the table. The cycles are suggested based upon validation of a single, wrapped, multiple insert case, within a properly maintained autoclave. It is critical that process parameters be validated for each facility's individual type of sterilization equipment and product load configuration.

The use of flash sterilization is not recommended.

Cycle Type	Minimum Temperature	Min Exposure Time / Dry Time [Wrapped]
Prevacuum (with 4 vacuum pulses)	132 – 134°C	10 min / 20 minute dry time

Storage Instructions

Store in a cool dry place and keep away from direct sunlight. Prior to use, inspect product packaging for signs of tampering, damage, or water contamination.

How Supplied

The required components and specialized instruments are supplied non-sterile in a container suitable for steam sterilization. All components and instruments may be purchased separately. The system includes the items shown below.

Nomenclature

DIAPHYSEAL PLATE

MULTIDIRECTIONAL SCREW (MD)

90° LOCKING SCREW (NL)

90° SET SCREW (NL-SS)

Order No.	Description
131213000	Diaphyseal Plate, 7 Holes
131213001	Diaphyseal Plate, 9 Holes
131213002	Diaphyseal Plate, 11 Holes
131213003	Diaphyseal Plate, 13 Holes
131213004	Diaphyseal Plate, 15 Holes
MD-Series*	3.8mm diameter Shoulder Multidirectional Screw, 20mm to 38mm in 2mm increments
NL-Series*	3.8mm diameter Shoulder 90° Lock Screw, 20mm to 38mm in 2 mm increments
NLSS	90° Set Screw
KW20	K-Wire, 2.0mm
DB28	Drill Bit, 2.8mm
FHDS	Hex Driver
SDI	Square Driver Insert
MQC	Mini Quick-Connect Handle
SBDG	Bone Depth Gauge
SSTG	Shoulder Soft Tissue Guide
KPE	Key Periosteal Elevator

* Order number ends with component nominal length in mm. Individual Components may be supplied separately.

Contraindications

If any of the following are suspected, tests are to be performed prior to implantation. Active or latent infection. Sepsis. Insufficient quantity or quality of bone and/or soft tissue. Material sensitivity. Patients who are unwilling or incapable of following post operative care instructions.

Precautions

An implant shall never be reused. Previous stresses may have created imperfections that can potentially lead to device failure. Instrumentation must be inspected for wear or damage prior to usage. Protect implant appliances against scratching or nicking. Such stress concentration can lead to failure.

Warnings

Although the surgeon is the learned intermediary between the company and the patient, the important information conveyed in this document should be conveyed to the patient. The patient must be cautioned about the use, limitations and possible adverse effects of these implants including the potential for these devices failing as a result of loose fixation and/or loosening, stress, excessive activity, load bearing particularly when the implants experience increased loads due to a delayed union, nonunion, or incomplete healing. The patient must be warned that failure to follow postoperative care instructions may cause the implant or treatment to fail.

Clinical History

These instructions are based upon current experience. The physician may wish to vary the procedure in accordance with clinical judgment.

Technique for the Diaphyseal Plate (PLT)

- Expose and debride fracture using the appropriate surgical technique.
- Use traction and direct manipulation to reduce the fracture.
- If appropriate, secure the reduction by using 2.0 mm k-wires inserted through the dual k-wire/suture holes on the plate.
- Position the spherical end of the tissue guide in an oblong hole of the plate. Using the 2.8mm drill bit, drill hole through both cortexes.
- Remove the soft tissue guide from oblong hole and insert the depth gauge and determine the cortical screw length required.
- Select the appropriate length Multi-Directional cortical screw and secure the plate to the bone using this screw.
- If utilized, remove the K-wire/s from plate and adjust plate position, if necessary, by loosening the cortical screw first then securing the screw when the optimal position of the plate is achieved.
- Evaluate reduction and placement using fluoroscopy and readjust as necessary. Use both AP and lateral views for this purpose.
- Using the 2.8 mm drill bit, drill holes for the remaining cortical screws that will be used to fix the plate to the diaphysis.
- Multidirectional screws or 90° Lock Screws can be used in the standard holes.
- The cylindrical end of the Soft Tissue Guide must be used to drill holes when using the 90° Lock Screws. 90° Lock Screws provide very strong fixed angle characteristics. Multidirectional screws can be placed in any direction along a wide cone in space.
- The 90° Lock Screws can be locked into fixed angle mode by the use of set screws. Apply the set screws and fully tighten them using the square driver.
- If necessary, secure the surrounding soft tissue to the dual k-wire/suture holes of the plate using suture.
- Obtain final radiographic views.
- Close the wound using appropriate surgical technique and use drains as necessary.

SYMBOLS:

MATL: Material
QTY: Quantity
SS, SST: Stainless Steel
MADE IN: Made in <<Country>>
NON-STERILE: Supplied non-sterile

Manufacturer

Authorized Representative in the European Community

Biomet Trauma
 56 East Bell Drive
 P.O. Box 587
 Warsaw, Indiana, 46581 USA
 Tel: 1-800-348-9500
 Fax: 574-372-3968

Biomet UK LTD.
 Waterton Industrial Estate
 Bridgend CF31 3XA
 UK
 Tel: +44 (0) 1656 655221
 Fax: +44 (0) 1656 645454

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