EBA 21

ⓁⒹ	Bedienungsanleitung	5
ⓁⒺ	Operating Instructions	17
ⓁⒻ	Mode d'emploi	29
ⓁⒼ	Istruzioni per l'uso	41


Fig. 1


Fig. 2 EBA 21

**EG-Konformitätserklärung
EC Conformity Declaration
Déclaration de conformité CE
Dichiarazione di conformità alle norme CEE**

Andreas Hettich GmbH & Co. KG • Föhrenstraße 12 • D-78532 Tuttlingen • Germany

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Contents

1	Use according to specification	19
2	Residual risks	19
3	Notes on safety	19
4	Explanations of symbols and references	20
4.1	Symbol used on the machine	20
4.2	Symbols used in the operating instructions	20
5	Delivery checklist	20
6	Initial operation	21
7	Opening and closing the lid	21
7.1	Opening the lid	21
7.2	Closing the lid	21
8	Installation and removal of the rotor	21
9	Loading the rotor	22
10	Aerosol tight sealing of angle rotors	22
11	Control and display elements	22
11.1	Control panel symbols	22
11.2	Control panel pushbuttons (keys)	22
11.3	Adjustment possibilities	23
12	Entering centrifugation parameter	23
13	Programming	23
13.1	Programme input/alteration	23
13.2	Programme recall	23
14	Centrifugation	24
14.1	Centrifugation with pre-set time	24
14.2	Continuous run	24
14.3	Short-term centrifugation	24
15	Emergency Stop	24
16	Acoustic Signal	25
17	Recall hours of operation	25
18	Relative centrifugal force (RCF)	25
19	Centrifugation of materials with higher density	25
20	Rotor Identification	26
21	Emergency release	26
22	Maintenance and servicing	26
22.1	Centrifuge	26
22.2	Rotors and Attachments	26
22.3	Autoclaving	27
22.4	Centrifuge containers	27
23	Faults	27
24	Change mains input fuse	27

25	Acceptance of the centrifuges for repair	28
26	Anhang / Appendix	52
26.1	Technische Daten / Technical specification	52
26.2	Rotoren und Zubehör / Rotors and accessories.....	53

1 Use according to specification

The machine on hand is a centrifuge that is exclusively conceived for the separation of substances or respectively substance mixtures with a density of max. 1.2 kg/dm³ and which is therefore only designated for this application. A different use or application over and above this is deemed not in accordance with the specifications. The company Andreas Hettich GmbH & Co. KG undertakes no liability for damages resulting therefrom.

Belonging to the application according to specification is also the observance of all references contained in the Instruction Manual and compliance with the inspection and maintenance works.

2 Residual risks

The machine is constructed according to the state of the art and the recognized technical safety regulations. However during your application of the device hazards towards life and limb of the user or third parties or respectively impairments to the device or other material assets can arise. The machine is only to be used for the specified applications and only in an impeccable technical safety condition.

Disturbances that can interfere with the safety are to be immediately rectified.

3 Notes on safety

- **Before the initial operation of your centrifuge you should read and pay attention to the operating instructions.**
Only personnel that has read and understood the operating instructions are allowed to operate the device.
- Along with the operating instructions and the legal regulations on accident prevention, you should also follow the recognised professional regulations for working in a safe and professional manner.
These operating instructions should be read in conjunction with any other instructions concerning accident prevention and environmental protection based on the national regulations of the country where the device is to be used.
- This centrifuge is a state-of-the-art piece of equipment which is extremely safe to operate.
 - However, it can lead to danger for users or others if used by untrained staff, in an inappropriate way or for a purpose other than that it was designed for.
- The centrifuge should be installed on a good, stable base.
- Before using the centrifuge absolutely check the rotor for firm placement.
- Ensure that no persons and / or dangerous materials are located within a safety zone of 300 mm around the equipment when a centrifuge is running.
- The centrifuge must not be moved or knocked during operation.
- In case of fault or emergency release, never touch the rotor before it has stopped turning.
- To avoid damage due to condensate, when changing from a cold to a warm room the centrifuge must either run hot in the cold room for 30 minutes, or warm up for 3 hours in the warm room, before connecting to the mains.
- Load centrifuge rotor evenly. All positions on rotor must be filled.
- Do not fill centrifuge containers inside the centrifuge.
- Centrifuge containers must not be filled beyond the capacity specified by the manufacturer.
- Standard centrifuge containers of glass will not stand RCF values exceeding 4000 (DIN 58970, pg. 2)
- Centrifuge containers may only be centrifuged with accessories (reducing adapters, frames, suspensions, etc.) authorised by the manufacturer (see section "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories").
- When centrifuging with maxim revolutions per minute the density of the materials or the material mixtures may not exceed 1.2 kg/dm³.
- The centrifuge may only be operated when the balance is within the bounds of acceptability.
- The centrifuge must not be operated in areas subject to danger of explosions.
- The centrifuge must not be used with:
 - inflammable or explosive materials
 - materials that react with one another producing a lot of energy.
- If users have to centrifuge hazardous materials or compounds contaminated with toxic, radioactive or pathogenic micro-organisms, they must take appropriate measures.
Without additional proceedings (like an additional bioseal between bucket and lid of bucket or angle rotor with a special bioseal between rotor and lid) a centrifuge is not a biosafety system in accordance to the regulation EN 61010-2-020. In the case of material belonging to risk group II (see the World Health Organisation's "Laboratory Biosafety Manual") they should employ a biosafety system. Under this system small drips and aerosols are prevented from escaping by a bioseal (packing ring) located between the hanger and the lid. Centrifuge containers with special screw caps, as obtainable through trade suppliers, can also be used for hazardous substances.
In the case of materials from the higher risk groups greater safety provision is required than the arrangements described above. In a biosafety system, centrifuge containers with special screw caps must be used.
- For further details of available biosafety systems see section "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories". If in doubt, you should obtain relevant information from the manufacturer.
- The centrifuge must not be operated with highly corrosive substances which could impair the mechanical integrity of rotors, hangers and accessories.

- Any rotors, hangers or accessories showing clear signs of corrosion or mechanical defects must not be used for centrifuging.
- Repairs must only be carried out by personnel authorised to do so by the manufacturer.
- Only original spare parts and original accessories licensed by the Hettich company are allowed to be utilised.
- This centrifuge is classified in Germany as a Group 3 device according to the *Medizinische Geräteverordnung MedGV* (the regulations on medical equipment).
- It conforms to safety regulations based on: IEC 61010-1, IEC 61010-2-020
- The safe operation and reliability of the centrifuge can only be guaranteed if:
 - the centrifuge is operated in accordance with the operating instructions,
 - the electrical installation on the site where the centrifuge is installed conforms to the demands of IEC stipulations,
 - prescribed tests to BGV A1, BGR 261 are carried out by an expert.
- With centrifuges for robotic use please pay attention to the notes of the key operated switch.

No claim under guarantee will be considered by the manufacturer unless the above instructions have been adhered to.

4 Explanations of symbols and references

4.1 Symbol used on the machine

Attention, general hazard area.

Before using the centrifuge implicitly read the operating instructions and pay attention to the safety relevant references!

4.2 Symbols used in the operating instructions

This symbol refers to important circumstances.

This symbol refers to safety relevant warnings and indicates possibly dangerous situations.

The non-adherence to these warnings can lead to material damage and injury to personal.

5 Delivery checklist

- 1 Connecting cable
- 1 Hex. pin driver
- 1 Operating instructions
- 1 Notes on moving the equipment safely

The rotor(s) and associated accessories are included in the delivery in the quantity ordered.

6 Initial operation

- According to the laboratory instrument standards EN 61010-2-020 an emergency switch to separate power supply in the event of a failure must be installed in the building electrical system. This switch has to be placed remote from the centrifuge, preferred outside of the room in which the centrifuge is installed or near by the exit of this room.
- Remove the transportation safety device from the bottom of the housing, see sheet "Transportation safety device"
- Position the centrifuge in a stable and level manner in a suitable place. When setting up the equipment, care should be taken to provide the required safety area of 300 mm around the centrifuge in accordance with IEC 61010-1.

No persons and / or dangerous materials may be located within the safety zone when the centrifuge is running.

- Do not place any object in front of the ventiduct. Keep a ventilation area of 300 mm around the ventiduct.
- Check whether the mains voltage tallies with the statement on the type plate.
- Connect the centrifuge with the connection cable to a standard mains socket. For connection ratings refer to Chapter "Appendix, Technical data/Technical specification".
- Turn on the mains switch. Switch position "I". The machine type and program version will be displayed and the LEDs light up. **⏏ OPEN ⏏ OEFFNEN** will be displayed after 8 seconds. The LED **⏏** lights up.
- Open the lid. The last used centrifuge data will be displayed.
- Remove the transportation safety device in the centrifuge, see sheet "Transportation safety device".

7 Opening and closing the lid

7.1 Opening the lid

- Rotate the turning handle on the front panel to the left. The **⏏** LED will now go off.
- Open the lid.

The lid can only be opened when the centrifuge is switched on and the rotor is at rest. If it cannot be opened under these circumstances, see the section on "Emergency release".


7.2 Closing the lid

- Place the lid and lightly press down the front edge of the lid.
- Turn the handle on the front panel to the right. The **⏏** LED will come on.

Do not bang the lid shut.

8 Installation and removal of the rotor


- Clean the motor shaft (B) and the collet chuck of the rotor (A) and then lightly grease the motor shaft. Particles of dirt between the motor shaft and the collet chuck prevent perfect positioning of the rotor and cause erratic running.
- Position the rotor vertically on the motor shaft and press down as far as the stop. The rotor must audibly engage in the collet chuck.
- Tighten the rotor tension nut with the supplied wrench by turning in a clockwise direction.
- Check the rotor for firm seating.
- Loosening the rotor: Loosen the tensioning nut by turning anti-clockwise. Turn the tensioning nut until the rotor is able to be lifted from the motor shaft.

9 Loading the rotor

- With swing-out rotors all rotor positions must be lined with **identical** hangers. Certain hangers are marked with the number of the rotor position. These hangers may only be used in the respective rotor position.
- The rotors and hangers may only be loaded symmetrically. For authorised combinations see Chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories".
- In containers with rubber inserts, the same number of rubber inserts must always be among the centrifuge containers.
- Always fill the centrifuge containers outside of the centrifuge.
- No liquid should be allowed to enter the centrifugal chamber during filling and swinging out of the hangers.
- In order to maintain the weight differences within the centrifuge container as marginal as possible, a consistent fill level in the containers is to be heeded.
- The weight of the permissible filling quantity is specified on each rotor. This weight may not be exceeded.

10 Aerosol tight sealing of angle rotors

To ensure aerosol sealing the lid of an aerosol sealed angular rotor must be firmly closed.

This is performed with the help of the supplied spanner that is inserted through the borehole in the rotary handle. If the rotary handle does not have a borehole then the lid must be firmly closed by turning in clockwise direction.

For available aerosol tight angle rotors, please see chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories".

11 Control and display elements

Please turn to the rear side of this cover sheet.

Fig. 2: Display and control panel

11.1 Control panel symbols

Unbalance display. The unbalance display lights up if the permissible weight difference within the rotor loading has been exceeded.

Rotation display. The rotation display lights up during the centrifugation run as long as the rotor is turning.

Lid display. The lid display lights up if the lid is closed and locked. During the centrifugation run the lid display is extinguished.

Operating errors or any occurring disturbances are symbolised in the display (see Chapter "Faults").

11.2 Control panel pushbuttons (keys)

Selection control key for selection of specific parameter.
The subsequent parameter is selected by every further keystroke.

By pushing the arrow keys the value of a previously selected parameter is increased or decreased. By keeping an arrow key depressed the value is increased or decreased with increasing speed.

Start the centrifugation run. The rotation display  lights up.
Acquisition of input and changes.

Stop the centrifugation run.
The rotor runs down with pre-selected brake step. The LED in the key lights up until the rotor stops.
Pushing the key twice triggers the EMERGENCY STOP.

Selecting the RCF display.
The LED in the key lights up.




Short-term centrifugation.
The centrifugation run is effected as long as the key is held down.




Storing and recalling programmes.
The subsequent programme is displayed by every further keystroke (1 - 2 - 3 - #).


11.3 Adjustment possibilities

t/min	Adjustable from 1 - 99 min. in 1-min. intervals.
t/sec	Adjustable from 1 - 59 sec. in 1-sec. intervals.
∞	Continuous run.
RPM	Revolutions per minute. A numerical value from 500 RPM up to the maximum speed of the rotor can be set. Maximum speed of the rotor, see Chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories". Adjustable from 500 RPM to 10000 RPM in steps of 10 and from 10000 RPM up to the maximum speed in steps of 100.
RAD/mm	Centrifugation radius. Input in mm. For centrifugation radius see Chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories". The input of the radius is only possible if the RCF display is selected (LED in the key  lights up).
RCF	Relative Centrifugal Acceleration. A numerical value can be set, which gives a speed between 500 RPM and the maximum speed of the rotor. Adjustable up to 100 in intervals of 1, and from 100 in intervals of 10. The RCF value is automatically rounded up or rounded down with regard to the RPM interval. The input of the RCF is only possible if the RCF display is selected (LED in the key  lights up).
	Starting steps 1 - 9. Step 9 = shortest starting time, Step 1 = longest starting time.
	Brake steps 0 - 9. Step 9 = shortest run-down time, Step 1 = long run-down time, Step 0 = longest run-down time (brakeless run-down).
PROG	Programme Number. 3 programmes can be stored (programme positions 1 - 2 - 3). The programme position # serves as temporary storage for altered adjustments.

12 Entering centrifugation parameter






 If no key is pressed for 8 seconds long after the selection or during the input of parameters, the previous values will be shown in the display. The input of parameter then has to be executed again.

- Select the RPM or RCF display with the key . If the RCF display is selected, the LED in the key  lights up.
- Select the desired parameter with the key , and adjust with the keys  .
- In order to adjust the symbol ∞, the parameters **t/min** and **t/sec** must be set to zero with the key .
- After input of all parameters, press the key  in order to store the adjustments on the programme position #. As confirmation, ***** ok ***** will be displayed for a short period.


 The data on the programme position # will be overwritten with every input of parameters and pressing of the key .

13 Programming

13.1 Programme input/alteration

- Select the RPM or RCF display with the key . If the RCF display is selected, the LED in the key  lights up.
- Select the desired parameter with the key , and adjust with the keys  .
- Select the parameter **PROG** with the key , and set the desired programme position with the keys  . The LED in the key  lights up.
- Press the key  in order to store the adjustments to the desired programme position. As confirmation, ***** ok ***** will be displayed for a short period.

13.2 Programme recall

- Select the desired programme position by pressing the key .
- The centrifugation data of the selected programme position will be displayed.
- The parameters can be checked by pressing the key .

14 Centrifugation

If the permissible weight difference within the rotor loading has been exceeded, the drive shuts down during the start-up, the unbalance display lights up, and **IMBALANCE** is displayed.

If the speed in the selected program is higher than the maximum speed of the rotor, it will not be possible to start a centrifugation run. **N > ROTOR MAX** will be displayed (see Chapter "Faults").

A centrifugation run can be stopped at any time by pushing the key **STOP**.

All parameters can be selected and altered during the centrifugation run (see Chapter "Entering centrifugation parameter").

You can switch-over at any time between the RPM and RCF display with the key **RCF**. The input of the centrifugation radius is necessary if you are working with the RCF display.

If **OPEN OEFFNEN** is displayed, a further operation of the centrifuge is only possible after opening the lid once.

If **R xx n-max xxxxx** is displayed, then no centrifugation run has taken place as the rotor was changed beforehand, refer to Chapter "Rotor Identification".

- Turn on the mains switch. Switch position I.
- Load the rotor and close the centrifuge lid.

14.1 Centrifugation with pre-set time

- Adjusting time or recall a programme with pre-set time (see Chapter "Entering centrifugation parameter" or "Programme recall").
- Press the key **START**. The rotation display  lights up as long as the rotor is turning.
- After expiration of the time or with truncation of the centrifugation run by pushing the key **STOP**, the run-down is effected with the selected brake step. The brake step is displayed.

During the centrifugation run the rotational speed of the rotor or the subsequently resulting RCF value, the sample temperature (only in centrifuges with cooling) and the remaining time will be displayed.

14.2 Continuous run

- Adjusting the symbol ∞ or recall a continuous run programme (see Chapter "Entering centrifugation parameter" or "Programme recall").
- Press the key **START**. The rotation display  lights up as long as the rotor is turning. The time metering begins at 00:00.
- Press the key **STOP** in order to stop the centrifugation run. The run-down is effected with the selected brake step. The brake step is displayed.

During the centrifugation run the rotational speed of the rotor or the subsequently resulting RCF value, the sample temperature (only in centrifuges with cooling) and the expired time will be displayed.

14.3 Short-term centrifugation

- Hold down the key **IMPULS**. The rotation display  lights up as long as the rotor is turning.
- Let go of the key **IMPULS** again in order to stop the centrifugation run. The run-down is effected with the selected brake step. The brake step is displayed.

During the centrifugation run the rotational speed of the rotor or the subsequently resulting RCF value, the sample temperature (only in centrifuges with cooling) and the expired time will be displayed.

15 Emergency Stop

- Press the key **STOP** twice.

With Emergency Stop the run-down is effected with brake step 9 (shortest run-down time). Brake step 9 is displayed. If brake step 0 was pre-selected, the run-down time is technically longer than with brake step 9.


16 Acoustic Signal

The acoustic signal sounds:

- Upon the appearance of a disturbance in 2 second intervals.
- After completion of a centrifugation run and rotor standstill in 30 second intervals.


The acoustic signal is stopped by opening the lid or pressing any key.

The signal after completion of the centrifugation run can be activated or deactivated in the following manner, if the rotor is at standstill:

- Hold down the key  for 8 seconds.
After 8 seconds, **SOUND / BELL** appears in the display.
- Set **OFF** or **ON** with the key  or .
- Press the key  in order to store the setting.
As confirmation, ***** ok ***** will be displayed for a short period.

17 Recall hours of operation

Recall hours of operation is only possible during rotor standstill.

- Hold down the key  for 8 seconds.
After 8 seconds, **SOUND / BELL** appears in the display.
- Press the key  once again.
The centrifuge's hours of operation (**CONTROL:**) are displayed.
- Press the key  or  to exit the hours of operation recall.

18 Relative centrifugal force (RCF)

The relative centrifugal force (RCF) is given as a multiple of the acceleration of gravity (g). It is a unit-free value and serves to compare the separation and sedimentation performance.

These values are calculated using the formula below:

$$\text{RCF} = \left(\frac{\text{RPM}}{1000} \right)^2 \times r \times 1,118 \quad \Rightarrow \quad \text{RPM} = \sqrt{\frac{\text{RCF}}{r \times 1,118}} \times 1000$$

RCF = relative centrifugal force

RPM = rotational speed (revolutions per minute)

r = centrifugal radius in mm = distance from the centre of the turning axis to the bottom of the centrifuge.
For more on the centrifugal radius see the chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories".

The relative centrifugal force (RCF) stands in relation to the revolutions per minute and the centrifugal radius.

19 Centrifugation of materials with higher density

The rotors are designed to centrifuge substances up to a maximum mean homogenous density of 1.2 kg/dm³ when rotating at the stated speed.

Denser substances must be centrifuged at lower speed.

The permissible speed can be calculated using the following formula:

$$\text{Reduced speed (n}_{\text{red}}) = \sqrt{\frac{1.2}{\text{Greater density}}} \times \text{Rated speed}$$

e.g.: RPM 4000, density 1.6 kg/dm³

$$n_{\text{red}} = \sqrt{\frac{1.2}{1.6}} \times 4000 = 3464 \text{ RPM}$$

If in doubt you should obtain clarification from the manufacturer.

20 Rotor Identification

After every start of a centrifugation run the rotor utilised is identified.

After a change of rotor the drive switches off and the rotor code (R xx) as well as the maximum rotational speed (n-max=xxxxx) of the rotor are displayed.

A further operation of the centrifuge is only possible after a single opening of the lid.

If, following a rotor change, the maximum speed of the rotor is less than the set speed, the speed is limited to the maximum speed of the rotor.

21 Emergency release

The lid cannot be opened in the event of a power failure. An emergency release has to be executed by hand.

For emergency release disconnect the centrifuge from the mains.
Open the lid only during rotor standstill.

Please turn to the rear side of the cover sheet.

- Position the centrifuge on the table sufficiently far forward so that the hole (Fig. 1, A) for the emergency release on the base of the casing can be reached.
- Insert the hexagonal pin driver into this hole, push it upwards and at the same time rotate the lid-lock handle to the left.
- Open the lid.

22 Maintenance and servicing

Pull the mains plug before cleaning.

Before any other cleaning or decontamination process other than that recommended by the manufacturer is applied, the user has to check with the manufacturer that the planned process does not damage the device.

- Cleaning agents and disinfectants which lie in the pH range 5 – 8 are to be utilised. Alkaline cleaning agents with a pH value > 8 are to be avoided.
- In order to prevent appearances of corrosion through cleaning agents or disinfectants, the application guide from the manufacturer of the cleaning agent or disinfectant are absolutely to be heeded.

22.1 Centrifuge

- Regularly clean the centrifuge housing and the centrifugal chamber; and if necessary, clean with soap or a mild cleaning agent and water. For one thing, this services purposes of hygiene, and it also prevents corrosion through adhering impurities.
- In the event of condensation water formation, dry the centrifugal chamber by wiping out with an absorbent cloth.
- If infectious materials penetrates into the centrifugal chamber this is to be disinfected immediately. For surface disinfection we recommend Bacillol Plus from the company Bode Chemicals in Hamburg. After use the disinfection medium Bacillol Plus must be thoroughly removed with water.
- Lightly grease the rubber seal of the centrifugal chamber after every cleaning.

22.2 Rotors and Attachments

- Rotors and accessory parts must be regularly cleaned with soap or a mild cleaning agent and water in order to prevent corrosion and changes of material. Cleaning is recommended at least once a week, even better after every usage.
- If the rotor or accessory parts are contaminated by pathogenic or radioactive material, a suitable cleaning has to be executed. For disinfection we recommend Helipur H plus N from the company B. Braun Melsungen. For the removal of radioactive material we recommend deco neutracon from the company Decon Laboratories Limited.
- The rotors and accessory parts must be dried immediately after cleaning.
- Angle rotors, container and hanger made of aluminium are to be lightly greased after drying using acid-free grease, e.g. vaseline.
- With aerosol-sealed rotors and bio safety systems (see Chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories") the sealing rings are to be checked and cleaned regularly (weekly). The sealing ring is to be replaced immediately upon indication of crack formation, embrittlement or abrasive wear.
- In order to prevent corrosion as a result of moisture between the rotor and the motor shaft, the rotor should be disassembled and cleaned at least once a month, and the motor shaft should be lightly greased.
- The rotors and accessory parts are to be checked on a monthly basis for corrosion damage.

Rotors and attachments may no longer be utilised upon indication of wear and tear or corrosion.

22.3 Autoclaving

Swing-out rotors, angle rotors made of aluminium and accessory parts can be autoclaved at 121°C (20 min).

	The lids of the rotors and containers must be removed prior to autoclaving. Autoclaving accelerates the ageing process of plastics. In addition, autoclaving may discolour plastics.
---	---

22.4 Centrifuge containers

- With leakiness or after the breakage of centrifuging containers broken container parts and leaked centrifugation material are to be completely removed.
- The rubber inserts as well as the plastic sleeves of the rotors are to be replaced after a glass breakage.

	Remaining glass splitters cause further glass breakage!
---	---

- If this concerns infectious material, a disinfection process is to be executed immediately.

23 Faults

If the fault cannot be rectified according to the faults table, Hettich customer services must be informed.


Please state the type of centrifuge and the factory serial number. Both values are visible on the centrifuge type plate.

- MAINS RESET:
- Mains switch OFF for longer than 10 secs.
 - Mains switch ON.

Message / fault		Cause	Remedy
No display	---	- No voltage. - Overvoltage protection tripped out.	- Check supply voltage. - Mains switch ON.
TACHO - ERROR	01	- Faulty speedometer.	- Open lid.
	02	- No rotor installed. - Defective motor, frequency converter or drive.	- Turn rotor manually. - MAINS-RESET when power is switched on, rotor should turn.
IMBALANCE	---	Imbalance about motor axis through weight differential in rotor assembly.	- Open lid. - Correct imbalance.
CONTROL - ERROR	04, 06 - 09	Error in lid locking or lid closure.	- Open lid. - MAINS-RESET
N > MAX	05	Rotation too fast	
N < MIN	13	Rotation too slow	
ROTORCODE	10	Incorrect rotor coding	
MAINS INTERRUPT	---	Power failure, centrifuging not properly completed	- Open lid. - Push (START) button.
VERSIONS-ERROR	12	Mismatch between electronic components	- Open lid. - MAINS-RESET
CONTROL-ERROR	21 - 27	Error / defective control unit	
SER I/O - ERROR	30 - 38	Error / defective interface	
FU / CCI - ERROR	60 - 83	Error / defective motor control	
N > ROTOR-MAX	---	Speed in the selected program greater than the maximum speed of the rotor.	- Check the set speed. - Reduce the set speed

24 Change mains input fuse

	Switch off the mains switch and separate the centrifuge from the mains!
---	---

The fuse holder (A) with the mains input fuses is located next to the mains switch.

- Remove the connecting cable from the machine plug socket.
- Press the snap-fit (B) against the fuse holder (A) and remove.
- Exchange defective mains input fuse.
- Reinsert the fuse holder until the snap-fit clicks shut.
- Reconnect the centrifuge to the mains supply.

25 Acceptance of the centrifuges for repair

If the centrifuge is returned to the manufacturer for repair, it must be decontaminated and cleaned to protect persons, environment and material.

We reserve the right to accept contaminated centrifuges.

Costs incurred for cleaning and disinfection are to be charged to the customer.

We ask for your understanding in this matter.

You will find further information on our homepage: <http://www.HettichLab.com>

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






1) Sample temp. rise during maximum speed and 1 hour running time






1416		Reduzierung / adapter							
 <p>Winkelrotor 12-fach / Angle rotor 12-times</p> <p>∠ 35°</p>									
	---	1056		1054					
	Röhrchen / tubes								
	0507		0509	Vacutainer	Vacutainer				
									
Kapazität / capacity	ml	15	5	15	1,6 – 5,0	4 - 7			
Maße / dimensions Ø x L	mm	17 x 100	13 x 75	17 x 120	13 x 75	13 x 100			
Anzahl p. Red./number p. adapter		---	---	---	---	---			
Anzahl p. Rotor / number p. rotor		12	12	6	12	12			
Drehzahl / speed	RPM	6000	6000	6000	6000	6000			
RZB / RCF	⁴⁾	4105	3421	4105	3300	4105			
Radius / radius	mm	102	85	102	82	102			
 9 (97%)	sec	18	18	18	18	18			
 9	sec	23	23	23	23	23			
 0	sec	391	391	391	391	391			
Probenerwärmung/Sample temp. rise	K ¹⁾	12	12	12	12	12			

1112		Reduzierung / adapter							
 <p>Winkelrotor 12-fach / Angle rotor 12-times</p> <p>∠ 45°</p>		2031	2023	2024					
									
	Röhrchen / tube ³⁾								
	0536	2078							
									
Kapazität / capacity	ml	2,0	1,5	0,5	0,8	0,4	0,2		
Maße / dimensions Ø x L	mm	11 x 38	8 x 30	8 x 30	8 x 45	6 x 45	6 x 18		
Anzahl p. Red./number p. adapter		---	---	---	---	---	---		
Anzahl p. Rotor / number p. rotor		12	12	12	12	12	12		
Drehzahl / speed	RPM	18000	18000	18000	18000	18000	18000		
RZB / RCF		23907	23907	23907	23907	23907	23907		
Radius / radius	mm	66	66	66	66	66	66		
 9 (97%)	sec	14	14	14	14	14	14		
 9	sec	18	18	18	18	18	18		
 0	sec	172	172	172	172	172	172		
Probenerwärmung/Sample temp. rise	K ¹⁾²⁾	15	15	15	15	15	15		

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

- 1) Sample temp. rise during maximum speed and 1 hour running time
- 2) If the closed rubber seal (order no. E612) is installed, the sample temperature rise will increase.
- 3) If the rubber seal with ventiduct (order no. E1877) is installed, only centrifuge containers made of plastic may be used.
- 4) Centrifuge containers of glass will not stand RCF values exceeding 4000



1095	Reduzierung / adapter							
		2031	2023	2024				
	Winkelrotor 24-fach / Angle rotor 24-times  45° aerosoldicht / aerosol proof	---						
Röhrchen / tube ³⁾								
								
Kapazität / capacity	ml	2,0	1,5	0,8	0,5	0,4	0,2	
Maße / dimensions	∅ x L	mm	11 x 38	11 x 38	8 x 45	8 x 30	6 x 45	6 x 18
Anzahl p. Red./number p. adapter		---	---	---	---	---	---	
Anzahl p. Rotor / number p. rotor		24	24	24	24	24	24	
Drehzahl / speed	RPM	15000	15000	15000	15000	15000	15000	
RZB / RCF		21885	21885	21480	21480	21480	21480	
Radius / radius	mm	87	87	87	87	87	87	
 9 (97%)	sec	32	32	32	32	32	32	
 9	sec	32	32	32	32	32	32	
 0	sec	690	690	690	690	690	690	
Probenerwärmung/Sample temp. rise	K ¹⁾²⁾	18	18	18	18	18	18	

1089	Reduzierung / adapter							
		2031	2023	2024				
	Winkelrotor 30-fach / Angle rotor 30-times  45° aerosoldicht / aerosol proof	---						
Röhrchen / tube ³⁾								
								
Kapazität / capacity	ml	2,0	1,5	0,8	0,5	0,4	0,2	
Maße / dimensions	∅ x L	mm	11 x 38	11 x 38	8 x 45	8 x 30	6 x 45	6 x 18
Anzahl p. Red./number p. adapter		---	---	---	---	---	---	
Anzahl p. Rotor / number p. rotor		30	30	30	30	30	30	
Drehzahl / speed	RPM	14000	14000	14000	14000	14000	14000	
RZB / RCF		21255	21255	20850	20850	20850	20850	
Radius / radius	mm	97	97	97	97	97	97	
 9 (97%)	sec	40	40	40	40	40	40	
 9	sec	40	40	40	40	40	40	
 0	sec	800	800	800	800	800	800	
Probenerwärmung/Sample temp. rise	K ¹⁾²⁾	19	19	19	19	19	19	

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






- 1) Sample temp. rise during maximum speed and 1 hour running time
- 2) If the closed rubber seal (order no. E612) is installed, the sample temperature rise will increase.
- 3) If the rubber seal with ventiduct (order no. E1877) is installed, only centrifuge containers made of plastic may be used.

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All rights reserved. No part of this publication may be reproduced without the written prior permission of the copyright owner.