Instructions for Use - Mode d'Emploi - Gebrauchsanleitung  
Istruzioni per l'Uso - Instrucciones de Uso

# Horizontal-Vertical Lumbar Valve Systems\*

\*Incorporating the original Hakim™ valve mechanism

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

**Sterile - Sterilized with ethylene oxide gas. Non pyrogenic. For single use only. Do not use open or damaged packages. Do not resterilize.**

**WARRANTY DISCLAIMER AND LIMITATION OF REMEDY**

**THERE IS NO EXPRESS OR IMPLIED WARRANTY OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, OR OTHER WARRANTY ON INTEGRA NEUROSCIENCES PRODUCTS DESCRIBED IN THESE INSTRUCTIONS FOR USE. UNDER NO CIRCUMSTANCES SHALL INTEGRA NEUROSCIENCES IMPLANTS BE LIABLE FOR MEDICAL EXPENSES OR ANY DIRECT, INDIRECT, INCIDENTAL, OR CONSEQUENTIAL DAMAGES OTHER THAN AS EXPRESSLY PROVIDED BY SPECIFIC LAW. NO PERSON HAS ANY AUTHORITY TO BIND INTEGRA NEUROSCIENCES IMPLANTS TO ANY REPRESENTATION OR WARRANTY CONCERNING INTEGRA NEUROSCIENCES IMPLANTS PRODUCTS EXCEPT AS SPECIFICALLY SET FORTH HEREIN.**

Descriptions and specifications appearing in Integra NeuroSciences printed matter, including this publication, are meant solely to generally describe the product at the time of manufacture and do not constitute any express warranties. Integra NeuroSciences hydrocephalus valve systems, components and accessories are designed to perform reliably over an extended period of time. However, any implanted component may have to be replaced at any time because of failure or for medical reasons, and the patient and/or the patient's family should be made aware of this possibility.

Note: After implantation, complete the patient identification card for the patient or family.

**Rx only.**

Product returns to Integra NeuroSciences: Contact your local Representative for instructions on return of products.

**DESCRIPTION**

The Horizontal-Vertical Lumbar Valve System is an implantable device that provides controlled intraventricular pressure and cerebrospinal fluid (CSF) drainage from the lumbar subarachnoid region to the peritoneal cavity. The system minimizes overdrainage by automatically compensating for changes in CSF hydrostatic pressure when the patient changes from a recumbent to an upright position.

**COMPONENTS**

The valve System consists of:

- a valve unit,
- a closed-tip lumbar catheter,
- a peritoneal catheter,
- a Tuohy needle,
- a stepdown connector,
- a straight connector,
- a suture clamp,
- an antechamber,
- a PTFE-coated guidewire in dispenser,
- a Luer-lock connector.

*Materials in contact with body tissues or body fluids are silicone elastomer with or without barium sulfate, stainless steel, ruby, polypropylene with barium sulfate, polycarbonate and PTFE.*

**Compatibility with Magnetic Resonance Imaging (MRI)**

Integra NeuroSciences tests indicate that MRI exposure does not affect pressure valve settings. No significant forces or temperature changes were noticed during MRI exposure.

**VALVE UNIT**

The Valve Unit (Fig. 1 - Page 2) contains two ball-in cone valve mechanisms: A spring-actuated (lower) pressure mechanism on the inlet side and a gravity-actuated (higher) pressure mechanism on the outlet side. The two valve modules are contained in a clear silicone elastomer housing with a clear inlet tubing and a white (radiopaque) outlet tubing. Four suture holes are provided on the valve housing. An arrow indicates the flow direction. Each unit has a color-coded band that represents a specific combination of horizontal and vertical pressure ranges. Each Valve Unit is identified by a serial number, this number should be recorded in the patient's file for future reference.

Each Horizontal-Vertical (H-V) Lumbar Valve System adjusts for changes in posture to maintain intraventricular pressure at appropriate levels in the standing and recumbent positions. When the Valve Unit is horizontal, the spring-actuated, ball-in-cone inlet valve alone controls the lower pressure limit. When the Valve Unit is vertical, the gravity-actuated, ball-in-cone outlet valve combines with the inlet valve to provide the higher pressure.

## PRESSURE RANGES

The horizontal pressure range should generally be selected according to the extent of the patient's ventricular enlargement. If enlargement is minor, higher horizontal pressures are generally applied. The vertical valve pressure should normally be selected on the basis of the patient's height from the base of the spine to the head. Six pressure ranges are available. Color codes, closing pressure ranges, and suggested applications of the system are given in the table.

**Table of Pressure Ranges**

<b>For Patients with Slight Ventricular Enlargement</b>				
	<b>Horizontal (Closing Pressure, mmH<sub>2</sub>O)</b>	<b>Vertical</b>	<b>Color code</b>	<b>Suggested Application</b>
903-325A	85 - 125	205 - 285	Blue	Small children
903-335A	85 - 125	265 - 365	Yellow	Tall children, adults
903-345A	85 - 125	325 - 445	Green	Tall adults
<b>For Patients with Extreme Ventricular Enlargement</b>				
	<b>Horizontal (Closing Pressure, mmH<sub>2</sub>O)</b>	<b>Vertical</b>	<b>Color code</b>	<b>Suggested Application</b>
903-320A	50 - 80	170 - 240	Red	Small children
903-330A	50 - 80	230 - 320	White	Tall children, adults
903-340A	50 - 80	290 - 400	Brown	Tall adults

**Note:** Each valve is individually tested for its horizontal and vertical closing pressures and color code accordingly. The closing pressure of each valve has been verified to be within the pressure range limits indicated on the package label at the time of manufacture.

## INDICATIONS

Horizontal-Vertical Lumbar Valve Systems are implantable devices used in the treatment of patients with communicating hydrocephalus to shunt cerebrospinal fluid (CSF) from the lumbar subarachnoid region to the peritoneal cavity. They provide controlled intraventricular pressure and CSF drainage in patients with hydrocephalus. The antechamber can be electively mounted in line with the Valve Unit to allow for CSF sampling or injections in the subarachnoid space.

## CONTRAINDICATIONS

The H-V Lumbar Valve System should not be used in patients with non-communicating hydrocephalus. Integra NeuroSciences hydrocephalus valve system components should not be used when infection exists along the shunt pathway (meningitis, ventriculitis, peritonitis, septicemia, bacteremia). Postponement of shunt or component implantation is advisable if infection is present anywhere in the body.

## ADVERSE REACTIONS

- In addition to risks associated with shunt system or component implantation, major complications include mechanical failure, shunt obstruction, tubing fracture, infection and CSF leakage along the shunt pathway. These complications may occur during the procedure or at any time thereafter. System disconnection can lead to catheter migration into the subarachnoid space or peritoneal cavity.
- System obstruction can lead to intracranial hypertension syndrome.
- Overdrainage may result in subdural hematomas, craniostenosis, intracranial hypotension syndrome, sunken fontanelles in infants or acquired tonsillar herniation.
- Low-grade bacterial colonization may cause recurrent fevers, anemia, splenomegaly, shunt nephritis or pulmonary hypertension.
- Possible complications of lumboperitoneal shunting include arachnoiditis, scoliosis, hyperlordosis, and back pain.
- Perforation of an abdominal viscus is a rare complication.

## PRECAUTIONS

This product is recommended for single use only. Do not implant explanted systems or components in another patient. *This product is sterilized with ethylene oxide. Do not use if the package is open or damaged. Use the device prior to the "Use Before" date on the package label.*

**Caution - Do not resterilize. Integra NeuroSciences will not be liable for any direct, indirect, incidental or consequential damages resulting from or related to resterilization.**

Observe patients closely for symptoms of shunt failure. Headache, irritability, vomiting, drowsiness or mental deterioration and other symptoms associated with elevated intracranial pressure may be signs of a nonfunctional shunt. Redness, tenderness or erosion of the skin along the shunt pathway may indicate an infection in the shunt system. These complications, along with CSF leakage or overdrainage, require prompt replacement of the shunt system or affected component(s).

- Have duplicates of each component available during implantation.
- Avoid nicking the tubing. Cover forceps and clamp jaws with silicone elastomer tubing.
- Do not lubricate valves. They will be lubricated adequately by the CSF.
- Avoid contaminants such as glove talc and contact with skin oils when handling tubing.
- Other than trimming catheters, modifications are not recommended.

## RECOMMENDED IMPLANTATION PROCEDURE

### Clinical Judgment

Incision size, size and type of nonabsorbable suture material, etc., should be selected according to the surgeon's experience and preference. All references "to solution" in these instructions, unless otherwise noted, refer to sterile, pyrogen-free saline or antibiotic solution. Procedure may be varied in accordance with the surgeon's clinical judgment.

### Modifications

The lumbar and drainage catheters can be trimmed to the desired length during implantation. However, the peritoneal drainage catheter should be trimmed only on the beveled (proximal) end. The distal end contains slits that aid in drainage and should not be trimmed. The lumbar catheter should be trimmed only at its distal end, the proximal end contains drainage holes necessary for proper function. **Other modifications are not recommended.**

### Intraoperative Testing - Patency Test

Integra NeuroSciences recommends testing of valve patency only. The peritoneal catheter and antechamber must be connected to the valve Unit in order to perform these tests. Before implantation, valve patency may be assessed using the following procedure:

- 1 - Place inlet tubing in solution.
- 2 - Attach syringe to distal tubing and aspirate solution to prefill the system.
- 3 - Hold proximal tubing end in solution, allowing distal end to hang freely. Patency is confirmed if solution drips from distal end.

## SUGGESTED PROCEDURE

### I. Preparation

Place the patient in the lateral position, left side up, flex the right (lower) leg and extend the left (upper) leg. Prepare and drape the lower back, hip, and abdomen. Administer local or general anesthesia. Make an incision above the iliac crest.

### II. Placement of the lumbar catheter

1. At the physician's option, make a 2 cm incision prior to Tuohy needle insertion.
2. Connect a syringe to the Luer-lock hub of the guidewire dispenser, and flush the guidewire with sterile, saline solution.

3. Insert the Luer-lock connector into the lumbar catheter. Connect a syringe to the connector and flush the catheter with sterile, saline solution.
4. Insert the guidewire, PTFE-coated tip first, into the catheter.
5. Insert the Tuohy needle, with obturator in place, into the subarachnoid space. Insert the needle bevel up for the cephalad placement of the catheter, bevel down for caudad placement.
6. To verify that the lumen of the needle has entered the subarachnoid space, withdraw the obturator slowly until CSF appears. If no CSF flows, reinsert the obturator and reposition the needle.
7. Verify that the needle is correctly positioned. Remove the obturator and insert the catheter and guidewire through the needle to desired lumbar level.
8. Remove the Luer-lock connector and withdraw the needle over the distal end of the catheter.  
**Caution:** To avoid catheter damage, do not withdraw the catheter after it has been inserted into the needle. If the catheter must be removed, withdraw needle and catheter simultaneously.
9. Remove the guidewire.  
**Caution:** To avoid catheter damage, do not withdraw the guidewire with the Tuohy needle still in place.
10. Clamp the lumbar catheter.
11. Anchor the lumbar catheter to the fascia using the butterfly suture clamp.
12. Using a tunneler or other suitable instrument, pass the lumbar catheter subcutaneously from the back to previous incision.
13. Reconnect the Luer-lock connector to the lumbar catheter.
14. Release the clamp on the lumbar catheter and remove air from the system by allowing CSF to flow through the catheter. Clamp the catheter again. **Do not reduce catheter length at this time.**

### III. Placement of the Valve Unit

**Caution:** To achieve correct compensation for postural pressure changes, the Valve Unit must be properly oriented and immobilized. To avoid misaligning the Valve Unit and twisting the tubings connect the Valve Unit to both catheters before positioning it.

1. Connect the antechamber integrated connector to the inlet side of the Valve Unit (clear tubing).
2. Using a blunt needle inserted in the tubing, fill the Valve Unit and the antechamber with sterile saline solution.
3. Trim the free end of the lumbar catheter to desired length. Connect the lumbar catheter to the tubing of the antechamber using the stepdown connector, place the small end in lumbar catheter. For easy connection, insert the guidewire through the stepdown connector and a short distance into the distal end of the catheter. Slide the catheter onto the connector and remove the guidewire.
4. Secure all connections with non absorbable sutures.
5. To check valve patency, release the clamp on the lumbar catheter. Hold the valve in the horizontal position until CSF flows through the valve. Reclamp the lumbar catheter.
- 6.

---

Placement of the Peritoneal Catheter: The peritoneum is opened using of any technique consistent with the surgeon's experience:

- Make a small incision in the peritoneum. Insert the peritoneal catheter square cut (not beveled) end into the incision. Advance the peritoneal catheter into the peritoneal cavity until the catheter reaches or extends beyond the midline. After the peritoneal catheter is positioned, close the peritoneal incision.
  - If there is no contraindication such as previous abdominal surgery or extreme obesity, a peritoneal trocar might be used to reduce the total duration of surgery.
7. Trim the free (beveled) end of the peritoneal catheter to desired length. Connect the peritoneal catheter to the outlet side (white tubing) of the Valve Unit using the straight connector.
  8. Secure the catheter and outlet tubing to the connector with two non absorbable sutures.
  9. Position the Valve Unit so that the arrow on the inlet valve points toward the patient's feet. Using permanent suture material, secure the Valve Unit to the fascia, passing sutures through the four holes in the valve housing.
  10. Secure the antechamber to the fascia using the four pre-cut suture holes on the reservoir base.
  11. Remove the clamp from the lumbar catheter. Close the superficial fascia and skin in a routine manner.  
**Caution:** When sampling or infusing solution using the antechamber, a 22 G or smaller needle should be used, after the skin has been cleansed with antiseptic. To prevent leakage, the antechamber should always be punctured with the needle held perpendicular to the dome.

### **Valve Unit Replacement in Integra NeuroSciences or other Competitive Valve Systems**

Make an incision next to the Valve Unit. Disconnect the lumbar catheter and verify CSF flow. If the catheter is patent, disconnect the peritoneal catheter and proceed to check for distal patency. If both catheters are patent, remove the Valve Unit and replace it.

**Caution:** Some components manufactured by other companies incorporate slit valve mechanisms that will raise the pressure of the valve system. In this case or in doubt, replace components with equivalent supplied with H-V Lumbar Valve System.

Fill Valve Unit with solution and insert Valve Unit (note flow direction arrow) between proximal and distal catheters using straight and stepdown connectors. Tie and remove clamp.

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The wave logo and Integra NeuroSciences are trademarks of Integra LifeSciences Corporation.