Manuale di Istruzione Uso e Manutenzione
User and Maintenance Manual
Manuel d'Utilisation et d'Entretien
Gebrauchs- und Wartungsanleitung
Manual de Instrucciones de Uso y Mantenimiento

Elettrobisturi
Electrosurgical knife
Bistouri électrique
Elektroskalpell
Electrobisturí

carlo de giorgi s.r.l.

Dati modello**Model data / Donnés du modèle / Modelldaten / Datos del modelo**

Modello Model Modèle Modell Modelo	Alimentazione Power supply Alimentation Betriebs-spannung Alimentación	Anno di costruzione Year of manufacture Année de fabrication Baujahr Año de fabricación
640/00	<input type="checkbox"/> 230 V 50 Hz <input type="checkbox"/> 120 V 60 Hz	<input type="checkbox"/> 2007 <input type="checkbox"/> 2008 <input type="checkbox"/> 2009 <input type="checkbox"/> 2010 <input type="checkbox"/> 2011

**carlo de giorgi** s.r.l.

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



CONTENTS

1	STANDARDS AND GENERAL PRECAUTIONS	3
1.1	INTRODUCTION.....	3
1.2	DESCRIPTION OF THE APPLIANCE	3
1.3	REFERENCE STANDARDS	4
1.3.1	OBLIGATORY COMPLIANCE	4
1.3.2	VOLUNTARY COMPLIANCE.....	4
1.4	RESPONSIBILITIES OF THE CUSTOMER.....	4
1.5	REQUESTS FOR SERVICE AND SPARE PARTS	5
2	TECHNICAL FEATURES	5
3	COMPONENTS	5
4	OPERATION AND USE	6
4.1	OPERATION	6
4.2	OPERATOR	7
4.3	INTENDED USE.....	7
4.4	SAFETY WARNINGS.....	7
5	OPERATOR INSTRUCTIONS	8
5.1	CONTROLS	8
5.2	INSTRUCTIONS FOR USE	9
5.2.1	SETTING INTO OPERATION.....	9
5.2.2	NORMAL DEACTIVATION AND EMERGENCY DEACTIVATION PROCEDURES	9
5.2.3	RECOMMENDATIONS FOR USE	10
5.2.4	WARNINGS.....	10
6	MAINTENANCE	11
6.1	CLEANING	11
6.2	STERILIZATION.....	11
7	DIAGNOSTICS	12
8	ACCESSORIES	13
8.1	ELECTRODES	13
8.2	VETERINARY NEUTRAL ELECTRODE.....	13
9	DISPOSAL	13
10	ANNEXES	14
10.1	OUTPUT POWER DIAGRAMS – NOMINAL LOAD SETTING (1200 Ω)	14
10.2	DIAGRAMS SHOWING OUTPUT POWER IN RELATION TO LOAD – FULL POWER AND HALF POWER.....	15
10.3	ELECTRODE LIST	16

1 - STANDARDS AND GENERAL PRECAUTIONS

1.1 - INTRODUCTION

This manual is the property of CARLO DE GIORGI S.r.l. The contents of this manual cannot be reproduced or disclosed to other parties. All rights reserved.

1.2 - DESCRIPTION OF THE APPLIANCE

The appliance is a diathermic coagulator (radio frequency ESK) designed for use in medical, aesthetic (electro-epilation) and veterinary applications for minor surgery of short duration (from 10 to 20 seconds), with intermittent duty. This device must be used exclusively by persons who are properly trained and qualified to perform surgery and who have adequate experience in the use of electrosurgical knives. Before using the electrosurgical knife it is mandatory to have read the contents of this manual in their entirety and to know and understand the position and function of all the relative commands, paying particular attention to all the safety warnings and precautions prescribed herein.

Failure to comply with this requirement or incorrect use of the appliance can result in injury to the operator, the patient and other persons in the vicinity. This appliance is covered by the terms of the Medical Devices Directive 93/42/EEC and is classified as type IIb in accordance with annex IX, regulation 9.

IT IS STRICTLY FORBIDDEN TO OPEN OR REPAIR THE APPLIANCE; IN THE EVENT OF MALFUNCTIONS ALWAYS CONSULT CARLO DE GIORGI.

1.3 - REFERENCE STANDARDS

1.3.1 - Obligatory compliance

- EEC directive n. 93/42 – Medical Devices (D.Lgs. n° 46/1997)
- EEC directive n.73/23 – Low voltage Directive (Italian legislation: Law n° 791/1977, Decree Law n° 626/1996, Decree Law n° 277/97)
- EEC directive n.89/336 relative to Electromagnetic Compatibility (EMC) - (Italian legislation: Decree Law n° 615/1996)
- DPR 27.4.1955 n.547, "Norme per la prevenzione degli infortuni sul lavoro" (work accident prevention regulations)
- DPR 27.4.1956 n.303, "Norme generali per l'igiene del lavoro" (general work hygiene regulations)
- Italian Decree Law 12 August 1991 n.277 implementing EEC directives n.80/605; n.83/477; 86/188; 88/642 concerning the protection of workers against risks deriving from exposure to chemical, physical and biological agents in the workplace, as defined by art. 7 of Law 30 July 1990, n.212.
- Legislative Decree No. 151 of 25th July 2005 implementing EC Directives 2002/95/EC, 2002/96/EC and 2003/108/EC relative to reduction of the use of hazardous substances in electrical and electronic equipment as well as waste disposal.

1.3.2 - Voluntary compliance

- CEI EN 60601-1 (1998) Electromedical appliances - Part 1 – General safety regulations.
- CEI EN 60601-2-2 (2000) Electromedical appliances - Part 1 - Special regulations for the safety of high frequency apparatuses for electrosurgery.
- CEI EN 60601-1-2 (2001) Electromedical appliances - Part 1 – General safety regulations - Supplementary standard: electromagnetic compatibility – Prescriptions and tests.

1.4 - RESPONSIBILITIES OF THE CUSTOMER

The user must install the appliance in a suitable place equipped with an electrical system that complies with established legislation.

The appliance must be installed in a dry environment that is illuminated in accordance with the requirements of established legislation.

NOTE: The expression "established legislation" refers to the laws in force in the user's country.

1.5 - REQUESTS FOR SERVICE AND SPARE PARTS

For all electrical maintenance requirements, contact CARLO DE GIORGI di Baranzate, Via Tonale, 1, Milan, Italy.

In the case of faults, notify CARLO DE GIORGI technical personnel of the exact nature of the problem, specifying also the lot number shown on the label affixed to the rear of the appliance.

2 - TECHNICAL FEATURES

		640/00	
Voltage	V	230/120	Class II B app.
Frequency	Hz	50/60	
Peak power consumption	W	175	Type BF app.
Power output at 750 kHz/1200Ω	VA	50 ± 20%	
Delayed acting fuse T	A	1.25 (230V)	
Mass	kg	2.5	

The appliance features double electrical insulation (Class II) and a dedicated system protecting against stray currents and ground faults, with particular regard for the ground electrode on the patient (BF type appliance).

3 - COMPONENTS

The package contains:

- a) power supply unit;
- b) handpiece;
- c) pedal;
- d) cable with neutral cylindrical electrode;
- e) case with n° 8 electrodes;
- f) user instructions.

The appliance can be transported by hand.

4 - OPERATION AND USE

4.1 - OPERATION

The power supply/control unit must be placed on a surface next to the operator in such a way that:

- it is in a stable position with no risk of falling;
- the handpiece lead is not subject to excessive traction, thus avoiding the risk of breakage of internal wires or dragging of the power supply unit;
- the functional status of the appliance can be checked easily.

Connect handpiece (1) to connector (2) located on the rear of the appliance (by inserting the handpiece connector and turning the locking ring through about 1/4 of a revolution until it stops - do not use excess force); Connect electrode (3) to cable (4), which must be inserted in the "plate" socket (5) on the front panel; the appliance can be activated using the button on the handpiece (1); if you wish to activate the appliance using the pedal control (6), connect the pedal cable to connector (7) on the rear of the appliance. Finally, insert the plug on the power cable into a mains socket outlet.


Fig. 1

Now unscrew and remove the head, insert the required type of electrode and then refit the head, tightening it appropriately to ensure the best possible electrical continuity between the parts.

If the electrode cannot be inserted into the collet, extract the collet from its location and push in the electrode with force; when the electrode has been inserted, refit the resulting assembly into the handpiece and then refit the head as described above.

Always sterilize the handpiece and electrodes before use; the handpiece is specifically designed for sterilization in an autoclave.

Always switch OFF the appliance after use by setting the main power switch to "0".

4.2 - OPERATOR

The operator starts the appliance by means of a voluntary command via the button on the handpiece or the control pedal, thereby assuming direct control of its operation.

4.3 - INTENDED USE

The electrosurgical knife is designed for use in medical, dental and veterinary surgeries for minor surgery and micro cauterization for brief and intermittent interventions (20 seconds max.).

The appliance is not designed for any uses other than those specified by the manufacturer.

4.4 - SAFETY WARNINGS

- The appliance must be used exclusively by skilled and properly qualified personnel who are aware of the potential hazards.
- Do not use the appliance in the vicinity of inflammable materials or gases.
- Do not open or disassemble the appliance.
- Patients with pacemakers or stimulatory electrodes may be subject to disturbances or damage to the relative systems caused by the electrosurgical knife; check the compatibility with such devices with the appliance (seek cardiology information).
- Wear suitable work clothing and observe the safety regulations in force in the workplace.
- Always use the neutral electrode, which must be as close as possible to the area of surgical intervention, and must be gripped tightly by the patient. This procedure must be performed before you start setting up the appliance.
- The patient must not come into contact with metal parts that are connected to ground or that have reduced ground impedance (e.g. operating tables, stands, etc.).
- Avoid all skin-skin contact (e.g. between arms and body); this can be achieved if the patient is clothed or by interposing dry gauzes.
- If the appliance is used in conjunction with monitoring equipment, all the monitoring equipment electrodes must be kept as far away as possible from the intervention site.
- The use of needle type monitoring electrodes should be avoided wherever possible.

- The appliance cables must be arranged so that they are not in contact with the patient or with other cables.
- Use the lowest possible output power to obtain the desired results.
- Pay attention to the use of inflammable liquids or gases such as disinfectants, anaesthetics, solvents, etc.; the appliance is not designed for use in atmosphere subject to explosion hazards so ensure that any such substances have been completely eliminated or evaporated and dispersed before activating the appliance (ensure there are no residual substances in the area).
- For the above reasons, also pay attention to the presence of endogenous gases and materials impregnated with oxygen (e.g. cotton wadding, gauzes, etc.).
- Periodically inspect the accessories, with particular reference to the condition of insulation of cables and electrodes.

5 - OPERATOR INSTRUCTIONS

5.1 - CONTROLS

The appliance is equipped with the following control devices (Fig. 2):

- illuminated main power ON/OFF switch (1) (green);
- red AL indicator LED (alarm) (2);
- plate electrode socket (PLATE) (3);
- 3-position selector (4) with three function LEDs:
 - CUT (yellow LED)
 - CUT-COAG (orange LED)
 - COAG (blue LED)
- output power control potentiometer (5).


Fig. 2

5.2 - INSTRUCTIONS FOR USE

5.2.1 - Setting into operation

After performing the installation procedures, the appliance can be set into operation (Fig. 2):

- 1) fit the required electrode on the handpiece as indicated in § 4.1;
- 2) select the required function using selector (4);
- 3) switch on the appliance by means of the main power switch (1);
- 4) select the required power level by means of potentiometer (5);
- 5) instruct the patient to grip the neutral electrode;
- 6) press the button on the handpiece or press the pedal control; this action will start the audible signal and deliver power to the electrode attached to the handpiece.

Note that the appliance will assume alarm status after 20 seconds of continuous use; this condition is signalled by illumination of the red AL indicator (2) and a simultaneous intermittent beeping to inform the operator that the current operation must be terminated as rapidly as possible without prejudice to the patient; to continue using the appliance without the alarm condition it must be switched off for approximately 10 seconds (by setting the main power switch to OFF).

The appliance can also be utilized without the ground electrode; this procedure is permitted only for performing epilation or for low power operations (less than 1/3 of the maximum output power).

5.2.2 - Normal deactivation and emergency deactivation procedures

The appliance can be deactivated as follows:

- by releasing the button on the handpiece or the pedal control;
- by setting the main power switch to OFF;
- by unplugging the power cable from the mains socket outlet.

5.2.3 - Recommendations for use

- Use the level of currently strictly necessary in relation to the operational conditions and the conditions of the patient's tissue, keeping the knife in constant movement without exerting pressure.
- Keep the knife electrode tip clean.
- The CUT setting delivers higher current to allow cleaner incisions; the COAG setting is suitable for working on haemorrhagic areas or bleeding tissue in general.
- If you find that cutting power is insufficient or if fragments of tissue tend to adhere to the knife electrode, increase output power by turning potentiometer (5) clockwise.
- Sparking in the area of the electrode indicates that output power is too high; in this case reduce power by turning potentiometer (5) counter-clockwise.

5.2.4 - Warnings

BEFORE USING THE APPLIANCE CHECK ITS COMPATIBILITY WITH THE OTHER EQUIPMENT/APPLIANCES IN USE IN THE AREA; THE ELECTROSURGICAL KNIFE GENERATES ELECTROMAGNETIC DISTURBANCES DURING NORMAL OPERATION WHICH COULD INTERFERE WITH THE CORRECT OPERATION OF OTHER DEVICES IN THE AREA.

ALWAYS SWITCH OFF THE APPLIANCE WHEN YOU FINISH WORKING WITH IT AND FOLLOWING ALL STOPPAGES, IRRESPECTIVE OF THE MOTIVE.

IF YOU NOTICE MALFUNCTIONS STOP THE APPLIANCE IMMEDIATELY AND DISCONNECT IT FROM THE ELECTRICAL POWER SUPPLY. BEFORE RECONNECTING THE POWER SUPPLY FIND THE CAUSE OF THE FAULT AND REMEDY IT. IN THE CASE OF PROBLEMS, CONSULT CARLO DE GIORGI.

6 - MAINTENANCE

6.1 - CLEANING

The appliance does not require particular maintenance with the exception of regular cleaning of the working area: before cleaning the appliance always disconnect it from the power supply.

Use non flammable cleaning products and liquids that do not conduct electricity.

6.2 - STERILIZATION

Pay particular attention to cleaning and sterilizing the handpiece and electrodes, which can be treated in an autoclave for 13 minutes at 127°C after bagging.

The cable and neutral electrode must be cleaned with alcohol or sterilizing products.

IMPORTANT: do not sterilize the handpiece with liquids that conduct electricity.

7 - DIAGNOSTICS

If the output power level is too low when the appliance is set up to deliver higher power levels, check that the neutral electrode is properly in contact with the patient and properly connected to the appliance (incorrect application or bad contact).

If the appliance fails to operate check the fuse at the rear (8, Fig. 1). If the fuse must be changed refer to the specifications in the table in chapter 2 of this manual.

Before checking the fuse, set the main power switch to "0" (OFF) and unplug the mains power cable; check the condition of the fuse and replace if necessary. Once the fuse cap has been refitted you can insert the power cable plug into the mains socket outlet.

If malfunctions and faulty operation persist, unplug the appliance and consult CARLO DE GIORGI.

WARNING!

ALL WORK ON THE APPLIANCE MUST BE PERFORMED WITH
THE POWER SUPPLY DISCONNECTED

8 - ACCESSORIES

8.1 - ELECTRODES

Refer to the annexes for a list of available electrodes (see par. 10.3).

8.2 - VETERINARY NEUTRAL ELECTRODE

The appliance can be equipped with a plate electrode specifically designed for positioning beneath the animal in veterinary applications.

9 - DISPOSAL

The crossed out wheellie bin symbol indicates that the apparatus must not be disposed of with other waste.

The differentiated disposal of this apparatus is organized and managed by the manufacturer. When the user wishes to throw this apparatus away he must contact and follow the system adopted by the manufacturer to dispose of it correctly.


Differentiated disposal to enable possible recycling or environmentally compatible elimination, helps to limit undesirable effects on health and the environment and promotes the reuse and/or recycling of the materials which make up the apparatus.

Sanctions, in accordance with the relative legislation in force in the country of use, will be imposed on the holder if the apparatus is disposed of illegally.


10 - ANNEXES


10.1 - OUTPUT POWER DIAGRAMS – NOMINAL LOAD SETTING (1200 Ω)

CUTTING


CAUTERIZING


COMBINED



10.2 - DIAGRAMS SHOWING OUTPUT POWER IN RELATION TO LOAD – full power and half power

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