realcheck™ *MICRO 2*

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1. INTRODUCTION

1.1. Features of the REALCHECK MICRO 2

This completely automated digital sphygmomanometer uses an oscillometric measurement procedure that permits fast and reliable measurements of the systolic and diastolic pressure and the heart rate.

REALCHECK MICRO 2 uses Fuzzy Logic technology to measure arterial pressure. As the cuff deflates, the inner system takes a pressure reading based on oscillations. Once the reading has been taken, the rapid deflating process commences and the display shows the systolic and diastolic pressure together with the heart rates. The high precision has been clinically tested and the system has been specially designed for the utmost ease of use.

REALCHECK MICRO 2 features a large LCD display that permits clear reading of the functional status and cuff pressure throughout the measurement process. Before using, read the contents of this manual carefully and then store it in a safe place. Contact a physician for further queries concerning blood pressure and how to measure it.

Attention!

1.2. Important information about self-measurement

- Do not forget: **self-measurement means Control**, not diagnosis or treatment. Unusual values must always be discussed with the doctor. Any change in dosage of medications prescribed by the doctor should be discussed with him/her personally.
- The pulse display is **not** suitable for checking the frequency of heart pacemakers!
- In case of cardiac irregularity (Arrhythmia), measurements made with instrument should only be evaluated after a consultation with the doctor.

2. INFORMATION REGARDING BLOOD PRESSURE AND HOW TO MEASURE IT CORRECTLY

2.1. Hypertension

Hypertension is one of the most important vascular risk factors and is directly linked to the main causes of death in the western world (heart attack and ictus). Checking the arterial pressure must then be considered an indication of the health of the cardiovascular system. Like all biological variables, blood pressure is changeable and is influenced by numerous parameters which can be internal, (hormones, stress, anxiety etc.) or external (temperature, smoke, coffee etc.), to our organism; it indicates the pressure that the blood, pumped by the heart, exerts on the arteries: *maximum or systolic pressure – minimum or diastolic pressure*.

It tends to increase with age and is influenced by the seasons, the climate and the altitude.

By measuring it (in millimetres of Mercury – mmHg -) normal conditions can be established, as well a situation of hypertension.


Fig. 1

Source: Guideline (1999) World Health Organization (WHO) International Society of Hypertension (ISH)

Variations in blood pressure in response to diverse stimuli during the day make it hard to judge on the basis of a single measuring. It is consequently advisable to measure it twice a day (once in the morning on waking and once in the late afternoon). Persistent hypertension should be treated with medications under the strict supervision of a doctor, and certain habits should be corrected (salt consumption, a sedentary way of life, obesity, smoke, alcohol, stress).

The measuring of blood pressure enables you to verify the effectiveness of the measures adopted to keep the situation under control. Monitoring your own blood pressure is a practical and easy-to-repeat process when reliable devices are used and has the added advantage of reducing the feeling of unease often experienced by patients having their blood pressure taken in doctors' surgeries or hospitals, which can result in altered values.

However it is important to underline that the option to measure your own blood pressure should not become a reason to avoid a visit to the doctor, and it certainly should not lead to 'do-it-yourself' treatment.

It is therefore essential that frequency of measuring, the results of the measuring and eventual therapies should be agreed with your doctor. It is always unwise to alter doses, or drugs prescribed by your doctor, on your own initiative.

**READ INSTRUCTION MANUAL CAREFULLY BEFORE USE
THE DEVICE MUST NOT BE DISASSEMBLED
FOR ANY PROBLEM PLEASE CONTACT THE TECHNICAL SERVICE OF DISTRIBUTOR FOR SPECIFIC INFORMATION ON YOUR OWN BLOOD PRESSURE
PLEASE CONTACT A DOCTOR**

3. IMPORTANT SAFETY RULES

1. Before using the device , always check the electrical requirements on manufacturer label and verify the type of battery to use itself.
2. Particular precautions must be made concerning electromagnetic compatibility. The medical device must be set-up and used according to information supplied with the present manual;
3. Respect the safety rules and particularly:
 - Children and/or disable persons may only use the device under the supervision of an adult in possession of his/her full faculties.
 - If the device is stored at low temperature, do not use immediately. Leave the device at room temperature for at least one hour.
 - Avoid extreme environmental conditions and the exposure to direct sunlight.
 - **NEVER DIP ANY PART OF THE DEVICE AND ACCESSORIES IN WATER, IT ISN'T WATER PROOF.**
 - The device must be used only by adult people.
 - Cleaning the device with soft and dry cloth. Do not use petrol, thinners or similar solvents. Spot on the cuff can be removed carefully with a damp cloth. Do not wash any part of the device.
 - If you do not use the device for a long time, please take the batteries off. **KEEP THE BATTERIES OUT OF THE REACH OF CHILDREN.**
 - It is recommended the use of alkaline batteries.
 - Do not wash or iron the cuff
 - Do not squeeze the cuff
 - Keep the cuff away from cutting edges and do not twist it
 - Avoid to drop the device or threat it roughly in any way.
4. The device must not be disassembled. For any problem please contact the Technical Service of the Distributor. **Failure to comply with the above invalidates the safety of the device and the guarantee claims.**
5. **This medical device must be used exclusively for the use for which it has been designed and as described in this manual.** Any different use must be considered incorrect and therefore dangerous; the manufacturer and the distributor cannot be considered liable for damage caused to people, things or animals by improper, incorrect and/or unreasonable use.
6. Measurement on people with serious peripheral circulation or renal problems, people who have had an ictus or if are unconscious and the case of serious arrhythmia may be difficult and bring extra error; it is recommended to consult the physician.
7. In case of malfunctions contact the technical service of the distributor. If by technical check the device will result opened by user, the guarantee will be

- expired.
8. The manufacturer and distributor cannot be considered liable for damages to people, things or animals, caused by, repairs or changes by user or people not authorized.

4. PUTTING THE BLOOD-PRESSURE MONITOR INTO OPERATION

4.1. Inserting the batteries

After you have unpacked your device, first insert the batteries. The battery compartment is located

on the side of the device (see illustration).

a) Remove cover as illustrated in Fig. 2

b) Insert the batteries (4 x size AA 1,5V), **thereby observing the indicated polarity** (Fig.3) and close the cover.

Caution, the life span of the batteries included can be shorter than expected in consequence of storage period.

Fig. 2


Fig. 3

4.2 Battery replacement

1. Remember, once the batteries are removed or used up, the data will be cleared out.
2. A use of new and good quality alkaline batteries, the device can make more than 300 measurements under room temperature of 22°C, inflation pressure 185 mmHg


3. When weak battery symbol appears on the display (Fig. 4), please replace with two new batteries at the same time. Do not leave exhaust batteries into device compartment.

Fig. 4

4. Do not use rechargeable batteries. Use only batteries in compliance with specification declared in this manual.
5. When install or replace the batteries, the device will automatically reset and display all symbols on the LCD for 0.5 seconds (Fig. 5)
6. Remove the batteries from the compartment if the device will be left unused for long periods.

4.3 Display Explanation

The Pulse mark has also a function of data interpretation in order to advise in case of abnormal results. When measuring result is abnormal (Hypertension or Hypotension range) the Pulse mark (Fig. 5) will flash automatically. These symbols are a quick, convenient way of gauging the condition of one's blood pressure, following the table below.

Figure	Range	Systolic (mmHg)	Diastolic (mmHg)
	Hypertension	> 140	> 90
	Normal-High	130~139	85~89
	Normal	< 129	< 84
	Hypotension	< 90	< 50

Data in accordance with:
 (1999) WHO (World Health Organization) / ISH (International Society of Hypertension) Definitions and Classification of Blood Pressure Level
 2003 ESH (European Society of Hypertension) / ESC (European Society of Cardiology) Guidelines

5. HOW TO MEASURE YOUR BLOOD PRESSURE

5.1. Before the measurement:

- Avoid eating, smoking as well as all forms of exercise before the measurement. All these factors influence the measurement result. Try and find time to relax by sitting in an armchair in a quite atmosphere for about ten minutes before the measurement.
- Remove any garment that fits closely to your wrist.
- Measure always on the same arm (normally left).
- Attempt to carry out the measurements regularly at the same time of day, since the blood-pressure changes during the course of the day.

5.2. Common sources of error:

Note: Comparable blood-pressure measurements always require the same conditions! These are normally always quiet conditions.

- All efforts by the patient to support the arm can increase the blood-pressure. Make sure you are in a comfortable, relaxed position and do not activate any of the muscles in the measurement arm during the measurement. Use a cushion for support if necessary.
- If the arm artery lies considerably lower (or higher) than the heart, an erroneously higher (or lower) blood-pressure will be measured!

- A loose cuff or bad position of instrument may cause false measurement values.
- Please avoid strong electromagnetism interference, (such as mobile telephone, microwave oven, etc). This can cause low accuracy in the measurement.
- Use only alkaline batteries.
- Do not judge the measuring results and adjust your treatment yourself, it may be dangerous. Please give measuring results to physician to make a diagnose. Please take your medications as by physician's instructions.
- This equipment may stop working if there is a strong interference of electrostatic discharge.
- If the device is stored at low temperature, do not use it immediately. Leave the device at room temperature for at least one hour before using it. The Arrhythmia may cause mistakes in the measurements, so is good to use with great caution the device on peoples that have arrhythmia symptoms.

5.3. How to place the cuff on the wrist

- a) Place your left arm on a table, with the palm face upward. Put the cuff on the left wrist and fasten the Velcro strap. The circumference of the wrist cuff is about between 13.5 and 21.5cm and you can adjust it to be fitted snugly.
- b) Caution! The distance between the edge of the wrist cuff and the hand must be 1cm max.

- c) Relax your body and don't bend your wrist.
- d) Be aware of the posture and adjust the cuff height. The height of the cuff should be the same as your heart during measurement (see below picture).

Note:

If it is not possible to fit the cuff to the left arm, it can also be placed on the right one. However all measurements should be made using the same arm.

5.4. Measuring procedure

1. Press I/O button to operate the device.
2. The display will switch on with the indication 00 blinking.

3. About 1-2 seconds later, the pressure will be reset, then the cuff begins to inflate and the pressure value will be displayed.
4. When the cuff pressure rises to 185 mmHg, the inflation movements stop and the indicated number drops down to start measuring.

5. If pressure of 185 mmHg is not enough to measure, the device will increase automatically the cuff until a mass pressure of 300 mmHg.
6. If the value of 300 mmHg is not enough to make measure, over that the "Er" sign will appear and the device will stop inflating and release the cuff pressure.

7. Along with the measurements, the heart beat rate is also recorded, at this time, the heart sign will appear on the display.

8. The cuff pressure will fall slowly during measuring. Do not shake the cuff or the monitor. Keep quiet.
9. After finishing measurement the cuff deflates rapidly and automatically and the values of systolic and diastolic pressures and pulse rate will be displayed on the display together with the heart symbol that must be interpreted as indicated at paragraph 4.4 of this instruction manual.

10. The measuring results will be shown until the I/O button is pushed again or till automatic switch off 2 minutes later.
11. When the results was shown, press “M1” or “M2” key to record data in the Memory Set.

5.5. FAST INTERRUPTION DURING A MEASUREMENTS

If it is necessary to interrupt a blood pressure measurement for any reason (e.g. the patient feels unwell), the I/O button can simply be pressed at any time. The device immediately reduces the cuff pressure.

5.6.Memory – displaying the last measurement

The device can store twin sets (total of 100sets) of recorded data, and each set data includes values of diastolic and systolic pressure and pulse rates. When “M1” or “M2” button is pressed, “01” will be displayed for 1 second and then the first stored data will be displayed on the display. If you want to read the second stored data, press “M1” or “M2” again and the second stored data will be displayed. With the same procedure you can read, backwards, all other stored data. When press “M1”, the display shows number 1. When press “M2”, the display shows number 2.

Further information

Measurements should not occur soon after each other, since otherwise, the results will be falsified. Wait therefore for several minutes (at least 15 min.) in a relaxed position, sitting or lying, before repeating a measurement.





6. CARE, STORAGE AND REPAIR

1. It's necessary to protect this device against high moisture and heat, direct sunlight, shock, solvents, alcohol.
2. Remove the batteries if the device will be unused for a long time, and keep the batteries out of the reach of children.
3. Keep the device away from sharp subjects and don't extend or twist the cuff.
4. Use only soft and dry cloth to clean the device.

WARNING: Under no circumstances any part of device and its accessories could be washed.

Keep the device and all of its accessories out of the reach of children.

7. ERROR MESSAGE

INDICATION	PROBLEM	CAUSE AND SOLUTION
	Inflation fault	<ol style="list-style-type: none"> 1. Cuff is not fastened well, please refasten it. 2. Cuff is broken, change it with new one. 3. Pump, venting or deflation valve fail, please refer to the distributor.
	Artificial interference and noise Abnormal reading	<ol style="list-style-type: none"> 1. Move or talk during measurement. 2. Posture, etc. was different for each measurement. 3. User is too nervous. 4. User has severe heart problem.
	Deflation too fast	<ol style="list-style-type: none"> 1. Deflation valve failure, refer to the distributor. 2. Cuff not fastened, re-fasten it.
	Over safety pressure	<ol style="list-style-type: none"> 1. Air pipe blocked. 2. Inflation controller failure. 3. Pressure sensor fault. Please check and repeat the measurement and/or ask to your distributor.
	Low Battery	Replace with two new batteries at the same time.

8. TROUBLESHOOTING

TROUBLE	CHECK POINT	SOLUTION
With the batteries installed, there is no indication on the display when the I/O (On/Off) button is pressed.	<ol style="list-style-type: none"> 1. The batteries are exhausted. 2. The polarity of batteries is wrong. 3. The contacts of batteries compartment is dirty. 	<ol style="list-style-type: none"> 1. Replace the batteries. 2. Install the batteries correctly. 3. Clean the battery terminals with dry cloth.
The inflation action cannot be done or the air pressure cannot rise.	<ol style="list-style-type: none"> 1. The cuff position is not correct. 	<ol style="list-style-type: none"> 1. Check the cuff position and repeat the measurement.

TROUBLE	CHECK POINT	SOLUTION
The weak battery sign appears on the display	Battery low	Replace the batteries
The blood pressure cannot be taken and on the display appears error message.	<ol style="list-style-type: none"> 1. The cuff is not well fastened. 2. User is nervous. 3. The posture is not correct. 	<ol style="list-style-type: none"> 1. Re-fasten the cuff. 2. Sit down and relax. Keep quiet during the measurement. 3. Keep the cuff and heart at the same height during the measurement. <p>If the patient has severe heart problem, then the blood pressure may not be read correctly.</p>
Under normal measurement circumstances, the reading at home is different from that of the clinics, or each measurement has various reading.	<ol style="list-style-type: none"> 1. The variation is due to the different environments. 2. The blood pressure changes with physiological or psychological status. 	Record everyday value trend and consult it with a physician.

9. TECHNICAL CHARACTERISTICS

Model	REALCHECK MICRO 2
Size	67(W) X 73(H) X 70(D)mm
Weight	Approximately 110g not including batteries
Measuring method	Oscillometry
Measuring range	20 to 280 mmHg (pressure)
	40 to 200 beats/minute (pulse rate)
Measuring accuracy	± 3 mmHg for systolic and diastolic pressure
	± 5% of the reading for the pulse rate
Inflation	Automatic
Rapid deflation	Automatic valve
Pressure release device	Automatic valve
Operating out	Automatic (after 2 minutes)
Batteries	2x1.5V LR03 alkaline batteries type AAA
Memory	Twin zones (50x2 memories)
Operation temperature and humidity	+10°C to +40°C, 30% to 85%RH
Storage temperature and humidity	-5°C to +50°C, 15% to 90%RH
Cuff size	Applicable for wrist with size between 13.5 to 21.5 cm

Notes: Specification are subject to change without further notice.

Guidance and manufacture's declaration – Electromagnetic immunity (as request by regulation EN 60601-1-2:2001)			
The model Automatic Blood Pressure Monitor is intended for use in the electromagnetic environment specified below. The customer or the user of the model Automatic Blood Pressure Monitor should assure that is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floor should be wood, concrete or ceramic tile. If floors are covered with synthetic material the relative humidity should be at least 30%.
Electrical fast transient/ burst IEC 61000-4-4	Not applicable	-	-
Surge IEC/EN 61000-4-5	Not applicable	-	-
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	Not applicable	-	-
Power frequency magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment

Guidance and manufacture's declaration – Electromagnetic immunity (as request by regulation EN 60601-1-2:2001)			
The model REALCHECK MICRO 2 Blood Pressure Monitor is intended for use in the electromagnetic environment specified below. The customer or the user of this device should assure that is used in such an environment.			
Immunity test	IEC / EN 60601 test level	Compliance level	Electromagnetic environment-guidance
Conducted RF IEC 61000-4-6	Not applicable	-	Portable and mobile RF communications equipment should be used no closer to any part of the REALCHECK MICRO 2 Blood Pressure Monitor, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
Radiated RF IEC 61000-4-3	3V/m 80MHz-2.5GHz	3V/m	

			$d = \left[\frac{3.5}{E_f} \right] \sqrt{P}$ 80MHz to 800MHz $d = \left[\frac{7}{E_f} \right] \sqrt{P}$ 800MHz to 2,5GHz Where P is the maximum out put power rating of the transmitter in watts (W) according to the transmitter manufacturer and is recommended separation distance in meters (m). Field strength from fixed RF transmitters, as determined by an electromagnetic site survey (a), should be less than the compliance level in each frequency range (b). Interference may occur in the vicinity of equipment marked with the following symbol: 
<p>Note 1 At 80MHz and 800MHz the higher frequency range applies.</p> <p>Note 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>			
	<p>(a) Field strength from fixed transmitters, such as base stations for radio telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the REALCHECK MICRO 2 Blood Pressure Monitor should be observed to verify normal operation. If abnormal performance is observed additional measures may be necessary, such as re-orienting or relocating the device.</p> <p>(b) Over the frequency range 150KHz to 80MHz, field strength should be less than 3V/m</p>		

10. SIMBOLOGY

	Type BF
CE 0123	CE marking in conformity with EC directive 93/42/EEC
	Warning, consult the instruction manual
	To keep in dry and cool place
	Storage temperature
	Direct Current
I	Start
O	Off

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“Actuation of European directives 2002/95/EC, 2002/96/EC and 2003/108/EC, for reduction in use of dangerous substances in the electric and electronic device and for garbage disposal”

The symbol as over applied on the device or its packaging means that at the end of its useful life the product must not be disposed of with domestic waste. At the end of device useful, the user will must deliver it to the able collecting centres for electric and electronic garbage, or give back to the retailer in the moment of equivalent new device purchasing, one against one.

Disposing of the product separately prevents possible negative consequences for the environment and for health, deriving from inadequate disposal. It also allows the recovery of materials of which it's made up in order to obtain an important saving of energy and resources and to avoid negative effects to the ambient and health.

In case of abusive disposal of device by user, will be applied administrative endorsements in compliance with current standard.

11. WARRANTY

Warranty conditions

- The unit is covered by a 2-year warranty against defects in materials and workmanship. The warranty does not extend to units that have undergone unauthorised repairs or been tampered with.
- The warranty period commences on the purchase date. It is mandatory to attach the sales receipt or invoice serving as proof of the purchase date to the warranty certificate below, which must be filled out in every part.
- Under the warranty, any components with manufacturing defects will be replaced and/or repaired free of charge. The Purchaser will bear all expenses for transporting, delivering and collecting the unit.

The warranty does not cover:

- The accessories supplied with the unit or parts subject to normal wear.
- Repair of unsubstantiated defects.
- Repairs of units that have failed or been rendered defective as a result of improper use, carelessness or negligence or if the damage is not attributable to the manufacturer (accidental falls, careless transport, etc.).

The warranty does not cover any direct or indirect damage or injury caused to persons, animals or property as a result of improper use of the product or during its period of inefficiency.

WHEN SUBMITTING A WARRANTY CLAIM, IN OBSERVANCE OF THE ABOVE CONDITIONS, SEND THE UNIT TOGETHER WITH THE CERTIFICATE BELOW AND PROOF OF PURCHASE TO:
SSL HEALTHCARE ITALIA S.P.A. VIA MARCO EMILIO LEPIDO, 178/5 - BOLOGNA (ITALY).

CERTIFICATE TO BE RETURNED IN CASE OF REPLACEMENT/REPAIRS

THE WARRANTY IS VALID ONLY IF THE SALES RECEIPT OR INVOICE IS ATTACHED

Unit Type/Model: _____

Lot: _____ Serial No: (SN) _____

Date of purchase: _____

PURCHASER'S DETAILS

Last Name and Name: _____

Address: _____

Telephone: _____

Description of the defect: _____

I hereby authorise the use of the above information in accordance with the Law on Privacy 196/03

Signature confirming acceptance of warranty conditions _____

CE 0123

Medical Device
Classe IIa
(DIR 93/42/CEE)

realcheck™

Fabbricante/Manufacturer:

CA.MI snc

Via Ugo La Malfa, 31
43010 Pilastro (PR) - Italy

Distributore/Distributor

SSL Healthcare Italia SPA

Via M.E.Lepido, 178/5
40132 Bologna - Italy