Peri-Strips Dry®

With Veritas® Collagen Matrix Circular Staple Line Reinforcement
Avec renforcement de ligne d'agrafes circulaire de la matrice de collagène Veritas®
Mit Veritas® Collagen-Matrix kreisförmiger Klemmnahtverstärkung
Con rinforzo circolare della linea di sutura Veritas® a matrice di collagene
con refuerzo de la línea de grapado circular con matriz de colágeno Veritas®
Met Veritas® Collageenmatrix Ronde nieterlijnversterking
Med Veritas® kollagen-matrix • cirkulær hæftesømsforstærkning
Med Veritas® kollagenmatris med cirkulär agraflinje förstärkning
Med Veritas® kollagenmatrise kantforsterkning med sirkulære stifter
Veritas®'li, Kollajen Matriks Dairesel Zimba Hattı Kuvvetlendirici
Με υλικό ενίσχυσης γραμμής κυκλικών συνδετήρων με μήτρα κολλαγόνου Veritas®

INSTRUCTIONS FOR USE.....	4
MODE D'EMPLOI.....	9
GEBRAUCHSANWEISUNG	14
ISTRUZIONI PER L'USO.....	19
INSTRUCCIONES DE USO	24
GEBRUIKSAANWIJZING	29
BRUGSANVISNING	34
BRUKSANVISNING	39
BRUKSANVISNING	44
KULLANMA TALÍMATI	49
ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ.....	54

Peri-Strips Dry® with Veritas® Collagen Matrix Components

Illustrations 1-4


Illustration 1


Illustration 2


Illustration 3


Illustration 4

Peri-Strips Dry® with Veritas® Collagen Matrix

Figures 1-7


FIGURE 1


FIGURE 2


FIGURE 3


FIGURE 4


FIGURE 5


FIGURE 6

FIGURE 7

A

PSD Gel
Gel PSD
PSD-Gel
Gel PSD
Gel PSD
PSD Gel
PSD-gel
PSD-gel
PSD Gel
PSD Jel
Γέλη PSD

B

Anvil Loading Unit
Unité de charge de l'enclume
Amboss-Ladevorrichtung
Unità di caricamento incudine
Unidad de carga del yunque
Aambeeld-ladeenheid
Ambolt-klargøringsenhet
Städladdningsenhet
Ambolt-ladeenhet
Örs Yükleme Ünitesi
Μονάδα φόρτωσης άκμωνα

C

Cartridge cone
Cône de cartouche
Kassettenkappe
Cono cartuccia
Cono del cartucho

Cassetteconus
Patronkonus
Kassettkon
Patronkjegle
Kartuş konisi
Κωνική φύστιγγα

D

Cartridge Loading Unit
Unité de charge de la cartouche
Kassetten-Ladevorrichtung
Unità di caricamento cartuccia
Unidad de carga del cartucho

E

Cassette-laadeenheid
Patronklargøringsenhet
Kasset-laddningsenhet
Patron-ladeenhet
Kartuş Yükleme Ünitesi
Μονάδα φόρτωσης φύστιγγας

F

Buttress
Renfort
Verstärkung
Materiale di rinforzo
Refuerzo
Steun
Forstærkning
Stöd
Forstøtning

G

Inner Plastic Ring
Anneau en plastique interne
Innerer Kunststoffring
Anello interno di plastica
Anillo de plástico interior
Plastic binnenring
Indvendig plastikring
Inre plastering
Indre plastering
İç Plastik Halka
Εσωτερικός πλαστικός δακτύλιος

H

Destek
Υλικό ενίσχυσης
Attachment tool
Outil de fixation
Befestigungsteil
Strumento di fissaggio
Herramienta de acoplamiento
Hulpgereedschap
Fastgørelsesværktøj
Fästverktyg
Festeverktøy
Ataşman aracı
Εργαλείο προσάρτησης
Innerer Kunststoffring
Anneau en plastique interne
Innerer Kunststoffring
Anello interno di plastica
Anillo de plástico interior
Plastic binnenring
Indvendig plastikring
Inre plastering
Indre plastering
İç Plastik Halka
Εσωτερικός πλαστικός δακτύλιος

I

C-Clip
Fermo a C
Pasador C
C-klem
C-klemme
C-klämna
C-klemme
C klips
Συνδετήρας "C"

J

Notches
Encoches
Kerben
Tacche
Muescas
Nokken
Riller
Skårer
Hakk
Çentikler
Έγκοπές
Inner Plastic Ring
Anneau en plastique interne
Innerer Kunststoffring
Anello interno di plastica
Anillo de plástico interior
Plastic binnenring
Indvendig plastikring
Inre plastering
Indre plastering
İç Plastik Halka
Εσωτερικός πλαστικός δακτύλιος

SYMBOL DEFINITIONS:

Store at a controlled room temperature (20°-25°C/68° -77° F). Do not use if exposed to heat. Check heat indicator.

Do not reuse

Consult *Instructions for Use*

STERILE

Sterilized

STERILE EO

Ethylene Oxide Sterilized

Steam Heat Sterilized

NaOH

Product treated with sodium hydroxide.

BOVINE

Product derived from USDA inspected cattle.

MADE IN THE USA

Rx Only

CAUTION: Federal (USA) law restricts this device to sale by, or on the order of, a physician.

Peel

REF

Catalog number

Use by Date

LOT

Lot number

Compatible stapler models

C

Cartridge

A

Anvil

Manufacturer

PN SSI part number. TN SSI tracking number. C SSI internal code. D SSI internal code.

DESCRIPTION:

Peri-Strips Dry® with Veritas® Collagen Matrix (PSD-V) Circular Staple Line Reinforcement is prepared from dehydrated bovine pericardium procured from cattle under 30 months of age in the United States.

The packaging contains two loading units w/buttress material, one for the anvil and one for the cartridge side of the stapler. The packaging also includes one cartridge cone to aid atraumatic advancement of the stapler into the surgical site, and one tube of PSD Gel with Tip to create a temporary bond between the buttress and the surgical stapler until the stapler is positioned and fired.

All components are packaged sterile (tray within a tray).

INDICATIONS FOR USE:

PSD-V is intended for use as a prosthesis for the surgical repair of soft tissue deficiencies using surgical staplers when staple line reinforcement is needed.

PSD-V can be used for reinforcement of staple lines during gastric, bariatric, small bowel, colon and colo-rectal procedures.

CONTRAINDICATIONS:

The use of PSD-V is contraindicated in patients with known sensitivity to bovine material.

ADVERSE REACTIONS:

As with any surgical procedure, adverse reactions are possible and include but are not limited to: infection, rejection, erosion, and allergic reaction.

WARNINGS:

Do not re-sterilize.

Discard all open unused components.

Do not use product if there is damage to the tray seals.

Ensure the staple line is completely covered with the buttress or inadequate coverage after firing may result.

PSD-V is not designed, sold, or intended for use except as indicated; doing so may result in surgical complications.

Synovis products differ; substitution of one product for another product may reduce product performance.

Do not use PSD-V if the product has been exposed (1) to solutions above room temperature (2) to chemicals antibiotics or other substances other than specifically addressed in the instructions.

CAUTIONS:

Do not get the anvil or cartridge buttress wet before applying PSD Gel, or the buttress may not adhere to the stapler properly.

Ensure the anvil buttress and cartridge buttress are on the corresponding stapler components or the buttress may not adhere to the stapler properly.

Cartridge cones may be used for entry into the surgical site on the cartridge side of the stapler during endoscopic surgery. Before completion of surgery be sure to remove the cartridge cone from the surgical site.

Final tissue compression, including Peri-Strips Dry with Veritas (PSD-V), must meet the range specified by the stapler manufacturer. PSD-V increases the total thickness of the area stapled by 0.4 mm - 1.2 mm (0.016" - 0.048").

Follow *Instructions for Use* supplied by the stapler manufacturer. Do not use PSD-V contrary to the stapler manufacturer's instructions.

The cartridge and anvil loading units differ; substitution of one loading unit for another may interfere with closing and firing of the stapler.

PERI-STRIPS DRY WITH VERITAS COMPONENTS:

Each PSD-V package contains the following components (See A-J):

Cartridge and anvil loading units (See Illus. 1): The cartridge loading unit includes a circular PSD-V buttress affixed on an inner plastic ring mounted onto attachment tool, the cartridge assembly for the 21 mm Autosuture™ stapler does not contain an inner plastic ring.

The anvil loading unit includes a circular buttress affixed on an inner plastic ring loaded onto attachment tool; however, the anvil assemblies for the 21 mm staplers do not contain an inner plastic ring.

The buttress is mounted on a plastic attachment tool designed to properly align and deliver the

buttress onto the stapler cartridge or anvil.

Cartridge cone (See Illus. 1): The cartridge cone is designed to protect the cartridge buttress during advancement of the stapler into the surgical site. There are two features to facilitate removal of the cone from the stapler and from the surgical site before completion of the surgery. The indentations at the top of the cone are designed to allow the use of a forceps/hemostat to remove the cone from the stapler and/or surgical site. The holes at the tip of the cartridge cone are designed to allow a suture to be placed in the cone to facilitate the removal of the cone from the surgical site.

PSD Gel: PSD Gel is used to create a temporary bond between the buttress and the stapler until the stapler is positioned and fired. The plastic tip is intended to help control the amount of gel applied.

C-clip (See Illus. 3): The anvil assemblies for the 21 mm staplers contain a C-clip attached to the narrow end of the attachment tool. The C-clip is designed to secure the buttress to the anvil.

Individually packaged Attachment Tool (Anvil only REF Attachment Tool 21E or REF Attachment Tool 21U): The individually packaged Attachment Tool is for use with the anvil assembly only and is for use in the event the attachment tool supplied with the anvil assembly (see Figure 1) malfunctions during the process of applying the buttress to the anvil side of the stapler. The individually packaged Attachment Tool is supplied for use with 21mm staplers only.

INSTRUCTIONS FOR USE:

Each model of PSD-V has been designed specifically for the stapler models indicated on the label; verify that the correct model of PSD-V has been selected.

Loading technique for PSD-V varies. Please follow appropriate technique as indicated below.

After firing the stapler, plastic alignment components (inner plastic ring and C-clip if applicable) will remain in the cut-out within the stapler head.

A. Instructions for Circular Staplers (Ethicon Endo-Surgery, Inc.):

- A1. Remove the Tyvek lid from the outer tray and aseptically remove the inner tray; the inner tray may be placed in the sterile field. Inspect the inner tray. Do not use if the inner tray is damaged or if the seals are not intact.
- A2. Remove the Tyvek lid and plastic tray insert from the inner tray.

- A3. Break off the metal tip of the PSD Gel tube.
- A4. Attach the plastic gel tip to the top of the PSD Gel tube.
- A5. Remove the cartridge cone from the tray. The indentation in the cone contains several small holes through which a suture might be placed. Loop and tie a suture at least the length of the stapler (approximately 15 inches) through two of the holes in the cartridge cone (See Figure 7).

Note: *The recommended suture is an absorbable or non-absorbable material containing no memory for shape, e.g. silk, attached to a 2-0 curved needle.*

- A6. Using the identification on the plastic storage tray identify the anvil (A) and cartridge (C) loading units.

Caution: The cartridge and anvil loading units differ; substitution of one unit for another may interfere with closing and firing of the stapler.

Preparing the anvil side:

- A7. Ensure the anvil side of the stapler is clean and dry.
- A8. Apply a continuous thin line of PSD Gel onto the anvil side of the circular stapler (See Figure 1).

Caution: Do not get the buttress or stapler components wet before applying PSD Gel, or the buttress may not adhere to the stapler properly.

- A9. Place the anvil loading unit (buttress side down) over the shaft of the anvil (See Figure 2).

For Ethicon staplers 25 mm or >: Align the notches on the inner plastic ring with the locking springs on the anvil shaft. If necessary rotate the loading unit until the inner plastic ring fits into the anvil (See Illus. 4).

- A10. Remove the attachment tool by pulling it away from the stapler (See Figure 3).
- A11. Using the flat side of the attachment tool or gloved fingers, secure the anvil buttress to the anvil. Ensure there are no bubbles or gaps between the buttress and anvil and that the buttress completely covers the staple line.

For Ethicon 21 mm staplers:

- A12. The narrow end of the attachment tool contains a C-clip to secure the anvil buttress to the anvil. Using the attachment tool, slide the C-clip onto the shaft of the anvil. (See Figure 4).

- A13. Remove attachment tool leaving the C-clip behind.

Note: *In the event the c-clip becomes dislodged from the shaft of the anvil during steps A12 or A13, discard the c-clip and open the individually packaged Attachment Tool. Proceed with preparation of the anvil side as instructed beginning with step A12.*

- A14. Using gloved fingers, slide the C-clip down into the head of the anvil towards the purse string notch (See Figure 5).
- A15. Allow the buttress to adhere for a minimum of 5 minutes before using.

Preparing the cartridge side:

- A16. Ensure the cartridge side of the stapler is clean and dry.
- A17. Apply a continuous thin line of PSD Gel onto the cartridge side of the circular stapler.
- Caution:** Do not get the buttress or stapler components wet before applying PSD Gel, or the buttress may not adhere to the stapler properly.
- A18. Place the cartridge loading unit (buttress side down) over the inner shaft of the cartridge head. If necessary rotate the loading unit until the inner plastic ring fits into the cartridge head (See Figure 6).
- A19. Remove the attachment tool by pulling it away from the stapler (See Figure 3).
- A20. Using the flat side of the attachment tool or gloved fingers, secure the cartridge buttress to the stapler. Ensure there are no bubbles or gaps between the buttress and stapler and that the buttress completely covers the staples.
- A21. Allow the buttress to adhere for a minimum of 5 minutes before using.
- A22. Lower the inner shaft of the stapler and firmly secure the cartridge cone over the end of the stapler (See Figure 7).

B. Instructions for using Peri-Strips Dry with Veritas Collagen Matrix (PSD-V) 21mm and 25mm Circular Staple Line Reinforcement on Autosuture™ DST Series EEA 21mm and 25mm Staplers

Select current Ethicon product PSD 21-E-V or PSD 25-E-V for corresponding size of Autosuture™

DST Series EEA stapler. Verify that the correct size of Peri-Strips Dry with Veritas has been selected.

B1. Remove the Tyvek lid from the outer tray and aseptically remove the inner tray; the inner tray may be placed in the sterile field. Inspect the inner tray. Do not use if the inner tray is damaged or if the seals are not intact.

B2. Remove the Tyvek lid and plastic tray insert from the inner tray.

B3. Attached to the PSD Gel tube is a plastic gel tip, snap off the top of the PSD Gel tube

B4. Attach the gel tip to the top of the plastic Gel tube (this gel tip helps control the amount of gel that is placed on the stapler).

Using the identification on the plastic storage tray identify the anvil (A) and cartridge (C) loading units.

The cartridge assembly includes a buttress affixed on an inner plastic ring. The cartridge assembly for Autosuture™ 21mm does not include an inner plastic ring. The anvil assembly for the 25mm staplers also includes a buttress affixed on an inner plastic ring (See Illus. 1).

B5. Remove the buttress material from the inner plastic ring and the attachment tool.

Preparing the anvil side (21 or 25 mm stapler):

Ensure the anvil side of the stapler is clean and dry.

B6. Apply a thin continuous line of PSD gel onto the anvil side of the circular stapler (See Figure 1).

B7. Place the anvil buttress over the shaft of the anvil.

B8. Using gloved fingers secure the anvil buttress to the anvil. Ensure there are no bubbles or gaps between the buttress and the anvil and that buttress completely cover the staple line. Using gloved fingers squeeze the excess gel out from underneath the buttress, smooth down the top of the buttress. Wipe off any remaining gel with finger.

B9. Allow the buttress to adhere for a minimum of 1 hour before using.

Preparing the OrVil device:

B10. Transect one cut in tissue ring.

B11. Apply a thin continuous line of PSD gel onto the anvil side of the circular stapler.

B12. Wrap anvil buttress around shaft from left to right, adjoining the two ends on the other side.

B13. Using gloved fingers secure the anvil buttress to the anvil. Ensure there are no bubbles or gaps between the buttress and the anvil and that buttress completely cover the staple line. Using gloved fingers squeeze the excess gel out from underneath the buttress, smooth down the top of the buttress. Wipe off any remaining gel with finger.

B14. Allow the buttress to adhere for a minimum of 1 hour before using.

Preparing the cartridge side (21 or 25mm stapler):

Ensure the cartridge side of the stapler is clean and dry.

B15. Apply a thin continuous line of PSD gel onto the cartridge side of the circular stapler (See Figure 1).

B16. Place the cartridge buttress over the shaft of the cartridge.

B17. Using gloved fingers secure the cartridge buttress to the cartridge. Ensure there are no bubbles or gaps between the buttress and the cartridge and that buttress completely cover the staple line. Using gloved fingers squeeze the excess gel out from underneath the buttress, smooth down the top of the buttress. Wipe off any remaining gel with finger.

B18. Allow the buttress to adhere for a minimum of 1 hour before using.

Implanting Peri-Strips Dry with Veritas:

After one hour, check the buttress to see if it is dry, if the tissue is still tacky or wet, too much gel was used and more time is needed for gel to dry.

B19. For optimal performance, limit the anvil dwell time between introduction into the surgical site and stapler firing to no more than 10 minutes. Longer dwell time may cause delamination.

DISCLAIMER OF WARRANTIES:

Synovis Surgical Innovations (SSI), a division of Synovis Life Technologies, Inc., warrants that reasonable care has been used in the manufacture of this device. This warranty is exclusive and in lieu of all other warranties whether expressed, implied, written, or oral, including, but not limited to, any implied warranties of merchantability or fitness. As a result of biological differences in individuals, no product is 100% effective under all circumstances. Because of this fact and since SSI has no control over the conditions under which the device is used, diagnosis of the patient, methods of administration or its handling after it leaves its possession, SSI does not warrant either a good effect or against an ill effect following its use. The manufacturer shall not be liable for any incidental or consequential loss, damage or expense arising directly or indirectly from the use of this device. SSI will replace any device, which is defective at the time of shipment. No representative of SSI may change any of the foregoing or assume any additional liability or responsibility in connection with this device.

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Peri-Strips Dry and Veritas are registered trademarks of Synovis Life Technologies, Inc., (US Patent #6,312,474; 6,652,594; 7,128,748; Additional Patents Pending)
AR-MED Limited (Company Number 3154742)
A Division of Synovis Life Technologies, Inc.