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*All the information given in this manual is what available at the time of printing.
The manufacturer reserves the right to make changes to the product at any time without notice.
It is advisable to verify if there are any updates.*

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All rights reserved.

*Forbidden **any form** of print, duplication, reproduction and publication of this instruction manual, or
part of it, without the written agreement of GIBERTINI ELETTRONICA Srl*

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1 UNPACKING

If the unit has been shipped by carrier, carefully unpack the cartons and stand the distillation unit on a bench. Identify the box containing the reagents, bottles and various connectors.

A visual check should be made of all items for evidence of any damage that, if found, must be reported to the carrier and the suppliers within 3 days.

1.1 STANDARD ACCESSORY KIT

Glassware

3	100 Volumetric flasks with plug
1	200 Volumetric flasks with plug
2	300 ml Conical flasks with plug
1	500 ml Dispenser
1	100 ml Dropper bottle
3	Dropper plastic pipettes

Reagents

1	500 ml Calcium oxide CaO suspension 12%
1	250 ml Tartaric acid 50 %
1	100 ml Concentrated silicone antifoam
1	500 ml Solution for liqueurs
1	500 ml solution to clean glass nozzle (solution B)

Fittings

1	Water inlet tubing and clips
1	Water drainage tubing and clip
1	Tube for phial
1	Tube for peristaltic pump
1	Tube for steam

Spare parts

1	Cartridge for water (column of deionisation)
---	--

NOTE

The flasks are calibrated and interchangeable in mass, not in volume.
Remember that our flasks explicitly do not indicate ml, but only 100 and 200.
We DO NOT withdraw 100 ml, but a repeated quantity. So accuracy in ml is NOT necessary, but it is imperative that the volume of wine and distillate is the same at the same temperature.

1.2 ELECTRICAL REQUIREMENTS

Power surge:

220V (-10/+15%) - 50 Hz

Power consumption:

3300 W max

2600 W normal


Figure 1

1.3 WATER REQUIREMENTS

WATER

Minimum pressure : 2 bar

Max temperature : 20 °C

Flow rate : 5 litres/min.

Warning – Verify the cleanness in relation to the water “hardness”.


Figure 2

The water of the cooling column is used for the cleaning of the drainage basin. If an other cleaning is necessary , select “Alcohol”. Put the volumetric flask and start the distillation with empty distillation chamber

1.4 DIMENSIONS AND WEIGHT

Dimensions: L30 x D48 x H102 cm

Net weight: 41 kg

2 SPECIFICATIONS

2.1 GENERAL DESCRIPTION

An automatic distillation unit designed for distillation of wines, spirits and beers for determination of alcohol content and volatile acidity.

The measured quantity of sample, (100 or 200 ml for alcohol content, 20 ml for volatile acidity), is placed in the distillation chamber, together with the reagents. Reagents are added to neutralise the sample and to make a conductive solution.

The distillate is collected in a volumetric flask mounted on a precision balance arm. The balance will stop the distillation when a preset volume has been collected.

For volatile acidity, the distillate is collected in a 300 ml conical flask.

The SUPER DEE can be used for determination of other analytical parameters of the wines (e.g. sorbic acid, total sulphur dioxide). For further informations see the book "Analysis methods of wines and alcoholic beverages" Gibertini – current Edition.

The SUPER DEE has an ecologic device for the waste water: the water passes through the cooling coil only during the distillation, and it cleans the drainage basin.

The SUPER DEE complies with the directives 2006/42 EEC and 85/375/EEC for the safety of the analysts. The instruments are also in accordance with the directive 2004/108 EEC (electromagnetic compatibility).

2.2 EXECUTION TIMES

Wines	100 ml	3 - 4 minutes
Wines	200 ml	6 - 7 minutes
Volatile acidity	250 ml	6 - 7 minutes
Spirits		4 - 8 minutes

2.3 COMPARISON WITH THE OFFICIAL METHODS

The instruments are in conformity with OIV methods as prescribed by Regulation n. 606/2009 EEC, both for the Alcoholic Strength by Volume (ASV) and for volatile acidity.

2.3.1 Alcoholic strength by volume (ASV)

For the wines, the current regulations permit the use of any model of apparatus for distillation or extraction with a steam current on condition that it satisfies the following test:

- Distill 5 times in succession a 10 % vol water-alcohol mixture. After the last distillation, the distillate should present an Alcoholic Strength by Volume of almost 9,9 % vol. In other words, during each distillation, the instrument does not lose a quantity of alcohol higher than 0,02 % vol.

For spirits, the following test must be satisfied:

- Distil 200 ml of water-alcohol solution, with well-known titre near 50 % vol. The instrument does not lose a quantity of alcohol higher than 0,1 % vol.

2.3.2 Volatile acidity

The ECC legislation does not describe a unique type of distillation apparatus, but it states that any apparatus or technique, which satisfies the following 3 requirements, constitutes an official international apparatus or technique:

- 1) The steam of water produced by the generator must be free from CO₂. Distill 20 ml water collecting 250 ml of distillate, adding 0,1 ml di NaOH 0,1 N and 2 drops phenolphthalein 1 % in alcohol: the colour must remain stable pink for at least 10 seconds.
- 2) 99,5 % of acetic acid contained in a 0,1 M solution must be found in the distillate. Distil 20 ml acetic acid 0,1 M. Collect 250 ml distillate and titrate with at least 19,9 ml NaOH 0,1 N.
- 3) A maximum of 0,5 % of a 1 M solution of lactic acid may distil. Distil 20 ml of a M solution of lactic acid, collect 250 ml of distillate and tritrate with maximum 1,0 ml NaOH 0,1 N.

3. OPERATIVE MENU

3.1 LANGUAGE SELECTION MENU

To select a language different than the one currently set, operate as follows:

- Switch off the distillator with the ON/OFF SWITCH.
- Switch on again pressing simultaneously the **SHIFT** key on the keyboard and keep pressed until the end of the initial test.
- The display shows the current language and the shifting arrows in correspondence of **A** and **C** keys. Select the language with these two keys.
- Confirm with the **MODE** key.

3.2 OPERATIVE MENU

See also Chapter 6 (ANALYSIS – Example of distillation)

The instrument display is composed by two lines of text. The first shows the operative mode and the second describes the function corresponding to **A**, **B** and **C** keys.

Mode: Alcohol Alcohol Volat.

Switch On the instrument. The display shows the menu with the analysis type selection (alcohol or volatile)

By selecting the Alcohol Mode, the instrument it is ready for a distillation with electrodes. By selecting Volatile Mode, the instrument is ready for a distillation with steam and it starts to warm the boilers.

For a correct distillation with steam, it is recommended to predispose the instrument in volatile mode some minutes before the distillation and if needed to perform a distillation of water only.

By pressing the **MODE** key, it is possible to proceed to the next menu,

Set : 80 80 160 240
--

The software has 3 definite set points : **80**, **160** and **240** grams, selectable by pressing the **A**, **B** or **C** key respectively

The precision balance shows the gram.

Set : 87 - +
--

It is possible to modify these values by pressing the **SHIFT** key with the **A** (-) or **B** (+) key

By pressing the **MODE** key, it is possible to proceed to the distillation start menu,

Standby Start

It is possible to start the distillation by pressing the **A** key, correspondent to the **START** function. The distillation does not start without a volumetric flask or a conical flask on the balance.

Once the distillation has started, the display shows the weight value collected and the total volume to be collected.

15 -> 80 Stop
--

When the distillate reaches the selected quantity (80), the instrument does not accept a new start if the flask is not removed.

It is possible, at any moment, to stop the distillation by pressing the **B** key (**Stop**).

If the distillation is running, only the **STOP** key (**B**) is active. If the instrument is not active, it is possible to go to the following menu by the **MODE** key, to select the quantity of distillate to collect.

4 DETERMINATION OF VOLATILE ACIDITY

Pour into the distillation chamber 20 ml wine, free of CO₂. It is advisable to add 1 ml tartaric acid 50 % (for example if the wine has a fixed acidity less than 6 g/l).

Lower the plexiglas shield, close the distillation chamber and insert a conical flask on the balance, otherwise the distillation DOES NOT START (see figure 3), and press START push button.

The distillation stops automatically as soon as about 250 ml of distilled liquid have been collected. The time of distillation is 6-7 minutes.

At the end of the distillation, add to the distillate a few drops of phenolphthaleine 1 % and titrate with NaOH 0,1 N until the distillate appears of a stable pink colour for at least 10 seconds.

Be **n** the millilitres used.

$$\text{Gross volatile acidity} = n \times 0,3$$

In order to obtain the correct value of volatile acidity, it is necessary to determine the amount of sulphur dioxide present in the distillate, whose quantity must be subtracted from the gross volatile acidity.

Figure 3

4.1 SUBTRACTION OF SULPHUR DIOXIDE SO₂

Add to the already titrated distillate a few drops of sulphuric acid solution diluted 1:4 and about 5 ml of starch paste 1%. The solution will become colourless.

Titrate the free sulphur dioxide with I₂ N/100 until the solution turns blue. Be **n'** the millilitres used.

Now make the solution alkaline by adding 20 ml of saturated sodium borate solution.

The solution will turn pink .

Titrate with I₂ N/100 until the appearance of the blue colour. **n''** is the millilitres used.

Subtract from ml of NaOH 0,1 N (**n**) a quantity equivalent to: $\frac{n'}{10} + \frac{n''}{20} = n'''$

Correcy volatile acidity = (**n** – **n'''**) x 0,3 (g/l acetic acid).

The automatic enological analyser GIBERTINI QUICK allows the titration of the distillate for the determination of volatile acidity, either gross or correct.

5 WINES AND SPIRITS DISTILLATION

5.1 THE DISTILLATION PROCEDURE FOR WINES

May the young fermenting wines or sparkling wines have any CO₂ removed by agitation or filtering through filter paper, or by short immersion in an ultrasonic bath.

Distillation procedure	Observations
1. The wine is poured up to the mark of the special flask ⁵	The temperature of the sample should be 20 °C. One can also operate at room temperature: from the sampling and the bringing back to volume after distillation a maximum deviation of ± 2 °C is allowed.
2. The content is poured into the distillation chamber of the SUPER DEE	
3. The flask is washed 3 time with a few ml of distilled water that is also poured into the chamber	
4. The sample is fully alkalized by the addition ⁶ of a suspension of milk of lime (2M), containing 120 g/litre CaO	The use of phenolphthalein indicator paper can be useful. However, the coulour change in the wine is obtained, visually, because of the alkalinity of the wine. This procedure is necessary to avoid the distillation of the other volatile substances (acetic acid, sulphurous anhydride, aldehydes, etc.), which will interfere with the density. In fact, the desired goal is to obtain a distillate similar to a hydroalcolohic solution, that is to say a pure alcohol-water mixture.
5. A few ml of diluted antifoam is added	This is an aqueous suspension of silicones or of oxyethylenic polymers. Foaming in the bulb slows the distillation and disrupts the passage of ethanol into the distillate. The presence of sugars in the sample enhances the foam formation.
6. Add approximately 2 ml distilled water to the same flask the sample was measured in. Collect approximately 75 ml of distillate therein.	The addition of water serves to avoid the loss of alcohol. The distillation is stopped automatically by a balance mechanism when the correct quantity has been distilled.
7. Add distilled water to the flask nearly up to volume. It is stoppered and brought to volume once temperature of 20 °C or the temperature of the sample before distillation is reached.	The use of a thermostatic bath is suggested
8. The flask is stoppered and mixed with care to avoid bubbles	The aim is to ensure a homogeneous distilate because, the determination of gravity will proceed from this.

⁵ The flasks are calibrated and interchangeable in mass, not in volume.

Remember that our flask explicitly DO NOT indicate ml, but only 100-200.

We DO NOT withdraw 100 ml, but a repeated quantity. So accuracy in ml is NOT necessary, but it is imperative that the volume of wine and distillate is the same at the same temperature.

⁶ The exact quantity depends on the sample.

5.2 THE DISTILLATION PROCEDURE FOR SPIRITS

Spirits present different characteristics from those of wines:

- Values of ASV up to 70%
- Variable composition

In fact spirits with very different composition exist: from cream whisky to fruit creams, from very sweet liqueurs such as vodka with fruit or limoncello to liqueurs with extracts of aromatic plants or roots, or even drinks with egg or coffee.

In addition, the higher alcoholic strength results in:

- a very slow distillation in order to make sure of the complete extraction of alcohol;
- collection of a greater volume of distillate (around 95 ml).

It is difficult to find a common distillation method for all the different kinds of spirits because the various ingredients interfere with the distillation in different ways, sometimes blocking the passage of the alcohol. The basic procedure for spirits is listed below. This procedure allows carrying out a good distillation in a short time, without excessive foam, and resulting in a clear distillate; the specific methods applicable to different alcohol products are indicated in the subsequent paragraph.

The distillation procedure	Observations
1. The sample is poured up to the mark of the special flask ⁷	It is recommended that the temperature of the sample is 20 °C. One can also operate at room temperature: a maximum deviation of ± 2 °C is allowed between the sampling and the bringing back to volume after distillation.
2. Pour the contents into the bulb of the distillatory	
3. Wash the flask 3 times with 10 ml of distilled water and add washings to distillatory bulb	
4. Add the “ solution for liqueurs ” ⁸ : solution 1:1 of NaCl 1% and [KAl(SO ₄) ₂] 12%	The aim is to increase the conductivity of the sample and to suppress the aromatic substances.
5. Add a few ml of diluted antifoam	Its addition avoids the formation of foam which would be an obstacle to the normal distillation of the alcohol.
6. Add about 2 ml cold distilled water to the same flask used for measurement of the sample. Collect \approx 95 ml of distillate therein.	
7. Add distilled water nearly up to volume, stopper and bring to volume when the temperature reaches 20 °C or the temperature of sample before distillation	USE OF A THERMOSTATIC BATH IS ADVISED
8. Stopper the flask and mix with care to avoid bubbles	

Regarding brandy, etc. the method is the same as the one indicated for wines, with the exception of the addition of a small quantity of lime milk: because of the low acidity of such products, 1-2 ml is sufficient. In case of alcoholic strength higher than 50 % vol, it is advisable to dilute to 1:2

⁷ See note 1

⁸ For the quantity see Table 1

5.3 COMPARISON BETWEEN THE SUPER DEE AND THE OIV GLASS DISTILLING UNIT

Tests of distillation with spirits with extremely varied compositions were carried out. Each liqueur was distilled :

- with the SUPER DEE distilling unit following the above procedure;
- with a glass distilling unit using the official procedures of OIV.

This latter distilling unit satisfies the official OIV condition.

The results obtained with the two apparatus are shown in the **TABLE 1**.

A specific method was developed to be used with the SUPER DEE distilling unit for each different sample type. Generally it is recommended:

- very sweet liqueurs: use small quantities of distilled water to wash the flask (less than 10 ml); in many cases this solution blocks the formation of emulsions which can result in an incorrect distillation;
- liqueurs containing milk or cream: the solution for liqueurs should not be added (it may curdle the milk creating a clot that is difficult to clean out of the flask).

The results obtained with the use of the SUPER DEE distilling unit are without doubt comparable with those obtained with the OIV glass distilling unit, with the advantage of much shorter distillation times. The differences found from the experimental data on some samples compared to the values indicated on the label fall within the tolerance expected by the current legislation (0,3 % vol on indicated ASV).

For further information see the book "Analysis methods of wines and alcoholic beverages", Gibertini – current Edition.

TABLE 1

SAMPLE	INGREDIENTS	DEE PROCEDURE	%vol DEE	Declared
Liqueur type Strega	Alcohol water, sugar, natural flavours, saffron	100 ml +1.5 ml antifoam +3 ml solution for liqueurs	40.37	40
Lemon cream	Water, sugar, lemon infusion, milk, cream, alcohol, glucose syrup	100 ml +12 ml antifoam	22.05	22
Rakee of raspberries	rakee, raspberries juice, sugar, natural flavours	100 ml +1.5 ml antifoam +3 ml solution for liqueurs	32.02	32
Limoncello	Alcohol, sugar, infusion of lemon rind, natural flavours	100 ml + 3 ml antifoam +10 ml milk of lime	30.29	30
Whisky cream	Dairy cream, sugar, Whisky, water, natural flavours	100 ml +15 ml antifoam	14.71	15
Rakee	Rakee	100 ml +5 drop antifoam +3 ml solution for liqueurs	40.12	40
Pina Colada	sugar, alcohol, cream, rum, milk, flavours	100 ml + 9 ml antifoam	16.80	17
China Martini	Sugar, alcohol, flavours, colorants: caramel; acidifier: citric acid	100 ml +15 ml antifoam + 10 ml milk of lime	31.00	31

(continue)

SAMPLE	INGREDIENTS	DEE PROCEDURE	%vol DEE	Declared
White Vermuth	Wine, sugar, alcohol, flavours	100 ml +1.5 ml antifoam + 5 ml milk of lime	16.06	16
Strong anisette	Alcohol, sugar, natural flavours	100 ml +1.5 ml antifoam +5 ml solution for liqueurs	40.02	40
Liqueur with coffee	sugar, alcohol, coffee, natural flavours	100 ml + 12 ml antifoam	32.89	33
Cherry	Alcohol, sugar, egriot juice, natural flavours	100 ml +1.5 ml antifoam +5 ml solution for liqueurs	30.03	30
Tequila	Distillate of agave techilana	100 ml +1.5 ml antifoam +2 ml solution for liqueurs	38.14	38
Pastis	Alcohol, sugar, 2 ⁰ / ₀₀ anetolo, infusion of vegetable substances, natural flavours,	100 ml +1.5 ml antifoam +5 ml solution for liqueurs	44.76	45
Banana Vodka	vodka, sugar, banana juice, flavours	100 ml +1.5 ml antifoam +3 ml solution for liqueurs	25.17	25
Brandy	Alcohol hydrate, Brandy, natural flavours, colorant: caramel	100 ml + 3 ml antifoam +2 ml solution for liqueurs	36.07	36
Liqueur with chocolate and mint	Milk, sugar, cream, glucose syrup, alcohol, milk protein, infusion of cacao, flavours, colorant: caramel	100 ml + 12 ml antifoam	16.06	16
Cachaca	Distillate of sugar-cane	100 ml + 3 ml antifoam +3 ml solution for liqueurs	43.66	43
Amaretto (alcoholized for exclusive use in confectionery)	Alcohol, caramel, flavours	25 ml + 75 ml distilled water (collect 100 ml) + 3 ml antifoam +3 ml solution for liqueurs	% vol value × 4 = 70.00	70
Sake	Distillate of rice	100 ml + 1 ml antifoam +10 ml solut. for liqueurs	15.85	16
Sambuca	Alcohol, sugars, anicette, natural flavours	100 ml + 5 ml antifoam +8 ml solution for liqueurs	39.86	40
Cointreau	Alcohol, sugars, orange flavours	100 ml + 5 ml antifoam +7 ml solution for liqueurs	39,81	40

NOTE

The samples are only indicative, so they must not be considered exhaustive.

6 ANALYSIS – Example of distillation

6.1 ANALYSIS TYPE SELECTION

Figure 4

80 / 160 for alcohol distillation
 240 for volatile

Push **SHIFT** to modify the data.
 The display shows + 0 -

TO SELECT THE KEY FOR USE :

80	=	100 ml	key A	} ALCOHOL
160	=	200 ml	key B	
240	=	240 ml	key C	} VOLATILE

6.2 SELECTION MODE - ALCOHOL / VOLATILE

To select ALCOHOL mode push the **A** key.

Confirm with the **MODE** key.

Figure 5

To select VOLATILE mode push the **B** key.

Confirm with the **MODE** key.

Figure 6

6.3 START DISTILLATION

Start distillation pushing the **A** key.

Figure 7

WARNING
Before starting the distillation check that the plug is indeed closed (turn clockwise the knob), the flask is installed on the balance and the protection is completely lowered. If there aren't these conditions the distillation don't start and an acoustic signal informs of the anomaly.

Example of distillation:
-plug and drainage closed
-flask on the balance
-protection down.

Figure 8

6.4 DISTILLATION CHAMBER CLEANING

After the distillation open the drainage and clean the distillation chamber with the shower supplied with SUPER DEE

Figure 9

When the SUPER DEE distilling unit is not in use leave the drainage open.

Figure 10

7 ERROR MESSAGES

There are indications for the main anomalies. The display shows the messages with the operative situations and the lights **WARNING** and **ALARM**.

Remember that, for a correct operation, the system must have a capacity of 5 litres/min. and that the water must have a temperature of 20°C max and a pressure of at least 2 bar.

1) In case of NO water or low water pressure:

Error !
Water press. low

THE INSTRUMENT DOES NOT START (Fig. 11)

2) After 10 seconds the display shows the message

No water

THE INSTRUMENT DOES NOT START

If the error situation continues the distillation will be cancelled and the program will wait to confirm the reading of the message by operator.

No water
Ok

CONFIRM THE READING OF THE MESSAGE.

3) In case of cooling water too hot the yellow light flashes. Push the START key to start a distillation. The instrument will wait 10 seconds and if the error situation remains the operator will select if to continue or cancel the distillation.

Water too hot
Go on Cancel

CONFIRM OR CANCEL THE DISTILLATION (Fig. 12)

ATTENTION!
The temperature of the cooling water upper 20 °C compromises the result of the distillation whether in alcohol or Volatile acidity mode.


Figure 11


Figure 12

Many situations of error have the indication of the error message until the restored normal situation.

4) In case of drainage open the following message will appear:

Attention!
Drainage open

*THE DISTILLATION DON'T START.
THE INSTRUMENT WAITS FOR THE CLOSURE OF THE DRAINAGE (Fig. 13)*

5) In case of open plexiglas the following message will appear:

Attention!
Plexiglass open

*THE DISTILLATION DOES NOT START.
THE INSTRUMENT WAIT FOR THE CLOSURE OF THE PLEXIGLASS (Fig. 14)*

Figure 13


Figure 14

8 MAINTENANCE

8.1 DISTILLATION CHAMBER

After each distillation the distillation chamber must be washed with water to prevent contamination of the successive sample.

Periodically clean the distillation chamber introducing approximately 500 ml of distilled or demineralised water and approximately 20 ml of sodium chloride NaCl 1% solution. Start the distillation in ALCOHOL MODE.

To remove deposits or incrustations that the samples can have left on the electrodes and/or on the glass can be done in the two following ways:

- 1- Pour into the distillation chamber a dilute solution of HCl 0,1N. Await some hour and eventually use a tube-brush to clean the chamber.
- 2- Introduce into the distillation chamber the solution B, supplied with the instrument, in sufficient quantity to cover the incrustations. Also in this case to await some hour and use eventually a tube-brush.

In any case rinse always accurately the distillation chamber before reuse the instrument.

WARNING! Don't press the START button when into the distillation chamber there is a cleaning solution. For more safety remove the power socket from the main.

8.2 ELECTRODES

The electrodes will gradually diminish in size and will have to be replaced at some stage. It is difficult to predict the frequency of replacement because it will depend on total usage and type of sample (approx. 2000 tests).

For the replacement call the SAT (Service) of GIBERTINI - Phone (+39) 023541434 or your supplier.

8.3 RECTIFICATION COLUMN

During the normal use of the SUPER DEE, the column doesn't require a particular cleanness. However, if the column is contaminated by acid residues, it is necessary to clean it in this way :

- fill the distillation chamber with about 500 ml of distilled water and 20 m of sodium chloride NaCl 1 % solution. Start distillation in ALCOHOL MODE.

It is not necessary to collect the distillate. Stop the distillation holding pressed toward the balance for some second.

8.4 EMPTYING FUNCTION

- Switch off the DISTILLATOR with the ON/OFF SWITCH.
- Switch on pressing simultaneously the **B** key on the keyboard.
- The display shows the message "**Take off water tube Wait for tank drain**"
- Empty the cooling coil using the shower.
- Wait for the complete emptying of the distilled water tank by means of the steam pump (the water finishes to go out by the nozzle and the pump becomes very noisy). Switch off the instrument.
- Incline the instrument backwards to help the exit of water.

8.5 DEIONISATION COLUMN (REFILL R3)


Figure 15

Replace the column when it is completely red (after about 400 volatiles). It is advisable however to replace the cartridge when it is black only in the last 3-5 cm.

Replacement of the cartridge:

- unscrew the nuts that keep the cartridge in position;
- incline the cartridge leaving the bottom end in position. Close the water inlet hole (on top) with a finger to prevent water from spilling over the components underneath;
- substitute the cartridge paying attention that the inlet and outlet gaskets of the cartridge are in the correct position;
- screw in the fixing nuts.


Verify that the water pipes do not have bents that may prevent water to flow freely.

8.6 BOILERS

If the display shows **CHECK BOILERS**, push the 2 red keys (Fig. 15).

NOTE

Arrange a possible periodic maintenance with the GIBERTINI SAT (Service) - tel. (+39) 02 3541434

customer care

info@gibertini.com

Tel. (+39) 023541434

ISO 9001:2008

SIT n. 94

GIBERTINI

ELETTRONICA s.r.l.

Via Bellini 37
20026 Novate Mil.se
Milano (Italy)

Tel. (+39) 02 3541434
Fax (+39) 02 3541438
www.gibertini.com
sales@gibertini.com