**LEVITAS®**

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---

**EXCEL SERIES  
INSTRUCTIONS MANUAL**

**MORETTI S.P.A.**  
Via Bruxelles 3 - Melegnano  
52022 Cavriglia (Arezzo)  
Tel. +39 055 96 21 11

Fax +39 055 96 21 200  
[www.morettispa.com](http://www.morettispa.com)  
[info@morettispa.com](mailto:info@morettispa.com)

## INDICE

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I Class Medical Device  
according to Directive 93/42/EEC  
and further modifications

## CODE

- LTM661** Kit mattress and pump without pressure adjustment (EXCEL 1000 MINI)  
**LTM662** Kit mattress and pump with pressure adjustment (EXCEL 2000 DELUXE)  
**LTM663** Kit mattress with interchangeable cells and pump with pressure adjustment (EXCEL 4000 DELUXE)  
**LTM603** Kit mattress and Junior pump with pressure adjustment (EXCEL 2000 JUNIOR)

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Thank you for purchasing a LEVITAS kit mattress and pump Excel by Moretti S.p.A. product range. This user's manual provides some suggestions as how to correctly use the product you have chosen and gives some valuable advice for your safety. Please read through the manual carefully before using the product. Should you have any queries, please contact your retailer for advice and assistance. In case of doubts, we recommend to contact the dealer who can help and advice you properly.

## INTENDED USE

The kit pump and mattress have been designed for the bed sore treatment and prevention, which may occur during long hospital and/or house confinement.

### CAUTION!

- Do not use the product for a purpose not indicated in this manual
- Only qualified personnel or trained for the treatment and prevention of bed sore, can operate with the device
- Moretti S.p.A declines all responsibilities for any consequences resulting from an incorrect use of this product and from unauthorized alteration to the frame of the product.
- The manufacturer reserves the right to change the information contained in this document without previous notice.

## EC DECLARATION OF CONFORMITY

Moretti S.p.A. hereby declares, under its own exclusive responsibility, that the products manufactured and sold by Moretti S.p.A. and belonging to the kit pump and mattress category are in compliance with the applicable provisions outlined in the Medical Devices Directive 93/42/EEC and further modifications.

To this end Moretti S.p.A. guarantees and declares as follows, under its own responsibility:

1. The devices in question satisfy the essential requirements of Annex I of directive 93/42/EEC as prescribed by Annex VII of the above directive
2. The devices in question ARE NOT MEASUREMENT DEVICES
3. The devices in question ARE NOT INTENDED FOR CLINICAL INVESTIGATION
4. The devices in question are sold in NON STERILE packing

5. The devices in question should be considered as belonging to class I
6. DO NOT use/install the devices in question for any other purpose than that expressly indicated by Moretti S.p.A.
7. Moretti S.p.A. maintains and submits the documentation certifying conformity with Directive 93/42/EEC to the proper authorities and further modifications

**STANDARDS**

The system has been tested and successfully approved to the following standards:

- EN 60601-1
- EN 60601-1-2

**GENERAL WARNINGS**

- DO NOT use this product or any available optional equipment without first completely reading and understanding this instruction manual. If you are unable to understand the warnings, cautions or instructions, contact a healthcare professional, dealer or technical personnel before attempting to use this equipment, otherwise injury or damage may occur
- Advise with physician or therapist to determinate the correct adjustment and the correct use of the device
- Keep the packed kit away from heat sources
- SERVICE LIFE- the device use limit is defined by the wear of the parts
- DO NOT allow children to play on or operate the pump

**ELECTROMAGNETIC COMPATIBILITY WARNINGS**

This equipment has been tested and it complies with medical devices limits belonging to EN 60601-1-2:2007. These limits have been designed to provide reasonable protection against harmful interference in a typical medical installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to other devices in the vicinity. However, there is no guaranteed that the interference will not occur in a particular installation. If this equipment does cause harmful interference to other devices, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference following one or more of the suggestions below:

- Reorient or relocate the receiving device
- Increase the separation between the equipment
- Connect the equipment into an outlet on a circuit different from that to which the other devices are connected
- Consult the manufacturer or field service technician for help

**SYMBOLGY**

 Code	 Warning	 Temperature limit	 "BF" type	 Follow instruction manual	 Class II isolation	
 Production date	 Disposal according to EC/2002/96	 ON	 OFF	 Manufacturer	 Lot number	<b>IP2X</b> The device is protected against foreign objects the size of 12,5mm or more

**KIT CONTENT**

The EXCEL series is composed by four kind of kit which can satisfy all the requirement about treatment and prevention of complications caused by long hospital and/or house confinement, in particular bedsores.

*LTM661 KIT EXCEL 1000 Mini - Bubble mattress and pump without pressure adjustment*

*LTM662 KIT EXCEL 2000 Deluxe - Bubble mattress and pump with pressure adjustment*

*LTM663 KIT EXCEL 4000 Deluxe - Mattress with interchangeable cells and pump with pressure adjustment*


*LTLTM603 KIT EXCEL 2000 Junior - bubble mattress and Junior pump with pressure adjustment*

UNIT CODE	DESCRIPTION	LTM661	LTM662	LTM663	LTM603
<b>LTM659</b>	Pump without pressure adjustment	•			
<b>LTM660</b>	Pump with pressure adjustment		•	•	
<b>LTM669</b>	Junior pump with pressure adjustment				•
<b>LTM655</b>	Bubble mattress with flaps	•	•		•
<b>LTM651</b>	Mattress with interchangeable cells 87x200cm			•	
<b>RV355</b>	PVC mattress mending kit	•	•		•
<b>LTM670</b>	Spare tubes for PVC mattress LTM655	•	•	•	•
-----	Instruction manual	•	•	•	•

## GENERAL DESCRIPTION

### Mattress and pump

**NOTE** Pictures are only for demonstration and do not refer to a specific product


1. Front panel
2. Mounting brackets
3. Air hose port
4. Power cord
5. Mattress system or bubble pad system
6. Pump unit
7. Main power switch
8. Pressure adjust knob

**INSTALLATION**

1. Place the mattress on a common mattress, with the air hose connectors oriented to the footboard

**NOTE** lean the mattress on a 5cm mattress



2. Hang the pump onto the bed rail (footboard side), or place it on a flat surface

3. Connect air hose connectors from air mattress to the pump unit. Ensure the air hoses are not kinked or tucked under mattress

4. Plug the power cord into the electrical outlet

5. Turn the switch ON

**WARNING!**

- Make sure the pump unit is suitable for the local power voltage
- Do not position the equipment so that is difficult disconnect the plug

**BEFORE USE**

Check the correct assembly, paying particular attention to air hoses connector

**USING WARNINGS****DANGER -to reduce risk of electrocution-**

1. Always unplug this product immediately after using
2. Do not use while bathing
3. Do not place or store this product where it can fall or be pulled into a tub or sink
4. Do not place in or drop into water or other liquid
5. Do not reach for a product that has fallen into water. Unplug immediately

**WARNING -to reduce risk of burns, electrocution, fine or injury to people**

1. Evaluate patients for entrapment risk according to facility protocol and monitor patient appropriately
2. This system is not intended to be used with patients who have spinal cord injury.
3. Close supervision is necessary when this product is used on or near children. Electrical burns or choking accident may occur if a child swallow a small part detached from the device
4. Use this product only for its intended use as described in this manual. Do not use attachment not recommended by the manufacturer
5. Never operate this product if it has a damaged cord or plug, if it is not working properly, or if it has been dropped, damaged, or dropped into water. Return the product to a service centre for examination and repair
6. Keep the cord away from heated surfaces
7. Never block the air openings of this product or place it on a soft surface, such as a bed or couch, where their openings may be blocked. Keep the air openings free of lint, hairs, and other similar particles
8. Never drop or insert any object into any opening or hose
9. Do not modify this equipment without authorization of the manufacturer
10. Mattress cover has passed skin sensitization and skin irritation test. However, if you suspect that you may have had or are having an allergic reaction, please consult a physician immediately
11. If there is a possibility of electro-magnetic interference with mobile phones, please increase the distance to 3,3m between devices or turn off the mobile phone

**HOW TO USE**

The EXCEL mattress and pump should be used on a common mattress. Referring to the interchangeable cells mattress, we suggest to lean on a 5cm mattress. Ensure the mattress is always clean and in good condition. For the mattress cleaning follow the appropriate chapter. Ensure the mattress is duly blown up before laying the patient on. Keep the pump turned on during all the confinement period.

The LTM662, LTM663 and LTM603 series are equipped with a pressure adjustment which must vary depending on the patient weight. Follow the table to choose the correct pressure.

ONLY FOR LTM662 AND LTM663 (because provided with a graduated scale on the pressure knob):

Patient wight	Position on the pressure knob
Up to 40 kg	2-3
from 40 to 60 kg	3-4
from 60 to 70 kg	4-5
from 70 to 80 kg	5-6
from 80 to 90 kg	6-7
from 90 to 100 kg	8

**ONLY FOR LTM603 (because is not provided with a graduated scale on the pressure knob):**

1. For the first usage it is recommend the pressure knob to be set to "FIRM" for a quickest inflation
2. Once the patient lay on the mattress verify the pressure is correct by putting an hand between the mattress and the patient. The contact must be minimum. If it is necessary decrease the pressure by rotating the knob to "SOFT"

**NOTE** In case the pressure is constantly low, check for any holes (on the interchangeable elements or on the connection tubes). If it is necessary substitute the damaged connection tubes or the interchangeable elements, otherwise contact your dealer.

If a rapid deflation is necessary unplug the connection tube from the mattress

### STORAGE

- Fold the mattress accurately to prevent creases
- Roll from the head end towards the foot end with the valve open

### MAINTENANCE

The LEVITAS devices by MORETTI are checked carefully and supplied with EC mark, once launched on the market.

For patient and physician's safety we recommend to let check the device by the manufacturer or by laboratory authorized every year.

In case of repair, use only original spare parts and accessories.

## CLEANING AND DISINFECTION

### WARNING!

Unplug the pump before starting the cleaning procedure

Before using the device with people follow this procedure:

- Clean the mattress with mild soap
- Do not immerse or soak pump unit
- After cleaning, dry the mattress without direct sunlight exposure.
- Check the connection tubes are not obstructed by dust or else. In this case use a soft swab to clean the tubes

**NOTE** Never use acids, alkalis or solvents such as acetone

## CONDITIONS OF DISPOSAL

### -General conditions of disposal

In case of disposal don't use the inserting container for municipal waste. We recommend to dispose the kit in the appropriate disposal areas for recycling.

### TREATMENT OF BATTERIES (Directive 2002/96/EC)

At the end of its life, the product must not be disposed of along with other domestic waste. The users must dispose of this equipment by bringing it to a specific recycling point for electric and electronic equipment or at retailers that provide this service. By ensuring these batteries are disposed of correctly, you will help prevent potentially negative consequences for the environment and human health which could otherwise be caused by inappropriate waste handling of the battery. To remark the need to dispose of electrical equipment separately, the products is marked with crossed mobile waste bin

## SPARE PARTS/ACCESSORIES

For spare parts and accessories refer to main catalogue

## TROUBLESHOOTING

If your questions can't be answered with above information, please contact your local agent directly. They might require a technician to take care of the problem

PROBLEM	SOLUTION
<b>Pump is not ON</b>	<ul style="list-style-type: none"> <li>• Check if the plug is connected to mains</li> <li>• Check if there is any blown fuse</li> </ul>
<b>The mattress doesn't inflate</b>	<ul style="list-style-type: none"> <li>• Check if the connection between mattress and pump unit is tightly secured</li> <li>• Check if there is any leakage on air cells</li> </ul>
<b>Mattress form is loose</b>	<ul style="list-style-type: none"> <li>• Check if all the snap buttons or straps of mattress are all securely fastened</li> <li>• Check if the mattress is fixed to the bed frame by straps</li> </ul>

## ELECTROMAGNETIC COMPATIBILITY DECLARATION

**Guidance and manufacturer's declaration- Electromagnetic emissions-**

This device is intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

EMISSION TEST	COMPLIANCE	ELECTROMAGNETIC ENVIRONMENT- GUIDANCE
Rf emissions CISPR 11	Group 1	The device uses RF Energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment
RF emissions CISPR11	Class B	The device is suitable for use in all establishment. Including domestic establishments and those directly connected to the public low-voltage power supply network.
Harmonic emissions IEC61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC61000-3-3	Complies	

**Guidance and manufacturer's declaration- Electromagnetic immunity**

The Excel devices are intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

Immunity Test	IEC60601 test level	Compliance	Electromagnetic Environment-Guidance
Electrostatic Discharge (ESD) IEC61000-4-2	± 6kV contact ± 8kV air	± 6kV contact ± 8kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/ burst IEC61000-4-4	± 2kV for power supply line ± 1kV for input/out line	± 2kV for power supply line ± 1kV for input/out line	Mains power quality should be that of atypical commercial or hospital environment
Surge IEC61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s)	Mains power quality should be that of atypical commercial or hospital environment.

<p>Voltage dips, short interruptions and voltage variations on power supply input lines IEC61000-4-11</p>	<p>&lt;5 % UT (&gt;95 % dip in UT) for 0,5 cycle  40 % UT (60 % dip in UT) for 5 cycles  70 % UT (30 % dip in UT) for 25 cycles  &lt;5 % UT (&gt;95 % dip in UT) for 5 sec</p>	<p>&lt;5 % UT (&gt;95 % dip in UT) for 0,5 cycle  40 % UT (60 % dip in UT) for 5 cycles  70 % UT (30 % dip in UT) for 25 cycles  &lt;5 % UT (&gt;95 % dip in UT) for 5 sec</p>	<p>Mains power quality should be that of atypical commercial or hospital environment. If the user of this device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.</p>
<p>Power frequency (50/60Hz) magnetic field IEC61000-4-8</p>	<p>3 A/m</p>	<p>3 A/m</p>	<p>Power frequency magnetic fields should be at levels characteristic of atypical location in a typical commercial or hospital environment.</p>
<p><b>NOTE</b> UT is the a.c. mains voltage prior to the application of the test level</p>			

**Guidance and Manufacturer's Declaration - Electromagnetic Immunity:**

This Excel devices are intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

Immunity Test	IEC60601 test level	Compliance	Electromagnetic Environment-Guidance
<p>Conducted RF IEC 61000-4-6</p>	<p>3Vrms150 kHz to 80 MHz outside ISM bandsa</p>	<p>3 Vrms</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of this device, including cables, than there commended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance  <math>d = 1.2\sqrt{P}</math> 150 kHz to 80 MHz  <math>d = 1.2\sqrt{P}</math> 150 kHz to 80 MHz  <math>d = 2.3\sqrt{P}</math> 80 MHz to 2.5G MHz                      Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).b                      Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey c, should be less than the compliance level in each frequency ranged. Interference may occur in the vicinity of equipment marked with the following symbol: </p>
<p>Radiated RF IEC 61000-4-3</p>	<p>3 V/m 80 MHz to 2.5 GHz</p>	<p>3 V/m</p>	<p>(This cell is merged with the one above for content purposes)</p>

**NOTE 1** At 80 MHz and 800 MHz, the higher frequency range applies.

- NOTE 2** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.
- NOTE 3** Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and landmobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device
- NOTE 4** Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

**Recommended separation distances between portable and mobile RF communications equipment and this device:**

This device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of this device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and this device as recommended below, according to the maximum output power of the communications equipment. For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 KHz to 80MHz $d=1.2\sqrt{P}$	80MHz to 800MHz $d=1.2\sqrt{P}$	800MHz to 2.5 GHz $d=2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

- NOTE 1** At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.
- NOTE 2** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

## TECHNICAL FEATURES

-Pump unit-

FEATURES	LTM661/662/663	LTM603
<b>Power supply</b>	AC 230V ~ 50 HZ 0,05A	AC 230V ~ 50 HZ 0,05A
<b>Power</b>	10 W	10 W
<b>Fuse rating</b>	T1AL 250V	T1AL 250V
<b>Pressure</b>	LTM661 50-105 mmHg LTM662/663 50-105 mmHg	30-90 mmHg
<b>Dimensions</b>	LTM661- 21,5x11x9,5 cm LTM662- 25x13x9,5 cm LTM663- 25x13x9,5 cm	LTM603 25x12,5 x8,5 cm
<b>Cycle</b>	6 minutes	9,5 minutes
<b>Weight</b>	LTM661- 1,35 KG LTM662/663- 1,6 KG	LTM603- 1,4 KG
<b>Environment temperature</b>	Operation: +10°C / +40°C Storage: -15°C/+50°C Shipping: -15°C/+70°C	Operation: +10°C / +40°C Storage: -15°C/+50°C Shipping: -15°C/+70°C
<b>Environment humidity</b>	Operation: 10%-90% non- condensing Storage: 10%-90% non- condensing Shipping: 10%-90% non- condensing	Operation: 10%-90% non- condensing Storage: 10%-90% non- condensing Shipping: 10%-90% non- condensing
<b>Atmospheric pressure</b>	Operation: 700-1013,25 hPa	Operation: 700-1013,25 hPa
<b>Classification</b>	Class II, Type BF, IP2X (keep dry) Applied part: air mattress Not suitable for use in the presence of a flammable anesthetic mixture (No AP or APG protection)	Class II, Type BF, IP2X (keep dry) Applied part: air mattress Not suitable for use in the presence of a flammable anesthetic mixture (No AP or APG protection)

**-Mattress-**

<b>FEATURES</b>	<b>BUBBLE MATTRESS (LTM655)</b>	<b>INTERCHANGEABLE ELEMENTS MATTRESS ( LTM651)</b>
<b>Dimensions (without flaps)</b>	Inflated 200 x 90 cm Deflated 240x105 cm	200X90X12,5 mm
<b>Cell height (inflated)</b>	6,4 mm	12,5 mm
<b>Thickness</b>	0,64 mm (+0,02%)	-----
<b>Single layer thickness</b>	0,32 mm	-----
<b>Material</b>	PVC EN 71	Nylon + PU
<b>Elements number</b>	-----	18
<b>Max weight</b>	100 kg	120 kg
<b>Weight</b>	2,4 kg	4,3 kg
<b>Pressure</b>	50-105 mmHg	50-105 mmHg

## WARRANTY

Moretti products are guaranteed from material or manufacturing faults for 2 years from the purchasing date, except possible exclusion or restriction as follows. The warranty shall not be applied in the possible damages caused by improper use, abuse or alteration, and the warranty will not be valid if the instructions for use are not strictly adhered to. The correct intended use is specified in this manual.

Moretti is not responsible for consequent damages, personal injuries or whatever caused by or in relation to wrong installation or improper use.

Moretti warranty does not cover damages resulting from: natural disaster, not authorized maintenance or repairs, faults caused by problems on electricity supply (when necessary), use of spare parts not covered by Moretti, improper use, not authorized alteration, shipment damages (different from original Moretti shipment), or in case of insufficient maintenance as indicated in the manual.

The warranty doesn't cover components subject to wear and tear during the correct use of the device.

## REPAIRING

### **-Warranty repair-**

If a Moretti item presents material or manufacturing faults during the warranty period, Moretti will confirm with customer if the fault can be covered from warranty. Moretti, at its unquestionable discretion, can repair or replace the item, by a Moretti dealer or to Moretti headquarters. Labor cost can be charged to Moretti if the repair is covered by warranty. A repair or a replacement doesn't extend the warranty.

### **-Repair not covered by warranty-**

A product out of warranty can be sent after Moretti authorization. The labor and shipping costs for good out of warranty are to be paid by the customer or by the dealer. The repairs are guaranteed for 6 months from the good received.

### **-Non-defective device-**

The customer will be informed if, after the device return and examination, Moretti declares that the device is not faulty. In this case the good will be sent back to customer, the shipping charge will be paid by the customer.

**SPARE PARTS**

The original Moretti spare parts are guaranteed for 6 months from the receiving date.

**EXEMPT CLAUSES**

Moretti does not offer any other declaration, explicit or implicit warranty or conditions, including possible declarations, warranties or conditions of merchantability, fitness for a specific purpose, non infringement and non interference, all but what expressly specified in this warranty. Moretti does not guarantee the nonstop and faultless usage. The duration of possible implicit warranties which can be imposed by the law is limited by the warranty period, in the limits of law. Some states or countries don't allow limiting the implicit warranty or the exclusion or the limitation for accidental damages. In such countries, some of those exclusions or limitation may not be applied to the user. The present warranty may be modified without prior notification.

**WARRANTY CERTIFICATE**

**Product** \_\_\_\_\_

**Purchased on (date)** \_\_\_\_\_

**Retailer** \_\_\_\_\_

**Address** \_\_\_\_\_ **Town/city** \_\_\_\_\_

**Sold to** \_\_\_\_\_

**Address** \_\_\_\_\_ **Town/city** \_\_\_\_\_

**MORETTI S.P.A.**

Via Bruxelles, 3 - Meleto 52022 Cavriglia (Arezzo) Tel. +39 055 96 21 11

[www.morettispa.com](http://www.morettispa.com) email: [info@morettispa.com](mailto:info@morettispa.com)

**MADE IN P.R.C.**

**Distribuito da / Distributed by:**

Via Bruxelles, 3 - Melegnano  
52022 Cavriglia (Arezzo)  
Tel. +39 055 96 21 11

Fax. +39 055 96 21 200  
[www.morettispa.com](http://www.morettispa.com)  
[info@morettispa.com](mailto:info@morettispa.com)