# StarPhoresis

## **Wide Mini Gel Horizontal Electrophoresis System**


N2023-1410

## **Long Gel Horizontal Electrophoresis Systems**

N2113-2510 N2120-2510

## **Wide Long Gel Horizontal Electrophoresis Systems**

N2123-2510 N2123-4010

User Manual • Manuel d'utilisation • Bedienungsanleitung • Manuale Utilizzatore

<b>English</b>	2 – 12
<b>Français</b>	13 – 23
<b>Deutsch</b>	24 – 35
<b>Italiano</b>	36 – 46

---

## ENGLISH

### StarPhoresis Wide Mini Gel Horizontal Systems

Catalogue No. N2023-1410

### StarPhoresis Long Gel Horizontal Systems

Catalogue Nos. N2113-2510, N2120-2510

### StarPhoresis Wide Long Gel Horizontal Systems

Catalogue Nos. N2123-2510, N2123-4010

**Please read the entire User Manual before operating this unit. These units are capable of delivering potentially lethal voltage when connected to a power supply and should be operated only by qualified and technically trained personnel.**


### Warranty

Please check that the unit has been received complete and undamaged. Refer to the illustration below and to *TABLE A* to check that all components are present. Save all packaging and documents until you have thoroughly inspected your shipment; if you find that your order is incorrect or damaged please call STARLAB for instructions for its return.

This warranty is valid for 36 months, but only if the product has been used and cared for according to this User Manual. No liability is accepted for loss or damage arising from incorrect use. The liability is limited to the repair or replacement of the unit, or a refund of the purchase price, at STARLAB's option. STARLAB is not liable for any consequential damages. We refer to our general terms for sale.

STARLAB reserves the right to alter the specifications of the Electrophoresis Systems without prior notice. This will enable us to implement improvements as soon as they are available.

### StarPhoresis Horizontal Electrophoresis System components:

**TABLE A:** System components included

	N2023-1410	N2113-2510	N2120-2510	N2123-2510	N2123-4010
<b>Buffer Chamber</b> (inc. levelling base)	1	1	1	1	1
<b>Cover</b> (w/power cords)	1	1	1	1	1
<b>Gel Tray</b> (gasketed)	1	1	1	1	1
<b>Combs</b>	4	3	3	4	4

**TABLE OF CONTENTS**

I.	ENVIRONMENTAL CONDITIONS FOR USE .....	3
II.	SAFETY PRECAUTIONS .....	4
III.	INTRODUCTION .....	4
IV.	GENERAL CARE AND CLEANING .....	5
V.	SPECIFICATIONS .....	6
VI.	CASTING THE GEL .....	6
VII.	OPERATION .....	7
VIII.	VISUALISATION .....	8
IX.	TROUBLESHOOTING .....	8
X.	COMB SPECIFICATIONS .....	9
XI.	ACCESSORIES .....	12

**I. ENVIRONMENTAL CONDITIONS FOR USE**

- This unit is intended for indoor use only.
- This unit can be operated safely at an altitude of 2,000m.
- The normal operating temperature range is between 4°C and 65°C.
- Maximum relative humidity 80% for temperatures up to 31°C, decreasing linearly to 50%, relative humidity at 40°C.

## II. SAFETY PRECAUTIONS

Please read the User Manual carefully before using the STARLAB Horizontal Electrophoresis Unit. This manual contains important operating and safety information. STARLAB's electrophoresis units are designed to perform flawlessly for years in the most demanding laboratories. Please take the time to read this manual so you understand the safety and operating instructions to ensure the successful use of the unit. Alterations could cause serious injury to the user or the system.

Power to the unit is supplied by an external power supply. The power supply must meet safety standards for IEC 1010-1 regulations and must be ground isolated and incorporate a no-load detecting circuit.

Power is supplied to the gel through the cover of the system providing a safety interlock to the user. Users should not attempt to operate this unit without the safety inter-locked cover in place.

- To avoid the risk of personal shock always disconnect the unit from the power supply before:  
a) removing cover , b) moving the unit or, c) adding running buffer.
- Running conditions should not exceed the maximum operating voltage or current.
- DO NOT fill the buffer chamber with running buffer above the maximum fill line.
- Use this apparatus only for its intended purpose and as described in this manual.
- Do not use the product if the power cords are damaged or if its surfaces are cracked.
- Do not operate electrophoresis units in metal trays.

**CAUTION:** During electrophoresis very low quantities of various gases are produced at the electrodes. The type of gas produced depends on the composition of the buffer employed. To disperse these gases make sure that the unit is run in a well-ventilated area.

## III. INTRODUCTION

Thank you for purchasing a StarPhoresis Long Gel Horizontal Electrophoresis System. Our systems have increased capacity and are ideally suited for the analysis of up to 600 sample sets of DNA or RNA, or labs that need to clearly resolve samples containing a large number of fragments. All units incorporate our tape-free gel casting to expediate gel preparation.

The UV transparent gel trays are available with a variety of comb slot configurations. Combs are available in 1.0mm and 1.5mm thicknesses, as well as a variety of tooth configurations to accommodate varying numbers of samples and sample volumes.

### **Outstanding features ensure trouble-free use**

- Robust, acrylic construction stands up to daily usage without breakage, warping or leakage.
- All of STARLAB's Long Horizontal Electrophoresis Systems are hand-made to provide unmatched durability for years of daily use.
- Systems meet or exceed all safety requirements for the IEC 1010-1 Standards.
- Gel trays accommodate gasketed end gates, allowing casting of gels in the unit without tape, dams or external casting stands.

## IV. GENERAL CARE AND CLEANING

**WARNING!** Acrylic is not resistant to aromatic or halogenated hydrocarbons, ketones or esters. Organic solvents cause acrylic to “craze” or crack. Do not use ethanol or other organic solvents to clean the unit.

- Do not autoclave, bake or microwave the unit.
- Before using, clean and dry the unit. Clean with DISTILLED WATER ONLY; dry parts with clean, lint-free lab wipes or, preferably, air dry. Use care when cleaning or drying the unit near the platinum wire. The connectors should be clean and dry before usage or storage.
- Do not use abrasive creams or scourers.
- Do not use cleaning brushes in the electrode area.
- A thorough rinse with distilled water is all that is generally required to clean the unit after use. A mild detergent may also be used. Acrylic can also be exposed to a mild bleach solution (10:1). In addition RNase removal products are also safe for acrylic.

### Gel tray

All STARLAB gel trays are made from UV transmissible acrylic, which allows a stained agarose gel to be viewed directly in the tray on a transilluminator.

**NOTE:** care should be taken to avoid pouring freshly melted agarose directly into the gel tray, as this may damage the tray. **Always allow agarose to cool to 60°C before pouring into the tray.**

### Combs

Combs are available in thicknesses of 1.0mm and 1.5mm. Comb spines are acrylic and the teeth are polycarbonate. Combs designated “Microtitre” are compatible for use with a multi-channel pipettor. Each comb is labelled with its thickness, number of teeth and catalogue number. Please see Chapter 10, *Comb Specifications*, for details of the combs available for each device.

### Buffer Tank

The buffer tank holds the buffer and gel tray. The larger systems are manufactured with a three point levelling base. The need to cast and run larger gels on a level surface is very important for consistent reproducible results. Current is passed from the power supply through the leads to the cover connectors, through the platinum wires and into the buffer and from there into the gel. The cover incorporates safety rings that prevent human contact with the electrode connections to prevent electrical shock if the power supply is not turned off prior to removing the lid. Users should not attempt to operate the system without the lid in place. The electrode banana plugs are colour-coded red and black to ensure that the gel is placed in the correct orientation in the tank **[almost all electrophoretic separations of biomolecules run from the black electrode (negative) to the red electrode (positive)]**. The connection between the banana plug and the electrode wire is sealed with a proprietary composite that eliminates corrosion of the connection. The electrode wires are 0.254mm platinum and are secured to the inner walls of the tank. Care should be taken not to damage the electrodes when cleaning the unit.

## V. SPECIFICATIONS

**TABLE B:** General specifications

	<b>N2023-1410</b>	<b>N2113-2510</b>	<b>N2120-2510</b>	<b>N2123-2510</b>	<b>N2123-4010</b>
Gel dimensions (w x l)	23.5 x 14cm	13 x 25cm	20 x 25cm	23.5 x 25cm	23.5 x 40cm
System dimensions (w x l x h)	26 x 29 x 10cm	21.5 x 37.5 x 10.5cm	29 x 37.5 x 10.5cm	32 x 37.5 x 10.5cm	32 x 53 x 10.5cm
Tray configuration	4 comb slots	6 comb slots	6 comb slots	14 comb slots	12 comb slots
UVT gel tray	Runs 1, 2, 3, or 4 equal length gels	Runs 1, 2, 3 or 4 equal length gels	Runs 1, 2, 3 or 4 equal length gels	Runs 1, 2, 3, 4, 5 or 10 equal length gels	Runs 1, 2, 3, 4, 6 or 12 equal length gels
Agarose capacity (0.5cm T)	161ml	163ml	250ml	294ml	460ml
Max. buffer volume	1000ml	1600ml	2300ml	3000ml	4500ml
Buffer exchange ports	optional	optional	optional	optional	yes
Max. sample capacity	200	112	168	500	600
Max. sample run, full length	50	28	42	50	50
Comb thickness	1.5mm, plus optional 1.0mm				
Comb specifications	4 x MT50 teeth, 1.5mm thickness	12, 16 & 20 teeth, 1.5mm thickness	16, 24 & 36 teeth, 1.5mm thickness	2 x MT25 & 2 x MT50 teeth, 1.5mm thickness	2 x MT25 & 2 x MT50 teeth, 1.5mm thickness
Safety lid with attached power cords	1	1	1	1	1

## VI. CASTING THE GEL

1. Remove the lid from the gel box by holding the front of the buffer chamber with one hand and pulling the lid off by holding the centre of the back of the lid. The cover is attached to the back of the unit at the connection of the power cords to the banana plugs located on the unit.
2. To cast gels using the StarPhoresis gel tray, place the gel tray into the chamber making sure the gel tray rests level and is centred on the platform of the buffer chamber. Slide the gasketed end gates into the outermost grooves on either side of the gel tray. The end gates should be inserted tightly into the grooves with the gasket side facing out.
3. The larger systems are manufactured with a three-point levelling base. The need to cast and run larger gels on a level surface is very important for consistent reproducible results. Level the unit using the two thumb screws on the front of the base by slowly turning one thumb screw at a time, lining up the bubble in the level with the centre circle.
4. Prepare agarose as described in your individual application.
5. **IMPORTANT:** to avoid damage to UVT tray, always cool agarose to 60°C before pouring into tray. Any bubbles that form while pouring may be popped with a small glass pipette or the edge of the comb. Place combs into the comb slots on the gel tray being careful not to create bubbles along the teeth of the comb. If the entire gel length is required for separation, place one comb in the first comb slot on the tray. To separate more samples on one gel, use both comb slots on the tray. Allow gel to solidify for 30 minutes.

## VII. OPERATION

1. Once the agarose has solidified, carefully remove the end gates to expose the open ends of the agarose to the buffer. Be sure the gel tray is placed in the proper orientation; the first comb should be closest to the cathode (black) electrode. The banana plugs have been sealed against corrosion with colour-coded sealant to ensure proper orientation.

**NOTE:** If a high percentage agarose gel is cast it may be difficult to remove the casting comb from the tray. If resistance is felt when removing the casting comb, place the gel tray into the buffer tank and add buffer before removing the casting comb.

2. Pour the buffer into the tank. The level of buffer should cover the gel by at least 2-3mm (see *Table B: General specifications* for buffer volume guidelines). A fill line is located on the outside wall of the tank for easy reference. Too little buffer may cause the gel to dry out during the run, while excess buffer may slow DNA migration in the gel.
3. Load your samples on the gel. To load samples, angle the pipette tip into the well and slowly underlay sample into the well. Take care not to push pipette tip through the bottom of the sample well as this will result in significant well-to-well leakage of samples
4. Slide the cover onto the electrophoresis buffer tank, insuring complete connection of the tank's banana plugs to the female connector ends of the power cords that are fixed to the cover. The cover acts as a safety inter-lock to prevent accidental shock during operation.
5. Attach the power cords to the power supply. Turn on the power supply and set to the appropriate voltage level and begin electrophoresis.

Recommended running conditions are 5 volts/cm of inter-electrode distance, see *TABLE C* below.

**TABLE C:** Recommended run conditions

Model	Distance between electrodes	Max. Voltage
N2023-1410	22cm (110V)	200V
N2113-2510	34.3cm (171.5V)	200V
N2120-2510	34.3cm (171.5V)	200V
N2123-2510	34.3cm (171.5V)	200V
N2123-4010	49cm (recommended 200V)	200V

6. Run the gel for the appropriate amount of time for the specific sample being analysed or until the dye front has migrated to approximately 5mm from the end of the gel. Turn off power supply and detach the power cords. Slide the cover back to remove. Lift the gel tray out of the buffer tank and proceed with visualisation.
7. It is important to rinse the tank with distilled water after every use in order to keep it clean. It is recommended to allow the tank to air dry, rather than drying with a wipe or towel, to avoid damage to the electrode wires.

## VIII. VISUALISATION

Stain the gel with ethidium bromide solution or other suitable stain for agarose gels.

**Note:** Ethidium bromide is a mutagen and can cause serious skin and eye irritation. STARLAB recommends wearing its Microflex® Freeform™ gloves, and safety glasses, when handling.

After the staining process gel may be viewed using a transilluminator. The gel tray is made of UV transmissible acrylic and, therefore, may be placed directly on the transilluminator to view the gel.

## IX. TROUBLESHOOTING

Many factors may affect the quality of gel preparations. For example, preparation of gel and sample buffers, gel casting and tank assembly and/or run conditions. Reading and following the instructions in this user manual can solve most problems. Some of those most commonly experienced problems, along with suggestions for solving them, are listed below:

Problem: **COMBS WELL GET DAMAGED WHEN COMB IS REMOVED**

- Always allow the gel to solidify completely before moving the gel tray or the unit, and before removing the comb. Gently rock the comb back and forth to loosen the comb — the rocking helps to avoid suction as the comb is removed.
- If a high percentage agarose gel is cast it may be difficult to remove the casting comb, place the gel tray into the buffer tank and add running buffer before removing the casting comb.

Problem: **SAMPLES DO NOT BAND SHARPLY & APPEAR DIFFUSED IN THE GEL BEFORE RUNNING**

- Gels should be allowed to solidify completely. For standard agarose this would be approx. 30 mins.
- If STARLAB low-melt agarose is used it may be necessary to completely solidify gels at a cooler temperature in the refrigerator or cold room.
- Gels should be submerged in 3–5mm of buffer to avoid gel drying out, but excess buffer  $\geq$  5mm can cause decreased DNA mobility and band distortion.

Problem: **BANDS ARE NOT SHARP, CLEAR AND EVEN**

- Always follow the proper procedure for preparing the agarose according to the manufacturer's instructions.
- When preparing the agarose, be sure all the agarose powder is in solution before heating. Swirl occasionally while heating to melt and mix completely.

Problem: **GEL SMILING**

- Gel may be uneven. Be sure to pour gel on a level surface.
- Voltage may be too high. Lower voltage setting on power supply. See *TABLE C*.

Problem: **NO MIGRATION**

- No power is contacting the gel. Be sure all contacts are made between the power supply and the electrophoresis tank.
- Check buffer preparation to be sure the proper ionic strength buffer was prepared.

- Check integrity of electrodes to be sure no breakage has occurred. Any breaks in the electrode will cause slow migration or none at all. If a break is found contact STARLAB.
- To test the unit, fill the tank with running buffer and attach to power supply without the gel tray and gel in the tank. The platinum wire electrodes will produce small bubbles as the current passes through. If a complete circuit does not exist there will be few or no bubbles.

Problem: **SLOW MIGRATION**

- Check buffer preparation to be sure the proper ionic strength buffer was prepared.
- Voltage too low, increase voltage, see Instructions for proper running conditions.
- Check integrity of electrodes to be sure no breakage has occurred. Any breaks in the electrode will cause slow migration or none at all. If a break is found contact STARLAB.

Problem: **GEL SEEMS TO RUN SLOWER THAN USUAL RUNNING CONDITIONS**

- The volume of running buffer used to submerge the gel should only be 3-5mm over the gel surface. Gel should be completely submerged to avoid the gel from drying out, which can smear the bands and possibly melt the gel due to overheating.
- If excess buffer is added the mobility of the DNA decreases and band distortion may result.

Problem: **EXCESSIVE HEAT**

- Ionic strength of buffer is too high. Check buffer preparation to be sure the proper ionic strength buffer was prepared.
- Voltage is too high, decrease the voltage; see *TABLE C* for proper running conditions.

Problem: **ONE SIDE OF GEL RUNNING SLOWER**

- Gel may be uneven. Be sure to pour gel on a level surface.
- Check integrity of electrodes to be sure no breakage has occurred. Any breaks in the electrode may cause uneven migration or none at all. If a break is found contact STARLAB.

Problem: **NO BANDS VISIBLE**

- Be sure gel was placed in the electrophoresis tank in the proper orientation. If orientation or polarity is reversed, the samples will migrate off the gel.
- If DNA marker was used and is present on the gel, check the integrity of the sample.
- Shorten run time. Gel should be stopped when the dye front has migrated to approximately 5mm from the end of the gel.

## **X. COMB SPECIFICATIONS**

**NOTE:** Well volumes are calculated assuming a 0.5cm thick gel. Since the comb placement leaves 0.2 cm of agarose under the sample well, a 0.5cm thick gel will have a sample well that is 0.3cm deep (a 1cm thick gel would leave a sample well depth of 0.8cm).

Loading volume is calculated as:

**Thickness of Tooth x Width of Tooth x Depth of Sample (gel thickness – 0.2cm) x 75%**

<b>Combs for N2113-2510 (13 x 25cm)</b>				<b>Well volume (µl)</b>		
<b>STARLAB Comb Cat No.</b>	<b>Number of Teeth</b>	<b>Thickness (mm)</b>	<b>Width (mm)</b>	<b>0.5cm</b>	<b>0.75cm</b>	<b>1.0cm</b>
<b>N2113-2538</b>	8	1.0	13.9	31	57	83
<b>N2113-2541</b>	12	1.0	8.7	20	36	52
<b>N2113-2535</b>	16	1.0	6.1	14	25	37
<b>N2113-2537</b>	20	1.0	4.5	10	19	27
<b>N2113-2539</b>	24	1.0	3.5	8	14	21
<b>N2113-2067</b> Microtitre comb	14	1.0	7.2	16	30	43
<b>N2113-2528</b>	8	1.5	13.9	47	86	125
<b>N2113-2532</b>	12	1.5	8.7	29	54	78
<b>N2113-2527</b>	16	1.5	6.1	21	38	55
<b>N2113-2530</b>	20	1.5	4.5	15	28	41
<b>N2113-2531</b>	24	1.5	3.5	12	22	32
<b>N2113-2590</b> Wall Comb	1	1.5	-	-	-	-
<b>N2113-2065</b> Microtitre Comb	14	1.5	7.2	24	45	65
<b>N2113-2050</b> Preparative Comb	2	1.5	117.3/4.7	396/17	725/29	1056/42

<b>Combs for N2023-1410 (23.5 x 14cm)</b>				<b>Well volume (µl)</b>		
<b>N2123-2510 (23.5 x 25cm) &amp; N2123-4010 (23.5 x 40cm)</b>				<b>0.5cm</b>	<b>0.75cm</b>	<b>1.0cm</b>
<b>STARLAB Comb Cat No.</b>	<b>Number of Teeth</b>	<b>Thickness (mm)</b>	<b>Width (mm)</b>	<b>0.5cm</b>	<b>0.75cm</b>	<b>1.0cm</b>
<b>N2023-1470</b> Microtitre comb	25	1.0	7.2	16	30	43
<b>N2023-1486</b> Microtitre comb	50	1.0	3.0	7	12	18
<b>N2023-1464</b> Microtitre comb	25	1.5	7.2	24	45	65
<b>N2023-1485</b> Microtitre comb	50	1.5	3.0	10	19	27
<b>N2023-1490</b> Wall comb	1	1.5	-	-	-	-

**Combs for N2120-2510 (20 x 25cm)**

				Well volume (µl)		
				0.5cm	0.75cm	1.0cm
STARLAB Comb Cat No.	Number of Teeth	Thickness (mm)	Width (mm)			
<b>N2120-2538</b>	8	1.0	22.7	51	94	136
<b>N2120-2541</b>	12	1.0	14.5	33	60	87
<b>N2120-2535</b>	16	1.0	10.4	23	43	62
<b>N2120-2537</b>	20	1.0	8.0	18	33	48
<b>N2120-2539</b>	24	1.0	6.4	14	26	38
<b>N2120-2546</b>	28	1.0	5.2	12	21	31
<b>N2120-2547</b>	32	1.0	4.3	10	18	26
<b>N2120-2545</b>	36	1.0	3.7	8	15	22
<b>N2120-2553</b> Microtitre Comb	18	1.0	7.2	16	30	43
<b>N2120-2554</b> Microtitre Comb	21	1.0	7.2	16	30	43
<b>N2120-2555</b> Microtitre Comb	42	1.0	2.7	6	11	16
<b>N2120-2528</b>	8	1.5	22.7	77	140	204
<b>N2120-2532</b>	12	1.5	14.5	49	90	131
<b>N2120-2527</b>	16	1.5	10.4	35	64	94
<b>N2120-2530</b>	20	1.5	8.0	27	50	72
<b>N2120-2531</b>	24	1.5	6.4	22	40	58
<b>N2120-2543</b>	28	1.5	5.2	18	32	47
<b>N2120-2544</b>	32	1.5	4.3	15	27	39
<b>N2120-2542</b>	36	1.5	3.7	13	23	33
<b>N2120-2590</b> Wall Comb	1	1.5	-	-	-	-
<b>N2120-2549</b> Microtitre Comb	18	1.5	7.2	24	45	65
<b>N2120-2551</b> Microtitre Comb	21	1.5	7.2	24	45	65
<b>N2120-2552</b> Microtitre Comb	42	1.5	2.7	9	17	24
<b>N2120-2550</b> Preparative Comb	2	1.5	187/4.7	631/16	1157/29	1683/42

## XI. ACCESSORIES

### Accessories for N2023-1410 (23.5 x 14cm)

<b>N2023-1416</b>	StarPhoresis Gel Tray 23.5 x 14cm gasketed, 4 comb slots, UVT
<b>N2113-2588</b>	Replacement Gaskets for Gel Tray N2023-1416, Pack of 2

### Accessories for N2113-2510 (13 x 25cm)

<b>N2113-2516</b>	StarPhoresis Gel Tray 13 x 25cm gasketed, 6 comb slots, UVT
<b>N2113-2588</b>	Replacement Gaskets for Gel Tray N2113-2516, Pack of 2

### Accessories for N2120-2510 (20 x 25cm)

<b>N2120-2516</b>	StarPhoresis Gel Tray 20 x 25cm gasketed, 6 comb slots, UVT
<b>N2120-2588</b>	Replacement Gaskets for Gel Tray N2120-2516, Pack of 2

### Accessories for N2123-2510 (23.5 x 25cm)

<b>N2123-2516</b>	StarPhoresis Gel Tray 23.5 x 25cm gasketed, 14 comb slots, UVT
<b>N2123-4088</b>	Replacement Gaskets for Gel Tray N2123-2516 and N2123-4016, Pack of 2

### Accessories for N2123-4010 (23.5 x 40cm)

<b>N2123-4016</b>	StarPhoresis Gel Tray 23.5 x 40cm gasketed, 12 comb slots, UVT
<b>N2123-4088</b>	Replacement Gaskets for Gel Tray N2123-2516 and N2123-4016, Pack of 2

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**STARLAB (UK), Ltd**
**STARLAB S.r.l.**