PS2A200

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Indice

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Power Supply Safety – English

- If this equipment is used in a manner not specified by Hoefer, Inc. the protection provided by the equipment may be impaired.
- This instrument is designed for indoor laboratory use only.
- Only accessories and parts approved or supplied by Hoefer, Inc. may be used for operating, maintaining, and servicing this product.
- Warning! Because this instrument can develop sufficient voltage and current to produce a lethal shock, care must be exercised in its operation.
- This instrument is designed in accordance with the EN61010-1:2001 electrical safety standard. Nevertheless, it should be used only by properly trained operators. Read this entire manual before using the instrument and use only according to the instructions.
- The instrument must always be used with the earth lead of the power cord correctly grounded to earth at the mains outlet.
- Use only undamaged electrical wire and equipment specific for the voltages you will use. All equipment connected to high voltage should be in accordance with EN61010-1:2001.
- Keep the instrument as dry and clean as possible. Wipe regularly with a soft, damp cloth. Let the instrument dry completely before use.
- Do not operate the instrument in extreme humidity

(above 80%). Avoid condensation by letting the unit equilibrate to ambient temperature when taking the instrument from a colder to a warmer environment.

- To permit sufficient cooling, ensure that the vents of the instrument are not covered.
- Because this instrument has high-current capability, check the rating of any equipment attached. Do not exceed its rating when setting the power supply limits.

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```
Hoefer PS2A200 Power Supply
Serial No. test-EPS4
Software Ver. B0
Protocol # 1 Protocol 1
***** Run conditions *****
Timed Hrs
Const Volt
Step #1 0:05 Hrs
50 Volts
Step #2 0:30 Hrs
150 Volts

***** Run started *****
At total run time 0:00 Hrs
50Volts 246mA
12Watts
Step:1 0:05 Hrs
0:00 Hrs Elap

***** End Time Step *****
At total run time 0:05 Hrs
50Volts 246mA
12Watts
Step:1 0:05 Hrs
0:05 Hrs Elap

At total run time 0:10 Hrs
150Volts 732mA
110Watts
Step:2 0:30 Hrs
0:05 Hrs Elap

At total run time 0:30 Hrs
150Volts 727mA
109Watts
Step:2 0:30 Hrs
0:25 Hrs Elap

***** End Time Step *****
At total run time 0:35 Hrs
150Volts 726mA
109Watts
Step:2 0:30 Hrs
0:30 Hrs Elap

***** Run ended *****
At total run time 0:35 Hrs
```

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© 2012 Hoefer, Inc.