P5 Newtron XS

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ENGLISH

CONTENTS

1-	INTRODUCTION	5	10-	ELECTROMAGNETIC COMPATIBILITY	
				Electromagnetic emissions	16
2-	WARNINGS	5-6		Electromagnetic immunity	17
				Electromagnetic immunity /	
3	DESCRIPTION			radiofrequency mobile equipment	17-18
	Physical description	6-7		Recommended separation distances	18
	Technical description	7-8		Cable lengths	19
4	INSTALLATION/FIRST USE		11-	DISPOSAL AND RECYCLING	19
	Unpacking the device	8			
	Recommendations	8	12-	MANUFACTURER'S LIABILITY	19
	Installation	9			
	First use	9	13-	REGULATIONS	20
5-	IRRIGATION SOLUTIONS	10	14-	SYMBOLS	20
6-	USING THE DEVICE	10	15-	CUSTOMER RELATIONS	
7-	SHUTTING DOWN THE DEVICE	11		Manufacturer identification	112
				Subsidiaries	112-113
8-	ROUTINE MAINTENANCE/STERILIZATION				
	Maintaining the irrigation circuit/ irrigation tank	11			
	Handpiece maintenance	11			
	Tip maintenance	12-13			
	Device maintenance	13			
9-	MAINTENANCE				
	Fuse replacement	13-14			
	Pump replacement	14			
	Operation faults	14-15			

The list of accessories referred to in this manual is non exhaustive. Consult your local dealer for further details.

I - INTRODUCTION

Congratulations! You have just taken possession of your P5 NEWTRON XS®.

Designed by SATELEC®, the P5 NEWTRON XS is a piezoelectric ultrasonic generator that uses ultrasonic tips for dentistry.

In order to get the most out of this high-technology device and ensure it has a long working life, please read this manual carefully before using or maintaining the device.

Points preceded by the symbol  should paid particular attention.

II - WARNINGS

CAUTION:

United States Federal Law restricts the use of this device solely to qualified, trained and competent dental health practitioners or under their supervision.

To reduce the risks of accidents, the precautions stipulated below must be taken:

Device users:

- Use of the P5 NEWTRON XS is restricted solely to qualified, trained and competent dental health practitioners in the normal context of their work.

- If you have received this device by error, please contact the supplier so that it can be removed.

Interactions:

- The P5 NEWTRON XS must not be used if the patient and/or the operator has a cardiac stimulator (pacemaker) or any other active implant (e.g. a cochlear implant).
- The device complies with applicable electromagnetic compatibility standards. The user should nevertheless ensure that any potential electromagnetic interference does not cause an additional risk (presence of radiofrequency emitters, electronic devices, etc.).
- The device is not designed to withstand shocks delivered by an electric defibrillator.

Electrical connection:

- Your device must be connected to the electric power supply by a certified dental installation technician.
- To avoid any risk of electric shock, this device must be connected to an electric supply with a protective ground.
- The electric supply to which the device is connected must comply with the standards in force in your country.

Using the device:

- Do not use the device if it appears to be damaged or faulty.
- Turn the device off before unplugging the power cord.
- To unplug the power cord, grip the cord plug and hold the wall socket.
- Never fill the irrigation solution tank while it is placed on the device.

- Do not fill the irrigation solution tank beyond the level indicated by the border of the tank.
- When the device is not to be used for a long period of time, rinse thoroughly, then unplug the device from the electric supply.
- Do not attempt to move the device during use.

Environment:

- Do not cover the device or obstruct the ventilation vents.
- Do not immerse the device in liquid, and do not use it outdoors.
- Do not tilt the device at an angle greater than 5°.
- Do not place the device near a heat source.
- Make sure that the cords are not in the traffic path.
- The device should be stored in its original packaging in a safe place.
- The device is not designed for use in the presence of anesthetic gases or any other flammable gas.
- Do not expose the device to water vapor, or splashes.
- The device is not designed to work near ionizing radiation.
- Do not insert metal objects into the device (risk of electric shock, short-circuit or emission of hazardous substances).

Maintenance:

- The device must be rinsed before and after each use.
- To avoid any reaction between different irrigation solutions which may be used (crystallization, precipitation, etc.), the irrigation solution tank and irrigation circuit must be rinsed with distilled or demineralized water each time the irrigation solution is changed.

Accessories:

- The device can only accept accessories designed by SATELEC for the particular use for which they are intended.
- Use of accessories from other manufacturers is a potential hazard for you or your patients.

Repair:

- Do not repair or modify the device without prior authorization from SATELEC.
- In the case of a fault, contact the supplier of your device. Do not use unauthorized repairers, who might make your device dangerous for you and your patients.

If you have any doubt, contact an approved dealer or the SATELEC customer support department:

- www.acteongroup.com
- E-mail: satelec@acteongroup.com.

III - DESCRIPTION

3.1 PHYSICAL DESCRIPTION

The P5 NEWTRON XS is composed of five elements:

- a control unit (Fig. 1-1),
- a footswitch (Fig. 1-2),
- a cable with connector for handpiece (Fig. 1-3),
- a tank for irrigation solution (Fig. 1-4),
- an irrigation solution tank cover (Fig. 1-5).

The P5 NEWTRON XS control unit has a control dial on the upper side marked OI (ON / OFF) and marks ranging from 1 to 20 (power level).

Located at the back of the unit are:

- a power cord socket with ground pin (Fig. 2-1),
- a non-disconnectable footswitch cord (Fig. 2-2).

The handpiece connect to the scaler cord (non-disconnectable) located on the right-hand side (Fig. 1-3).

The handpiece connector can be placed on the holder (Fig. 3-4).

The irrigation flow adjustment and Flush/Prime button is located on the left-hand side (Fig. 2-3).

The device has ventilation vents on the lower side.

The device has four rubber feet to allow proper air flow.

The accessories: SUPRASSON NEWTRON handpiece and tips have their own individual instructions manuals.

3. 2 TECHNICAL DESCRIPTION

a) Control unit upper panel

Dial (Fig. 3-1) on the top is used to:

- Switch the machine on and off: OI (Fig. 3-2).
- Adjust operating power: 1 to 20.

The display is separated into four colored sections, each with five positions (fig. 3-2)

The "ON" indicator (Fig. 3-3) lights up with the color and number corresponding to the chosen power level:

- **Green: 1 to 9:** very low to low power, mainly used for periodontics.
- **Yellow: 6 to 14:** medium power, mainly used for endodontics.
- **Blue: 12 to 18:** high power, mainly used for scaling.
- **Orange: 17 to 20:** very high power, used mainly for prosthesis removal.

The technology of this device allows the practitioner to work at constant amplitude depending on the power selected with compensation according to the effort applied.

Note:

The tip usage power level is selected according to the directions in the TipBook™ and the color-coding system (CCS).

b) Control unit back panel

The P5 NEWTRON XS power receptacle with its grounding pin is connected to the power network using the disconnectable power cable.

The practitioner switches on the P5 NEWTRON XS using the ON/OFF footswitch.

Pressing the footswitch automatically switches on the ultrasonic handpiece and the irrigation function (if this is not set to "0").

The weight of the device and its anti-slip pad guarantee good stability.

c) Right-hand side of the control unit

The cord connecting the device and the handpiece connector governs:

- Fluid circulation for the spray.
- The electric connection between the device and the handpiece.

d) Left-hand side of the control unit

The knob (Fig. 2-3) governs:

- Switching off the irrigation function when turned completely counterclockwise.
- Irrigation flow rate adjustment: from "min" to "max".

FLUSH/PRIME function, when pressed. Press this button to start draining with a delay of 8 min (time required to empty a full tank); press again to stop the draining and delay.

Note:

- Knob in position *MIN* corresponds to a flow of around 5ml/min.
- Knob in position *MAX* corresponds to a flow of around 40ml/min.
- *FLUSH/PRIME* function corresponds to a flow of around 40ml/min.

e) Technical characteristics

Device name: P5 NEWTRON XS

Electrical power supply:

Voltage: 100VAC / 115VAC / 230VAC (factory set)

Frequency: 50Hz/60Hz

Rated power: 50VA @ 230VAC

Operating:

Intermittent service: 10min ON / 5min OFF

Output characteristics:

Vibration min. frequency: 28kHz

Irrigation flow rate: 0 to 40ml/min (nominal value)

Protection:

Electrical safety: Class 1, type BF

2 fuses (power receptacle):

5 mm x 20 mm / 315mAT @ 230 V rating 35A

5 mm x 20 mm / 400mATT @ 115 VAC rating 35A

5 mm x 20 mm / 400mATT @ 100 VAC rating 35A

1 FU3 fuse : CMS / 2AT

Overall dimensions:

Control unit:

Width: 158mm

Height: 88mm

Height with tank: 125mm

Depth: 290mm

Weight: 1700 g

Footswitch:

Width: 70mm

Height: 30mm

Depth: 95mm

Weight: 150g

Handpiece cable: 2.000mm

Tank capacity: 300ml

Temperatures:

Operating: +10°C to +40°C

Storage: -20°C to +70°C

Humidity:

Operating: 30% to 75%

Storage: 10% to 100% including condensation

IV - INSTALLATION / FIRST USE

4. 1 UNPACKING THE DEVICE

Upon receipt of the device, check for any damage caused in transit.

Contact your supplier if necessary.

4. 2 RECOMMENDATIONS

Have your device connected to the electric supply by a certified dental installation technician.

P5 NEWTRON XS electrical connections must be in compliance with current standards. A grounded power socket must be used.

4. 3 INSTALLATION

 **Important:**

Never rotate the handpiece connector on its cord, otherwise you may damage the device.

Do not place the P5 NEWTRON XS close to or on top of another device.

Do not place the power cord and the footswitch cord in a wire cover or in a cable gland.

- Place the control unit in position on a fixed, horizontal surface, or with a slope of no more than 5°.
- Check that the power control dial is on OI (stop).
- Connect the power cord to the power receptacle of the device.
- Connect the power cord to a socket with a ground pin.
- Position the footswitch so that it is easily accessible for your feet.
- Check that there is no humidity in the handpiece connections. If there is, dry it out by wiping, then blowing with the multifunction syringe.
- Connect the handpiece by simply plugging in the connector cord and lining up the electrical contacts.
- Place the handpiece on its rest.
- Install the irrigation solution tank on its support, after removing the protection cap on the device's receiving connector (Fig. 2-5).
- Place the cover on the irrigation tank.

4. 4 FIRST USE

Before using the P5 NEWTRON XS for the first time you are advised to carry out the following operations:

- Remove the tank from the device.
- Remove the cap from the tank.
- Fill the irrigation solution tank with disinfectant, without exceeding the level indicated by the border (Fig. 4-1).
- Place the cover on the irrigation solution tank.
- Replace the tank on the device.
- Connect the handpiece to the handpiece cord.
- Set the device to minimum power (1).
- Press the irrigation flow setting button (Fig. 2-3) to rinse the irrigation circuit on the device.
- Completely empty the disinfectant solution from the tank.
- Remove the tank from the device.
- Remove the cover from the irrigation solution tank.
- Fill the irrigation solution tank with distilled water or demineralized water.
- Place the cover on the irrigation solution tank.
- Replace the tank on the device.
- Press the irrigation flow setting button (Fig. 2-3) to rinse the irrigation circuit on the device.
- Completely empty the tank.
- Disconnect the handpiece.
- Remove the cover from the irrigation solution tank.
- Remove the irrigation solution tank.
- Perform maintenance/sterilization on all the material according to the instructions in chapter 8.

V - IRRIGATION SOLUTIONS

The P5 NEWTRON XS is not designed to administer drugs.

However it can perform irrigation circuit rinsing and operating field irrigation with medical products such as:

- Physiological saline.
- Hypochlorite.
- Chlorhexidine.
- Dakin®.
- Eludril®.
- E.D.T.A. Salvizol®.

VI - USING THE DEVICE

Important:

- **Do not disconnect the handpiece from the cord while the machine is plugged in and the footswitch is pressed.**
 - **For your safety and that of your patient, the P5 NEWTRON XS must not be used with accessories other than those supplied by SATELEC.**
 - **For correct tip vibration the tip must be correctly screwed in, and not forced past its stop point. To avoid tip self-blockage, it must always be removed after use.**
 - **If it is screwed too far, the tip or handpiece can fracture.**
 - **The irrigation tank cover must always be in place while the device is in use.**
 - **Check the integrity of the device and its accessories before and after each use in order to detect any problem.**
 - **If necessary, replace all faulty parts.**
- Always remove the irrigation tank from the device before filling it.
- The P5 NEWTRON XS is started as follows:
- Make sure that the device is connected correctly and properly insulated.
 - Check that there is no humidity in the handpiece connections or the handpiece cord connector. Otherwise, dry it out by wiping and blowing with a multifunction syringe.
 - Plug the handpiece into the handpiece cord connector.
 - Remove the tank from the device.
 - Remove the irrigation tank cover.
 - Fill the tank with irrigation solution without exceeding the level indicated by the border (Fig. 4-1).
 - Replace the cover on the irrigation tank.
 - Install the filled tank on its support.
 - Set the device to minimum power (1) and maximum flow.
 - Start up the irrigation current by pressing the footswitch.
 - When the irrigation solution reaches the end of the handpiece, release the footswitch.
 - Screw in the desired tip.
 - Tighten the tip with moderate force (using the wrench) to ensure optimal ultrasonic operation.
 - Set the device to the required power (according to color and number) using the control dial.
 - If necessary, refer to the TipBook for ideal SATELEC tip power settings.
 - Adjust the irrigation solution flow using the flow adjustment knob.
 - Press the footswitch to use the device.

VII - SHUTTING DOWN THE DEVICE

At the end of the working day, you must:

- Empty the irrigation solution out of the irrigation circuit and the irrigation solution tanks, by pressing the irrigation flow setting button (Fig. 2, item 3).
- Rinse the circuit and irrigation tank (see chapter IV, 4.4) with distilled or demineralized water.
- Disconnect the handpiece.
- Switch the device off (O).

VIII - ROUTINE MAINTENANCE/STERILIZATION

Important:

To maintain the sterility or asepsis of the accessories (handpiece, tips, etc.), store them in sealed pouches or leak-proof containers suitable for medical practice. The handpiece cord and tank are not sterilizable.

The following maintenance/sterilization instructions should be followed before using for the first time any of the accessories delivered with the device, and before each subsequent use.

Avoid using cleaning agents and disinfectants containing inflammable substances.

If this is not possible, make sure that all of the product has evaporated and that there is no combustible matter on the device or its accessories before it is switched on.

8.1 MAINTAINING THE IRRIGATION CIRCUIT / IRRIGATION TANK

In preparation for a long period of disuse of your P5 NEWTRON XS, the tank and irrigation circuit of the device must be rinsed (see 4.4).

After a period of disuse of your P5 NEWTRON XS, bacteria may develop on the inside of the irrigation solution tank and irrigation solution circuit.

You are therefore strongly advised to disinfect the irrigation solution tank and irrigation circuit of the device (see 4.4) before using it again.

8.2 SUPRASSON NEWTRON HANDPIECE MAINTENANCE

After each use, the whole of the handpiece irrigation circuit and the tip must be rinsed with distilled or demineralized water for 20 to 30 seconds.

For cleaning, disinfecting and sterilizing the handpiece, see its specific instructions manual.

8. 3 TIP MAINTENANCE

Important:

During sterilization, items made of different metals must not come into contact. Any contact would establish galvanic couples and cause local damage. To avoid this happening, place the objects in a sterilization pouch.

a) Lifetime

Tip shape and weight are the determining characteristics for obtaining maximal efficiency of the ultrasonic generator.

Careful monitoring of these two characteristics by the user will ensure that optimum performance of the device is maintained.

Consequently, you are strongly advised not to alter the structure of the tips by filing or twisting them or modifying them in any other way.

Similarly, tip aging leads to alteration of its characteristics through normal wear.

Always replace a tip that has been damaged by wear or by accidental impact (dropped, deformed, etc.).

Use the TIPCARD™ to check scaler tip wear. Replace routinely-used tips at least one to three times a year.

b) Pre-disinfection/Chemical Cleaning

To be carried out immediately after treatment, wearing thick gloves.

Immerse in an ultrasonic bath (enzymatic, quaternary ammonium base detergent / disinfectant solution), carefully respecting the concentration and contact time recommended by the manufacturer of the solution.

Use a solution that bears the CE marking or that is in compliance with any standard required by the national regulations.

Rinse under running water for at least 30 seconds.

c) Mechanical/Chemical Cleaning

To be carried out immediately after pre-disinfection / chemical cleaning, wearing thick gloves.

Scrub the products in a new bath of enzymatic, quaternary ammonium base detergent / disinfectant solution, for at least 30 seconds, and in any case until all remaining traces of contamination have been removed.

Rinse again under running water for at least 30 seconds.

d) Drying

Dry with a clean, non-woven single-use cloth, so as to remove any traces of liquid.

Keep in single-use sterilization sleeves or pouches in compliance with the specifications defined in standard EN ISO 11607-1 or in any equivalent standard required by the national regulations.

e) Sterilization

The tips must be sterilized individually in an autoclave, according to the following parameters:

- autoclave: type B, compliant with standard EN 13060;
- sterilization temperature: 134°C;
- sterilization stabilization time: 18 minutes;
- pressure: 2 Bars min.

f) Storage

Then, store the sterilized products in a dry place protected from dust.

Before re-use, if the integrity of the packaging is not conformable, recondition the re-sterilize in accordance with the defined protocol.

If contamination is visible inside the bag, destroy the product.

g) Disposal of the product

Dispose of the product in receptacles for waste materials of healthcare activities involving infectious risks.

CAUTION: Not all autoclaves can reach 134° C. Not all autoclaves draw a pre-vacuum. Please refer to your autoclave manufacturer for specific sterilization instructions.

8. 4 DEVICE MAINTENANCE

Important:

- Do not use the device if it appears to be faulty.
- Do not use an abrasive cleaning agent on the device.

a) Maintenance

The P5 NEWTRON XS control unit, the handpiece cord and the footswitch must be cleaned and disinfected on a daily basis (using alcohol or disinfectants, dental practice wipes, such as SEPTOL™ wipes).

b) Monitoring

The accessories and cords should be regularly checked for faulty insulation. Replace them if necessary.

It is important to keep the control unit ventilation vents clean in order to avoid abnormal heating.

Check that there is no humidity in the handpiece connector. If there is, dry it out by wiping, then blowing with the multifunction syringe.

Note:

In the case of a fault, you are advised to contact the supplier of your P5 NEWTRON XS. Do not use unauthorized repairers, who might make your device dangerous for you and your patients.

IX - MAINTENANCE

9. 1 FUSE REPLACEMENT

Your P5 NEWTRON XS is protected by two fuses located in the mains receptacle (Fig. 2-1).

The replacement procedure is as follows:

- Switch the device off (position O).
- Unplug the power cord from the electric supply.
- Unplug the power cord from the power receptacle (Fig. 2-1).
- Insert the tip of a flat screwdriver into the notch above the fuse drawer to release it (Fig. 2-4).
- Remove the old fuses.
- Replace with new fuses of the same type and rating.
- Push the fuse drawer back into its housing until it clicks to indicate the correct position.
- Plug the power cord into the receptacle (Fig. 2-1).
- Plug the power cord into the electric supply.

Note:

The device also has an internal fuse (ref. FU3 on the printed circuit) which is not accessible to the user.

Contact the SATELEC after-sales service (refer to chapter 2 - Repairs).

All the information necessary for the repair of faulty components is made available by SATELEC to the technical personnel of the SATELEC-approved dealer network who carry out the work.

9. 2 PUMP REPLACEMENT

**Important:**

- It is essential to replace a worn pump with a pump supplied by SATELEC.
- In case of difficulty, you can contact the SATELEC after-sales service (see Chapter 2 - Repairs).

The P5 NEWTRON XS pump is a "wearing" mechanical part.

To maintain your device in perfect working order, SATELEC advises you to change the pump every year.

Please refer to the accompanying documents supplied with the replacement pump.

9. 3 OPERATION FAULTS

Please refer to the following table.

Fault detected	Possible causes	Solutions
Device does not operate (LED indicator off)	Power cord not plugged in properly.	- Check the power socket. - Return to SATELEC after-sales.
	Power power switch in position O.	- Place the power switch in position I.
	Power receptacle fuse(s) blown.	Replace the power receptacle fuses.
	Internal fuse blown.	Return to SATELEC after-sales.
No spray	Irrigation solution tank empty.	Fill the irrigation solution tank.
	Faulty pump.	Replace the pump.
	Irrigation deactivated	Return to SATELEC after-sales
	Blocked irrigation circuit	Return to SATELEC after-sales.
	Blocked tip or file.	Unblock the tip or file.
	Wrong choice of tip.	Check the tip.
Power fault/ low vibrations	Tip worn or distorted.	Replace the tip.
	Incorrect power adjustment.	Refer to the TipBook for the necessary information.
	Incorrect use: incorrect angle of attack or inadequate pressure on the tooth..	Refer to the TipBook for the necessary information.
	Liquid or moisture present between handpiece and cord.	Dry the electrical contacts thoroughly.
No ultrasonic vibration	Incorrect tip tightening.	Tighten the tip with the wrench.
	Faulty connector contact.	Clean the connector contacts.
	Handpiece cord wire cut.	Return to SATELEC after-sales to replace the cord.
Leak between the handpiece cord and the connector on the device	Worn handpiece 1.15x1 seal.	Replace the seal (kit F12304).
Leak of liquid under the device	Tube rupture.	Return to SATELEC after-sales.
Leak of liquid at the tank connector level.	Worn seals.	- Replace the two seals (E15121). - Return to SATELEC after-sales.

X - ELECTROMAGNETIC COMPATIBILITY

Important:

The power cable, the handpiece cable and the footswitch cable must be kept apart.

The P5 NEWTRON XS requires special precautions to be taken with regard to electromagnetic compatibility.

It must be installed and prepared for use as described in chapter 4.

Certain types of mobile telecommunication devices such as mobile telephones are likely to interfere with the P5 NEWTRON XS.

The recommended separation distances in this paragraph must therefore be complied with.

The P5 NEWTRON XS must not be used near or on top of another device.

If this cannot be avoided, its operation under the conditions of use must be checked beforehand.

The use of accessories other than those specified or sold by SATELEC as replacement parts may have the consequence of increasing the emissions or decreasing the immunity of the P5 NEWTRON XS.

10. 1 ELECTROMAGNETIC EMISSIONS

The P5 NEWTRON XS is intended for use in the electromagnetic environment specified in the table below.

The user and/or installer must ensure that the P5 NEWTRON XS is used in such an environment.

Emission test	Compliance	Electromagnetic environment - guidance
RF emission - CISPR 11.	Group 1	P5 NEWTRON XS uses RF energy for internal operation. Therefore, its radiofrequency emissions are very low and are not likely to cause any interference in nearby equipment.
RF emission - CISPR 11.	Class B	P5 NEWTRON XS is suitable for use in all establishments, including domestic establishments and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
Harmonic current emission IEC 61000-3-2.	Class A	
Voltage fluctuation and flicker IEC 61000-3-3.	Complies	

10. 2 ELECTROMAGNETIC IMMUNITY

The P5 NEWTRON XS is intended for use in the electromagnetic environment specified in the table below.

The user and/or installer must ensure that the device is used in such an electromagnetic environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2.	± 6KV contact ± 8KV air	± 6KV contact ± 8KV air	Floors must be wood, concrete, cement or tiled. If floors are covered with synthetic material (carpet, etc.), the relative humidity must be at least 30%.
Electrical fast transients IEC 61000-4-4.	± 2KV for power supply lines	± 2KV for power supply lines	Power quality should be that of a typical commercial or hospital environment (hospital, clinic).
Surges IEC 61000-4-5.	± 1KV differential mode ± 2KV common mode	± 1KV differential mode ± 2KV common mode	Power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations IEC 61000-4-11.	<5% U _T (>95% dip in U _T) for 0.5 cycles. 40% U _T (60% dip in U _T) for 5 cycles 70% U _T (30% dip in U _T) for 25 cycles <5% U _T (>95% dip in U _T) for 250 cycles	<5% U _T (>95% dip in U _T) for 0.5 cycles. 40% U _T (60% dip in U _T) for 5 cycles 70% U _T (30% dip in U _T) for 25 cycles <5% U _T (>95% dip in U _T) for 250 cycles	Power quality should be that of a typical commercial or hospital environment. If the use of the P5 NEWTRON XS requires continued operation during a power cut, it is recommended that the product be powered from a separate power supply (UPS, etc.).

10. 3 ELECTROMAGNETIC IMMUNITY / MOBILE RADIOFREQUENCY EQUIPMENT

The P5 NEWTRON XS is intended for use in the electromagnetic environment specified in the table below.

The user and/or installer must ensure that the P5 NEWTRON XS is used in such an electromagnetic environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Portable and mobile radiofrequency communications devices must not be used near the P5 NEWTRON XS (including its cables) at a distance less than that recommended and calculated according to the frequency and power of the emitter.			
Conducted disturbance, radiofrequency fields. IEC61000-4-6	3 V/m 150 KHz to 80 MHz	3 V/m	Recommended separation distance: $d = 1.2 \sqrt{P}$
Radiated radiofrequency electromagnetic field. IEC61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = 1.2 \sqrt{P}$ 80MHz à 800MHz. $d = 2.3 \sqrt{P}$ 800MHz à 2.5GHz. Where P is the maximum power rating of the emitter in watts (W) according to the manufacturer's specifications and d is the recommended minimum separation distance in meters (m).
The electromagnetic field strengths of fixed radiofrequency emitters, as determined by an electromagnetic environment measurement (a), must be less than the compliance level in each frequency range (b). Interference may occur near equipment marked with the symbol below: 			

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These specifications may not be applicable in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and persons.

(a): The electromagnetic field strengths of fixed radiofrequency emitters, such as base stations for mobile telephones (cellular/cordless), mobile radios, amateur radios, AM/FM radio broadcasts and TV broadcasts cannot be determined exactly by theory. To assess the electromagnetic environment due to fixed radiofrequency emitters, an electromagnetic environment measurement must be made. If the measured radiofrequency field strength in the immediate environment where the product is used exceeds the compliance level specified above, the performance of the product must be tested to verify whether it conforms to the specifications. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the product.

(b): In the 150 kHz to 80 MHz frequency range, the electromagnetic field strengths must be less than 3 V/m.

10. 4 RECOMMENDED SEPARATION DISTANCES

The P5 NEWTRON XS is intended for use in an electromagnetic environment in which radiated radiofrequency disturbances are controlled.

The P5 NEWTRON XS user and/or installer can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile radiofrequency communications equipment (emitters) and the P5 NEWTRON XS, according to the maximum output power of the equipment, as recommended in the table below.

Rated max. power of the emitter (W)	Separation distance in meters (m) according to emitter frequency		
	150 KHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2,5 GHz $d = 2.3 \sqrt{P}$
0,01	0.12 m	0.12 m	0.23 m
0,1	0.38 m	0.38 m	0.73 m
1	1.2 m	1.2 m	2.3 m
10	3.8 m	3.8 m	7.3 m
100	12 m	12 m	23 m

For emitters rated at max. power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the emitter, where P is the max. power rating of the emitter in watts (W) according the manufacturer.

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These specifications may not be applicable in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and persons.

10. 5 CABLE LENGTHS

Cables and accessories	Maximum length	Complies with:
Handpiece cord Control footswitch Power cord	<3m	RF emission, CISPR 1 - Class B/Group 1 Harmonic current emission: IEC61000-3-2 Voltage fluctuation: IEC61000-3-3 Immunity to electrostatic discharge: IEC61000-4-2 Immunity to electrical fast transients/bursts: IEC61000-4-4 Immunity to surges: IEC61000-4-5 Immunity to voltage dips, short interruptions and voltage variations: IEC61000-4-11 Immunity to conducted disturbances induced by radiofrequency fields: IEC61000-4-6 Immunity to radiated radiofrequency electromagnetic fields: IEC61000-4-3

XI - DISPOSAL AND RECYCLING

As electrical and electronic equipment, the device must be disposed of according to a specialized procedure for collection, pick-up and recycling or destruction (in particular on the European market, with reference to Directive 2002/96/EC of 23/01/2003).

When your device reaches the end of its life, we consequently recommend that you contact your dental equipment dealer (or, failing this, the nearest ACTEON GROUP office, the list of which is given in chapter 15), for information on how to proceed.

XII - LIABILITY

The manufacturer is not liable if:

- the manufacturer's installation recommendations have not been followed (supply voltage, electromagnetic environment, etc.),
- repairs have been performed by persons not authorized by the manufacturer,
- the device has been used in an electrical installation which does not comply with current standards,
- the device has been used in a way which is not stipulated in this manual,
- accessories (tips, handpiece, pump, etc.) other than those supplied by SATELEC have been used,
- the instructions in this document have not been followed.

Note:

The manufacturer reserves the right to modify the device and/or the user manual without notice.

XII - REGULATIONS

This medical device is classified as class IIa according to European Directive 93/42/EEC.

This equipment is manufactured in compliance with the current IEC 60601-1 standard.

This equipment has been designed and manufactured according to an ISO 13485-certified quality assurance system.

XIV - SYMBOLS

	Alternating current
	BF device
	Warning, please refer to the accompanying documentation
	Not to be disposed of as household waste
	CE Marking
	Irrigation

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