GE Healthcare

# ÄKTApilot

Safety Handbook  
Sicherheitshandbuch  
Manual de seguridad  
Manual de segurança  
Manuel de sécurité  
Manuale sulla sicurezza  
Säkerhetshandbok

ÄKTA™

## Important user information

---

All users must read this entire manual to fully understand the safe use of ÄKTApilot.

### WARNING!

---

The WARNING! sign highlights instructions that must be followed to avoid personal injury. Do not proceed until all stated conditions are clearly understood and met.

### CAUTION!

---

The CAUTION! sign highlights instructions that must be followed to avoid damage to the product or other equipment. Do not proceed until all stated conditions are met and clearly understood.

### Note

---

The Note sign is used to indicate information important for trouble-free and optimal use of the product.

### CE Certifying

---

This product meets all requirements of applicable CE-directives. A copy of the corresponding Declaration of Conformity is available on request.

The **CE** symbol and corresponding declaration of conformity, is valid for the instrument when it is:

- used as a stand-alone unit, or
- connected to other CE-marked GE Healthcare instruments, or
- connected to other products recommended or described in this manual, and
- used in the same state as it was delivered from GE Healthcare except for alterations described in this manual.

### Recycling

---

This symbol indicates that the waste of electrical and electronic equipment must not be disposed as unsorted municipal waste and must be collected separately. Please contact  an authorized representative of the manufacturer for information concerning the decommissioning of equipment.

### WARNING!

---

This is a Class A product. In a domestic environment, it might cause radio interference, in which case the user might be required to take appropriate measures.

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---

## 1.1 Introduction

The purpose of this manual is to describe safety precautions and emergency and disposal procedures for ÄKTA™pilot as well as instructions for recycling. A presentation of all safety labels is also included.

**IMPORTANT!** All users must read this entire manual to fully understand the safe use of ÄKTApilot.

**IMPORTANT!** The separation unit and the PC should be installed and prepared by GE Healthcare personnel.

**IMPORTANT!** To avoid risk of injury, operation and user maintenance should be performed by properly trained personnel only, and in accordance with the instructions.

### 1.1.1 ÄKTApilot system

ÄKTApilot is a high performance, automated liquid chromatography system. The system consists of the separation unit and a Windowsbased computer with UNICORN control system version 4.11 or higher.

The system is designed for process development, process scale-up and scale-down, and small scale production.

Fig 1-1. ÄKTApilot system.

### 1.1.2 Associated documentation

The *ÄKTApilot User Manual* describes the operation of ÄKTApilot. It provides instructions for preparing the system and performing a run using UNICORN control system.

The *ÄKTApilot Instrument Handbook* contains technical information, maintenance schedules and instructions for troubleshooting and user maintenance.

UNICORN control system is described in three manuals:

- *Getting Started*
- *User Reference Manual (2 pcs.)*
- *Administration and Technical Manual*

### 1.1.3 Warning symbols

**WARNING!** The Warning symbol highlights instructions that must be strictly followed to avoid personal injury. Be sure not to proceed until the instructions are clearly understood and all stated conditions are met.

### 1.1.4 Caution notices

**CAUTION!** The Caution sign highlights instructions that must be followed to avoid damage to the product or other equipment. Be sure not to proceed until the instructions are clearly understood and all stated conditions are met.

### 1.1.5 Note signs

*Note:* Indicates information important for trouble-free and optimal use of the product.

## 1.2 Safety precautions

- Read the user manuals before using ÄKTApilot.
- The separation unit and the PC should be installed and prepared by GE Healthcare personnel.
- This system is designed for indoor use only.
- To avoid risk of injury, operation and user maintenance should be performed by properly trained personnel only, and in accordance with the instructions.
- Any equipment connected to the system should meet the requirements of the EN 61 010-1 or other international safety standard. Within EU, connected equipment must be CE-labelled.
- Do not operate ÄKTApilot hardware or UNICORN control system in any other way than described in ÄKTApilot and UNICORN user manuals.

1 Safety instructions  
1.2 Safety precautions

**1.2.1 Warnings when using the system**

**WARNING!** The instrument must not be opened by the user. It contains electrical circuits which can give a lethal electric shock.

**WARNING!** Always disconnect the mains supply, unless stated otherwise, before attempting to replace any item on the system during maintenance.

**WARNING!** The instrument must always be used with the protective earth lead of the power cord correctly grounded to earth at the mains outlet.

**WARNING! HIGH PRESSURE.** The system operates under high pressure. Protective glasses must be worn at all times.

**WARNING! HAZARDOUS CHEMICALS.** When using hazardous chemicals, all suitable protective measures, such as wearing protective glasses, must be taken.

**WARNING! OVER-PRESSURE.** Use columns that withstand expected pressures. If not, the columns might rupture, resulting in injury.

**WARNING! OVER-PRESSURE.** The flow rate may under no circumstances exceed the specified column maximum flow rate. High flows might affect the packed medium, causing the pressure to exceed the specified column maximum pressure.

**WARNING!** Use ONLY tubings supplied by Amersham Biosciences to ensure that the pressure specifications of the tubings are fulfilled.

**WARNING! OVER-PRESSURE.** Never block the outlet tubing with, for instance, stop plugs, since this will create over-pressure and might result in injury.

**WARNING! HAZARDOUS CHEMICALS.** Incorrectly fitted tubing might loosen, causing a jet of liquid to spray out. This is especially dangerous if hazardous chemicals are in use. Connect the tubing by first inserting the tubing fully, then tightening the connector fingertight. Do not tighten ETFE tubing further as this will damage the end of the tubing. If a leakage occurs, immediately switch off the system before rectifying the leakage.

**WARNING! HAZARDOUS CHEMICALS.** When using hazardous chemicals, ensure that the entire system tubing has been flushed thoroughly with distilled water before service and maintenance.

**WARNING! CORROSIVE CHEMICALS.** NaOH is corrosive and therefore dangerous to health. Avoid spillage and wear protective glasses.

**WARNING!** This is a Class A product. In a domestic environment, this product may cause radio interference, in which case the user may be required to take adequate measures.

**WARNING! HEAVY OBJECT!** The separation unit weighs 114 kg. Use suitable lifting equipment when moving the system.

**WARNING! HIGH VOLTAGE.** The mains cable should only be connected by authorized service personnel. Faulty connection might result in live system parts that can give a lethal electric shock.

**WARNING! HAZARDOUS CHEMICALS.** Sulfuric acid included in the test solutions is corrosive and therefore dangerous to health. Avoid spillage and use eye protection.

**WARNING! NO SERVICEABLE PARTS INSIDE.** Do not open covers. Service and planned maintenance should be performed by personnel authorized by GE Healthcare.

1 Safety instructions  
1.2 Safety precautions

**1.2.2 Cautions when using the system**

**CAUTION!** Always make sure that the chromatography media, columns and system components are compatible with sodium hydroxide at the concentration, time, and temperatures used.

**CAUTION!** Never leave the pH electrode in the cell holder when the system is not used, since this might cause the glass membrane of the electrode to dry out. Remove the pH electrode from the cell holder and fit the end cover filled with a 1:1 mixture of pH 4 buffer and 2 M  $\text{KNO}_3$ . Do NOT store in water only.

**CAUTION!** Do not allow solutions containing dissolved salts, proteins or other solid solutes to dry out in the UV cell. Do not allow particles to enter the cell as damage to the cell might occur.

**CAUTION!** Do not touch the tips of the optical fiber with anything other than lens paper.

**CAUTION!** Only use spare parts supplied or specified by GE Healthcare

**CAUTION!** Check valves have precision matched components and should only be disassembled by a trained service engineer. If the problem cannot be corrected, the check valve should be replaced completely.

**CAUTION!** Handle the check valves with care when they have been removed from the pump heads to prevent loss of any internal components.

**CAUTION!** Do not disassemble the pump head unless there is good reason to believe that there is an internal leakage. Always make sure that sufficient spare components are available before attempting to replace a spare part.

**CAUTION! REPLACING SPARE PARTS!** Read the instructions carefully. For example, some individual parts of the pump head can be assembled incorrectly. Check the orientation of each part before continuing with the next instruction.

**CAUTION!** NaOH might cause the outer protection cylinder on the air trap to rupture. Avoid spillage.

**CAUTION!** Over-tightening might damage threads. Use a torque wrench to tighten the components.

- 1 Safety instructions
- 1.3 Labels and indicators

### 1.3 Labels and indicators

#### 1.3.1 Rating labels

The identification and rating labels (see Fig 1-3) are located on the lower part of the rear panel on the separation unit.

Voltage	100–240 V~
Frequency	50–60 Hz
Power, max	800 VA

System No. 123456	Code No. 18-1164-95
----------------------	------------------------

Fig 1-2. Layout of identification and rating labels.

#### 1.3.2 Safety label

The safety label (see Fig 1-3) is attached to the lower part of the rear panel on the separation unit. The safety label warns the user of risk for personal injury. Do not proceed until the instructions are clearly understood and all stated conditions are met.



Fig 1-3. Safety label on ÄKTApilot.

### 1.3.3 Controls and indicators

The separation unit is equipped with the following indicators, switch and push-buttons:

Indicator/Switch/ Push-button	Color	Description
<b>POWER</b> (switch)	–	Switches on/off power to the system. <b>Located on the left panel.</b>
<b>POWER</b> (indicator)	Green	Flashes rapidly for a few seconds during the internal communication test in the separation unit when switching on power to the unit. Flashes slowly when the communication test is finished. Steady light when UNICORN is connected to the separation unit.
<b>CONTINUE/RUN</b> (push-button/ indicator)	Green	Pressing this button with the system in Pause mode, causes the pump to continue at the set flow rate and gradient values. The button is lit in Run mode.
<b>PAUSE</b> (push-button/ indicator)	Yellow	Pressing this button, stops the pump but retains the set flow rate and gradient values. All inlet and outlet valves are closed. The button is lit in Pause mode.
Alarm buzzer	–	Indicates an alarm in UNICORN. <b>Located behind the front panel.</b>

**Table 1-1.** Indicators, switches and push-buttons on the separation unit.

Using the **CONTINUE/RUN** and **PAUSE** push-buttons to set the system in "Service mode" (the valve membranes move to a position between open and closed):

- 1 Disconnect the UniNet-1 cable between the separation unit and the computer. Wait for the **POWER** indicator to start flashing slowly (after approx. 1 min).
- 2 Press **CONTINUE/RUN** and **PAUSE** simultaneously for 3 seconds.

- 1 Safety instructions
- 1.4 Emergency procedures

## **1.4 Emergency procedures**

### **1.4.1 Emergency shutdown**

In a situation where there is a risk of injury, switch off the mains power by turning the mains power switch (located on the left side panel) to position **0**.

### **1.4.2 Power failure routine**

In the event of a power failure, the run is interrupted in an undefined state.

### **1.4.3 Restart procedure**

- 1 Switch on the separation unit by turning the mains power switch (located on the left side panel) to position **1**.

The **POWER** indicator on the separation unit flashes rapidly for a few seconds during the internal communication test. After the test, the indicator flashes slowly.

- 2 Connect UNICORN to the separation unit. When the communication is established, the **POWER** indicator on the separation unit should show a steady light.

## 1.5 Recycling

ÄKTApilot system contains the following materials:

- Stainless steel
- Aluminum
- Glass
- Titanium
- Elgiloy
- Ruby/sapphire
- Plastic:
  - Polypropylene (PP)
  - Polyetheretherketone (PEEK)
  - Polyethylene (PE)
  - Polytetrafluoroethylene (PTFE)
  - Perfluoroethylenepropylene copolymer (FEP)
  - Ethylenetetrafluoroethylene (ETFE)
  - Chlorotrifluoroethylene (CTFE)
- Perfluororubber (FFKM)
- Electronic components

When taking ÄKTApilot out of service, the different materials must be separated and recycled according to local regulations.

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[www.gehealthcare.com/akta](http://www.gehealthcare.com/akta)

GE Healthcare Bio-Sciences AB  
Björkgatan 30  
751 84 Uppsala  
Sweden

GE, imagination at work and GE monogram are trademarks of General Electric Company.

Drop Design, UNICORN and ÄKTA are trademarks of GE Healthcare companies.

All third party trademarks are the property of their respective owners.

All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare which supplies them. A copy of these terms and conditions is available on request. Contact your local GE Healthcare representative for the most current information

© 2003- 2007 General Electric Company—All rights reserved.  
First published Feb. 2003

GE Healthcare Europe GmbH  
Munzinger Strasse 5, D-79111 Freiburg, Germany

GE Healthcare UK Ltd  
Amersham Place, Little Chalfont, Buckinghamshire, HP7 9NA, UK

GE Healthcare Bio-Sciences Corp  
800 Centennial Avenue, P.O. Box 1327, Piscataway, NJ 08855-1327, USA

GE Healthcare Bio-Sciences KK  
Sanken Bldg. 3-25-1, Hyakunincho, Shinjuku-ku, Tokyo 169-0073, Japan

Asia Pacific T +85 65 62751830 F +85 65 62751829 • Australasia T +61 2 8820 8299 F +61 2 8820 8200 • Austria T 01 /57606 1613 F 01 /57606 1614 • Belgium T 0800 73 890 F 02 416 8206  
• Canada T 1 800 463 5800 F 1 800 567 1008 • Central & East Europe T +43 1 972 720 F +43 1 972 722 750 • Denmark T +45 70 25 24 50 F +45 45 16 2424 • Eire T 1 800 709992 F +44 1494  
542010 • Finland & Baltics T +358 9 512 3940 F +358 9 512 39439 • France T 01 69 35 67 00 F 01 69 41 98 77 • Germany T 0800 9080 711 F 0800 9080 712 • Greater China T +852 2100  
6300 F +852 2100 6338 • Italy T 02 26001 320 F 02 26001 399 • Japan T 81 3 5331 9336 F 81 3 5331 9370 • Korea T 82 2 6201 3700 F 82 2 6201 3803 • Latin America T +55 11 3933 7300  
F +55 11 3933 7304 • Middle East & Africa T +30 210 96 00 687 F +30 210 96 00 693 • Netherlands T 0800-82 82 82 1 F 0800-82 82 82 4 • Norway T +47 815 65 777 F +47 815 65 666 •  
Portugal T 21 417 7035 F 21 417 3184 • Russia, CIS & NIS T +7 495 956 5177 F +7 495 956 5176 • Spain T 902 11 72 65 F 935 94 49 65 • Sweden T 018 612 1900 F 018 612 1910 • Switzerland  
T 0848 8028 10 F 0848 8028 11 • UK T 0800 515 313 F 0800 616 927 • USA T +1 800 526 3593 F +1 877 295 8102