Piezotome

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ENGLISH

C O N T E N T S

I - INTRODUCTION	5	VI - MAINTENANCE	
		Generator unit and accessories	18
II - WARNINGS	5-6	Maintenance of handpieces and tips	18
		Taking the pump apart	18-19
III - DESCRIPTION		Changing the fuses	19
Items delivered	6-7	Fan	19
Components	7-8	Inspection	19
Technical specification	8-12		
IV - INSTALLATION		VII - MAINTENANCE	
Recommendations	12	Troubleshooting	19-22
Unpacking the device	12-13	Replacement of faulty components	22
Installation	13-14		
V - USE		VIII - DISPOSAL AND RECYCLING	22
Using for the first time	14		
Booting	14-15	IX - LIABILITY	22-23
Switching the active connection	15		
PIEZOTOME function	15-16	X - MANUFACTURER	23
NEWTRON function	16-17		
Changing tips	17	XI - SYMBOLS	23
		AFTER-SALES SERVICE	131

I - INTRODUCTION

Congratulations! You have just taken possession of your PIEZOTOME™.

Designed by SATELEC®, the PIEZOTOME is an ultrasound generator for use in dental surgery.

- The PIEZOTOME function can be used for surgical procedures, including osteotomy, osteoplasty, periodontal surgery and implantation.
- NEWTRON function can be used for mechanical ultrasound treatment in prophylaxis, periodontics or endodontics.

This device is fitted with two handpiece cord connectors. Depending on the application, the practitioner may connect a PIEZOTOME handpiece or a NEWTRON handpiece.

For greater efficiency, two handpieces of the same type can be connected at the same time.

The PIEZOTOME function offers four utilization modes at pre-set ultrasound power settings.

In NEWTRON function, the power of the ultrasound can be finely adjusted by the user.

In both cases, the user regulates the flow rate of the irrigation fluid.

The practitioner controls the device using a keyboard and a multi-function footswitch.

The irrigation fluid flow rate and the ultrasound power are monitored on a screen.

For enhanced efficacy, the last settings are memorized by the machine.

The PIEZOTOME makes it possible to work in optimal sterile conditions.

In order to get the most out of this high-technology device for as long a time as possible—as well as cutting down the risk of accident—please read this manual carefully before using the device for the first time.

II - WARNINGS

2. 1 INSTALLATION PRECAUTIONS

- The device should only be connected to the main electric supply by a qualified dental installation engineer who distributes SATELEC products.
- The electrical supply to which the device is connected should comply with current standards. Use of an earthed mains plug is compulsory.
- On reception of the device, check for damage caused in transit. If there is any problem, contact the SATELEC After-Sales Service (see Section 10.2).
- Follow the installation instructions in the correct order (see Chapter IV).

2. 2 ENVIRONMENTAL PRECAUTIONS

- Do not cover the device or obstruct the air vents.
- Do not immerse in liquid.
- Operating temperature: +10 to +40°C at a maximum of 95% relative humidity without condensation.
- Storage temperature: -20 to +70°C at a maximum of 95% relative humidity without condensation.
- Do not use the device in an inflammable atmosphere or in the presence of explosive gas (e.g. volatile anesthetics).

- Position the device on a flat surface, tilted to an angle of no more than 10°.
- Make sure that the cords do not get in people's way.
- The PIEZOTOME complies with current Electromagnetic Compatibility standards (CEI 601-1-2). However, the user should ensure that any possible electromagnetic interference does not present an additional risk of interfering with or damaging the device or any other electronic devices in the vicinity.

The device should not be used if either the practitioner or the patient is carrying an active implanted device (e.g. a pacemaker). If in doubt, seek advice from a qualified source.

- The mechanical ultrasound waves produced by the tips may damage prostheses.
- The device and its accessories are not designed to withstand shocks delivered by an electric defibrillator.
- The device should be stored in its original packaging in an appropriate place, where it does not pose any danger to people.

2. 3 USAGE PRECAUTIONS

- The PIEZOTOME should only ever be used by trained dental practitioners in the context of their regular activity.
- Before using the device, check the patient's condition (for fragile tissues, incomplete anesthesia, if the patient is in a particularly emotional or nervous state, etc.).
- The device is designed to be used exclusively with SATELEC accessories.
- Do not use the device if it or any of its accessories appear damaged or faulty.
- Only ever use a sterile irrigation fluid such as normal saline. The irrigation line is not designed for the administration of drugs.

- In the event of any problem, you are strongly recommended to contact the SATELEC After-Sales Service (see Section 10.2) rather than having it repaired by an unauthorized person who may make the device dangerous for patient and/or practitioner.
- Turn the device off before unplugging from the mains.
- Before use, make sure that any disinfectant products (which might contain inflammable substances) have completely evaporated.

2. 4 MAINTENANCE

- Follow maintenance instructions to the letter (see Chapter VI).
- The user must never destroy the device or any of its accessories (see Chapter VIII).

III - DESCRIPTION

3. 1 ITEMS DELIVERED

The components which make up the PIEZOTOME are delivered in separate packages:

- One "Generator Unit and Accessories" package
- One or more "Handpiece Packs" and/or kit(s) (according to the option chosen)

A detailed list of the items delivered is given below.

Refer to Figure 1 for the purposes of identification.

a) Generator Unit and Accessories

- A. 1 PIEZOTOME Generator Unit
- B. 1 mains cord
- C. 1 multifunction footswitch
- D. 2 handpiece supports
- E. 5 sterile irrigation lines with integrated cassette
- F. 1 bracket + 3 hooks
- G. 1 pump dismantling wrench
- H. 1 manual

b) Packs and instrument kits (depending on the option chosen)

PIEZOTOME BONE SURGERY Pack

- I. 1 sterilization box
- J. 1 PIEZOTOME handpiece and its cord
- K. 1 set of BONE SURGERY tips (BS1, BS2L, BS2R, BS4, BS5, BS6)
- L. 1 torque wrench
- M. 1 manual
- P. 5 fixing clips

PIEZOTOME SINUS LIFT Pack

- I. 1 sterilization box
- J. 1 PIEZOTOME handpiece and its cord
- K. 1 set of SINUS LIFT tips (SL1, SL2, SL3, SL4, SL5)
- L. 1 torque wrench + 1 flat universal wrench
- M. 1 manual
- P. 5 fixing clips

No-tip PIEZOTOME Pack

- I. 1 sterilization box
- J. 1 PIEZOTOME handpiece and its cord
- L. 1 torque wrench
- M. 1 notice
- P. 5 fixing clips

NEWTRON Periodontics Pack

- I. 1 sterilization box
- N. 1 NEWTRON handpiece
- O. 1 NEWTRON handpiece cord

- K. 1 set of Perio tips (H1, H2L, H2R, H3, H4L, H4R)
- L. 1 torque wrench
- M. 1 notice
- P. 5 fixing clips

NEWTRON Retro Pack

- I. 1 sterilization box
- N. 1 NEWTRON handpiece
- O. 1 NEWTRON handpiece cord
- K. 1 set of Retro tips (S12-70D, S12-90ND, S13LD, S13RD)
- L. 1 torque wrench
- M. 1 notice
- P. 5 fixing clips

No-tip NEWTRON Pack

- I. 1 sterilization box
- N. 1 NEWTRON handpiece
- O. 1 NEWTRON handpiece cord
- L. 1 torque wrench
- M. 1 notice
- P. 5 fixing clips

3. 2 COMPONENTS

a) PIEZOTOME Unit: Front panel

Refer to Figure 2 to identify the items listed below.

1. Control screen (for details, see Section 3.2.f)
2. Keyboard (for details, see Section 3.2.g)
3. Left handpiece cord connector
4. Right handpiece cord connector
5. Left pump unit
6. Right pump unit

b) PIEZOTOME Unit: Back panel

Refer to Figure 3 to identify the items listed below.

1. Bracket support
2. Bracket

3. Fuse box
4. Mains cable socket
5. Air vents
6. Power supply switch
7. Footswitch connector

c) PIEZOTOME Unit: Pumps

Refer to Figure 4 to identify the items listed below.

1. Cover (open)
2. Irrigation line cassette compartment
3. Pump cartridges

d) Multifunction footswitch

Refer to Figure 6 to identify the items listed below.

1. Change footswitch:
 - mode (in PIEZOTOME function)
 - range (in NEWTRON function)
2. Switch-pedal: starts/stops the output of ultrasound
3. Switch-button: starts/stops irrigation
4. Switch-button: Right connection active/Left connection active
5. Flush irrigation circuit button
6. Connector cord for the PIEZOTOME Unit

e) Handpieces

Refer to the User's manual delivered with the equipment.

f) Control screen

Refer to Figure 7 to identify the items listed below.

1. LCD display (shows the fluid flow rate or ultrasound power)
2. Irrigation On light
3. Light indicating the type of information being displayed on the LCD
 - On: à fluid flow rate (in ml/mn) displayed
 - Off: à ultrasound power displayed

4. Lights indicating mode selected in PIEZOTOME functioning (always off in NEWTRON function)

g) Keyboard

Refer to Figure 8 to identify the items listed below.

1. Left connection active light
2. Right connection active light
3. Activate left connection key
4. Activate right connection key
5. Select ultrasound power range key (NEWTRON function)
6. Fine ultrasound power tuning key (within the selected range) (NEWTRON function)
7. Select keys for PIEZOTOME function
8. Set irrigation flow rate keys
9. Switch-key: Starts/stops irrigation
10. Flush irrigation circuit key

h) Taking the pump apart

Refer to Figure 9 to identify the items listed below.

1. Pump dismantling wrench
2. Pump cover hinge
3. Pump cover
4. Pump cartridge

3. 3 TECHNICAL SPECIFICATIONS

a) Standards

This equipment complies with the applicable regulations and standards.

European Directive no. 93/42/CEE (Class IIa Medical Device).

Norms CEI 60601-1, CEI 601-1-2 and related norms.

PIEZOTOME needs special precautions regarding electromagnetic compatibility.

It must be installed and put into service as described in the "Installation" chapter.

Certain types of mobile telecommunication devices such as mobile telephones are likely to interfere with PIEZOTOME.

The separation distances recommended in the "Electromagnetic Compatibility" chapter must be complied with.

PIEZOTOME should not be used adjacent to or stacked with other equipment.

If adjacent or stacked use is necessary, the Piezotome should be observed to verify normal operation in the configuration in which it will be used.

The use of accessories other than those specified or sold by Satelec as replacement parts may increase emissions or reduce the immunity of PIEZOTOME.

Electromagnetic emissions

PIEZOTOME is intended for use in the electromagnetic environment specified in the table below.

The user and/or the installer must ensure that PIEZOTOME is used in such an environment.

Emission test	Compliance	Electromagnetic environment - guidance
RF emission CISPR 11.	Group 1	PIEZOTOME uses RF energy only for its internal function. Therefore, its radiofrequency emissions are very low and are not likely to cause any interference in nearby equipment. PIEZOTOME is suitable for use in all establishments, including domestic establishments and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
RF emission CISPR 11.	Class B	
Harmonic current emission IEC 61000-3-2.	Class A	
Voltage fluctuation and flicker IEC 61000-3-3.	Complies	

Electromagnetic immunity

PIEZOTOME is intended for use in the electromagnetic environment specified in the table below.

The user and/or installer must ensure that PIEZOTOME is used in such an electromagnetic environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2.	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors must be wood, concrete, cement or tiled. If the floors are covered with synthetic material, the relative humidity must be at least 30%.
Electrical fast transients IEC 61000-4-4.	± 2 kV for power supply lines	± 2 kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surges IEC 61000-4-5.	± 1 kV differential mode ± 2 kV common mode	± 1 kV differential mode ± 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations IEC 61000-4-11.	<5% U _T (>95% dip in U _T) for 0.5 cycles 40% U _T (60% dip in U _T) for 5 cycles 70% U _T (30% dip in U _T) for 25 cycles <5% U _T (>95% dip in U _T) for 250 cycles	<5% U _T (>95% dip in U _T) for 0.5 cycles 40% U _T (60% dip in U _T) for 5 cycles 70% U _T (30% dip in U _T) for 25 cycles <5% U _T (>95% dip in U _T) for 250 cycles	Mains power quality should be that of a typical commercial or hospital environment. If the use of Piezotome requires continued operation during power mains interruptions, it is recommended that the product be powered from a separate power supply (UPS, etc.).

Electromagnetic immunity/mobile radiofrequency equipment

PIEZOTOME is intended for use in the electromagnetic environment specified in the table below. The user and/or installer must ensure that PIEZOTOME is used in such an electromagnetic environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Portable and mobile radiofrequency communications devices must not be used near PIEZOTOME (including its cables) at a distance less than that recommended and calculated according to the frequency and power of the emitter.			
Conducted disturbance, radiofrequency fields. IEC 61000-4-6	3 V/m 150 KHz to 80 MHz	3 V/m	Recommended separation distance: $d = 1.2 \sqrt{P}$
Radiated radiofrequency electromagnetic field. IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ 80 MHz to 800 MHz Where P is the maximum power rating of the emitter in watts (W) according to the manufacturer's specification and d is the recommended minimum separation distance in metres (m).
The electromagnetic field strengths of fixed radiofrequency emitters, as determined by an electromagnetic environment measurement (a), must be less than the compliance level in each frequency range (b). Interference may occur near equipment marked with the symbol below: 			

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These specifications may not be applicable in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and persons.

(a): The electromagnetic field strengths of fixed emitters, such as base stations for mobile telephones (cellular/cordless), mobile radios, amateur radio, AM/FM radio broadcasts and TV broadcasts cannot be determined exactly by theory. To assess the electromagnetic environment due to fixed radiofrequency emitters, an electromagnetic environment measurement must be made. If the measured radiofrequency field strength in the immediate environment where the product is used exceeds the radiofrequency compliance level specified above, the performance of the product must be tested to verify whether it conforms to the specification. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the product.

(b): In the 150 kHz to 80 MHz frequency range, the electromagnetic field strengths must be less than 3 V/m.

Recommended separation distances

PIEZOTOME is intended for use in an electromagnetic environment in which radiated radiofrequency disturbances are controlled.

The PIEZOTOME user and/or installer can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile radiofrequency communications equipment (transmitters) and PIEZOTOME, according to the maximum output power of the equipment, as recommended in the table below.

Rated max. power of the transmitter (W)	Separation distance in metres (m) according to transmitter frequency		
	150 KHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$
0.01	0.12 m	0.12 m	0.23 m
0.1	0.38 m	0.38 m	0.73 m
1	1.2 m	1.2 m	2.3 m
10	3.8 m	3.8 m	7.3 m
100	12 m	12 m	23 m
For transmitters rated at max. power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the max. power rating of the emitter in watts (W) according to the manufacturer.			

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These specifications may not be applicable in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and persons.

PIEZOTOME cable lengths

Cables and accessories	Maximum length	Complies with
Handpiece cord Control footswitch	2 m	RF emission, CISPR 1 - Class B/Group 1 Harmonic current emission - IEC 61000-3-2 Voltage fluctuation - IEC 61000-3-3 Immunity to electrostatic discharge - IEC 61000-4-2 Immunity to electrical fast transients/bursts - IEC 61000-4-4 Immunity to surges - IEC 61000-4-5 Immunity to Voltage dips, short interruptions and voltage variations - IEC 1000-4-11 Immunity to conducted disturbances induced by radiofrequency fields - IEC 61000-4-6 Immunity to radiated radiofrequency electromagnetic fields - IEC 61000-4-3

b) PIEZOTOME Unit

Dimensions/weight

Width: 350 mm (13.79 in)

Height: 110 mm (4.33 in) without bracket
490 mm (19.31 in) with bracket

Depth: 240 mm (9.46 in)

Weight: 2.6 kg (5.73 lb) (without any accessories)

Electrical specifications

Class: Class I, type BF

Operating mode: Intermittent 10 min/5 min

Power supply: 115 V AC, 50/60 Hz
230 V AC, 50/60 Hz

Power rating: 55 VA

Fuses: 230 V: 2 x 315 mA T (5 x 20 mm)
115 V: 2 x 630 mA T (5 x 20 mm)

Internal protector fuses: 2 x 1.5 AT (not accessible to the user)

Stored memory retention: 10 years at +25°C

Peristaltic pumps

Flow rate in irrigation mode: 5-80 ml/min $\pm$ 10%

Flow rate in flush mode: 80 ml/min $\pm$ 10%

Each motor is enclosed in a pump housing mounted on rubber shock absorbers.

Ultrasound

Output frequency: 28-36 KHz

c) Multifunction footswitch

Width: 185 mm (7.29 in)

Height: 40 mm (1.58 in)

Depth: 145 mm (5.71 in)

Number of footswitches: 2

Number of sealed buttons: 3

Connecting cord: about 2,000 mm (78.80 in)

d) Handpiece supports

Sterilizables

e) Irrigation line

Single-use sterile irrigation line with integrated pump cartridge

f) Bracket

Height: 473 mm (18.64 in)

Number of attachment points: 3

Number of hooks: 3

Can accommodate: Bags (not provided)

Bottles (not provided)

g) Cords

Earthed mains cord:

Length: 2,000 mm (78.80 in)

PIEZOTOME cord:

Length: 2000 mm (78,80 in)

Ultrasound cord:

Length: 2000 mm (78,80 in)

h) PIEZOTOME and NEWTRON Packs

 Refer to the User's manual delivered with the equipment.

NB: A handpiece compatibility system makes it impossible to screw a PIEZOTOME surgery tip (Bone Surgery, Sinus Lift or ligament cutting) to a NEWTRON handpiece or a CCS tip (periodontics, endodontics, etc.) to a PIEZOTOME handpiece.

IV - INSTALLATION

4. 1 RECOMMENDATIONS

See Section 2.1 "Installation precautions".

4. 2 UNPACKING THE DEVICE

On reception of the device, check for any damage caused in transit.

If you observe a problem:

- 1) Keep the original packaging.
- 2) Note the device's serial number.

- 3) Contact the SATELEC After-Sales Service (see Section 10.2).

4. 3 INSTALLATION

It is important to follow the instructions below in order.

a) Unit connections

- 1) Position the Unit (Figure 1, #A) in easy reach on a flat surface which is tilted to an angle of no more than 10°.
- 2) Insert the two rods of the bracket (Figure 3, #2) into the bracket supports (Figure 3, #1) so that the bracket is suspended above the Unit.
- 3) Plug the multifunction footswitch cord connector (Figure 5, #6) into the Unit (Figure 3, #7).
- 4) Place the footswitch on the ground so it is easy to operate with the foot.
- 5) Plug the mains cord (Figure 1, #B) into the power supply socket (Figure 3, #4).
- 6) Plug the PIEZOTOME handpiece and its cord (Figure 1, #J) into one of the unit cord connector (Figure 2, #3 or #4), taking care to align the red marks on each of the two components.
- 7) If the NEWTRON function is chosen, plug the NEWTRON cord (Figure 1, #O) into the unoccupied cord connector (Figure 2, #3 or #4).
- 8) Connect the NEWTRON handpiece (Figure 1, #N) to its cord (Figure 1, #O).

 Refer to the User's manual delivered with the equipment.

Nota:

- In the same way, it is possible to connect a second PIEZOTOME handpiece and its cord (Figure 1, #J) or a second NEWTRON cord to the unoccupied cord connector.

- When it is switched on, the device automatically identifies which type of cord is connected (PIEZOTOME or NEWTRON) (see Section 5.2.c).

b) Fitting the tip to the handpiece

- 1) Screw a tip (Figure 1, #K) into the handpiece (Figure 1, #J ou #N) using the wrench (Figure 1, #L).

 Refer to the User's manual delivered with the equipment.

- 2) Lay the handpiece on its support (Figure 1, #D).

NB: A handpiece compatibility system makes it impossible to screw a PIEZOTOME surgery tip (Bone Surgery, Sinus Lift or ligament cutting) to a NEWTRON handpiece or a CCS tip (periodontics, endodontics, etc.) to a PIEZOTOME handpiece.

c) Installing the irrigation line

The following instructions describe how to install an irrigation line. You will have to repeat the process to install a second irrigation line.

- 1) Open the sterile pouch and remove the irrigation line.
- 2) Connect the end of the silicone tube to the end of the cord connector.
- 3) Lift the pump cover on the same side of the device (Figure 4, #1).
- 4) Introduce the pump cartridge (Figure 4, #3) into its special compartment (Figure 4, #2).
- 5) Close the pump cover (Figure 4, #1), taking care not to trap any of the tubes.
- 6) Push the pump cartridge pin into the cap of the irrigation solution bottle/bag and open the air vent (if there is one).
- 7) If there is only one bottle/bag: hang the hook from the frame's central attachment point (Figure 1, #F-a).

- If there are two bottles/bags: hang the hooks from the lateral attachment points (Figure 1, #F-b).
- 8) Hang the bottle(s)/bag(s) off the bracket (Figure 1, #F).
 - 9) Position the fixing clips of the irrigation line (provided with the unit) on the handpiece cord (Figure 5) and push gently.
- d) Switching on**
- 1) Check that the power supply switch (Figure 3, #6) is Off (O).
 - 2) Plug the mains cord into the power supply.
 - 3) Switch the power supply switch (Figure 3, #6) On (I).
- 1) The 6 indicator lights on the control screen (Figure 7, #2, #3 and #4) and the 2 active connection indicator lights on the keyboard (Figure 8, #1 and #2) will light up.
 - 2) The LCD display on the control screen (Figure 7, #1) will show ".8.8" (the points and segments should always be lit up).
 - 3) The background of the LCD display (Figure 7, #1) goes from green to yellow to blue to orange to white. Check that all the colors appear in the right order.
 - 4) The 6 control screen indicator lights (Figure 7, #2, #3 and #4) and the 2 keyboard active connection indicator lights (Figure 8, #1 and #2) will go out.

V - USE

5. 1 USING FOR THE FIRST TIME

Accessories and the unit must always be sterilized before the device is used for the first time (see Section 6.1).

5. 2 BOOTING

Every time it is switched on, the device automatically boots up, a process which lasts a few seconds. The device should not be used until this process is complete. Booting up involves a series of different steps, as described below:

a) Step 1: Testing

This step checks that the LCD display is working properly. The user should check that it proceeds correctly:

b) Step 2: Default active connection

The device contains a single ultrasound generator with two outlets (Left and Right) (Figure 2, #3 and #4).

Only one of these can be used at a time. The connection being used is referred to as the "active connection".

- 1) The device detects the type of cord present at each connector (Left and Right).
- 2) It automatically decides which connection will be active on the following basis:
 - a PIEZOTOME connection takes priority over a NEWTRON connection.
 - the right connection takes priority over the left connection.
- 3) The light corresponding to the active connection lights up (Figure 8, #1 or #2).

NB: In the course of use, the operator can of course switch from one connection to the other (see Section 5.3).

c) Step 3: Switching to PIEZOTOME function or NEWTRON function

The device identifies the type of connection which is active (PIEZOTOME or NEWTRON).

- If the active connection has been identified as PIEZOTOME, the device automatically goes into PIEZOTOME function:

The PIEZOTOME block of the keyboard is activated (Figure 8, #7).

The NEWTRON block of the keyboard is inactivated (Figure 8, #5 and #6).

- If the active connection has been identified as NEWTRON, the device automatically goes into NEWTRON function:

The NEWTRON block of the keyboard is activated (Figure 8, #5 and #6).

The PIEZOTOME block of the keyboard is inactivated (Figure 8, #7).

d) Step 4: Loading settings

In this step, the device reloads the settings (ultrasound power and irrigation flow rate) stored at the end of the last session.

The settings are displayed on the control screen:

- In PIEZOTOME function:

A light indicates the current mode (Figure 7, #4).

The fluid flow rate is shown on the LCD display (Figure 7, #1).

- In NEWTRON function:

The background color of the LCD display indicates the current range (Figure 7, #1).

The fluid flow rate is shown on the LCD display (Figure 7, #1).

NB: the first time the device is started up, factory default settings are used.

e) Step 5: Activation of irrigation

At the end of the booting up process, the device automatically starts irrigation. The corresponding indicator (Figure 7, #2) lights up.

5. 3 SWITCHING THE ACTIVE CONNECTION (LEFT/RIGHT)

The user can switch from one connection to the other (Left/Right) during use:

- From the keyboard:

The left connection key (Figure 8, #3) activates the left connection.

The right connection key (Figure 8, #4) activates the right connection.

- Using the footswitch switch-button (Figure 6, #4).

Every time the button is stepped on, the active connection (Left/Right) switches.

The device immediately identifies the type of cord at the active connector (PIEZOTOME or NEWTRON) and accordingly goes into either PIEZOTOME function or NEWTRON function; the settings associated with the identified connection (side + type) are loaded (see Section 5.2.d).

NB: the keys on the keyboard and the footswitch button will be inactive if:

- the device is emitting ultrasound.
- there is no cord plugged into the connector.

5. 4 PIEZOTOME FUNCTION

a) Mode selection

In PIEZOTOME function, the PIEZOTOME block on the keyboard is active (Figure 8, #7).

The device has 4 operating modes with pre-set ultrasound power levels.

Mode	Density	Cortical	Bone
1	Very High	Thick	
2	High	Thick	
3	Medium	Thin	Dense Spongy
4	Low	Thin or none	Spongy

The practitioner has to choose the mode which is most appropriate to the procedure.

The mode can be selected either:

- From the keyboard (Figure 8, #7): The user should press the key corresponding to the mode to be selected.
- Using the left footswitch (Figure 6, #1): Every time the left footswitch is stepped on, the device goes on to the next mode (1, 2, 3, 4, 1...)

An indicator light on the control screen confirms the selection (Figure 7, #4).

The fluid flow rate associated with the selected mode is shown on the LCD display (Figure 7, #1).

NB: the keys on the keyboard and the footswitch button will be inactive if the device is emitting ultrasound.

b) Irrigation settings

In PIEZOTOME function, the control screen LCD display (Figure 7, #1) still shows the fluid flow rate in ml/mn (on a white background). An indicator light (Figure 7, #3) confirms the rate.

The user can adjust the fluid flow rate using the +/- keys of the irrigation block on the keyboard (Figure 8, #8).

The rate can be varied from 10-80 ml/mn in steps of 5 ml/mn.

Irrigation can be activated either:

- Using the keyboard switch-key (Figure 8, #9)
- Using the footswitch switch-button (Figure 6, #3)

If irrigation is On, the indicator (Figure 7, #2) lights up.

If irrigation is Off, the indicator (Figure 7, #2) goes out.

NB:

- the pump only operates while ultrasound is being emitted.
- the keys on the keyboard and the footswitch button are inactivated as long as the device is emitting ultrasound.

c) Turning the ultrasound on

Ultrasound is turned On or Off using the footswitch (Figure 6, #2). Every time the footswitch is stepped on, the generator is switched (On/Off).

- If irrigation is On, the indicator (Figure 7, #2) lights up and the pump switches on at the same time as the generator.
- If irrigation is Off, the indicator (Figure 7, #2) goes out and the pump does not switch on at the same time as the generator.

5. 5 NEWTRON FUNCTION

a) Power range

Range selection

In NEWTRON function, the NEWTRON block on the keyboard is active (Figure 8, #5 and #6).

The device has 4 ultrasound power ranges. The practitioner has to choose the range which is most appropriate to the procedure.

The range can be selected:

- From the keyboard (Figure 8, #5): The user should press the key corresponding to the range to be selected.
- Using the footswitch button (Figure 6, #1): Every time the button is stepped on, the device goes on to the next range (G, Y, B, O, G, ...).

Range	Ultrasound Power	Procedure
Green	Low	Periodontics mainly
Yellow	Medium	Endodontics mainly
Blue	High	Prophylaxis mainly
Orange	Very High	Specific treatment modalities

The background color of the LCD display (Figure 7, #1) confirms the selected range.

NB: The keyboard keys and the footswitch button will be inactive if the device is emitting ultrasound.

Fine adjustment

Within each ultrasound range, the user can fine-tune the power.

In NEWTRON function, the fine-tuning of the ultrasound within the selected range is shown by default on the LCD display (Figure 7, #1) (see Section 5.2.d).

The indicator light (Figure 7, #3) being out confirms that the LCD display is showing the ultrasound power.

The fine-tuning includes 10 levels which can be selected using the +/- keys of the NEWTRON block on the keyboard (Figure 8, #6).

b) Modulating irrigation

In NEWTRON function, the LCD display (Figure 7, #1) shows the ultrasound level by default. The screen's background color corresponds to the range selected (green, yellow, blue or orange).

The indicator light (Figure 7, #3) being out confirms that the LCD display is showing the ultrasound power.

The user can nevertheless modulate the irrigation flow rate using the +/- keys of the irrigation block on the keyboard (Figure 8, #8).

The rate can be varied from 5-40 ml/min in steps of 1 ml/min.

During adjustment, the LCD display (Figure 7, #1) shows the current fluid flow rate.

The indicator light (Figure 7, #3) being on confirms that the LCD display is showing the flow rate.

Three seconds after the flow rate has been set, the display switches back to the ultrasound level.

Irrigation is turned on by:

- Using the keyboard switch-key (Figure 8, #9)
- Using the footswitch switch-button (Figure 6, #3)

If irrigation is On, the indicator light (Figure 7, #2) lights up.

If irrigation is Off, the indicator light (Figure 7, #2) goes out.

NB:

- the pump only turns on when ultrasound is being emitted.
- the keyboard keys and the footswitch button are inactive if the device is emitting ultrasound.

c) Turning the ultrasound on

The ultrasound is turned On or Off using the footswitch (Figure 6, #2). Every time the footswitch is stepped on, the generator switches.

- If irrigation is On (and the indicator light Figure 7, #2) is lit up, the pump turns on at the same time as the generator.
- If irrigation is not On (and the indicator light Figure 7, #2) is out, the pump does not turn on at the same time as the generator.

5. 6 CHANGING TIPS

Refer to the User's manual delivered with the equipment.

NB: A handpiece compatibility system makes it impossible to screw a PIEZOTOME surgery tip (Bone Surgery, Sinus Lift or ligament cutting) to a NEWTRON handpiece or a CCS tip (periodontics, endodontics, etc.) to a PIEZOTOME handpiece.

VI - MAINTENANCE

Before every procedure, check that all the instruments are sterile and install a new sterile irrigation line with the handpiece cord.

6. 1 UNIT AND ACCESSORIES

(not including handpieces and tips)

A preliminary maintenance procedure should be carried out on the unit and its accessories every time the device is to be used. This maintenance procedure can be broken down into 3 phases:

- Cleaning (unit and accessories)
- Decontamination (unit and accessories)
- Sterilization (accessories)

Recommended maintenance procedure:

- 1) Clean and disinfect the unit and accessories using a suitable product or with special disinfectant wipes for dental applications. Avoid using cleaning agents and disinfectants containing inflammable substances.
- 2) Remove any fluid that might have penetrated into openings in the course of decontamination.
- 3) Put the decontaminated accessories into the sterilization box (Figure 1, #1).
Components made of different types of metal should not be allowed to come into contact with one another during sterilization (because of the risk of damage if electrolytic couples are formed).
- 4) Wrap the sterilization box in cloth or put it into a sterilization pouch.
- 5) Sterilize in a Type B autoclave (18 minutes steam-sterilization at 134°C and a pressure of 2 bar (200hPa).

6. 2 MAINTENANCE OF HANDPIECES AND TIPS

 Refer to the User's Manual for the specific piece of equipment.

During this procedure, tips should be taken off the handpieces.

6. 3 TAKING THE PUMP APART

(basic maintenance)

PIEZOTOME pumps can be taken apart to facilitate cleaning of internal components.

Proceed as follows:

- 1) Set the power supply switch (Figure 3, #6) to "0".
- 2) Pull the mains cord out of its socket (Figure 3, #4).
- 3) Position the lugs of the pump wrench (Figure 9, #1) under the claws of the pump cover hinge (Figure 9, #2).
- 4) Push the pump wrench towards the pump until it can advance no further.
- 5) Release the pump cover/pump cartridge compartment unit (Figure 9, #3 and #4) by gently lifting it upwards.
- 6) Remove the pump wrench (Figure , #1).
- 7) Remove the pump cartridge compartment (Figure 9, #4) from the pump cover (Figure 9, #3).
- 8) Clean the pump and components using a suitable product, making sure that no fluids run down the motor axle (getting fluid in the motor may damage the device).
- 9) Put the pump cartridge compartment (Figure 9, #4) and the cover (Figure 9, #3) back together.
- 10) Replace the unit on the pump, taking care to introduce each of the three metal rods (on the pump side) into the pump cartridge compartment guides. A click will be heard when they are in the correct position.

- 11) Put the mains cord back into its socket (Figure 3, #4).

6. 4 CHANGING THE FUSES

- 1) Set the power supply switch (Figure 3, #6) to "0".
- 2) Pull the mains cord out of its socket (Figure 3, #4).
- 3) Introduce the point of a flat screwdriver into the notch above the fuse drawer to release it (Figure 3, #3).
- 4) Remove old fuses.
- 5) Replace with new fuses of the same type and rating.
- 6) Push the fuse drawer back into its compartment until a click is heard to indicate that it is in the correct position.
- 7) Put the mains cord back into its socket (Figure 3, #4).

NB: the device also contains two internal fuses which the user cannot access. Contact the SATELEC After-Sales Service (see Section 10.2).

6. 5 FAN

Regularly check the fan for dust.

Proceed as follows:

- 1) Set the power supply switch (Figure 3, #6) to "0".
- 2) Pull the mains cord out of its socket (Figure 3, #4).
- 3) Clean off the fan by gently blowing through the air vents (Figure 3, #5).
- 4) Put the mains cord back into its socket (Figure 3, #4).

6. 6 INSPECTION

The accessories and cords should be regularly checked for faulty insulation. Replace, if necessary.

If any problem is observed, you are recommended to contact the SATELEC After-Sales Service (see Section 10.2) rather than some other repair service.

VII - MAINTENANCE

7. 1 TROUBLESHOOTING

NB: the user should never attempt to repair or modify the device without prior permission from SATELEC.

FAULT	POSSIBLE CAUSE	REMEDY
The device does not turn on	Mains cord disconnected or faulty	Check and, if necessary, replace mains cord
	Fuses damaged	Change the fuses (see Section 6.4) NB: the device also contains two internal fuses which the user cannot access. Contact the SATELEC After-Sales Service (see Section 10.2)
	Mains supply socket faulty	Contact an electrician
	Power supply switch (Figure 3, #6) set to Off (O)	Set the power supply switch (Figure 3, #6) to On (I)
Ultrasound is working but there is no irrigation	Irrigation not on	Switch irrigation on using the keyboard key-switch (Figure 8, #9) or the footswitch switch-button (Figure 6, #3)
	Irrigation line not connected to the handpiece	Connect the irrigation line (see Section 4.3.c)
	Flow rate incorrectly set	Adjust the flow rate setting (see Section 5.4.b or 5.5.b)
	Handpiece and/or tip clogged up	Disconnect the handpiece and check to see if fluid is coming out of the cord: - If not, contact the SATELEC After-Sales Service (see Section 10.2). - If it is, unscrew the tip from the handpiece, blow air into the tip and handpiece. If the problem persists, contact the SATELEC After-Sales Service of (see Section 10.2)
	Wrong tip or the selected ultrasound power is too weak	Change the tip or the mode (in Piezotome function) or power range (In Newtron function) (see Section 5.4.a or 5.5.a)
No atomization	Tip does not screw in properly	Screw the tip back in again Refer to the manual provided with the equipment
	Tip deformed or worn	Replace the tip Refer to the manual provided with the equipment
	Handpiece and/or tip not being used properly	Refer to the manual provided with the equipment
	Poor contact between the handpiece and its connector	Push the handpiece in without turning it
	Fluid present in the NEWTRON handpiece connector	Dry out the connector. If there is a black residue around the connector pins and/or the handpiece terminals, contact the SATELEC After-Sales Service (see Section 10.2)
	Fluid flow rate poorly adjusted	Adjust the flow rate setting (see Section 5.4.b or 5.5.b)

FAULT	POSSIBLE CAUSE	REMEDY
The control screen is lighting up but neither the ultrasound nor the irrigation systems are working	Footswitch disconnected or faulty	Check that the footswitch is properly connected Contact the After-Sales Service of SATELEC (see Section 10.2)
Ultrasound too strong or too weak	Tip not screwed in properly	Screw the tip back in. Refer to the manual provided with the equipment
	Tip deformed or worn	Replace the tip. Refer to the manual provided with the equipment
	Handpiece and/or tip not being used properly	Refer to the manual provided with the equipment
	Fluid present in the NEWTRON handpiece connector	Dry out the connector. If there is a black residue around the connector pins and/or the handpiece terminals, contact the SATELEC After-Sales Service (see Section 10.2)
No ultrasound output	Tip not screwed in properly	Screw the tip back in. Refer to the manual provided with the equipment
	Handpiece faulty	Replace the handpiece. Refer to the manual provided with the equipment
	Oxidation and/or a black residue around the NEWTRON connector pins	Clean the contacts, replace the cord or contact the SATELEC After-Sales Service (see Section 10.2)
	Cord cut	Replace the cord
	Handpiece not suitable	Refer to the manual et provided with the equipment
	Fluid present in the NEWTRON handpiece connector	Dry out the connector. If there is a black residue around the connector pins and/or the handpiece terminals, contact the SATELEC After-Sales Service (see Section 10.2)
Fluid leakage in the NEWTRON handpiece connector	Small water-proof washer in the NEWTRON handpiece worn out	Replace the washer. Refer to the manual provided with the equipment
Fluid leakage at the end of the tip	Disposable pump cartridge damaged	Change the pump cartridge (see Section 4.3.c)
TR OR TO DISPLAYED and no ultrasound output or irrigation	Internal problem with the device's functioning	Turn the device off and then turn it on again. If the problem persists, contact the SATELEC After-Sales Service (see Section 10.2)

FAULT	POSSIBLE CAUSE	REMEDY
Mode (in PIEZOTOME function) or power range (In NEWTRON function) cannot be changed	The device is emitting ultrasound	Stop ultrasound by stepping on the footswitch (Figure 6, #2)
	There is no cord plugged into the connector	Check that the handpiece cord is correctly connected to the connector
The irrigation circuit cannot be flushed	The device is emitting ultrasound	Stop ultrasound by stepping on the footswitch (Figure 6, #2)
	There is no line plugged into the handpiece	Check that the irrigation line is correctly connected to the handpiece
Active connection (left/right) cannot be switched	The device is emitting ultrasound	Stop ultrasound by stepping on the footswitch (Figure 6, #2)
	There is no cord plugged into the connector	Check that the handpiece cord is correctly connected to the connector

7. 2 REPLACEMENT OF FAULTY COMPONENTS

For the technical staff of agents in its distribution network, SATELEC will provide all the information necessary for the repair of faulty components that SATELEC considers as reparable.

- Used irrigation lines should be disposed of in a special contaminated waste container.
- Used tips should be disposed of in special safety containers.

IX - LIABILITY

VIII - DISPOSAL AND RECYCLING

In no circumstances should worn out equipment be disposed of by the user.

The presence of this symbol  means that this product should not be treated as household waste but must be processed separately, in line with Directive 2002/96/CE.

By making sure that this product is disposed of correctly, you are helping to protect the environment.

For more information on the recycling of this product, please contact your local authority, an official dump or the outlet where you bought the product.

The manufacturer is not liable if:

- the manufacturer's recommendations have not been followed (see Chapter 2),
- repairs have been performed by anybody other than an engineer authorized by the manufacturer
- the device has been used in an electrical installation which does not comply with current standards,
- the device has been used in a way which is not stipulated in this manual,
- the instructions in this document have not been followed.

NB: the manufacturer reserves the right to change the device and/or the User's Manual without prior notice.

X - MANUFACTURER

Note the device's serial number before contacting your distributor.

SATELEC
17, avenue Gustave Eiffel
BP 30216
33708 MERIGNAC - France
Tel : +33 (0) 556 34 06 07
Fax : +33 (0) 556 34 92 92
E.mail : satelec@acteongroup.com
www.acteongroup.com

XI - SYMBOLOGY

	Alternating current
	BF device
	Warning, please refer to the accompanying documentation
	CE Marking
	Not to be disposed of as household waste
	Steam-sterilization (autoclave)
	On/Off

Note: refer to the label on the back of the device.

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