Vari**Lase*[®]

Laser Fiber and Procedure Kits

VARI-LASE[®] BRIGHT TIP[™] Endovenous Laser Procedure Kit and Components

Instructions For Use

USA CAUTION

Federal (USA) law restricts this device to sale by or on the order of a physician.

CAUTION

The VARI-LASE Bright Tip procedure should be performed by physicians with adequate training in the use of the device.

DEVICE DESCRIPTION

The VARI-LASE Bright Tip endovenous laser procedure kit and components are disposable items used in conjunction with a solid state diode laser console operating at wavelengths of 810nm, 940nm, 980nm, or 1470nm and a maximum power of 14W.

Each package contains one or more of the following components:

- 600µm core laser fiber, 3.5m, with or without markings, with ceramic distal tip
- Fiber lock
- Introducer sheath with or without marks (multiple lengths)
- Vari-Lase Flex[™] catheter with radiopaque marker at distal tip (multiple lengths)
- 0.035" (0.89mm) or 0.018" (0.45mm) guidewires
- 19 Gauge or 21 Gauge percutaneous entry needle
- Micro-introducer kit or micro-introducer components consisting of:
 - 21 Gauge percutaneous entry needle
 - Micro-access introducer with dilator
 - 0.018" (0.45mm) guidewire

INDICATIONS

The VARI-LASE Bright Tip kit is indicated for the treatment of varicose veins and varicosities associated with superficial reflux of the Great Saphenous Vein, and for the treatment of incompetence and reflux of superficial veins in the lower extremity.

CONTRAINDICATIONS

The VARI-LASE Bright Tip procedure is contraindicated in patients with an aneurysmal section in the vein segment to be treated.

The VARI-LASE Bright Tip procedure is contraindicated in patients with severe peripheral vascular disease, as evidenced by an ankle-brachial index of <0.5.

The VARI-LASE Bright Tip procedure is contraindicated in patients with thrombus in the vein segment to be treated.

WARNINGS

The VARI-LASE Bright Tip components are supplied sterile for single use only. Do not resterilize and/or reuse any components.

Treatment of a vein located close to the skin surface may result in a skin burn.

Paresthesia may occur from thermal damage to adjacent sensory nerves.

When treating the Great Saphenous Vein, the distal tip of the VARI-LASE Bright Tip laser fiber must not be positioned closer than 2cm distal to the saphenofemoral junction, as verified by ultrasound visualization or fluoroscopy, to protect from venous thrombus embolization.

Use of excessive energy has been reported to compromise the integrity of the laser fiber and/or sheath and may result in embolization of device components.

Appropriate eye protection must be worn by the patient and all operating personnel to protect from damage by direct or reflected laser energy.

The marks on the introducer sheath, catheter and VARI-LASE Bright Tip laser fiber have been calibrated to the VARI-LASE Bright Tip components. Do not use alternative introducer sheaths or laser fibers, or incorrect positioning, as vessel damage may occur.

COMPLICATIONS

As with all medical procedures, complications may occur. For this procedure, these may include:

- vessel perforation
- thrombosis
- pulmonary embolism
- phlebitis
- hematoma
- infection
- paresthesia
- skin burns

PRECAUTIONS

The VARI-LASE Bright Tip procedure should be performed by physicians thoroughly trained in percutaneous endovenous techniques and procedures.

Do not use the VARI-LASE Bright Tip product if the packaging has been damaged.

Inspect the VARI-LASE Bright Tip laser fiber and introducer sheath prior to use for any bends, kinks or breaks. Do not use a damaged VARI-LASE Bright Tip laser fiber or introducer sheath, as vessel damage and/or inability to advance or withdraw the VARI-LASE Bright Tip laser fiber may occur.

Exercise care in handling of the VARI-LASE Bright Tip laser fiber during a procedure to reduce the possibility of accidental breakage, bending or kinking.

When removing the VARI-LASE Bright Tip laser fiber from the packaging, first remove the tip gently from the foam protection. Then, lift the remaining laser fiber from the retention grooves of the laser fiber tray. Do not slide the Bright Tip laser fiber through the retention grooves in order to reduce the possibility of accidental breakage of the laser fiber or the ceramic tip.

Avoid touching the needle tip to the sheath and laser fiber when injecting anesthesia, as this could damage the protective cladding on the laser fiber.

Connect the VARI-LASE Bright Tip laser fiber to the laser console immediately after the protective cap is removed to prevent contamination of the laser fiber optic surface.

When the VARI-LASE Bright Tip laser fiber is in the body, it should be manipulated only under ultrasound visualization. Do not attempt to move the VARI-LASE Bright Tip laser fiber without observing the resultant tip response.

When treating shorter vein segments, the distal tip of the VARI-LASE Bright Tip laser fiber must maintain a safe distance from the deep venous system as confirmed by ultrasound visualization or fluoroscopy to protect from venous thrombus embolization.

Never advance or withdraw an endovascular device against resistance until the cause of the resistance is determined. Movement of the VARI-LASE Bright Tip laser fiber against resistance may result in damage to the laser fiber or vessel perforation.

English/Instructions for Use1

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Vascular Solutions, Inc.
6464 Sycamore Court North
Minneapolis, MN 55369
USA
(888) 240-6001 USA
 (763) 656-4300
 (763) 656-4250
www.vasc.com

The VARI-LASE Flex catheter is not rated for injection of fluid at more than 200psi.

CLINICAL PROCEDURE

The following instructions provide technical direction but do not obviate the necessity of formal training in the use of the VARI-LASE Bright Tip procedure. The techniques and procedures described do not represent ALL medically acceptable protocols, nor are they intended as a substitute for the physician's experience and judgment in treating any specific patient.

Preparations for Use

1. Prepare the laser console in accordance with its operating instructions. Verify that the connection of the VARI-LASE Bright Tip laser fiber to the console is appropriate.
2. Map the vessel treatment area using duplex ultrasound, being careful to mark the vessel location on the patient's skin for a guide in treatment.
3. Prep and drape the limb in sterile fashion and wrap the ultrasound transducer with a sterile cover.
4. Using sterile technique, open the VARI-LASE Bright Tip pouches and transfer the components into the sterile field.
5. Remove the components from the packaging and inspect for any defects or missing components.

Access Procedure

1. Cannulate the vein using the 19 Gauge or 21 Gauge access needle.
2. Insert the corresponding guidewire into the vessel and remove the needle over the wire. Insert the introducer sheath and dilator over the wire and into the vessel.
3. A Y-adapter guide tube is provided to facilitate the use of a J-tip guidewire with the VARI-LASE Flex catheter.
4. Position the introducer sheath at the desired treatment site.

NOTE: Marked introducer sheaths or catheters have marks located at intervals of 1cm to aid in location during placement and removal.
5. Remove the wire and dilator and flush the introducer sheath or catheter using standard technique.

Treatment Procedure

WARNING: For the remainder of the procedure, laser safety glasses must be worn by the patient and all operating personnel.

1. If using a VARI-LASE Bright Tip sheath-specific marked laser fiber with a pre-positioned fiber lock, pass the connector end of the VARI-LASE Bright Tip laser fiber outside of the sterile field.
2. Outside of the sterile field, remove the cap from the VARI-LASE Bright Tip laser fiber connector and connect the VARI-LASE Bright Tip laser fiber to the laser console.
3. Insert the VARI-LASE Bright Tip laser fiber following one of the two following procedures:
 - If the VARI-LASE Bright Tip laser fiber is supplied without a fiber lock, insert the VARI-LASE Bright Tip laser fiber into the introducer sheath and advance until approximately 2.5cm of the Bright Tip laser fiber is exposed beyond the distal tip of the introducer sheath, using ultrasound for guidance. Note: The ceramic tip of the VARI-LASE Bright Tip
4. Confirm the position of the VARI-LASE Bright Tip laser fiber using ultrasound or fluoroscopy, if appropriate. If treating the Great Saphenous Vein, the VARI-LASE Bright Tip laser fiber should be exposed from the introducer sheath and located at least 2cm distal to the saphenofemoral junction. The VARI-LASE Bright Tip laser fiber should be exposed from the catheter ~2.5cm, and the red aiming beam should be visible through the skin.
5. The ceramic tip of the VARI-LASE Bright Tip laser fiber prevents the vessel wall from contacting the laser fiber transmission core. When using the VARI-LASE Bright Tip laser fiber, anesthetic delivery should be adjusted accordingly. Bathe the surrounding tissue with dilute anesthetic to provide thermal protection without causing vessel compression or vein wall apposition. Due to the VARI-LASE Bright Tip protection of the vessel wall, lower power settings, slower pull-back rates and significantly less anesthesia are recommended, all as compared to bare tip laser fibers.
6. Place the laser console in the "Ready" mode.
7. Set power level of the console to the settings recommended in the operating instructions of the laser console.

laser fiber is highly echogenic and is approximately 10mm in length.

- If the VARI-LASE Bright Tip laser fiber is supplied with a fiber lock, the lock will be pre-positioned proximal to the most proximal mark of the Bright Tip laser fiber. If it has moved during shipping, position the fiber lock appropriately and tighten into place. The laser fiber should then be advanced into the sheath until the distal mark reaches the hub of the sheath until the distal mark reaches the hub of the sheath or catheter. Then, pull back the introducer sheath and snap the fiber lock into the end of the sheath. The laser fiber tip will be exposed approximately 2.5cm from the end of the sheath.
- When using the Flex catheter, advance the Bright Tip laser fiber until the distal mark reaches the hemostasis valve, then pull the Flex catheter back to the proximal mark to reveal 2.5cm of laser fiber. Tighten the valve to secure the laser fiber positioning.
- When using the VARI-LASE Short kit with a Vascular Solutions micro-introducer sheath, the fiber lock will be pre-positioned proximal to the mark on the laser fiber. If it has moved during shipping, position the fiber lock appropriately and tighten in place. The laser fiber should then be advanced through the micro-introducer sheath until the single mark reaches the hub of the micro-introducer sheath. Then, pull back the micro-introducer sheath and twist to secure the fiber lock to the micro-introducer sheath and expose 1cm of laser fiber from the end of the sheath.

WARNING: When treating the Great Saphenous Vein, the distal tip of the VARI-LASE Bright Tip laser fiber must be positioned not closer than 2cm distal to the saphenofemoral junction as verified by ultrasound visualization.

4. Turn down the overhead lights and check that the patient and all personnel are wearing laser safety glasses.
5. Hold the introducer sheath or catheter by the hub. Press down the foot pedal to fire the laser and simultaneously withdraw the introducer sheath and the VARI-LASE Bright Tip laser fiber while delivering approximately 50-70 joules/cm. Do not compress or attempt to place the laser fiber in contact with vein wall.
6. When treating shorter vein segments, the power setting of the Vari-Lase console can be lowered and speed of withdrawal slowed to deliver approximately 110 joules/cm.
7. After completing the procedure along the entire treatment length, turn the laser to standby, remove the introducer sheath and the VARI-LASE Bright Tip laser fiber from the vessel, and hold compression until bleeding stops.
8. After removing the VARI-LASE Bright Tip laser fiber from the vessel, inspect the laser fiber for visual damage.

Post-Procedure Care

1. Apply a hemostatic bandage over the percutaneous treatment site.
2. Apply a compression stocking along the entire treatment site.
3. The patient may ambulate immediately after hemostasis is achieved. Evaluate the patient after ambulation for pain, bleeding or swelling.
4. Discharge with post-procedure care instructions.
5. Follow-up exams are recommended at intervals determined by the health care provider.

PACKAGING & STORAGE

All VARI-LASE Bright Tip kits and components have been sterilized with ethylene oxide.

STERILE EO LIMITED WARRANTY

Vascular Solutions, Inc. warrants that the VARI-LASE Bright Tip endovenous laser components are free from defects in workmanship and materials prior to the stated expiration date. Liability under this warranty is limited to refund or replacement of any product which has been found by Vascular Solutions, Inc. to be defective in workmanship or materials. Vascular Solutions, Inc. shall not be liable for any incidental, special or consequential damages arising from the use of the VARI-LASE Bright Tip endovenous laser components. Damage to the product through misuse, alteration, improper storage or improper handling shall void this limited warranty.

No employee, agent or distributor of Vascular Solutions, Inc. has any authority to alter or amend this limited warranty in any respect. Any purported alteration or amendment shall not be enforceable against Vascular Solutions, Inc.

THIS WARRANTY IS EXPRESSLY IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE OR ANY OTHER OBLIGATION OF VASCULAR SOLUTIONS, INC.

PATENTS AND TRADEMARKS

May be covered by one or more U.S. or international patents.

See: www.vasc.com/products/patents

Vari-Lase® is a registered trademark of Vascular Solutions, Inc.

Vari-Lase Flex™ and Vari-Lase Bright Tip™ are trademarks of Vascular Solutions, Inc.
See the International Symbols Glossary on page 8.

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Vari-Lase® es una marca registrada de Vascular Solutions, Inc.
Vari-Lase Flex™ y Vari-Lase Bright Tip™ son marcas de Vascular Solutions, Inc.